








General instructions


Titan abutments polishing (Raw-Abutment®)

Required tools and materials

-  Laboratory analogue
-  Diamond polishers/Tungsten-carbide burs
-  Steam cleaner
-  Brush for high-gloss polishing (ø 14 mm, ~15.000 revolutions)
-  Microscope
-  Polishing paste (Diamond Polishing Paste, article code ZBAA9821)
-  Cotton polishing brush
-  Separating disc
-  Product for high-gloss polishing

! * ATTENTION! Please follow respective instructions of the product manufacturer! Differences are possible due to the simplified process description.

Titan abutments preparation



1 Fix the abutment on the laboratory analogue

The laboratory analogue protects the connection with the implant: do not polish the connection!


Hold the abutment by the analogue, by the Raw-Abutment®, with forceps.




2 Apply the polishing paste on the brush for high-gloss polishing

Conventional cementation requires a rougher surface than friction: the steps described here must then be adapted to the technique used.

A smooth surface reduces bacteria deposition: high-gloss polish emergence profiles!



3 Polish the abutment uniformly



4 Separate the abutment from the block with a separating disc

For certain restorations, it could be easier to separate the abutment first.



5 Neaten up the abutment with a diamond or tungsten-carbide bur



6 Polish the abutment uniformly (optional)

! **WARNING:** burs and products for high-gloss polishing must be binder-free, otherwise they might alloy with the metal of the abutment.

! * ATTENTION! Please follow respective instructions of the product manufacturer! Differences are possible due to the simplified process description.

Abutments polishing under the microscope (optional) and finishing



1 Apply the polishing paste on the brush for high-gloss polishing




2 Polish the abutment uniformly

Conventional cementation requires a rougher surface than friction: the steps described here must then be adapted to the technique used.

A smooth surface reduces bacterial growth: high-gloss polish emergence profiles!



3 Polish with the cotton brush and the polishing product for a high-gloss polish*



4 Steam clean the abutment