

Zirkonzahn®

Human Zirconium Technology




100% PRETTAU ZIRCONIA

“Zirconia needs heroes” Enrico Steger

Prettau Zirconia

This case presentation features a Prettau Bridge with occlusal screw retention.

Difficult implant placement necessitated the construction of a separate secondary part which was integrated and retained on the primary frame via transversal screws. In this way frame strength was maintained and aesthetics were optimized. New Screw-Tec-Tools by Zirkonzahn facilitated the transversal screw design. The possibility of actually cutting a thread into zirconia opens the way for a multitude of new applications. Our patient was particularly impressed with the perfect fit and natural appearance of the new bridge.



The Prettau Bridge

“It has become a passion of mine to create smiles for people by using fully anatomically contoured zirconia restorations!”

Jörg Iowa - Tallinn, Estonia




Zirkonzahn®



HUMAN ZIRCONIUM TECHNOLOGY

ZIRKONZAHN World Wide - Tel +39 0474 066 680 - info@zirkonzahn.com - www.zirkonzahn.com



Milling without limits



WEAA5419=