

Zirkonzahn[®]

Human Zirconium Technology




NEAT*

CAD/CAM for Fixed and Removable Restorations made of Prettau® Zirconia

A 68 year old male patient is no longer satisfied with his dental esthetics. The challenge is to restore the full upper and lower jaw without orthodontic and only minimal surgical pretreatment.

To disguise the inappropriate tooth positions and to reestablish esthetic length/width ratios of the crowns, the restoration is planned to include mucosal aspects of pink ceramic. Conventional diagnostic setups are transformed by the CAD/CAM technique into acrylic provisionals. After a clinical adjustment and testing period, the provisionals are exactly copied by double scan into a fixed zirconia restoration for the upper and into a removable zirconia restoration for the lower jaw.



The Prettau® Bridge

Dr. Nico Bühler

Mauricio Cuéllar de la Torre – Zirkonzahn Education Center



Zirkonzahn®



HUMAN ZIRCONIUM TECHNOLOGY

ZIRKONZAHN Worldwide – Tel +39 0474 066 680 – info@zirkonzahn.com – www.zirkonzahn.com



WEAA5432=

Milled with Zirkonzahn CAD/CAM System 5-TEC