

ALPHA-BIO TEC® CONICAL HEX CONNECTION IMPLANTS

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE AB-CH-ABCH (Raw-Abutment® customizable hex and screw)	BSAF11381	20			10-12
SET ZZ BASE AB-CH-CCTB (Conical cemented titan base non hex and screw)	BSAD11331	20			10-12
SET ZZ BASE AB-CH-ICOTH (Impression coping open tray hex and Screw and Driver)	BSAQ11341	15			10-12
SET ZZ BASE AB-CH-NTB (Narrow titan base non hex and screw)	BSAP11331	20			10-12
SET ZZ BASE AB-CH-NTBH (Narrow titan base hex and screw)	BSAP11381	20			10-12
SET ZZ BASE AB-CH-PCTBH (Parallel cemented titan base hex and screw)	BSAB11381	20			10-12
SET ZZ BASE AB-CH-SA (Scananalog and screw)	ZBDE11331	15			10-12
SET ZZ BASE AB-CH-SM (Scanmarker and screw)	ZBAD11331	20			10-12
ZZ BASE AB-CH-ICSC (Impression coping Screw)	BSFE11301	15			10-12
ZZ BASE AB-CH-SC (Screw)	BSAE11301	20			10-12
ZZ BASE AB-CH-WMSM-L12 (White Metal Scanmarker L12 with Screw)	ZBED11331L12	15			10-12
ZZ BASE AB-CH-ZSC (Zirconia screw)	BSAE11331	20			10-12

ALPHA-BIO TEC® SPI/DFI/ATID

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE AB-SP-ABCH 3.75 (Raw-Abutment® customizable hex and screw)	BSAF1881	20			10-12
SET ZZ BASE AB-SP-ABCH-D14 3.75 (Raw-Abutment® customizable hex D14 and screw)	BSAF1881D14	20			10-12
SET ZZ BASE AB-SP-ASCTBH-K80 3.75 (Angled screw-channel titan base hex K80 and screw)	BSAT1881K80	20			10-12
SET ZZ BASE AB-SP-ASCTB-K80 3.75 (Angled screw-channel titan base non hex K80 and screw)	BSAT1831K80	20			10-12
SET ZZ BASE AB-SP-CCTB 3.75 (Conical cemented titan base non hex and screw)	BSAD1831	20			10-12
SET ZZ BASE AB-SP-CCTB TC (Conical cemented titan base non hex and screw)	BSAD1832	25			10-12
SET ZZ BASE AB-SP-CCTB-K85 3.75 (Conical cemented titan base non hex K85 and screw)	BSAD1831K85	20			10-12
SET ZZ BASE AB-SP-CCTB-K85 TC (Conical cemented titan base non hex K85 and screw)	BSAD1832K85	25			10-12
SET ZZ BASE AB-SP-ICOT TC (Impression coping open tray non hex and Screw and Driver)	BSAQ1862	15			10-12
SET ZZ BASE AB-SP-ICOTH 3.75 (Impression coping open tray hex and Screw and Driver)	BSAQ1841	15			10-12
SET ZZ BASE AB-SP-PCTBH 3.75 (Parallel cemented titan base hex and screw)	BSAB1881	20			10-12
SET ZZ BASE AB-SP-PCTBH-K85 3.75 (Parallel cemented titan base hex K85 and screw)	BSAB1881K85	20			10-12
SET ZZ BASE AB-SP-SA 3.75 (Scananalog and screw)	ZBDE1831	15			10-12
SET ZZ BASE AB-SP-SA TC (Scananalog and screw)	ZBDE1832	15			10-12
SET ZZ BASE AB-SP-SM 3.75 (Scanmarker and screw)	ZBAD1831	20			10-12
SET ZZ BASE AB-SP-SM TC (Scanmarker and screw)	ZBAD1832	25			10-12
SET ZZ BASE AB-SP-WSM 3.75 (White Scanmarker and screw)	ZBCD1831	20			10-12
SET ZZ BASE AB-SP-WSM TC (White Scanmarker and screw)	ZBCD1832	25			10-12
ZZ BASE AB-SP-ASCSC 3.75 (Angled screw-channel Screw)	BSKE1801	20			10-12
ZZ BASE AB-SP-ICSC 3.75 (Impression coping Screw)	BSFE1801	15			10-12
ZZ BASE AB-SP-ICSC TC (Impression coping Screw)	BSFE1802	15			10-12
ZZ BASE AB-SP-MI 3.75 (Model implant)	MIAG1801				10-12
ZZ BASE AB-SP-MI TC (Model implant)	MIAG1802				10-12
ZZ BASE AB-SP-SC 3.75 (Screw)	BSAE1801	20			10-12
ZZ BASE AB-SP-SC TC (Screw)	BSAE1802	25			10-12
ZZ BASE AB-SP-SC-G 3.75 (Screw Gold)	BSCE1801	20			10-12
ZZ BASE AB-SP-WMSM-L12 3.75 (White Metal Scanmarker L12 with Screw)	ZBED1831L12	15			10-12
ZZ BASE AB-SP-WMSM-L12 TC (White Metal Scanmarker L12 with Screw)	ZBED1832L12	15			10-12
ZZ BASE AB-SP-ZSC 3.75 (Zirconia screw)	BSAE1831	20			10-12

ALPHA-BIO TEC® SPI/DFI/ATID

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE AB-SP-ZSC TC (Zirconia screw)	BSAE1832	25			10-12

ANTHOGYR® AXIOM®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE AG-AX-ABCH REG/PX (Raw-Abutment® customizable hex and screw)	BSAF8381	25			10-12
SET ZZ BASE AG-AX-ABCH-D14 REG/PX (Raw-Abutment® customizable hex D14 and screw)	BSAF8381D14	25			10-12
SET ZZ BASE AG-AX-ASCTBH-K80 REG/PX (Angled screw-channel titan base hex K80 and screw)	BSAT8381K80	20			10-12
SET ZZ BASE AG-AX-ASCTB-K80 REG/PX (Angled screw-channel titan base non hex K80 and screw)	BSAT8331K80	20			10-12
SET ZZ BASE AG-AX-CCTB ACA (Conical cemented titan base non hex and screw)	BSAD8333	25			10-12
SET ZZ BASE AG-AX-CCTB MUA 4.0 (Conical cemented titan base non hex and screw)	BSAD8334	15			10-12
SET ZZ BASE AG-AX-CCTB MUA 4.8 (Conical cemented titan base non hex and screw)	BSAD8335	15			10-12
SET ZZ BASE AG-AX-CCTB REG/PX (Conical cemented titan base non hex and screw)	BSAD8331	25			10-12
SET ZZ BASE AG-AX-CCTB SCA (Conical cemented titan base non hex and screw)	BSAD8332	25			10-12
SET ZZ BASE AG-AX-CCTB-K85 MUA 4.0 (Conical cemented titan base non hex K85 and screw)	BSAD8334K85	15			10-12
SET ZZ BASE AG-AX-CCTB-K85 MUA 4.8 (Conical cemented titan base non hex K85 and screw)	BSAD8335K85	15			10-12
SET ZZ BASE AG-AX-CCTB-K85 REG/PX (Conical cemented titan base non hex K85 and screw)	BSAD8331K85	25			10-12
SET ZZ BASE AG-AX-CCTB-L10 REG/PX (Conical cemented titan base non hex L10 and screw)	BSAD8331L10	25			10-12
SET ZZ BASE AG-AX-CCTB-L20 REG/PX (Conical cemented titan base non hex L20 and screw)	BSAD8331L20	25			10-12
SET ZZ BASE AG-AX-CCTB-L30 REG/PX (Conical cemented titan base non hex L30 and screw)	BSAD8331L30	25			10-12
SET ZZ BASE AG-AX-CCTB-L40 REG/PX (Conical cemented titan base non hex L40 and screw)	BSAD8331L40	25			10-12
SET ZZ BASE AG-AX-ICOT MUA 4.0 (Impression coping open tray non hex and Screw and Driver)	BSAQ8364	15			10-12
SET ZZ BASE AG-AX-ICOT MUA 4.8 (Impression coping open tray non hex and Screw and Driver)	BSAQ8365	15			10-12
SET ZZ BASE AG-AX-ICOTH REG/PX (Impression coping open tray hex and Screw and Driver)	BSAQ8341	15			10-12
SET ZZ BASE AG-AX-NTB REG/PX (Narrow titan base non hex and screw)	BSAP8331	25			10-12
SET ZZ BASE AG-AX-NTBH REG/PX (Narrow titan base hex and screw)	BSAP8381	25			10-12
SET ZZ BASE AG-AX-PCTBH REG/PX (Parallel cemented titan base hex and screw)	BSAB8381	25			10-12
SET ZZ BASE AG-AX-PCTBH-K85 REG/PX (Parallel cemented titan base hex K85 and screw)	BSAB8381K85	25			10-12
SET ZZ BASE AG-AX-PCTBH-L10 REG/PX (Parallel cemented titan base hex L10 and screw)	BSAB8391L10	25			10-12
SET ZZ BASE AG-AX-PCTBH-L20 REG/PX (Parallel cemented titan base hex L20 and screw)	BSAB8391L20	25			10-12
SET ZZ BASE AG-AX-PCTBH-L30 REG/PX (Parallel cemented titan base hex L30 and screw)	BSAB8391L30	25			10-12
SET ZZ BASE AG-AX-PCTBH-L40 REG/PX (Parallel cemented titan base hex L40 and screw)	BSAB8391L40	25			10-12
SET ZZ BASE AG-AX-SA ACA (Scananalog and screw)	ZBDE8333	15			10-12
SET ZZ BASE AG-AX-SA MUA 4.0 (Scananalog and screw)	ZBDE8334	15			10-12

ANTHOGYR® AXIOM®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE AG-AX-SA MUA 4.8 (Scanalog and screw)	ZBDE8335	15			10-12
SET ZZ BASE AG-AX-SA REG/PX (Scanalog and screw)	ZBDE8331	15			10-12
SET ZZ BASE AG-AX-SA SCA (Scanalog and screw)	ZBDE8332	15			10-12
SET ZZ BASE AG-AX-SM ACA (Scanmarker and screw)	ZBAD8333	25			10-12
SET ZZ BASE AG-AX-SM MUA 4.0 (Scanmarker and screw)	ZBAD8334	15			10-12
SET ZZ BASE AG-AX-SM MUA 4.8 (Scanmarker and screw)	ZBAD8335	15			10-12
SET ZZ BASE AG-AX-SM REG/PX (Scanmarker and screw)	ZBAD8331	25			10-12
SET ZZ BASE AG-AX-SM SCA (Scanmarker and screw)	ZBAD8332	25			10-12
SET ZZ BASE AG-AX-WSM ACA (White Scanmarker and screw)	ZBCD8333	25			10-12
SET ZZ BASE AG-AX-WSM MUA 4.0 (White Scanmarker and screw)	ZBCD8334	15			10-12
SET ZZ BASE AG-AX-WSM MUA 4.8 (White Scanmarker and screw)	ZBCD8335	15			10-12
SET ZZ BASE AG-AX-WSM REG/PX (White Scanmarker and screw)	ZBCD8331	25			10-12
SET ZZ BASE AG-AX-WSM SCA (White Scanmarker and screw)	ZBCD8332	25			10-12
ZZ BASE AG-AX-ASCSC REG/PX (Angled screw-channel Screw)	BSKE8301	20			10-12
ZZ BASE AG-AX-ICSC MUA 4.0 (Impression coping Screw)	BSFE8304	15			10-12
ZZ BASE AG-AX-ICSC MUA 4.8 (Impression coping Screw)	BSFE8305	15			10-12
ZZ BASE AG-AX-ICSC REG/PX (Impression coping Screw)	BSFE8301	15			10-12
ZZ BASE AG-AX-MI REG/PX (Model implant)	MIAG8301				10-12
ZZ BASE AG-AX-SC ACA (Screw)	BSAE8303	25			10-12
ZZ BASE AG-AX-SC MUA 4.0 (Screw)	BSAE8304	15			10-12
ZZ BASE AG-AX-SC MUA 4.8 (Screw)	BSAE8305	15			10-12
ZZ BASE AG-AX-SC REG/PX (Screw)	BSAE8301	25			10-12
ZZ BASE AG-AX-SC SCA (Screw)	BSAE8302	25			10-12
ZZ BASE AG-AX-WMSM-L12 MUA 4.0 (White Metal Scanmarker L12 with Screw)	ZBED8334L12	15			10-12
ZZ BASE AG-AX-WMSM-L12 MUA 4.8 (White Metal Scanmarker L12 with Screw)	ZBED8335L12	15			10-12
ZZ BASE AG-AX-WMSM-L12 REG/PX (White Metal Scanmarker L12 with Screw)	ZBED8331L12	15			10-12
ZZ BASE AG-AX-ZSC ACA (Zirconia Screw)	BSAE8333	25			10-12
ZZ BASE AG-AX-ZSC MUA 4.0 (Zirconia Screw)	BSAE8334	15			10-12
ZZ BASE AG-AX-ZSC MUA 4.8 (Zirconia Screw)	BSAE8335	15			10-12

ANTHOGR® AXIOM®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE AG-AX-ZSC REG/PX (Zirconia Screw)	BSAE8331	25			10-12
ZZ BASE AG-AX-ZSC SCA (Zirconia Screw)	BSAE8332	25			10-12

ASTRA TECH IMPLANT SYSTEM EV/UNIABUTMENT® EV

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE AS-EV-ABCH 3.0 (Raw-Abutment® customizable hex and screw)	BSAF7681	25			10-12
SET ZZ BASE AS-EV-ABCH 3.6 (Raw-Abutment® customizable hex and screw)	BSAF7682	25			10-12
SET ZZ BASE AS-EV-ABCH 4.2 (Raw-Abutment® customizable hex and screw)	BSAF7683	25			10-12
SET ZZ BASE AS-EV-ABCH 4.8 (Raw-Abutment® customizable hex and screw)	BSAF7684	25			10-12
SET ZZ BASE AS-EV-ABCH 5.4 (Raw-Abutment® customizable hex and screw)	BSAF7685	25			10-12
SET ZZ BASE AS-EV-ABCH-D14 3.0 (Raw-Abutment® customizable hex D14 and screw)	BSAF7681D14	25			10-12
SET ZZ BASE AS-EV-ABCH-D14 3.6 (Raw-Abutment® customizable hex D14 and screw)	BSAF7682D14	25			10-12
SET ZZ BASE AS-EV-ABCH-D14 4.2 (Raw-Abutment® customizable hex D14 and screw)	BSAF7683D14	25			10-12
SET ZZ BASE AS-EV-ABCH-D14 4.8 (Raw-Abutment® customizable hex D14 and screw)	BSAF7684D14	25			10-12
SET ZZ BASE AS-EV-ABCH-D14 5.4 (Raw-Abutment® customizable hex D14 and screw)	BSAF7685D14	25			10-12
SET ZZ BASE AS-EV-ASCTBH-K80 3.0 (Angled screw-channel titan base hex K80 and screw)	BSAT7681K80	20			10-12
SET ZZ BASE AS-EV-ASCTBH-K80 3.6 (Angled screw-channel titan base hex K80 and screw)	BSAT7682K80	20			10-12
SET ZZ BASE AS-EV-ASCTBH-K80 4.2 (Angled screw-channel titan base hex K80 and screw)	BSAT7683K80	20			10-12
SET ZZ BASE AS-EV-ASCTBH-K80 Profile4.2 (Angled screw-channel titan base hex K80 and screw)	BSAT7687K80	20			10-12
SET ZZ BASE AS-EV-ASCTB-K80 3.0 (Angled screw-channel titan base non hex K80 and screw)	BSAT7631K80	20			10-12
SET ZZ BASE AS-EV-ASCTB-K80 3.6 (Angled screw-channel titan base non hex K80 and screw)	BSAT7632K80	20			10-12
SET ZZ BASE AS-EV-ASCTB-K80 4.2 (Angled screw-channel titan base non hex K80 and screw)	BSAT7633K80	20			10-12
SET ZZ BASE AS-EV-ASCTB-K80 P4.2 (Angled screw-channel titan base non hex K80 and screw)	BSAT7637K80	20			10-12
SET ZZ BASE AS-EV-CCTB 3.0 (Conical cemented titan base non hex and screw)	BSAD7631	25			10-12
SET ZZ BASE AS-EV-CCTB 3.6 (Conical cemented titan base non hex and screw)	BSAD7632	25			10-12
SET ZZ BASE AS-EV-CCTB 4.2 (Conical cemented titan base non hex and screw)	BSAD7633	25			10-12
SET ZZ BASE AS-EV-CCTB 4.8 (Conical cemented titan base non hex and screw)	BSAD7634	25			10-12
SET ZZ BASE AS-EV-CCTB 5.4 (Conical cemented titan base non hex and screw)	BSAD7635	25			10-12
SET ZZ BASE AS-EV-CCTB Profile4.2 (Conical cemented titan base non hex and screw)	BSAD7637	25			10-12
SET ZZ BASE AS-EV-CCTB Profile4.8 (Conical cemented titan base non hex and screw)	BSAD7638	25			10-12
SET ZZ BASE AS-EV-CCTB UA (Conical cemented titan base non hex and screw)	BSAD7636	15			10-12
SET ZZ BASE AS-EV-CCTB-G 3.0 (Conical cemented titan base non hex and screw Gold)	BSCD7631	25			10-12
SET ZZ BASE AS-EV-CCTB-G 3.6 (Conical cemented titan base non hex and screw Gold)	BSCD7632	25			10-12
SET ZZ BASE AS-EV-CCTB-G 4.2 (Conical cemented titan base non hex and screw Gold)	BSCD7633	25			10-12

ASTRA TECH IMPLANT SYSTEM EV/UNIABUTMENT® EV

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE AS-EV-CCTB-G Profile4.2 (Conical cemented titan base non hex and screw Gold)	BSCD7637	25			10-12
SET ZZ BASE AS-EV-CCTB-G Profile4.8 (Conical cemented titan base non hex and screw Gold)	BSCD7638	25			10-12
SET ZZ BASE AS-EV-CCTB-G UA (Conical cemented titan base non hex and screw Gold)	BSCD7636	15			10-12
SET ZZ BASE AS-EV-CCTB-K85 3.0 (Conical cemented titan base non hex K85 and screw)	BSAD7631K85	25			10-12
SET ZZ BASE AS-EV-CCTB-K85 3.6 (Conical cemented titan base non hex K85 and screw)	BSAD7632K85	25			10-12
SET ZZ BASE AS-EV-CCTB-K85 4.2 (Conical cemented titan base non hex K85 and screw)	BSAD7633K85	25			10-12
SET ZZ BASE AS-EV-CCTB-K85 4.8 (Conical cemented titan base non hex K85 and screw)	BSAD7634K85	25			10-12
SET ZZ BASE AS-EV-CCTB-K85 5.4 (Conical cemented titan base non hex K85 and screw)	BSAD7635K85	25			10-12
SET ZZ BASE AS-EV-CCTB-K85 Profile4.2 (Conical cemented titan base non hex K85 and screw)	BSAD7637K85	25			10-12
SET ZZ BASE AS-EV-CCTB-K85 Profile4.8 (Conical cemented titan base non hex K85 and screw)	BSAD7638K85	25			10-12
SET ZZ BASE AS-EV-CCTB-K85 UA (Conical cemented titan base non hex K85 and screw)	BSAD7636K85	15			10-12
SET ZZ BASE AS-EV-CCTB-L15 3.6 (Conical cemented titan base non hex L15 and screw)	BSAD7632L15	25			10-12
SET ZZ BASE AS-EV-CCTB-L15 4.2 (Conical cemented titan base non hex L15 and screw)	BSAD7633L15	25			10-12
SET ZZ BASE AS-EV-CCTB-L15 4.8 (Conical cemented titan base non hex L15 and screw)	BSAD7634L15	25			10-12
SET ZZ BASE AS-EV-CCTB-L15 5.4 (Conical cemented titan base non hex L15 and screw)	BSAD7635L15	25			10-12
SET ZZ BASE AS-EV-CCTB-L15-G 3.6 (Conical cemented titan base non hex L15 and screw Gold)	BSCD7632L15	25			10-12
SET ZZ BASE AS-EV-CCTB-L20 3.0 (Conical cemented titan base non hex L20 and screw)	BSAD7631L20	25			10-12
SET ZZ BASE AS-EV-CCTB-L20 3.6 (Conical cemented titan base non hex L20 and screw)	BSAD7632L20	25			10-12
SET ZZ BASE AS-EV-CCTB-L20 4.2 (Conical cemented titan base non hex L20 and screw)	BSAD7633L20	25			10-12
SET ZZ BASE AS-EV-CCTB-L20 4.8 (Conical cemented titan base non hex L20 and screw)	BSAD7634L20	25			10-12
SET ZZ BASE AS-EV-CCTB-L20 5.4 (Conical cemented titan base non hex L20 and screw)	BSAD7635L20	25			10-12
SET ZZ BASE AS-EV-CCTB-L20-G 3.0 (Conical cemented titan base non hex L20 and screw Gold)	BSCD7631L20	25			10-12
SET ZZ BASE AS-EV-CCTB-L20-G 3.6 (Conical cemented titan base non hex L20 and screw Gold)	BSCD7632L20	25			10-12
SET ZZ BASE AS-EV-CCTB-L20-G 4.2 (Conical cemented titan base non hex L20 and screw Gold)	BSCD7633L20	25			10-12
SET ZZ BASE AS-EV-CCTB-L20-G 4.8 (Conical cemented titan base non hex L20 and screw Gold)	BSCD7634L20	25			10-12
SET ZZ BASE AS-EV-CCTB-L30 3.0 (Conical cemented titan base non hex L30 and screw)	BSAD7631L30	25			10-12
SET ZZ BASE AS-EV-CCTB-L30 3.6 (Conical cemented titan base non hex L30 and screw)	BSAD7632L30	25			10-12
SET ZZ BASE AS-EV-CCTB-L30 4.2 (Conical cemented titan base non hex L30 and screw)	BSAD7633L30	25			10-12
SET ZZ BASE AS-EV-CCTB-L30 4.8 (Conical cemented titan base non hex L30 and screw)	BSAD7634L30	25			10-12

ASTRA TECH IMPLANT SYSTEM EV/UNIABUTMENT® EV

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE AS-EV-CCTB-L30 5.4 (Conical cemented titan base non hex L30 and screw)	BSAD7635L30	25			10-12
SET ZZ BASE AS-EV-CCTB-L30-G 3.0 (Conical cemented titan base non hex L30 and screw Gold)	BSCD7631L30	25			10-12
SET ZZ BASE AS-EV-CCTB-L30-G 3.6 (Conical cemented titan base non hex L30 and screw Gold)	BSCD7632L30	25			10-12
SET ZZ BASE AS-EV-CCTB-L30-G 4.2 (Conical cemented titan base non hex L30 and screw Gold)	BSCD7633L30	25			10-12
SET ZZ BASE AS-EV-CCTB-L30-G 4.8 (Conical cemented titan base non hex L30 and screw Gold)	BSCD7634L30	25			10-12
SET ZZ BASE AS-EV-CCTB-L40 3.0 (Conical cemented titan base non hex L40 and screw)	BSAD7631L40	25			10-12
SET ZZ BASE AS-EV-CCTB-L40 3.6 (Conical cemented titan base non hex L40 and screw)	BSAD7632L40	25			10-12
SET ZZ BASE AS-EV-CCTB-L40 4.2 (Conical cemented titan base non hex L40 and screw)	BSAD7633L40	25			10-12
SET ZZ BASE AS-EV-CCTB-L40 4.8 (Conical cemented titan base non hex L40 and screw)	BSAD7634L40	25			10-12
SET ZZ BASE AS-EV-CCTB-L40 5.4 (Conical cemented titan base non hex L40 and screw)	BSAD7635L40	25			10-12
SET ZZ BASE AS-EV-CCTB-L40-G 3.0 (Conical cemented titan base non hex L40 and screw Gold)	BSCD7631L40	25			10-12
SET ZZ BASE AS-EV-CCTB-L40-G 3.6 (Conical cemented titan base non hex L40 and screw Gold)	BSCD7632L40	25			10-12
SET ZZ BASE AS-EV-CCTB-L40-G 4.2 (Conical cemented titan base non hex L40 and screw Gold)	BSCD7633L40	25			10-12
SET ZZ BASE AS-EV-ICOT 3.0 (Impression coping open tray non hex and Screw and Driver)	BSAQ7661	15			10-12
SET ZZ BASE AS-EV-ICOT 3.6 (Impression coping open tray non hex and Screw and Driver)	BSAQ7662	15			10-12
SET ZZ BASE AS-EV-ICOT 4.2 (Impression coping open tray non hex and Screw and Driver)	BSAQ7663	15			10-12
SET ZZ BASE AS-EV-ICOT 4.8 (Impression coping open tray non hex and Screw and Driver)	BSAQ7664	15			10-12
SET ZZ BASE AS-EV-ICOT 5.4 (Impression coping open tray non hex and Screw and Driver)	BSAQ7665	15			10-12
SET ZZ BASE AS-EV-ICOT Profile4.2 (Impression coping open tray non hex and Screw and Driver)	BSAQ7667	15			10-12
SET ZZ BASE AS-EV-ICOT Profile4.8 (Impression coping open tray non hex and Screw and Driver)	BSAQ7668	15			10-12
SET ZZ BASE AS-EV-ICOT UA (Impression coping open tray non hex and Screw and Driver)	BSAQ7666	15			10-12
SET ZZ BASE AS-EV-ICOTH 3.0 (Impression coping open tray hex and Screw and Driver)	BSAQ7641	15			10-12
SET ZZ BASE AS-EV-ICOTH 3.6 (Impression coping open tray hex and Screw and Driver)	BSAQ7642	15			10-12
SET ZZ BASE AS-EV-ICOTH 4.2 (Impression coping open tray hex and Screw and Driver)	BSAQ7643	15			10-12
SET ZZ BASE AS-EV-ICOTH 4.8 (Impression coping open tray hex and Screw and Driver)	BSAQ7644	15			10-12
SET ZZ BASE AS-EV-ICOTH 5.4 (Impression coping open tray hex and Screw and Driver)	BSAQ7645	15			10-12
SET ZZ BASE AS-EV-ICOTH Profile4.2 (Impression coping open tray hex and Screw and Driver)	BSAQ7647	15			10-12
SET ZZ BASE AS-EV-ICOTH Profile4.8 (Impression coping open tray hex and Screw and Driver)	BSAQ7648	15			10-12
SET ZZ BASE AS-EV-NTB 3.0 (Narrow titan base non hex and screw)	BSAP7631	25			10-12

All information is subject to change. Errors and omissions excepted. Version: 21.09.2023

ASTRA TECH IMPLANT SYSTEM EV/UNIABUTMENT® EV

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE AS-EV-NTB 3.6 (Narrow titan base non hex and screw)	BSAP7632	25			10-12
SET ZZ BASE AS-EV-NTB 4.2 (Narrow titan base non hex and screw)	BSAP7633	25			10-12
SET ZZ BASE AS-EV-NTB 4.8 (Narrow titan base non hex and screw)	BSAP7634	25			10-12
SET ZZ BASE AS-EV-NTBH 3.0 (Narrow titan base hex and screw)	BSAP7681	25			10-12
SET ZZ BASE AS-EV-NTBH 3.6 (Narrow titan base hex and screw)	BSAP7682	25			10-12
SET ZZ BASE AS-EV-NTBH 4.2 (Narrow titan base hex and screw)	BSAP7683	25			10-12
SET ZZ BASE AS-EV-NTBH 4.8 (Narrow titan base hex and screw)	BSAP7684	25			10-12
SET ZZ BASE AS-EV-PCTBH 3.0 (Parallel cemented titan base hex and screw)	BSAB7681	25			10-12
SET ZZ BASE AS-EV-PCTBH 3.6 (Parallel cemented titan base hex and screw)	BSAB7682	25			10-12
SET ZZ BASE AS-EV-PCTBH 4.2 (Parallel cemented titan base hex and screw)	BSAB7683	25			10-12
SET ZZ BASE AS-EV-PCTBH 4.8 (Parallel cemented titan base hex and screw)	BSAB7684	25			10-12
SET ZZ BASE AS-EV-PCTBH 5.4 (Parallel cemented titan base hex and screw)	BSAB7685	25			10-12
SET ZZ BASE AS-EV-PCTBH Profile4.2 (Parallel cemented titan base hex and screw)	BSAB7687	25			10-12
SET ZZ BASE AS-EV-PCTBH Profile4.8 (Parallel cemented titan base hex and screw)	BSAB7688	25			10-12
SET ZZ BASE AS-EV-PCTBH-G 3.0 (Parallel cemented titan base hex and screw Gold)	BSCB7681	25			10-12
SET ZZ BASE AS-EV-PCTBH-G 3.6 (Parallel cemented titan base hex and screw Gold)	BSCB7682	25			10-12
SET ZZ BASE AS-EV-PCTBH-G 4.2 (Parallel cemented titan base hex and screw Gold)	BSCB7683	25			10-12
SET ZZ BASE AS-EV-PCTBH-G 4.8 (Parallel cemented titan base hex and screw Gold)	BSCB7684	25			10-12
SET ZZ BASE AS-EV-PCTBH-G Profile4.2 (Parallel cemented titan base hex and screw Gold)	BSCB7687	25			10-12
SET ZZ BASE AS-EV-PCTBH-G Profile4.8 (Parallel cemented titan base hex and screw Gold)	BSCB7688	25			10-12
SET ZZ BASE AS-EV-PCTBH-K85 3.0 (Parallel cemented titan base hex K85 and screw)	BSAB7681K85	25			10-12
SET ZZ BASE AS-EV-PCTBH-K85 3.6 (Parallel cemented titan base hex K85 and screw)	BSAB7682K85	25			10-12
SET ZZ BASE AS-EV-PCTBH-K85 4.2 (Parallel cemented titan base hex K85 and screw)	BSAB7683K85	25			10-12
SET ZZ BASE AS-EV-PCTBH-K85 4.8 (Parallel cemented titan base hex K85 and screw)	BSAB7684K85	25			10-12
SET ZZ BASE AS-EV-PCTBH-K85 5.4 (Parallel cemented titan base hex K85 and screw)	BSAB7685K85	25			10-12
SET ZZ BASE AS-EV-PCTBH-K85 Profile4.2 (Parallel cemented titan base hex K85 and screw)	BSAB7687K85	25			10-12
SET ZZ BASE AS-EV-PCTBH-K85 Profile4.8 (Parallel cemented titan base hex K85 and screw)	BSAB7688K85	25			10-12
SET ZZ BASE AS-EV-PCTBH-L15 3.0 (Parallel cemented titan base hex L15 and screw)	BSAB7691L15	25			10-12
SET ZZ BASE AS-EV-PCTBH-L15 3.6 (Parallel cemented titan base hex L15 and screw)	BSAB7692L15	25			10-12

ASTRA TECH IMPLANT SYSTEM EV/UNIABUTMENT® EV

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE AS-EV-PCTBH-L15 4.2 (Parallel cemented titan base hex L15 and screw)	BSAB7693L15	25			10-12
SET ZZ BASE AS-EV-PCTBH-L15 4.8 (Parallel cemented titan base hex L15 and screw)	BSAB7694L15	25			10-12
SET ZZ BASE AS-EV-PCTBH-L15 5.4 (Parallel cemented titan base hex L15 and screw)	BSAB7695L15	25			10-12
SET ZZ BASE AS-EV-PCTBH-L15-G 3.0 (Parallel cemented titan base hex L15 and screw Gold)	BSCB7691L15	25			10-12
SET ZZ BASE AS-EV-PCTBH-L15-G 3.6 (Parallel cemented titan base hex L15 and screw Gold)	BSCB7692L15	25			10-12
SET ZZ BASE AS-EV-PCTBH-L15-G 4.2 (Parallel cemented titan base hex L15 and screw Gold)	BSCB7693L15	25			10-12
SET ZZ BASE AS-EV-PCTBH-L15-G 4.8 (Parallel cemented titan base hex L15 and screw Gold)	BSCB7694L15	25			10-12
SET ZZ BASE AS-EV-PCTBH-L20 3.0 (Parallel cemented titan base hex L20 and screw)	BSAB7691L20	25			10-12
SET ZZ BASE AS-EV-PCTBH-L20 3.6 (Parallel cemented titan base hex L20 and screw)	BSAB7692L20	25			10-12
SET ZZ BASE AS-EV-PCTBH-L20 4.2 (Parallel cemented titan base hex L20 and screw)	BSAB7693L20	25			10-12
SET ZZ BASE AS-EV-PCTBH-L20 4.8 (Parallel cemented titan base hex L20 and screw)	BSAB7694L20	25			10-12
SET ZZ BASE AS-EV-PCTBH-L20 5.4 (Parallel cemented titan base hex L20 and screw)	BSAB7695L20	25			10-12
SET ZZ BASE AS-EV-PCTBH-L20-G 3.0 (Parallel cemented titan base hex L20 and screw Gold)	BSCB7691L20	25			10-12
SET ZZ BASE AS-EV-PCTBH-L20-G 3.6 (Parallel cemented titan base hex L20 and screw Gold)	BSCB7692L20	25			10-12
SET ZZ BASE AS-EV-PCTBH-L20-G 4.2 (Parallel cemented titan base hex L20 and screw Gold)	BSCB7693L20	25			10-12
SET ZZ BASE AS-EV-PCTBH-L30 3.0 (Parallel cemented titan base hex L30 and screw)	BSAB7691L30	25			10-12
SET ZZ BASE AS-EV-PCTBH-L30 3.6 (Parallel cemented titan base hex L30 and screw)	BSAB7692L30	25			10-12
SET ZZ BASE AS-EV-PCTBH-L30 4.2 (Parallel cemented titan base hex L30 and screw)	BSAB7693L30	25			10-12
SET ZZ BASE AS-EV-PCTBH-L30 4.8 (Parallel cemented titan base hex L30 and screw)	BSAB7694L30	25			10-12
SET ZZ BASE AS-EV-PCTBH-L30 5.4 (Parallel cemented titan base hex L30 and screw)	BSAB7695L30	25			10-12
SET ZZ BASE AS-EV-PCTBH-L30-G 3.0 (Parallel cemented titan base hex L30 and screw Gold)	BSCB7691L30	25			10-12
SET ZZ BASE AS-EV-PCTBH-L30-G 3.6 (Parallel cemented titan base hex L30 and screw Gold)	BSCB7692L30	25			10-12
SET ZZ BASE AS-EV-PCTBH-L30-G 4.2 (Parallel cemented titan base hex L30 and screw Gold)	BSCB7693L30	25			10-12
SET ZZ BASE AS-EV-PCTBH-L40 3.0 (Parallel cemented titan base hex L40 and screw)	BSAB7691L40	25			10-12
SET ZZ BASE AS-EV-PCTBH-L40 3.6 (Parallel cemented titan base hex L40 and screw)	BSAB7692L40	25			10-12
SET ZZ BASE AS-EV-PCTBH-L40 4.2 (Parallel cemented titan base hex L40 and screw)	BSAB7693L40	25			10-12
SET ZZ BASE AS-EV-PCTBH-L40 4.8 (Parallel cemented titan base hex L40 and screw)	BSAB7694L40	25			10-12
SET ZZ BASE AS-EV-PCTBH-L40 5.4 (Parallel cemented titan base hex L40 and screw)	BSAB7695L40	25			10-12
SET ZZ BASE AS-EV-PCTBH-L40-G 3.0 (Parallel cemented titan base hex L40 and screw Gold)	BSCB7691L40	25			10-12

ASTRA TECH IMPLANT SYSTEM EV/UNIABUTMENT® EV

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
SET ZZ BASE AS-EV-PCTBH-L40-G 3.6 (Parallel cemented titan base hex L40 and screw Gold)	BSCB7692L40	25			10-12
SET ZZ BASE AS-EV-PCTBH-L40-G 4.2 (Parallel cemented titan base hex L40 and screw Gold)	BSCB7693L40	25			10-12
SET ZZ BASE AS-EV-SA 3.0 (Scananalog and screw)	ZBDE7631	15			10-12
SET ZZ BASE AS-EV-SA 3.6 (Scananalog and screw)	ZBDE7632	15			10-12
SET ZZ BASE AS-EV-SA 4.2 (Scananalog and screw)	ZBDE7633	15			10-12
SET ZZ BASE AS-EV-SA 4.8 (Scananalog and screw)	ZBDE7634	15			10-12
SET ZZ BASE AS-EV-SA 5.4 (Scananalog and screw)	ZBDE7635	15			10-12
SET ZZ BASE AS-EV-SA Profile4.2 (Scananalog and screw)	ZBDE7637	15			10-12
SET ZZ BASE AS-EV-SA Profile4.8 (Scananalog and screw)	ZBDE7638	15			10-12
SET ZZ BASE AS-EV-SA UA (Scananalog and screw)	ZBDE7636	15			10-12
SET ZZ BASE AS-EV-SM 3.0 (Scanmarker and screw)	ZBAD7631	25			10-12
SET ZZ BASE AS-EV-SM 3.6 (Scanmarker and screw)	ZBAD7632	25			10-12
SET ZZ BASE AS-EV-SM 4.2 (Scanmarker and screw)	ZBAD7633	25			10-12
SET ZZ BASE AS-EV-SM 4.8 (Scanmarker and screw)	ZBAD7634	25			10-12
SET ZZ BASE AS-EV-SM 5.4 (Scanmarker and screw)	ZBAD7635	25			10-12
SET ZZ BASE AS-EV-SM Profile4.2 (Scanmarker and screw)	ZBAD7637	25			10-12
SET ZZ BASE AS-EV-SM Profile4.8 (Scanmarker and screw)	ZBAD7638	25			10-12
SET ZZ BASE AS-EV-SM UA (Scanmarker and screw)	ZBAD7636	15			10-12
SET ZZ BASE AS-EV-WSM 3.0 (White Scanmarker and screw)	ZBCD7631	25			10-12
SET ZZ BASE AS-EV-WSM 3.6 (White Scanmarker and screw)	ZBCD7632	25			10-12
SET ZZ BASE AS-EV-WSM 4.8 (White Scanmarker and screw)	ZBCD7634	25			10-12
SET ZZ BASE AS-EV-WSM 5.4 (White Scanmarker and screw)	ZBCD7635	25			10-12
SET ZZ BASE AS-EV-WSM UA (White Scanmarker and screw)	ZBCD7636	15			10-12
SET ZZ BASE AS-EV-ZZMUA17 3.0 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and implk)	BSAI7631	15			10-12
SET ZZ BASE AS-EV-ZZMUA17 3.0 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and implk)	BSAI7641	15			10-12
SET ZZ BASE AS-EV-ZZMUA17 3.6 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and implk)	BSAI7632	20			10-12
SET ZZ BASE AS-EV-ZZMUA17 3.6 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and implk)	BSAI7642	20			10-12
SET ZZ BASE AS-EV-ZZMUA17 4.2 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and implk)	BSAI7633	25			10-12
SET ZZ BASE AS-EV-ZZMUA17 4.2 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and implk)	BSAI7643	25			10-12

ASTRA TECH IMPLANT SYSTEM EV/UNIABUTMENT® EV

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
SET ZZ BASE AS-EV-ZZMUA17 4.8 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and impl)	BSAI7634	30			10-12
SET ZZ BASE AS-EV-ZZMUA17 4.8 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and impl)	BSAI7644	30			10-12
SET ZZ BASE AS-EV-ZZMUA17 5.4 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and impl)	BSAI7635	30			10-12
SET ZZ BASE AS-EV-ZZMUA17 5.4 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and impl)	BSAI7645	30			10-12
SET ZZ BASE AS-EV-ZZMUA17-L10 3.0 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10)	BSAI7631L10	15			10-12
SET ZZ BASE AS-EV-ZZMUA17-L10 3.0 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 L10)	BSAI7641L10	15			10-12
SET ZZ BASE AS-EV-ZZMUA17-L10 3.6 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10)	BSAI7632L10	20			10-12
SET ZZ BASE AS-EV-ZZMUA17-L10 3.6 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 L10)	BSAI7642L10	20			10-12
SET ZZ BASE AS-EV-ZZMUA17-L10 4.2 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10)	BSAI7633L10	25			10-12
SET ZZ BASE AS-EV-ZZMUA17-L10 4.2 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 L10)	BSAI7643L10	25			10-12
SET ZZ BASE AS-EV-ZZMUA17-L10 4.8 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10)	BSAI7634L10	30			10-12
SET ZZ BASE AS-EV-ZZMUA17-L10 4.8 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 L10)	BSAI7644L10	30			10-12
SET ZZ BASE AS-EV-ZZMUA17-L10 5.4 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10)	BSAI7635L10	30			10-12
SET ZZ BASE AS-EV-ZZMUA17-L10 5.4 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 L10)	BSAI7645L10	30			10-12
SET ZZ BASE AS-EV-ZZMUA-L05 5.4 (Zirkonzahn Multi Unit Abutment non hex L05 and Screw)	BSAH7635L05	25	25		10-12
SET ZZ BASE AS-EV-ZZMUA-L05 Profile4.8 (Zirkonzahn Multi Unit Abutment non hex L05 and Sc	BSAH7638L05	25	25		10-12
SET ZZ BASE AS-EV-ZZMUA-L15 3.0 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH7631L15	25	25		10-12
SET ZZ BASE AS-EV-ZZMUA-L15 3.6 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH7632L15	25	25		10-12
SET ZZ BASE AS-EV-ZZMUA-L15 4.2 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH7633L15	25	25		10-12
SET ZZ BASE AS-EV-ZZMUA-L15 4.8 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH7634L15	25	25		10-12
SET ZZ BASE AS-EV-ZZMUA-L15 5.4 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH7635L15	25	25		10-12
SET ZZ BASE AS-EV-ZZMUA-L15 Profile4.2 (Zirkonzahn Multi Unit Abutment non hex L15 and Sc	BSAH7637L15	25	25		10-12
SET ZZ BASE AS-EV-ZZMUA-L15 Profile4.8 (Zirkonzahn Multi Unit Abutment non hex L15 and Sc	BSAH7638L15	25	25		10-12
SET ZZ BASE AS-EV-ZZMUA-L25 3.0 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH7631L25	25	25		10-12
SET ZZ BASE AS-EV-ZZMUA-L25 3.6 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH7632L25	25	25		10-12
SET ZZ BASE AS-EV-ZZMUA-L25 4.2 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH7633L25	25	25		10-12
SET ZZ BASE AS-EV-ZZMUA-L25 4.8 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH7634L25	25	25		10-12
SET ZZ BASE AS-EV-ZZMUA-L25 5.4 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH7635L25	25	25		10-12
SET ZZ BASE AS-EV-ZZMUA-L25 Profile4.2 (Zirkonzahn Multi Unit Abutment non hex L25 and Sc	BSAH7637L25	25	25		10-12

ASTRA TECH IMPLANT SYSTEM EV/UNIABUTMENT® EV

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE AS-EV-ZZMUA-L25 Profile4.8 (Zirkonzahn Multi Unit Abutment non hex L25 and Sc	BSAH7638L25	25	25		10-12
SET ZZ BASE AS-EV-ZZMUA-L35 3.0 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH7631L35	25	25		10-12
SET ZZ BASE AS-EV-ZZMUA-L35 3.6 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH7632L35	25	25		10-12
SET ZZ BASE AS-EV-ZZMUA-L35 4.2 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH7633L35	25	25		10-12
SET ZZ BASE AS-EV-ZZMUA-L35 4.8 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH7634L35	25	25		10-12
SET ZZ BASE AS-EV-ZZMUA-L35 5.4 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH7635L35	25	25		10-12
SET ZZ BASE AS-EV-ZZMUA-L35 Profile4.2 (Zirkonzahn Multi Unit Abutment non hex L35 and Sc	BSAH7637L35	25	25		10-12
SET ZZ BASE AS-EV-ZZMUA-L35 Profile4.8 (Zirkonzahn Multi Unit Abutment non hex L35 and Sc	BSAH7638L35	25	25		10-12
SET ZZ BASE AS-EV-ZZMUA-L45 3.0 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH7631L45	25	25		10-12
SET ZZ BASE AS-EV-ZZMUA-L45 3.6 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH7632L45	25	25		10-12
SET ZZ BASE AS-EV-ZZMUA-L45 4.2 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH7633L45	25	25		10-12
SET ZZ BASE AS-EV-ZZMUA-L45 4.8 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH7634L45	25	25		10-12
SET ZZ BASE AS-EV-ZZMUA-L45 5.4 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH7635L45	25	25		10-12
SET ZZ BASE AS-EV-ZZMUA-L45 Profile4.2 (Zirkonzahn Multi Unit Abutment non hex L45 and Sc	BSAH7637L45	25	25		10-12
SET ZZ BASE AS-EV-ZZMUA-L45 Profile4.8 (Zirkonzahn Multi Unit Abutment non hex L45 and Sc	BSAH7638L45	25	25		10-12
ZZ BASE AS-EV-ASCSC 3.0 (Angled screw-channel Screw)	BSKE7601	20			10-12
ZZ BASE AS-EV-ASCSC 3.6 (Angled screw-channel Screw)	BSKE7602	20			10-12
ZZ BASE AS-EV-ASCSC 4.2 (Angled screw-channel Screw)	BSKE7603	20			10-12
ZZ BASE AS-EV-ASCSC 4.8 (Angled screw-channel Screw)	BSKE7604	20			10-12
ZZ BASE AS-EV-ASCSC Profile4.2 (Angled screw-channel Screw)	BSKE7607	20			10-12
ZZ BASE AS-EV-ASCSC UA (Angled screw-channel Screw)	BSKE7606	15			10-12
ZZ BASE AS-EV-BSC 3.0 (Black Screw)	BSIE7601	25			10-12
ZZ BASE AS-EV-BSC 3.6 (Black Screw)	BSIE7602	25			10-12
ZZ BASE AS-EV-BSC 4.2 (Black Screw)	BSIE7603	25			10-12
ZZ BASE AS-EV-BSC 4.8 (Black Screw)	BSIE7604	25			10-12
ZZ BASE AS-EV-BSC 5.4 (Black Screw)	BSIE7605	25			10-12
ZZ BASE AS-EV-BSC Profile4.2 (Black Screw)	BSIE7607	25			10-12
ZZ BASE AS-EV-BSC Profile4.8 (Black Screw)	BSIE7608	25			10-12
ZZ BASE AS-EV-BSC UA (Black Screw)	BSIE7606	15			10-12

ASTRA TECH IMPLANT SYSTEM EV/UNIABUTMENT® EV

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE AS-EV-BZSC 3.0 (Black Zirconia screw)	BSIE7631	25			10-12
ZZ BASE AS-EV-BZSC 3.6 (Black Zirconia screw)	BSIE7632	25			10-12
ZZ BASE AS-EV-BZSC 4.2 (Black Zirconia screw)	BSIE7633	25			10-12
ZZ BASE AS-EV-BZSC 4.8 (Black Zirconia screw)	BSIE7634	25			10-12
ZZ BASE AS-EV-BZSC 5.4 (Black Zirconia screw)	BSIE7635	25			10-12
ZZ BASE AS-EV-BZSC Profile4.2 (Black Zirconia screw)	BSIE7637	25			10-12
ZZ BASE AS-EV-BZSC Profile4.8 (Black Zirconia screw)	BSIE7638	25			10-12
ZZ BASE AS-EV-BZSC UA (Black Zirconia screw)	BSIE7636	15			10-12
ZZ BASE AS-EV-ICSC 3.0 (Impression coping Screw)	BSFE7601	15			10-12
ZZ BASE AS-EV-ICSC 3.6 (Impression coping Screw)	BSFE7602	15			10-12
ZZ BASE AS-EV-ICSC 4.2 (Impression coping Screw)	BSFE7603	15			10-12
ZZ BASE AS-EV-ICSC 4.8 (Impression coping Screw)	BSFE7604	15			10-12
ZZ BASE AS-EV-ICSC 5.4 (Impression coping Screw)	BSFE7605	15			10-12
ZZ BASE AS-EV-ICSC Profile4.2 (Impression coping Screw)	BSFE7607	15			10-12
ZZ BASE AS-EV-ICSC Profile4.8 (Impression coping Screw)	BSFE7608	15			10-12
ZZ BASE AS-EV-ICSC UA (Impression coping Screw)	BSFE7606	15			10-12
ZZ BASE AS-EV-LBSC 3.0 (Lab Screw)	BSHE7601	25			10-12
ZZ BASE AS-EV-LBSC 3.6 (Lab Screw)	BSHE7602	25			10-12
ZZ BASE AS-EV-LBSC 4.2 (Lab Screw)	BSHE7603	25			10-12
ZZ BASE AS-EV-LBSC 4.8 (Lab Screw)	BSHE7604	25			10-12
ZZ BASE AS-EV-LBSC 5.4 (Lab Screw)	BSHE7605	25			10-12
ZZ BASE AS-EV-LBSC Profile4.2 (Lab Screw)	BSHE7607	25			10-12
ZZ BASE AS-EV-LBSC Profile4.8 (Lab Screw)	BSHE7608	25			10-12
ZZ BASE AS-EV-LBSC UA (Lab Screw)	BSHE7606	15			10-12
ZZ BASE AS-EV-LBZSC 3.0 (Lab Zirconia screw)	BSHE7631	25			10-12
ZZ BASE AS-EV-LBZSC 3.6 (Lab Zirconia screw)	BSHE7632	25			10-12
ZZ BASE AS-EV-LBZSC 4.2 (Lab Zirconia screw)	BSHE7633	25			10-12
ZZ BASE AS-EV-LBZSC 4.8 (Lab Zirconia screw)	BSHE7634	25			10-12
ZZ BASE AS-EV-LBZSC 5.4 (Lab Zirconia screw)	BSHE7635	25			10-12

ASTRA TECH IMPLANT SYSTEM EV/UNIABUTMENT® EV

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
ZZ BASE AS-EV-LBZSC Profile4.2 (Lab Zirconia screw)	BSHE7637	25			10-12
ZZ BASE AS-EV-LBZSC Profile4.8 (Lab Zirconia screw)	BSHE7638	25			10-12
ZZ BASE AS-EV-LBZSC UA (Lab Zirconia screw)	BSHE7636	15			10-12
ZZ BASE AS-EV-MI 3.0 (Model implant)	MIAG7601				10-12
ZZ BASE AS-EV-MI 3.6 (Model implant)	MIAG7602				10-12
ZZ BASE AS-EV-MI 4.2 (Model implant)	MIAG7603				10-12
ZZ BASE AS-EV-MI 4.8 (Model implant)	MIAG7604				10-12
ZZ BASE AS-EV-MI 5.4 (Model implant)	MIAG7605				10-12
ZZ BASE AS-EV-MI Profile4.2 (Model implant)	MIAG7607				10-12
ZZ BASE AS-EV-MI Profile4.8 (Model implant)	MIAG7608				10-12
ZZ BASE AS-EV-MI UA (Model implant)	MIAG7606				10-12
ZZ BASE AS-EV-SC 3.0 (Screw)	BSAE7601	25			10-12
ZZ BASE AS-EV-SC 3.6 (Screw)	BSAE7602	25			10-12
ZZ BASE AS-EV-SC 4.2 (Screw)	BSAE7603	25			10-12
ZZ BASE AS-EV-SC 4.8 (Screw)	BSAE7604	25			10-12
ZZ BASE AS-EV-SC 5.4 (Screw)	BSAE7605	25			10-12
ZZ BASE AS-EV-SC Profile4.2 (Screw)	BSAE7607	25			10-12
ZZ BASE AS-EV-SC Profile4.8 (Screw)	BSAE7608	25			10-12
ZZ BASE AS-EV-SC UA (Screw)	BSAE7606	15			10-12
ZZ BASE AS-EV-SC-G 3.0 (Screw Gold)	BSCE7601	25			10-12
ZZ BASE AS-EV-SC-G 3.6 (Screw Gold)	BSCE7602	25			10-12
ZZ BASE AS-EV-SC-G 4.8 (Screw Gold)	BSCE7604	25			10-12
ZZ BASE AS-EV-SC-G 5.4 (Screw Gold)	BSCE7605	25			10-12
ZZ BASE AS-EV-SC-G Profile4.8 (Screw Gold)	BSCE7608	25			10-12
ZZ BASE AS-EV-SC-G UA (Screw Gold)	BSCE7606	15			10-12
ZZ BASE AS-EV-WMSM-L10 3.0 (White Metal Scanmarker L10 with Screw)	ZBED7631L10	15			10-12
ZZ BASE AS-EV-WMSM-L10 3.6 (White Metal Scanmarker L10 with Screw)	ZBED7632L10	15			10-12
ZZ BASE AS-EV-WMSM-L10 4.2 (White Metal Scanmarker L10 with Screw)	ZBED7633L10	15			10-12
ZZ BASE AS-EV-WMSM-L10 4.8 (White Metal Scanmarker L10 with Screw)	ZBED7634L10	15			10-12

All information is subject to change. Errors and omissions excepted. Version: 21.09.2023

ASTRA TECH IMPLANT SYSTEM EV/UNIABUTMENT® EV

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE AS-EV-WMSM-L10 5.4 (White Metal Scanmarker L10 with Screw)	ZBED7635L10	15			10-12
ZZ BASE AS-EV-WMSM-L10 Profile4.2 (White Metal Scanmarker L10 with Screw)	ZBED7637L10	15			10-12
ZZ BASE AS-EV-WMSM-L10 Profile4.8 (White Metal Scanmarker L10 with Screw)	ZBED7638L10	15			10-12
ZZ BASE AS-EV-WMSM-L10 UA (White Metal Scanmarker L10 with Screw)	ZBED7636L10	15			10-12
ZZ BASE AS-EV-WMSM-L12 3.0 (White Metal Scanmarker L12 with Screw)	ZBED7631L12	15			10-12
ZZ BASE AS-EV-WMSM-L12 3.6 (White Metal Scanmarker L12 with Screw)	ZBED7632L12	15			10-12
ZZ BASE AS-EV-WMSM-L12 4.2 (White Metal Scanmarker L12 with Screw)	ZBED7633L12	15			10-12
ZZ BASE AS-EV-WMSM-L12 4.8 (White Metal Scanmarker L12 with Screw)	ZBED7634L12	15			10-12
ZZ BASE AS-EV-WMSM-L12 5.4 (White Metal Scanmarker L12 with Screw)	ZBED7635L12	15			10-12
ZZ BASE AS-EV-WMSM-L12 Profile4.2 (White Metal Scanmarker L12 with Screw)	ZBED7637L12	15			10-12
ZZ BASE AS-EV-WMSM-L12 Profile4.8 (White Metal Scanmarker L12 with Screw)	ZBED7638L12	15			10-12
ZZ BASE AS-EV-WMSM-L12 UA (White Metal Scanmarker L12 with Screw)	ZBED7636L12	15			10-12
ZZ BASE AS-EV-WMSM-L15 3.0 (White Metal Scanmarker L15 with Screw)	ZBED7631L15	15			10-12
ZZ BASE AS-EV-WMSM-L15 3.6 (White Metal Scanmarker L15 with Screw)	ZBED7632L15	15			10-12
ZZ BASE AS-EV-WMSM-L15 4.2 (White Metal Scanmarker L15 with Screw)	ZBED7633L15	15			10-12
ZZ BASE AS-EV-WMSM-L15 4.8 (White Metal Scanmarker L15 with Screw)	ZBED7634L15	15			10-12
ZZ BASE AS-EV-WMSM-L15 5.4 (White Metal Scanmarker L15 with Screw)	ZBED7635L15	15			10-12
ZZ BASE AS-EV-WMSM-L15 Profile4.2 (White Metal Scanmarker L15 with Screw)	ZBED7637L15	15			10-12
ZZ BASE AS-EV-WMSM-L15 Profile4.8 (White Metal Scanmarker L15 with Screw)	ZBED7638L15	15			10-12
ZZ BASE AS-EV-WMSM-L15 UA (White Metal Scanmarker L15 with Screw)	ZBED7636L15	15			10-12
ZZ BASE AS-EV-ZSC 3.0 (Zirconia Screw)	BSAE7631	25			10-12
ZZ BASE AS-EV-ZSC 3.6 (Zirconia Screw)	BSAE7632	25			10-12
ZZ BASE AS-EV-ZSC 4.2 (Zirconia Screw)	BSAE7633	25			10-12
ZZ BASE AS-EV-ZSC 4.8 (Zirconia Screw)	BSAE7634	25			10-12
ZZ BASE AS-EV-ZSC 5.4 (Zirconia Screw)	BSAE7635	25			10-12
ZZ BASE AS-EV-ZSC Profile4.2 (Zirconia Screw)	BSAE7637	25			10-12
ZZ BASE AS-EV-ZSC Profile4.8 (Zirconia Screw)	BSAE7638	25			10-12
ZZ BASE AS-EV-ZSC UA (Zirconia Screw)	BSAE7636	15			10-12
ZZ BASE AS-EV-ZSC-G 3.0 (Zirconia Screw Gold)	BSCE7631	25			10-12

ASTRA TECH IMPLANT SYSTEM EV/UNIABUTMENT® EV

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE AS-EV-ZSC-G 3.6 (Zirconia Screw Gold)	BSCE7632	25			10-12
ZZ BASE AS-EV-ZSC-G 4.2 (Zirconia Screw Gold)	BSCE7633	25			10-12
ZZ BASE AS-EV-ZSC-G 4.8 (Zirconia Screw Gold)	BSCE7634	25			10-12
ZZ BASE AS-EV-ZSC-G 5.4 (Zirconia Screw Gold)	BSCE7635	25			10-12
ZZ BASE AS-EV-ZSC-G Profile4.2 (Zirconia Screw Gold)	BSCE7637	25			10-12
ZZ BASE AS-EV-ZSC-G Profile4.8 (Zirconia Screw Gold)	BSCE7638	25			10-12
ZZ BASE AS-EV-ZSC-G UA (Zirconia Screw Gold)	BSCE7636	15			10-12
ZZ BASE AS-EV-ZZMUASC 3.0 (Zirkonzahn Multi Unit Abutment Screw)	BSEE7601	15			10-12
ZZ BASE AS-EV-ZZMUASC 3.6 (Zirkonzahn Multi Unit Abutment Screw)	BSEE7602	20			10-12
ZZ BASE AS-EV-ZZMUASC 4.2 (Zirkonzahn Multi Unit Abutment Screw)	BSEE7603	25			10-12
ZZ BASE AS-EV-ZZMUASC 4.8 (Zirkonzahn Multi Unit Abutment Screw)	BSEE7604	30			10-12
ZZ BASE AS-EV-ZZMUASC 5.4 (Zirkonzahn Multi Unit Abutment Screw)	BSEE7605	30			10-12

ASTRA TECH IMPLANT SYSTEM OSSEOSPEED™

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE AS-OS-ABCH 3.0 (Raw-Abutment® customizable hex and screw)	BSAF0181	15			10-12
SET ZZ BASE AS-OS-ABCH 3.5/4.0 (Raw-Abutment® customizable hex and screw)	BSAF0182	20			10-12
SET ZZ BASE AS-OS-ABCH 4.5/5.0 (Raw-Abutment® customizable hex and screw)	BSAF0183	25			10-12
SET ZZ BASE AS-OS-ABCH-D14 3.0 (Raw-Abutment® customizable hex D14 and screw)	BSAF0181D14	15			10-12
SET ZZ BASE AS-OS-ABCH-D14 3.5/4.0 (Raw-Abutment® customizable hex D14 and screw)	BSAF0182D14	20			10-12
SET ZZ BASE AS-OS-ABCH-D14 4.5/5.0 (Raw-Abutment® customizable hex D14 and screw)	BSAF0183D14	25			10-12
SET ZZ BASE AS-OS-ASCTBH-K80 3.0 (Angled screw-channel titan base hex K80 and screw)	BSAT0181K80	15			10-12
SET ZZ BASE AS-OS-ASCTBH-K80 3.5/4.0 (Angled screw-channel titan base hex K80 and screw)	BSAT0182K80	20			10-12
SET ZZ BASE AS-OS-ASCTBH-K80 4.5/5.0 (Angled screw-channel titan base hex K80 and screw)	BSAT0183K80	20			10-12
SET ZZ BASE AS-OS-ASCTB-K80 3.0 (Angled screw-channel titan base non hex K80 and screw)	BSAT0131K80	15			10-12
SET ZZ BASE AS-OS-ASCTB-K80 3.5/4.0 (Angled screw-channel titan base non hex K80 and screw)	BSAT0132K80	20			10-12
SET ZZ BASE AS-OS-ASCTB-K80 4.5/5.0 (Angled screw-channel titan base non hex K80 and screw)	BSAT0133K80	20			10-12
SET ZZ BASE AS-OS-CCTB 3.0 (Conical cemented titan base non hex and screw)	BSAD0131	15			10-12
SET ZZ BASE AS-OS-CCTB 3.5/4.0 (Conical cemented titan base non hex and screw)	BSAD0132	20			10-12
SET ZZ BASE AS-OS-CCTB 4.5/5.0 (Conical cemented titan base non hex and screw)	BSAD0133	25			10-12
SET ZZ BASE AS-OS-CCTB-G 3.0 (Conical cemented titan base non hex and screw Gold)	BSCD0131	15			10-12
SET ZZ BASE AS-OS-CCTB-G 3.5/4.0 (Conical cemented titan base non hex and screw Gold)	BSCD0132	20			10-12
SET ZZ BASE AS-OS-CCTB-G 4.5/5.0 (Conical cemented titan base non hex and screw Gold)	BSCD0133	25			10-12
SET ZZ BASE AS-OS-CCTB-K85 3.0 (Conical cemented titan base non hex K85 and screw)	BSAD0131K85	15			10-12
SET ZZ BASE AS-OS-CCTB-K85 3.5/4.0 (Conical cemented titan base non hex K85 and screw)	BSAD0132K85	20			10-12
SET ZZ BASE AS-OS-CCTB-K85 4.5/5.0 (Conical cemented titan base non hex K85 and screw)	BSAD0133K85	25			10-12
SET ZZ BASE AS-OS-CCTB-L10 3.0 (Conical cemented titan base non hex L10 and screw)	BSAD0131L10	15			10-12
SET ZZ BASE AS-OS-CCTB-L10 3.5/4.0 (Conical cemented titan base non hex L10 and screw)	BSAD0132L10	20			10-12
SET ZZ BASE AS-OS-CCTB-L10 4.5/5.0 (Conical cemented titan base non hex L10 and screw)	BSAD0133L10	25			10-12
SET ZZ BASE AS-OS-CCTB-L10-G 3.0 (Conical cemented titan base non hex L10 and screw Gold)	BSCD0131L10	15			10-12
SET ZZ BASE AS-OS-CCTB-L10-G 3.5/4.0 (Conical cemented titan base non hex L10 and screw Gold)	BSCD0132L10	20			10-12
SET ZZ BASE AS-OS-CCTB-L10-G 4.5/5.0 (Conical cemented titan base non hex L10 and screw Gold)	BSCD0133L10	25			10-12
SET ZZ BASE AS-OS-CCTB-L20 3.0 (Conical cemented titan base non hex L20 and screw)	BSAD0131L20	15			10-12
SET ZZ BASE AS-OS-CCTB-L20 3.5/4.0 (Conical cemented titan base non hex L20 and screw)	BSAD0132L20	20			10-12

ASTRA TECH IMPLANT SYSTEM OSSEOSPEED™

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
SET ZZ BASE AS-OS-CCTB-L20 4.5/5.0 (Conical cemented titan base non hex L20 and screw)	BSAD0133L20	25			10-12
SET ZZ BASE AS-OS-CCTB-L20-G 3.0 (Conical cemented titan base non hex L20 and screw Gold)	BSCD0131L20	15			10-12
SET ZZ BASE AS-OS-CCTB-L20-G 3.5/4.0 (Conical cemented titan base non hex L20 and screw Gold)	BSCD0132L20	20			10-12
SET ZZ BASE AS-OS-CCTB-L20-G 4.5/5.0 (Conical cemented titan base non hex L20 and screw Gold)	BSCD0133L20	25			10-12
SET ZZ BASE AS-OS-CCTB-L30 3.0 (Conical cemented titan base non hex L30 and screw)	BSAD0131L30	15			10-12
SET ZZ BASE AS-OS-CCTB-L30 3.5/4.0 (Conical cemented titan base non hex L30 and screw)	BSAD0132L30	20			10-12
SET ZZ BASE AS-OS-CCTB-L30 4.5/5.0 (Conical cemented titan base non hex L30 and screw)	BSAD0133L30	25			10-12
SET ZZ BASE AS-OS-CCTB-L30-G 3.0 (Conical cemented titan base non hex L30 and screw Gold)	BSCD0131L30	15			10-12
SET ZZ BASE AS-OS-CCTB-L30-G 3.5/4.0 (Conical cemented titan base non hex L30 and screw Gold)	BSCD0132L30	20			10-12
SET ZZ BASE AS-OS-CCTB-L30-G 4.5/5.0 (Conical cemented titan base non hex L30 and screw Gold)	BSCD0133L30	25			10-12
SET ZZ BASE AS-OS-CCTB-L40 3.5/4.0 (Conical cemented titan base non hex L40 and screw)	BSAD0132L40	20			10-12
SET ZZ BASE AS-OS-CCTB-L40 4.5/5.0 (Conical cemented titan base non hex L40 and screw)	BSAD0133L40	25			10-12
SET ZZ BASE AS-OS-CCTB-L40-G 3.5/4.0 (Conical cemented titan base non hex L40 and screw Gold)	BSCD0132L40	20			10-12
SET ZZ BASE AS-OS-CCTB-L40-G 4.5/5.0 (Conical cemented titan base non hex L40 and screw Gold)	BSCD0133L40	25			10-12
SET ZZ BASE AS-OS-ICOT 3.0 (Impression coping open tray non hex and Screw and Driver)	BSAQ0161	15			10-12
SET ZZ BASE AS-OS-ICOT 3.5/4.0 (Impression coping open tray non hex and Screw and Driver)	BSAQ0162	15			10-12
SET ZZ BASE AS-OS-ICOT 4.5/5.0 (Impression coping open tray non hex and Screw and Driver)	BSAQ0163	15			10-12
SET ZZ BASE AS-OS-ICOTH 3.0 (Impression coping open tray hex and Screw and Driver)	BSAQ0141	15			10-12
SET ZZ BASE AS-OS-ICOTH 3.5/4.0 (Impression coping open tray hex and Screw and Driver)	BSAQ0142	15			10-12
SET ZZ BASE AS-OS-ICOTH 4.5/5.0 (Impression coping open tray hex and Screw and Driver)	BSAQ0143	15			10-12
SET ZZ BASE AS-OS-NTB 3.0 (Narrow titan base non hex and screw)	BSAP0131	15			10-12
SET ZZ BASE AS-OS-NTB 3.5/4.0 (Narrow titan base non hex and screw)	BSAP0132	20			10-12
SET ZZ BASE AS-OS-NTBH 3.0 (Narrow titan base hex and screw)	BSAP0181	15			10-12
SET ZZ BASE AS-OS-NTBH 3.5/4.0 (Narrow titan base hex and screw)	BSAP0182	20			10-12
SET ZZ BASE AS-OS-PCTBH 3.0 (Parallel cemented titan base hex and screw)	BSAB0181	15			10-12
SET ZZ BASE AS-OS-PCTBH 3.0 (Parallel cemented titan base hex and screw)	BSAB0191	15			10-12
SET ZZ BASE AS-OS-PCTBH 3.5/4.0 (Parallel cemented titan base hex and screw)	BSAB0192	20			10-12
SET ZZ BASE AS-OS-PCTBH 4.5/5.0 (Parallel cemented titan base hex and screw)	BSAB0183	25			10-12
SET ZZ BASE AS-OS-PCTBH 4.5/5.0 (Parallel cemented titan base hex and screw)	BSAB0193	25			10-12

ASTRA TECH IMPLANT SYSTEM OSSEOSPEED™

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE AS-OS-PCTBH-G 3.0 (Parallel cemented titan base hex and screw Gold)	BSCB0181	15			10-12
SET ZZ BASE AS-OS-PCTBH-G 3.5/4.0 (Parallel cemented titan base hex and screw Gold)	BSCB0192	20			10-12
SET ZZ BASE AS-OS-PCTBH-G 4.5/5.0 (Parallel cemented titan base hex and screw Gold)	BSCB0193	25			10-12
SET ZZ BASE AS-OS-PCTBH-K85 3.0 (Parallel cemented titan base hex K85 and screw)	BSAB0181K85	15			10-12
SET ZZ BASE AS-OS-PCTBH-K85 3.5/4.0 (Parallel cemented titan base hex K85 and screw)	BSAB0182K85	20			10-12
SET ZZ BASE AS-OS-PCTBH-K85 4.5/5.0 (Parallel cemented titan base hex K85 and screw)	BSAB0183K85	25			10-12
SET ZZ BASE AS-OS-PCTBH-L20 3.0 (Parallel cemented titan base hex L20 and screw)	BSAB0181L20	15			10-12
SET ZZ BASE AS-OS-PCTBH-L20 3.0 (Parallel cemented titan base hex L20 and screw)	BSAB0191L20	15			10-12
SET ZZ BASE AS-OS-PCTBH-L20 3.5/4.0 (Parallel cemented titan base hex L20 and screw)	BSAB0182L20	20			10-12
SET ZZ BASE AS-OS-PCTBH-L20 4.5/5.0 (Parallel cemented titan base hex L20 and screw)	BSAB0183L20	25			10-12
SET ZZ BASE AS-OS-PCTBH-L20 4.5/5.0 (Parallel cemented titan base hex L20 and screw)	BSAB0193L20	25			10-12
SET ZZ BASE AS-OS-PCTBH-L20-G 3.0 (Parallel cemented titan base hex L20 and screw Gold)	BSCB0181L20	15			10-12
SET ZZ BASE AS-OS-PCTBH-L20-G 3.5/4.0 (Parallel cemented titan base hex L20 and screw Gol)	BSCB0182L20	20			10-12
SET ZZ BASE AS-OS-PCTBH-L20-G 4.5/5.0 (Parallel cemented titan base hex L20 and screw Gol)	BSCB0183L20	25			10-12
SET ZZ BASE AS-OS-PCTBH-L30 3.0 (Parallel cemented titan base hex L30 and screw)	BSAB0181L30	15			10-12
SET ZZ BASE AS-OS-PCTBH-L30 3.0 (Parallel cemented titan base hex L30 and screw)	BSAB0191L30	15			10-12
SET ZZ BASE AS-OS-PCTBH-L30 3.5/4.0 (Parallel cemented titan base hex L30 and screw)	BSAB0182L30	20			10-12
SET ZZ BASE AS-OS-PCTBH-L30 4.5/5.0 (Parallel cemented titan base hex L30 and screw)	BSAB0183L30	25			10-12
SET ZZ BASE AS-OS-PCTBH-L30 4.5/5.0 (Parallel cemented titan base hex L30 and screw)	BSAB0193L30	25			10-12
SET ZZ BASE AS-OS-PCTBH-L30-G 3.0 (Parallel cemented titan base hex L30 and screw Gold)	BSCB0181L30	15			10-12
SET ZZ BASE AS-OS-PCTBH-L30-G 3.5/4.0 (Parallel cemented titan base hex L30 and screw Gol)	BSCB0182L30	20			10-12
SET ZZ BASE AS-OS-PCTBH-L30-G 4.5/5.0 (Parallel cemented titan base hex L30 and screw Gol)	BSCB0183L30	25			10-12
SET ZZ BASE AS-OS-PCTBH-L40 3.5/4.0 (Parallel cemented titan base hex L40 and screw)	BSAB0182L40	20			10-12
SET ZZ BASE AS-OS-PCTBH-L40 4.5/5.0 (Parallel cemented titan base hex L40 and screw)	BSAB0183L40	25			10-12
SET ZZ BASE AS-OS-PCTBH-L40 4.5/5.0 (Parallel cemented titan base hex L40 and screw)	BSAB0193L40	25			10-12
SET ZZ BASE AS-OS-PCTBH-L40-G 3.5/4.0 (Parallel cemented titan base hex L40 and screw Gol)	BSCB0182L40	20			10-12
SET ZZ BASE AS-OS-PCTBH-L40-G 4.5/5.0 (Parallel cemented titan base hex L40 and screw Gol)	BSCB0183L40	25			10-12
SET ZZ BASE AS-OS-SA 3.0 (Scananalog and screw)	ZBDE0131	15			10-12
SET ZZ BASE AS-OS-SA 3.5/4.0 (Scananalog and screw)	ZBDE0132	15			10-12

ASTRA TECH IMPLANT SYSTEM OSSEOSPEED™

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE AS-OS-SA 4.5/5.0 (Scananalog and screw)	ZBDE0133	15			10-12
SET ZZ BASE AS-OS-SM 3.0 (Scanmarker and screw)	ZBAD0131	15			10-12
SET ZZ BASE AS-OS-SM 3.5/4.0 (Scanmarker and screw)	ZBAD0132	20			10-12
SET ZZ BASE AS-OS-SM 4.5/5.0 (Scanmarker and screw)	ZBAD0133	25			10-12
SET ZZ BASE AS-OS-WSM 3.0 (White Scanmarker and screw)	ZBCD0131	15			10-12
SET ZZ BASE AS-OS-WSM 3.5/4.0 (White Scanmarker and screw)	ZBCD0132	20			10-12
SET ZZ BASE AS-OS-WSM 4.5/5.0 (White Scanmarker and screw)	ZBCD0133	25			10-12
SET ZZ BASE AS-OS-ZZMUA17 3.0 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and impli	BSAI0131	15			10-12
SET ZZ BASE AS-OS-ZZMUA17 3.5/4.0 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and i	BSAI0132	20			10-12
SET ZZ BASE AS-OS-ZZMUA17 4.5/5.0 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and i	BSAI0133	30			10-12
SET ZZ BASE AS-OS-ZZMUA17-L10 3.0 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10	BSAI0131L10	15			10-12
SET ZZ BASE AS-OS-ZZMUA17-L10 3.5/4.0 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L	BSAI0132L10	20			10-12
SET ZZ BASE AS-OS-ZZMUA17-L10 4.5/5.0 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L	BSAI0133L10	30			10-12
SET ZZ BASE AS-OS-ZZMUA-L15 3.0 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH0131L15	25	15		10-12
SET ZZ BASE AS-OS-ZZMUA-L15 3.5/4.0 (Zirkonzahn Multi Unit Abutment non hex L15 and Scre	BSAH0132L15	25	20		10-12
SET ZZ BASE AS-OS-ZZMUA-L15 4.5/5.0 (Zirkonzahn Multi Unit Abutment non hex L15 and Scre	BSAH0133L15	25	25		10-12
SET ZZ BASE AS-OS-ZZMUA-L25 3.0 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH0131L25	25	15		10-12
SET ZZ BASE AS-OS-ZZMUA-L25 3.5/4.0 (Zirkonzahn Multi Unit Abutment non hex L25 and Scre	BSAH0132L25	25	20		10-12
SET ZZ BASE AS-OS-ZZMUA-L25 4.5/5.0 (Zirkonzahn Multi Unit Abutment non hex L25 and Scre	BSAH0133L25	25	25		10-12
SET ZZ BASE AS-OS-ZZMUA-L35 3.0 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH0131L35	25	15		10-12
SET ZZ BASE AS-OS-ZZMUA-L35 3.5/4.0 (Zirkonzahn Multi Unit Abutment non hex L35 and Scre	BSAH0132L35	25	20		10-12
SET ZZ BASE AS-OS-ZZMUA-L35 4.5/5.0 (Zirkonzahn Multi Unit Abutment non hex L35 and Scre	BSAH0133L35	25	25		10-12
SET ZZ BASE AS-OS-ZZMUA-L45 3.0 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH0131L45	25	15		10-12
SET ZZ BASE AS-OS-ZZMUA-L45 3.5/4.0 (Zirkonzahn Multi Unit Abutment non hex L45 and Scre	BSAH0132L45	25	20		10-12
SET ZZ BASE AS-OS-ZZMUA-L45 4.5/5.0 (Zirkonzahn Multi Unit Abutment non hex L45 and Scre	BSAH0133L45	25	25		10-12
ZZ BASE AS-OS-ASCSC 3.0 (Angled screw-channel Screw)	BSKE0101	15			10-12
ZZ BASE AS-OS-ASCSC 3.5/4.0 (Angled screw-channel Screw)	BSKE0102	20			10-12
ZZ BASE AS-OS-ASCSC 4.5/5.0 (Angled screw-channel Screw)	BSKE0103	20			10-12
ZZ BASE AS-OS-BSC 3.0 (Black Screw)	BSIE0101	15			10-12

ASTRA TECH IMPLANT SYSTEM OSSEOSPEED™

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
ZZ BASE AS-OS-BSC 3.5/4.0 (Black Screw)	BSIE0102	20			10-12
ZZ BASE AS-OS-BSC 4.5/5.0 (Black Screw)	BSIE0103	25			10-12
ZZ BASE AS-OS-BZSC 3.5/4.0 (Black Zirconia screw)	BSIE0132	20			10-12
ZZ BASE AS-OS-BZSC 4.5/5.0 (Black Zirconia screw)	BSIE0133	25			10-12
ZZ BASE AS-OS-ICSC 3.0 (Impression coping Screw)	BSFE0101	15			10-12
ZZ BASE AS-OS-ICSC 3.5/4.0 (Impression coping Screw)	BSFE0102	15			10-12
ZZ BASE AS-OS-ICSC 4.5/5.0 (Impression coping Screw)	BSFE0103	15			10-12
ZZ BASE AS-OS-LBSC 3.0 (Lab Screw)	BSHE0101	15			10-12
ZZ BASE AS-OS-LBSC 3.5/4.0 (Lab Screw)	BSHE0102	20			10-12
ZZ BASE AS-OS-LBSC 4.5/5.0 (Lab Screw)	BSHE0103	25			10-12
ZZ BASE AS-OS-LBZSC 3.5/4.0 (Lab Zirconia screw)	BSHE0132	20			10-12
ZZ BASE AS-OS-LBZSC 4.5/5.0 (Lab Zirconia screw)	BSHE0133	25			10-12
ZZ BASE AS-OS-MI 3.0 (Model implant)	MIAG0101				10-12
ZZ BASE AS-OS-MI 3.5/4.0 (Model implant)	MIAG0102				10-12
ZZ BASE AS-OS-MI 4.5/5.0 (Model implant)	MIAG0103				10-12
ZZ BASE AS-OS-SC 3.0 (Screw)	BSAE0101	15			10-12
ZZ BASE AS-OS-SC 3.5/4.0 (Screw)	BSAE0102	20			10-12
ZZ BASE AS-OS-SC 4.5/5.0 (Screw)	BSAE0103	25			10-12
ZZ BASE AS-OS-SC-G 3.0 (Screw Gold)	BSCE0101	15			10-12
ZZ BASE AS-OS-SC-G 3.5/4.0 (Screw Gold)	BSCE0102	20			10-12
ZZ BASE AS-OS-SC-G 4.5/5.0 (Screw Gold)	BSCE0103	25			10-12
ZZ BASE AS-OS-WMSM-L10 3.0 (White Metal Scanmarker L10 with Screw)	ZBED0131L10	15			10-12
ZZ BASE AS-OS-WMSM-L10 3.5/4.0 (White Metal Scanmarker L10 with Screw)	ZBED0132L10	15			10-12
ZZ BASE AS-OS-WMSM-L10 4.5/5.0 (White Metal Scanmarker L10 with Screw)	ZBED0133L10	15			10-12
ZZ BASE AS-OS-WMSM-L12 3.0 (White Metal Scanmarker L12 with Screw)	ZBED0131L12	15			10-12
ZZ BASE AS-OS-WMSM-L12 3.5/4.0 (White Metal Scanmarker L12 with Screw)	ZBED0132L12	15			10-12
ZZ BASE AS-OS-WMSM-L12 4.5/5.0 (White Metal Scanmarker L12 with Screw)	ZBED0133L12	15			10-12
ZZ BASE AS-OS-WMSM-L15 3.0 (White Metal Scanmarker L15 with Screw)	ZBED0131L15	15			10-12
ZZ BASE AS-OS-WMSM-L15 3.5/4.0 (White Metal Scanmarker L15 with Screw)	ZBED0132L15	15			10-12

ASTRA TECH IMPLANT SYSTEM OSSEOSPEED™

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE AS-OS-WMSM-L15 4.5/5.0 (White Metal Scanmarker L15 with Screw)	ZBED0133L15	15			10-12
ZZ BASE AS-OS-ZSC 3.5/4.0 (Zirconia screw)	BSAE0132	20			10-12
ZZ BASE AS-OS-ZSC 4.5/5.0 (Zirconia screw)	BSAE0133	25			10-12
ZZ BASE AS-OS-ZSC-G 3.5/4.0 (Zirconia screw Gold)	BSCE0132	20			10-12
ZZ BASE AS-OS-ZSC-G 4.5/5.0 (Zirconia screw Gold)	BSCE0133	25			10-12
ZZ BASE AS-OS-ZZMUASC 3.0 (Zirkonzahn Multi Unit Abutment Screw)	BSEE0101	15			10-12
ZZ BASE AS-OS-ZZMUASC 3.5/4.0 (Zirkonzahn Multi Unit Abutment Screw)	BSEE0102	20			10-12
ZZ BASE AS-OS-ZZMUASC 4.5/5.0 (Zirkonzahn Multi Unit Abutment Screw)	BSEE0103	30			10-12

ASTRA TECH UNIABUTMENT®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE AS-UA-CCTB 20° (Conical cemented titan base non hex and screw)	BSAD0731	15			10-12
SET ZZ BASE AS-UA-CCTB 45° (Conical cemented titan base non hex and screw)	BSAD0732	15			10-12
SET ZZ BASE AS-UA-CCTB-G 20° (Conical cemented titan base non hex and screw Gold)	BSCD0731	15			10-12
SET ZZ BASE AS-UA-CCTB-G 45° (Conical cemented titan base non hex and screw Gold)	BSCD0732	15			10-12
SET ZZ BASE AS-UA-CCTB-K85 20° (Conical cemented titan base non hex K85 and screw)	BSAD0731K85	15			10-12
SET ZZ BASE AS-UA-CCTB-K85 45° (Conical cemented titan base non hex K85 and screw)	BSAD0732K85	15			10-12
SET ZZ BASE AS-UA-ICOT 20° (Impression coping open tray non hex and Screw and Driver)	BSAQ0761	15			10-12
SET ZZ BASE AS-UA-ICOT 45° (Impression coping open tray non hex and Screw and Driver)	BSAQ0762	15			10-12
SET ZZ BASE AS-UA-SA 20° (Scananalog and screw)	ZBDE0731	15			10-12
SET ZZ BASE AS-UA-SA 45° (Scananalog and screw)	ZBDE0732	15			10-12
SET ZZ BASE AS-UA-SM 20° (Scanmarker and screw)	ZBAD0731	15			10-12
SET ZZ BASE AS-UA-SM 45° (Scanmarker and screw)	ZBAD0732	15			10-12
SET ZZ BASE AS-UA-WSM 20° (White Scanmarker and screw)	ZBCD0731	15			10-12
SET ZZ BASE AS-UA-WSM 45° (White Scanmarker and screw)	ZBCD0732	15			10-12
ZZ BASE AS-MUA-ZSC-G 20° (Zirconia screw Gold)	BSCE0731	15			10-12
ZZ BASE AS-MUA-ZSC-G 45° (Zirconia screw Gold)	BSCE0732	15			10-12
ZZ BASE AS-UA-BSC 20° (Black Screw)	BSIE0701	15			10-12
ZZ BASE AS-UA-BSC 45° (Black Screw)	BSIE0702	15			10-12
ZZ BASE AS-UA-BZSC 20° (Black Zirconia screw)	BSIE0731	15			10-12
ZZ BASE AS-UA-BZSC 45° (Black Zirconia screw)	BSIE0732	15			10-12
ZZ BASE AS-UA-ICSC 20° (Impression coping Screw)	BSFE0701	15			10-12
ZZ BASE AS-UA-ICSC 45° (Impression coping Screw)	BSFE0702	15			10-12
ZZ BASE AS-UA-LBSC 20° (Lab Screw)	BSHE0701	15			10-12
ZZ BASE AS-UA-LBSC 45° (Lab Screw)	BSHE0702	15			10-12
ZZ BASE AS-UA-LBZSC 20° (Lab Zirconia screw)	BSHE0731	15			10-12
ZZ BASE AS-UA-LBZSC 45° (Lab Zirconia screw)	BSHE0732	15			10-12
ZZ BASE AS-UA-MI 20° (Model implant)	MIAG0701				10-12
ZZ BASE AS-UA-MI 45° (Model implant)	MIAG0702				10-12
ZZ BASE AS-UA-SC 20° (Screw)	BSAE0701	15			10-12

ASTRA TECH UNIABUTMENT®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE AS-UA-SC 45° (Screw)	BSAE0702	15			10-12
ZZ BASE AS-UA-SC-G 20° (Screw Gold)	BSCE0701	15			10-12
ZZ BASE AS-UA-SC-G 45° (Screw Gold)	BSCE0702	15			10-12
ZZ BASE AS-UA-WMSM-L10 20° (White Metal Scanmarker L10 with Screw)	ZBED0731L10	15			10-12
ZZ BASE AS-UA-WMSM-L10 45° (White Metal Scanmarker L10 with Screw)	ZBED0732L10	15			10-12
ZZ BASE AS-UA-WMSM-L12 20° (White Metal Scanmarker L12 with Screw)	ZBED0731L12	15			10-12
ZZ BASE AS-UA-WMSM-L12 45° (White Metal Scanmarker L12 with Screw)	ZBED0732L12	15			10-12
ZZ BASE AS-UA-WMSM-L15 20° (White Metal Scanmarker L15 with Screw)	ZBED0731L15	15			10-12
ZZ BASE AS-UA-WMSM-L15 45° (White Metal Scanmarker L15 with Screw)	ZBED0732L15	15			10-12
ZZ BASE AS-UA-ZSC 20° (Zirconia screw)	BSAE0731	15			10-12
ZZ BASE AS-UA-ZSC 45° (Zirconia screw)	BSAE0732	15			10-12

A-Z IMPLANT® MC

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE AZ-MC-CCTB A (Conical cemented titan base non hex and screw)	BSAD10531	25			10-12
SET ZZ BASE AZ-MC-SM A (Scanmarker and screw)	ZBAD10531	25			10-12
ZZ BASE AZ-MC-MI A (Model implant)	MIAG10501				10-12
ZZ BASE AZ-MC-SC A (Screw)	BSAE10501	25			10-12
ZZ BASE AZ-MC-ZSC A (Zirconia screw)	BSAE10531	25			10-12

A-Z IMPLANT® VL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE AZ-VL-ABCH 3.2 (Raw-Abutment® customizable hex and screw)	BSAF5694	20-25			10-12
SET ZZ BASE AZ-VL-ABCH 3.8 (Raw-Abutment® customizable hex and screw)	BSAF5691	20-25			10-12
SET ZZ BASE AZ-VL-ABCH 4.4 (Raw-Abutment® customizable hex and screw)	BSAF5692	20-25			10-12
SET ZZ BASE AZ-VL-ABCH 5.5 (Raw-Abutment® customizable hex and screw)	BSAF5693	20-25			10-12
SET ZZ BASE AZ-VL-CCTB 3.2 (Conical cemented titan base non hex and screw)	BSAD5634	20-25			10-12
SET ZZ BASE AZ-VL-CCTB 3.8 (Conical cemented titan base non hex and screw)	BSAD5631	20-25			10-12
SET ZZ BASE AZ-VL-CCTB 4.4 (Conical cemented titan base non hex and screw)	BSAD5632	20-25			10-12
SET ZZ BASE AZ-VL-CCTB 5.5 (Conical cemented titan base non hex and screw)	BSAD5633	20-25			10-12
SET ZZ BASE AZ-VL-CCTB-K85 3.2 (Conical cemented titan base non hex K85 and screw)	BSAD5634K85	20-25			10-12
SET ZZ BASE AZ-VL-CCTB-K85 3.8 (Conical cemented titan base non hex K85 and screw)	BSAD5631K85	20-25			10-12
SET ZZ BASE AZ-VL-CCTB-K85 4.4 (Conical cemented titan base non hex K85 and screw)	BSAD5632K85	20-25			10-12
SET ZZ BASE AZ-VL-CCTB-K85 5.5 (Conical cemented titan base non hex K85 and screw)	BSAD5633K85	20-25			10-12
SET ZZ BASE AZ-VL-ICOTH 3.2 (Impression coping open tray hex and Screw and Driver)	BSAQ5644	15			10-12
SET ZZ BASE AZ-VL-ICOTH 3.8 (Impression coping open tray hex and Screw and Driver)	BSAQ5641	15			10-12
SET ZZ BASE AZ-VL-ICOTH 4.4 (Impression coping open tray hex and Screw and Driver)	BSAQ5642	15			10-12
SET ZZ BASE AZ-VL-ICOTH 5.5 (Impression coping open tray hex and Screw and Driver)	BSAQ5643	15			10-12
SET ZZ BASE AZ-VL-NTB 3.2 (Narrow titan base non hex and screw)	BSAP5634	20-25			10-12
SET ZZ BASE AZ-VL-NTB 3.8 (Narrow titan base non hex and screw)	BSAP5631	20-25			10-12
SET ZZ BASE AZ-VL-NTBH 3.2 (Narrow titan base hex and screw)	BSAP5684	20-25			10-12
SET ZZ BASE AZ-VL-NTBH 3.8 (Narrow titan base hex and screw)	BSAP5681	20-25			10-12
SET ZZ BASE AZ-VL-PCTBH 3.2 (Parallel cemented titan base hex and screw)	BSAB5684	20-25			10-12
SET ZZ BASE AZ-VL-PCTBH 3.8 (Parallel cemented titan base hex and screw)	BSAB5681	20-25			10-12
SET ZZ BASE AZ-VL-PCTBH 4.4 (Parallel cemented titan base hex and screw)	BSAB5682	20-25			10-12
SET ZZ BASE AZ-VL-PCTBH 5.5 (Parallel cemented titan base hex and screw)	BSAB5683	20-25			10-12
SET ZZ BASE AZ-VL-PCTBH-K85 3.2 (Parallel cemented titan base hex K85 and screw)	BSAB5684K85	20-25			10-12
SET ZZ BASE AZ-VL-PCTBH-K85 3.8 (Parallel cemented titan base hex K85 and screw)	BSAB5681K85	20-25			10-12
SET ZZ BASE AZ-VL-PCTBH-K85 4.4 (Parallel cemented titan base hex K85 and screw)	BSAB5682K85	20-25			10-12
SET ZZ BASE AZ-VL-PCTBH-K85 5.5 (Parallel cemented titan base hex K85 and screw)	BSAB5683K85	20-25			10-12
SET ZZ BASE AZ-VL-SA 3.2 (Scananalog and screw)	ZBDE5634	15			10-12

A-Z IMPLANT® VL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE AZ-VL-SA 3.8 (Scananalog and screw)	ZBDE5631	15			10-12
SET ZZ BASE AZ-VL-SA 4.4 (Scananalog and screw)	ZBDE5632	15			10-12
SET ZZ BASE AZ-VL-SA 5.5 (Scananalog and screw)	ZBDE5633	15			10-12
SET ZZ BASE AZ-VL-SM 3.2 (Scanmarker and screw)	ZBAD5634	20-25			10-12
SET ZZ BASE AZ-VL-SM 3.8 (Scanmarker and screw)	ZBAD5631	20-25			10-12
SET ZZ BASE AZ-VL-SM 4.4 (Scanmarker and screw)	ZBAD5632	20-25			10-12
SET ZZ BASE AZ-VL-SM 5.5 (Scanmarker and screw)	ZBAD5633	20-25			10-12
SET ZZ BASE AZ-VL-WSM 3.2 (White Scanmarker and screw)	ZBCD5634	20-25			10-12
SET ZZ BASE AZ-VL-WSM 3.8 (White Scanmarker and screw)	ZBCD5631	20-25			10-12
SET ZZ BASE AZ-VL-WSM 4.4 (White Scanmarker and screw)	ZBCD5632	20-25			10-12
SET ZZ BASE AZ-VL-WSM 5.5 (White Scanmarker and screw)	ZBCD5633	20-25			10-12
ZZ BASE AZ-VL-ASCSC 3.2 (Angled screw-channel Screw)	BSKE5604	20			10-12
ZZ BASE AZ-VL-ASCSC 3.8 (Angled screw-channel Screw)	BSKE5601	20			10-12
ZZ BASE AZ-VL-ASCSC 4.4 (Angled screw-channel Screw)	BSKE5602	20			10-12
ZZ BASE AZ-VL-ICSC 3.2 (Impression coping Screw)	BSFE5604	15			10-12
ZZ BASE AZ-VL-ICSC 3.8 (Impression coping Screw)	BSFE5601	15			10-12
ZZ BASE AZ-VL-ICSC 4.4 (Impression coping Screw)	BSFE5602	15			10-12
ZZ BASE AZ-VL-ICSC 5.5 (Impression coping Screw)	BSFE5603	15			10-12
ZZ BASE AZ-VL-SC 3.2 (Screw)	BSAE5604	20-25			10-12
ZZ BASE AZ-VL-SC 3.8 (Screw)	BSAE5601	20-25			10-12
ZZ BASE AZ-VL-SC 4.4 (Screw)	BSAE5602	20-25			10-12
ZZ BASE AZ-VL-SC 5.5 (Screw)	BSAE5603	20-25			10-12
ZZ BASE AZ-VL-WMSM-L12 3.2 (White Metal Scanmarker L12 with Screw)	ZBED5634L12	15			10-12
ZZ BASE AZ-VL-WMSM-L12 3.8 (White Metal Scanmarker L12 with Screw)	ZBED5631L12	15			10-12
ZZ BASE AZ-VL-WMSM-L12 4.4 (White Metal Scanmarker L12 with Screw)	ZBED5632L12	15			10-12
ZZ BASE AZ-VL-WMSM-L12 5.5 (White Metal Scanmarker L12 with Screw)	ZBED5633L12	15			10-12
ZZ BASE AZ-VL-ZSC 3.2 (Zirconia screw)	BSAE5634	20-25			10-12
ZZ BASE AZ-VL-ZSC 3.8 (Zirconia screw)	BSAE5631	20-25			10-12
ZZ BASE AZ-VL-ZSC 4.4 (Zirconia screw)	BSAE5632	20-25			10-12

A-Z IMPLANT® VL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE AZ-VL-ZSC 5.5 (Zirconia screw)	BSAE5633	20-25			10-12

BEGO SEMADOS® MINI

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE BE-MI-ABCH 2.7-3.1 (Raw-Abutment® customizable hex and screw)	BSAF7591	25			10-12
SET ZZ BASE BE-MI-CCTB 2.7-3.1 (Conical cemented titan base non hex and screw)	BSAD7531	25			10-12
SET ZZ BASE BE-MI-CCTB-K85 2.7-3.1 (Conical cemented titan base non hex K85 and screw)	BSAD7531K85	25			10-12
SET ZZ BASE BE-MI-ICOTH 2.7-3.1 (Impression coping open tray hex and Screw and Driver)	BSAQ7541	15			10-12
SET ZZ BASE BE-MI-PCTBH 2.7-3.1 (Parallel cemented titan base hex and screw)	BSAB7581	25			10-12
SET ZZ BASE BE-MI-SM 2.7-3.1 (Scanmarker and screw)	ZBAD7531	25			10-12
SET ZZ BASE BE-MI-WSM 2.7-3.1 (White Scanmarker and screw)	ZBCD7531	25			10-12
ZZ BASE BE-MI-ICSC 2.7-3.1 (Impression coping Screw)	BSFE7501	15			10-12
ZZ BASE BE-MI-SC 2.7-3.1 (Screw)	BSAE7501	25			10-12
ZZ BASE BE-MI-WMSM-L12 2.7-3.0 (White Metal Scanmarker L12 with Screw)	ZBED7531L12	15			10-12
ZZ BASE BE-MI-ZSC 2.7-3.1 (Zirconia Screw)	BSAE7531	25			10-12

BEGO SEMADOS® S/RI/SC/SCX/RS/RX-LINE / MULTIPLUS SYSTEM

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
SET ZZ BASE BE-SE-S/RI-ABCH 3.0 (Raw-Abutment® customizable hex and screw)	BSAF5087	20			10-12
SET ZZ BASE BE-SE-S/RI-ABCH 3.25 (Raw-Abutment® customizable hex and screw)	BSAF5081	30			10-12
SET ZZ BASE BE-SE-S/RI-ABCH 3.75 (Raw-Abutment® customizable hex and screw)	BSAF5082	30			10-12
SET ZZ BASE BE-SE-S/RI-ABCH 4.1 (Raw-Abutment® customizable hex and screw)	BSAF5083	30			10-12
SET ZZ BASE BE-SE-S/RI-ABCH 4.5 (Raw-Abutment® customizable hex and screw)	BSAF5084	30			10-12
SET ZZ BASE BE-SE-S/RI-ABCH 5.5 (Raw-Abutment® customizable hex and screw)	BSAF5085	30			10-12
SET ZZ BASE BE-SE-S/RI-ASCTBH-K80 3.0 (Angled screw-channel titan base hex K80 and screw)	BSAT5087K80	20			10-12
SET ZZ BASE BE-SE-S/RI-ASCTBH-K80 3.25 (Angled screw-channel titan base hex K80 and screw)	BSAT5081K80	20			10-12
SET ZZ BASE BE-SE-S/RI-ASCTBH-K80 3.75 (Angled screw-channel titan base hex K80 and screw)	BSAT5082K80	20			10-12
SET ZZ BASE BE-SE-S/RI-ASCTBH-K80 4.1 (Angled screw-channel titan base hex K80 and screw)	BSAT5083K80	20			10-12
SET ZZ BASE BE-SE-S/RI-ASCTBH-K80 4.5 (Angled screw-channel titan base hex K80 and screw)	BSAT5084K80	20			10-12
SET ZZ BASE BE-SE-S/RI-ASCTB-K80 3.0 (Angled screw-channel titan base non hex K80 and screw)	BSAT5037K80	20			10-12
SET ZZ BASE BE-SE-S/RI-ASCTB-K80 3.25 (Angled screw-channel titan base non hex K80 and screw)	BSAT5031K80	20			10-12
SET ZZ BASE BE-SE-S/RI-ASCTB-K80 3.75 (Angled screw-channel titan base non hex K80 and screw)	BSAT5032K80	20			10-12
SET ZZ BASE BE-SE-S/RI-ASCTB-K80 4.1 (Angled screw-channel titan base non hex K80 and screw)	BSAT5033K80	20			10-12
SET ZZ BASE BE-SE-S/RI-ASCTB-K80 4.5 (Angled screw-channel titan base non hex K80 and screw)	BSAT5034K80	20			10-12
SET ZZ BASE BE-SE-S/RI-CCTB 3.0 (Conical cemented titan base non hex and screw)	BSAD5037	20			10-12
SET ZZ BASE BE-SE-S/RI-CCTB 3.25 (Conical cemented titan base non hex and screw)	BSAD5031	30			10-12
SET ZZ BASE BE-SE-S/RI-CCTB 3.75 (Conical cemented titan base non hex and screw)	BSAD5032	30			10-12
SET ZZ BASE BE-SE-S/RI-CCTB 4.1 (Conical cemented titan base non hex and screw)	BSAD5033	30			10-12
SET ZZ BASE BE-SE-S/RI-CCTB 4.5 (Conical cemented titan base non hex and screw)	BSAD5034	30			10-12
SET ZZ BASE BE-SE-S/RI-CCTB 5.5 (Conical cemented titan base non hex and screw)	BSAD5035	30			10-12
SET ZZ BASE BE-SE-S/RI-CCTB MP (Conical cemented titan base non hex and screw)	BSAD5036	30			10-12
SET ZZ BASE BE-SE-S/RI-CCTB MP (Conical cemented titan base non hex and screw)	BSAD5046	30			10-12
SET ZZ BASE BE-SE-S/RI-CCTB-K85 3.0 (Conical cemented titan base non hex K85 and screw)	BSAD5037K85	20			10-12
SET ZZ BASE BE-SE-S/RI-CCTB-K85 3.25 (Conical cemented titan base non hex K85 and screw)	BSAD5031K85	30			10-12
SET ZZ BASE BE-SE-S/RI-CCTB-K85 3.75 (Conical cemented titan base non hex K85 and screw)	BSAD5032K85	30			10-12
SET ZZ BASE BE-SE-S/RI-CCTB-K85 4.1 (Conical cemented titan base non hex K85 and screw)	BSAD5033K85	30			10-12
SET ZZ BASE BE-SE-S/RI-CCTB-K85 4.5 (Conical cemented titan base non hex K85 and screw)	BSAD5034K85	30			10-12

All information is subject to change. Errors and omissions excepted. Version: 21.09.2023

BEGO SEMADOS® S/RI/SC/SCX/RS/RX-LINE / MULTIPLUS SYSTEM

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE BE-SE-S/RI-CCTB-K85 5.5 (Conical cemented titan base non hex K85 and screw)	BSAD5035K85	30			10-12
SET ZZ BASE BE-SE-S/RI-CCTB-K85 MP (Conical cemented titan base non hex K85 and screw)	BSAD5036K85	30			10-12
SET ZZ BASE BE-SE-S/RI-CCTB-L10 4.5 (Conical cemented titan base non hex L10 and screw)	BSAD5034L10	30			10-12
SET ZZ BASE BE-SE-S/RI-CCTB-L20 3.25 (Conical cemented titan base non hex L20 and screw)	BSAD5031L20	30			10-12
SET ZZ BASE BE-SE-S/RI-CCTB-L20 3.75 (Conical cemented titan base non hex L20 and screw)	BSAD5032L20	30			10-12
SET ZZ BASE BE-SE-S/RI-CCTB-L20 4.1 (Conical cemented titan base non hex L20 and screw)	BSAD5033L20	30			10-12
SET ZZ BASE BE-SE-S/RI-CCTB-L20 4.5 (Conical cemented titan base non hex L20 and screw)	BSAD5034L20	30			10-12
SET ZZ BASE BE-SE-S/RI-CCTB-L20 5.5 (Conical cemented titan base non hex L20 and screw)	BSAD5035L20	30			10-12
SET ZZ BASE BE-SE-S/RI-CCTB-L30 3.25 (Conical cemented titan base non hex L30 and screw)	BSAD5031L30	30			10-12
SET ZZ BASE BE-SE-S/RI-CCTB-L30 3.75 (Conical cemented titan base non hex L30 and screw)	BSAD5032L30	30			10-12
SET ZZ BASE BE-SE-S/RI-CCTB-L30 4.1 (Conical cemented titan base non hex L30 and screw)	BSAD5033L30	30			10-12
SET ZZ BASE BE-SE-S/RI-CCTB-L30 4.5 (Conical cemented titan base non hex L30 and screw)	BSAD5034L30	30			10-12
SET ZZ BASE BE-SE-S/RI-CCTB-L30 5.5 (Conical cemented titan base non hex L30 and screw)	BSAD5035L30	30			10-12
SET ZZ BASE BE-SE-S/RI-CCTB-L40 3.25 (Conical cemented titan base non hex L40 and screw)	BSAD5031L40	30			10-12
SET ZZ BASE BE-SE-S/RI-CCTB-L40 3.75 (Conical cemented titan base non hex L40 and screw)	BSAD5032L40	30			10-12
SET ZZ BASE BE-SE-S/RI-CCTB-L40 4.1 (Conical cemented titan base non hex L40 and screw)	BSAD5033L40	30			10-12
SET ZZ BASE BE-SE-S/RI-CCTB-L40 4.5 (Conical cemented titan base non hex L40 and screw)	BSAD5034L40	30			10-12
SET ZZ BASE BE-SE-S/RI-CCTB-L40 5.5 (Conical cemented titan base non hex L40 and screw)	BSAD5035L40	30			10-12
SET ZZ BASE BE-SE-S/RI-ICOT 3.0 (Impression coping open tray non hex and Screw and Driver)	BSAQ5067	15			10-12
SET ZZ BASE BE-SE-S/RI-ICOT 3.25 (Impression coping open tray non hex and Screw and Driver)	BSAQ5061	15			10-12
SET ZZ BASE BE-SE-S/RI-ICOT 3.75 (Impression coping open tray non hex and Screw and Driver)	BSAQ5062	15			10-12
SET ZZ BASE BE-SE-S/RI-ICOT 4.1 (Impression coping open tray non hex and Screw and Driver)	BSAQ5063	15			10-12
SET ZZ BASE BE-SE-S/RI-ICOT 4.5 (Impression coping open tray non hex and Screw and Driver)	BSAQ5064	15			10-12
SET ZZ BASE BE-SE-S/RI-ICOT 5.5 (Impression coping open tray non hex and Screw and Driver)	BSAQ5065	15			10-12
SET ZZ BASE BE-SE-S/RI-ICOT MP (Impression coping open tray non hex and Screw and Driver)	BSAQ5066	15			10-12
SET ZZ BASE BE-SE-S/RI-ICOTH 3.0 (Impression coping open tray hex and Screw and Driver)	BSAQ5047	15			10-12
SET ZZ BASE BE-SE-S/RI-ICOTH 3.25 (Impression coping open tray hex and Screw and Driver)	BSAQ5041	15			10-12
SET ZZ BASE BE-SE-S/RI-ICOTH 3.75 (Impression coping open tray hex and Screw and Driver)	BSAQ5042	15			10-12
SET ZZ BASE BE-SE-S/RI-ICOTH 4.1 (Impression coping open tray hex and Screw and Driver)	BSAQ5043	15			10-12

All information is subject to change. Errors and omissions excepted. Version: 21.09.2023

BEGO SEMADOS® S/RI/SC/SCX/RS/RX-LINE / MULTIPLUS SYSTEM

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE BE-SE-S/RI-ICOTH 4.5 (Impression coping open tray hex and Screw and Driver)	BSAQ5044	15			10-12
SET ZZ BASE BE-SE-S/RI-ICOTH 5.5 (Impression coping open tray hex and Screw and Driver)	BSAQ5045	15			10-12
SET ZZ BASE BE-SE-S/RI-NTB 3.0 (Narrow titan base non hex and screw)	BSAP5037	20			10-12
SET ZZ BASE BE-SE-S/RI-NTB 3.25 (Narrow titan base non hex and screw)	BSAP5031	30			10-12
SET ZZ BASE BE-SE-S/RI-NTB 3.75 (Narrow titan base non hex and screw)	BSAP5032	30			10-12
SET ZZ BASE BE-SE-S/RI-NTB 4.1 (Narrow titan base non hex and screw)	BSAP5033	30			10-12
SET ZZ BASE BE-SE-S/RI-NTB 4.5 (Narrow titan base non hex and screw)	BSAP5034	30			10-12
SET ZZ BASE BE-SE-S/RI-NTBH 3.0 (Narrow titan base hex and screw)	BSAP5087	20			10-12
SET ZZ BASE BE-SE-S/RI-NTBH 3.25 (Narrow titan base hex and screw)	BSAP5081	30			10-12
SET ZZ BASE BE-SE-S/RI-NTBH 3.75 (Narrow titan base hex and screw)	BSAP5082	30			10-12
SET ZZ BASE BE-SE-S/RI-NTBH 4.1 (Narrow titan base hex and screw)	BSAP5083	30			10-12
SET ZZ BASE BE-SE-S/RI-NTBH 4.5 (Narrow titan base hex and screw)	BSAP5084	30			10-12
SET ZZ BASE BE-SE-S/RI-PCTBH 3.0 (Parallel cemented titan base hex and screw)	BSAB5087	20			10-12
SET ZZ BASE BE-SE-S/RI-PCTBH 3.25 (Parallel cemented titan base hex and screw)	BSAB5081	30			10-12
SET ZZ BASE BE-SE-S/RI-PCTBH 3.75 (Parallel cemented titan base hex and screw)	BSAB5082	30			10-12
SET ZZ BASE BE-SE-S/RI-PCTBH 4.1 (Parallel cemented titan base hex and screw)	BSAB5083	30			10-12
SET ZZ BASE BE-SE-S/RI-PCTBH 4.5 (Parallel cemented titan base hex and screw)	BSAB5084	30			10-12
SET ZZ BASE BE-SE-S/RI-PCTBH 5.5 (Parallel cemented titan base hex and screw)	BSAB5085	30			10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-K85 3.0 (Parallel cemented titan base hex K85 and screw)	BSAB5087K85	20			10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-K85 3.25 (Parallel cemented titan base hex K85 and screw)	BSAB5081K85	30			10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-K85 3.75 (Parallel cemented titan base hex K85 and screw)	BSAB5082K85	30			10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-K85 4.1 (Parallel cemented titan base hex K85 and screw)	BSAB5083K85	30			10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-K85 4.5 (Parallel cemented titan base hex K85 and screw)	BSAB5084K85	30			10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-K85 5.5 (Parallel cemented titan base hex K85 and screw)	BSAB5085K85	30			10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-L20 3.25 (Parallel cemented titan base hex L20 and screw)	BSAB5091L20	30			10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-L20 3.75 (Parallel cemented titan base hex L20 and screw)	BSAB5092L20	30			10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-L20 4.1 (Parallel cemented titan base hex L20 and screw)	BSAB5093L20	30			10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-L20 4.5 (Parallel cemented titan base hex L20 and screw)	BSAB5094L20	30			10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-L20 5.5 (Parallel cemented titan base hex L20 and screw)	BSAB5095L20	30			10-12

BEGO SEMADOS® S/RI/SC/SCX/RS/RX-LINE / MULTIPLUS SYSTEM

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE BE-SE-S/RI-PCTBH-L30 3.25 (Parallel cemented titan base hex L30 and screw)	BSAB5091L30	30			10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-L30 3.75 (Parallel cemented titan base hex L30 and screw)	BSAB5092L30	30			10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-L30 4.1 (Parallel cemented titan base hex L30 and screw)	BSAB5093L30	30			10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-L30 4.5 (Parallel cemented titan base hex L30 and screw)	BSAB5094L30	30			10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-L30 5.5 (Parallel cemented titan base hex L30 and screw)	BSAB5095L30	30			10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-L40 3.25 (Parallel cemented titan base hex L40 and screw)	BSAB5091L40	30			10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-L40 3.75 (Parallel cemented titan base hex L40 and screw)	BSAB5092L40	30			10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-L40 4.1 (Parallel cemented titan base hex L40 and screw)	BSAB5093L40	30			10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-L40 4.5 (Parallel cemented titan base hex L40 and screw)	BSAB5094L40	30			10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-L40 5.5 (Parallel cemented titan base hex L40 and screw)	BSAB5095L40	30			10-12
SET ZZ BASE BE-SE-S/RI-SA 3.0 (Scananalog and screw)	ZBDE5037	15			10-12
SET ZZ BASE BE-SE-S/RI-SA 3.25 (Scananalog and screw)	ZBDE5031	15			10-12
SET ZZ BASE BE-SE-S/RI-SA 3.75 (Scananalog and screw)	ZBDE5032	15			10-12
SET ZZ BASE BE-SE-S/RI-SA 4.1 (Scananalog and screw)	ZBDE5033	15			10-12
SET ZZ BASE BE-SE-S/RI-SA 4.5 (Scananalog and screw)	ZBDE5034	15			10-12
SET ZZ BASE BE-SE-S/RI-SA 5.5 (Scananalog and screw)	ZBDE5035	15			10-12
SET ZZ BASE BE-SE-S/RI-SA MP (Scananalog and screw)	ZBDE5036	15			10-12
SET ZZ BASE BE-SE-S/RI-SM 3.0 (Scanmarker and screw)	ZBAD5037	30			10-12
SET ZZ BASE BE-SE-S/RI-SM 3.25 (Scanmarker and screw)	ZBAD5031	30			10-12
SET ZZ BASE BE-SE-S/RI-SM 3.75 (Scanmarker and screw)	ZBAD5032	30			10-12
SET ZZ BASE BE-SE-S/RI-SM 4.1 (Scanmarker and screw)	ZBAD5033	30			10-12
SET ZZ BASE BE-SE-S/RI-SM 4.5 (Scanmarker and screw)	ZBAD5034	30			10-12
SET ZZ BASE BE-SE-S/RI-SM 5.5 (Scanmarker and screw)	ZBAD5035	30			10-12
SET ZZ BASE BE-SE-S/RI-SM MP (Scanmarker and screw)	ZBAD5036	30			10-12
SET ZZ BASE BE-SE-S/RI-WSM 3.25 (White Scanmarker and screw)	ZBCD5031	30			10-12
SET ZZ BASE BE-SE-S/RI-WSM 3.75 (White Scanmarker and screw)	ZBCD5032	30			10-12
SET ZZ BASE BE-SE-S/RI-WSM 4.1 (White Scanmarker and screw)	ZBCD5033	30			10-12
SET ZZ BASE BE-SE-S/RI-WSM 4.5 (White Scanmarker and screw)	ZBCD5034	30			10-12
SET ZZ BASE BE-SE-S/RI-WSM 5.5 (White Scanmarker and screw)	ZBCD5035	30			10-12

BEGO SEMADOS® S/RI/SC/SCX/RS/RXS-LINE / MULTIPLUS SYSTEM

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE BE-SE-S/RI-WSM MP (White Scanmarker and screw)	ZBCD5036	30			10-12
ZZ BASE BE-SE-S/RI-ASCSC 3.0 (Angled screw-channel Screw)	BSKE5007	20			10-12
ZZ BASE BE-SE-S/RI-ASCSC 3.25 (Angled screw-channel Screw)	BSKE5001	20			10-12
ZZ BASE BE-SE-S/RI-ASCSC 3.75 (Angled screw-channel Screw)	BSKE5002	20			10-12
ZZ BASE BE-SE-S/RI-ASCSC 4.1 (Angled screw-channel Screw)	BSKE5003	20			10-12
ZZ BASE BE-SE-S/RI-ASCSC 4.5 (Angled screw-channel Screw)	BSKE5004	20			10-12
ZZ BASE BE-SE-S/RI-ICSC 3.0 (Impression coping Screw)	BSFE5007	15			10-12
ZZ BASE BE-SE-S/RI-ICSC 3.25 (Impression coping Screw)	BSFE5001	15			10-12
ZZ BASE BE-SE-S/RI-ICSC 3.75 (Impression coping Screw)	BSFE5002	15			10-12
ZZ BASE BE-SE-S/RI-ICSC 4.1 (Impression coping Screw)	BSFE5003	15			10-12
ZZ BASE BE-SE-S/RI-ICSC 4.5 (Impression coping Screw)	BSFE5004	15			10-12
ZZ BASE BE-SE-S/RI-ICSC 5.5 (Impression coping Screw)	BSFE5005	15			10-12
ZZ BASE BE-SE-S/RI-ICSC MP (Impression coping Screw)	BSFE5006	15			10-12
ZZ BASE BE-SE-S/RI-MI 3.0 (Model implant)	MIAG5007				10-12
ZZ BASE BE-SE-S/RI-MI 3.25 (Model implant)	MIAG5001				10-12
ZZ BASE BE-SE-S/RI-MI 3.75 (Model implant)	MIAG5002				10-12
ZZ BASE BE-SE-S/RI-MI 4.1 (Model implant)	MIAG5003				10-12
ZZ BASE BE-SE-S/RI-MI 4.5 (Model implant)	MIAG5004				10-12
ZZ BASE BE-SE-S/RI-MI 5.5 (Model implant)	MIAG5005				10-12
ZZ BASE BE-SE-S/RI-MI MP (Model implant)	MIAG5006				10-12
ZZ BASE BE-SE-S/RI-SC 3.0 (Screw)	BSAE5007	20			10-12
ZZ BASE BE-SE-S/RI-SC 3.25 (Screw)	BSAE5001	30			10-12
ZZ BASE BE-SE-S/RI-SC 3.75 (Screw)	BSAE5002	30			10-12
ZZ BASE BE-SE-S/RI-SC 4.1 (Screw)	BSAE5003	30			10-12
ZZ BASE BE-SE-S/RI-SC 4.5 (Screw)	BSAE5004	30			10-12
ZZ BASE BE-SE-S/RI-SC 5.5 (Screw)	BSAE5005	30			10-12
ZZ BASE BE-SE-S/RI-SC MP (Screw)	BSAE5006	30			10-12
ZZ BASE BE-SE-S/RI-WMSM-L12 3.0 (White Metal Scanmarker L12 with Screw)	ZBED5037L12	15			10-12
ZZ BASE BE-SE-S/RI-WMSM-L12 3.25 (White Metal Scanmarker L12 with Screw)	ZBED5031L12	15			10-12

BEGO SEMADOS® S/RI/SC/SCX/RS/R SX-LINE / MULTIPLUS SYSTEM

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE BE-SE-S/RI-WMSM-L12 3.75 (White Metal Scanmarker L12 with Screw)	ZBED5032L12	15			10-12
ZZ BASE BE-SE-S/RI-WMSM-L12 4.1 (White Metal Scanmarker L12 with Screw)	ZBED5033L12	15			10-12
ZZ BASE BE-SE-S/RI-WMSM-L12 4.5 (White Metal Scanmarker L12 with Screw)	ZBED5034L12	15			10-12
ZZ BASE BE-SE-S/RI-WMSM-L12 5.5 (White Metal Scanmarker L12 with Screw)	ZBED5035L12	15			10-12
ZZ BASE BE-SE-S/RI-WMSM-L12 MP (White Metal Scanmarker L12 with Screw)	ZBED5036L12	15			10-12
ZZ BASE BE-SE-S/RI-ZSC 3.0 (Zirconia screw)	BSAE5037	20			10-12
ZZ BASE BE-SE-S/RI-ZSC 3.25 (Zirconia screw)	BSAE5031	30			10-12
ZZ BASE BE-SE-S/RI-ZSC 3.75 (Zirconia screw)	BSAE5032	30			10-12
ZZ BASE BE-SE-S/RI-ZSC 4.1 (Zirconia screw)	BSAE5033	30			10-12
ZZ BASE BE-SE-S/RI-ZSC 4.5 (Zirconia screw)	BSAE5034	30			10-12
ZZ BASE BE-SE-S/RI-ZSC 5.5 (Zirconia screw)	BSAE5035	30			10-12
ZZ BASE BE-SE-S/RI-ZSC MP (Zirconia screw)	BSAE5036	30			10-12

BIOHORIZONS® EXTERNAL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE BH-EX-ABCH 3.5 (Raw-Abutment® customizable hex and screw)	BSAF3891	30			10-12
SET ZZ BASE BH-EX-ABCH 4.0 (Raw-Abutment® customizable hex and screw)	BSAF3892	30			10-12
SET ZZ BASE BH-EX-ABCH 5.0 (Raw-Abutment® customizable hex and screw)	BSAF3893	30			10-12
SET ZZ BASE BH-EX-ABCH 6.0 (Raw-Abutment® customizable hex and screw)	BSAF3894	30			10-12
SET ZZ BASE BH-EX-CCTB 3.5 (Conical cemented titan base non hex and screw)	BSAD3831	30			10-12
SET ZZ BASE BH-EX-CCTB 4.0 (Conical cemented titan base non hex and screw)	BSAD3832	30			10-12
SET ZZ BASE BH-EX-CCTB 5.0 (Conical cemented titan base non hex and screw)	BSAD3833	30			10-12
SET ZZ BASE BH-EX-CCTB 6.0 (Conical cemented titan base non hex and screw)	BSAD3834	30			10-12
SET ZZ BASE BH-EX-ICOTH 3.5 (Impression coping open tray hex and Screw and Driver)	BSAQ3841	15			10-12
SET ZZ BASE BH-EX-ICOTH 4.0 (Impression coping open tray hex and Screw and Driver)	BSAQ3842	15			10-12
SET ZZ BASE BH-EX-ICOTH 5.0 (Impression coping open tray hex and Screw and Driver)	BSAQ3843	15			10-12
SET ZZ BASE BH-EX-ICOTH 6.0 (Impression coping open tray hex and Screw and Driver)	BSAQ3844	15			10-12
SET ZZ BASE BH-EX-NTB 3.5 (Narrow titan base non hex and screw)	BSAP3831	30			10-12
SET ZZ BASE BH-EX-NTBH 3.5 (Narrow titan base hex and screw)	BSAP3881	30			10-12
SET ZZ BASE BH-EX-PCTBH 3.5 (Parallel cemented titan base hex and screw)	BSAB3881	30			10-12
SET ZZ BASE BH-EX-PCTBH 4.0 (Parallel cemented titan base hex and screw)	BSAB3882	30			10-12
SET ZZ BASE BH-EX-PCTBH 5.0 (Parallel cemented titan base hex and screw)	BSAB3883	30			10-12
SET ZZ BASE BH-EX-PCTBH 6.0 (Parallel cemented titan base hex and screw)	BSAB3884	30			10-12
SET ZZ BASE BH-EX-SA 3.5 (Scananalog and screw)	ZBDE3831	15			10-12
SET ZZ BASE BH-EX-SA 4.0 (Scananalog and screw)	ZBDE3832	15			10-12
SET ZZ BASE BH-EX-SA 5.0 (Scananalog and screw)	ZBDE3833	15			10-12
SET ZZ BASE BH-EX-SA 6.0 (Scananalog and screw)	ZBDE3834	15			10-12
SET ZZ BASE BH-EX-SM 3.5 (Scanmarker and screw)	ZBAD3831	30			10-12
SET ZZ BASE BH-EX-SM 4.0 (Scanmarker and screw)	ZBAD3832	30			10-12
SET ZZ BASE BH-EX-SM 5.0 (Scanmarker and screw)	ZBAD3833	30			10-12
SET ZZ BASE BH-EX-SM 6.0 (Scanmarker and screw)	ZBAD3834	30			10-12
SET ZZ BASE BH-EX-WSM 3.5 (White Scanmarker and screw)	ZBCD3831	30			10-12
SET ZZ BASE BH-EX-WSM 4.0 (White Scanmarker and screw)	ZBCD3832	30			10-12
SET ZZ BASE BH-EX-WSM 5.0 (White Scanmarker and screw)	ZBCD3833	30			10-12

BIOHORIZONS® EXTERNAL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE BH-EX-WSM 6.0 (White Scanmarker and screw)	ZBCD3834	30			10-12
SET ZZ BASE BH-EX-ZZMUA-L15 3.5 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH3831L15	25	30		10-12
SET ZZ BASE BH-EX-ZZMUA-L15 4.0 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH3832L15	25	30		10-12
SET ZZ BASE BH-EX-ZZMUA-L15 5.0 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH3833L15	25	30		10-12
SET ZZ BASE BH-EX-ZZMUA-L25 3.5 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH3831L25	25	30		10-12
SET ZZ BASE BH-EX-ZZMUA-L25 4.0 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH3832L25	25	30		10-12
SET ZZ BASE BH-EX-ZZMUA-L25 5.0 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH3833L25	25	30		10-12
SET ZZ BASE BH-EX-ZZMUA-L35 3.5 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH3831L35	25	30		10-12
SET ZZ BASE BH-EX-ZZMUA-L35 4.0 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH3832L35	25	30		10-12
SET ZZ BASE BH-EX-ZZMUA-L35 5.0 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH3833L35	25	30		10-12
SET ZZ BASE BH-EX-ZZMUA-L45 3.5 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH3831L45	25	30		10-12
SET ZZ BASE BH-EX-ZZMUA-L45 4.0 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH3832L45	25	30		10-12
SET ZZ BASE BH-EX-ZZMUA-L45 5.0 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH3833L45	25	30		10-12
ZZ BASE BH-EX-ICSC 3.5 (Impression coping Screw)	BSFE3801	15			10-12
ZZ BASE BH-EX-ICSC 4.0 (Impression coping Screw)	BSFE3802	15			10-12
ZZ BASE BH-EX-ICSC 5.0 (Impression coping Screw)	BSFE3803	15			10-12
ZZ BASE BH-EX-ICSC 6.0 (Impression coping Screw)	BSFE3804	15			10-12
ZZ BASE BH-EX-MI 3.5 (Model implant)	MIAG3801				10-12
ZZ BASE BH-EX-MI 5.0 (Model implant)	MIAG3803				10-12
ZZ BASE BH-EX-MI 6.0 (Model implant)	MIAG3804				10-12
ZZ BASE BH-EX-SC 3.5 (Screw)	BSAE3801	30			10-12
ZZ BASE BH-EX-SC 4.0 (Screw)	BSAE3802	30			10-12
ZZ BASE BH-EX-SC 5.0 (Screw)	BSAE3803	30			10-12
ZZ BASE BH-EX-SC 6.0 (Screw)	BSAE3804	30			10-12
ZZ BASE BH-EX-WMSM-L12 3.5 (White Metal Scanmarker L12 with Screw)	ZBED3831L12	15			10-12
ZZ BASE BH-EX-WMSM-L12 4.0 (White Metal Scanmarker L12 with Screw)	ZBED3832L12	15			10-12
ZZ BASE BH-EX-WMSM-L12 5.0 (White Metal Scanmarker L12 with Screw)	ZBED3833L12	15			10-12
ZZ BASE BH-EX-WMSM-L12 6.0 (White Metal Scanmarker L12 with Screw)	ZBED3834L12	15			10-12
ZZ BASE BH-EX-ZSC 3.5 (Zirconia screw)	BSAE3831	30			10-12

BIOHORIZONS® EXTERNAL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
ZZ BASE BH-EX-ZSC 4.0 (Zirconia screw)	BSAE3832	30			10-12
ZZ BASE BH-EX-ZSC 5.0 (Zirconia screw)	BSAE3833	30			10-12
ZZ BASE BH-EX-ZSC 6.0 (Zirconia screw)	BSAE3834	30			10-12

BIOHORIZONS® INTERNAL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE BH-INT-ABCH 3.0 (Raw-Abutment® customizable hex and screw)	BSAF3781	30			10-12
SET ZZ BASE BH-INT-ABCH 3.5 (Raw-Abutment® customizable hex and screw)	BSAF3782	30			10-12
SET ZZ BASE BH-INT-ABCH 4.5 (Raw-Abutment® customizable hex and screw)	BSAF3783	30			10-12
SET ZZ BASE BH-INT-ABCH 5.7 (Raw-Abutment® customizable hex and screw)	BSAF3784	30			10-12
SET ZZ BASE BH-INT-ASCTBH-K80 3.0 (Angled screw-channel titan base hex K80 and screw)	BSAT3781K80	20			10-12
SET ZZ BASE BH-INT-ASCTBH-K80 3.5 (Angled screw-channel titan base hex K80 and screw)	BSAT3782K80	20			10-12
SET ZZ BASE BH-INT-ASCTBH-K80 4.5 (Angled screw-channel titan base hex K80 and screw)	BSAT3783K80	20			10-12
SET ZZ BASE BH-INT-ASCTB-K80 3.0 (Angled screw-channel titan base non hex K80 and screw)	BSAT3731K80	20			10-12
SET ZZ BASE BH-INT-ASCTB-K80 3.5 (Angled screw-channel titan base non hex K80 and screw)	BSAT3732K80	20			10-12
SET ZZ BASE BH-INT-ASCTB-K80 4.5 (Angled screw-channel titan base non hex K80 and screw)	BSAT3733K80	20			10-12
SET ZZ BASE BH-INT-CCTB 3.0 (Conical cemented titan base non hex and screw)	BSAD3731	30			10-12
SET ZZ BASE BH-INT-CCTB 3.5 (Conical cemented titan base non hex and screw)	BSAD3732	30			10-12
SET ZZ BASE BH-INT-CCTB 4.5 (Conical cemented titan base non hex and screw)	BSAD3733	30			10-12
SET ZZ BASE BH-INT-CCTB 5.7 (Conical cemented titan base non hex and screw)	BSAD3734	30			10-12
SET ZZ BASE BH-INT-CCTB-K85 3.0 (Conical cemented titan base non hex K85 and screw)	BSAD3731K85	30			10-12
SET ZZ BASE BH-INT-CCTB-K85 3.5 (Conical cemented titan base non hex K85 and screw)	BSAD3732K85	30			10-12
SET ZZ BASE BH-INT-CCTB-K85 4.5 (Conical cemented titan base non hex K85 and screw)	BSAD3733K85	30			10-12
SET ZZ BASE BH-INT-CCTB-K85 5.7 (Conical cemented titan base non hex K85 and screw)	BSAD3734K85	30			10-12
SET ZZ BASE BH-INT-CCTB-L10 4.5 (Conical cemented titan base non hex L10 and screw)	BSAD3733L10	30			10-12
SET ZZ BASE BH-INT-CCTB-L10 5.7 (Conical cemented titan base non hex L10 and screw)	BSAD3734L10	30			10-12
SET ZZ BASE BH-INT-CCTB-L20 3.0 (Conical cemented titan base non hex L20 and screw)	BSAD3731L20	30			10-12
SET ZZ BASE BH-INT-CCTB-L20 3.5 (Conical cemented titan base non hex L20 and screw)	BSAD3732L20	30			10-12
SET ZZ BASE BH-INT-CCTB-L20 4.5 (Conical cemented titan base non hex L20 and screw)	BSAD3733L20	30			10-12
SET ZZ BASE BH-INT-CCTB-L20 5.7 (Conical cemented titan base non hex L20 and screw)	BSAD3734L20	30			10-12
SET ZZ BASE BH-INT-CCTB-L30 3.0 (Conical cemented titan base non hex L30 and screw)	BSAD3731L30	30			10-12
SET ZZ BASE BH-INT-CCTB-L30 3.5 (Conical cemented titan base non hex L30 and screw)	BSAD3732L30	30			10-12
SET ZZ BASE BH-INT-CCTB-L30 4.5 (Conical cemented titan base non hex L30 and screw)	BSAD3733L30	30			10-12
SET ZZ BASE BH-INT-CCTB-L30 5.7 (Conical cemented titan base non hex L30 and screw)	BSAD3734L30	30			10-12
SET ZZ BASE BH-INT-CCTB-L40 3.0 (Conical cemented titan base non hex L40 and screw)	BSAD3731L40	30			10-12

BIOHORIZONS® INTERNAL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE BH-INT-CCTB-L40 3.5 (Conical cemented titan base non hex L40 and screw)	BSAD3732L40	30			10-12
SET ZZ BASE BH-INT-CCTB-L40 4.5 (Conical cemented titan base non hex L40 and screw)	BSAD3733L40	30			10-12
SET ZZ BASE BH-INT-CCTB-L40 5.7 (Conical cemented titan base non hex L40 and screw)	BSAD3734L40	30			10-12
SET ZZ BASE BH-INT-ICOTH 3.0 (Impression coping open tray hex and Screw and Driver)	BSAQ3741	15			10-12
SET ZZ BASE BH-INT-ICOTH 3.5 (Impression coping open tray hex and Screw and Driver)	BSAQ3742	15			10-12
SET ZZ BASE BH-INT-ICOTH 4.5 (Impression coping open tray hex and Screw and Driver)	BSAQ3743	15			10-12
SET ZZ BASE BH-INT-ICOTH 5.7 (Impression coping open tray hex and Screw and Driver)	BSAQ3744	15			10-12
SET ZZ BASE BH-INT-NTB 3.0 (Narrow titan base non hex and screw)	BSAP3731	30			10-12
SET ZZ BASE BH-INT-NTB 3.5 (Narrow titan base non hex and screw)	BSAP3732	30			10-12
SET ZZ BASE BH-INT-NTBH 3.0 (Narrow titan base hex and screw)	BSAP3781	30			10-12
SET ZZ BASE BH-INT-NTBH 3.5 (Narrow titan base hex and screw)	BSAP3782	30			10-12
SET ZZ BASE BH-INT-PCTBH 3.0 (Parallel cemented titan base hex and screw)	BSAB3781	30			10-12
SET ZZ BASE BH-INT-PCTBH 3.5 (Parallel cemented titan base hex and screw)	BSAB3782	30			10-12
SET ZZ BASE BH-INT-PCTBH 4.5 (Parallel cemented titan base hex and screw)	BSAB3783	30			10-12
SET ZZ BASE BH-INT-PCTBH 5.7 (Parallel cemented titan base hex and screw)	BSAB3784	30			10-12
SET ZZ BASE BH-INT-PCTBH-K85 3.0 (Parallel cemented titan base hex K85 and screw)	BSAB3781K85	30			10-12
SET ZZ BASE BH-INT-PCTBH-K85 3.5 (Parallel cemented titan base hex K85 and screw)	BSAB3782K85	30			10-12
SET ZZ BASE BH-INT-PCTBH-K85 4.5 (Parallel cemented titan base hex K85 and screw)	BSAB3783K85	30			10-12
SET ZZ BASE BH-INT-PCTBH-K85 5.7 (Parallel cemented titan base hex K85 and screw)	BSAB3784K85	30			10-12
SET ZZ BASE BH-INT-PCTBH-L10 4.5 (Parallel cemented titan base hex L10 and screw)	BSAB3783L10	30			10-12
SET ZZ BASE BH-INT-PCTBH-L10 5.7 (Parallel cemented titan base hex L10 and screw)	BSAB3784L10	30			10-12
SET ZZ BASE BH-INT-PCTBH-L20 3.0 (Parallel cemented titan base hex L20 and screw)	BSAB3781L20	30			10-12
SET ZZ BASE BH-INT-PCTBH-L20 3.5 (Parallel cemented titan base hex L20 and screw)	BSAB3782L20	30			10-12
SET ZZ BASE BH-INT-PCTBH-L20 4.5 (Parallel cemented titan base hex L20 and screw)	BSAB3783L20	30			10-12
SET ZZ BASE BH-INT-PCTBH-L20 5.7 (Parallel cemented titan base hex L20 and screw)	BSAB3784L20	30			10-12
SET ZZ BASE BH-INT-PCTBH-L30 3.0 (Parallel cemented titan base hex L30 and screw)	BSAB3781L30	30			10-12
SET ZZ BASE BH-INT-PCTBH-L30 3.5 (Parallel cemented titan base hex L30 and screw)	BSAB3782L30	30			10-12
SET ZZ BASE BH-INT-PCTBH-L30 4.5 (Parallel cemented titan base hex L30 and screw)	BSAB3783L30	30			10-12
SET ZZ BASE BH-INT-PCTBH-L30 5.7 (Parallel cemented titan base hex L30 and screw)	BSAB3784L30	30			10-12

BIOHORIZONS® INTERNAL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE BH-INT-PCTBH-L40 3.0 (Parallel cemented titan base hex L40 and screw)	BSAB3781L40	30			10-12
SET ZZ BASE BH-INT-PCTBH-L40 3.5 (Parallel cemented titan base hex L40 and screw)	BSAB3782L40	30			10-12
SET ZZ BASE BH-INT-PCTBH-L40 4.5 (Parallel cemented titan base hex L40 and screw)	BSAB3783L40	30			10-12
SET ZZ BASE BH-INT-PCTBH-L40 5.7 (Parallel cemented titan base hex L40 and screw)	BSAB3784L40	30			10-12
SET ZZ BASE BH-INT-SA 3.0 (Scananalog and screw)	ZBDE3731	15			10-12
SET ZZ BASE BH-INT-SA 3.5 (Scananalog and screw)	ZBDE3732	15			10-12
SET ZZ BASE BH-INT-SA 4.5 (Scananalog and screw)	ZBDE3733	15			10-12
SET ZZ BASE BH-INT-SA 5.7 (Scananalog and screw)	ZBDE3734	15			10-12
SET ZZ BASE BH-INT-SM 3.0 (Scanmarker and screw)	ZBAD3731	30			10-12
SET ZZ BASE BH-INT-SM 3.5 (Scanmarker and screw)	ZBAD3732	30			10-12
SET ZZ BASE BH-INT-SM 4.5 (Scanmarker and screw)	ZBAD3733	30			10-12
SET ZZ BASE BH-INT-SM 5.7 (Scanmarker and screw)	ZBAD3734	30			10-12
SET ZZ BASE BH-INT-WSM 3.0 (White Scanmarker and screw)	ZBCD3731	30			10-12
SET ZZ BASE BH-INT-WSM 3.5 (White Scanmarker and screw)	ZBCD3732	30			10-12
SET ZZ BASE BH-INT-WSM 4.5 (White Scanmarker and screw)	ZBCD3733	30			10-12
SET ZZ BASE BH-INT-WSM 5.7 (White Scanmarker and screw)	ZBCD3734	30			10-12
SET ZZ BASE BH-INT-ZZMUA-L15 3.0 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH3731L15	25	30		10-12
SET ZZ BASE BH-INT-ZZMUA-L15 3.5 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH3732L15	25	30		10-12
SET ZZ BASE BH-INT-ZZMUA-L15 4.5 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH3733L15	25	30		10-12
SET ZZ BASE BH-INT-ZZMUA-L25 3.0 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH3731L25	25	30		10-12
SET ZZ BASE BH-INT-ZZMUA-L25 3.5 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH3732L25	25	30		10-12
SET ZZ BASE BH-INT-ZZMUA-L25 4.5 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH3733L25	25	30		10-12
SET ZZ BASE BH-INT-ZZMUA-L35 3.0 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH3731L35	25	30		10-12
SET ZZ BASE BH-INT-ZZMUA-L35 3.5 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH3732L35	25	30		10-12
SET ZZ BASE BH-INT-ZZMUA-L35 4.5 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH3733L35	25	30		10-12
SET ZZ BASE BH-INT-ZZMUA-L45 3.0 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH3731L45	25	30		10-12
SET ZZ BASE BH-INT-ZZMUA-L45 3.5 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH3732L45	25	30		10-12
SET ZZ BASE BH-INT-ZZMUA-L45 4.5 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH3733L45	25	30		10-12
ZZ BASE BH-INT-ASCSC 3.0 (Angled screw-channel Screw)	BSKE3701	20			10-12

BIOHORIZONS® INTERNAL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE BH-INT-ASCSC 3.5 (Angled screw-channel Screw)	BSKE3702	20			10-12
ZZ BASE BH-INT-ASCSC 4.5 (Angled screw-channel Screw)	BSKE3703	20			10-12
ZZ BASE BH-INT-ICSC 3.0 (Impression coping Screw)	BSFE3701	15			10-12
ZZ BASE BH-INT-ICSC 3.5 (Impression coping Screw)	BSFE3702	15			10-12
ZZ BASE BH-INT-ICSC 4.5 (Impression coping Screw)	BSFE3703	15			10-12
ZZ BASE BH-INT-ICSC 5.7 (Impression coping Screw)	BSFE3704	15			10-12
ZZ BASE BH-INT-MI 3.0 (Model implant)	MIAG3701				10-12
ZZ BASE BH-INT-MI 3.5 (Model implant)	MIAG3702				10-12
ZZ BASE BH-INT-MI 4.5 (Model implant)	MIAG3703				10-12
ZZ BASE BH-INT-MI 5.7 (Model implant)	MIAG3704				10-12
ZZ BASE BH-INT-SC 3.0 (Screw)	BSAE3701	30			10-12
ZZ BASE BH-INT-SC 3.5 (Screw)	BSAE3702	30			10-12
ZZ BASE BH-INT-SC 4.5 (Screw)	BSAE3703	30			10-12
ZZ BASE BH-INT-SC 5.7 (Screw)	BSAE3704	30			10-12
ZZ BASE BH-INT-WMSM-L12 3.0 (White Metal Scanmarker L12 with Screw)	ZBED3731L12	15			10-12
ZZ BASE BH-INT-WMSM-L12 3.5 (White Metal Scanmarker L12 with Screw)	ZBED3732L12	15			10-12
ZZ BASE BH-INT-WMSM-L12 4.5 (White Metal Scanmarker L12 with Screw)	ZBED3733L12	15			10-12
ZZ BASE BH-INT-WMSM-L12 5.7 (White Metal Scanmarker L12 with Screw)	ZBED3734L12	15			10-12
ZZ BASE BH-INT-ZSC 3.0 (Zirconia screw)	BSAE3731	30			10-12
ZZ BASE BH-INT-ZSC 3.5 (Zirconia screw)	BSAE3732	30			10-12
ZZ BASE BH-INT-ZSC 4.5 (Zirconia screw)	BSAE3733	30			10-12
ZZ BASE BH-INT-ZSC 5.7 (Zirconia screw)	BSAE3734	30			10-12

BIOHORIZONS® MULTI-UNIT ABUTMENT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE BH-MUA-CCTB (Conical cemented titan base non hex and screw)	BSAD7731	15			10-12
SET ZZ BASE BH-MUA-CCTB-K85 (Conical cemented titan base non hex K85 and screw)	BSAD7731K85	15			10-12
SET ZZ BASE BH-MUA-ICOT (Impression coping open tray non hex and Screw and Driver)	BSAQ7761	15			10-12
SET ZZ BASE BH-MUA-SA (Scananalog and screw)	ZBDE7731	15			10-12
SET ZZ BASE BH-MUA-SM (Scanmarker and screw)	ZBAD7731	15			10-12
SET ZZ BASE BH-MUA-WSM (White Scanmarker and screw)	ZBCD7731	15			10-12
ZZ BASE BH-MUA-ICSC (Impression coping Screw)	BSFE7701	15			10-12
ZZ BASE BH-MUA-SC (Screw)	BSAE7701	15			10-12
ZZ BASE BH-MUA-WMSM-L12 (White Metal Scanmarker L12 with Screw)	ZBED7731L12	15			10-12
ZZ BASE BH-MUA-ZSC (Zirconia screw)	BSAE7731	15			10-12

BIOMET 3I™ CERTAIN® IMPLANT SYSTEM

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE B-C-ABCH 3.4 (Raw-Abutment® customizable hex and screw)	BSAF0891	20			10-12
SET ZZ BASE B-C-ABCH 4.1 (Raw-Abutment® customizable hex and screw)	BSAF0892	20			10-12
SET ZZ BASE B-C-ABCH 5.0 (Raw-Abutment® customizable hex and screw)	BSAF0893	20			10-12
SET ZZ BASE B-C-ABCH 6.0 (Raw-Abutment® customizable hex and screw)	BSAF0894	20			10-12
SET ZZ BASE B-C-ABCH-D14 3.4 (Raw-Abutment® customizable hex D14 and screw)	BSAF0891D14	20			10-12
SET ZZ BASE B-C-ABCH-D14 4.1 (Raw-Abutment® customizable hex D14 and screw)	BSAF0892D14	20			10-12
SET ZZ BASE B-C-ABCH-D14 5.0 (Raw-Abutment® customizable hex D14 and screw)	BSAF0893D14	20			10-12
SET ZZ BASE B-C-ABCH-D14 6.0 (Raw-Abutment® customizable hex D14 and screw)	BSAF0894D14	20			10-12
SET ZZ BASE B-C-ASCTBH-K80 3.4 (Angled screw-channel titan base hex K80 and screw)	BSAT0881K80	20			10-12
SET ZZ BASE B-C-ASCTBH-K80 4.1 (Angled screw-channel titan base hex K80 and screw)	BSAT0882K80	20			10-12
SET ZZ BASE B-C-ASCTB-K80 3.4 (Angled screw-channel titan base non hex K80 and thick screw)	BSAT0831K80	20			10-12
SET ZZ BASE B-C-ASCTB-K80 4.1 (Angled screw-channel titan base non hex K80 and thick screw)	BSAT0832K80	20			10-12
SET ZZ BASE B-C-CCTB 3.4 (Conical cemented titan base non hex and screw)	BSAD0831	20			10-12
SET ZZ BASE B-C-CCTB 3.4 (Conical cemented titan base non hex and thick screw)	BSAD0841	20			10-12
SET ZZ BASE B-C-CCTB 4.1 (Conical cemented titan base non hex and screw)	BSAD0832	20			10-12
SET ZZ BASE B-C-CCTB 4.1 (Conical cemented titan base non hex and thick screw)	BSAD0842	20			10-12
SET ZZ BASE B-C-CCTB 5.0 (Conical cemented titan base non hex and screw)	BSAD0833	20			10-12
SET ZZ BASE B-C-CCTB 5.0 (Conical cemented titan base non hex and thick screw)	BSAD0843	20			10-12
SET ZZ BASE B-C-CCTB 6.0 (Conical cemented titan base non hex and screw)	BSAD0834	20			10-12
SET ZZ BASE B-C-CCTB 6.0 (Conical cemented titan base non hex and thick screw)	BSAD0844	20			10-12
SET ZZ BASE B-C-CCTB-G 3.4 (Conical cemented titan base non hex and screw Gold)	BSCD0831	20			10-12
SET ZZ BASE B-C-CCTB-G 3.4 (Conical cemented titan base non hex and screw Gold)	BSCD0841	20			10-12
SET ZZ BASE B-C-CCTB-G 4.1 (Conical cemented titan base non hex and screw Gold)	BSCD0832	20			10-12
SET ZZ BASE B-C-CCTB-G 4.1 (Conical cemented titan base non hex and screw Gold)	BSCD0842	20			10-12
SET ZZ BASE B-C-CCTB-G 5.0 (Conical cemented titan base non hex and screw Gold)	BSCD0833	20			10-12
SET ZZ BASE B-C-CCTB-G 5.0 (Conical cemented titan base non hex and screw Gold)	BSCD0843	20			10-12
SET ZZ BASE B-C-CCTB-G 6.0 (Conical cemented titan base non hex and screw Gold)	BSCD0834	20			10-12
SET ZZ BASE B-C-CCTB-G 6.0 (Conical cemented titan base non hex and screw Gold)	BSCD0844	20			10-12
SET ZZ BASE B-C-CCTB-K85 3.4 (Conical cemented titan base non hex K85 and thick screw)	BSAD0831K85	20			10-12

BIOMET 3I™ CERTAIN® IMPLANT SYSTEM

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
SET ZZ BASE B-C-CCTB-K85 4.1 (Conical cemented titan base non hex K85 and thick screw)	BSAD0832K85	20			10-12
SET ZZ BASE B-C-CCTB-K85 5.0 (Conical cemented titan base non hex K85 and thick screw)	BSAD0833K85	20			10-12
SET ZZ BASE B-C-CCTB-K85 6.0 (Conical cemented titan base non hex K85 and thick screw)	BSAD0834K85	20			10-12
SET ZZ BASE B-C-CCTB-L10 3.4 (Conical cemented titan base non hex L10 and thick screw)	BSAD0831L10	20			10-12
SET ZZ BASE B-C-CCTB-L10 4.1 (Conical cemented titan base non hex L10 and thick screw)	BSAD0832L10	20			10-12
SET ZZ BASE B-C-CCTB-L10 5.0 (Conical cemented titan base non hex L10 and thick screw)	BSAD0833L10	20			10-12
SET ZZ BASE B-C-CCTB-L10 6.0 (Conical cemented titan base non hex L10 and thick screw)	BSAD0834L10	20			10-12
SET ZZ BASE B-C-CCTB-L10-G 3.4 (Conical cemented titan base non hex L10 and screw Gold)	BSCD0831L10	20			10-12
SET ZZ BASE B-C-CCTB-L10-G 4.1 (Conical cemented titan base non hex L10 and screw Gold)	BSCD0832L10	20			10-12
SET ZZ BASE B-C-CCTB-L10-G 5.0 (Conical cemented titan base non hex L10 and screw Gold)	BSCD0833L10	20			10-12
SET ZZ BASE B-C-CCTB-L10-G 6.0 (Conical cemented titan base non hex L10 and screw Gold)	BSCD0834L10	20			10-12
SET ZZ BASE B-C-CCTB-L10-K85 3.4 (Conical cemented titan base non hex L10 K85 and thick screw)	BSAD0831L10K85	20			10-12
SET ZZ BASE B-C-CCTB-L20 3.4 (Conical cemented titan base non hex L20 and thick screw)	BSAD0831L20	20			10-12
SET ZZ BASE B-C-CCTB-L20 4.1 (Conical cemented titan base non hex L20 and thick screw)	BSAD0832L20	20			10-12
SET ZZ BASE B-C-CCTB-L20 5.0 (Conical cemented titan base non hex L20 and thick screw)	BSAD0833L20	20			10-12
SET ZZ BASE B-C-CCTB-L20 6.0 (Conical cemented titan base non hex L20 and thick screw)	BSAD0834L20	20			10-12
SET ZZ BASE B-C-CCTB-L20-G 3.4 (Conical cemented titan base non hex L20 and screw Gold)	BSCD0831L20	20			10-12
SET ZZ BASE B-C-CCTB-L20-G 4.1 (Conical cemented titan base non hex L20 and screw Gold)	BSCD0832L20	20			10-12
SET ZZ BASE B-C-CCTB-L20-G 5.0 (Conical cemented titan base non hex L20 and screw Gold)	BSCD0833L20	20			10-12
SET ZZ BASE B-C-CCTB-L20-G 6.0 (Conical cemented titan base non hex L20 and screw Gold)	BSCD0834L20	20			10-12
SET ZZ BASE B-C-CCTB-L20-K85 3.4 (Conical cemented titan base non hex L20 K85 and thick screw)	BSAD0831L20K85	20			10-12
SET ZZ BASE B-C-CCTB-L30 3.4 (Conical cemented titan base non hex L30 and thick screw)	BSAD0831L30	20			10-12
SET ZZ BASE B-C-CCTB-L30 4.1 (Conical cemented titan base non hex L30 and thick screw)	BSAD0832L30	20			10-12
SET ZZ BASE B-C-CCTB-L30 5.0 (Conical cemented titan base non hex L30 and thick screw)	BSAD0833L30	20			10-12
SET ZZ BASE B-C-CCTB-L30 6.0 (Conical cemented titan base non hex L30 and thick screw)	BSAD0834L30	20			10-12
SET ZZ BASE B-C-CCTB-L30-G 3.4 (Conical cemented titan base non hex L30 and screw Gold)	BSCD0831L30	20			10-12
SET ZZ BASE B-C-CCTB-L30-G 4.1 (Conical cemented titan base non hex L30 and screw Gold)	BSCD0832L30	20			10-12
SET ZZ BASE B-C-CCTB-L30-G 5.0 (Conical cemented titan base non hex L30 and screw Gold)	BSCD0833L30	20			10-12
SET ZZ BASE B-C-CCTB-L30-G 6.0 (Conical cemented titan base non hex L30 and screw Gold)	BSCD0834L30	20			10-12

All information is subject to change. Errors and omissions excepted. Version: 21.09.2023

BIOMET 3I™ CERTAIN® IMPLANT SYSTEM

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
SET ZZ BASE B-C-CCTB-L30-K85 3.4 (Conical cemented titan base non hex L30 K85 and thick screw)	BSAD0831L30K85	20			10-12
SET ZZ BASE B-C-CCTB-L40 3.4 (Conical cemented titan base non hex L40 and thick screw)	BSAD0831L40	20			10-12
SET ZZ BASE B-C-CCTB-L40 4.1 (Conical cemented titan base non hex L40 and thick screw)	BSAD0832L40	20			10-12
SET ZZ BASE B-C-CCTB-L40 5.0 (Conical cemented titan base non hex L40 and thick screw)	BSAD0833L40	20			10-12
SET ZZ BASE B-C-CCTB-L40 6.0 (Conical cemented titan base non hex L40 and thick screw)	BSAD0834L40	20			10-12
SET ZZ BASE B-C-CCTB-L40-G 3.4 (Conical cemented titan base non hex L40 and screw Gold)	BSCD0831L40	20			10-12
SET ZZ BASE B-C-CCTB-L40-G 4.1 (Conical cemented titan base non hex L40 and screw Gold)	BSCD0832L40	20			10-12
SET ZZ BASE B-C-CCTB-L40-G 5.0 (Conical cemented titan base non hex L40 and screw Gold)	BSCD0833L40	20			10-12
SET ZZ BASE B-C-CCTB-L40-G 6.0 (Conical cemented titan base non hex L40 and screw Gold)	BSCD0834L40	20			10-12
SET ZZ BASE B-C-CCTB-L40-K85 3.4 (Conical cemented titan base non hex L40 K85 and thick screw)	BSAD0831L40K85	20			10-12
SET ZZ BASE B-C-ICOT 3.4 (Impression coping open tray non hex and Screw and Driver)	BSAQ0861	15			10-12
SET ZZ BASE B-C-ICOT 4.1 (Impression coping open tray non hex and Screw and Driver)	BSAQ0862	15			10-12
SET ZZ BASE B-C-ICOT 5.0 (Impression coping open tray non hex and Screw and Driver)	BSAQ0863	15			10-12
SET ZZ BASE B-C-ICOT 6.0 (Impression coping open tray non hex and Screw and Driver)	BSAQ0864	15			10-12
SET ZZ BASE B-C-ICOTH 3.4 (Impression coping open tray hex and Screw and Driver)	BSAQ0841	15			10-12
SET ZZ BASE B-C-ICOTH 4.1 (Impression coping open tray hex and Screw and Driver)	BSAQ0842	15			10-12
SET ZZ BASE B-C-ICOTH 5.0 (Impression coping open tray hex and Screw and Driver)	BSAQ0843	15			10-12
SET ZZ BASE B-C-ICOTH 6.0 (Impression coping open tray hex and Screw and Driver)	BSAQ0844	15			10-12
SET ZZ BASE B-C-NTB 3.4 (Narrow titan base non hex and thick screw)	BSAP0831	20			10-12
SET ZZ BASE B-C-NTB 4.1 (Narrow titan base non hex and thick screw)	BSAP0832	20			10-12
SET ZZ BASE B-C-NTBH 3.4 (Narrow titan base hex and screw)	BSAP0881	20			10-12
SET ZZ BASE B-C-NTBH 4.1 (Narrow titan base hex and screw)	BSAP0882	20			10-12
SET ZZ BASE B-C-PCTBH 3.4 (Parallel cemented titan base hex and screw)	BSAB0881	20			10-12
SET ZZ BASE B-C-PCTBH 4.1 (Parallel cemented titan base hex and screw)	BSAB0882	20			10-12
SET ZZ BASE B-C-PCTBH 5.0 (Parallel cemented titan base hex and screw)	BSAB0883	20			10-12
SET ZZ BASE B-C-PCTBH 6.0 (Parallel cemented titan base hex and screw)	BSAB0884	20			10-12
SET ZZ BASE B-C-PCTBH-G 3.4 (Parallel cemented titan base hex and screw Gold)	BSCB0881	20			10-12
SET ZZ BASE B-C-PCTBH-G 4.1 (Parallel cemented titan base hex and screw Gold)	BSCB0882	20			10-12
SET ZZ BASE B-C-PCTBH-G 5.0 (Parallel cemented titan base hex and screw Gold)	BSCB0883	20			10-12

All information is subject to change. Errors and omissions excepted. Version: 21.09.2023

BIOMET 3I™ CERTAIN® IMPLANT SYSTEM

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE B-C-PCTBH-G 6.0 (Parallel cemented titan base hex and screw Gold)	BSCB0884	20			10-12
SET ZZ BASE B-C-PCTBH-K85 3.4 (Parallel cemented titan base hex K85 and screw)	BSAB0881K85	20			10-12
SET ZZ BASE B-C-PCTBH-K85 4.1 (Parallel cemented titan base hex K85 and screw)	BSAB0882K85	20			10-12
SET ZZ BASE B-C-PCTBH-K85 5.0 (Parallel cemented titan base hex K85 and screw)	BSAB0883K85	20			10-12
SET ZZ BASE B-C-PCTBH-K85 6.0 (Parallel cemented titan base hex K85 and screw)	BSAB0884K85	20			10-12
SET ZZ BASE B-C-PCTBH-L10 3.4 (Parallel cemented titan base hex L10 and screw)	BSAB0891L10	20			10-12
SET ZZ BASE B-C-PCTBH-L10 4.1 (Parallel cemented titan base hex L10 and screw)	BSAB0892L10	20			10-12
SET ZZ BASE B-C-PCTBH-L10 5.0 (Parallel cemented titan base hex L10 and screw)	BSAB0893L10	20			10-12
SET ZZ BASE B-C-PCTBH-L10 6.0 (Parallel cemented titan base hex L10 and screw)	BSAB0894L10	20			10-12
SET ZZ BASE B-C-PCTBH-L10-G 3.4 (Parallel cemented titan base hex L10 and screw Gold)	BSCB0891L10	20			10-12
SET ZZ BASE B-C-PCTBH-L10-G 4.1 (Parallel cemented titan base hex L10 and screw Gold)	BSCB0892L10	20			10-12
SET ZZ BASE B-C-PCTBH-L10-G 5.0 (Parallel cemented titan base hex L10 and screw Gold)	BSCB0893L10	20			10-12
SET ZZ BASE B-C-PCTBH-L10-G 6.0 (Parallel cemented titan base hex L10 and screw Gold)	BSCB0894L10	20			10-12
SET ZZ BASE B-C-PCTBH-L20 3.4 (Parallel cemented titan base hex L20 and screw)	BSAB0891L20	20			10-12
SET ZZ BASE B-C-PCTBH-L20 4.1 (Parallel cemented titan base hex L20 and screw)	BSAB0892L20	20			10-12
SET ZZ BASE B-C-PCTBH-L20 5.0 (Parallel cemented titan base hex L20 and screw)	BSAB0893L20	20			10-12
SET ZZ BASE B-C-PCTBH-L20 6.0 (Parallel cemented titan base hex L20 and screw)	BSAB0894L20	20			10-12
SET ZZ BASE B-C-PCTBH-L20-G 3.4 (Parallel cemented titan base hex L20 and screw Gold)	BSCB0891L20	20			10-12
SET ZZ BASE B-C-PCTBH-L20-G 4.1 (Parallel cemented titan base hex L20 and screw Gold)	BSCB0892L20	20			10-12
SET ZZ BASE B-C-PCTBH-L20-G 5.0 (Parallel cemented titan base hex L20 and screw Gold)	BSCB0893L20	20			10-12
SET ZZ BASE B-C-PCTBH-L20-G 6.0 (Parallel cemented titan base hex L20 and screw Gold)	BSCB0894L20	20			10-12
SET ZZ BASE B-C-PCTBH-L30 3.4 (Parallel cemented titan base hex L30 and screw)	BSAB0891L30	20			10-12
SET ZZ BASE B-C-PCTBH-L30 4.1 (Parallel cemented titan base hex L30 and screw)	BSAB0892L30	20			10-12
SET ZZ BASE B-C-PCTBH-L30 5.0 (Parallel cemented titan base hex L30 and screw)	BSAB0893L30	20			10-12
SET ZZ BASE B-C-PCTBH-L30 6.0 (Parallel cemented titan base hex L30 and screw)	BSAB0894L30	20			10-12
SET ZZ BASE B-C-PCTBH-L30-G 3.4 (Parallel cemented titan base hex L30 and screw Gold)	BSCB0891L30	20			10-12
SET ZZ BASE B-C-PCTBH-L30-G 4.1 (Parallel cemented titan base hex L30 and screw Gold)	BSCB0892L30	20			10-12
SET ZZ BASE B-C-PCTBH-L30-G 5.0 (Parallel cemented titan base hex L30 and screw Gold)	BSCB0893L30	20			10-12
SET ZZ BASE B-C-PCTBH-L30-G 6.0 (Parallel cemented titan base hex L30 and screw Gold)	BSCB0894L30	20			10-12

BIOMET 3I™ CERTAIN® IMPLANT SYSTEM

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE B-C-PCTBH-L40 3.4 (Parallel cemented titan base hex L40 and screw)	BSAB0891L40	20			10-12
SET ZZ BASE B-C-PCTBH-L40 4.1 (Parallel cemented titan base hex L40 and screw)	BSAB0892L40	20			10-12
SET ZZ BASE B-C-PCTBH-L40 5.0 (Parallel cemented titan base hex L40 and screw)	BSAB0893L40	20			10-12
SET ZZ BASE B-C-PCTBH-L40 6.0 (Parallel cemented titan base hex L40 and screw)	BSAB0894L40	20			10-12
SET ZZ BASE B-C-PCTBH-L40-G 3.4 (Parallel cemented titan base hex L40 and screw Gold)	BSCB0891L40	20			10-12
SET ZZ BASE B-C-PCTBH-L40-G 4.1 (Parallel cemented titan base hex L40 and screw Gold)	BSCB0892L40	20			10-12
SET ZZ BASE B-C-PCTBH-L40-G 5.0 (Parallel cemented titan base hex L40 and screw Gold)	BSCB0893L40	20			10-12
SET ZZ BASE B-C-PCTBH-L40-G 6.0 (Parallel cemented titan base hex L40 and screw Gold)	BSCB0894L40	20			10-12
SET ZZ BASE B-C-SA 3.4 (Scananalog and screw)	ZBDE0831	15			10-12
SET ZZ BASE B-C-SA 4.1 (Scananalog and screw)	ZBDE0832	15			10-12
SET ZZ BASE B-C-SA 5.0 (Scananalog and screw)	ZBDE0833	15			10-12
SET ZZ BASE B-C-SA 6.0 (Scananalog and screw)	ZBDE0834	15			10-12
SET ZZ BASE B-C-SM 3.4 (Scanmarker and screw)	ZBAD0831	20			10-12
SET ZZ BASE B-C-SM 4.1 (Scanmarker and screw)	ZBAD0832	20			10-12
SET ZZ BASE B-C-SM 5.0 (Scanmarker and screw)	ZBAD0833	20			10-12
SET ZZ BASE B-C-SM 6.0 (Scanmarker and screw)	ZBAD0834	20			10-12
SET ZZ BASE B-C-WSM 5.0 (White Scanmarker and screw)	ZBCD0833	20			10-12
SET ZZ BASE B-C-WSM 6.0 (White Scanmarker and screw)	ZBCD0834	20			10-12
SET ZZ BASE B-C-ZZMUA17 3.4 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and implant)	BSAI0831	20			10-12
SET ZZ BASE B-C-ZZMUA17 4.1 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and implant)	BSAI0832	20			10-12
SET ZZ BASE B-C-ZZMUA17 4.1 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and implant)	BSAI0842	20			10-12
SET ZZ BASE B-C-ZZMUA17 5.0 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and implant)	BSAI0833	20			10-12
SET ZZ BASE B-C-ZZMUA17 5.0 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and implant)	BSAI0843	20			10-12
SET ZZ BASE B-C-ZZMUA17-L10 3.4 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10 and implant)	BSAI0831L10	20			10-12
SET ZZ BASE B-C-ZZMUA17-L10 4.1 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10 and implant)	BSAI0832L10	20			10-12
SET ZZ BASE B-C-ZZMUA17-L10 4.1 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 L10 and implant)	BSAI0842L10	20			10-12
SET ZZ BASE B-C-ZZMUA17-L10 5.0 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10 and implant)	BSAI0833L10	20			10-12
SET ZZ BASE B-C-ZZMUA17-L10 5.0 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 L10 and implant)	BSAI0843L10	20			10-12
SET ZZ BASE B-C-ZZMUA-L15 3.4 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH0831L15	25	20		10-12

BIOMET 3I™ CERTAIN® IMPLANT SYSTEM

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE B-C-ZZMUA-L15 4.1 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH0832L15	25	20		10-12
SET ZZ BASE B-C-ZZMUA-L15 5.0 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH0833L15	25	20		10-12
SET ZZ BASE B-C-ZZMUA-L15 6.0 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH0834L15	25	20		10-12
SET ZZ BASE B-C-ZZMUA-L25 3.4 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH0831L25	25	20		10-12
SET ZZ BASE B-C-ZZMUA-L25 4.1 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH0832L25	25	20		10-12
SET ZZ BASE B-C-ZZMUA-L25 5.0 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH0833L25	25	20		10-12
SET ZZ BASE B-C-ZZMUA-L25 6.0 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH0834L25	25	20		10-12
SET ZZ BASE B-C-ZZMUA-L35 3.4 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH0831L35	25	20		10-12
SET ZZ BASE B-C-ZZMUA-L35 4.1 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH0832L35	25	20		10-12
SET ZZ BASE B-C-ZZMUA-L35 5.0 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH0833L35	25	20		10-12
SET ZZ BASE B-C-ZZMUA-L35 6.0 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH0834L35	25	20		10-12
SET ZZ BASE B-C-ZZMUA-L45 3.4 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH0831L45	25	20		10-12
SET ZZ BASE B-C-ZZMUA-L45 4.1 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH0832L45	25	20		10-12
SET ZZ BASE B-C-ZZMUA-L45 5.0 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH0833L45	25	20		10-12
SET ZZ BASE B-C-ZZMUA-L45 6.0 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH0834L45	25	20		10-12
ZZ BASE B-C-ASCSC 3.4 (Angled screw-channel Screw)	BSKE0801	20			10-12
ZZ BASE B-C-ASCSC 4.1 (Angled screw-channel Screw)	BSKE0802	20			10-12
ZZ BASE B-C-ASCTSC 3.4 (Angled screw-channel Thick Screw)	BSKE0811	20			10-12
ZZ BASE B-C-ASCTSC 4.1 (Angled screw-channel Thick Screw)	BSKE0812	20			10-12
ZZ BASE B-C-BSC 3.4 (Black Screw)	BSIE0801	20			10-12
ZZ BASE B-C-BSC 4.1 (Black Screw)	BSIE0802	20			10-12
ZZ BASE B-C-BSC 5.0 (Black Screw)	BSIE0803	20			10-12
ZZ BASE B-C-BSC 6.0 (Black Screw)	BSIE0804	20			10-12
ZZ BASE B-C-BTSC 3.4 (Black Thick screw)	BSIE0811	20			10-12
ZZ BASE B-C-BTSC 4.1 (Black Thick screw)	BSIE0812	20			10-12
ZZ BASE B-C-BTSC 5.0 (Black Thick screw)	BSIE0813	20			10-12
ZZ BASE B-C-BTSC 6.0 (Black Thick screw)	BSIE0814	20			10-12
ZZ BASE B-C-BZSC 3.4 (Black Zirconia screw)	BSIE0831	20			10-12
ZZ BASE B-C-BZSC 4.1 (Black Zirconia screw)	BSIE0832	20			10-12

BIOMET 3I™ CERTAIN® IMPLANT SYSTEM

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE B-C-BZSC 5.0 (Black Zirconia screw)	BSIE0833	20			10-12
ZZ BASE B-C-BZSC 6.0 (Black Zirconia screw)	BSIE0834	20			10-12
ZZ BASE B-C-ICSC 3.4 (Impression coping Screw)	BSFE0801	15			10-12
ZZ BASE B-C-ICSC 4.1 (Impression coping Screw)	BSFE0802	15			10-12
ZZ BASE B-C-ICSC 5.0 (Impression coping Screw)	BSFE0803	15			10-12
ZZ BASE B-C-ICSC 6.0 (Impression coping Screw)	BSFE0804	15			10-12
ZZ BASE B-C-LBSC 3.4 (Lab Screw)	BSHE0801	20			10-12
ZZ BASE B-C-LBSC 4.1 (Lab Screw)	BSHE0802	20			10-12
ZZ BASE B-C-LBSC 5.0 (Lab Screw)	BSHE0803	20			10-12
ZZ BASE B-C-LBSC 6.0 (Lab Screw)	BSHE0804	20			10-12
ZZ BASE B-C-LBTSC 3.4 (Lab Thick screw)	BSHE0811	20			10-12
ZZ BASE B-C-LBTSC 4.1 (Lab Thick screw)	BSHE0812	20			10-12
ZZ BASE B-C-LBTSC 5.0 (Lab Thick screw)	BSHE0813	20			10-12
ZZ BASE B-C-LBTSC 6.0 (Lab Thick screw)	BSHE0814	20			10-12
ZZ BASE B-C-LBZSC 3.4 (Lab Zirconia screw)	BSHE0831	20			10-12
ZZ BASE B-C-LBZSC 4.1 (Lab Zirconia screw)	BSHE0832	20			10-12
ZZ BASE B-C-LBZSC 5.0 (Lab Zirconia screw)	BSHE0833	20			10-12
ZZ BASE B-C-LBZSC 6.0 (Lab Zirconia screw)	BSHE0834	20			10-12
ZZ BASE B-C-MI 3.4 (Model implant)	MIAG0801				10-12
ZZ BASE B-C-MI 4.1 (Model implant)	MIAG0802				10-12
ZZ BASE B-C-MI 5.0 (Model implant)	MIAG0803				10-12
ZZ BASE B-C-MI 6.0 (Model implant)	MIAG0804				10-12
ZZ BASE B-C-SC 3.4 (Screw)	BSAE0801	20			10-12
ZZ BASE B-C-SC 4.1 (Screw)	BSAE0802	20			10-12
ZZ BASE B-C-SC 5.0 (Screw)	BSAE0803	20			10-12
ZZ BASE B-C-SC 6.0 (Screw)	BSAE0804	20			10-12
ZZ BASE B-C-SC-G 3.4 (Screw Gold)	BSCE0801	20			10-12
ZZ BASE B-C-SC-G 4.1 (Screw Gold)	BSCE0802	20			10-12
ZZ BASE B-C-SC-G 5.0 (Screw Gold)	BSCE0803	20			10-12

All information is subject to change. Errors and omissions excepted. Version: 21.09.2023

BIOMET 3I™ CERTAIN® IMPLANT SYSTEM

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE B-C-SC-G 6.0 (Screw Gold)	BSCE0804	20			10-12
ZZ BASE B-C-TSC 3.4 (Thick screw)	BSAE0811	20			10-12
ZZ BASE B-C-TSC 4.1 (Thick screw)	BSAE0812	20			10-12
ZZ BASE B-C-TSC 5.0 (Thick screw)	BSAE0813	20			10-12
ZZ BASE B-C-TSC 6.0 (Thick screw)	BSAE0814	20			10-12
ZZ BASE B-C-TSC-G 3.4 (Thick screw Gold)	BSCE0811	20			10-12
ZZ BASE B-C-TSC-G 4.1 (Thick screw Gold)	BSCE0812	20			10-12
ZZ BASE B-C-TSC-G 5.0 (Thick screw Gold)	BSCE0813	20			10-12
ZZ BASE B-C-TSC-G 6.0 (Thick screw Gold)	BSCE0814	20			10-12
ZZ BASE B-C-WMSM-L10 3.4 (White Metal Scanmarker L10 with Screw)	ZBED0831L10	15			10-12
ZZ BASE B-C-WMSM-L10 4.1 (White Metal Scanmarker L10 with Screw)	ZBED0832L10	15			10-12
ZZ BASE B-C-WMSM-L10 5.0 (White Metal Scanmarker L10 with Screw)	ZBED0833L10	15			10-12
ZZ BASE B-C-WMSM-L10 6.0 (White Metal Scanmarker L10 with Screw)	ZBED0834L10	15			10-12
ZZ BASE B-C-WMSM-L12 3.4 (White Metal Scanmarker L12 with Screw)	ZBED0831L12	15			10-12
ZZ BASE B-C-WMSM-L12 4.1 (White Metal Scanmarker L12 with Screw)	ZBED0832L12	15			10-12
ZZ BASE B-C-WMSM-L12 5.0 (White Metal Scanmarker L12 with Screw)	ZBED0833L12	15			10-12
ZZ BASE B-C-WMSM-L12 6.0 (White Metal Scanmarker L12 with Screw)	ZBED0834L12	15			10-12
ZZ BASE B-C-WMSM-L15 3.4 (White Metal Scanmarker L15 with Screw)	ZBED0831L15	15			10-12
ZZ BASE B-C-WMSM-L15 4.1 (White Metal Scanmarker L15 with Screw)	ZBED0832L15	15			10-12
ZZ BASE B-C-WMSM-L15 5.0 (White Metal Scanmarker L15 with Screw)	ZBED0833L15	15			10-12
ZZ BASE B-C-WMSM-L15 6.0 (White Metal Scanmarker L15 with Screw)	ZBED0834L15	15			10-12
ZZ BASE B-C-ZSC 3.4 (Zirconia screw)	BSAE0831	20			10-12
ZZ BASE B-C-ZSC 4.1 (Zirconia screw)	BSAE0832	20			10-12
ZZ BASE B-C-ZSC 5.0 (Zirconia screw)	BSAE0833	20			10-12
ZZ BASE B-C-ZSC 6.0 (Zirconia screw)	BSAE0834	20			10-12
ZZ BASE B-C-ZSC-G 3.4 (Zirconia screw Gold)	BSCE0831	20			10-12
ZZ BASE B-C-ZSC-G 4.1 (Zirconia screw Gold)	BSCE0832	20			10-12
ZZ BASE B-C-ZSC-G 5.0 (Zirconia screw Gold)	BSCE0833	20			10-12

BIOMET 3I™ EXTERNAL CONNECTION IMPLANT SYSTEM (OSSEOTITE®)

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE B-O-ABCH 3.4 (Raw-Abutment® customizable hex and screw)	BSAF0991	35			10-12
SET ZZ BASE B-O-ABCH 4.1 (Raw-Abutment® customizable hex and screw)	BSAF0992	35			10-12
SET ZZ BASE B-O-ABCH 5.0 (Raw-Abutment® customizable hex and screw)	BSAF0993	35			10-12
SET ZZ BASE B-O-ABCH 6.0 (Raw-Abutment® customizable hex and screw)	BSAF0994	35			10-12
SET ZZ BASE B-O-ASCTBH-K80 3.4 (Angled screw-channel titan base hex K80 and screw)	BSAT0981K80	20			10-12
SET ZZ BASE B-O-ASCTBH-K80 4.1 (Angled screw-channel titan base hex K80 and screw)	BSAT0982K80	20			10-12
SET ZZ BASE B-O-ASCTB-K80 3.4 (Angled screw-channel titan base non hex K80 and screw)	BSAT0931K80	20			10-12
SET ZZ BASE B-O-ASCTB-K80 4.1 (Angled screw-channel titan base non hex K80 and screw)	BSAT0932K80	20			10-12
SET ZZ BASE B-O-CCTB 3.4 (Conical cemented titan base non hex and screw)	BSAD0931	35			10-12
SET ZZ BASE B-O-CCTB 4.1 (Conical cemented titan base non hex and screw)	BSAD0932	35			10-12
SET ZZ BASE B-O-CCTB 5.0 (Conical cemented titan base non hex and screw)	BSAD0933	35			10-12
SET ZZ BASE B-O-CCTB 6.0 (Conical cemented titan base non hex and screw)	BSAD0934	35			10-12
SET ZZ BASE B-O-CCTB-K85 3.4 (Conical cemented titan base non hex K85 and screw)	BSAD0931K85	35			10-12
SET ZZ BASE B-O-CCTB-K85 4.1 (Conical cemented titan base non hex K85 and screw)	BSAD0932K85	35			10-12
SET ZZ BASE B-O-CCTB-K85 5.0 (Conical cemented titan base non hex K85 and screw)	BSAD0933K85	35			10-12
SET ZZ BASE B-O-CCTB-K85 6.0 (Conical cemented titan base non hex K85 and screw)	BSAD0934K85	35			10-12
SET ZZ BASE B-O-CCTB-L10 3.4 (Conical cemented titan base non hex L10 and screw)	BSAD0931L10	35			10-12
SET ZZ BASE B-O-CCTB-L10 4.1 (Conical cemented titan base non hex L10 and screw)	BSAD0932L10	35			10-12
SET ZZ BASE B-O-CCTB-L10 5.0 (Conical cemented titan base non hex L10 and screw)	BSAD0933L10	35			10-12
SET ZZ BASE B-O-CCTB-L10 6.0 (Conical cemented titan base non hex L10 and screw)	BSAD0934L10	35			10-12
SET ZZ BASE B-O-CCTB-L20 3.4 (Conical cemented titan base non hex L20 and screw)	BSAD0931L20	35			10-12
SET ZZ BASE B-O-CCTB-L20 4.1 (Conical cemented titan base non hex L20 and screw)	BSAD0932L20	35			10-12
SET ZZ BASE B-O-CCTB-L20 5.0 (Conical cemented titan base non hex L20 and screw)	BSAD0933L20	35			10-12
SET ZZ BASE B-O-CCTB-L20 6.0 (Conical cemented titan base non hex L20 and screw)	BSAD0934L20	35			10-12
SET ZZ BASE B-O-CCTB-L30 3.4 (Conical cemented titan base non hex L30 and screw)	BSAD0931L30	35			10-12
SET ZZ BASE B-O-CCTB-L30 4.1 (Conical cemented titan base non hex L30 and screw)	BSAD0932L30	35			10-12
SET ZZ BASE B-O-CCTB-L30 5.0 (Conical cemented titan base non hex L30 and screw)	BSAD0933L30	35			10-12
SET ZZ BASE B-O-CCTB-L30 6.0 (Conical cemented titan base non hex L30 and screw)	BSAD0934L30	35			10-12
SET ZZ BASE B-O-CCTB-L40 3.4 (Conical cemented titan base non hex L40 and screw)	BSAD0931L40	35			10-12

BIOMET 3I™ EXTERNAL CONNECTION IMPLANT SYSTEM (OSSEOTITE®)

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
SET ZZ BASE B-O-CCTB-L40 4.1 (Conical cemented titan base non hex L40 and screw)	BSAD0932L40	35			10-12
SET ZZ BASE B-O-CCTB-L40 5.0 (Conical cemented titan base non hex L40 and screw)	BSAD0933L40	35			10-12
SET ZZ BASE B-O-CCTB-L40 6.0 (Conical cemented titan base non hex L40 and screw)	BSAD0934L40	35			10-12
SET ZZ BASE B-O-ICOT 3.4 (Impression coping open tray non hex and Screw and Driver)	BSAQ0961	15			10-12
SET ZZ BASE B-O-ICOT 4.1 (Impression coping open tray non hex and Screw and Driver)	BSAQ0962	15			10-12
SET ZZ BASE B-O-ICOT 5.0 (Impression coping open tray non hex and Screw and Driver)	BSAQ0963	15			10-12
SET ZZ BASE B-O-ICOT 6.0 (Impression coping open tray non hex and Screw and Driver)	BSAQ0964	15			10-12
SET ZZ BASE B-O-ICOTH 3.4 (Impression coping open tray hex and Screw and Driver)	BSAQ0941	15			10-12
SET ZZ BASE B-O-ICOTH 4.1 (Impression coping open tray hex and Screw and Driver)	BSAQ0942	15			10-12
SET ZZ BASE B-O-ICOTH 5.0 (Impression coping open tray hex and Screw and Driver)	BSAQ0943	15			10-12
SET ZZ BASE B-O-ICOTH 6.0 (Impression coping open tray hex and Screw and Driver)	BSAQ0944	15			10-12
SET ZZ BASE B-O-NTB 3.4 (Narrow titan base non hex and screw)	BSAP0931	35			10-12
SET ZZ BASE B-O-NTB 4.1 (Narrow titan base non hex and screw)	BSAP0932	35			10-12
SET ZZ BASE B-O-NTBH 3.4 (Narrow titan base hex and screw)	BSAP0981	35			10-12
SET ZZ BASE B-O-NTBH 4.1 (Narrow titan base hex and screw)	BSAP0982	35			10-12
SET ZZ BASE B-O-PCTBH 3.4 (Parallel cemented titan base hex and screw)	BSAB0981	35			10-12
SET ZZ BASE B-O-PCTBH 4.1 (Parallel cemented titan base hex and screw)	BSAB0982	35			10-12
SET ZZ BASE B-O-PCTBH 5.0 (Parallel cemented titan base hex and screw)	BSAB0983	35			10-12
SET ZZ BASE B-O-PCTBH 6.0 (Parallel cemented titan base hex and screw)	BSAB0984	35			10-12
SET ZZ BASE B-O-PCTBH-K85 3.4 (Parallel cemented titan base hex K85 and screw)	BSAB0981K85	35			10-12
SET ZZ BASE B-O-PCTBH-K85 4.1 (Parallel cemented titan base hex K85 and screw)	BSAB0982K85	35			10-12
SET ZZ BASE B-O-PCTBH-K85 5.0 (Parallel cemented titan base hex K85 and screw)	BSAB0983K85	35			10-12
SET ZZ BASE B-O-PCTBH-K85 6.0 (Parallel cemented titan base hex K85 and screw)	BSAB0984K85	35			10-12
SET ZZ BASE B-O-PCTBH-L10 3.4 (Parallel cemented titan base hex L10 and screw)	BSAB0991L10	35			10-12
SET ZZ BASE B-O-PCTBH-L10 4.1 (Parallel cemented titan base hex L10 and screw)	BSAB0992L10	35			10-12
SET ZZ BASE B-O-PCTBH-L10 5.0 (Parallel cemented titan base hex L10 and screw)	BSAB0993L10	35			10-12
SET ZZ BASE B-O-PCTBH-L10 6.0 (Parallel cemented titan base hex L10 and screw)	BSAB0994L10	35			10-12
SET ZZ BASE B-O-PCTBH-L20 3.4 (Parallel cemented titan base hex L20 and screw)	BSAB0991L20	35			10-12
SET ZZ BASE B-O-PCTBH-L20 4.1 (Parallel cemented titan base hex L20 and screw)	BSAB0992L20	35			10-12

All information is subject to change. Errors and omissions excepted. Version: 21.09.2023

BIOMET 3I™ EXTERNAL CONNECTION IMPLANT SYSTEM (OSSEOTITE®)

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE B-O-PCTBH-L20 5.0 (Parallel cemented titan base hex L20 and screw)	BSAB0993L20	35			10-12
SET ZZ BASE B-O-PCTBH-L20 6.0 (Parallel cemented titan base hex L20 and screw)	BSAB0994L20	35			10-12
SET ZZ BASE B-O-PCTBH-L30 3.4 (Parallel cemented titan base hex L30 and screw)	BSAB0991L30	35			10-12
SET ZZ BASE B-O-PCTBH-L30 4.1 (Parallel cemented titan base hex L30 and screw)	BSAB0992L30	35			10-12
SET ZZ BASE B-O-PCTBH-L30 5.0 (Parallel cemented titan base hex L30 and screw)	BSAB0993L30	35			10-12
SET ZZ BASE B-O-PCTBH-L30 6.0 (Parallel cemented titan base hex L30 and screw)	BSAB0994L30	35			10-12
SET ZZ BASE B-O-PCTBH-L40 3.4 (Parallel cemented titan base hex L40 and screw)	BSAB0991L40	35			10-12
SET ZZ BASE B-O-PCTBH-L40 4.1 (Parallel cemented titan base hex L40 and screw)	BSAB0992L40	35			10-12
SET ZZ BASE B-O-PCTBH-L40 5.0 (Parallel cemented titan base hex L40 and screw)	BSAB0993L40	35			10-12
SET ZZ BASE B-O-PCTBH-L40 6.0 (Parallel cemented titan base hex L40 and screw)	BSAB0994L40	35			10-12
SET ZZ BASE B-O-SA 3.4 (Scananalog and screw)	ZBDE0931	15			10-12
SET ZZ BASE B-O-SA 4.1 (Scananalog and screw)	ZBDE0932	15			10-12
SET ZZ BASE B-O-SA 5.0 (Scananalog and screw)	ZBDE0933	15			10-12
SET ZZ BASE B-O-SA 6.0 (Scananalog and screw)	ZBDE0934	15			10-12
SET ZZ BASE B-O-SM 3.4 (Scanmarker and screw)	ZBAD0931	35			10-12
SET ZZ BASE B-O-SM 4.1 (Scanmarker and screw)	ZBAD0932	35			10-12
SET ZZ BASE B-O-SM 5.0 (Scanmarker and screw)	ZBAD0933	35			10-12
SET ZZ BASE B-O-SM 6.0 (Scanmarker and screw)	ZBAD0934	35			10-12
SET ZZ BASE B-O-WSM 3.4 (White Scanmarker and screw)	ZBCD0931	35			10-12
SET ZZ BASE B-O-WSM 4.1 (White Scanmarker and screw)	ZBCD0932	35			10-12
SET ZZ BASE B-O-WSM 5.0 (White Scanmarker and screw)	ZBCD0933	35			10-12
SET ZZ BASE B-O-WSM 6.0 (White Scanmarker and screw)	ZBCD0934	35			10-12
ZZ BASE B-O-ASCSC 3.4 (Angled screw-channel Screw)	BSKE0901	20			10-12
ZZ BASE B-O-ASCSC 4.1 (Angled screw-channel Screw)	BSKE0902	20			10-12
ZZ BASE B-O-ICSC 3.4 (Impression coping Screw)	BSFE0901	15			10-12
ZZ BASE B-O-ICSC 4.1 (Impression coping Screw)	BSFE0902	15			10-12
ZZ BASE B-O-ICSC 5.0 (Impression coping Screw)	BSFE0903	15			10-12
ZZ BASE B-O-ICSC 6.0 (Impression coping Screw)	BSFE0904	15			10-12
ZZ BASE B-O-MI 3.4 (Model implant)	MIAG0901				10-12

BIOMET 3I™ EXTERNAL CONNECTION IMPLANT SYSTEM (OSSEOTITE®)

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
ZZ BASE B-O-MI 4.1 (Model implant)	MIAG0902				10-12
ZZ BASE B-O-MI 5.0 (Model implant)	MIAG0903				10-12
ZZ BASE B-O-MI 6.0 (Model implant)	MIAG0904				10-12
ZZ BASE B-O-SC 3.4 (Screw)	BSAE0901	35			10-12
ZZ BASE B-O-SC 4.1 (Screw)	BSAE0902	35			10-12
ZZ BASE B-O-SC 5.0 (Screw)	BSAE0903	35			10-12
ZZ BASE B-O-SC 6.0 (Screw)	BSAE0904	35			10-12
ZZ BASE B-O-WMSM-L12 3.4 (White Metal Scanmarker L12 with Screw)	ZBED0931L12	15			10-12
ZZ BASE B-O-WMSM-L12 4.1 (White Metal Scanmarker L12 with Screw)	ZBED0932L12	15			10-12
ZZ BASE B-O-WMSM-L12 5.0 (White Metal Scanmarker L12 with Screw)	ZBED0933L12	15			10-12
ZZ BASE B-O-WMSM-L12 6.0 (White Metal Scanmarker L12 with Screw)	ZBED0934L12	15			10-12
ZZ BASE B-O-ZSC 3.4 (Zirconia screw)	BSAE0931	35			10-12
ZZ BASE B-O-ZSC 4.1 (Zirconia screw)	BSAE0932	35			10-12
ZZ BASE B-O-ZSC 5.0 (Zirconia screw)	BSAE0933	35			10-12
ZZ BASE B-O-ZSC 6.0 (Zirconia screw)	BSAE0934	35			10-12

BIOMET 3I™ LOW PROFILE ABUTMENT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
SET ZZ BASE B-LPA-ASCTB-K80 (Angled screw-channel titan base non hex K80 and screw)	BSAT6431K80	20			10-12
SET ZZ BASE B-LPA-CCTB (Conical cemented titan base non hex and screw)	BSAD6431	20			10-12
SET ZZ BASE B-LPA-CCTB (Conical cemented titan base non hex and screw)	BSAD6441	20			10-12
SET ZZ BASE B-LPA-CCTB-K85 (Conical cemented titan base non hex K85 and screw)	BSAD6431K85	20			10-12
SET ZZ BASE B-LPA-CCTB-L20 (Conical cemented titan base non hex L20 and screw)	BSAD6431L20	20			10-12
SET ZZ BASE B-LPA-CCTB-L30 (Conical cemented titan base non hex L30 and screw)	BSAD6431L30	20			10-12
SET ZZ BASE B-LPA-CCTB-L40 (Conical cemented titan base non hex L40 and screw)	BSAD6431L40	20			10-12
SET ZZ BASE B-LPA-ICOT (Impression coping open tray non hex and Screw and Driver)	BSAQ6461	15			10-12
SET ZZ BASE B-LPA-SA (Scananalog and screw)	ZBDE6431	15			10-12
SET ZZ BASE B-LPA-SM (Scanmarker and screw)	ZBAD6431	20			10-12
SET ZZ BASE B-LPA-WSM (White Scanmarker and screw)	ZBCD6431	20			10-12
ZZ BASE B-LPA-ICSC (Impression coping Screw)	BSFE6401	15			10-12
ZZ BASE B-LPA-MI (Model implant)	MIAG6401				10-12
ZZ BASE B-LPA-SC (Screw)	BSAE6401	20			10-12
ZZ BASE B-LPA-WMSM-L12 (White Metal Scanmarker L12 with Screw)	ZBED6431L12	15			10-12
ZZ BASE B-LPA-ZSC (Zirconia Screw)	BSAE6431	20			10-12

BIOTECH DENTAL KONTAKT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
BI-KO-EX-A15 3.6/4.2/4.8/5.4 (Titanium Base Extractor A15)	EXAB6202A15				10-12
SET ZZ BASE BI-KO-ABCH 3.0 (Raw-Abutment® customizable hex and screw)	BSAF6281	20			10-12
SET ZZ BASE BI-KO-ABCH 3.6/4.2/4.8/5.4 (Raw-Abutment® customizable hex and screw)	BSAF6282	20			10-12
SET ZZ BASE BI-KO-ABCH-D14 3.6/4.2/4.8/5.4 (Raw-Abutment® customizable hex D14 and screw)	BSAF6282D14	20			10-12
SET ZZ BASE BI-KO-ASCTBH-K80 3.0 (Angled screw-channel titan base hex K80 and screw)	BSAT6281K80	20			10-12
SET ZZ BASE BI-KO-ASCTBH-K80 3.6/4.2/4.8/5.4 (Angled screw-channel titan base hex K80 and screw)	BSAT6282K80	20			10-12
SET ZZ BASE BI-KO-ASCTB-K80 3.0 (Angled screw-channel titan base non hex K80 and screw)	BSAT6231K80	20			10-12
SET ZZ BASE BI-KO-ASCTB-K80 3.6/4.2/4.8/5.4 (Angled screw-channel titan base non hex K80 and screw)	BSAT6232K80	20			10-12
SET ZZ BASE BI-KO-CCTB 3.0 (Conical cemented titan base non hex and screw)	BSAD6231	20			10-12
SET ZZ BASE BI-KO-CCTB 3.6/4.2/4.8/5.4 (Conical cemented titan base non hex and screw)	BSAD6232	20			10-12
SET ZZ BASE BI-KO-CCTB 4.0 (Conical cemented titan base non hex and screw)	BSAD6233	20			10-12
SET ZZ BASE BI-KO-CCTB 4.9 (Conical cemented titan base non hex and screw)	BSAD6234	20			10-12
SET ZZ BASE BI-KO-CCTB-K85 3.0 (Conical cemented titan base non hex K85 and screw)	BSAD6231K85	20			10-12
SET ZZ BASE BI-KO-CCTB-K85 3.6/4.2/4.8/5.4 (Conical cemented titan base non hex K85 and screw)	BSAD6232K85	20			10-12
SET ZZ BASE BI-KO-CCTB-L15 3.6/4.2/4.8/5.4 (Conical cemented titan base non hex L15 and screw)	BSAD6232L15	20			10-12
SET ZZ BASE BI-KO-CCTB-L20 3.0 (Conical cemented titan base non hex L20 and screw)	BSAD6231L20	20			10-12
SET ZZ BASE BI-KO-CCTB-L20 3.6/4.2/4.8/5.4 (Conical cemented titan base non hex L20 and screw)	BSAD6232L20	20			10-12
SET ZZ BASE BI-KO-CCTB-L30 3.0 (Conical cemented titan base non hex L30 and screw)	BSAD6231L30	20			10-12
SET ZZ BASE BI-KO-CCTB-L30 3.6/4.2/4.8/5.4 (Conical cemented titan base non hex L30 and screw)	BSAD6232L30	20			10-12
SET ZZ BASE BI-KO-CCTB-L40 3.0 (Conical cemented titan base non hex L40 and screw)	BSAD6231L40	20			10-12
SET ZZ BASE BI-KO-CCTB-L40 3.6/4.2/4.8/5.4 (Conical cemented titan base non hex L40 and screw)	BSAD6232L40	20			10-12
SET ZZ BASE BI-KO-ICOT 3.0 (Impression coping open tray non hex and Screw and Driver)	BSAQ6261	15			10-12
SET ZZ BASE BI-KO-ICOT 3.6/4.2/4.8/5.4 (Impression coping open tray non hex and Screw and Driver)	BSAQ6262	15			10-12
SET ZZ BASE BI-KO-ICOTH 3.0 (Impression coping open tray hex and Screw and Driver)	BSAQ6241	15			10-12
SET ZZ BASE BI-KO-ICOTH 3.6/4.2/4.8/5.4 (Impression coping open tray hex and Screw and Driver)	BSAQ6242	15			10-12
SET ZZ BASE BI-KO-NTB 3.0 (Narrow titan base non hex and screw)	BSAP6231	20			10-12
SET ZZ BASE BI-KO-NTB 3.6/4.2/4.8/5.4 (Narrow titan base non hex and screw)	BSAP6232	20			10-12
SET ZZ BASE BI-KO-NTBH 3.0 (Narrow titan base hex and screw)	BSAP6281	20			10-12
SET ZZ BASE BI-KO-NTBH 3.6/4.2/4.8/5.4 (Narrow titan base hex and screw)	BSAP6282	20			10-12

All information is subject to change. Errors and omissions excepted. Version: 21.09.2023

BIOTECH DENTAL KONTAKT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE BI-KO-PCTBH 3.0 (Parallel cemented titan base hex and screw)	BSAB6281	20			10-12
SET ZZ BASE BI-KO-PCTBH 3.6/4.2/4.8/5.4 (Parallel cemented titan base hex and screw)	BSAB6282	20			10-12
SET ZZ BASE BI-KO-PCTBH-K85 3.0 (Parallel cemented titan base hex K85 and screw)	BSAB6281K85	20			10-12
SET ZZ BASE BI-KO-PCTBH-K85 3.6/4.2/4.8/5.4 (Parallel cemented titan base hex K85 and screw)	BSAB6282K85	20			10-12
SET ZZ BASE BI-KO-PCTBH-L15 3.6/4.2/4.8/5.4 (Parallel cemented titan base hex L15 and screw)	BSAB6292L15	20			10-12
SET ZZ BASE BI-KO-PCTBH-L20 3.0 (Parallel cemented titan base hex L20 and screw)	BSAB6291L20	20			10-12
SET ZZ BASE BI-KO-PCTBH-L20 3.6/4.2/4.8/5.4 (Parallel cemented titan base hex L20 and screw)	BSAB6292L20	20			10-12
SET ZZ BASE BI-KO-PCTBH-L30 3.0 (Parallel cemented titan base hex L30 and screw)	BSAB6291L30	20			10-12
SET ZZ BASE BI-KO-PCTBH-L30 3.6/4.2/4.8/5.4 (Parallel cemented titan base hex L30 and screw)	BSAB6292L30	20			10-12
SET ZZ BASE BI-KO-PCTBH-L40 3.0 (Parallel cemented titan base hex L40 and screw)	BSAB6291L40	20			10-12
SET ZZ BASE BI-KO-PCTBH-L40 3.6/4.2/4.8/5.4 (Parallel cemented titan base hex L40 and screw)	BSAB6292L40	20			10-12
SET ZZ BASE BI-KO-SA 3.0 (Scanalog and screw)	ZBDE6231	15			10-12
SET ZZ BASE BI-KO-SA 3.6/4.2/4.8/5.4 (Scanalog and screw)	ZBDE6232	15			10-12
SET ZZ BASE BI-KO-SA 4.0 (Scanalog and screw)	ZBDE6233	15			10-12
SET ZZ BASE BI-KO-SA 4.9 (Scanalog and screw)	ZBDE6234	15			10-12
SET ZZ BASE BI-KO-SM 3.0 (Scanmarker and screw)	ZBAD6231	20			10-12
SET ZZ BASE BI-KO-SM 3.6/4.2/4.8/5.4 (Scanmarker and screw)	ZBAD6232	20			10-12
SET ZZ BASE BI-KO-SM 4.0 (Scanmarker and screw)	ZBAD6233	20			10-12
SET ZZ BASE BI-KO-SM 4.9 (Scanmarker and screw)	ZBAD6234	20			10-12
SET ZZ BASE BI-KO-WSM 3.0 (White Scanmarker and screw)	ZBCD6231	20			10-12
SET ZZ BASE BI-KO-WSM 3.6/4.2/4.8/5.4 (White Scanmarker and screw)	ZBCD6232	20			10-12
SET ZZ BASE BI-KO-WSM 4.0 (White Scanmarker and screw)	ZBCD6233	20			10-12
SET ZZ BASE BI-KO-WSM 4.9 (White Scanmarker and screw)	ZBCD6234	20			10-12
SET ZZ BASE BI-KO-ZZMUA-L15 3.0 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH6231L15	25	20		10-12
SET ZZ BASE BI-KO-ZZMUA-L15 3.6/4.2/4.8/5.4 (Zirkonzahn Multi Unit Abutment non hex L15 an	BSAH6232L15	25	20		10-12
SET ZZ BASE BI-KO-ZZMUA-L25 3.0 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH6231L25	25	20		10-12
SET ZZ BASE BI-KO-ZZMUA-L25 3.6/4.2/4.8/5.4 (Zirkonzahn Multi Unit Abutment non hex L25 an	BSAH6232L25	25	20		10-12
SET ZZ BASE BI-KO-ZZMUA-L35 3.0 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH6231L35	25	20		10-12
SET ZZ BASE BI-KO-ZZMUA-L35 3.6/4.2/4.8/5.4 (Zirkonzahn Multi Unit Abutment non hex L35 an	BSAH6232L35	25	20		10-12

BIOTECH DENTAL KONTAKT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE BI-KO-ZZMUA-L45 3.0 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH6231L45	25	20		10-12
SET ZZ BASE BI-KO-ZZMUA-L45 3.6/4.2/4.8/5.4 (Zirkonzahn Multi Unit Abutment non hex L45 an	BSAH6232L45	25	20		10-12
ZZ BASE BI-KO-ASCSC 3.0 (Angled screw-channel Screw)	BSKE6201	20			10-12
ZZ BASE BI-KO-ASCSC 3.6/4.2/4.8/5.4 (Angled screw-channel Screw)	BSKE6202	20			10-12
ZZ BASE BI-KO-ICSC 3.0 (Impression coping Screw)	BSFE6201	15			10-12
ZZ BASE BI-KO-ICSC 3.6/4.2/4.8/5.4 (Impression coping Screw)	BSFE6202	15			10-12
ZZ BASE BI-KO-MI 3.0 (Model implant)	MIAG6201				10-12
ZZ BASE BI-KO-MI 3.6/4.2/4.8/5.4 (Model implant)	MIAG6202				10-12
ZZ BASE BI-KO-MI 4.0 (Model implant)	MIAG6203				10-12
ZZ BASE BI-KO-MI 4.9 (Model implant)	MIAG6204				10-12
ZZ BASE BI-KO-SC 3.0 (Screw)	BSAE6201	20			10-12
ZZ BASE BI-KO-SC 3.6/4.2/4.8/5.4 (Screw)	BSAE6202	20			10-12
ZZ BASE BI-KO-SC 4.0 (Screw)	BSAE6203	20			10-12
ZZ BASE BI-KO-SC 4.9 (Screw)	BSAE6204	20			10-12
ZZ BASE BI-KO-WMSM-L10 3.0 (White Metal Scanmarker L10 with Screw)	ZBED6231L10	15			10-12
ZZ BASE BI-KO-WMSM-L10 3.6/4.2/4.8/5.4 (White Metal Scanmarker L10 with Screw)	ZBED6232L10	15			10-12
ZZ BASE BI-KO-WMSM-L12 3.0 (White Metal Scanmarker L12 with Screw)	ZBED6231L12	15			10-12
ZZ BASE BI-KO-WMSM-L12 3.6/4.2/4.8/5.4 (White Metal Scanmarker L12 with Screw)	ZBED6232L12	15			10-12
ZZ BASE BI-KO-WMSM-L15 3.0 (White Metal Scanmarker L15 with Screw)	ZBED6231L15	15			10-12
ZZ BASE BI-KO-WMSM-L15 3.6/4.2/4.8/5.4 (White Metal Scanmarker L15 with Screw)	ZBED6232L15	15			10-12
ZZ BASE BI-KO-ZSC 3.0 (Zirconia screw)	BSAE6231	20			10-12
ZZ BASE BI-KO-ZSC 3.6/4.2/4.8/5.4 (Zirconia screw)	BSAE6232	20			10-12
ZZ BASE BI-KO-ZSC 4.0 (Zirconia screw)	BSAE6233	20			10-12
ZZ BASE BI-KO-ZSC 4.9 (Zirconia screw)	BSAE6234	20			10-12

BRAINBASE CORPORATION MYTIS ARROW IMPLANT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE BB-MA-ABCH B (Raw-Abutment® customizable hex and screw)	BSAF9481	30			10-12
SET ZZ BASE BB-MA-ABCH C (Raw-Abutment® customizable hex and screw)	BSAF9482	30			10-12
SET ZZ BASE BB-MA-ABCH E (Raw-Abutment® customizable hex and screw)	BSAF9483	30			10-12
SET ZZ BASE BB-MA-ABCH EW/BW (Raw-Abutment® customizable hex and screw)	BSAF9484	30			10-12
SET ZZ BASE BB-MA-CCTB B (Conical cemented titan base non hex and screw)	BSAD9431	30			10-12
SET ZZ BASE BB-MA-CCTB C (Conical cemented titan base non hex and screw)	BSAD9432	30			10-12
SET ZZ BASE BB-MA-CCTB E (Conical cemented titan base non hex and screw)	BSAD9433	30			10-12
SET ZZ BASE BB-MA-CCTB EW/BW (Conical cemented titan base non hex and screw)	BSAD9434	30			10-12
SET ZZ BASE BB-MA-CCTB-K85 B (Conical cemented titan base non hex K85 and screw)	BSAD9431K85	30			10-12
SET ZZ BASE BB-MA-CCTB-K85 C (Conical cemented titan base non hex K85 and screw)	BSAD9432K85	30			10-12
SET ZZ BASE BB-MA-CCTB-K85 E (Conical cemented titan base non hex K85 and screw)	BSAD9433K85	30			10-12
SET ZZ BASE BB-MA-ICOTH B (Impression coping open tray hex and Screw and Driver)	BSAQ9441	15			10-12
SET ZZ BASE BB-MA-ICOTH C (Impression coping open tray hex and Screw and Driver)	BSAQ9442	15			10-12
SET ZZ BASE BB-MA-ICOTH E (Impression coping open tray hex and Screw and Driver)	BSAQ9443	15			10-12
SET ZZ BASE BB-MA-ICOTH EW/BW (Impression coping open tray hex and Screw and Driver)	BSAQ9444	15			10-12
SET ZZ BASE BB-MA-NTB C (Narrow titan base non hex and screw)	BSAP9432	30			10-12
SET ZZ BASE BB-MA-NTB E (Narrow titan base non hex and screw)	BSAP9433	30			10-12
SET ZZ BASE BB-MA-NTBH C (Narrow titan base hex and screw)	BSAP9482	30			10-12
SET ZZ BASE BB-MA-NTBH E (Narrow titan base hex and screw)	BSAP9483	30			10-12
SET ZZ BASE BB-MA-PCTBH B (Parallel cemented titan base hex and screw)	BSAB9481	30			10-12
SET ZZ BASE BB-MA-PCTBH C (Parallel cemented titan base hex and screw)	BSAB9482	30			10-12
SET ZZ BASE BB-MA-PCTBH E (Parallel cemented titan base hex and screw)	BSAB9483	30			10-12
SET ZZ BASE BB-MA-PCTBH EW/BW (Parallel cemented titan base hex and screw)	BSAB9484	30			10-12
SET ZZ BASE BB-MA-PCTBH-K85 B (Parallel cemented titan base hex K85 and screw)	BSAB9481K85	30			10-12
SET ZZ BASE BB-MA-PCTBH-K85 C (Parallel cemented titan base hex K85 and screw)	BSAB9482K85	30			10-12
SET ZZ BASE BB-MA-PCTBH-K85 E (Parallel cemented titan base hex K85 and screw)	BSAB9483K85	30			10-12
SET ZZ BASE BB-MA-SA B (Scananalog and screw)	ZBDE9431	15			10-12
SET ZZ BASE BB-MA-SA C (Scananalog and screw)	ZBDE9432	15			10-12
SET ZZ BASE BB-MA-SA E (Scananalog and screw)	ZBDE9433	15			10-12

BRAINBASE CORPORATION MYTIS ARROW IMPLANT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE BB-MA-SA EW/BW (Scanalog and screw)	ZBDE9434	15			10-12
SET ZZ BASE BB-MA-SM B (Scanmarker and screw)	ZBAD9431	30			10-12
SET ZZ BASE BB-MA-SM C (Scanmarker and screw)	ZBAD9432	30			10-12
SET ZZ BASE BB-MA-SM E (Scanmarker and screw)	ZBAD9433	30			10-12
SET ZZ BASE BB-MA-SM EW/BW (Scanmarker and screw)	ZBAD9434	30			10-12
SET ZZ BASE BB-MA-WSM B (White Scanmarker and screw)	ZBCD9431	30			10-12
SET ZZ BASE BB-MA-WSM C (White Scanmarker and screw)	ZBCD9432	30			10-12
SET ZZ BASE BB-MA-WSM E (White Scanmarker and screw)	ZBCD9433	30			10-12
SET ZZ BASE BB-MA-WSM EW/BW (White Scanmarker and screw)	ZBCD9434	30			10-12
ZZ BASE BB-MA-ASCSC B (Angled screw-channel Screw)	BSKE9401	20			10-12
ZZ BASE BB-MA-ASCSC C (Angled screw-channel Screw)	BSKE9402	20			10-12
ZZ BASE BB-MA-ASCSC E (Angled screw-channel Screw)	BSKE9403	20			10-12
ZZ BASE BB-MA-ICSC B (Impression coping Screw)	BSFE9401	15			10-12
ZZ BASE BB-MA-ICSC C (Impression coping Screw)	BSFE9402	15			10-12
ZZ BASE BB-MA-ICSC E (Impression coping Screw)	BSFE9403	15			10-12
ZZ BASE BB-MA-ICSC EW/BW (Impression coping Screw)	BSFE9404	15			10-12
ZZ BASE BB-MA-MI B (Model implant)	MIAG9401				10-12
ZZ BASE BB-MA-MI C (Model implant)	MIAG9402				10-12
ZZ BASE BB-MA-MI E (Model implant)	MIAG9403				10-12
ZZ BASE BB-MA-MI EW/BW (Model implant)	MIAG9404				10-12
ZZ BASE BB-MA-SC B (Screw)	BSAE9401	30			10-12
ZZ BASE BB-MA-SC C (Screw)	BSAE9402	30			10-12
ZZ BASE BB-MA-SC E (Screw)	BSAE9403	30			10-12
ZZ BASE BB-MA-SC EW/BW (Screw)	BSAE9404	30			10-12
ZZ BASE BB-MA-WMSM-L12 B (White Metal Scanmarker L12 with Screw)	ZBED9431L12	15			10-12
ZZ BASE BB-MA-WMSM-L12 C (White Metal Scanmarker L12 with Screw)	ZBED9432L12	15			10-12
ZZ BASE BB-MA-WMSM-L12 E (White Metal Scanmarker L12 with Screw)	ZBED9433L12	15			10-12
ZZ BASE BB-MA-WMSM-L12 EW/BW (White Metal Scanmarker L12 with Screw)	ZBED9434L12	15			10-12
ZZ BASE BB-MA-ZSC B (Zirconia Screw)	BSAE9431	30			10-12

BRAINBASE CORPORATION MYTIS ARROW IMPLANT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
ZZ BASE BB-MA-ZSC C (Zirconia Screw)	BSAE9432	30			10-12
ZZ BASE BB-MA-ZSC E (Zirconia Screw)	BSAE9433	30			10-12
ZZ BASE BB-MA-ZSC EW/BW (Zirconia Screw)	BSAE9434	30			10-12

BREDENT SKY® CLASSIC / BLUESKY

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE SY-CB-ABCH 3.5 (Raw-Abutment® customizable hex and screw)	BSAF2292	25			10-12
SET ZZ BASE SY-CB-ABCH 4.0/4.5/5.5 (Raw-Abutment® customizable hex and screw)	BSAF2291	25			10-12
SET ZZ BASE SY-CB-ASCTBH-K80 3.5 (Angled screw-channel titan base hex K80 and screw)	BSAT2282K80	20			10-12
SET ZZ BASE SY-CB-ASCTBH-K80 4.0/4.5/5.5 (Angled screw-channel titan base hex K80 and screw)	BSAT2281K80	20			10-12
SET ZZ BASE SY-CB-ASCTB-K80 3.5 (Angled screw-channel titan base non hex K80 and screw)	BSAT2232K80	20			10-12
SET ZZ BASE SY-CB-ASCTB-K80 4.0/4.5/5.5 (Angled screw-channel titan base non hex K80 and screw)	BSAT2231K80	20			10-12
SET ZZ BASE SY-CB-CCTB 3.5 (Conical cemented titan base non hex and screw)	BSAD2232	25			10-12
SET ZZ BASE SY-CB-CCTB 4.0/4.5/5.5 (Conical cemented titan base non hex and screw)	BSAD2231	25			10-12
SET ZZ BASE SY-CB-CCTB-K85 3.5 (Conical cemented titan base non hex K85 and screw)	BSAD2232K85	25			10-12
SET ZZ BASE SY-CB-CCTB-K85 4.0/4.5/5.5 (Conical cemented titan base non hex K85 and screw)	BSAD2231K85	25			10-12
SET ZZ BASE SY-CB-CCTB-L10 3.5 (Conical cemented titan base non hex L10 and screw)	BSAD2232L10	25			10-12
SET ZZ BASE SY-CB-CCTB-L10 4.0/4.5/5.5 (Conical cemented titan base non hex L10 and screw)	BSAD2231L10	25			10-12
SET ZZ BASE SY-CB-CCTB-L20 3.5 (Conical cemented titan base non hex L20 and screw)	BSAD2232L20	25			10-12
SET ZZ BASE SY-CB-CCTB-L20 4.0/4.5/5.5 (Conical cemented titan base non hex L20 and screw)	BSAD2231L20	25			10-12
SET ZZ BASE SY-CB-CCTB-L30 3.5 (Conical cemented titan base non hex L30 and screw)	BSAD2232L30	25			10-12
SET ZZ BASE SY-CB-CCTB-L30 4.0/4.5/5.5 (Conical cemented titan base non hex L30 and screw)	BSAD2231L30	25			10-12
SET ZZ BASE SY-CB-CCTB-L40 3.5 (Conical cemented titan base non hex L40 and screw)	BSAD2232L40	25			10-12
SET ZZ BASE SY-CB-CCTB-L40 4.0/4.5/5.5 (Conical cemented titan base non hex L40 and screw)	BSAD2231L40	25			10-12
SET ZZ BASE SY-CB-ICOTH 3.5 (Impression coping open tray hex and Screw and Driver)	BSAQ2242	15			10-12
SET ZZ BASE SY-CB-ICOTH 4.0/4.5/5.5 (Impression coping open tray hex and Screw and Driver)	BSAQ2241	15			10-12
SET ZZ BASE SY-CB-NTB 3.5 (Narrow titan base non hex and screw)	BSAP2232	25			10-12
SET ZZ BASE SY-CB-NTBH 3.5 (Narrow titan base hex and screw)	BSAP2282	25			10-12
SET ZZ BASE SY-CB-PCTBH 3.5 (Parallel cemented titan base hex and screw)	BSAB2282	25			10-12
SET ZZ BASE SY-CB-PCTBH 4.0/4.5/5.5 (Parallel cemented titan base hex and screw)	BSAB2281	25			10-12
SET ZZ BASE SY-CB-PCTBH-K85 3.5 (Parallel cemented titan base hex K85 and screw)	BSAB2282K85	25			10-12
SET ZZ BASE SY-CB-PCTBH-K85 4.0/4.5/5.5 (Parallel cemented titan base hex K85 and screw)	BSAB2281K85	25			10-12
SET ZZ BASE SY-CB-PCTBH-L10 3.5 (Parallel cemented titan base hex L10 and screw)	BSAB2282L10	25			10-12
SET ZZ BASE SY-CB-PCTBH-L10 4.0/4.5/5.5 (Parallel cemented titan base hex L10 and screw)	BSAB2281L10	25			10-12
SET ZZ BASE SY-CB-PCTBH-L20 3.5 (Parallel cemented titan base hex L20 and screw)	BSAB2282L20	25			10-12

BREDED SKY® CLASSIC / BLUESKY

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE SY-CB-PCTBH-L20 4.0/4.5/5.5 (Parallel cemented titan base hex L20 and screw)	BSAB2281L20	25			10-12
SET ZZ BASE SY-CB-PCTBH-L30 3.5 (Parallel cemented titan base hex L30 and screw)	BSAB2282L30	25			10-12
SET ZZ BASE SY-CB-PCTBH-L30 4.0/4.5/5.5 (Parallel cemented titan base hex L30 and screw)	BSAB2281L30	25			10-12
SET ZZ BASE SY-CB-PCTBH-L40 3.5 (Parallel cemented titan base hex L40 and screw)	BSAB2282L40	25			10-12
SET ZZ BASE SY-CB-PCTBH-L40 4.0/4.5/5.5 (Parallel cemented titan base hex L40 and screw)	BSAB2281L40	25			10-12
SET ZZ BASE SY-CB-SA 3.5 (Scananalog and screw)	ZBDE2232	15			10-12
SET ZZ BASE SY-CB-SA 4.0/4.5/5.5 (Scananalog and screw)	ZBDE2231	15			10-12
SET ZZ BASE SY-CB-SM 3.5 (Scanmarker and screw)	ZBAD2232	25			10-12
SET ZZ BASE SY-CB-SM 4.0/4.5/5.5 (Scanmarker and screw)	ZBAD2231	25			10-12
SET ZZ BASE SY-CB-WSM 3.5 (White Scanmarker and screw)	ZBCD2232	25			10-12
SET ZZ BASE SY-CB-WSM 4.0/4.5/5.5 (White Scanmarker and screw)	ZBCD2231	25			10-12
SET ZZ BASE SY-CB-ZZMUA-L15 3.5 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH2232L15	25	25		10-12
SET ZZ BASE SY-CB-ZZMUA-L15 4.0/4.5/5.5 (Zirkonzahn Multi Unit Abutment non hex L15 and S	BSAH2231L15	25	25		10-12
SET ZZ BASE SY-CB-ZZMUA-L25 3.5 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH2232L25	25	25		10-12
SET ZZ BASE SY-CB-ZZMUA-L25 4.0/4.5/5.5 (Zirkonzahn Multi Unit Abutment non hex L25 and S	BSAH2231L25	25	25		10-12
SET ZZ BASE SY-CB-ZZMUA-L35 3.5 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH2232L35	25	25		10-12
SET ZZ BASE SY-CB-ZZMUA-L35 4.0/4.5/5.5 (Zirkonzahn Multi Unit Abutment non hex L35 and S	BSAH2231L35	25	25		10-12
SET ZZ BASE SY-CB-ZZMUA-L45 3.5 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH2232L45	25	25		10-12
SET ZZ BASE SY-CB-ZZMUA-L45 4.0/4.5/5.5 (Zirkonzahn Multi Unit Abutment non hex L45 and S	BSAH2231L45	25	25		10-12
ZZ BASE SY-CB-ASCSC 3.5 (Angled screw-channel Screw)	BSKE2202	20			10-12
ZZ BASE SY-CB-ASCSC 4.0/4.5/5.5 (Angled screw-channel Screw)	BSKE2201	20			10-12
ZZ BASE SY-CB-ICSC 3.5 (Impression coping Screw)	BSFE2202	15			10-12
ZZ BASE SY-CB-ICSC 4.0/4.5/5.5 (Impression coping Screw)	BSFE2201	15			10-12
ZZ BASE SY-CB-MI 3.5 (Model implant)	MIAG2202				10-12
ZZ BASE SY-CB-MI 4.0/4.5/5.5 (Model implant)	MIAG2201				10-12
ZZ BASE SY-CB-SC 3.5 (Screw)	BSAE2202	25			10-12
ZZ BASE SY-CB-SC 4.0/4.5/5.5 (Screw)	BSAE2201	25			10-12
ZZ BASE SY-CB-WMSM-L12 3.5 (White Metal Scanmarker L12 with Screw)	ZBED2232L12	15			10-12
ZZ BASE SY-CB-WMSM-L12 4.0/4.5/5.5 (White Metal Scanmarker L12 with Screw)	ZBED2231L12	15			10-12

BREDENT SKY® CLASSIC / BLUESKY

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
ZZ BASE SY-CB-ZSC 3.5 (Zirconia screw)	BSAE2232	25			10-12
ZZ BASE SY-CB-ZSC 4.0/4.5/5.5 (Zirconia screw)	BSAE2231	25			10-12
ZZ BASE SY-CB-ZSC-G 4.0/4.5/5.5 (Zirconia Screw Gold)	BSCE2231	25			10-12

BREDENT SKY® FAST & FIXED

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE SY-FF-CCTB (Conical cemented titan base non hex and screw)	BSAD2331	20			10-12
SET ZZ BASE SY-FF-ICOT (Impression coping open tray non hex and Screw and Driver)	BSAQ2361	15			10-12
SET ZZ BASE SY-FF-SA (Scanalog and screw)	ZBDE2331	15			10-12
SET ZZ BASE SY-FF-SM (Scanmarker and screw)	ZBAD2331	20			10-12
ZZ BASE SY-FF-ICSC (Impression coping Screw)	BSFE2301	15			10-12
ZZ BASE SY-FF-MI (Model implant)	MIAG2301				10-12
ZZ BASE SY-FF-SC (Screw)	BSAE2301	20			10-12
ZZ BASE SY-FF-WMSM-L12 (White Metal Scanmarker L12 with Screw)	ZBED2331L12	15			10-12
ZZ BASE SY-FF-WMSM-L15 (White Metal Scanmarker L15 with Screw)	ZBED2331L15	15			10-12
ZZ BASE SY-FF-ZSC (Zirconia screw)	BSAE2331	20			10-12

BREDENT SKY® UNI.CONE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE SY-UC-CCTB (Conical cemented titan base non hex and screw)	BSAD8931	18			10-12
SET ZZ BASE SY-UC-CCTB-K85 (Conical cemented titan base non hex K85 and screw)	BSAD8931K85	18			10-12
SET ZZ BASE SY-UC-SA (Scananalog and screw)	ZBDE8931	15			10-12
SET ZZ BASE SY-UC-SM (Scanmarker and screw)	ZBAD8931	18			10-12
SET ZZ BASE SY-UC-WSM (White Scanmarker and screw)	ZBCD8931	18			10-12
ZZ BASE SY-UC-MI (Model implant)	MIAG8901				10-12
ZZ BASE SY-UC-SC (Screw)	BSAE8901	18			10-12
ZZ BASE SY-UC-ZSC (Zirconia screw)	BSAE8931	18			10-12

BTI® CONICAL SPACER

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE BT-MCS-CCTB 4.1/5.5 (Conical cemented titan base non hex and screw)	BSAD3333	10			10-12
SET ZZ BASE BT-MCS-SA 4.1/5.5 (Scanalog and screw)	ZBDE3333	15			10-12
SET ZZ BASE BT-MCS-SM 4.1/5.5 (Scanmarker and screw)	ZBAD3333	10			10-12
SET ZZ BASE BT-MCS-WSM 4.1/5.5 (White Scanmarker and screw)	ZBCD3333	10			10-12
ZZ BASE BT-MCS-SC 4.1/5.5 (Screw)	BSAE3303	10			10-12
ZZ BASE BT-MCS-ZSC 4.1/5.5 (Zirconia screw)	BSAE3333	10			10-12

BTI® EXTERNA®					
Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE BT-E-ABCH 4.1 (Raw-Abutment® customizable hex and screw)	BSAF3291	20			10-12
SET ZZ BASE BT-E-ABCH 5.5 (Raw-Abutment® customizable hex and screw)	BSAF3292	20			10-12
SET ZZ BASE BT-E-CCTB 4.1 (Conical cemented titan base non hex and screw)	BSAD3231	20			10-12
SET ZZ BASE BT-E-CCTB 5.5 (Conical cemented titan base non hex and screw)	BSAD3232	20			10-12
SET ZZ BASE BT-E-CCTB-K85 4.1 (Conical cemented titan base non hex K85 and screw)	BSAD3231K85	20			10-12
SET ZZ BASE BT-E-CCTB-K85 5.5 (Conical cemented titan base non hex K85 and screw)	BSAD3232K85	20			10-12
SET ZZ BASE BT-E-ICOTH 4.1 (Impression coping open tray hex and Screw and Driver)	BSAQ3241	15			10-12
SET ZZ BASE BT-E-ICOTH 5.5 (Impression coping open tray hex and Screw and Driver)	BSAQ3242	15			10-12
SET ZZ BASE BT-E-PCTBH 4.1 (Parallel cemented titan base hex and screw)	BSAB3281	20			10-12
SET ZZ BASE BT-E-PCTBH 5.5 (Parallel cemented titan base hex and screw)	BSAB3282	20			10-12
SET ZZ BASE BT-E-PCTBH-K85 4.1 (Parallel cemented titan base hex K85 and screw)	BSAB3281K85	20			10-12
SET ZZ BASE BT-E-PCTBH-K85 5.5 (Parallel cemented titan base hex K85 and screw)	BSAB3282K85	20			10-12
SET ZZ BASE BT-E-SA 4.1 (Scanalog and screw)	ZBDE3231	15			10-12
SET ZZ BASE BT-E-SA 5.5 (Scanalog and screw)	ZBDE3232	15			10-12
SET ZZ BASE BT-E-SM 4.1 (Scanmarker and screw)	ZBAD3231	20			10-12
SET ZZ BASE BT-E-SM 5.5 (Scanmarker and screw)	ZBAD3232	20			10-12
SET ZZ BASE BT-E-WSM 4.1 (White Scanmarker and screw)	ZBCD3231	20			10-12
SET ZZ BASE BT-E-WSM 5.5 (White Scanmarker and screw)	ZBCD3232	20			10-12
ZZ BASE BT-E-ASCSC 4.1 (Angled screw-channel Screw)	BSKE3201	20			10-12
ZZ BASE BT-E-ICSC 4.1 (Impression coping Screw)	BSFE3201	15			10-12
ZZ BASE BT-E-ICSC 5.5 (Impression coping Screw)	BSFE3202	15			10-12
ZZ BASE BT-E-SC 4.1 (Screw)	BSAE3201	20			10-12
ZZ BASE BT-E-SC 5.5 (Screw)	BSAE3202	20			10-12
ZZ BASE BT-E-SC-G 4.1 (Screw Gold)	BSCE3201	20			10-12
ZZ BASE BT-E-SC-G 5.5 (Screw Gold)	BSCE3202	20			10-12
ZZ BASE BT-E-WMSM-L12 4.1 (White Metal Scanmarker L12 with Screw)	ZBED3231L12	15			10-12
ZZ BASE BT-E-WMSM-L12 5.5 (White Metal Scanmarker L12 with Screw)	ZBED3232L12	15			10-12
ZZ BASE BT-E-ZSC 4.1 (Zirconia screw)	BSAE3231	20			10-12
ZZ BASE BT-E-ZSC 5.5 (Zirconia screw)	BSAE3232	20			10-12

BTI® EXTERNA®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
ZZ BASE BT-E-ZSC-G 4.1 (Zirconia Screw Gold)	BSCE3231	20			10-12
ZZ BASE BT-E-ZSC-G 5.5 (Zirconia Screw Gold)	BSCE3232	20			10-12

BTI® INTERNA®					
Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE BT-I-ABCH 3.0 (Raw-Abutment® customizable hex and screw)	BSAF3693	35			10-12
SET ZZ BASE BT-I-ABCH 3.5 (Raw-Abutment® customizable hex and screw)	BSAF3694	35			10-12
SET ZZ BASE BT-I-ABCH 4.1 (Raw-Abutment® customizable hex and screw)	BSAF3691	20			10-12
SET ZZ BASE BT-I-ABCH 5.5 (Raw-Abutment® customizable hex and screw)	BSAF3692	20			10-12
SET ZZ BASE BT-I-ASCTBH-K80 3.0 (Angled screw-channel titan base hex K80 and screw)	BSAT3683K80	20			10-12
SET ZZ BASE BT-I-ASCTBH-K80 3.5 (Angled screw-channel titan base hex K80 and screw)	BSAT3684K80	20			10-12
SET ZZ BASE BT-I-ASCTBH-K80 4.1 (Angled screw-channel titan base hex K80 and screw)	BSAT3681K80	20			10-12
SET ZZ BASE BT-I-ASCTB-K80 3.5 (Angled screw-channel titan base non hex K80 and screw)	BSAT3634K80	20			10-12
SET ZZ BASE BT-I-ASCTB-K80 4.1 (Angled screw-channel titan base non hex K80 and screw)	BSAT3631K80	20			10-12
SET ZZ BASE BT-I-CCTB 3.5 (Conical cemented titan base non hex and screw)	BSAD3634	35			10-12
SET ZZ BASE BT-I-CCTB 4.1 (Conical cemented titan base non hex and screw)	BSAD3631	20			10-12
SET ZZ BASE BT-I-CCTB 5.5 (Conical cemented titan base non hex and screw)	BSAD3632	20			10-12
SET ZZ BASE BT-I-CCTB-K85 3.5 (Conical cemented titan base non hex K85 and screw)	BSAD3634K85	35			10-12
SET ZZ BASE BT-I-CCTB-K85 4.1 (Conical cemented titan base non hex K85 and screw)	BSAD3631K85	20			10-12
SET ZZ BASE BT-I-CCTB-K85 5.5 (Conical cemented titan base non hex K85 and screw)	BSAD3632K85	20			10-12
SET ZZ BASE BT-I-ICOTH 3.0 (Impression coping open tray hex and Screw and Driver)	BSAQ3643	15			10-12
SET ZZ BASE BT-I-ICOTH 3.5 (Impression coping open tray hex and Screw and Driver)	BSAQ3644	15			10-12
SET ZZ BASE BT-I-ICOTH 4.1 (Impression coping open tray hex and Screw and Driver)	BSAQ3641	15			10-12
SET ZZ BASE BT-I-ICOTH 5.5 (Impression coping open tray hex and Screw and Driver)	BSAQ3642	15			10-12
SET ZZ BASE BT-I-NTB 3.5 (Narrow titan base non hex and screw)	BSAP3634	35			10-12
SET ZZ BASE BT-I-NTB 4.1 (Narrow titan base non hex and screw)	BSAP3631	20			10-12
SET ZZ BASE BT-I-NTBH 3.0 (Narrow titan base hex and screw)	BSAP3683	35			10-12
SET ZZ BASE BT-I-NTBH 3.5 (Narrow titan base hex and screw)	BSAP3684	35			10-12
SET ZZ BASE BT-I-NTBH 4.1 (Narrow titan base hex and screw)	BSAP3681	20			10-12
SET ZZ BASE BT-I-PCTBH 3.0 (Parallel cemented titan base hex and screw)	BSAB3683	35			10-12
SET ZZ BASE BT-I-PCTBH 3.5 (Parallel cemented titan base hex and screw)	BSAB3684	35			10-12
SET ZZ BASE BT-I-PCTBH 4.1 (Parallel cemented titan base hex and screw)	BSAB3681	20			10-12
SET ZZ BASE BT-I-PCTBH 5.5 (Parallel cemented titan base hex and screw)	BSAB3682	20			10-12
SET ZZ BASE BT-I-PCTBH-K85 3.0 (Parallel cemented titan base hex K85 and screw)	BSAB3683K85	35			10-12

BTI® INTERNA®					
Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE BT-I-PCTBH-K85 3.5 (Parallel cemented titan base hex K85 and screw)	BSAB3684K85	35			10-12
SET ZZ BASE BT-I-PCTBH-K85 4.1 (Parallel cemented titan base hex K85 and screw)	BSAB3681K85	20			10-12
SET ZZ BASE BT-I-PCTBH-K85 5.5 (Parallel cemented titan base hex K85 and screw)	BSAB3682K85	20			10-12
SET ZZ BASE BT-I-SA 3.5 (Scananalog and screw)	ZBDE3634	15			10-12
SET ZZ BASE BT-I-SA 4.1 (Scananalog and screw)	ZBDE3631	15			10-12
SET ZZ BASE BT-I-SA 5.5 (Scananalog and screw)	ZBDE3632	15			10-12
SET ZZ BASE BT-I-SM 3.0 (Scanmarker and screw)	ZBAD3633	35			10-12
SET ZZ BASE BT-I-SM 3.5 (Scanmarker and screw)	ZBAD3634	35			10-12
SET ZZ BASE BT-I-SM 4.1 (Scanmarker and screw)	ZBAD3631	20			10-12
SET ZZ BASE BT-I-SM 5.5 (Scanmarker and screw)	ZBAD3632	20			10-12
SET ZZ BASE BT-I-WSM 4.1 (White Scanmarker and screw)	ZBCD3631	20			10-12
SET ZZ BASE BT-I-WSM 5.5 (White Scanmarker and screw)	ZBCD3632	20			10-12
SET ZZ BASE BT-I-ZZMUA-L15 3.0 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH3633L15	25	35		10-12
SET ZZ BASE BT-I-ZZMUA-L25 3.0 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH3633L25	25	35		10-12
SET ZZ BASE BT-I-ZZMUA-L35 3.0 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH3633L35	25	35		10-12
SET ZZ BASE BT-I-ZZMUA-L45 3.0 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH3633L45	25	35		10-12
ZZ BASE BT-I-ASCSC 3.0 (Angled screw-channel Screw)	BSKE3603	20			10-12
ZZ BASE BT-I-ASCSC 3.5 (Angled screw-channel Screw)	BSKE3604	20			10-12
ZZ BASE BT-I-ASCSC 4.1 (Angled screw-channel Screw)	BSKE3601	20			10-12
ZZ BASE BT-I-ICSC 3.0 (Impression coping Screw)	BSFE3603	15			10-12
ZZ BASE BT-I-ICSC 3.5 (Impression coping Screw)	BSFE3604	15			10-12
ZZ BASE BT-I-ICSC 4.1 (Impression coping Screw)	BSFE3601	15			10-12
ZZ BASE BT-I-ICSC 5.5 (Impression coping Screw)	BSFE3602	15			10-12
ZZ BASE BT-I-MI 3.0 (Model implant)	MIAG3603				10-12
ZZ BASE BT-I-MI 3.5 (Model implant)	MIAG3604				10-12
ZZ BASE BT-I-MI 4.1 (Model implant)	MIAG3601				10-12
ZZ BASE BT-I-MI 5.5 (Model implant)	MIAG3602				10-12
ZZ BASE BT-I-SC 3.0 (Screw)	BSAE3603	35			10-12
ZZ BASE BT-I-SC 3.5 (Screw)	BSAE3604	35			10-12

BTI® INTERNA®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE BT-I-SC 4.1 (Screw)	BSAE3601	20			10-12
ZZ BASE BT-I-SC 5.5 (Screw)	BSAE3602	20			10-12
ZZ BASE BT-I-WMSM-L12 3.0 (White Metal Scanmarker L12 with Screw)	ZBED3633L12	15			10-12
ZZ BASE BT-I-WMSM-L12 3.5 (White Metal Scanmarker L12 with Screw)	ZBED3634L12	15			10-12
ZZ BASE BT-I-WMSM-L12 4.1 (White Metal Scanmarker L12 with Screw)	ZBED3631L12	15			10-12
ZZ BASE BT-I-WMSM-L12 5.5 (White Metal Scanmarker L12 with Screw)	ZBED3632L12	15			10-12
ZZ BASE BT-I-ZSC 3.0 (Zirconia screw)	BSAE3633	35			10-12
ZZ BASE BT-I-ZSC 3.5 (Zirconia screw)	BSAE3634	35			10-12
ZZ BASE BT-I-ZSC 4.1 (Zirconia screw)	BSAE3631	20			10-12
ZZ BASE BT-I-ZSC 5.5 (Zirconia screw)	BSAE3632	20			10-12

BTI® MULTI-IM®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE BT-MIM-CCTB 3.5/4.1 (Conical cemented titan base non hex and screw)	BSAD3331	10			10-12
SET ZZ BASE BT-MIM-CCTB 5.5 (Conical cemented titan base non hex and screw)	BSAD3332	10			10-12
SET ZZ BASE BT-MIM-SA 3.5/4.1 (Scanalog and screw)	ZBDE3331	15			10-12
SET ZZ BASE BT-MIM-SA 5.5 (Scanalog and screw)	ZBDE3332	15			10-12
SET ZZ BASE BT-MIM-SM 3.5/4.1 (Scanmarker and screw)	ZBAD3331	10			10-12
SET ZZ BASE BT-MIM-SM 5.5 (Scanmarker and screw)	ZBAD3332	10			10-12
SET ZZ BASE BT-MIM-WSM 3.5/4.1 (White Scanmarker and screw)	ZBCD3331	10			10-12
SET ZZ BASE BT-MIM-WSM 5.5 (White Scanmarker and screw)	ZBCD3332	10			10-12
ZZ BASE BT-MIM-SC 3.5/4.1 (Screw)	BSAE3301	10			10-12
ZZ BASE BT-MIM-SC 5.5 (Screw)	BSAE3302	10			10-12
ZZ BASE BT-MIM-ZSC 3.5/4.1 (Zirconia screw)	BSAE3331	10			10-12
ZZ BASE BT-MIM-ZSC 5.5 (Zirconia screw)	BSAE3332	10			10-12

BTI® MULTI-IM® ANGLED®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE BT-MIMA-CCTB 3.5/4.1 (Conical cemented titan base non hex and screw)	BSAD3334	10			10-12
SET ZZ BASE BT-MIMA-SA 3.5/4.1 (Scananalog and screw)	ZBDE3334	15			10-12
SET ZZ BASE BT-MIMA-SM 3.5/4.1 (Scanmarker and screw)	ZBAD3334	10			10-12
SET ZZ BASE BT-MIMA-WSM 3.5/4.1 (White Scanmarker and screw)	ZBCD3334	10			10-12
ZZ BASE BT-MIMA-SC 3.5/4.1 (Screw)	BSAE3304	10			10-12
ZZ BASE BT-MIMA-ZSC 3.5/4.1 (Zirconia screw)	BSAE3334	10			10-12

BTI® TINY®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
SET ZZ BASE BT-T-ABCH 3.5 (Raw-Abutment® customizable hex and screw)	BSAF3191	10			10-12
SET ZZ BASE BT-T-CCTB 3.5 (Conical cemented titan base non hex and screw)	BSAD3131	10			10-12
SET ZZ BASE BT-T-ICOTH 3.5 (Impression coping open tray hex and Screw and Driver)	BSAQ3141	15			10-12
SET ZZ BASE BT-T-NTB 3.5 (Narrow titan base non hex and screw)	BSAP3131	10			10-12
SET ZZ BASE BT-T-NTBH 3.5 (Narrow titan base hex and screw)	BSAP3181	10			10-12
SET ZZ BASE BT-T-PCTBH 3.5 (Parallel cemented titan base hex and screw)	BSAB3181	10			10-12
SET ZZ BASE BT-T-SA 3.5 (Scananalog and screw)	ZBDE3131	15			10-12
SET ZZ BASE BT-T-SM 3.5 (Scanmarker and screw)	ZBAD3131	10			10-12
SET ZZ BASE BT-T-WSM 3.5 (White Scanmarker and screw)	ZBCD3131	10			10-12
ZZ BASE BT-T-ICSC 3.5 (Impression coping Screw)	BSFE3101	15			10-12
ZZ BASE BT-T-SC 3.5 (Screw)	BSAE3101	10			10-12
ZZ BASE BT-T-WMSM-L12 3.5 (White Metal Scanmarker L12 with Screw)	ZBED3131L12	15			10-12
ZZ BASE BT-T-ZSC 3.5 (Zirconia screw)	BSAE3131	10			10-12

BTK THE SMILE SYSTEM® - BT-ISKONE/ BT-SAFE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE BK-IS-ABCH IR (Raw-Abutment® customizable hex and screw)	BSAF13381	20			10-12
SET ZZ BASE BK-IS-ABCH-D14 IR (Raw-Abutment® customizable hex D14 and screw)	BSAF13381D14	20			10-12
SET ZZ BASE BK-IS-CCTB IR (Conical cemented titan base non hex and screw)	BSAD13331	20			10-12
SET ZZ BASE BK-IS-CCTB-K85 IR (Conical cemented titan base non hex K85 and screw)	BSAD13331K85	20			10-12
SET ZZ BASE BK-IS-CCTB-L10 IR (Conical cemented titan base non hex L10 and screw)	BSAD13331L10	20			10-12
SET ZZ BASE BK-IS-CCTB-L20 IR (Conical cemented titan base non hex L20 and screw)	BSAD13331L20	20			10-12
SET ZZ BASE BK-IS-CCTB-L30 IR (Conical cemented titan base non hex L30 and screw)	BSAD13331L30	20			10-12
SET ZZ BASE BK-IS-CCTB-L40 IR (Conical cemented titan base non hex L40 and screw)	BSAD13331L40	20			10-12
SET ZZ BASE BK-IS-NTB IR (Narrow titan base non hex and screw)	BSAP13331	20			10-12
SET ZZ BASE BK-IS-NTBH IR (Narrow titan base hex and screw)	BSAP13381	20			10-12
SET ZZ BASE BK-IS-PCTBH IR (Parallel cemented titan base hex and screw)	BSAB13391	20			10-12
SET ZZ BASE BK-IS-PCTBH-K85 IR (Parallel cemented titan base hex K85 and screw)	BSAB13381K85	20			10-12
SET ZZ BASE BK-IS-PCTBH-L20 IR (Parallel cemented titan base hex L20 and screw)	BSAB13391L20	20			10-12
SET ZZ BASE BK-IS-PCTBH-L30 IR (Parallel cemented titan base hex L30 and screw)	BSAB13391L30	20			10-12
SET ZZ BASE BK-IS-PCTBH-L40 IR (Parallel cemented titan base hex L40 and screw)	BSAB13391L40	20			10-12
SET ZZ BASE BK-IS-SA IR (Scanalog and screw)	ZBDE13331	15			10-12
SET ZZ BASE BK-IS-SM IR (Scanmarker and screw)	ZBAD13331	20			10-12
SET ZZ BASE BK-IS-WSM IR (White Scanmarker and screw)	ZBCD13331	20			10-12
SET ZZ BASE BK-IS-ZZMUA17 IR Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and implant)	BSAI13331	20			10-12
SET ZZ BASE BK-IS-ZZMUA-L15 IR (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH13331L15	25	20		10-12
SET ZZ BASE BK-IS-ZZMUA-L25 IR (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH13331L25	25	20		10-12
SET ZZ BASE BK-IS-ZZMUA-L35 IR (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH13331L35	25	20		10-12
SET ZZ BASE BK-IS-ZZMUA-L45 IR (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH13331L45	25	20		10-12
ZZ BASE BK-IS-MI IR (Model implant)	MIAG13301				10-12
ZZ BASE BK-IS-SC IR (Screw)	BSAE13301	20			10-12
ZZ BASE BK-IS-SC-G IR (Screw Gold)	BSCE13301	20			10-12
ZZ BASE BK-IS-WMSM-L12 IR (White Metal Scanmarker L12 with Screw)	ZBED13331L12	15			10-12
ZZ BASE BK-IS-ZSC IR (Zirconia screw)	BSAE13331	20			10-12
ZZ BASE BK-IS-ZZMUASC IR (Zirkonzahn Multi Unit Abutment Screw)	BSEE13301	20			10-12

BTK THE SMILE SYSTEM® - BT-KLASSIC

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE BK-KL-ABCH EN (Raw-Abutment® customizable hex and screw)	BSAF13191	25			10-12
SET ZZ BASE BK-KL-ABCH ER (Raw-Abutment® customizable hex and screw)	BSAF13192	35			10-12
SET ZZ BASE BK-KL-ABCH-D14 EN (Raw-Abutment® customizable hex D14 and screw)	BSAF13191D14	25			10-12
SET ZZ BASE BK-KL-ABCH-D14 ER (Raw-Abutment® customizable hex D14 and screw)	BSAF13192D14	35			10-12
SET ZZ BASE BK-KL-CCTB EN (Conical cemented titan base non hex and screw)	BSAD13131	25			10-12
SET ZZ BASE BK-KL-CCTB ER (Conical cemented titan base non hex and screw)	BSAD13132	35			10-12
SET ZZ BASE BK-KL-CCTB-K85 EN (Conical cemented titan base non hex K85 and screw)	BSAD13131K85	25			10-12
SET ZZ BASE BK-KL-CCTB-K85 ER (Conical cemented titan base non hex K85 and screw)	BSAD13132K85	35			10-12
SET ZZ BASE BK-KL-CCTB-L10 EN (Conical cemented titan base non hex L10 and screw)	BSAD13131L10	25			10-12
SET ZZ BASE BK-KL-CCTB-L10 ER (Conical cemented titan base non hex L10 and screw)	BSAD13132L10	35			10-12
SET ZZ BASE BK-KL-CCTB-L10-G ER (Conical cemented titan base non hex L10 and screw Gold)	BSCD13132L10	35			10-12
SET ZZ BASE BK-KL-CCTB-L20 EN (Conical cemented titan base non hex L20 and screw)	BSAD13131L20	25			10-12
SET ZZ BASE BK-KL-CCTB-L20 ER (Conical cemented titan base non hex L20 and screw)	BSAD13132L20	35			10-12
SET ZZ BASE BK-KL-CCTB-L30 EN (Conical cemented titan base non hex L30 and screw)	BSAD13131L30	25			10-12
SET ZZ BASE BK-KL-CCTB-L30 ER (Conical cemented titan base non hex L30 and screw)	BSAD13132L30	35			10-12
SET ZZ BASE BK-KL-ICOTH EN (Impression coping open tray hex and Screw and Driver)	BSAQ13141	15			10-12
SET ZZ BASE BK-KL-ICOTH ER (Impression coping open tray hex and Screw and Driver)	BSAQ13142	15			10-12
SET ZZ BASE BK-KL-NTB EN (Narrow titan base non hex and screw)	BSAP13131	25			10-12
SET ZZ BASE BK-KL-NTB ER (Narrow titan base non hex and screw)	BSAP13132	35			10-12
SET ZZ BASE BK-KL-NTBH EN (Narrow titan base hex and screw)	BSAP13181	25			10-12
SET ZZ BASE BK-KL-NTBH ER (Narrow titan base hex and screw)	BSAP13182	35			10-12
SET ZZ BASE BK-KL-PCTBH EN (Parallel cemented titan base hex and screw)	BSAB13181	25			10-12
SET ZZ BASE BK-KL-PCTBH ER (Parallel cemented titan base hex and screw)	BSAB13182	35			10-12
SET ZZ BASE BK-KL-PCTBH-K85 EN (Parallel cemented titan base hex K85 and screw)	BSAB13181K85	25			10-12
SET ZZ BASE BK-KL-PCTBH-K85 ER (Parallel cemented titan base hex K85 and screw)	BSAB13182K85	35			10-12
SET ZZ BASE BK-KL-PCTBH-L10 EN (Parallel cemented titan base hex L10 and screw)	BSAB13191L10	25			10-12
SET ZZ BASE BK-KL-PCTBH-L10 ER (Parallel cemented titan base hex L10 and screw)	BSAB13192L10	35			10-12
SET ZZ BASE BK-KL-PCTBH-L10-G ER (Parallel cemented titan base hex L10 and screw Gold)	BSCB13192L10	35			10-12
SET ZZ BASE BK-KL-PCTBH-L20 EN (Parallel cemented titan base hex L20 and screw)	BSAB13191L20	25			10-12

BTK THE SMILE SYSTEM® - BT-KLASSIC

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE BK-KL-PCTBH-L20 ER (Parallel cemented titan base hex L20 and screw)	BSAB13192L20	35			10-12
SET ZZ BASE BK-KL-PCTBH-L30 EN (Parallel cemented titan base hex L30 and screw)	BSAB13191L30	25			10-12
SET ZZ BASE BK-KL-PCTBH-L30 ER (Parallel cemented titan base hex L30 and screw)	BSAB13192L30	35			10-12
SET ZZ BASE BK-KL-PCTBH-L40 EN (Parallel cemented titan base hex L40 and screw)	BSAB13191L40	25			10-12
SET ZZ BASE BK-KL-PCTBH-L40 ER (Parallel cemented titan base hex L40 and screw)	BSAB13192L40	35			10-12
SET ZZ BASE BK-KL-SA EN (Scananalog and screw)	ZBDE13131	15			10-12
SET ZZ BASE BK-KL-SA ER (Scananalog and screw)	ZBDE13132	15			10-12
SET ZZ BASE BK-KL-SM EN (Scanmarker and screw)	ZBAD13131	25			10-12
SET ZZ BASE BK-KL-SM ER (Scanmarker and screw)	ZBAD13132	35			10-12
SET ZZ BASE BK-KL-WSM EN (White Scanmarker and screw)	ZBCD13131	25			10-12
SET ZZ BASE BK-KL-WSM ER (White Scanmarker and screw)	ZBCD13132	35			10-12
SET ZZ BASE BK-KL-ZZMUA17 EN Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and imple	BSAI13131	20			10-12
SET ZZ BASE BK-KL-ZZMUA17 ER Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and imple	BSAI13132	30			10-12
SET ZZ BASE BK-KL-ZZMUA17-L10 EN Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10 ;	BSAI13131L10	20			10-12
SET ZZ BASE BK-KL-ZZMUA17-L10 ER Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10 ;	BSAI13132L10	30			10-12
SET ZZ BASE BK-KL-ZZMUA-L15 EN (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH13131L15	25	25		10-12
SET ZZ BASE BK-KL-ZZMUA-L15 ER (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH13132L15	25	35		10-12
SET ZZ BASE BK-KL-ZZMUA-L25 EN (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH13131L25	25	25		10-12
SET ZZ BASE BK-KL-ZZMUA-L25 ER (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH13132L25	25	35		10-12
SET ZZ BASE BK-KL-ZZMUA-L35 EN (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH13131L35	25	25		10-12
SET ZZ BASE BK-KL-ZZMUA-L35 ER (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH13132L35	25	35		10-12
SET ZZ BASE BK-KL-ZZMUA-L45 EN (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH13131L45	25	25		10-12
SET ZZ BASE BK-KL-ZZMUA-L45 ER (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH13132L45	25	35		10-12
ZZ BASE BK-KL-ICSC EN (Impression coping Screw)	BSFE13101	15			10-12
ZZ BASE BK-KL-ICSC ER (Impression coping Screw)	BSFE13102	15			10-12
ZZ BASE BK-KL-MI EN (Model implant)	MIAG13101				10-12
ZZ BASE BK-KL-MI ER (Model implant)	MIAG13102				10-12
ZZ BASE BK-KL-SC EN (Screw)	BSAE13101	25			10-12
ZZ BASE BK-KL-SC ER (Screw)	BSAE13102	35			10-12

BTK THE SMILE SYSTEM® - BT-KLASSIC

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE BK-KL-SC-G ER (Screw Gold)	BSCE13102	35			10-12
ZZ BASE BK-KL-WMSM-L12 EN (White Metal Scanmarker L12 with Screw)	ZBED13131L12	15			10-12
ZZ BASE BK-KL-WMSM-L12 ER (White Metal Scanmarker L12 with Screw)	ZBED13132L12	15			10-12
ZZ BASE BK-KL-ZSC EN (Zirconia screw)	BSAE13131	25			10-12
ZZ BASE BK-KL-ZSC ER (Zirconia screw)	BSAE13132	35			10-12
ZZ BASE BK-KL-ZSC-G EN (Zirconia screw Gold)	BSCE13131	25			10-12
ZZ BASE BK-KL-ZSC-G ER (Zirconia screw Gold)	BSCE13132	35			10-12
ZZ BASE BK-KL-ZZMUASC EN (Zirkonzahn Multi Unit Abutment Screw)	BSEE13101	20			10-12
ZZ BASE BK-KL-ZZMUASC ER (Zirkonzahn Multi Unit Abutment Screw)	BSEE13102	30			10-12

BTK THE SMILE SYSTEM® - BT-KONIC

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE BK-KO-ABCH IR (Raw-Abutment® customizable hex and screw)	BSAF13281	30			10-12
SET ZZ BASE BK-KO-ABCH IW (Raw-Abutment® customizable hex and screw)	BSAF13282	30			10-12
SET ZZ BASE BK-KO-ABCH-D14 IR (Raw-Abutment® customizable hex D14 and screw)	BSAF13281D14	30			10-12
SET ZZ BASE BK-KO-ABCH-D14 IW (Raw-Abutment® customizable hex D14 and screw)	BSAF13282D14	30			10-12
SET ZZ BASE BK-KO-CCTB 5.7 (Conical cemented titan base non hex and screw)	BSAD13233	30			10-12
SET ZZ BASE BK-KO-CCTB IR (Conical cemented titan base non hex and screw)	BSAD13231	30			10-12
SET ZZ BASE BK-KO-CCTB IW (Conical cemented titan base non hex and screw)	BSAD13232	30			10-12
SET ZZ BASE BK-KO-CCTB-K85 IR (Conical cemented titan base non hex K85 and screw)	BSAD13231K85	30			10-12
SET ZZ BASE BK-KO-CCTB-K85 IW (Conical cemented titan base non hex K85 and screw)	BSAD13232K85	30			10-12
SET ZZ BASE BK-KO-CCTB-L10 5.7 (Conical cemented titan base non hex L10 and screw)	BSAD13233L10	30			10-12
SET ZZ BASE BK-KO-CCTB-L10 IR (Conical cemented titan base non hex L10 and screw)	BSAD13231L10	30			10-12
SET ZZ BASE BK-KO-CCTB-L10 IW (Conical cemented titan base non hex L10 and screw)	BSAD13232L10	30			10-12
SET ZZ BASE BK-KO-CCTB-L20 5.7 (Conical cemented titan base non hex L20 and screw)	BSAD13233L20	30			10-12
SET ZZ BASE BK-KO-CCTB-L20 IR (Conical cemented titan base non hex L20 and screw)	BSAD13231L20	30			10-12
SET ZZ BASE BK-KO-CCTB-L20 IW (Conical cemented titan base non hex L20 and screw)	BSAD13232L20	30			10-12
SET ZZ BASE BK-KO-CCTB-L30 5.7 (Conical cemented titan base non hex L30 and screw)	BSAD13233L30	30			10-12
SET ZZ BASE BK-KO-CCTB-L30 IR (Conical cemented titan base non hex L30 and screw)	BSAD13231L30	30			10-12
SET ZZ BASE BK-KO-CCTB-L30 IW (Conical cemented titan base non hex L30 and screw)	BSAD13232L30	30			10-12
SET ZZ BASE BK-KO-CCTB-L40 5.7 (Conical cemented titan base non hex L40 and screw)	BSAD13233L40	30			10-12
SET ZZ BASE BK-KO-CCTB-L40 IR (Conical cemented titan base non hex L40 and screw)	BSAD13231L40	30			10-12
SET ZZ BASE BK-KO-CCTB-L40 IW (Conical cemented titan base non hex L40 and screw)	BSAD13232L40	30			10-12
SET ZZ BASE BK-KO-NTB IR (Narrow titan base non hex and screw)	BSAP13231	30			10-12
SET ZZ BASE BK-KO-NTBH IR (Narrow titan base hex and screw)	BSAP13281	30			10-12
SET ZZ BASE BK-KO-PCTBH 5.7 (Parallel cemented titan base hex and screw)	BSAB13283	30			10-12
SET ZZ BASE BK-KO-PCTBH IR (Parallel cemented titan base hex and screw)	BSAB13281	30			10-12
SET ZZ BASE BK-KO-PCTBH IW (Parallel cemented titan base hex and screw)	BSAB13282	30			10-12
SET ZZ BASE BK-KO-PCTBH-K85 IR (Parallel cemented titan base hex K85 and screw)	BSAB13281K85	30			10-12
SET ZZ BASE BK-KO-PCTBH-K85 IW (Parallel cemented titan base hex K85 and screw)	BSAB13282K85	30			10-12
SET ZZ BASE BK-KO-PCTBH-L10 5.7 (Parallel cemented titan base hex L10 and screw)	BSAB13293L10	30			10-12

BTK THE SMILE SYSTEM® - BT-KONIC

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE BK-KO-PCTBH-L10 IR (Parallel cemented titan base hex L10 and screw)	BSAB13291L10	30			10-12
SET ZZ BASE BK-KO-PCTBH-L10 IW (Parallel cemented titan base hex L10 and screw)	BSAB13292L10	30			10-12
SET ZZ BASE BK-KO-PCTBH-L20 5.7 (Parallel cemented titan base hex L20 and screw)	BSAB13293L20	30			10-12
SET ZZ BASE BK-KO-PCTBH-L20 IR (Parallel cemented titan base hex L20 and screw)	BSAB13291L20	30			10-12
SET ZZ BASE BK-KO-PCTBH-L20 IW (Parallel cemented titan base hex L20 and screw)	BSAB13292L20	30			10-12
SET ZZ BASE BK-KO-PCTBH-L30 5.7 (Parallel cemented titan base hex L30 and screw)	BSAB13293L30	30			10-12
SET ZZ BASE BK-KO-PCTBH-L30 IR (Parallel cemented titan base hex L30 and screw)	BSAB13291L30	30			10-12
SET ZZ BASE BK-KO-PCTBH-L30 IW (Parallel cemented titan base hex L30 and screw)	BSAB13292L30	30			10-12
SET ZZ BASE BK-KO-PCTBH-L40 5.7 (Parallel cemented titan base hex L40 and screw)	BSAB13293L40	30			10-12
SET ZZ BASE BK-KO-PCTBH-L40 IR (Parallel cemented titan base hex L40 and screw)	BSAB13291L40	30			10-12
SET ZZ BASE BK-KO-PCTBH-L40 IW (Parallel cemented titan base hex L40 and screw)	BSAB13292L40	30			10-12
SET ZZ BASE BK-KO-SA IR (Scananalog and screw)	ZBDE13231	15			10-12
SET ZZ BASE BK-KO-SA IW (Scananalog and screw)	ZBDE13232	15			10-12
SET ZZ BASE BK-KO-SM IR (Scanmarker and screw)	ZBAD13231	30			10-12
SET ZZ BASE BK-KO-SM IW (Scanmarker and screw)	ZBAD13232	30			10-12
SET ZZ BASE BK-KO-WSM IW (White Scanmarker and screw)	ZBCD13232	30			10-12
SET ZZ BASE BK-KO-ZZMUA-L05 IR (Zirkonzahn Multi Unit Abutment non hex L05 and Screw)	BSAH13231L05	25	30		10-12
SET ZZ BASE BK-KO-ZZMUA-L15 IR (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH13231L15	25	30		10-12
SET ZZ BASE BK-KO-ZZMUA-L15 IW (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH13232L15	25	30		10-12
SET ZZ BASE BK-KO-ZZMUA-L25 IR (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH13231L25	25	30		10-12
SET ZZ BASE BK-KO-ZZMUA-L25 IW (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH13232L25	25	30		10-12
SET ZZ BASE BK-KO-ZZMUA-L35 IR (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH13231L35	25	30		10-12
SET ZZ BASE BK-KO-ZZMUA-L35 IW (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH13232L35	25	30		10-12
SET ZZ BASE BK-KO-ZZMUA-L45 IR (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH13231L45	25	30		10-12
SET ZZ BASE BK-KO-ZZMUA-L45 IW (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH13232L45	25	30		10-12
ZZ BASE BK-KO-MI IR (Model implant)	MIAG13201				10-12
ZZ BASE BK-KO-MI IW (Model implant)	MIAG13202				10-12
ZZ BASE BK-KO-SC 5.7 (Screw)	BSAE13203	30			10-12
ZZ BASE BK-KO-SC IR (Screw)	BSAE13201	30			10-12

BTK THE SMILE SYSTEM® - BT-KONIC

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE BK-KO-SC IW (Screw)	BSAE13202	30			10-12
ZZ BASE BK-KO-WMSM-L12 IR (White Metal Scanmarker L12 with Screw)	ZBED13231L12	15			10-12
ZZ BASE BK-KO-ZSC 5.7 (Zirconia screw)	BSAE13233	30			10-12
ZZ BASE BK-KO-ZSC IR (Zirconia screw)	BSAE13231	30			10-12
ZZ BASE BK-KO-ZSC IW (Zirconia screw)	BSAE13232	30			10-12
ZZ BASE BK-KO-ZSC-G IR (Zirconia Screw Gold)	BSCE13231	30			10-12
ZZ BASE BK-KO-ZSC-G IW (Zirconia Screw Gold)	BSCE13232	30			10-12

BTK THE SMILE SYSTEM® - BT-SAFE TL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE BK-SA-ABCH KW (Raw-Abutment® customizable hex and screw)	BSAF13481	35			10-12
SET ZZ BASE BK-SA-ABCH-D14 KW (Raw-Abutment® customizable hex D14 and screw)	BSAF13481D14	35			10-12
SET ZZ BASE BK-SA-CCTB KW (Conical cemented titan base non hex and screw)	BSAD13431	35			10-12
SET ZZ BASE BK-SA-CCTB-K85 KW (Conical cemented titan base non hex K85 and screw)	BSAD13431K85	35			10-12
SET ZZ BASE BK-SA-CCTB-L15 KW (Conical cemented titan base non hex L15 and screw)	BSAD13431L15	35			10-12
SET ZZ BASE BK-SA-PCTBH KW (Parallel cemented titan base hex and screw)	BSAB13471	35			10-12
SET ZZ BASE BK-SA-PCTBH KW (Parallel cemented titan base hex and screw)	BSAB13481	35			10-12
SET ZZ BASE BK-SA-PCTBH-K85 KW (Parallel cemented titan base hex K85 and screw)	BSAB13481K85	35			10-12
SET ZZ BASE BK-SA-PCTBH-L15 KW (Parallel cemented titan base hex L15 and screw)	BSAB13491L15	35			10-12
SET ZZ BASE BK-SA-SA KW (Scananalog and screw)	ZBDE13431	15			10-12
SET ZZ BASE BK-SA-SM KW (Scanmarker and screw)	ZBAD13431	35			10-12
SET ZZ BASE BK-SA-WSM KW (White Scanmarker and screw)	ZBCD13431	35			10-12
SET ZZ BASE BK-SA-ZZMUA17 KW Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and impl)	BSAI13431	30			10-12
SET ZZ BASE BK-SA-ZZMUA17 KW Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and impl)	BSAI13441	30			10-12
SET ZZ BASE BK-SA-ZZMUA-L15 KW (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH13431L15	25	35		10-12
SET ZZ BASE BK-SA-ZZMUA-L25 KW (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH13431L25	25	35		10-12
SET ZZ BASE BK-SA-ZZMUA-L35 KW (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH13431L35	25	35		10-12
SET ZZ BASE BK-SA-ZZMUA-L45 KW (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH13431L45	25	35		10-12
ZZ BASE BK-SA-BLZSC KW (Black Long Zirconia Screw)	BSIE13481	35			10-12
ZZ BASE BK-SA-BSC KW (Black Screw)	BSIE13401	35			10-12
ZZ BASE BK-SA-BZSC KW (Black Zirconia screw)	BSIE13431	35			10-12
ZZ BASE BK-SA-LBLZSC KW (Lab Long Zirconia Screw)	BSHE13481	35			10-12
ZZ BASE BK-SA-LBSC KW (Lab Screw)	BSHE13401	35			10-12
ZZ BASE BK-SA-LBZSC KW (Lab Zirconia screw)	BSHE13431	35			10-12
ZZ BASE BK-SA-LZSC KW (Long zirconia screw)	BSAE13481	35			10-12
ZZ BASE BK-SA-LZSC-G KW (Long Zirconia Screw Gold)	BSCE13481	35			10-12
ZZ BASE BK-SA-MI KW (Model implant)	MIAG13401				10-12
ZZ BASE BK-SA-SC KW (Screw)	BSAE13401	35			10-12
ZZ BASE BK-SA-WMSM-L12 KW (White Metal Scanmarker L12 with Screw)	ZBED13431L12	15			10-12

BTK THE SMILE SYSTEM® - BT-SAFE TL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
ZZ BASE BK-SA-WMSM-L15 KW (White Metal Scanmarker L15 with Screw)	ZBED13431L15	15			10-12
ZZ BASE BK-SA-ZSC KW (Zirconia screw)	BSAE13431	35			10-12
ZZ BASE BK-SA-ZSC-G KW (Zirconia Screw Gold)	BSCE13431	35			10-12
ZZ BASE BK-SA-ZZMUASC KW (Zirkonzahn Multi Unit Abutment Screw)	BSEE13401	30			10-12

CAMLOG® BAR ABUTMENTS (COMFOUR®)

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
SET ZZ BASE C-BA-CCTB 3.3/3.8/4.3 (Conical cemented titan base non hex and screw)	BSAD2131	20			10-12
SET ZZ BASE C-BA-CCTB 5.0/6.0 (Conical cemented titan base non hex and screw)	BSAD2132	35			10-12
SET ZZ BASE C-BA-CCTB-K85 3.3/3.8/4.3 (Conical cemented titan base non hex K85 and screw)	BSAD2131K85	20			10-12
SET ZZ BASE C-BA-CCTB-K85 5.0/6.0 (Conical cemented titan base non hex K85 and screw)	BSAD2132K85	35			10-12
SET ZZ BASE C-BA-ICOT 3.3/3.8/4.3 (Impression coping open tray non hex and Screw and Driver)	BSAQ2161	15			10-12
SET ZZ BASE C-BA-ICOT 5.0/6.0 (Impression coping open tray non hex and Screw and Driver)	BSAQ2162	15			10-12
SET ZZ BASE C-BA-SA 3.3/3.8/4.3 (Scanalog and screw)	ZBDE2131	15			10-12
SET ZZ BASE C-BA-SA 5.0/6.0 (Scanalog and screw)	ZBDE2132	15			10-12
SET ZZ BASE C-BA-SM 3.3/3.8/4.3 (Scanmarker and screw)	ZBAD2131	20			10-12
SET ZZ BASE C-BA-SM 5.0/6.0 (Scanmarker and screw)	ZBAD2132	35			10-12
SET ZZ BASE C-BA-WSM 3.3/3.8/4.3 (White Scanmarker and screw)	ZBCD2131	20			10-12
SET ZZ BASE C-BA-WSM 5.0/6.0 (White Scanmarker and screw)	ZBCD2132	35			10-12
ZZ BASE C-BA-ICSC 3.3/3.8/4.3 (Impression coping Screw)	BSFE2101	15			10-12
ZZ BASE C-BA-ICSC 5.0/6.0 (Impression coping Screw)	BSFE2102	15			10-12
ZZ BASE C-BA-SC 3.3/3.8/4.3 (Screw)	BSAE2101	20			10-12
ZZ BASE C-BA-SC 5.0/6.0 (Screw)	BSAE2102	35			10-12
ZZ BASE C-BA-WMSM-L12 3.3/3.8/4.3 (White Metal Scanmarker L12 with Screw)	ZBED2131L12	15			10-12
ZZ BASE C-BA-WMSM-L12 5.0/6.0 (White Metal Scanmarker L12 with Screw)	ZBED2132L12	15			10-12
ZZ BASE C-BA-ZSC 3.3/3.8/4.3 (Zirconia screw)	BSAE2131	20			10-12
ZZ BASE C-BA-ZSC 5.0/6.0 (Zirconia screw)	BSAE2132	35			10-12

CAMLOG® CERALOG®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE C-CE-ICOTH Hexalobe (Impression coping open tray hex and Screw and Driver)	BSAQ12741	15			10-12
SET ZZ BASE C-CE-SM Hexalobe (Scanmarker and screw)	ZBAD12731	25			10-12
SET ZZ BASE C-CE-WSM Hexalobe (White Scanmarker and screw)	ZBCD12731	25			10-12
ZZ BASE C-CE-ICSC Hexalobe (Impression coping Screw)	BSFE12701	15			10-12
ZZ BASE C-CE-MI Hexalobe (Model implant)	MIAG12701				10-12
ZZ BASE C-CE-SC Hexalobe (Screw)	BSAE12701	25			10-12
ZZ BASE C-CE-WMSM-L12 Hexalobe (White Metal Scanmarker L12 with Screw)	ZBED12731L12	15			10-12
ZZ BASE C-CE-ZSC Hexalobe (Zirconia Screw)	BSAE12731	25			10-12

CAMLOG® CONELOG®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
C-CL-EX-A15 3.3/3.8/4.3 (Titanium Base Extractor A15)	EXAB4701A15				10-12
C-CL-EX-A15 5.0 (Titanium Base Extractor A15)	EXAB4704A15				10-12
SET ZZ BASE C-CL-ABCH 3.3 (Raw-Abutment® customizable hex and screw)	BSAF4781	20			10-12
SET ZZ BASE C-CL-ABCH 3.8 (Raw-Abutment® customizable hex and screw)	BSAF4782	20			10-12
SET ZZ BASE C-CL-ABCH 4.3 (Raw-Abutment® customizable hex and screw)	BSAF4783	20			10-12
SET ZZ BASE C-CL-ABCH 5.0 (Raw-Abutment® customizable hex and screw)	BSAF4784	20			10-12
SET ZZ BASE C-CL-ASCTBH-K80 3.3 (Angled screw-channel titan base hex K80 and screw)	BSAT4781K80	20			10-12
SET ZZ BASE C-CL-ASCTBH-K80 3.3 Type 2 (Angled screw-channel titan base hex K80 Type 2 a	BSAT4791K80	20			10-12
SET ZZ BASE C-CL-ASCTBH-K80 3.8 (Angled screw-channel titan base hex K80 and screw)	BSAT4782K80	20			10-12
SET ZZ BASE C-CL-ASCTBH-K80 3.8 Type 2 (Angled screw-channel titan base hex K80 Type 2 a	BSAT4792K80	20			10-12
SET ZZ BASE C-CL-ASCTBH-K80 4.3 (Angled screw-channel titan base hex K80 and screw)	BSAT4783K80	20			10-12
SET ZZ BASE C-CL-ASCTBH-K80 4.3 Type 2 (Angled screw-channel titan base hex K80 Type 2 a	BSAT4793K80	20			10-12
SET ZZ BASE C-CL-ASCTBH-K80 5.0 (Angled screw-channel titan base hex K80 and screw)	BSAT4784K80	20			10-12
SET ZZ BASE C-CL-ASCTB-K80 3.3 (Angled screw-channel titan base non hex K80 and screw)	BSAT4731K80	20			10-12
SET ZZ BASE C-CL-ASCTB-K80 3.8 (Angled screw-channel titan base non hex K80 and screw)	BSAT4732K80	20			10-12
SET ZZ BASE C-CL-ASCTB-K80 4.3 (Angled screw-channel titan base non hex K80 and screw)	BSAT4733K80	20			10-12
SET ZZ BASE C-CL-ASCTB-K80 5.0 (Angled screw-channel titan base non hex K80 and screw)	BSAT4734K80	20			10-12
SET ZZ BASE C-CL-CCTB 3.3 (Conical cemented titan base non hex and screw)	BSAD4731	20			10-12
SET ZZ BASE C-CL-CCTB 3.8 (Conical cemented titan base non hex and screw)	BSAD4732	20			10-12
SET ZZ BASE C-CL-CCTB 4.3 (Conical cemented titan base non hex and screw)	BSAD4733	20			10-12
SET ZZ BASE C-CL-CCTB 5.0 (Conical cemented titan base non hex and screw)	BSAD4734	20			10-12
SET ZZ BASE C-CL-CCTB-K85 3.3 (Conical cemented titan base non hex K85 and screw)	BSAD4731K85	20			10-12
SET ZZ BASE C-CL-CCTB-K85 3.8 (Conical cemented titan base non hex K85 and screw)	BSAD4732K85	20			10-12
SET ZZ BASE C-CL-CCTB-K85 4.3 (Conical cemented titan base non hex K85 and screw)	BSAD4733K85	20			10-12
SET ZZ BASE C-CL-CCTB-K85 5.0 (Conical cemented titan base non hex K85 and screw)	BSAD4734K85	20			10-12
SET ZZ BASE C-CL-CCTB-L15 3.3 (Conical cemented titan base non hex L15 and screw)	BSAD4731L15	20			10-12
SET ZZ BASE C-CL-CCTB-L15 3.8 (Conical cemented titan base non hex L15 and screw)	BSAD4732L15	20			10-12
SET ZZ BASE C-CL-CCTB-L15 4.3 (Conical cemented titan base non hex L15 and screw)	BSAD4733L15	20			10-12
SET ZZ BASE C-CL-CCTB-L15 5.0 (Conical cemented titan base non hex L15 and screw)	BSAD4734L15	20			10-12

CAMLOG® CONELOG®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE C-CL-CCTB-L20 3.3 (Conical cemented titan base non hex L20 and screw)	BSAD4731L20	20			10-12
SET ZZ BASE C-CL-CCTB-L20 3.8 (Conical cemented titan base non hex L20 and screw)	BSAD4732L20	20			10-12
SET ZZ BASE C-CL-CCTB-L20 4.3 (Conical cemented titan base non hex L20 and screw)	BSAD4733L20	20			10-12
SET ZZ BASE C-CL-CCTB-L20 5.0 (Conical cemented titan base non hex L20 and screw)	BSAD4734L20	20			10-12
SET ZZ BASE C-CL-CCTB-L30 3.3 (Conical cemented titan base non hex L30 and screw)	BSAD4731L30	20			10-12
SET ZZ BASE C-CL-CCTB-L30 3.8 (Conical cemented titan base non hex L30 and screw)	BSAD4732L30	20			10-12
SET ZZ BASE C-CL-CCTB-L30 4.3 (Conical cemented titan base non hex L30 and screw)	BSAD4733L30	20			10-12
SET ZZ BASE C-CL-CCTB-L30 5.0 (Conical cemented titan base non hex L30 and screw)	BSAD4734L30	20			10-12
SET ZZ BASE C-CL-CCTB-L40 3.3 (Conical cemented titan base non hex L40 and screw)	BSAD4731L40	20			10-12
SET ZZ BASE C-CL-CCTB-L40 3.8 (Conical cemented titan base non hex L40 and screw)	BSAD4732L40	20			10-12
SET ZZ BASE C-CL-CCTB-L40 4.3 (Conical cemented titan base non hex L40 and screw)	BSAD4733L40	20			10-12
SET ZZ BASE C-CL-CCTB-L40 5.0 (Conical cemented titan base non hex L40 and screw)	BSAD4734L40	20			10-12
SET ZZ BASE C-CL-ICOTH 3.3 (Impression coping open tray hex and Screw and Driver)	BSAQ4741	15			10-12
SET ZZ BASE C-CL-ICOTH 3.8 (Impression coping open tray hex and Screw and Driver)	BSAQ4742	15			10-12
SET ZZ BASE C-CL-ICOTH 4.3 (Impression coping open tray hex and Screw and Driver)	BSAQ4743	15			10-12
SET ZZ BASE C-CL-ICOTH 5.0 (Impression coping open tray hex and Screw and Driver)	BSAQ4744	15			10-12
SET ZZ BASE C-CL-NTB 3.3 (Narrow titan base non hex and screw)	BSAP4731	20			10-12
SET ZZ BASE C-CL-NTB 3.8 (Narrow titan base non hex and screw)	BSAP4732	20			10-12
SET ZZ BASE C-CL-NTB 4.3 (Narrow titan base non hex and screw)	BSAP4733	20			10-12
SET ZZ BASE C-CL-NTBH 3.3 (Narrow titan base hex and screw)	BSAP4781	20			10-12
SET ZZ BASE C-CL-NTBH 3.8 (Narrow titan base hex and screw)	BSAP4782	20			10-12
SET ZZ BASE C-CL-NTBH 4.3 (Narrow titan base hex and screw)	BSAP4783	20			10-12
SET ZZ BASE C-CL-PCTBH 3.3 (Parallel cemented titan base hex and screw)	BSAB4791	20			10-12
SET ZZ BASE C-CL-PCTBH 3.8 (Parallel cemented titan base hex and screw)	BSAB4782	20			10-12
SET ZZ BASE C-CL-PCTBH 3.8 (Parallel cemented titan base hex and screw)	BSAB4792	20			10-12
SET ZZ BASE C-CL-PCTBH 4.3 (Parallel cemented titan base hex and screw)	BSAB4783	20			10-12
SET ZZ BASE C-CL-PCTBH 4.3 (Parallel cemented titan base hex and screw)	BSAB4793	20			10-12
SET ZZ BASE C-CL-PCTBH 5.0 (Parallel cemented titan base hex and screw)	BSAB4784	20			10-12
SET ZZ BASE C-CL-PCTBH 5.0 (Parallel cemented titan base hex and screw)	BSAB4794	20			10-12

CAMLOG® CONELOG®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE C-CL-PCTBH-K85 3.3 (Parallel cemented titan base hex K85 and screw)	BSAB4781K85	20			10-12
SET ZZ BASE C-CL-PCTBH-K85 3.8 (Parallel cemented titan base hex K85 and screw)	BSAB4782K85	20			10-12
SET ZZ BASE C-CL-PCTBH-K85 4.3 (Parallel cemented titan base hex K85 and screw)	BSAB4783K85	20			10-12
SET ZZ BASE C-CL-PCTBH-K85 5.0 (Parallel cemented titan base hex K85 and screw)	BSAB4784K85	20			10-12
SET ZZ BASE C-CL-PCTBH-L15 3.3 (Parallel cemented titan base hex L15 and screw)	BSAB4791L15	20			10-12
SET ZZ BASE C-CL-PCTBH-L15 3.8 (Parallel cemented titan base hex L15 and screw)	BSAB4792L15	20			10-12
SET ZZ BASE C-CL-PCTBH-L15 4.3 (Parallel cemented titan base hex L15 and screw)	BSAB4793L15	20			10-12
SET ZZ BASE C-CL-PCTBH-L15 5.0 (Parallel cemented titan base hex L15 and screw)	BSAB4794L15	20			10-12
SET ZZ BASE C-CL-PCTBH-L20 3.3 (Parallel cemented titan base hex L20 and screw)	BSAB4791L20	20			10-12
SET ZZ BASE C-CL-PCTBH-L20 3.8 (Parallel cemented titan base hex L20 and screw)	BSAB4792L20	20			10-12
SET ZZ BASE C-CL-PCTBH-L20 4.3 (Parallel cemented titan base hex L20 and screw)	BSAB4793L20	20			10-12
SET ZZ BASE C-CL-PCTBH-L20 5.0 (Parallel cemented titan base hex L20 and screw)	BSAB4794L20	20			10-12
SET ZZ BASE C-CL-PCTBH-L30 3.3 (Parallel cemented titan base hex L30 and screw)	BSAB4791L30	20			10-12
SET ZZ BASE C-CL-PCTBH-L30 3.8 (Parallel cemented titan base hex L30 and screw)	BSAB4792L30	20			10-12
SET ZZ BASE C-CL-PCTBH-L30 4.3 (Parallel cemented titan base hex L30 and screw)	BSAB4793L30	20			10-12
SET ZZ BASE C-CL-PCTBH-L30 5.0 (Parallel cemented titan base hex L30 and screw)	BSAB4794L30	20			10-12
SET ZZ BASE C-CL-PCTBH-L40 3.3 (Parallel cemented titan base hex L40 and screw)	BSAB4791L40	20			10-12
SET ZZ BASE C-CL-PCTBH-L40 3.8 (Parallel cemented titan base hex L40 and screw)	BSAB4792L40	20			10-12
SET ZZ BASE C-CL-PCTBH-L40 4.3 (Parallel cemented titan base hex L40 and screw)	BSAB4793L40	20			10-12
SET ZZ BASE C-CL-PCTBH-L40 5.0 (Parallel cemented titan base hex L40 and screw)	BSAB4794L40	20			10-12
SET ZZ BASE C-CL-SA 3.3 (Scanalog and screw)	ZBDE4731	15			10-12
SET ZZ BASE C-CL-SA 3.8 (Scanalog and screw)	ZBDE4732	15			10-12
SET ZZ BASE C-CL-SA 4.3 (Scanalog and screw)	ZBDE4733	15			10-12
SET ZZ BASE C-CL-SA 5.0 (Scanalog and screw)	ZBDE4734	15			10-12
SET ZZ BASE C-CL-SM 3.3 (Scanmarker and screw)	ZBAD4731	20			10-12
SET ZZ BASE C-CL-SM 3.8 (Scanmarker and screw)	ZBAD4732	20			10-12
SET ZZ BASE C-CL-SM 4.3 (Scanmarker and screw)	ZBAD4733	20			10-12
SET ZZ BASE C-CL-SM 5.0 (Scanmarker and screw)	ZBAD4734	20			10-12
SET ZZ BASE C-CL-WSM 3.3 (White Scanmarker and screw)	ZBCD4731	20			10-12

CAMLOG® CONELOG®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE C-CL-WSM 3.8 (White Scanmarker and screw)	ZBCD4732	20			10-12
SET ZZ BASE C-CL-WSM 4.3 (White Scanmarker and screw)	ZBCD4733	20			10-12
SET ZZ BASE C-CL-WSM 5.0 (White Scanmarker and screw)	ZBCD4734	20			10-12
SET ZZ BASE C-CL-ZZMUA-L15 3.3 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH4731L15	25	30		10-12
SET ZZ BASE C-CL-ZZMUA-L15 3.8 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH4732L15	25	30		10-12
SET ZZ BASE C-CL-ZZMUA-L15 4.3 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH4733L15	25	30		10-12
SET ZZ BASE C-CL-ZZMUA-L15 5.0 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH4734L15	25	30		10-12
SET ZZ BASE C-CL-ZZMUA-L25 3.3 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH4731L25	25	30		10-12
SET ZZ BASE C-CL-ZZMUA-L25 3.8 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH4732L25	25	30		10-12
SET ZZ BASE C-CL-ZZMUA-L25 4.3 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH4733L25	25	30		10-12
SET ZZ BASE C-CL-ZZMUA-L25 5.0 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH4734L25	25	30		10-12
SET ZZ BASE C-CL-ZZMUA-L35 3.3 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH4731L35	25	30		10-12
SET ZZ BASE C-CL-ZZMUA-L35 3.8 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH4732L35	25	30		10-12
SET ZZ BASE C-CL-ZZMUA-L35 4.3 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH4733L35	25	30		10-12
SET ZZ BASE C-CL-ZZMUA-L35 5.0 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH4734L35	25	30		10-12
SET ZZ BASE C-CL-ZZMUA-L45 3.3 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH4731L45	25	30		10-12
SET ZZ BASE C-CL-ZZMUA-L45 3.8 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH4732L45	25	30		10-12
SET ZZ BASE C-CL-ZZMUA-L45 4.3 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH4733L45	25	30		10-12
SET ZZ BASE C-CL-ZZMUA-L45 5.0 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH4734L45	25	30		10-12
ZZ BASE C-CL-ASCSC 3.3 (Angled screw-channel Screw)	BSKE4701	20			10-12
ZZ BASE C-CL-ASCSC 3.8 (Angled screw-channel Screw)	BSKE4702	20			10-12
ZZ BASE C-CL-ASCSC 4.3 (Angled screw-channel Screw)	BSKE4703	20			10-12
ZZ BASE C-CL-ASCSC 5.0 (Angled screw-channel Screw)	BSKE4704	20			10-12
ZZ BASE C-CL-ICSC 3.3 (Impression coping Screw)	BSFE4701	15			10-12
ZZ BASE C-CL-ICSC 3.8 (Impression coping Screw)	BSFE4702	15			10-12
ZZ BASE C-CL-ICSC 4.3 (Impression coping Screw)	BSFE4703	15			10-12
ZZ BASE C-CL-ICSC 5.0 (Impression coping Screw)	BSFE4704	15			10-12
ZZ BASE C-CL-MI 3.3 (Model implant)	MIAG4701				10-12
ZZ BASE C-CL-MI 3.8 (Model implant)	MIAG4702				10-12

All information is subject to change. Errors and omissions excepted. Version: 21.09.2023

CAMLOG® CONELOG®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE C-CL-MI 4.3 (Model implant)	MIAG4703				10-12
ZZ BASE C-CL-MI 5.0 (Model implant)	MIAG4704				10-12
ZZ BASE C-CL-SC 3.3 (Screw)	BSAE4701	20			10-12
ZZ BASE C-CL-SC 3.8 (Screw)	BSAE4702	20			10-12
ZZ BASE C-CL-SC 4.3 (Screw)	BSAE4703	20			10-12
ZZ BASE C-CL-SC 5.0 (Screw)	BSAE4704	20			10-12
ZZ BASE C-CL-WMSM-L10 3.3 (White Metal Scanmarker L10 with Screw)	ZBED4731L10	15			10-12
ZZ BASE C-CL-WMSM-L10 3.8 (White Metal Scanmarker L10 with Screw)	ZBED4732L10	15			10-12
ZZ BASE C-CL-WMSM-L10 4.3 (White Metal Scanmarker L10 with Screw)	ZBED4733L10	15			10-12
ZZ BASE C-CL-WMSM-L10 5.0 (White Metal Scanmarker L10 with Screw)	ZBED4734L10	15			10-12
ZZ BASE C-CL-WMSM-L12 3.3 (White Metal Scanmarker L12 with Screw)	ZBED4731L12	15			10-12
ZZ BASE C-CL-WMSM-L12 3.8 (White Metal Scanmarker L12 with Screw)	ZBED4732L12	15			10-12
ZZ BASE C-CL-WMSM-L12 4.3 (White Metal Scanmarker L12 with Screw)	ZBED4733L12	15			10-12
ZZ BASE C-CL-WMSM-L12 5.0 (White Metal Scanmarker L12 with Screw)	ZBED4734L12	15			10-12
ZZ BASE C-CL-WMSM-L15 3.3 (White Metal Scanmarker L15 with Screw)	ZBED4731L15	15			10-12
ZZ BASE C-CL-WMSM-L15 3.8 (White Metal Scanmarker L15 with Screw)	ZBED4732L15	15			10-12
ZZ BASE C-CL-WMSM-L15 4.3 (White Metal Scanmarker L15 with Screw)	ZBED4733L15	15			10-12
ZZ BASE C-CL-WMSM-L15 5.0 (White Metal Scanmarker L15 with Screw)	ZBED4734L15	15			10-12
ZZ BASE C-CL-ZSC 3.3 (Zirconia screw)	BSAE4731	20			10-12
ZZ BASE C-CL-ZSC 3.8 (Zirconia screw)	BSAE4732	20			10-12
ZZ BASE C-CL-ZSC 4.3 (Zirconia screw)	BSAE4733	20			10-12
ZZ BASE C-CL-ZSC 5.0 (Zirconia screw)	BSAE4734	20			10-12

CAMLOG® J-LINE/K-LINE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE C-J-SM 3.3 (Scanmarker and screw)	ZBAD2031	20			10-12
SET ZZ BASE C-J-SM 3.8 (Scanmarker and screw)	ZBAD2032	20			10-12
SET ZZ BASE C-J-SM 4.3 (Scanmarker and screw)	ZBAD2033	20			10-12
SET ZZ BASE C-J-SM 5. 0 (Scanmarker and screw)	ZBAD2034	20			10-12
SET ZZ BASE C-J-ABCH 3.3 (Raw-Abutment® customizable hex and screw)	BSAF2091	20			10-12
SET ZZ BASE C-J-ABCH 3.8 (Raw-Abutment® customizable hex and screw)	BSAF2092	20			10-12
SET ZZ BASE C-J-ABCH 4.3 (Raw-Abutment® customizable hex and screw)	BSAF2093	20			10-12
SET ZZ BASE C-J-ABCH 5.0 (Raw-Abutment® customizable hex and screw)	BSAF2094	20			10-12
SET ZZ BASE C-J-ABCH 6.0 (Raw-Abutment® customizable hex and screw)	BSAF2095	20			10-12
SET ZZ BASE C-J-ABCH-D14 3.3 (Raw-Abutment® customizable hex D14 and screw)	BSAF2091D14	20			10-12
SET ZZ BASE C-J-ABCH-D14 3.8 (Raw-Abutment® customizable hex D14 and screw)	BSAF2092D14	20			10-12
SET ZZ BASE C-J-ABCH-D14 4.3 (Raw-Abutment® customizable hex D14 and screw)	BSAF2093D14	20			10-12
SET ZZ BASE C-J-ABCH-D14 5.0 (Raw-Abutment® customizable hex D14 and screw)	BSAF2094D14	20			10-12
SET ZZ BASE C-J-ABCH-D14 6.0 (Raw-Abutment® customizable hex D14 and screw)	BSAF2095D14	20			10-12
SET ZZ BASE C-J-ASCTBH-K80 3.3 (Angled screw-channel titan base hex K80 and screw)	BSAT2081K80	15			10-12
SET ZZ BASE C-J-ASCTBH-K80 3.3 Type 2 (Angled screw-channel titan base hex K80 Type 2 and screw)	BSAT2091K80	15			10-12
SET ZZ BASE C-J-ASCTBH-K80 3.8 (Angled screw-channel titan base hex K80 and screw)	BSAT2082K80	15			10-12
SET ZZ BASE C-J-ASCTBH-K80 3.8 Type 2 (Angled screw-channel titan base hex K80 Type 2 and screw)	BSAT2092K80	15			10-12
SET ZZ BASE C-J-ASCTBH-K80 4.3 (Angled screw-channel titan base hex K80 and screw)	BSAT2083K80	15			10-12
SET ZZ BASE C-J-ASCTBH-K80 4.3 Type 2 (Angled screw-channel titan base hex K80 Type 2 and screw)	BSAT2093K80	15			10-12
SET ZZ BASE C-J-ASCTB-K80 3.3 (Angled screw-channel titan base non hex K80 and screw)	BSAT2031K80	15			10-12
SET ZZ BASE C-J-ASCTB-K80 3.8 (Angled screw-channel titan base non hex K80 and screw)	BSAT2032K80	15			10-12
SET ZZ BASE C-J-ASCTB-K80 4.3 (Angled screw-channel titan base non hex K80 and screw)	BSAT2033K80	15			10-12
SET ZZ BASE C-J-CCTB 3.3 (Conical cemented titan base non hex and screw)	BSAD2031	20			10-12
SET ZZ BASE C-J-CCTB 3.8 (Conical cemented titan base non hex and screw)	BSAD2032	20			10-12
SET ZZ BASE C-J-CCTB 4.3 (Conical cemented titan base non hex and screw)	BSAD2033	20			10-12
SET ZZ BASE C-J-CCTB 5.0 (Conical cemented titan base non hex and screw)	BSAD2034	20			10-12
SET ZZ BASE C-J-CCTB 6.0 (Conical cemented titan base non hex and screw)	BSAD2035	20			10-12
SET ZZ BASE C-J-CCTB-G 4.3 (Conical cemented titan base non hex and screw Gold)	BSCD2033	20			10-12

CAMLOG® J-LINE/K-LINE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE C-J-CCTB-G 5.0 (Conical cemented titan base non hex and screw Gold)	BSCD2034	20			10-12
SET ZZ BASE C-J-CCTB-G 6.0 (Conical cemented titan base non hex and screw Gold)	BSCD2035	20			10-12
SET ZZ BASE C-J-CCTB-K85 3.3 (Conical cemented titan base non hex K85 and screw)	BSAD2031K85	20			10-12
SET ZZ BASE C-J-CCTB-K85 3.8 (Conical cemented titan base non hex K85 and screw)	BSAD2032K85	20			10-12
SET ZZ BASE C-J-CCTB-K85 4.3 (Conical cemented titan base non hex K85 and screw)	BSAD2033K85	20			10-12
SET ZZ BASE C-J-CCTB-K85 5.0 (Conical cemented titan base non hex K85 and screw)	BSAD2034K85	20			10-12
SET ZZ BASE C-J-CCTB-K85 6.0 (Conical cemented titan base non hex K85 and screw)	BSAD2035K85	20			10-12
SET ZZ BASE C-J-CCTB-L10 3.3 (Conical cemented titan base non hex L10 and screw)	BSAD2031L10	20			10-12
SET ZZ BASE C-J-CCTB-L10 3.8 (Conical cemented titan base non hex L10 and screw)	BSAD2032L10	20			10-12
SET ZZ BASE C-J-CCTB-L10 4.3 (Conical cemented titan base non hex L10 and screw)	BSAD2033L10	20			10-12
SET ZZ BASE C-J-CCTB-L10 5.0 (Conical cemented titan base non hex L10 and screw)	BSAD2034L10	20			10-12
SET ZZ BASE C-J-CCTB-L10 6.0 (Conical cemented titan base non hex L10 and screw)	BSAD2035L10	20			10-12
SET ZZ BASE C-J-CCTB-L10-G 3.8 (Conical cemented titan base non hex L10 and screw Gold)	BSCD2032L10	20			10-12
SET ZZ BASE C-J-CCTB-L10-G 4.3 (Conical cemented titan base non hex L10 and screw Gold)	BSCD2033L10	20			10-12
SET ZZ BASE C-J-CCTB-L10-G 6.0 (Conical cemented titan base non hex L10 and screw Gold)	BSCD2035L10	20			10-12
SET ZZ BASE C-J-CCTB-L20 3.3 (Conical cemented titan base non hex L20 and screw)	BSAD2031L20	20			10-12
SET ZZ BASE C-J-CCTB-L20 3.8 (Conical cemented titan base non hex L20 and screw)	BSAD2032L20	20			10-12
SET ZZ BASE C-J-CCTB-L20 4.3 (Conical cemented titan base non hex L20 and screw)	BSAD2033L20	20			10-12
SET ZZ BASE C-J-CCTB-L20 5.0 (Conical cemented titan base non hex L20 and screw)	BSAD2034L20	20			10-12
SET ZZ BASE C-J-CCTB-L20 6.0 (Conical cemented titan base non hex L20 and screw)	BSAD2035L20	20			10-12
SET ZZ BASE C-J-CCTB-L20-G 3.8 (Conical cemented titan base non hex L20 and screw Gold)	BSCD2032L20	20			10-12
SET ZZ BASE C-J-CCTB-L20-G 4.3 (Conical cemented titan base non hex L20 and screw Gold)	BSCD2033L20	20			10-12
SET ZZ BASE C-J-CCTB-L20-G 6.0 (Conical cemented titan base non hex L20 and screw Gold)	BSCD2035L20	20			10-12
SET ZZ BASE C-J-CCTB-L30 3.3 (Conical cemented titan base non hex L30 and screw)	BSAD2031L30	20			10-12
SET ZZ BASE C-J-CCTB-L30 3.8 (Conical cemented titan base non hex L30 and screw)	BSAD2032L30	20			10-12
SET ZZ BASE C-J-CCTB-L30 4.3 (Conical cemented titan base non hex L30 and screw)	BSAD2033L30	20			10-12
SET ZZ BASE C-J-CCTB-L30 5.0 (Conical cemented titan base non hex L30 and screw)	BSAD2034L30	20			10-12
SET ZZ BASE C-J-CCTB-L30 6.0 (Conical cemented titan base non hex L30 and screw)	BSAD2035L30	20			10-12
SET ZZ BASE C-J-CCTB-L30-G 3.3 (Conical cemented titan base non hex L30 and screw Gold)	BSCD2031L30	20			10-12

CAMLOG® J-LINE/K-LINE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE C-J-CCTB-L30-G 3.8 (Conical cemented titan base non hex L30 and screw Gold)	BSCD2032L30	20			10-12
SET ZZ BASE C-J-CCTB-L30-G 4.3 (Conical cemented titan base non hex L30 and screw Gold)	BSCD2033L30	20			10-12
SET ZZ BASE C-J-CCTB-L30-G 5.0 (Conical cemented titan base non hex L30 and screw Gold)	BSCD2034L30	20			10-12
SET ZZ BASE C-J-CCTB-L30-G 6.0 (Conical cemented titan base non hex L30 and screw Gold)	BSCD2035L30	20			10-12
SET ZZ BASE C-J-ICOT 3.3 (Impression coping open tray non hex and Screw and Driver)	BSAQ2061	15			10-12
SET ZZ BASE C-J-ICOT 3.8 (Impression coping open tray non hex and Screw and Driver)	BSAQ2062	15			10-12
SET ZZ BASE C-J-ICOT 4.3 (Impression coping open tray non hex and Screw and Driver)	BSAQ2063	15			10-12
SET ZZ BASE C-J-ICOT 5.0 (Impression coping open tray non hex and Screw and Driver)	BSAQ2064	15			10-12
SET ZZ BASE C-J-ICOT 6.0 (Impression coping open tray non hex and Screw and Driver)	BSAQ2065	15			10-12
SET ZZ BASE C-J-ICOTH 3.3 (Impression coping open tray hex and Screw and Driver)	BSAQ2041	15			10-12
SET ZZ BASE C-J-ICOTH 3.8 (Impression coping open tray hex and Screw and Driver)	BSAQ2042	15			10-12
SET ZZ BASE C-J-ICOTH 4.3 (Impression coping open tray hex and Screw and Driver)	BSAQ2043	15			10-12
SET ZZ BASE C-J-ICOTH 5.0 (Impression coping open tray hex and Screw and Driver)	BSAQ2044	15			10-12
SET ZZ BASE C-J-ICOTH 6.0 (Impression coping open tray hex and Screw and Driver)	BSAQ2045	15			10-12
SET ZZ BASE C-J-NTB 3.3 (Narrow titan base non hex and screw)	BSAP2031	20			10-12
SET ZZ BASE C-J-NTB 3.8 (Narrow titan base non hex and screw)	BSAP2032	20			10-12
SET ZZ BASE C-J-NTBH 3.3 (Narrow titan base hex and screw)	BSAP2081	20			10-12
SET ZZ BASE C-J-NTBH 3.8 (Narrow titan base hex and screw)	BSAP2082	20			10-12
SET ZZ BASE C-J-PCTBH 3.3 (Parallel cemented titan base hex and screw)	BSAB2081	20			10-12
SET ZZ BASE C-J-PCTBH 3.8 (Parallel cemented titan base hex and screw)	BSAB2082	20			10-12
SET ZZ BASE C-J-PCTBH 4.3 (Parallel cemented titan base hex and screw)	BSAB2083	20			10-12
SET ZZ BASE C-J-PCTBH 5.0 (Parallel cemented titan base hex and screw)	BSAB2084	20			10-12
SET ZZ BASE C-J-PCTBH 6.0 (Parallel cemented titan base hex and screw)	BSAB2085	20			10-12
SET ZZ BASE C-J-PCTBH-G 3.3 (Parallel cemented titan base hex and screw Gold)	BSCB2081	20			10-12
SET ZZ BASE C-J-PCTBH-G 3.8 (Parallel cemented titan base hex and screw Gold)	BSCB2082	20			10-12
SET ZZ BASE C-J-PCTBH-G 5.0 (Parallel cemented titan base hex and screw Gold)	BSCB2084	30			10-12
SET ZZ BASE C-J-PCTBH-G 6.0 (Parallel cemented titan base hex and screw Gold)	BSCB2085	20			10-12
SET ZZ BASE C-J-PCTBH-K85 3.3 (Parallel cemented titan base hex K85 and screw)	BSAB2081K85	20			10-12
SET ZZ BASE C-J-PCTBH-K85 3.8 (Parallel cemented titan base hex K85 and screw)	BSAB2082K85	20			10-12

CAMLOG® J-LINE/K-LINE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE C-J-PCTBH-K85 4.3 (Parallel cemented titan base hex K85 and screw)	BSAB2083K85	20			10-12
SET ZZ BASE C-J-PCTBH-K85 5.0 (Parallel cemented titan base hex K85 and screw)	BSAB2084K85	20			10-12
SET ZZ BASE C-J-PCTBH-K85 6.0 (Parallel cemented titan base hex K85 and screw)	BSAB2085K85	20			10-12
SET ZZ BASE C-J-PCTBH-L10 3.3 (Parallel cemented titan base hex L10 and screw)	BSAB2091L10	20			10-12
SET ZZ BASE C-J-PCTBH-L10 3.8 (Parallel cemented titan base hex L10 and screw)	BSAB2092L10	20			10-12
SET ZZ BASE C-J-PCTBH-L10 4.3 (Parallel cemented titan base hex L10 and screw)	BSAB2093L10	20			10-12
SET ZZ BASE C-J-PCTBH-L10 5.0 (Parallel cemented titan base hex L10 and screw)	BSAB2094L10	20			10-12
SET ZZ BASE C-J-PCTBH-L10 6.0 (Parallel cemented titan base hex L10 and screw)	BSAB2095L10	20			10-12
SET ZZ BASE C-J-PCTBH-L10-G 3.3 (Parallel cemented titan base hex L10 and screw Gold)	BSCB2091L10	20			10-12
SET ZZ BASE C-J-PCTBH-L10-G 3.8 (Parallel cemented titan base hex L10 and screw Gold)	BSCB2092L10	20			10-12
SET ZZ BASE C-J-PCTBH-L10-G 4.3 (Parallel cemented titan base hex L10 and screw Gold)	BSCB2093L10	20			10-12
SET ZZ BASE C-J-PCTBH-L10-G 5.0 (Parallel cemented titan base hex L10 and screw Gold)	BSCB2094L10	20			10-12
SET ZZ BASE C-J-PCTBH-L10-G 6.0 (Parallel cemented titan base hex L10 and screw Gold)	BSCB2095L10	20			10-12
SET ZZ BASE C-J-PCTBH-L20 3.3 (Parallel cemented titan base hex L20 and screw)	BSAB2091L20	20			10-12
SET ZZ BASE C-J-PCTBH-L20 3.8 (Parallel cemented titan base hex L20 and screw)	BSAB2092L20	20			10-12
SET ZZ BASE C-J-PCTBH-L20 4.3 (Parallel cemented titan base hex L20 and screw)	BSAB2093L20	20			10-12
SET ZZ BASE C-J-PCTBH-L20 5.0 (Parallel cemented titan base hex L20 and screw)	BSAB2094L20	20			10-12
SET ZZ BASE C-J-PCTBH-L20 6.0 (Parallel cemented titan base hex L20 and screw)	BSAB2095L20	20			10-12
SET ZZ BASE C-J-PCTBH-L20-G 3.3 (Parallel cemented titan base hex L20 and screw Gold)	BSCB2091L20	20			10-12
SET ZZ BASE C-J-PCTBH-L20-G 3.8 (Parallel cemented titan base hex L20 and screw Gold)	BSCB2092L20	20			10-12
SET ZZ BASE C-J-PCTBH-L20-G 4.3 (Parallel cemented titan base hex L20 and screw Gold)	BSCB2093L20	20			10-12
SET ZZ BASE C-J-PCTBH-L20-G 5.0 (Parallel cemented titan base hex L20 and screw Gold)	BSCB2094L20	20			10-12
SET ZZ BASE C-J-PCTBH-L20-G 6.0 (Parallel cemented titan base hex L20 and screw Gold)	BSCB2095L20	20			10-12
SET ZZ BASE C-J-PCTBH-L30 3.3 (Parallel cemented titan base hex L30 and screw)	BSAB2091L30	20			10-12
SET ZZ BASE C-J-PCTBH-L30 3.8 (Parallel cemented titan base hex L30 and screw)	BSAB2092L30	20			10-12
SET ZZ BASE C-J-PCTBH-L30 4.3 (Parallel cemented titan base hex L30 and screw)	BSAB2093L30	20			10-12
SET ZZ BASE C-J-PCTBH-L30 5.0 (Parallel cemented titan base hex L30 and screw)	BSAB2094L30	20			10-12
SET ZZ BASE C-J-PCTBH-L30 6.0 (Parallel cemented titan base hex L30 and screw)	BSAB2095L30	20			10-12
SET ZZ BASE C-J-PCTBH-L30-G 3.3 (Parallel cemented titan base hex L30 and screw Gold)	BSCB2091L30	20			10-12

CAMLOG® J-LINE/K-LINE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE C-J-PCTBH-L30-G 3.8 (Parallel cemented titan base hex L30 and screw Gold)	BSCB2092L30	20			10-12
SET ZZ BASE C-J-PCTBH-L30-G 4.3 (Parallel cemented titan base hex L30 and screw Gold)	BSCB2093L30	20			10-12
SET ZZ BASE C-J-PCTBH-L30-G 5.0 (Parallel cemented titan base hex L30 and screw Gold)	BSCB2094L30	20			10-12
SET ZZ BASE C-J-PCTBH-L30-G 6.0 (Parallel cemented titan base hex L30 and screw Gold)	BSCB2095L30	20			10-12
SET ZZ BASE C-J-SA 3.3 (Scananalog and screw)	ZBDE2031	15			10-12
SET ZZ BASE C-J-SA 3.8 (Scananalog and screw)	ZBDE2032	15			10-12
SET ZZ BASE C-J-SA 4.3 (Scananalog and screw)	ZBDE2033	15			10-12
SET ZZ BASE C-J-SA 5.0 (Scananalog and screw)	ZBDE2034	15			10-12
SET ZZ BASE C-J-SA 6.0 (Scananalog and screw)	ZBDE2035	15			10-12
SET ZZ BASE C-J-SM 6.0 (Scanmarker and screw)	ZBAD2035	20			10-12
SET ZZ BASE C-J-WSM 3.3 (White Scanmarker and screw)	ZBCD2031	20			10-12
SET ZZ BASE C-J-WSM 3.8 (White Scanmarker and screw)	ZBCD2032	20			10-12
SET ZZ BASE C-J-WSM 5.0 (White Scanmarker and screw)	ZBCD2034	20			10-12
SET ZZ BASE C-J-WSM 6.0 (White Scanmarker and screw)	ZBCD2035	20			10-12
SET ZZ BASE C-J-ZZMUA17 3.3 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and implant)	BSAI2031	20			10-12
SET ZZ BASE C-J-ZZMUA17 3.3 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and implant)	BSAI2041	20			10-12
SET ZZ BASE C-J-ZZMUA17 3.8 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and implant)	BSAI2032	20			10-12
SET ZZ BASE C-J-ZZMUA17 3.8 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and implant)	BSAI2042	20			10-12
SET ZZ BASE C-J-ZZMUA17 4.3 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and implant)	BSAI2033	20			10-12
SET ZZ BASE C-J-ZZMUA17 4.3 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and implant)	BSAI2043	20			10-12
SET ZZ BASE C-J-ZZMUA17 5.0 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and implant)	BSAI2034	30			10-12
SET ZZ BASE C-J-ZZMUA17 5.0 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and implant)	BSAI2044	30			10-12
SET ZZ BASE C-J-ZZMUA17-L10 3.3 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10 anc	BSAI2031L10	20			10-12
SET ZZ BASE C-J-ZZMUA17-L10 3.3 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 L10 anc	BSAI2041L10	20			10-12
SET ZZ BASE C-J-ZZMUA17-L10 3.8 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10 anc	BSAI2032L10	20			10-12
SET ZZ BASE C-J-ZZMUA17-L10 3.8 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 L10 anc	BSAI2042L10	20			10-12
SET ZZ BASE C-J-ZZMUA17-L10 4.3 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10 anc	BSAI2033L10	20			10-12
SET ZZ BASE C-J-ZZMUA17-L10 4.3 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 L10 anc	BSAI2043L10	20			10-12
SET ZZ BASE C-J-ZZMUA17-L10 5.0 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10 anc	BSAI2034L10	30			10-12

CAMLOG® J-LINE/K-LINE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE C-J-ZZMUA17-L10 5.0 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 L10 and Screw)	BSAI2044L10	30			10-12
SET ZZ BASE C-J-ZZMUA-L05 3.8 (Zirkonzahn Multi Unit Abutment non hex L05 and Screw)	BSAH2032L05	25	20		10-12
SET ZZ BASE C-J-ZZMUA-L05 4.3 (Zirkonzahn Multi Unit Abutment non hex L05 and Screw)	BSAH2033L05	25	20		10-12
SET ZZ BASE C-J-ZZMUA-L05 5.0 (Zirkonzahn Multi Unit Abutment non hex L05 and Screw)	BSAH2034L05	25	30		10-12
SET ZZ BASE C-J-ZZMUA-L05 6.0 (Zirkonzahn Multi Unit Abutment non hex L05 and Screw)	BSAH2035L05	25	30		10-12
SET ZZ BASE C-J-ZZMUA-L15 3.3 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH2031L15	25	20		10-12
SET ZZ BASE C-J-ZZMUA-L15 3.8 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH2032L15	25	20		10-12
SET ZZ BASE C-J-ZZMUA-L15 4.3 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH2033L15	25	20		10-12
SET ZZ BASE C-J-ZZMUA-L15 5.0 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH2034L15	25	30		10-12
SET ZZ BASE C-J-ZZMUA-L15 6.0 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH2035L15	25	30		10-12
SET ZZ BASE C-J-ZZMUA-L25 3.3 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH2031L25	25	20		10-12
SET ZZ BASE C-J-ZZMUA-L25 3.8 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH2032L25	25	20		10-12
SET ZZ BASE C-J-ZZMUA-L25 4.3 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH2033L25	25	20		10-12
SET ZZ BASE C-J-ZZMUA-L25 5.0 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH2034L25	25	30		10-12
SET ZZ BASE C-J-ZZMUA-L25 6.0 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH2035L25	25	30		10-12
SET ZZ BASE C-J-ZZMUA-L35 3.3 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH2031L35	25	20		10-12
SET ZZ BASE C-J-ZZMUA-L35 3.8 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH2032L35	25	20		10-12
SET ZZ BASE C-J-ZZMUA-L35 4.3 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH2033L35	25	20		10-12
SET ZZ BASE C-J-ZZMUA-L35 5.0 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH2034L35	25	30		10-12
SET ZZ BASE C-J-ZZMUA-L35 6.0 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH2035L35	25	30		10-12
SET ZZ BASE C-J-ZZMUA-L45 3.3 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH2031L45	25	20		10-12
SET ZZ BASE C-J-ZZMUA-L45 3.8 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH2032L45	25	20		10-12
SET ZZ BASE C-J-ZZMUA-L45 4.3 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH2033L45	25	20		10-12
SET ZZ BASE C-J-ZZMUA-L45 5.0 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH2034L45	25	30		10-12
SET ZZ BASE C-J-ZZMUA-L45 6.0 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH2035L45	25	30		10-12
ZZ BASE C-J-ASCSC 3.3 (Angled screw-channel Screw)	BSKE2001	15			10-12
ZZ BASE C-J-ASCSC 3.8 (Angled screw-channel Screw)	BSKE2002	15			10-12
ZZ BASE C-J-ASCSC 4.3 (Angled screw-channel Screw)	BSKE2003	15			10-12
ZZ BASE C-J-BSC 3.3 (Black Screw)	BSIE2001	20			10-12

CAMLOG® J-LINE/K-LINE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE C-J-BSC 3.8 (Black Screw)	BSIE2002	20			10-12
ZZ BASE C-J-BSC 4.3 (Black Screw)	BSIE2003	20			10-12
ZZ BASE C-J-BSC 5.0 (Black Screw)	BSIE2004	20			10-12
ZZ BASE C-J-BSC 6.0 (Black Screw)	BSIE2005	20			10-12
ZZ BASE C-J-BZSC 3.3 (Black Zirconia screw)	BSIE2031	20			10-12
ZZ BASE C-J-BZSC 3.8 (Black Zirconia screw)	BSIE2032	20			10-12
ZZ BASE C-J-BZSC 4.3 (Black Zirconia screw)	BSIE2033	20			10-12
ZZ BASE C-J-BZSC 5.0 (Black Zirconia screw)	BSIE2034	20			10-12
ZZ BASE C-J-BZSC 6.0 (Black Zirconia screw)	BSIE2035	30			10-12
ZZ BASE C-J-ICSC 3.3 (Impression coping Screw)	BSFE2001	15			10-12
ZZ BASE C-J-ICSC 3.8 (Impression coping Screw)	BSFE2002	15			10-12
ZZ BASE C-J-ICSC 4.3 (Impression coping Screw)	BSFE2003	15			10-12
ZZ BASE C-J-ICSC 5.0 (Impression coping Screw)	BSFE2004	15			10-12
ZZ BASE C-J-ICSC 6.0 (Impression coping Screw)	BSFE2005	15			10-12
ZZ BASE C-J-LBSC 3.3 (Lab Screw)	BSHE2001	20			10-12
ZZ BASE C-J-LBSC 3.8 (Lab Screw)	BSHE2002	20			10-12
ZZ BASE C-J-LBSC 4.3 (Lab Screw)	BSHE2003	20			10-12
ZZ BASE C-J-LBSC 5.0 (Lab Screw)	BSHE2004	20			10-12
ZZ BASE C-J-LBSC 6.0 (Lab Screw)	BSHE2005	20			10-12
ZZ BASE C-J-LBZSC 3.3 (Lab Zirconia screw)	BSHE2031	20			10-12
ZZ BASE C-J-LBZSC 3.8 (Lab Zirconia screw)	BSHE2032	20			10-12
ZZ BASE C-J-LBZSC 4.3 (Lab Zirconia screw)	BSHE2033	20			10-12
ZZ BASE C-J-LBZSC 5.0 (Lab Zirconia screw)	BSHE2034	20			10-12
ZZ BASE C-J-LBZSC 6.0 (Lab Zirconia screw)	BSHE2035	20			10-12
ZZ BASE C-J-MI 3.3 (Model implant)	MIAG2001				10-12
ZZ BASE C-J-MI 3.8 (Model implant)	MIAG2002				10-12
ZZ BASE C-J-MI 4.3 (Model implant)	MIAG2003				10-12
ZZ BASE C-J-MI 5.0 (Model implant)	MIAG2004				10-12
ZZ BASE C-J-MI 6.0 (Model implant)	MIAG2005				10-12

CAMLOG® J-LINE/K-LINE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE C-J-SC 3.3 (Screw)	BSAE2001	20			10-12
ZZ BASE C-J-SC 3.8 (Screw)	BSAE2002	20			10-12
ZZ BASE C-J-SC 4.3 (Screw)	BSAE2003	20			10-12
ZZ BASE C-J-SC 5.0 (Screw)	BSAE2004	20			10-12
ZZ BASE C-J-SC 6.0 (Screw)	BSAE2005	20			10-12
ZZ BASE C-J-SC-G 3.3 (Screw Gold)	BSCE2001	20			10-12
ZZ BASE C-J-SC-G 4.3 (Screw Gold)	BSCE2003	20			10-12
ZZ BASE C-J-SC-G 5.0 (Screw Gold)	BSCE2004	20			10-12
ZZ BASE C-J-SC-G 6.0 (Screw Gold)	BSCE2005	20			10-12
ZZ BASE C-J-WMSM-L10 3.3 (White Metal Scanmarker L10 with Screw)	ZBED2031L10	15			10-12
ZZ BASE C-J-WMSM-L10 3.8 (White Metal Scanmarker L10 with Screw)	ZBED2032L10	15			10-12
ZZ BASE C-J-WMSM-L10 4.3 (White Metal Scanmarker L10 with Screw)	ZBED2033L10	15			10-12
ZZ BASE C-J-WMSM-L10 5.0 (White Metal Scanmarker L10 with Screw)	ZBED2034L10	15			10-12
ZZ BASE C-J-WMSM-L10 6.0 (White Metal Scanmarker L10 with Screw)	ZBED2035L10	15			10-12
ZZ BASE C-J-WMSM-L12 3.3 (White Metal Scanmarker L12 with Screw)	ZBED2031L12	15			10-12
ZZ BASE C-J-WMSM-L12 3.8 (White Metal Scanmarker L12 with Screw)	ZBED2032L12	15			10-12
ZZ BASE C-J-WMSM-L12 4.3 (White Metal Scanmarker L12 with Screw)	ZBED2033L12	15			10-12
ZZ BASE C-J-WMSM-L12 5.0 (White Metal Scanmarker L12 with Screw)	ZBED2034L12	15			10-12
ZZ BASE C-J-WMSM-L12 6.0 (White Metal Scanmarker L12 with Screw)	ZBED2035L12	15			10-12
ZZ BASE C-J-WMSM-L15 3.3 (White Metal Scanmarker L15 with Screw)	ZBED2031L15	15			10-12
ZZ BASE C-J-WMSM-L15 3.8 (White Metal Scanmarker L15 with Screw)	ZBED2032L15	15			10-12
ZZ BASE C-J-WMSM-L15 4.3 (White Metal Scanmarker L15 with Screw)	ZBED2033L15	15			10-12
ZZ BASE C-J-WMSM-L15 5.0 (White Metal Scanmarker L15 with Screw)	ZBED2034L15	15			10-12
ZZ BASE C-J-WMSM-L15 6.0 (White Metal Scanmarker L15 with Screw)	ZBED2035L15	15			10-12
ZZ BASE C-J-ZSC 3.3 (Zirconia screw)	BSAE2031	20			10-12
ZZ BASE C-J-ZSC 3.8 (Zirconia screw)	BSAE2032	20			10-12
ZZ BASE C-J-ZSC 4.3 (Zirconia screw)	BSAE2033	20			10-12
ZZ BASE C-J-ZSC 5.0 (Zirconia screw)	BSAE2034	20			10-12
ZZ BASE C-J-ZSC 6.0 (Zirconia screw)	BSAE2035	20			10-12

CAMLOG® J-LINE/K-LINE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE C-J-ZSC-G 3.3 (Zirconia Screw Gold)	BSCE2031	20			10-12
ZZ BASE C-J-ZSC-G 3.8 (Zirconia Screw Gold)	BSCE2032	20			10-12
ZZ BASE C-J-ZSC-G 4.3 (Zirconia Screw Gold)	BSCE2033	20			10-12
ZZ BASE C-J-ZSC-G 5.0 (Zirconia Screw Gold)	BSCE2034	20			10-12
ZZ BASE C-J-ZSC-G 6.0 (Zirconia Screw Gold)	BSCE2035	20			10-12
ZZ BASE C-J-ZZMUASC 3.3 (Zirkonzahn Multi Unit Abutment Screw)	BSEE2001	20			10-12
ZZ BASE C-J-ZZMUASC 3.8 (Zirkonzahn Multi Unit Abutment Screw)	BSEE2002	20			10-12
ZZ BASE C-J-ZZMUASC 4.3 (Zirkonzahn Multi Unit Abutment Screw)	BSEE2003	20			10-12
ZZ BASE C-J-ZZMUASC 5.0 (Zirkonzahn Multi Unit Abutment Screw)	BSEE2004	30			10-12

CAMLOG® VARIO SR

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE C-VSR-CCTB 3.8/4.3 (Conical cemented titan base non hex and screw)	BSAD5431	20			10-12
SET ZZ BASE C-VSR-CCTB 5.0/6.0 (Conical cemented titan base non hex and screw)	BSAD5432	20			10-12
SET ZZ BASE C-VSR-SA 3.8/4.3 (Scananalog and screw)	ZBDE5431	15			10-12
SET ZZ BASE C-VSR-SA 5.0/6.0 (Scananalog and screw)	ZBDE5432	15			10-12
SET ZZ BASE C-VSR-SM 3.8/4.3 (Scanmarker and screw)	ZBAD5431	20			10-12
SET ZZ BASE C-VSR-SM 5.0/6.0 (Scanmarker and screw)	ZBAD5432	20			10-12
ZZ BASE C-VSR-SC 3.8/4.3 (Screw)	BSAE5401	20			10-12
ZZ BASE C-VSR-SC 5.0/6.0 (Screw)	BSAE5402	20			10-12
ZZ BASE C-VSR-ZSC 3.8/4.3 (Zirconia screw)	BSAE5431	20			10-12
ZZ BASE C-VSR-ZSC 5.0/6.0 (Zirconia screw)	BSAE5432	20			10-12

CHAMPIONS® IMPLANTS (R)EVOLUTION

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
SET ZZ BASE CI-RE-ABCH 3.5/4.0/4.5/5.5 (Raw-Abutment® customizable hex and screw)	BSAF10781	30			10-12
SET ZZ BASE CI-RE-CCTB 3.5/4.0/4.5/5.5 (Conical cemented titan base non hex and screw)	BSAD10731	30			10-12
SET ZZ BASE CI-RE-CCTB MUA (Conical cemented titan base non hex and screw)	BSAD10732	20			10-12
SET ZZ BASE CI-RE-CCTB-K85 3.5/4.0/4.5/5.5 (Conical cemented titan base non hex K85 and screw)	BSAD10731K85	30			10-12
SET ZZ BASE CI-RE-CCTB-K85 MUA (Conical cemented titan base non hex K85 and screw)	BSAD10732K85	20			10-12
SET ZZ BASE CI-RE-ICOT MUA (Impression coping open tray non hex and Screw and Driver)	BSAQ10762	15			10-12
SET ZZ BASE CI-RE-ICOTH 3.5/4.0/4.5/5.5 (Impression coping open tray hex and Screw and Driver)	BSAQ10741	15			10-12
SET ZZ BASE CI-RE-NTB 3.5/4.0/4.5/5.5 (Narrow titan base non hex and screw)	BSAP10731	30			10-12
SET ZZ BASE CI-RE-NTBH 3.5/4.0/4.5/5.5 (Narrow titan base hex and screw)	BSAP10781	30			10-12
SET ZZ BASE CI-RE-PCTBH 3.5/4.0/4.5/5.5 (Parallel cemented titan base hex and screw)	BSAB10781	30			10-12
SET ZZ BASE CI-RE-PCTBH-K85 3.5/4.0/4.5/5.5 (Parallel cemented titan base hex K85 and screw)	BSAB10781K85	30			10-12
SET ZZ BASE CI-RE-SA 3.5/4.0/4.5/5.5 (Scanalog and screw)	ZBDE10731	15			10-12
SET ZZ BASE CI-RE-SA MUA (Scanalog and screw)	ZBDE10732	15			10-12
SET ZZ BASE CI-RE-SM 3.5/4.0/4.5/5.5 (Scanmarker and screw)	ZBAD10731	30			10-12
SET ZZ BASE CI-RE-SM MUA (Scanmarker and screw)	ZBAD10732	20			10-12
ZZ BASE CI-RE-ICSC 3.5/4.0/4.5/5.5 (Impression coping Screw)	BSFE10701	15			10-12
ZZ BASE CI-RE-ICSC MUA (Impression coping Screw)	BSFE10702	15			10-12
ZZ BASE CI-RE-MI 3.5/4.0/4.5/5.5 (Model implant)	MIAG10701				10-12
ZZ BASE CI-RE-MI MUA (Model implant)	MIAG10702				10-12
ZZ BASE CI-RE-SC 3.5/4.0/4.5/5.5 (Screw)	BSAE10701	30			10-12
ZZ BASE CI-RE-SC MUA (Screw)	BSAE10702	20			10-12
ZZ BASE CI-RE-WMSM-L12 3.5/4.0/4.5/5.5 (White Metal Scanmarker L12 with Screw)	ZBED10731L12	15			10-12
ZZ BASE CI-RE-WMSM-L12 MUA (White Metal Scanmarker L12 with Screw)	ZBED10732L12	15			10-12
ZZ BASE CI-RE-ZSC 3.5/4.0/4.5/5.5 (Zirconia screw)	BSAE10731	30			10-12
ZZ BASE CI-RE-ZSC MUA (Zirconia screw)	BSAE10732	20			10-12

CONMET® HEX

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
SET ZZ BASE CO-H-CCTB HEX 2.2 (Conical cemented titan base non hex and screw)	BSAD7831	20			10-12
SET ZZ BASE CO-H-CCTB HEX 2.7 (Conical cemented titan base non hex and screw)	BSAD7832	20			10-12
SET ZZ BASE CO-H-ICOTH HEX 2.2 (Impression coping open tray hex and Screw and Driver)	BSAQ7841	15			10-12
SET ZZ BASE CO-H-ICOTH HEX 2.7 (Impression coping open tray hex and Screw and Driver)	BSAQ7842	15			10-12
SET ZZ BASE CO-H-PCTBH HEX 2.2 (Parallel cemented titan base hex and screw)	BSAB7881	20			10-12
SET ZZ BASE CO-H-PCTBH HEX 2.7 (Parallel cemented titan base hex and screw)	BSAB7882	20			10-12
SET ZZ BASE CO-H-SM HEX 2.2 (Scanmarker and screw)	ZBAD7831	20			10-12
SET ZZ BASE CO-H-SM HEX 2.7 (Scanmarker and screw)	ZBAD7832	20			10-12
ZZ BASE CO-H-ASCSC HEX 2.2 (Angled screw-channel Screw)	BSKE7801	20			10-12
ZZ BASE CO-H-ASCSC HEX 2.7 (Angled screw-channel Screw)	BSKE7802	20			10-12
ZZ BASE CO-H-ICSC HEX 2.2 (Impression coping Screw)	BSFE7801	15			10-12
ZZ BASE CO-H-ICSC HEX 2.7 (Impression coping Screw)	BSFE7802	15			10-12
ZZ BASE CO-H-SC HEX 2.2 (Screw)	BSAE7801	20			10-12
ZZ BASE CO-H-SC HEX 2.7 (Screw)	BSAE7802	20			10-12
ZZ BASE CO-H-WMSM-L12 HEX 2.2 (White Metal Scanmarker L12 with Screw)	ZBED7831L12	15			10-12
ZZ BASE CO-H-WMSM-L12 HEX 2.7 (White Metal Scanmarker L12 with Screw)	ZBED7832L12	15			10-12
ZZ BASE CO-H-ZSC HEX 2.2 (Zirconia Screw)	BSAE7831	20			10-12
ZZ BASE CO-H-ZSC HEX 2.7 (Zirconia Screw)	BSAE7832	20			10-12

CONNECT®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
SET ZZ BASE CN-MC-CCTB (Conical cemented titan base non hex and screw)	BSAD8831	20			10-12
SET ZZ BASE CN-MC-ICOTH (Impression coping open tray hex and Screw and Driver)	BSAQ8841	15			10-12
SET ZZ BASE CN-MC-PCTBH (Parallel cemented titan base hex and screw)	BSAB8881	20			10-12
SET ZZ BASE CN-MC-SA (Scananalog and screw)	ZBDE8831	15			10-12
SET ZZ BASE CN-MC-SM (Scanmarker and screw)	ZBAD8831	20			10-12
SET ZZ BASE CN-ST-CCTB (Conical cemented titan base non hex and screw)	BSAD8832	35			10-12
SET ZZ BASE CN-ST-ICOTH (Impression coping open tray hex and Screw and Driver)	BSAQ8842	15			10-12
SET ZZ BASE CN-ST-PCTBH (Parallel cemented titan base hex and screw)	BSAB8882	35			10-12
SET ZZ BASE CN-ST-SA (Scananalog and screw)	ZBDE8832	15			10-12
SET ZZ BASE CN-ST-SM (Scanmarker and screw)	ZBAD8832	35			10-12
ZZ BASE CN-MC-ICSC (Impression coping Screw)	BSFE8801	15			10-12
ZZ BASE CN-MC-MI (Model implant)	MIAG8801				10-12
ZZ BASE CN-MC-SC (Screw)	BSAE8801	20			10-12
ZZ BASE CN-MC-WMSM-L12 (White Metal Scanmarker L12 with Screw)	ZBED8831L12	15			10-12
ZZ BASE CN-MC-ZSC (Zirconia Screw)	BSAE8831	20			10-12
ZZ BASE CN-ST-ICSC (Impression coping Screw)	BSFE8802	15			10-12
ZZ BASE CN-ST-MI (Model implant)	MIAG8802				10-12
ZZ BASE CN-ST-SC (Screw)	BSAE8802	35			10-12
ZZ BASE CN-ST-WMSM-L12 (White Metal Scanmarker L12 with Screw)	ZBED8832L12	15			10-12
ZZ BASE CN-ST-ZSC (Zirconia Screw)	BSAE8832	35			10-12

COWELLMEDI INNO INTERNAL IMPLANT SYSTEM™

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
SET ZZ BASE CM-IN-ABCH 3.5-6.0 (Raw-Abutment® customizable hex and screw)	BSAF6981	30			10-12
SET ZZ BASE CM-IN-SA 3.5-6.0 (Scananalog and screw)	ZBDE6931	15			10-12
SET ZZ BASE CM-IN-SM 3.5-6.0 (Scanmarker and screw)	ZBAD6931	30			10-12
ZZ BASE CM-IN-SC 3.5-6.0 (Screw)	BSAE6901	30			10-12

CUMDENTE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE CU-AS/AS-V-ABCH 3.5/4.0 (Raw-Abutment® customizable hex and screw)	BSAF10183	25			10-12
SET ZZ BASE CU-AS/AS-V-ABCH 4.7/5.0/5.5 (Raw-Abutment® customizable hex and screw)	BSAF10184	25			10-12
SET ZZ BASE CU-AS/AS-V-CCTB 3.5/4.0 (Conical cemented titan base non hex and screw)	BSAD10133	25			10-12
SET ZZ BASE CU-AS/AS-V-CCTB 4.7/5.0/5.5 (Conical cemented titan base non hex and screw)	BSAD10134	25			10-12
SET ZZ BASE CU-AS/AS-V-CCTB-K85 3.5/4.0 (Conical cemented titan base non hex K85 and scr	BSAD10133K85	25			10-12
SET ZZ BASE CU-AS/AS-V-CCTB-K85 4.7/5.0/5.5 (Conical cemented titan base non hex K85 and	BSAD10134K85	25			10-12
SET ZZ BASE CU-AS/AS-V-ICOTH 3.5/4.0 (Impression coping open tray hex and Screw and Drive	BSAQ10143	15			10-12
SET ZZ BASE CU-AS/AS-V-ICOTH 4.7/5.0/5.5 (Impression coping open tray hex and Screw and C	BSAQ10144	15			10-12
SET ZZ BASE CU-AS/AS-V-NTB 3.5/4.0 (Narrow titan base non hex and screw)	BSAP10133	25			10-12
SET ZZ BASE CU-AS/AS-V-NTBH 3.5/4.0 (Narrow titan base hex and screw)	BSAP10183	25			10-12
SET ZZ BASE CU-AS/AS-V-PCTBH 3.5/4.0 (Parallel cemented titan base hex and screw)	BSAB10183	25			10-12
SET ZZ BASE CU-AS/AS-V-PCTBH 4.7/5.0/5.5 (Parallel cemented titan base hex and screw)	BSAB10184	25			10-12
SET ZZ BASE CU-AS/AS-V-PCTBH-K85 3.5/4.0 (Parallel cemented titan base hex K85 and screw	BSAB10183K85	25			10-12
SET ZZ BASE CU-AS/AS-V-PCTBH-K85 4.7/5.0/5.5 (Parallel cemented titan base hex K85 and sc	BSAB10184K85	25			10-12
SET ZZ BASE CU-AS/AS-V-SA 3.5/4.0 (Scananalog and screw)	ZBDE10133	15			10-12
SET ZZ BASE CU-AS/AS-V-SA 4.7/5.0/5.5 (Scananalog and screw)	ZBDE10134	15			10-12
SET ZZ BASE CU-AS/AS-V-SM 3.5/4.0 (Scanmarker and screw)	ZBAD10133	25			10-12
SET ZZ BASE CU-AS/AS-V-SM 4.7/5.0/5.5 (Scanmarker and screw)	ZBAD10134	25			10-12
SET ZZ BASE CU-AS/AS-V-WSM 3.5/4.0 (White Scanmarker and screw)	ZBCD10133	25			10-12
SET ZZ BASE CU-AS-ABCH 3.0 (Raw-Abutment® customizable hex and screw)	BSAF10182	25			10-12
SET ZZ BASE CU-AS-CCTB 3.0 (Conical cemented titan base non hex and screw)	BSAD10132	25			10-12
SET ZZ BASE CU-AS-CCTB-K85 3.0 (Conical cemented titan base non hex K85 and screw)	BSAD10132K85	25			10-12
SET ZZ BASE CU-AS-ICOTH 3.0 (Impression coping open tray hex and Screw and Driver)	BSAQ10142	15			10-12
SET ZZ BASE CU-AS-NTB 3.0 (Narrow titan base non hex and screw)	BSAP10132	25			10-12
SET ZZ BASE CU-AS-NTBH 3.0 (Narrow titan base hex and screw)	BSAP10182	25			10-12
SET ZZ BASE CU-AS-PCTBH 3.0 (Parallel cemented titan base hex and screw)	BSAB10182	25			10-12
SET ZZ BASE CU-AS-PCTBH-K85 3.0 (Parallel cemented titan base hex K85 and screw)	BSAB10182K85	25			10-12
SET ZZ BASE CU-AS-SA 3.0 (Scananalog and screw)	ZBDE10132	15			10-12
SET ZZ BASE CU-AS-SM 3.0 (Scanmarker and screw)	ZBAD10132	25			10-12

CUMDENTE					
Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE CU-MVX-ABCH (Raw-Abutment® customizable hex and screw)	BSAF10191	25			10-12
SET ZZ BASE CU-MVX-CCTB (Conical cemented titan base non hex and screw)	BSAD10131	25			10-12
SET ZZ BASE CU-MVX-CCTB-K85 (Conical cemented titan base non hex K85 and screw)	BSAD10131K85	25			10-12
SET ZZ BASE CU-MVX-ICOTH (Impression coping open tray hex and Screw and Driver)	BSAQ10141	15			10-12
SET ZZ BASE CU-MVX-NTB (Narrow titan base non hex and screw)	BSAP10131	25			10-12
SET ZZ BASE CU-MVX-NTBH (Narrow titan base hex and screw)	BSAP10181	25			10-12
SET ZZ BASE CU-MVX-PCTBH (Parallel cemented titan base hex and screw)	BSAB10181	25			10-12
SET ZZ BASE CU-MVX-PCTBH-K85 (Parallel cemented titan base hex K85 and screw)	BSAB10181K85	25			10-12
SET ZZ BASE CU-MVX-SA (Scananalog and screw)	ZBDE10131	15			10-12
SET ZZ BASE CU-MVX-SM (Scanmarker and screw)	ZBAD10131	25			10-12
SET ZZ BASE CU-MVX-WSM (White Scanmarker and screw)	ZBCD10131	25			10-12
ZZ BASE CU-AS/AS-V-ASCSC 3.5/4.0 (Angled screw-channel Screw)	BSKE10103	20			10-12
ZZ BASE CU-AS/AS-V-ICSC 3.5/4.0 (Impression coping Screw)	BSFE10103	15			10-12
ZZ BASE CU-AS/AS-V-ICSC 4.7/5.0/5.5 (Impression coping Screw)	BSFE10104	15			10-12
ZZ BASE CU-AS/AS-V-MI 3.5/4.0 (Model implant)	MIAG10103				10-12
ZZ BASE CU-AS/AS-V-MI 4.7/5.0/5.5 (Model implant)	MIAG10104				10-12
ZZ BASE CU-AS/AS-V-SC 3.5/4.0 (Screw)	BSAE10103	25			10-12
ZZ BASE CU-AS/AS-V-SC 4.7/5.0/5.5 (Screw)	BSAE10104	25			10-12
ZZ BASE CU-AS/AS-V-WMSM-L12 3.5/4.0 (White Metal Scanmarker L12 with Screw)	ZBED10133L12	15			10-12
ZZ BASE CU-AS/AS-V-WMSM-L12 4.7/5.0/5.5 (White Metal Scanmarker L12 with Screw)	ZBED10134L12	15			10-12
ZZ BASE CU-AS/AS-V-ZSC 3.5/4.0 (Zirconia screw)	BSAE10133	25			10-12
ZZ BASE CU-AS/AS-V-ZSC 4.7/5.0/5.5 (Zirconia screw)	BSAE10134	25			10-12
ZZ BASE CU-AS-ICSC 3.0 (Impression coping Screw)	BSFE10102	15			10-12
ZZ BASE CU-AS-MI 3.0 (Model implant)	MIAG10102				10-12
ZZ BASE CU-AS-SC 3.0 (Screw)	BSAE10102	25			10-12
ZZ BASE CU-AS-WMSM-L12 3.0 (White Metal Scanmarker L12 with Screw)	ZBED10132L12	15			10-12
ZZ BASE CU-AS-ZSC 3.0 (Zirconia screw)	BSAE10132	25			10-12
ZZ BASE CU-MVX-ICSC (Impression coping Screw)	BSFE10101	15			10-12
ZZ BASE CU-MVX-MI (Model implant)	MIAG10101				10-12

CUMDENTE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
ZZ BASE CU-MVX-SC (Screw)	BSAE10101	25			10-12
ZZ BASE CU-MVX-WMSM-L12 (White Metal Scanmarker L12 with Screw)	ZBED10131L12	15			10-12
ZZ BASE CU-MVX-ZSC (Zirconia screw)	BSAE10131	25			10-12

DENTAL RATIO® OKTAGON® BONE LEVEL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE DO-B-ABCH NC (Raw-Abutment® customizable hex and screw)	BSAF13781	35			10-12
SET ZZ BASE DO-B-ABCH RC (Raw-Abutment® customizable hex and screw)	BSAF13782	35			10-12
SET ZZ BASE DO-B-ABCH-D14 NC (Raw-Abutment® customizable hex D14 and screw)	BSAF13781D14	35			10-12
SET ZZ BASE DO-B-ABCH-D14 RC (Raw-Abutment® customizable hex D14 and screw)	BSAF13782D14	35			10-12
SET ZZ BASE DO-B-CCTB NC (Conical cemented titan base non hex and screw)	BSAD13731	35			10-12
SET ZZ BASE DO-B-CCTB RC (Conical cemented titan base non hex and screw)	BSAD13732	35			10-12
SET ZZ BASE DO-B-CCTB-K85 NC (Conical cemented titan base non hex K85 and screw)	BSAD13731K85	35			10-12
SET ZZ BASE DO-B-CCTB-K85 RC (Conical cemented titan base non hex K85 and screw)	BSAD13732K85	35			10-12
SET ZZ BASE DO-B-CCTB-L10 NC (Conical cemented titan base non hex L10 and screw)	BSAD13731L10	35			10-12
SET ZZ BASE DO-B-CCTB-L10 RC (Conical cemented titan base non hex L10 and screw)	BSAD13732L10	35			10-12
SET ZZ BASE DO-B-CCTB-L20 NC (Conical cemented titan base non hex L20 and screw)	BSAD13731L20	35			10-12
SET ZZ BASE DO-B-CCTB-L20 RC (Conical cemented titan base non hex L20 and screw)	BSAD13732L20	35			10-12
SET ZZ BASE DO-B-CCTB-L30 NC (Conical cemented titan base non hex L30 and screw)	BSAD13731L30	35			10-12
SET ZZ BASE DO-B-CCTB-L30 RC (Conical cemented titan base non hex L30 and screw)	BSAD13732L30	35			10-12
SET ZZ BASE DO-B-ICOTH NC (Impression coping open tray hex and Screw and Driver)	BSAQ13741	15			10-12
SET ZZ BASE DO-B-ICOTH RC (Impression coping open tray hex and Screw and Driver)	BSAQ13742	15			10-12
SET ZZ BASE DO-B-NTB NC (Narrow titan base non hex and screw)	BSAP13731	35			10-12
SET ZZ BASE DO-B-NTB RC (Narrow titan base non hex and screw)	BSAP13732	35			10-12
SET ZZ BASE DO-B-NTBH NC (Narrow titan base hex and screw)	BSAP13781	35			10-12
SET ZZ BASE DO-B-NTBH RC (Narrow titan base hex and screw)	BSAP13782	35			10-12
SET ZZ BASE DO-B-PCTBH NC (Parallel cemented titan base hex and screw)	BSAB13791	35			10-12
SET ZZ BASE DO-B-PCTBH RC (Parallel cemented titan base hex and screw)	BSAB13792	35			10-12
SET ZZ BASE DO-B-PCTBH-K85 NC (Parallel cemented titan base hex K85 and screw)	BSAB13781K85	35			10-12
SET ZZ BASE DO-B-PCTBH-K85 RC (Parallel cemented titan base hex K85 and screw)	BSAB13782K85	35			10-12
SET ZZ BASE DO-B-PCTBH-L20 NC (Parallel cemented titan base hex L20 and screw)	BSAB13791L20	35			10-12
SET ZZ BASE DO-B-PCTBH-L20 RC (Parallel cemented titan base hex L20 and screw)	BSAB13792L20	35			10-12
SET ZZ BASE DO-B-PCTBH-L30 NC (Parallel cemented titan base hex L30 and screw)	BSAB13791L30	35			10-12
SET ZZ BASE DO-B-PCTBH-L30 RC (Parallel cemented titan base hex L30 and screw)	BSAB13792L30	35			10-12
SET ZZ BASE DO-B-PCTBH-L40 NC (Parallel cemented titan base hex L40 and screw)	BSAB13791L40	35			10-12

DENTAL RATIO® OKTAGON® BONE LEVEL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE DO-B-PCTBH-L40 RC (Parallel cemented titan base hex L40 and screw)	BSAB13792L40	35			10-12
SET ZZ BASE DO-B-SA NC (Scananalogue and screw)	ZBDE13731	15			10-12
SET ZZ BASE DO-B-SA RC (Scananalogue and screw)	ZBDE13732	15			10-12
SET ZZ BASE DO-B-SM NC (Scanmarker and screw)	ZBAD13731	35			10-12
SET ZZ BASE DO-B-SM RC (Scanmarker and screw)	ZBAD13732	35			10-12
SET ZZ BASE DO-B-WSM NC (White Scanmarker and screw)	ZBCD13731	35			10-12
SET ZZ BASE DO-B-WSM RC (White Scanmarker and screw)	ZBCD13732	35			10-12
SET ZZ BASE DO-B-ZZMUA17 NC Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and implai	BSAI13731	20			10-12
SET ZZ BASE DO-B-ZZMUA17 NC Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and implai	BSAI13741	20			10-12
SET ZZ BASE DO-B-ZZMUA17 RC Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and implai	BSAI13732	20			10-12
SET ZZ BASE DO-B-ZZMUA17 RC Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and implai	BSAI13742	20			10-12
SET ZZ BASE DO-B-ZZMUA17-L10 NC Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10 a	BSAI13731L10	20			10-12
SET ZZ BASE DO-B-ZZMUA17-L10 NC Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 L10 a	BSAI13741L10	20			10-12
SET ZZ BASE DO-B-ZZMUA17-L10 RC Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10 a	BSAI13732L10	20			10-12
SET ZZ BASE DO-B-ZZMUA17-L10 RC Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 L10 a	BSAI13742L10	20			10-12
SET ZZ BASE DO-B-ZZMUA-L15 NC (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH13731L15	25	30		10-12
SET ZZ BASE DO-B-ZZMUA-L15 RC (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH13732L15	25	30		10-12
SET ZZ BASE DO-B-ZZMUA-L25 NC (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH13731L25	25	30		10-12
SET ZZ BASE DO-B-ZZMUA-L25 RC (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH13732L25	25	30		10-12
SET ZZ BASE DO-B-ZZMUA-L35 NC (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH13731L35	25	30		10-12
SET ZZ BASE DO-B-ZZMUA-L35 RC (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH13732L35	25	30		10-12
SET ZZ BASE DO-B-ZZMUA-L45 NC (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH13731L45	25	30		10-12
SET ZZ BASE DO-B-ZZMUA-L45 RC (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH13732L45	25	30		10-12
ZZ BASE DO-B-BSC NC (Black Screw)	BSIE13701	35			10-12
ZZ BASE DO-B-BSC RC (Black Screw)	BSIE13702	35			10-12
ZZ BASE DO-B-BZSC NC (Black Zirconia screw)	BSIE13731	35			10-12
ZZ BASE DO-B-BZSC RC (Black Zirconia screw)	BSIE13732	35			10-12
ZZ BASE DO-B-ICSC NC (Impression coping Screw)	BSFE13701	15			10-12
ZZ BASE DO-B-ICSC RC (Impression coping Screw)	BSFE13702	15			10-12

DENTAL RATIO® OKTAGON® BONE LEVEL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE DO-B-LBSC NC (Lab Screw)	BSHE13701	35			10-12
ZZ BASE DO-B-LBSC RC (Lab Screw)	BSHE13702	35			10-12
ZZ BASE DO-B-LBZSC NC (Lab Zirconia screw)	BSHE13731	35			10-12
ZZ BASE DO-B-LBZSC RC (Lab Zirconia screw)	BSHE13732	35			10-12
ZZ BASE DO-B-MI NC (Model implant)	MIAG13701				10-12
ZZ BASE DO-B-MI RC (Model implant)	MIAG13702				10-12
ZZ BASE DO-B-SC NC (Screw)	BSAE13701	35			10-12
ZZ BASE DO-B-SC RC (Screw)	BSAE13702	35			10-12
ZZ BASE DO-B-SC-G NC (Screw Gold)	BSCE13701	35			10-12
ZZ BASE DO-B-SC-G RC (Screw Gold)	BSCE13702	35			10-12
ZZ BASE DO-B-WMSM-L12 NC (White Metal Scanmarker L12 with Screw)	ZBED13731L12	15			10-12
ZZ BASE DO-B-WMSM-L12 RC (White Metal Scanmarker L12 with Screw)	ZBED13732L12	15			10-12
ZZ BASE DO-B-ZSC NC (Zirconia screw)	BSAE13731	35			10-12
ZZ BASE DO-B-ZSC RC (Zirconia screw)	BSAE13732	35			10-12
ZZ BASE DO-B-ZSC-G NC (Zirconia Screw Gold)	BSCE13731	35			10-12
ZZ BASE DO-B-ZSC-G RC (Zirconia Screw Gold)	BSCE13732	35			10-12
ZZ BASE DO-B-ZZMUASC NC (Zirkonzahn Multi Unit Abutment Screw)	BSEE13701	20			10-12
ZZ BASE DO-B-ZZMUASC RC (Zirkonzahn Multi Unit Abutment Screw)	BSEE13702	20			10-12

DENTAL RATIO® OKTAGON® MULTI UNITS ABUTMENT - NOBEL SCREWDRI

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE DO-M-CCTB NC/RC (Conical cemented titan base non hex and screw)	BSAD13831	15			10-12
SET ZZ BASE DO-M-CCTB NC/RC (Conical cemented titan base non hex and screw)	BSAD13841	15			10-12
SET ZZ BASE DO-M-CCTB-K85 NC/RC (Conical cemented titan base non hex K85 and screw)	BSAD13831K85	15			10-12
SET ZZ BASE DO-M-ICOT NC/RC (Impression coping open tray non hex and Screw and Driver)	BSAQ13861	15			10-12
SET ZZ BASE DO-M-SA NC/RC (Scanalog and screw)	ZBDE13831	15			10-12
SET ZZ BASE DO-M-SM NC/RC (Scanmarker and screw)	ZBAD13831	15			10-12
ZZ BASE DO-M-BSC NC/RC (Black Screw)	BSIE13801	15			10-12
ZZ BASE DO-M-BZSC NC/RC (Black Zirconia screw)	BSIE13831	15			10-12
ZZ BASE DO-M-ICSC NC/RC (Impression coping Screw)	BSFE13801	15			10-12
ZZ BASE DO-M-LBSC NC/RC (Lab Screw)	BSHE13801	15			10-12
ZZ BASE DO-M-LBZSC NC/RC (Lab Zirconia screw)	BSHE13831	15			10-12
ZZ BASE DO-M-MI NC/RC (Model implant)	MIAG13801				10-12
ZZ BASE DO-M-SC NC/RC (Screw)	BSAE13801	15			10-12
ZZ BASE DO-M-SC-G NC/RC (Screw Gold)	BSCE13801	15			10-12
ZZ BASE DO-M-WMSM-L12 NC/RC (White Metal Scanmarker L12 with Screw)	ZBED13831L12	15			10-12
ZZ BASE DO-M-ZSC NC/RC (Zirconia screw)	BSAE13831	15			10-12

DENTAL RATIO® OKTAGON® TISSUE LEVEL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE DO-T-ABCH RP (Raw-Abutment® customizable hex and screw)	BSAF13981	35			10-12
SET ZZ BASE DO-T-ABCH WP (Raw-Abutment® customizable hex and screw)	BSAF13982	35			10-12
SET ZZ BASE DO-T-ABCH-D14 RP (Raw-Abutment® customizable hex D14 and screw)	BSAF13981D14	35			10-12
SET ZZ BASE DO-T-ABCH-D14 WP (Raw-Abutment® customizable hex D14 and screw)	BSAF13982D14	35			10-12
SET ZZ BASE DO-T-CCTB RP (Conical cemented titan base non hex and screw)	BSAD13931	35			10-12
SET ZZ BASE DO-T-CCTB WP (Conical cemented titan base non hex and screw)	BSAD13932	35			10-12
SET ZZ BASE DO-T-CCTB-G WP (Conical cemented titan base non hex and screw Gold)	BSCD13932	35			10-12
SET ZZ BASE DO-T-CCTB-K85 RP (Conical cemented titan base non hex K85 and screw)	BSAD13931K85	35			10-12
SET ZZ BASE DO-T-CCTB-K85 WP (Conical cemented titan base non hex K85 and screw)	BSAD13932K85	35			10-12
SET ZZ BASE DO-T-CCTB-L15 RP (Conical cemented titan base non hex L15 and screw)	BSAD13931L15	35			10-12
SET ZZ BASE DO-T-CCTB-L15 WP (Conical cemented titan base non hex L15 and screw)	BSAD13932L15	35			10-12
SET ZZ BASE DO-T-PCTBH RP (Parallel cemented titan base hex and screw)	BSAB13971	35			10-12
SET ZZ BASE DO-T-PCTBH RP (Parallel cemented titan base hex and screw)	BSAB13981	35			10-12
SET ZZ BASE DO-T-PCTBH WP (Parallel cemented titan base hex and screw)	BSAB13982	35			10-12
SET ZZ BASE DO-T-PCTBH-K85 RP (Parallel cemented titan base hex K85 and screw)	BSAB13981K85	35			10-12
SET ZZ BASE DO-T-PCTBH-K85 WP (Parallel cemented titan base hex K85 and screw)	BSAB13982K85	35			10-12
SET ZZ BASE DO-T-PCTBH-L15 RP (Parallel cemented titan base hex L15 and screw)	BSAB13991L15	35			10-12
SET ZZ BASE DO-T-PCTBH-L15 WP (Parallel cemented titan base hex L15 and screw)	BSAB13992L15	35			10-12
SET ZZ BASE DO-T-SA RP (Scanalog and screw)	ZBDE13931	15			10-12
SET ZZ BASE DO-T-SA WP (Scanalog and screw)	ZBDE13932	15			10-12
SET ZZ BASE DO-T-SM RP (Scanmarker and screw)	ZBAD13931	35			10-12
SET ZZ BASE DO-T-SM WP (Scanmarker and screw)	ZBAD13932	35			10-12
SET ZZ BASE DO-T-WSM RP (White Scanmarker and screw)	ZBCD13931	35			10-12
SET ZZ BASE DO-T-WSM WP (White Scanmarker and screw)	ZBCD13932	35			10-12
SET ZZ BASE DO-T-ZZMUA17 RP Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and implar	BSAI13931	30			10-12
SET ZZ BASE DO-T-ZZMUA17 RP Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and implar	BSAI13941	30			10-12
SET ZZ BASE DO-T-ZZMUA-L05 WP (Zirkonzahn Multi Unit Abutment non hex L05 and Screw)	BSAH13932L05	25	35		10-12
SET ZZ BASE DO-T-ZZMUA-L15 RP (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH13931L15	25	35		10-12
SET ZZ BASE DO-T-ZZMUA-L15 WP (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH13932L15	25	35		10-12

DENTAL RATIO® OKTAGON® TISSUE LEVEL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE DO-T-ZZMUA-L25 RP (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH13931L25	25	35		10-12
SET ZZ BASE DO-T-ZZMUA-L35 RP (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH13931L35	25	35		10-12
SET ZZ BASE DO-T-ZZMUA-L45 RP (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH13931L45	25	35		10-12
ZZ BASE DO-T-BLZSC RP (Black Long Zirconia Screw)	BSIE13981	35			10-12
ZZ BASE DO-T-BLZSC WP (Black Long Zirconia Screw)	BSIE13982	35			10-12
ZZ BASE DO-T-BSC RP (Black Screw)	BSIE13901	35			10-12
ZZ BASE DO-T-BSC WP (Black Screw)	BSIE13902	35			10-12
ZZ BASE DO-T-BZSC RP (Black Zirconia screw)	BSIE13931	35			10-12
ZZ BASE DO-T-BZSC WP (Black Zirconia screw)	BSIE13932	35			10-12
ZZ BASE DO-T-LBLZSC RP (Lab Long Zirconia Screw)	BSHE13981	35			10-12
ZZ BASE DO-T-LBLZSC WP (Lab Long Zirconia Screw)	BSHE13982	35			10-12
ZZ BASE DO-T-LBSC RP (Lab Screw)	BSHE13901	35			10-12
ZZ BASE DO-T-LBSC WP (Lab Screw)	BSHE13902	35			10-12
ZZ BASE DO-T-LBZSC RP (Lab Zirconia screw)	BSHE13931	35			10-12
ZZ BASE DO-T-LBZSC WP (Lab Zirconia screw)	BSHE13932	35			10-12
ZZ BASE DO-T-LZSC RP (Long zirconia screw)	BSAE13981	35			10-12
ZZ BASE DO-T-LZSC WP (Long zirconia screw)	BSAE13982	35			10-12
ZZ BASE DO-T-LZSC-G RP (Long Zirconia Screw Gold)	BSCE13981	35			10-12
ZZ BASE DO-T-LZSC-G WP (Long Zirconia Screw Gold)	BSCE13982	35			10-12
ZZ BASE DO-T-MI RP (Model implant)	MIAG13901				10-12
ZZ BASE DO-T-MI WP (Model implant)	MIAG13902				10-12
ZZ BASE DO-T-SC RP (Screw)	BSAE13901	35			10-12
ZZ BASE DO-T-SC WP (Screw)	BSAE13902	35			10-12
ZZ BASE DO-T-SC-G RP (Screw Gold)	BSCE13901	35			10-12
ZZ BASE DO-T-SC-G WP (Screw Gold)	BSCE13902	35			10-12
ZZ BASE DO-T-WMSM-L12 RP (White Metal Scanmarker L12 with Screw)	ZBED13931L12	15			10-12
ZZ BASE DO-T-ZSC RP (Zirconia screw)	BSAE13931	35			10-12
ZZ BASE DO-T-ZSC WP (Zirconia screw)	BSAE13932	35			10-12
ZZ BASE DO-T-ZSC-G RP (Zirconia Screw Gold)	BSCE13931	35			10-12

DENTAL RATIO® OKTAGON® TISSUE LEVEL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
ZZ BASE DO-T-ZSC-G WP (Zirconia Screw Gold)	BSCE13932	35			10-12
ZZ BASE DO-T-ZZMUASC RP (Zirkonzahn Multi Unit Abutment Screw)	BSEE13901	30			10-12

DENTALPOINT AG ZERAMEX® P6

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
ZZ BASE ZX-P6-SM RN (Scanmarker)	ZBAD10232				10-12
ZZ BASE ZX-P6-SM SN (Scanmarker)	ZBAD10231				10-12
ZZ BASE ZX-P6-WSM RN (White Scanmarker)	ZBCD10232				10-12
ZZ BASE ZX-P6-WSM SN (White Scanmarker)	ZBCD10231				10-12

DENTALPOINT AG ZERAMEX® XT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE ZX-XT-SA RB (Scanalog and screw)	ZBDE12931	15			10-12
SET ZZ BASE ZX-XT-SA WB (Scanalog and screw)	ZBDE12932	15			10-12
ZZ BASE ZX-XT-SM RB (Scanmarker)	ZBAD12931				10-12
ZZ BASE ZX-XT-SM WB (Scanmarker)	ZBAD12932				10-12
ZZ BASE ZX-XT-WSM RB (White Scanmarker)	ZBCD12931				10-12

DENTALTECH IMPLASSIC®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE DI-IL-ABCH HX2 (Raw-Abutment® customizable hex and screw)	BSAF11891	20			10-12
SET ZZ BASE DI-IL-ABCH TR (Raw-Abutment® customizable hex and screw)	BSAF11781	20			10-12
SET ZZ BASE DI-IL-ABCH TRW (Raw-Abutment® customizable hex and screw)	BSAF11782	20			10-12
SET ZZ BASE DI-IL-CCTB HX2 (Conical cemented titan base non hex and screw)	BSAD11831	20			10-12
SET ZZ BASE DI-IL-CCTB Pilastro PDM/PPM (Conical cemented titan base non hex and screw)	BSAD11832	20			10-12
SET ZZ BASE DI-IL-CCTB TR (Conical cemented titan base non hex and screw)	BSAD11731	20			10-12
SET ZZ BASE DI-IL-CCTB TRW (Conical cemented titan base non hex and screw)	BSAD11732	20			10-12
SET ZZ BASE DI-IL-PCTBH HX2 (Parallel cemented titan base hex and screw)	BSAB11881	20			10-12
SET ZZ BASE DI-IL-PCTBH TR (Parallel cemented titan base hex and screw)	BSAB11781	20			10-12
SET ZZ BASE DI-IL-PCTBH TRW (Parallel cemented titan base hex and screw)	BSAB11782	20			10-12
SET ZZ BASE DI-IL-SA HX2 (Scananalog and screw)	ZBDE11831	15			10-12
SET ZZ BASE DI-IL-SA Pilastro PDM/PPM (Scananalog and screw)	ZBDE11832	15			10-12
SET ZZ BASE DI-IL-SA TR (Scananalog and screw)	ZBDE11731	15			10-12
SET ZZ BASE DI-IL-SA TRW (Scananalog and screw)	ZBDE11732	15			10-12
SET ZZ BASE DI-IL-SM HX2 (Scanmarker and screw)	ZBAD11831	20			10-12
SET ZZ BASE DI-IL-SM Pilastro PDM/PPM (Scanmarker and screw)	ZBAD11832	20			10-12
SET ZZ BASE DI-IL-SM TR (Scanmarker and screw)	ZBAD11731	20			10-12
SET ZZ BASE DI-IL-SM TRW (Scanmarker and screw)	ZBAD11732	20			10-12
ZZ BASE DI-IL-MI HX2 (Model implant)	MIAG11801				10-12
ZZ BASE DI-IL-MI Pilastro PDM/PPM (Model implant)	MIAG11802				10-12
ZZ BASE DI-IL-MI TR (Model implant)	MIAG11701				10-12
ZZ BASE DI-IL-MI TRW (Model implant)	MIAG11702				10-12
ZZ BASE DI-IL-SC HX2 (Screw)	BSAE11801	20			10-12
ZZ BASE DI-IL-SC Pilastro PDM/PPM (Screw)	BSAE11802	20			10-12
ZZ BASE DI-IL-SC TR (Screw)	BSAE11701	20			10-12
ZZ BASE DI-IL-SC TRW (Screw)	BSAE11702	20			10-12
ZZ BASE DI-IL-ZSC HX2 (Zirconia screw)	BSAE11831	20			10-12
ZZ BASE DI-IL-ZSC Pilastro PDM/PPM (Zirconia screw)	BSAE11832	20			10-12
ZZ BASE DI-IL-ZSC TR (Zirconia screw)	BSAE11731	20			10-12

DENTALTECH IMPLASSIC®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE DI-IL-ZSC TRW (Zirconia screw)	BSAE11732	20			10-12

DENTALTECH IMPLASSIC®/IMPLOGIC®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE DI-IL-ABCH 3.25 (Raw-Abutment® customizable hex and screw)	BSAF11691	20			10-12
SET ZZ BASE DI-IL-ABCH 3.75 (Raw-Abutment® customizable hex and screw)	BSAF11682	20			10-12
SET ZZ BASE DI-IL-CCTB 3.25 (Conical cemented titan base non hex and screw)	BSAD11631	20			10-12
SET ZZ BASE DI-IL-CCTB 3.75 (Conical cemented titan base non hex and screw)	BSAD11632	20			10-12
SET ZZ BASE DI-IL-CCTB Pilaströ E.M. 3.25 (Conical cemented titan base non hex and screw)	BSAD11633	20			10-12
SET ZZ BASE DI-IL-CCTB Pilaströ E.M. 3.75 (Conical cemented titan base non hex and screw)	BSAD11634	20			10-12
SET ZZ BASE DI-IL-ICOT Pilaströ E.M. 3.25 (Impression coping open tray non hex and Screw and Driver)	BSAQ11663	15			10-12
SET ZZ BASE DI-IL-ICOT Pilaströ E.M. 3.75 (Impression coping open tray non hex and Screw and Driver)	BSAQ11664	15			10-12
SET ZZ BASE DI-IL-ICOTH 3.25 (Impression coping open tray hex and Screw and Driver)	BSAQ11641	15			10-12
SET ZZ BASE DI-IL-ICOTH 3.75 (Impression coping open tray hex and Screw and Driver)	BSAQ11642	15			10-12
SET ZZ BASE DI-IL-NTB 3.25 (Narrow titan base non hex and screw)	BSAP11631	20			10-12
SET ZZ BASE DI-IL-NTB 3.75 (Narrow titan base non hex and screw)	BSAP11632	20			10-12
SET ZZ BASE DI-IL-NTBH 3.25 (Narrow titan base hex and screw)	BSAP11681	20			10-12
SET ZZ BASE DI-IL-NTBH 3.75 (Narrow titan base hex and screw)	BSAP11682	20			10-12
SET ZZ BASE DI-IL-PCTBH 3.25 (Parallel cemented titan base hex and screw)	BSAB11681	20			10-12
SET ZZ BASE DI-IL-PCTBH 3.75 (Parallel cemented titan base hex and screw)	BSAB11682	20			10-12
SET ZZ BASE DI-IL-SA 3.25 (Scanalog and screw)	ZBDE11631	15			10-12
SET ZZ BASE DI-IL-SA 3.75 (Scanalog and screw)	ZBDE11632	15			10-12
SET ZZ BASE DI-IL-SA Pilaströ E.M. 3.25 (Scanalog and screw)	ZBDE11633	15			10-12
SET ZZ BASE DI-IL-SA Pilaströ E.M. 3.75 (Scanalog and screw)	ZBDE11634	15			10-12
SET ZZ BASE DI-IL-SM 3.25 (Scanmarker and screw)	ZBAD11631	20			10-12
SET ZZ BASE DI-IL-SM 3.75 (Scanmarker and screw)	ZBAD11632	20			10-12
SET ZZ BASE DI-IL-SM Pilaströ E.M. 3.25 (Scanmarker and screw)	ZBAD11633	20			10-12
SET ZZ BASE DI-IL-SM Pilaströ E.M. 3.75 (Scanmarker and screw)	ZBAD11634	20			10-12
ZZ BASE DI-IL-ICSC 3.25 (Impression coping Screw)	BSFE11601	15			10-12
ZZ BASE DI-IL-ICSC 3.75 (Impression coping Screw)	BSFE11602	15			10-12
ZZ BASE DI-IL-ICSC Pilaströ E.M. 3.25 (Impression coping Screw)	BSFE11603	15			10-12
ZZ BASE DI-IL-ICSC Pilaströ E.M. 3.75 (Impression coping Screw)	BSFE11604	15			10-12
ZZ BASE DI-IL-MI 3.25 (Model implant)	MIAG11601				10-12

DENTALTECH IMPLASSIC®/IMPLAGIC®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
ZZ BASE DI-IL-MI 3.75 (Model implant)	MIAG11602				10-12
ZZ BASE DI-IL-MI Pilaströ E.M. 3.25 (Model implant)	MIAG11603				10-12
ZZ BASE DI-IL-MI Pilaströ E.M. 3.75 (Model implant)	MIAG11604				10-12
ZZ BASE DI-IL-SC 3.25 (Screw)	BSAE11601	20			10-12
ZZ BASE DI-IL-SC 3.75 (Screw)	BSAE11602	20			10-12
ZZ BASE DI-IL-SC Pilaströ E.M. 3.25 (Screw)	BSAE11603	20			10-12
ZZ BASE DI-IL-SC Pilaströ E.M. 3.75 (Screw)	BSAE11604	20			10-12
ZZ BASE DI-IL-WMSM-L12 3.25 (White Metal Scanmarker L12 with Screw)	ZBED11631L12	15			10-12
ZZ BASE DI-IL-WMSM-L12 3.75 (White Metal Scanmarker L12 with Screw)	ZBED11632L12	15			10-12
ZZ BASE DI-IL-WMSM-L12 Pilaströ E.M. 3.25 (White Metal Scanmarker L12 with Screw)	ZBED11633L12	15			10-12
ZZ BASE DI-IL-WMSM-L12 Pilaströ E.M. 3.75 (White Metal Scanmarker L12 with Screw)	ZBED11634L12	15			10-12
ZZ BASE DI-IL-ZSC 3.25 (Zirconia screw)	BSAE11631	20			10-12
ZZ BASE DI-IL-ZSC 3.75 (Zirconia screw)	BSAE11632	20			10-12
ZZ BASE DI-IL-ZSC Pilaströ E.M. 3.25 (Zirconia screw)	BSAE11633	20			10-12
ZZ BASE DI-IL-ZSC Pilaströ E.M. 3.75 (Zirconia screw)	BSAE11634	20			10-12

DENTIUM IMPLANTIUM / SUPERLINE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE D-I-ABCH 3.6 (3.4) (Raw-Abutment® customizable hex and screw)	BSAF4881	30			10-12
SET ZZ BASE D-I-ABCH 4.0 (3.8) (Raw-Abutment® customizable hex and screw)	BSAF4882	30			10-12
SET ZZ BASE D-I-ABCH 4.5 (4.3) (Raw-Abutment® customizable hex and screw)	BSAF4883	30			10-12
SET ZZ BASE D-I-ABCH 5.0 (4.8) (Raw-Abutment® customizable hex and screw)	BSAF4884	30			10-12
SET ZZ BASE D-I-ABCH-D14 3.6 (3.4) (Raw-Abutment® customizable hex D14 and screw)	BSAF4881D14	30			10-12
SET ZZ BASE D-I-ABCH-D14 4.0 (3.8) (Raw-Abutment® customizable hex D14 and screw)	BSAF4882D14	30			10-12
SET ZZ BASE D-I-ABCH-D14 4.5 (4.3) (Raw-Abutment® customizable hex D14 and screw)	BSAF4883D14	30			10-12
SET ZZ BASE D-I-ABCH-D14 5.0 (4.8) (Raw-Abutment® customizable hex D14 and screw)	BSAF4884D14	30			10-12
SET ZZ BASE D-I-ASCTBH-K80 3.6 (3.4) (Angled screw-channel titan base hex K80 and screw)	BSAT4881K80	20			10-12
SET ZZ BASE D-I-ASCTBH-K80 4.0 (3.8) (Angled screw-channel titan base hex K80 and screw)	BSAT4882K80	20			10-12
SET ZZ BASE D-I-ASCTBH-K80 4.5 (4.3) (Angled screw-channel titan base hex K80 and screw)	BSAT4883K80	20			10-12
SET ZZ BASE D-I-ASCTBH-K80 5.0 (4.8) (Angled screw-channel titan base hex K80 and screw)	BSAT4884K80	20			10-12
SET ZZ BASE D-I-ASCTB-K80 3.6 (3.4) (Angled screw-channel titan base non hex K80 and screw)	BSAT4831K80	20			10-12
SET ZZ BASE D-I-ASCTB-K80 4.0 (3.8) (Angled screw-channel titan base non hex K80 and screw)	BSAT4832K80	20			10-12
SET ZZ BASE D-I-ASCTB-K80 4.5 (4.3) (Angled screw-channel titan base non hex K80 and screw)	BSAT4833K80	20			10-12
SET ZZ BASE D-I-ASCTB-K80 5.0 (4.8) (Angled screw-channel titan base non hex K80 and screw)	BSAT4834K80	20			10-12
SET ZZ BASE D-I-CCTB 3.6 (3.4) (Conical cemented titan base non hex and screw)	BSAD4831	30			10-12
SET ZZ BASE D-I-CCTB 4.0 (3.8) (Conical cemented titan base non hex and screw)	BSAD4832	30			10-12
SET ZZ BASE D-I-CCTB 4.5 (4.3) (Conical cemented titan base non hex and screw)	BSAD4833	30			10-12
SET ZZ BASE D-I-CCTB 5.0 (4.8) (Conical cemented titan base non hex and screw)	BSAD4834	30			10-12
SET ZZ BASE D-I-CCTB-G 4.5 (4.3) (Conical cemented titan base non hex and screw Gold)	BSCD4833	30			10-12
SET ZZ BASE D-I-CCTB-K85 3.6 (3.4) (Conical cemented titan base non hex K85 and screw)	BSAD4831K85	30			10-12
SET ZZ BASE D-I-CCTB-K85 4.0 (3.8) (Conical cemented titan base non hex K85 and screw)	BSAD4832K85	30			10-12
SET ZZ BASE D-I-CCTB-K85 4.5 (4.3) (Conical cemented titan base non hex K85 and screw)	BSAD4833K85	30			10-12
SET ZZ BASE D-I-CCTB-K85 5.0 (4.8) (Conical cemented titan base non hex K85 and screw)	BSAD4834K85	30			10-12
SET ZZ BASE D-I-CCTB-L10 3.6 (3.4) (Conical cemented titan base non hex L10 and screw)	BSAD4831L10	30			10-12
SET ZZ BASE D-I-CCTB-L10 4.0 (3.8) (Conical cemented titan base non hex L10 and screw)	BSAD4832L10	30			10-12
SET ZZ BASE D-I-CCTB-L10 4.5 (4.3) (Conical cemented titan base non hex L10 and screw)	BSAD4833L10	30			10-12
SET ZZ BASE D-I-CCTB-L10 5.0 (4.8) (Conical cemented titan base non hex L10 and screw)	BSAD4834L10	30			10-12

DENTIUM IMPLANTIUM / SUPERLINE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE D-I-CCTB-L20 3.6 (3.4) (Conical cemented titan base non hex L20 and screw)	BSAD4831L20	30			10-12
SET ZZ BASE D-I-CCTB-L20 4.0 (3.8) (Conical cemented titan base non hex L20 and screw)	BSAD4832L20	30			10-12
SET ZZ BASE D-I-CCTB-L20 4.5 (4.3) (Conical cemented titan base non hex L20 and screw)	BSAD4833L20	30			10-12
SET ZZ BASE D-I-CCTB-L20 5.0 (4.8) (Conical cemented titan base non hex L20 and screw)	BSAD4834L20	30			10-12
SET ZZ BASE D-I-CCTB-L20-G 4.5 (4.3) (Conical cemented titan base non hex L20 and screw Go	BSCD4833L20	30			10-12
SET ZZ BASE D-I-CCTB-L30 3.6 (3.4) (Conical cemented titan base non hex L30 and screw)	BSAD4831L30	30			10-12
SET ZZ BASE D-I-CCTB-L30 4.0 (3.8) (Conical cemented titan base non hex L30 and screw)	BSAD4832L30	30			10-12
SET ZZ BASE D-I-CCTB-L30 4.5 (4.3) (Conical cemented titan base non hex L30 and screw)	BSAD4833L30	30			10-12
SET ZZ BASE D-I-CCTB-L30 5.0 (4.8) (Conical cemented titan base non hex L30 and screw)	BSAD4834L30	30			10-12
SET ZZ BASE D-I-CCTB-L40 3.6 (3.4) (Conical cemented titan base non hex L40 and screw)	BSAD4831L40	30			10-12
SET ZZ BASE D-I-CCTB-L40 4.0 (3.8) (Conical cemented titan base non hex L40 and screw)	BSAD4832L40	30			10-12
SET ZZ BASE D-I-CCTB-L40 4.5 (4.3) (Conical cemented titan base non hex L40 and screw)	BSAD4833L40	30			10-12
SET ZZ BASE D-I-CCTB-L40 5.0 (4.8) (Conical cemented titan base non hex L40 and screw)	BSAD4834L40	30			10-12
SET ZZ BASE D-I-ICOTH 3.6 (3.4) (Impression coping open tray hex and Screw and Driver)	BSAQ4841	15			10-12
SET ZZ BASE D-I-ICOTH 4.0 (3.8) (Impression coping open tray hex and Screw and Driver)	BSAQ4842	15			10-12
SET ZZ BASE D-I-ICOTH 4.5 (4.3) (Impression coping open tray hex and Screw and Driver)	BSAQ4843	15			10-12
SET ZZ BASE D-I-ICOTH 5.0 (4.8) (Impression coping open tray hex and Screw and Driver)	BSAQ4844	15			10-12
SET ZZ BASE D-I-NTB 3.6 (3.4) (Narrow titan base non hex and screw)	BSAP4831	30			10-12
SET ZZ BASE D-I-NTB 4.0 (3.8) (Narrow titan base non hex and screw)	BSAP4832	30			10-12
SET ZZ BASE D-I-NTB 4.5 (4.3) (Narrow titan base non hex and screw)	BSAP4833	30			10-12
SET ZZ BASE D-I-NTB 5.0 (4.8) (Narrow titan base non hex and screw)	BSAP4834	30			10-12
SET ZZ BASE D-I-NTBH 3.6 (3.4) (Narrow titan base hex and screw)	BSAP4881	30			10-12
SET ZZ BASE D-I-NTBH 4.0 (3.8) (Narrow titan base hex and screw)	BSAP4882	30			10-12
SET ZZ BASE D-I-NTBH 4.5 (4.3) (Narrow titan base hex and screw)	BSAP4883	30			10-12
SET ZZ BASE D-I-NTBH 5.0 (4.8) (Narrow titan base hex and screw)	BSAP4884	30			10-12
SET ZZ BASE D-I-PCTBH 3.6 (3.4) (Parallel cemented titan base hex and screw)	BSAB4881	30			10-12
SET ZZ BASE D-I-PCTBH 4.0 (3.8) (Parallel cemented titan base hex and screw)	BSAB4882	30			10-12
SET ZZ BASE D-I-PCTBH 4.5 (4.3) (Parallel cemented titan base hex and screw)	BSAB4883	30			10-12
SET ZZ BASE D-I-PCTBH 5.0 (4.8) (Parallel cemented titan base hex and screw)	BSAB4884	30			10-12

DENTIUM IMPLANTUM / SUPERLINE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE D-I-PCTBH-G 3.6 (3.4) (Parallel cemented titan base hex and screw Gold)	BSCB4881	30			10-12
SET ZZ BASE D-I-PCTBH-G 4.5 (4.3) (Parallel cemented titan base hex and screw Gold)	BSCB4883	30			10-12
SET ZZ BASE D-I-PCTBH-G 5.0 (4.8) (Parallel cemented titan base hex and screw Gold)	BSCB4884	30			10-12
SET ZZ BASE D-I-PCTBH-K85 3.6 (3.4) (Parallel cemented titan base hex K85 and screw)	BSAB4881K85	30			10-12
SET ZZ BASE D-I-PCTBH-K85 4.0 (3.8) (Parallel cemented titan base hex K85 and screw)	BSAB4882K85	30			10-12
SET ZZ BASE D-I-PCTBH-K85 4.5 (4.3) (Parallel cemented titan base hex K85 and screw)	BSAB4883K85	30			10-12
SET ZZ BASE D-I-PCTBH-K85 5.0 (4.8) (Parallel cemented titan base hex K85 and screw)	BSAB4884K85	30			10-12
SET ZZ BASE D-I-PCTBH-L10 3.6 (3.4) (Parallel cemented titan base hex L10 and screw)	BSAB4891L10	30			10-12
SET ZZ BASE D-I-PCTBH-L10 4.0 (3.8) (Parallel cemented titan base hex L10 and screw)	BSAB4892L10	30			10-12
SET ZZ BASE D-I-PCTBH-L10 4.5 (4.3) (Parallel cemented titan base hex L10 and screw)	BSAB4893L10	30			10-12
SET ZZ BASE D-I-PCTBH-L10 5.0 (4.8) (Parallel cemented titan base hex L10 and screw)	BSAB4894L10	30			10-12
SET ZZ BASE D-I-PCTBH-L20 3.6 (3.4) (Parallel cemented titan base hex L20 and screw)	BSAB4891L20	30			10-12
SET ZZ BASE D-I-PCTBH-L20 4.0 (3.8) (Parallel cemented titan base hex L20 and screw)	BSAB4892L20	30			10-12
SET ZZ BASE D-I-PCTBH-L20 4.5 (4.3) (Parallel cemented titan base hex L20 and screw)	BSAB4893L20	30			10-12
SET ZZ BASE D-I-PCTBH-L20 5.0 (4.8) (Parallel cemented titan base hex L20 and screw)	BSAB4894L20	30			10-12
SET ZZ BASE D-I-PCTBH-L30 3.6 (3.4) (Parallel cemented titan base hex L30 and screw)	BSAB4891L30	30			10-12
SET ZZ BASE D-I-PCTBH-L30 4.0 (3.8) (Parallel cemented titan base hex L30 and screw)	BSAB4892L30	30			10-12
SET ZZ BASE D-I-PCTBH-L30 4.5 (4.3) (Parallel cemented titan base hex L30 and screw)	BSAB4893L30	30			10-12
SET ZZ BASE D-I-PCTBH-L30 5.0 (4.8) (Parallel cemented titan base hex L30 and screw)	BSAB4894L30	30			10-12
SET ZZ BASE D-I-PCTBH-L40 3.6 (3.4) (Parallel cemented titan base hex L40 and screw)	BSAB4891L40	30			10-12
SET ZZ BASE D-I-PCTBH-L40 4.0 (3.8) (Parallel cemented titan base hex L40 and screw)	BSAB4892L40	30			10-12
SET ZZ BASE D-I-PCTBH-L40 4.5 (4.3) (Parallel cemented titan base hex L40 and screw)	BSAB4893L40	30			10-12
SET ZZ BASE D-I-PCTBH-L40 5.0 (4.8) (Parallel cemented titan base hex L40 and screw)	BSAB4894L40	30			10-12
SET ZZ BASE D-I-SA 3.6 (3.4) (Scanalog and screw)	ZBDE4831	15			10-12
SET ZZ BASE D-I-SA 4.0 (3.8) (Scanalog and screw)	ZBDE4832	15			10-12
SET ZZ BASE D-I-SA 4.5 (4.3) (Scanalog and screw)	ZBDE4833	15			10-12
SET ZZ BASE D-I-SA 5.0 (4.8) (Scanalog and screw)	ZBDE4834	15			10-12
SET ZZ BASE D-I-SM 3.6 (3.4) (Scanmarker and screw)	ZBAD4831	30			10-12
SET ZZ BASE D-I-SM 4.0 (3.8) (Scanmarker and screw)	ZBAD4832	30			10-12

DENTIUM IMPLANTUM / SUPERLINE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE D-I-SM 4.5 (4.3) (Scanmarker and screw)	ZBAD4833	30			10-12
SET ZZ BASE D-I-SM 5.0 (4.8) (Scanmarker and screw)	ZBAD4834	30			10-12
SET ZZ BASE D-I-WSM 3.6 (3.4) (White Scanmarker and screw)	ZBCD4831	30			10-12
SET ZZ BASE D-I-WSM 4.0 (3.8) (White Scanmarker and screw)	ZBCD4832	30			10-12
SET ZZ BASE D-I-WSM 4.5 (4.3) (White Scanmarker and screw)	ZBCD4833	30			10-12
SET ZZ BASE D-I-WSM 5.0 (4.8) (White Scanmarker and screw)	ZBCD4834	30			10-12
ZZ BASE D-I-ASCSC 3.6 (3.4) (Angled screw-channel Screw)	BSKE4801	20			10-12
ZZ BASE D-I-ASCSC 4.0 (3.8) (Angled screw-channel Screw)	BSKE4802	20			10-12
ZZ BASE D-I-ASCSC 4.5 (4.3) (Angled screw-channel Screw)	BSKE4803	20			10-12
ZZ BASE D-I-ASCSC 5.0 (4.8) (Angled screw-channel Screw)	BSKE4804	20			10-12
ZZ BASE D-I-BSC 3.6 (3.4) (Black Screw)	BSIE4801	30			10-12
ZZ BASE D-I-BSC 4.0 (3.8) (Black Screw)	BSIE4802	30			10-12
ZZ BASE D-I-BSC 4.5 (4.3) (Black Screw)	BSIE4803	30			10-12
ZZ BASE D-I-BSC 5.0 (4.8) (Black Screw)	BSIE4804	30			10-12
ZZ BASE D-I-BZSC 3.6 (3.4) (Black Zirconia screw)	BSIE4831	30			10-12
ZZ BASE D-I-BZSC 4.0 (3.8) (Black Zirconia screw)	BSIE4832	30			10-12
ZZ BASE D-I-BZSC 4.5 (4.3) (Black Zirconia screw)	BSIE4833	30			10-12
ZZ BASE D-I-BZSC 5.0 (4.8) (Black Zirconia screw)	BSIE4834	30			10-12
ZZ BASE D-I-ICSC 3.6 (3.4) (Impression coping Screw)	BSFE4801	15			10-12
ZZ BASE D-I-ICSC 4.0 (3.8) (Impression coping Screw)	BSFE4802	15			10-12
ZZ BASE D-I-ICSC 4.5 (4.3) (Impression coping Screw)	BSFE4803	15			10-12
ZZ BASE D-I-ICSC 5.0 (4.8) (Impression coping Screw)	BSFE4804	15			10-12
ZZ BASE D-I-LBSC 3.6 (3.4) (Lab Screw)	BSHE4801	30			10-12
ZZ BASE D-I-LBSC 4.0 (3.8) (Lab Screw)	BSHE4802	30			10-12
ZZ BASE D-I-LBSC 4.5 (4.3) (Lab Screw)	BSHE4803	30			10-12
ZZ BASE D-I-LBSC 5.0 (4.8) (Lab Screw)	BSHE4804	30			10-12
ZZ BASE D-I-LBZSC 3.6 (3.4) (Lab Zirconia screw)	BSHE4831	30			10-12
ZZ BASE D-I-LBZSC 4.0 (3.8) (Lab Zirconia screw)	BSHE4832	30			10-12
ZZ BASE D-I-LBZSC 4.5 (4.3) (Lab Zirconia screw)	BSHE4833	30			10-12

DENTIUM IMPLANTUM / SUPERLINE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE D-I-LBZSC 5.0 (4.8) (Lab Zirconia screw)	BSHE4834	30			10-12
ZZ BASE D-I-MI 3.6 (3.4) (Model implant)	MIAG4801				10-12
ZZ BASE D-I-MI 4.0 (3.8) (Model implant)	MIAG4802				10-12
ZZ BASE D-I-MI 4.5 (4.3) (Model implant)	MIAG4803				10-12
ZZ BASE D-I-MI 5.0 (4.8) (Model implant)	MIAG4804				10-12
ZZ BASE D-I-SC 3.6 (3.4) (Screw)	BSAE4801	30			10-12
ZZ BASE D-I-SC 4.0 (3.8) (Screw)	BSAE4802	30			10-12
ZZ BASE D-I-SC 4.5 (4.3) (Screw)	BSAE4803	30			10-12
ZZ BASE D-I-SC 5.0 (4.8) (Screw)	BSAE4804	30			10-12
ZZ BASE D-I-SC-G 3.6 (3.4) (Screw Gold)	BSCE4801	30			10-12
ZZ BASE D-I-SC-G 4.0 (3.8) (Screw Gold)	BSCE4802	30			10-12
ZZ BASE D-I-WMSM-L12 3.6 (3.4) (White Metal Scanmarker L12 with Screw)	ZBED4831L12	15			10-12
ZZ BASE D-I-WMSM-L12 4.0 (3.8) (White Metal Scanmarker L12 with Screw)	ZBED4832L12	15			10-12
ZZ BASE D-I-WMSM-L12 4.5 (4.3) (White Metal Scanmarker L12 with Screw)	ZBED4833L12	15			10-12
ZZ BASE D-I-WMSM-L12 5.0 (4.8) (White Metal Scanmarker L12 with Screw)	ZBED4834L12	15			10-12
ZZ BASE D-I-ZSC 3.6 (3.4) (Zirconia screw)	BSAE4831	30			10-12
ZZ BASE D-I-ZSC 4.0 (3.8) (Zirconia screw)	BSAE4832	30			10-12
ZZ BASE D-I-ZSC 4.5 (4.3) (Zirconia screw)	BSAE4833	30			10-12
ZZ BASE D-I-ZSC 5.0 (4.8) (Zirconia screw)	BSAE4834	30			10-12
ZZ BASE D-I-ZSC-G 3.6 (3.4) (Zirconia Screw Gold)	BSCE4831	30			10-12
ZZ BASE D-I-ZSC-G 4.0 (3.8) (Zirconia Screw Gold)	BSCE4832	30			10-12
ZZ BASE D-I-ZSC-G 4.5 (4.3) (Zirconia Screw Gold)	BSCE4833	30			10-12

DENTIUM SCREW ABUTMENT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE D-SA-CCTB 4.5 (Conical cemented titan base non hex and screw)	BSAD4931	30			10-12
SET ZZ BASE D-SA-CCTB 5.5 (Conical cemented titan base non hex and screw)	BSAD4932	30			10-12
SET ZZ BASE D-SA-CCTB-K85 4.5 (Conical cemented titan base non hex K85 and screw)	BSAD4931K85	30			10-12
SET ZZ BASE D-SA-SA 4.5 (Scananalog and screw)	ZBDE4931	15			10-12
SET ZZ BASE D-SA-SA 5.5 (Scananalog and screw)	ZBDE4932	15			10-12
SET ZZ BASE D-SA-SM 4.5 (Scanmarker and screw)	ZBAD4931	30			10-12
SET ZZ BASE D-SA-SM 5.5 (Scanmarker and screw)	ZBAD4932	30			10-12
ZZ BASE D-SA-SC 4.5 (Screw)	BSAE4901	30			10-12
ZZ BASE D-SA-SC 5.5 (Screw)	BSAE4902	30			10-12
ZZ BASE D-SA-ZSC 4.5 (Zirconia screw)	BSAE4931	30			10-12
ZZ BASE D-SA-ZSC 5.5 (Zirconia screw)	BSAE4932	30			10-12

DENTSPLY SIRONA® ANKYLOS® (FRIADENT)

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE F-DA-ABC 3.5 (Raw-Abutment® customizable non hex and screw)	BSAF3031	15			10-12
SET ZZ BASE F-DA-ABC 4.5 (Raw-Abutment® customizable non hex and screw)	BSAF3032	15			10-12
SET ZZ BASE F-DA-ABC 5.5 (Raw-Abutment® customizable non hex and screw)	BSAF3033	15			10-12
SET ZZ BASE F-DA-ABC 7.0 (Raw-Abutment® customizable non hex and screw)	BSAF3034	15			10-12
SET ZZ BASE F-DA-ABCH 3.5 (Raw-Abutment® customizable hex and screw)	BSAF3091	15			10-12
SET ZZ BASE F-DA-ABCH 4.5 (Raw-Abutment® customizable hex and screw)	BSAF3092	15			10-12
SET ZZ BASE F-DA-ABCH 5.5 (Raw-Abutment® customizable hex and screw)	BSAF3093	15			10-12
SET ZZ BASE F-DA-ABCH 7.0 (Raw-Abutment® customizable hex and screw)	BSAF3094	15			10-12
SET ZZ BASE F-DA-ABCH-D14 3.5 (Raw-Abutment® customizable hex D14 and screw)	BSAF3091D14	15			10-12
SET ZZ BASE F-DA-ABCH-D14 4.5 (Raw-Abutment® customizable hex D14 and screw)	BSAF3092D14	15			10-12
SET ZZ BASE F-DA-ABCH-D14 5.5 (Raw-Abutment® customizable hex D14 and screw)	BSAF3093D14	15			10-12
SET ZZ BASE F-DA-ABCH-D14 7.0 (Raw-Abutment® customizable hex D14 and screw)	BSAF3094D14	15			10-12
SET ZZ BASE F-DA-ASCTB-K80 3.5 (Angled screw-channel titan base non hex K80 and screw)	BSAT3031K80	15			10-12
SET ZZ BASE F-DA-ASCTB-K80 4.5 (Angled screw-channel titan base non hex K80 and screw)	BSAT3032K80	15			10-12
SET ZZ BASE F-DA-CCTB 3.5 (Conical cemented titan base non hex and screw)	BSAD3031	15			10-12
SET ZZ BASE F-DA-CCTB 4.5 (Conical cemented titan base non hex and screw)	BSAD3032	15			10-12
SET ZZ BASE F-DA-CCTB 5.5 (Conical cemented titan base non hex and screw)	BSAD3033	15			10-12
SET ZZ BASE F-DA-CCTB 7.0 (Conical cemented titan base non hex and screw)	BSAD3034	15			10-12
SET ZZ BASE F-DA-CCTB-K85 3.5 (Conical cemented titan base non hex K85 and screw)	BSAD3031K85	15			10-12
SET ZZ BASE F-DA-CCTB-K85 4.5 (Conical cemented titan base non hex K85 and screw)	BSAD3032K85	15			10-12
SET ZZ BASE F-DA-CCTB-K85 5.5 (Conical cemented titan base non hex K85 and screw)	BSAD3033K85	15			10-12
SET ZZ BASE F-DA-CCTB-K85 7.0 (Conical cemented titan base non hex K85 and screw)	BSAD3034K85	15			10-12
SET ZZ BASE F-DA-ICOT 3.5 (Impression coping open tray non hex and Screw and Driver)	BSAQ3061	15			10-12
SET ZZ BASE F-DA-ICOT 4.5 (Impression coping open tray non hex and Screw and Driver)	BSAQ3062	15			10-12
SET ZZ BASE F-DA-ICOT 5.5 (Impression coping open tray non hex and Screw and Driver)	BSAQ3063	15			10-12
SET ZZ BASE F-DA-ICOT 7.0 (Impression coping open tray non hex and Screw and Driver)	BSAQ3064	15			10-12
SET ZZ BASE F-DA-NTB 3.5 (Narrow titan base non hex and screw)	BSAP3031	15			10-12
SET ZZ BASE F-DA-NTB 4.5 (Narrow titan base non hex and screw)	BSAP3032	15			10-12
SET ZZ BASE F-DA-PCTBH 3.5 (Parallel cemented titan base hex and screw)	BSAB3091	15			10-12

DENTSPLY SIRONA® ANKYLOS® (FRIADENT)

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE F-DA-PCTBH 4.5 (Parallel cemented titan base hex and screw)	BSAB3092	15			10-12
SET ZZ BASE F-DA-PCTBH 5.5 (Parallel cemented titan base hex and screw)	BSAB3093	15			10-12
SET ZZ BASE F-DA-PCTBH 7.0 (Parallel cemented titan base hex and screw)	BSAB3094	15			10-12
SET ZZ BASE F-DA-SA 3.5 (Scananalog and screw)	ZBDE3031	15			10-12
SET ZZ BASE F-DA-SA 4.5 (Scananalog and screw)	ZBDE3032	15			10-12
SET ZZ BASE F-DA-SA 5.5 (Scananalog and screw)	ZBDE3033	15			10-12
SET ZZ BASE F-DA-SA 7.0 (Scananalog and screw)	ZBDE3034	15			10-12
SET ZZ BASE F-DA-SM 3.5 (Scanmarker and screw)	ZBAD3031	15			10-12
SET ZZ BASE F-DA-SM 4.5 (Scanmarker and screw)	ZBAD3032	15			10-12
SET ZZ BASE F-DA-SM 5.5 (Scanmarker and screw)	ZBAD3033	15			10-12
SET ZZ BASE F-DA-SM 7.0 (Scanmarker and screw)	ZBAD3034	15			10-12
SET ZZ BASE F-DA-WSM 5.5 (White Scanmarker and screw)	ZBCD3033	15			10-12
SET ZZ BASE F-DA-WSM 7.0 (White Scanmarker and screw)	ZBCD3034	15			10-12
SET ZZ BASE F-DA-ZZMUA-L15 3.5 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH3031L15	25	35		10-12
SET ZZ BASE F-DA-ZZMUA-L15 4.5 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH3032L15	25	35		10-12
SET ZZ BASE F-DA-ZZMUA-L15 5.5 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH3033L15	25	35		10-12
SET ZZ BASE F-DA-ZZMUA-L15 7.0 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH3034L15	25	35		10-12
SET ZZ BASE F-DA-ZZMUA-L25 3.5 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH3031L25	25	35		10-12
SET ZZ BASE F-DA-ZZMUA-L25 4.5 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH3032L25	25	35		10-12
SET ZZ BASE F-DA-ZZMUA-L25 5.5 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH3033L25	25	35		10-12
SET ZZ BASE F-DA-ZZMUA-L25 7.0 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH3034L25	25	35		10-12
SET ZZ BASE F-DA-ZZMUA-L35 3.5 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH3031L35	25	35		10-12
SET ZZ BASE F-DA-ZZMUA-L35 4.5 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH3032L35	25	35		10-12
SET ZZ BASE F-DA-ZZMUA-L35 5.5 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH3033L35	25	35		10-12
SET ZZ BASE F-DA-ZZMUA-L35 7.0 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH3034L35	25	35		10-12
SET ZZ BASE F-DA-ZZMUA-L45 3.5 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH3031L45	25	35		10-12
SET ZZ BASE F-DA-ZZMUA-L45 4.5 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH3032L45	25	35		10-12
SET ZZ BASE F-DA-ZZMUA-L45 5.5 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH3033L45	25	35		10-12
SET ZZ BASE F-DA-ZZMUA-L45 7.0 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH3034L45	25	35		10-12

All information is subject to change. Errors and omissions excepted. Version: 21.09.2023

DENTSPLY SIRONA® ANKYLOS® (FRIADENT)

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE F-DA-ASCSC 3.5 (Angled screw-channel Screw)	BSKE3001	15			10-12
ZZ BASE F-DA-ASCSC 4.5 (Angled screw-channel Screw)	BSKE3002	15			10-12
ZZ BASE F-DA-ICSC 3.5 (Impression coping Screw)	BSFE3001	15			10-12
ZZ BASE F-DA-ICSC 4.5 (Impression coping Screw)	BSFE3002	15			10-12
ZZ BASE F-DA-ICSC 5.5 (Impression coping Screw)	BSFE3003	15			10-12
ZZ BASE F-DA-ICSC 7.0 (Impression coping Screw)	BSFE3004	15			10-12
ZZ BASE F-DA-LSC 3.5 (Long screw)	BSAE3081	15			10-12
ZZ BASE F-DA-LSC 4.5 (Long screw)	BSAE3082	15			10-12
ZZ BASE F-DA-LSC 5.5 (Long screw)	BSAE3083	15			10-12
ZZ BASE F-DA-LSC 7.0 (Long screw)	BSAE3084	15			10-12
ZZ BASE F-DA-SC 3.5 (Screw)	BSAE3001	15			10-12
ZZ BASE F-DA-SC 4.5 (Screw)	BSAE3002	15			10-12
ZZ BASE F-DA-SC 5.5 (Screw)	BSAE3003	15			10-12
ZZ BASE F-DA-SC 7.0 (Screw)	BSAE3004	15			10-12
ZZ BASE F-DA-SC-G 3.5 (Screw Gold)	BSCE3001	15			10-12
ZZ BASE F-DA-WMSM-L10 Hex 3.5 (White Metal Scanmarker L10 with Screw)	ZBED3031L10	15			10-12
ZZ BASE F-DA-WMSM-L10 Hex 4.5 (White Metal Scanmarker L10 with Screw)	ZBED3032L10	15			10-12
ZZ BASE F-DA-WMSM-L10 Hex 5.5 (White Metal Scanmarker L10 with Screw)	ZBED3033L10	15			10-12
ZZ BASE F-DA-WMSM-L10 Hex 7.0 (White Metal Scanmarker L10 with Screw)	ZBED3034L10	15			10-12
ZZ BASE F-DA-WMSM-L12 Hex 3.5 (White Metal Scanmarker L12 with Screw)	ZBED3031L12	15			10-12
ZZ BASE F-DA-WMSM-L12 Hex 4.5 (White Metal Scanmarker L12 with Screw)	ZBED3032L12	15			10-12
ZZ BASE F-DA-WMSM-L12 Hex 5.5 (White Metal Scanmarker L12 with Screw)	ZBED3033L12	15			10-12
ZZ BASE F-DA-WMSM-L12 Hex 7.0 (White Metal Scanmarker L12 with Screw)	ZBED3034L12	15			10-12
ZZ BASE F-DA-WMSM-L15 Hex 3.5 (White Metal Scanmarker L15 with Screw)	ZBED3031L15	15			10-12
ZZ BASE F-DA-WMSM-L15 Hex 4.5 (White Metal Scanmarker L15 with Screw)	ZBED3032L15	15			10-12
ZZ BASE F-DA-WMSM-L15 Hex 5.5 (White Metal Scanmarker L15 with Screw)	ZBED3033L15	15			10-12
ZZ BASE F-DA-WMSM-L15 Hex 7.0 (White Metal Scanmarker L15 with Screw)	ZBED3034L15	15			10-12

DENTSPLY SIRONA® ANKYLOS® (FRIADENT) / BALANCE BASE ABUTMENT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE F-DA-CCTB BBAN (Conical cemented titan base non hex and screw)	BSAD8231	10			10-12
SET ZZ BASE F-DA-CCTB-K85 BBAN (Conical cemented titan base non hex K85 and screw)	BSAD8231K85	10			10-12
SET ZZ BASE F-DA-ICOT BBAN (Impression coping open tray non hex and Screw and Driver)	BSAQ8261	15			10-12
SET ZZ BASE F-DA-SA BBAN (Scanalog and screw)	ZBDE8231	15			10-12
SET ZZ BASE F-DA-SM BBAN (Scanmarker and screw)	ZBAD8231	10			10-12
SET ZZ BASE F-DA-WSM BBAN (White Scanmarker and screw)	ZBCD8231	10			10-12
ZZ BASE F-DA-ICSC BBAN (Impression coping Screw)	BSFE8201	15			10-12
ZZ BASE F-DA-SC BBAN (Screw)	BSAE8201	10			10-12
ZZ BASE F-DA-WMSM-L12 BBAN (White Metal Scanmarker L12 with Screw)	ZBED8231L12	15			10-12
ZZ BASE F-DA-ZSC BBAN (Zirconia Screw)	BSAE8231	10			10-12

DENTSPLY SIRONA® XIVE® MP/TG

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE F-MP/TG-CCTB 3.4 (Conical cemented titan base non hex and screw)	BSAD4031	24			10-12
SET ZZ BASE F-MP/TG-CCTB 3.8 (Conical cemented titan base non hex and screw)	BSAD4032	24			10-12
SET ZZ BASE F-MP/TG-CCTB 4.5 (Conical cemented titan base non hex and screw)	BSAD4033	24			10-12
SET ZZ BASE F-MP/TG-CCTB 5.5 (Conical cemented titan base non hex and screw)	BSAD4034	24			10-12
SET ZZ BASE F-MP/TG-CCTB-K85 3.4 (Conical cemented titan base non hex K85 and screw)	BSAD4031K85	24			10-12
SET ZZ BASE F-MP/TG-CCTB-K85 3.8 (Conical cemented titan base non hex K85 and screw)	BSAD4032K85	24			10-12
SET ZZ BASE F-MP/TG-CCTB-K85 4.5 (Conical cemented titan base non hex K85 and screw)	BSAD4033K85	24			10-12
SET ZZ BASE F-MP/TG-CCTB-K85 5.5 (Conical cemented titan base non hex K85 and screw)	BSAD4034K85	24			10-12
SET ZZ BASE F-MP/TG-ICOT 3.4 (Impression coping open tray non hex and Screw and Driver)	BSAQ4061	15			10-12
SET ZZ BASE F-MP/TG-ICOT 3.8 (Impression coping open tray non hex and Screw and Driver)	BSAQ4062	15			10-12
SET ZZ BASE F-MP/TG-ICOT 4.5 (Impression coping open tray non hex and Screw and Driver)	BSAQ4063	15			10-12
SET ZZ BASE F-MP/TG-ICOT 5.5 (Impression coping open tray non hex and Screw and Driver)	BSAQ4064	15			10-12
SET ZZ BASE F-MP/TG-SA 3.4 (Scananalog and screw)	ZBDE4031	15			10-12
SET ZZ BASE F-MP/TG-SA 3.8 (Scananalog and screw)	ZBDE4032	15			10-12
SET ZZ BASE F-MP/TG-SA 4.5 (Scananalog and screw)	ZBDE4033	15			10-12
SET ZZ BASE F-MP/TG-SA 5.5 (Scananalog and screw)	ZBDE4034	15			10-12
SET ZZ BASE F-MP/TG-SM 3.4 (Scanmarker and screw)	ZBAD4031	24			10-12
SET ZZ BASE F-MP/TG-SM 3.8 (Scanmarker and screw)	ZBAD4032	24			10-12
SET ZZ BASE F-MP/TG-SM 4.5 (Scanmarker and screw)	ZBAD4033	24			10-12
SET ZZ BASE F-MP/TG-SM 5.5 (Scanmarker and screw)	ZBAD4034	24			10-12
ZZ BASE F-MP/TG-ICSC 3.4 (Impression coping Screw)	BSFE4001	15			10-12
ZZ BASE F-MP/TG-ICSC 3.8 (Impression coping Screw)	BSFE4002	15			10-12
ZZ BASE F-MP/TG-ICSC 4.5 (Impression coping Screw)	BSFE4003	15			10-12
ZZ BASE F-MP/TG-ICSC 5.5 (Impression coping Screw)	BSFE4004	15			10-12
ZZ BASE F-MP/TG-SC 3.4 (Screw)	BSAE4001	24			10-12
ZZ BASE F-MP/TG-SC 3.8 (Screw)	BSAE4002	24			10-12
ZZ BASE F-MP/TG-SC 4.5 (Screw)	BSAE4003	24			10-12
ZZ BASE F-MP/TG-SC 5.5 (Screw)	BSAE4004	24			10-12
ZZ BASE F-MP/TG-WMSM-L12 3.4 (White Metal Scanmarker L12 with Screw)	ZBED4031L12	15			10-12

DENTSPLY SIRONA® XIVE® MP/TG

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE F-MP/TG-WMSM-L12 3.8 (White Metal Scanmarker L12 with Screw)	ZBED4032L12	15			10-12
ZZ BASE F-MP/TG-WMSM-L12 4.5 (White Metal Scanmarker L12 with Screw)	ZBED4033L12	15			10-12
ZZ BASE F-MP/TG-WMSM-L12 5.5 (White Metal Scanmarker L12 with Screw)	ZBED4034L12	15			10-12
ZZ BASE F-MP/TG-ZSC 3.4 (Zirconia screw)	BSAE4031	24			10-12
ZZ BASE F-MP/TG-ZSC 3.8 (Zirconia screw)	BSAE4032	24			10-12
ZZ BASE F-MP/TG-ZSC 4.5 (Zirconia screw)	BSAE4033	24			10-12
ZZ BASE F-MP/TG-ZSC 5.5 (Zirconia screw)	BSAE4034	24			10-12

DENTSPLY SIRONA® XIVE®/FRIALIT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE F-DX-ABCH 3.0 (Raw-Abutment® customizable hex and screw)	BSAF1095	24			10-12
SET ZZ BASE F-DX-ABCH 3.4 (Raw-Abutment® customizable hex and screw)	BSAF1091	25			10-12
SET ZZ BASE F-DX-ABCH 3.8 (Raw-Abutment® customizable hex and screw)	BSAF1092	25			10-12
SET ZZ BASE F-DX-ABCH 4.5 (Raw-Abutment® customizable hex and screw)	BSAF1093	25			10-12
SET ZZ BASE F-DX-ABCH 5.5 (Raw-Abutment® customizable hex and screw)	BSAF1094	25			10-12
SET ZZ BASE F-DX-ABCH-D14 3.8 (Raw-Abutment® customizable hex D14 and screw)	BSAF1092D14	25			10-12
SET ZZ BASE F-DX-ASCTBH-K80 3.0 (Angled screw-channel titan base hex K80 and screw)	BSAT1085K80	20			10-12
SET ZZ BASE F-DX-ASCTBH-K80 3.4 (Angled screw-channel titan base hex K80 and screw)	BSAT1081K80	20			10-12
SET ZZ BASE F-DX-ASCTBH-K80 3.8 (Angled screw-channel titan base hex K80 and screw)	BSAT1082K80	20			10-12
SET ZZ BASE F-DX-ASCTBH-K80 4.5 (Angled screw-channel titan base hex K80 and screw)	BSAT1083K80	20			10-12
SET ZZ BASE F-DX-ASCTB-K80 3.0 (Angled screw-channel titan base non hex K80 and screw)	BSAT1035K80	20			10-12
SET ZZ BASE F-DX-ASCTB-K80 3.4 (Angled screw-channel titan base non hex K80 and screw)	BSAT1031K80	20			10-12
SET ZZ BASE F-DX-ASCTB-K80 3.8 (Angled screw-channel titan base non hex K80 and screw)	BSAT1032K80	20			10-12
SET ZZ BASE F-DX-ASCTB-K80 4.5 (Angled screw-channel titan base non hex K80 and screw)	BSAT1033K80	20			10-12
SET ZZ BASE F-DX-CCTB 3.0 (Conical cemented titan base non hex and screw)	BSAD1035	24			10-12
SET ZZ BASE F-DX-CCTB 3.4 (Conical cemented titan base non hex and screw)	BSAD1031	25			10-12
SET ZZ BASE F-DX-CCTB 3.8 (Conical cemented titan base non hex and screw)	BSAD1032	25			10-12
SET ZZ BASE F-DX-CCTB 4.5 (Conical cemented titan base non hex and screw)	BSAD1033	25			10-12
SET ZZ BASE F-DX-CCTB 5.5 (Conical cemented titan base non hex and screw)	BSAD1034	25			10-12
SET ZZ BASE F-DX-CCTB-G 3.0 (Conical cemented titan base non hex and screw Gold)	BSCD1035	24			10-12
SET ZZ BASE F-DX-CCTB-G 3.4 (Conical cemented titan base non hex and screw Gold)	BSCD1031	25			10-12
SET ZZ BASE F-DX-CCTB-G 3.8 (Conical cemented titan base non hex and screw Gold)	BSCD1032	25			10-12
SET ZZ BASE F-DX-CCTB-G 4.5 (Conical cemented titan base non hex and screw Gold)	BSCD1033	25			10-12
SET ZZ BASE F-DX-CCTB-G 5.5 (Conical cemented titan base non hex and screw Gold)	BSCD1034	25			10-12
SET ZZ BASE F-DX-CCTB-K85 3.0 (Conical cemented titan base non hex K85 and screw)	BSAD1035K85	24			10-12
SET ZZ BASE F-DX-CCTB-K85 3.4 (Conical cemented titan base non hex K85 and screw)	BSAD1031K85	25			10-12
SET ZZ BASE F-DX-CCTB-K85 3.8 (Conical cemented titan base non hex K85 and screw)	BSAD1032K85	25			10-12
SET ZZ BASE F-DX-CCTB-K85 4.5 (Conical cemented titan base non hex K85 and screw)	BSAD1033K85	25			10-12
SET ZZ BASE F-DX-CCTB-K85 5.5 (Conical cemented titan base non hex K85 and screw)	BSAD1034K85	25			10-12

DENTSPLY SIRONA® XIVE®/FRIALIT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE F-DX-CCTB-L10 3.0 (Conical cemented titan base non hex L10 and screw)	BSAD1035L10	24			10-12
SET ZZ BASE F-DX-CCTB-L10 3.4 (Conical cemented titan base non hex L10 and screw)	BSAD1031L10	25			10-12
SET ZZ BASE F-DX-CCTB-L10 3.8 (Conical cemented titan base non hex L10 and screw)	BSAD1032L10	25			10-12
SET ZZ BASE F-DX-CCTB-L10 4.5 (Conical cemented titan base non hex L10 and screw)	BSAD1033L10	25			10-12
SET ZZ BASE F-DX-CCTB-L10 5.5 (Conical cemented titan base non hex L10 and screw)	BSAD1034L10	25			10-12
SET ZZ BASE F-DX-CCTB-L10-G 3.8 (Conical cemented titan base non hex L10 and screw Gold)	BSCD1032L10	25			10-12
SET ZZ BASE F-DX-CCTB-L10-G 4.5 (Conical cemented titan base non hex L10 and screw Gold)	BSCD1033L10	25			10-12
SET ZZ BASE F-DX-CCTB-L20 3.0 (Conical cemented titan base non hex L20 and screw)	BSAD1035L20	24			10-12
SET ZZ BASE F-DX-CCTB-L20 3.4 (Conical cemented titan base non hex L20 and screw)	BSAD1031L20	25			10-12
SET ZZ BASE F-DX-CCTB-L20 3.8 (Conical cemented titan base non hex L20 and screw)	BSAD1032L20	25			10-12
SET ZZ BASE F-DX-CCTB-L20 4.5 (Conical cemented titan base non hex L20 and screw)	BSAD1033L20	25			10-12
SET ZZ BASE F-DX-CCTB-L20 5.5 (Conical cemented titan base non hex L20 and screw)	BSAD1034L20	25			10-12
SET ZZ BASE F-DX-CCTB-L20-G 3.4 (Conical cemented titan base non hex L20 and screw Gold)	BSCD1031L20	25			10-12
SET ZZ BASE F-DX-CCTB-L20-G 3.8 (Conical cemented titan base non hex L20 and screw Gold)	BSCD1032L20	25			10-12
SET ZZ BASE F-DX-CCTB-L30 3.0 (Conical cemented titan base non hex L30 and screw)	BSAD1035L30	24			10-12
SET ZZ BASE F-DX-CCTB-L30 3.4 (Conical cemented titan base non hex L30 and screw)	BSAD1031L30	25			10-12
SET ZZ BASE F-DX-CCTB-L30 3.8 (Conical cemented titan base non hex L30 and screw)	BSAD1032L30	25			10-12
SET ZZ BASE F-DX-CCTB-L30 4.5 (Conical cemented titan base non hex L30 and screw)	BSAD1033L30	25			10-12
SET ZZ BASE F-DX-CCTB-L30 5.5 (Conical cemented titan base non hex L30 and screw)	BSAD1034L30	25			10-12
SET ZZ BASE F-DX-CCTB-L40 3.0 (Conical cemented titan base non hex L40 and screw)	BSAD1035L40	24			10-12
SET ZZ BASE F-DX-CCTB-L40 3.4 (Conical cemented titan base non hex L40 and screw)	BSAD1031L40	25			10-12
SET ZZ BASE F-DX-CCTB-L40 3.8 (Conical cemented titan base non hex L40 and screw)	BSAD1032L40	25			10-12
SET ZZ BASE F-DX-CCTB-L40 4.5 (Conical cemented titan base non hex L40 and screw)	BSAD1033L40	25			10-12
SET ZZ BASE F-DX-CCTB-L40 5.5 (Conical cemented titan base non hex L40 and screw)	BSAD1034L40	25			10-12
SET ZZ BASE F-DX-ICOT 3.0 (Impression coping open tray non hex and Screw and Driver)	BSAQ1065	15			10-12
SET ZZ BASE F-DX-ICOT 3.4 (Impression coping open tray non hex and Screw and Driver)	BSAQ1061	15			10-12
SET ZZ BASE F-DX-ICOT 3.8 (Impression coping open tray non hex and Screw and Driver)	BSAQ1062	15			10-12
SET ZZ BASE F-DX-ICOT 4.5 (Impression coping open tray non hex and Screw and Driver)	BSAQ1063	15			10-12
SET ZZ BASE F-DX-ICOT 5.5 (Impression coping open tray non hex and Screw and Driver)	BSAQ1064	15			10-12

DENTSPLY SIRONA® XIVE®/FRIALIT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE F-DX-ICOTH 3.0 (Impression coping open tray hex and Screw and Driver)	BSAQ1045	15			10-12
SET ZZ BASE F-DX-ICOTH 3.4 (Impression coping open tray hex and Screw and Driver)	BSAQ1041	15			10-12
SET ZZ BASE F-DX-ICOTH 3.8 (Impression coping open tray hex and Screw and Driver)	BSAQ1042	15			10-12
SET ZZ BASE F-DX-ICOTH 4.5 (Impression coping open tray hex and Screw and Driver)	BSAQ1043	15			10-12
SET ZZ BASE F-DX-ICOTH 5.5 (Impression coping open tray hex and Screw and Driver)	BSAQ1044	15			10-12
SET ZZ BASE F-DX-NTB 3.0 (Narrow titan base non hex and screw)	BSAP1035	24			10-12
SET ZZ BASE F-DX-NTB 3.4 (Narrow titan base non hex and screw)	BSAP1031	25			10-12
SET ZZ BASE F-DX-NTB 3.8 (Narrow titan base non hex and screw)	BSAP1032	25			10-12
SET ZZ BASE F-DX-NTBH 3.0 (Narrow titan base hex and screw)	BSAP1085	24			10-12
SET ZZ BASE F-DX-NTBH 3.4 (Narrow titan base hex and screw)	BSAP1081	25			10-12
SET ZZ BASE F-DX-NTBH 3.8 (Narrow titan base hex and screw)	BSAP1082	25			10-12
SET ZZ BASE F-DX-PCTBH 3.0 (Parallel cemented titan base hex and screw)	BSAB1085	24			10-12
SET ZZ BASE F-DX-PCTBH 3.4 (Parallel cemented titan base hex and screw)	BSAB1081	25			10-12
SET ZZ BASE F-DX-PCTBH 3.8 (Parallel cemented titan base hex and screw)	BSAB1082	25			10-12
SET ZZ BASE F-DX-PCTBH 4.5 (Parallel cemented titan base hex and screw)	BSAB1083	25			10-12
SET ZZ BASE F-DX-PCTBH 5.5 (Parallel cemented titan base hex and screw)	BSAB1084	25			10-12
SET ZZ BASE F-DX-PCTBH-G 3.0 (Parallel cemented titan base hex and screw Gold)	BSCB1085	24			10-12
SET ZZ BASE F-DX-PCTBH-G 3.4 (Parallel cemented titan base hex and screw Gold)	BSCB1081	25			10-12
SET ZZ BASE F-DX-PCTBH-G 3.8 (Parallel cemented titan base hex and screw Gold)	BSCB1082	25			10-12
SET ZZ BASE F-DX-PCTBH-G 4.5 (Parallel cemented titan base hex and screw Gold)	BSCB1083	25			10-12
SET ZZ BASE F-DX-PCTBH-G 5.5 (Parallel cemented titan base hex and screw Gold)	BSCB1084	25			10-12
SET ZZ BASE F-DX-PCTBH-K85 3.0 (Parallel cemented titan base hex K85 and screw)	BSAB1085K85	24			10-12
SET ZZ BASE F-DX-PCTBH-K85 3.4 (Parallel cemented titan base hex K85 and screw)	BSAB1081K85	25			10-12
SET ZZ BASE F-DX-PCTBH-K85 3.8 (Parallel cemented titan base hex K85 and screw)	BSAB1082K85	25			10-12
SET ZZ BASE F-DX-PCTBH-K85 4.5 (Parallel cemented titan base hex K85 and screw)	BSAB1083K85	25			10-12
SET ZZ BASE F-DX-PCTBH-K85 5.5 (Parallel cemented titan base hex K85 and screw)	BSAB1084K85	25			10-12
SET ZZ BASE F-DX-PCTBH-L10 3.0 (Parallel cemented titan base hex L10 and screw)	BSAB1095L10	24			10-12
SET ZZ BASE F-DX-PCTBH-L10 3.4 (Parallel cemented titan base hex L10 and screw)	BSAB1091L10	25			10-12
SET ZZ BASE F-DX-PCTBH-L10 3.8 (Parallel cemented titan base hex L10 and screw)	BSAB1092L10	25			10-12

DENTSPLY SIRONA® XIVE®/FRIALIT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE F-DX-PCTBH-L10 4.5 (Parallel cemented titan base hex L10 and screw)	BSAB1093L10	25			10-12
SET ZZ BASE F-DX-PCTBH-L10 5.5 (Parallel cemented titan base hex L10 and screw)	BSAB1094L10	25			10-12
SET ZZ BASE F-DX-PCTBH-L20 3.0 (Parallel cemented titan base hex L20 and screw)	BSAB1095L20	24			10-12
SET ZZ BASE F-DX-PCTBH-L20 3.4 (Parallel cemented titan base hex L20 and screw)	BSAB1091L20	25			10-12
SET ZZ BASE F-DX-PCTBH-L20 3.8 (Parallel cemented titan base hex L20 and screw)	BSAB1092L20	25			10-12
SET ZZ BASE F-DX-PCTBH-L20 4.5 (Parallel cemented titan base hex L20 and screw)	BSAB1093L20	25			10-12
SET ZZ BASE F-DX-PCTBH-L20 5.5 (Parallel cemented titan base hex L20 and screw)	BSAB1094L20	25			10-12
SET ZZ BASE F-DX-PCTBH-L30 3.0 (Parallel cemented titan base hex L30 and screw)	BSAB1095L30	24			10-12
SET ZZ BASE F-DX-PCTBH-L30 3.4 (Parallel cemented titan base hex L30 and screw)	BSAB1091L30	25			10-12
SET ZZ BASE F-DX-PCTBH-L30 3.8 (Parallel cemented titan base hex L30 and screw)	BSAB1092L30	25			10-12
SET ZZ BASE F-DX-PCTBH-L30 4.5 (Parallel cemented titan base hex L30 and screw)	BSAB1093L30	25			10-12
SET ZZ BASE F-DX-PCTBH-L30 5.5 (Parallel cemented titan base hex L30 and screw)	BSAB1094L30	25			10-12
SET ZZ BASE F-DX-PCTBH-L40 3.0 (Parallel cemented titan base hex L40 and screw)	BSAB1095L40	24			10-12
SET ZZ BASE F-DX-PCTBH-L40 3.4 (Parallel cemented titan base hex L40 and screw)	BSAB1091L40	25			10-12
SET ZZ BASE F-DX-PCTBH-L40 3.8 (Parallel cemented titan base hex L40 and screw)	BSAB1092L40	25			10-12
SET ZZ BASE F-DX-PCTBH-L40 4.5 (Parallel cemented titan base hex L40 and screw)	BSAB1093L40	25			10-12
SET ZZ BASE F-DX-PCTBH-L40 5.5 (Parallel cemented titan base hex L40 and screw)	BSAB1094L40	25			10-12
SET ZZ BASE F-DX-SA 3.0 (Scanalog and screw)	ZBDE1035	15			10-12
SET ZZ BASE F-DX-SA 3.4 (Scanalog and screw)	ZBDE1031	15			10-12
SET ZZ BASE F-DX-SA 3.8 (Scanalog and screw)	ZBDE1032	15			10-12
SET ZZ BASE F-DX-SA 4.5 (Scanalog and screw)	ZBDE1033	15			10-12
SET ZZ BASE F-DX-SA 5.5 (Scanalog and screw)	ZBDE1034	15			10-12
SET ZZ BASE F-DX-SM 3.0 (Scanmarker and screw)	ZBAD1035	24			10-12
SET ZZ BASE F-DX-SM 3.4 (Scanmarker and screw)	ZBAD1031	25			10-12
SET ZZ BASE F-DX-SM 3.8 (Scanmarker and screw)	ZBAD1032	25			10-12
SET ZZ BASE F-DX-SM 4.5 (Scanmarker and screw)	ZBAD1033	25			10-12
SET ZZ BASE F-DX-SM 5.5 (Scanmarker and screw)	ZBAD1034	25			10-12
SET ZZ BASE F-DX-WSM 3.0 (White Scanmarker and screw)	ZBCD1035	24			10-12
SET ZZ BASE F-DX-WSM 3.4 (White Scanmarker and screw)	ZBCD1031	25			10-12

All information is subject to change. Errors and omissions excepted. Version: 21.09.2023

DENTSPLY SIRONA® XIVE®/FRIALIT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE F-DX-WSM 3.8 (White Scanmarker and screw)	ZBCD1032	25			10-12
SET ZZ BASE F-DX-WSM 4.5 (White Scanmarker and screw)	ZBCD1033	25			10-12
SET ZZ BASE F-DX-WSM 5.5 (White Scanmarker and screw)	ZBCD1034	25			10-12
SET ZZ BASE F-DX-ZZMUA17 3.0 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and implar	BSAI1035	15			10-12
SET ZZ BASE F-DX-ZZMUA17 3.0 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and implar	BSAI1045	15			10-12
SET ZZ BASE F-DX-ZZMUA17 3.4 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and implar	BSAI1031	20			10-12
SET ZZ BASE F-DX-ZZMUA17 3.4 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and implar	BSAI1041	20			10-12
SET ZZ BASE F-DX-ZZMUA17 3.8 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and implar	BSAI1032	20			10-12
SET ZZ BASE F-DX-ZZMUA17 3.8 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and implar	BSAI1042	20			10-12
SET ZZ BASE F-DX-ZZMUA17 4.5 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and implar	BSAI1033	20			10-12
SET ZZ BASE F-DX-ZZMUA17 4.5 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and implar	BSAI1043	20			10-12
SET ZZ BASE F-DX-ZZMUA17-L10 3.0 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10 ai	BSAI1035L10	15			10-12
SET ZZ BASE F-DX-ZZMUA17-L10 3.0 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 L10 ai	BSAI1045L10	15			10-12
SET ZZ BASE F-DX-ZZMUA17-L10 3.4 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10 ai	BSAI1031L10	20			10-12
SET ZZ BASE F-DX-ZZMUA17-L10 3.4 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 L10 ai	BSAI1041L10	20			10-12
SET ZZ BASE F-DX-ZZMUA17-L10 3.8 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10 ai	BSAI1032L10	20			10-12
SET ZZ BASE F-DX-ZZMUA17-L10 3.8 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 L10 ai	BSAI1042L10	20			10-12
SET ZZ BASE F-DX-ZZMUA17-L10 4.5 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10 ai	BSAI1033L10	20			10-12
SET ZZ BASE F-DX-ZZMUA17-L10 4.5 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 L10 ai	BSAI1043L10	20			10-12
SET ZZ BASE F-DX-ZZMUA-L05 4.5 (Zirkonzahn Multi Unit Abutment non hex L05 and Screw)	BSAH1033L05	25	24		10-12
SET ZZ BASE F-DX-ZZMUA-L15 3.0 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH1035L15	25	20		10-12
SET ZZ BASE F-DX-ZZMUA-L15 3.4 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH1031L15	25	24		10-12
SET ZZ BASE F-DX-ZZMUA-L15 3.8 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH1032L15	25	24		10-12
SET ZZ BASE F-DX-ZZMUA-L15 4.5 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH1033L15	25	24		10-12
SET ZZ BASE F-DX-ZZMUA-L15 5.5 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH1034L15	25	24		10-12
SET ZZ BASE F-DX-ZZMUA-L25 3.0 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH1035L25	25	20		10-12
SET ZZ BASE F-DX-ZZMUA-L25 3.4 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH1031L25	25	24		10-12
SET ZZ BASE F-DX-ZZMUA-L25 3.8 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH1032L25	25	24		10-12
SET ZZ BASE F-DX-ZZMUA-L25 4.5 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH1033L25	25	24		10-12

DENTSPLY SIRONA® XIVE®/FRIALIT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE F-DX-ZZMUA-L25 5.5 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH1034L25	25	24		10-12
SET ZZ BASE F-DX-ZZMUA-L35 3.0 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH1035L35	25	20		10-12
SET ZZ BASE F-DX-ZZMUA-L35 3.4 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH1031L35	25	24		10-12
SET ZZ BASE F-DX-ZZMUA-L35 3.8 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH1032L35	25	24		10-12
SET ZZ BASE F-DX-ZZMUA-L35 4.5 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH1033L35	25	24		10-12
SET ZZ BASE F-DX-ZZMUA-L35 5.5 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH1034L35	25	24		10-12
SET ZZ BASE F-DX-ZZMUA-L45 3.0 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH1035L45	25	20		10-12
SET ZZ BASE F-DX-ZZMUA-L45 3.4 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH1031L45	25	24		10-12
SET ZZ BASE F-DX-ZZMUA-L45 3.8 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH1032L45	25	24		10-12
SET ZZ BASE F-DX-ZZMUA-L45 4.5 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH1033L45	25	24		10-12
SET ZZ BASE F-DX-ZZMUA-L45 5.5 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH1034L45	25	24		10-12
ZZ BASE F-DX-ASCSC 3.0 (Angled screw-channel Screw)	BSKE1005	20			10-12
ZZ BASE F-DX-ASCSC 3.0 Conical (Angled Screw Channel Screw Conical)	BSKE1055	20			10-12
ZZ BASE F-DX-ASCSC 3.4 (Angled screw-channel Screw)	BSKE1001	20			10-12
ZZ BASE F-DX-ASCSC 3.4 Conical (Angled Screw Channel Screw Conical)	BSKE1051	20			10-12
ZZ BASE F-DX-ASCSC 3.8 (Angled screw-channel Screw)	BSKE1002	20			10-12
ZZ BASE F-DX-ASCSC 3.8 Conical (Angled Screw Channel Screw Conical)	BSKE1052	20			10-12
ZZ BASE F-DX-ASCSC 4.5 (Angled screw-channel Screw)	BSKE1003	20			10-12
ZZ BASE F-DX-ASCSC 4.5 Conical (Angled Screw Channel Screw Conical)	BSKE1053	20			10-12
ZZ BASE F-DX-ASCSC 5.5 Conical (Angled Screw Channel Screw Conical)	BSKE1054	20			10-12
ZZ BASE F-DX-BSC 3.0 (Black Screw)	BSIE1005	24			10-12
ZZ BASE F-DX-BSC 3.4 (Black Screw)	BSIE1001	25			10-12
ZZ BASE F-DX-BSC 3.8 (Black Screw)	BSIE1002	25			10-12
ZZ BASE F-DX-BSC 4.5 (Black Screw)	BSIE1003	25			10-12
ZZ BASE F-DX-BSC 5.5 (Black Screw)	BSIE1004	25			10-12
ZZ BASE F-DX-BZSC 3.0 (Black Zirconia screw)	BSIE1035	24			10-12
ZZ BASE F-DX-BZSC 3.4 (Black Zirconia screw)	BSIE1031	25			10-12
ZZ BASE F-DX-BZSC 3.8 (Black Zirconia screw)	BSIE1032	25			10-12
ZZ BASE F-DX-BZSC 4.5 (Black Zirconia screw)	BSIE1033	25			10-12

DENTSPLY SIRONA® XIVE®/FRIALIT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE F-DX-BZSC 5.5 (Black Zirconia screw)	BSIE1034	25			10-12
ZZ BASE F-DX-ICSC 3.0 (Impression coping Screw)	BSFE1005	15			10-12
ZZ BASE F-DX-ICSC 3.4 (Impression coping Screw)	BSFE1001	15			10-12
ZZ BASE F-DX-ICSC 3.8 (Impression coping Screw)	BSFE1002	15			10-12
ZZ BASE F-DX-ICSC 4.5 (Impression coping Screw)	BSFE1003	15			10-12
ZZ BASE F-DX-ICSC 5.5 (Impression coping Screw)	BSFE1004	15			10-12
ZZ BASE F-DX-LBSC 3.0 (Lab Screw)	BSHE1005	24			10-12
ZZ BASE F-DX-LBSC 3.4 (Lab Screw)	BSHE1001	25			10-12
ZZ BASE F-DX-LBSC 3.8 (Lab Screw)	BSHE1002	25			10-12
ZZ BASE F-DX-LBSC 4.5 (Lab Screw)	BSHE1003	25			10-12
ZZ BASE F-DX-LBSC 5.5 (Lab Screw)	BSHE1004	25			10-12
ZZ BASE F-DX-LBZSC 3.0 (Lab Zirconia screw)	BSHE1035	24			10-12
ZZ BASE F-DX-LBZSC 3.4 (Lab Zirconia screw)	BSHE1031	25			10-12
ZZ BASE F-DX-LBZSC 3.8 (Lab Zirconia screw)	BSHE1032	25			10-12
ZZ BASE F-DX-LBZSC 4.5 (Lab Zirconia screw)	BSHE1033	25			10-12
ZZ BASE F-DX-LBZSC 5.5 (Lab Zirconia screw)	BSHE1034	25			10-12
ZZ BASE F-DX-MI 3.0 (Model implant)	MIAG1005				10-12
ZZ BASE F-DX-MI 3.4 (Model implant)	MIAG1001				10-12
ZZ BASE F-DX-MI 3.8 (Model implant)	MIAG1002				10-12
ZZ BASE F-DX-MI 4.5 (Model implant)	MIAG1003				10-12
ZZ BASE F-DX-MI 5.5 (Model implant)	MIAG1004				10-12
ZZ BASE F-DX-SC 3.0 (Screw)	BSAE1005	24			10-12
ZZ BASE F-DX-SC 3.4 (Screw)	BSAE1001	25			10-12
ZZ BASE F-DX-SC 3.8 (Screw)	BSAE1002	25			10-12
ZZ BASE F-DX-SC 4.5 (Screw)	BSAE1003	25			10-12
ZZ BASE F-DX-SC 5.5 (Screw)	BSAE1004	25			10-12
ZZ BASE F-DX-SC-G 3.0 (Screw Gold)	BSCE1005	24			10-12
ZZ BASE F-DX-SC-G 3.4 (Screw Gold)	BSCE1001	25			10-12
ZZ BASE F-DX-SC-G 3.8 (Screw Gold)	BSCE1002	25			10-12

All information is subject to change. Errors and omissions excepted. Version: 21.09.2023

DENTSPLY SIRONA® XIVE®/FRIALIT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE F-DX-SC-G 4.5 (Screw Gold)	BSCE1003	25			10-12
ZZ BASE F-DX-SC-G 5.5 (Screw Gold)	BSCE1004	25			10-12
ZZ BASE F-DX-WMSM-L10 3.0 (White Metal Scanmarker L10 with Screw)	ZBED1035L10	15			10-12
ZZ BASE F-DX-WMSM-L10 3.4 (White Metal Scanmarker L10 with Screw)	ZBED1031L10	15			10-12
ZZ BASE F-DX-WMSM-L10 3.8 (White Metal Scanmarker L10 with Screw)	ZBED1032L10	15			10-12
ZZ BASE F-DX-WMSM-L10 4.5 (White Metal Scanmarker L10 with Screw)	ZBED1033L10	15			10-12
ZZ BASE F-DX-WMSM-L10 5.5 (White Metal Scanmarker L10 with Screw)	ZBED1034L10	15			10-12
ZZ BASE F-DX-WMSM-L12 3.0 (White Metal Scanmarker L12 with Screw)	ZBED1035L12	15			10-12
ZZ BASE F-DX-WMSM-L12 3.4 (White Metal Scanmarker L12 with Screw)	ZBED1031L12	15			10-12
ZZ BASE F-DX-WMSM-L12 3.8 (White Metal Scanmarker L12 with Screw)	ZBED1032L12	15			10-12
ZZ BASE F-DX-WMSM-L12 4.5 (White Metal Scanmarker L12 with Screw)	ZBED1033L12	15			10-12
ZZ BASE F-DX-WMSM-L12 5.5 (White Metal Scanmarker L12 with Screw)	ZBED1034L12	15			10-12
ZZ BASE F-DX-WMSM-L15 3.0 (White Metal Scanmarker L15 with Screw)	ZBED1035L15	15			10-12
ZZ BASE F-DX-WMSM-L15 3.4 (White Metal Scanmarker L15 with Screw)	ZBED1031L15	15			10-12
ZZ BASE F-DX-WMSM-L15 3.8 (White Metal Scanmarker L15 with Screw)	ZBED1032L15	15			10-12
ZZ BASE F-DX-WMSM-L15 4.5 (White Metal Scanmarker L15 with Screw)	ZBED1033L15	15			10-12
ZZ BASE F-DX-WMSM-L15 5.5 (White Metal Scanmarker L15 with Screw)	ZBED1034L15	15			10-12
ZZ BASE F-DX-ZSC 3.0 (Zirconia screw)	BSAE1035	24			10-12
ZZ BASE F-DX-ZSC 3.4 (Zirconia screw)	BSAE1031	25			10-12
ZZ BASE F-DX-ZSC 3.8 (Zirconia screw)	BSAE1032	25			10-12
ZZ BASE F-DX-ZSC 4.5 (Zirconia screw)	BSAE1033	25			10-12
ZZ BASE F-DX-ZSC 5.5 (Zirconia screw)	BSAE1034	25			10-12
ZZ BASE F-DX-ZSC-G 3.0 (Zirconia screw Gold)	BSCE1035	24			10-12
ZZ BASE F-DX-ZSC-G 3.4 (Zirconia screw Gold)	BSCE1031	25			10-12
ZZ BASE F-DX-ZSC-G 3.8 (Zirconia screw Gold)	BSCE1032	25			10-12
ZZ BASE F-DX-ZSC-G 4.5 (Zirconia screw Gold)	BSCE1033	25			10-12
ZZ BASE F-DX-ZSC-G 5.5 (Zirconia screw Gold)	BSCE1034	25			10-12
ZZ BASE F-DX-ZZMUASC 3.0 (Zirkonzahn Multi Unit Abutment Screw)	BSEE1005	15			10-12
ZZ BASE F-DX-ZZMUASC 3.4 (Zirkonzahn Multi Unit Abutment Screw)	BSEE1001	20			10-12

DYNA DENTAL® OCTALOCK/HELIX

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE DYN-OL-CCTB 3.6-4.0 (Conical cemented titan base non hex and screw)	BSAD4131	35			10-12
SET ZZ BASE DYN-OL-CCTB-K85 3.6-4.0 (Conical cemented titan base non hex K85 and screw)	BSAD4131K85	35			10-12
SET ZZ BASE DYN-OL-ICOTH 3.6-4.0 (Impression coping open tray hex and Screw and Driver)	BSAQ4141	15			10-12
SET ZZ BASE DYN-OL-PCTBH 3.6-4.0 (Parallel cemented titan base hex and screw)	BSAB4181	35			10-12
SET ZZ BASE DYN-OL-PCTBH-K85 3.6-4.0 (Parallel cemented titan base hex K85 and screw)	BSAB4181K85	35			10-12
SET ZZ BASE DYN-OL-SA 3.6-4.0 (Scanalog and screw)	ZBDE4131	15			10-12
SET ZZ BASE DYN-OL-SM 3.6-4.0 (Scanmarker and screw)	ZBAD4131	35			10-12
ZZ BASE DYN-OL-ASCSC 3.6-4.0 (Angled screw-channel Screw)	BSKE4101	20			10-12
ZZ BASE DYN-OL-ICSC 3.6-4.0 (Impression coping Screw)	BSFE4101	15			10-12
ZZ BASE DYN-OL-SC 3.6-4.0 (Screw)	BSAE4101	35			10-12
ZZ BASE DYN-OL-WMSM-L12 3.6-4.0 (White Metal Scanmarker L12 with Screw)	ZBED4131L12	15			10-12
ZZ BASE DYN-OL-ZSC 3.6-4.0 (Zirconia screw)	BSAE4131	35			10-12

FAIRIMPLANT FAIRTWO™

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE FI-FT-ABCH L (Raw-Abutment® customizable hex and screw)	BSAF9582	20			10-12
SET ZZ BASE FI-FT-ABCH S (Raw-Abutment® customizable hex and screw)	BSAF9581	20			10-12
SET ZZ BASE FI-FT-CCTB L (Conical cemented titan base non hex and screw)	BSAD9532	20			10-12
SET ZZ BASE FI-FT-CCTB S (Conical cemented titan base non hex and screw)	BSAD9531	20			10-12
SET ZZ BASE FI-FT-CCTB-K85 L (Conical cemented titan base non hex K85 and screw)	BSAD9532K85	20			10-12
SET ZZ BASE FI-FT-CCTB-K85 S (Conical cemented titan base non hex K85 and screw)	BSAD9531K85	20			10-12
SET ZZ BASE FI-FT-PCTBH L (Parallel cemented titan base hex and screw)	BSAB9582	20			10-12
SET ZZ BASE FI-FT-PCTBH S (Parallel cemented titan base hex and screw)	BSAB9581	20			10-12
SET ZZ BASE FI-FT-PCTBH-K85 L (Parallel cemented titan base hex K85 and screw)	BSAB9582K85	20			10-12
SET ZZ BASE FI-FT-PCTBH-K85 S (Parallel cemented titan base hex K85 and screw)	BSAB9581K85	20			10-12
SET ZZ BASE FI-FT-SA L (Scananalog and screw)	ZBDE9532	15			10-12
SET ZZ BASE FI-FT-SA S (Scananalog and screw)	ZBDE9531	15			10-12
SET ZZ BASE FI-FT-SM L (Scanmarker and screw)	ZBAD9532	20			10-12
SET ZZ BASE FI-FT-SM S (Scanmarker and screw)	ZBAD9531	20			10-12
SET ZZ BASE FI-FT-WSM L (White Scanmarker and screw)	ZBCD9532	20			10-12
SET ZZ BASE FI-FT-WSM S (White Scanmarker and screw)	ZBCD9531	20			10-12
ZZ BASE FI-FT-MI L (Model implant)	MIAG9502				10-12
ZZ BASE FI-FT-MI S (Model implant)	MIAG9501				10-12
ZZ BASE FI-FT-SC L (Screw)	BSAE9502	20			10-12
ZZ BASE FI-FT-SC S (Screw)	BSAE9501	20			10-12
ZZ BASE FI-FT-ZSC L (Zirconia Screw)	BSAE9532	20			10-12
ZZ BASE FI-FT-ZSC S (Zirconia Screw)	BSAE9531	20			10-12

GC TECH. AADVA™

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE GC-AA-ABCH N (Raw-Abutment® customizable hex and screw)	BSAF10681	20			10-12
SET ZZ BASE GC-AA-ABCH R (Raw-Abutment® customizable hex and screw)	BSAF10682	20			10-12
SET ZZ BASE GC-AA-ABCH W (Raw-Abutment® customizable hex and screw)	BSAF10683	20			10-12
SET ZZ BASE GC-AA-ABCH-D14 N (Raw-Abutment® customizable hex D14 and screw)	BSAF10681D14	20			10-12
SET ZZ BASE GC-AA-ABCH-D14 R (Raw-Abutment® customizable hex D14 and screw)	BSAF10682D14	20			10-12
SET ZZ BASE GC-AA-ABCH-D14 W (Raw-Abutment® customizable hex D14 and screw)	BSAF10683D14	20			10-12
SET ZZ BASE GC-AA-CCTB N (Conical cemented titan base non hex and screw)	BSAD10631	20			10-12
SET ZZ BASE GC-AA-CCTB R (Conical cemented titan base non hex and screw)	BSAD10632	20			10-12
SET ZZ BASE GC-AA-CCTB W (Conical cemented titan base non hex and screw)	BSAD10633	20			10-12
SET ZZ BASE GC-AA-CCTB-K85 N (Conical cemented titan base non hex K85 and screw)	BSAD10631K85	20			10-12
SET ZZ BASE GC-AA-CCTB-L20 N (Conical cemented titan base non hex L20 and screw)	BSAD10631L20	20			10-12
SET ZZ BASE GC-AA-CCTB-L20 R (Conical cemented titan base non hex L20 and screw)	BSAD10632L20	20			10-12
SET ZZ BASE GC-AA-CCTB-L20 W (Conical cemented titan base non hex L20 and screw)	BSAD10633L20	20			10-12
SET ZZ BASE GC-AA-CCTB-L30 N (Conical cemented titan base non hex L30 and screw)	BSAD10631L30	20			10-12
SET ZZ BASE GC-AA-CCTB-L30 R (Conical cemented titan base non hex L30 and screw)	BSAD10632L30	20			10-12
SET ZZ BASE GC-AA-CCTB-L30 W (Conical cemented titan base non hex L30 and screw)	BSAD10633L30	20			10-12
SET ZZ BASE GC-AA-CCTB-L40 N (Conical cemented titan base non hex L40 and screw)	BSAD10631L40	20			10-12
SET ZZ BASE GC-AA-CCTB-L40 R (Conical cemented titan base non hex L40 and screw)	BSAD10632L40	20			10-12
SET ZZ BASE GC-AA-CCTB-L40 W (Conical cemented titan base non hex L40 and screw)	BSAD10633L40	20			10-12
SET ZZ BASE GC-AA-ICOTH N (Impression coping open tray hex and Screw and Driver)	BSAQ10641	15			10-12
SET ZZ BASE GC-AA-ICOTH R (Impression coping open tray hex and Screw and Driver)	BSAQ10642	15			10-12
SET ZZ BASE GC-AA-ICOTH W (Impression coping open tray hex and Screw and Driver)	BSAQ10643	15			10-12
SET ZZ BASE GC-AA-NTB N (Narrow titan base non hex and screw)	BSAP10631	20			10-12
SET ZZ BASE GC-AA-NTB R (Narrow titan base non hex and screw)	BSAP10632	20			10-12
SET ZZ BASE GC-AA-NTB W (Narrow titan base non hex and screw)	BSAP10633	20			10-12
SET ZZ BASE GC-AA-NTBH N (Narrow titan base hex and screw)	BSAP10681	20			10-12
SET ZZ BASE GC-AA-NTBH R (Narrow titan base hex and screw)	BSAP10682	20			10-12
SET ZZ BASE GC-AA-NTBH W (Narrow titan base hex and screw)	BSAP10683	20			10-12
SET ZZ BASE GC-AA-PCTBH N (Parallel cemented titan base hex and screw)	BSAB10681	20			10-12

GC TECH. AADVA™

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE GC-AA-PCTBH R (Parallel cemented titan base hex and screw)	BSAB10682	20			10-12
SET ZZ BASE GC-AA-PCTBH W (Parallel cemented titan base hex and screw)	BSAB10683	20			10-12
SET ZZ BASE GC-AA-PCTBH-K85 N (Parallel cemented titan base hex K85 and screw)	BSAB10681K85	20			10-12
SET ZZ BASE GC-AA-PCTBH-K85 R (Parallel cemented titan base hex K85 and screw)	BSAB10682K85	20			10-12
SET ZZ BASE GC-AA-PCTBH-K85 W (Parallel cemented titan base hex K85 and screw)	BSAB10683K85	20			10-12
SET ZZ BASE GC-AA-PCTBH-L20 N (Parallel cemented titan base hex L20 and screw)	BSAB10691L20	20			10-12
SET ZZ BASE GC-AA-PCTBH-L20 R (Parallel cemented titan base hex L20 and screw)	BSAB10692L20	20			10-12
SET ZZ BASE GC-AA-PCTBH-L20 W (Parallel cemented titan base hex L20 and screw)	BSAB10693L20	20			10-12
SET ZZ BASE GC-AA-PCTBH-L30 N (Parallel cemented titan base hex L30 and screw)	BSAB10691L30	20			10-12
SET ZZ BASE GC-AA-PCTBH-L30 R (Parallel cemented titan base hex L30 and screw)	BSAB10692L30	20			10-12
SET ZZ BASE GC-AA-PCTBH-L30 W (Parallel cemented titan base hex L30 and screw)	BSAB10693L30	20			10-12
SET ZZ BASE GC-AA-PCTBH-L40 N (Parallel cemented titan base hex L40 and screw)	BSAB10691L40	20			10-12
SET ZZ BASE GC-AA-PCTBH-L40 R (Parallel cemented titan base hex L40 and screw)	BSAB10692L40	20			10-12
SET ZZ BASE GC-AA-PCTBH-L40 W (Parallel cemented titan base hex L40 and screw)	BSAB10693L40	20			10-12
SET ZZ BASE GC-AA-SA N (Scananalog and screw)	ZBDE10631	15			10-12
SET ZZ BASE GC-AA-SA R (Scananalog and screw)	ZBDE10632	15			10-12
SET ZZ BASE GC-AA-SA W (Scananalog and screw)	ZBDE10633	15			10-12
SET ZZ BASE GC-AA-SM N (Scanmarker and screw)	ZBAD10631	20			10-12
SET ZZ BASE GC-AA-SM R (Scanmarker and screw)	ZBAD10632	20			10-12
SET ZZ BASE GC-AA-SM W (Scanmarker and screw)	ZBAD10633	20			10-12
ZZ BASE GC-AA-ASCSC N (Angled screw-channel Screw)	BSKE10601	20			10-12
ZZ BASE GC-AA-ASCSC R (Angled screw-channel Screw)	BSKE10602	20			10-12
ZZ BASE GC-AA-ASCSC W (Angled screw-channel Screw)	BSKE10603	20			10-12
ZZ BASE GC-AA-ICSC N (Impression coping Screw)	BSFE10601	15			10-12
ZZ BASE GC-AA-ICSC R (Impression coping Screw)	BSFE10602	15			10-12
ZZ BASE GC-AA-ICSC W (Impression coping Screw)	BSFE10603	15			10-12
ZZ BASE GC-AA-MI N (Model implant)	MIAG10601				10-12
ZZ BASE GC-AA-MI R (Model implant)	MIAG10602				10-12
ZZ BASE GC-AA-MI W (Model implant)	MIAG10603				10-12

GC TECH. AADVA™

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE GC-AA-SC N (Screw)	BSAE10601	20			10-12
ZZ BASE GC-AA-SC R (Screw)	BSAE10602	20			10-12
ZZ BASE GC-AA-SC W (Screw)	BSAE10603	20			10-12
ZZ BASE GC-AA-WMSM-L12 N (White Metal Scanmarker L12 with Screw)	ZBED10631L12	15			10-12
ZZ BASE GC-AA-WMSM-L12 R (White Metal Scanmarker L12 with Screw)	ZBED10632L12	15			10-12
ZZ BASE GC-AA-WMSM-L12 W (White Metal Scanmarker L12 with Screw)	ZBED10633L12	15			10-12
ZZ BASE GC-AA-ZSC N (Zirconia screw)	BSAE10631	20			10-12
ZZ BASE GC-AA-ZSC R (Zirconia screw)	BSAE10632	20			10-12
ZZ BASE GC-AA-ZSC W (Zirconia screw)	BSAE10633	20			10-12

GLOBAL D TEKKA® IN-KONE®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE T-IK-ABCH 3.5/4.0/4.5/5.5 (Raw-Abutment® customizable hex and screw)	BSAF9081	15			10-12
SET ZZ BASE T-IK-ABCH-D14 3.5/4.0/4.5/5.5 (Raw-Abutment® customizable hex D14 and screw)	BSAF9081D14	15			10-12
SET ZZ BASE T-IK-ASCTBH-K80 3.5/4.0/4.5/5.5 (Angled screw-channel titan base hex K80 and screw)	BSAT9081K80	15			10-12
SET ZZ BASE T-IK-ASCTB-K80 3.5/4.0/4.5/5.5 (Angled screw-channel titan base non hex K80 and screw)	BSAT9031K80	15			10-12
SET ZZ BASE T-IK-CCTB 0° (Conical cemented titan base non hex and screw)	BSAD9032	10			10-12
SET ZZ BASE T-IK-CCTB 17°/30° (Conical cemented titan base non hex and screw)	BSAD9033	10			10-12
SET ZZ BASE T-IK-CCTB 3.5/4.0/4.5/5.5 (Conical cemented titan base non hex and screw)	BSAD9031	15			10-12
SET ZZ BASE T-IK-CCTB-K85 0° (Conical cemented titan base non hex K85 and screw)	BSAD9032K85	10			10-12
SET ZZ BASE T-IK-CCTB-K85 17°/30° (Conical cemented titan base non hex K85 and screw)	BSAD9033K85	10			10-12
SET ZZ BASE T-IK-CCTB-K85 3.5/4.0/4.5/5.5 (Conical cemented titan base non hex K85 and screw)	BSAD9031K85	15			10-12
SET ZZ BASE T-IK-CCTB-L20 3.5/4.0/4.5/5.5 (Conical cemented titan base non hex L20 and screw)	BSAD9031L20	15			10-12
SET ZZ BASE T-IK-CCTB-L30 3.5/4.0/4.5/5.5 (Conical cemented titan base non hex L30 and screw)	BSAD9031L30	15			10-12
SET ZZ BASE T-IK-CCTB-L40 3.5/4.0/4.5/5.5 (Conical cemented titan base non hex L40 and screw)	BSAD9031L40	15			10-12
SET ZZ BASE T-IK-ICOT 0° (Impression coping open tray non hex and Screw and Driver)	BSAQ9062	15			10-12
SET ZZ BASE T-IK-ICOT 17°/30° (Impression coping open tray non hex and Screw and Driver)	BSAQ9063	15			10-12
SET ZZ BASE T-IK-ICOTH 3.5/4.0/4.5/5.5 (Impression coping open tray hex and Screw and Driver)	BSAQ9041	15			10-12
SET ZZ BASE T-IK-NTB 3.5/4.0/4.5/5.5 (Narrow titan base non hex and screw)	BSAP9031	15			10-12
SET ZZ BASE T-IK-NTBH 3.5/4.0/4.5/5.5 (Narrow titan base hex and screw)	BSAP9081	15			10-12
SET ZZ BASE T-IK-PCTBH 3.5/4.0/4.5/5.5 (Parallel cemented titan base hex and screw)	BSAB9081	15			10-12
SET ZZ BASE T-IK-PCTBH-K85 3.5/4.0/4.5/5.5 (Parallel cemented titan base hex K85 and screw)	BSAB9081K85	15			10-12
SET ZZ BASE T-IK-PCTBH-L20 3.5/4.0/4.5/5.5 (Parallel cemented titan base hex L20 and screw)	BSAB9091L20	15			10-12
SET ZZ BASE T-IK-PCTBH-L30 3.5/4.0/4.5/5.5 (Parallel cemented titan base hex L30 and screw)	BSAB9091L30	15			10-12
SET ZZ BASE T-IK-PCTBH-L40 3.5/4.0/4.5/5.5 (Parallel cemented titan base hex L40 and screw)	BSAB9091L40	15			10-12
SET ZZ BASE T-IK-SA 0° (Scanalog and screw)	ZBDE9032	15			10-12
SET ZZ BASE T-IK-SA 17°/30° (Scanalog and screw)	ZBDE9033	15			10-12
SET ZZ BASE T-IK-SA 3.5/4.0/4.5/5.5 (Scanalog and screw)	ZBDE9031	15			10-12
SET ZZ BASE T-IK-SM 0° (Scanmarker and screw)	ZBAD9032	10			10-12
SET ZZ BASE T-IK-SM 17°/30° (Scanmarker and screw)	ZBAD9033	10			10-12
SET ZZ BASE T-IK-SM 3.5/4.0/4.5/5.5 (Scanmarker and screw)	ZBAD9031	15			10-12

GLOBAL D TEKKA® IN-KONE®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE T-IK-WSM 0° (White Scanmarker and screw)	ZBCD9032	10			10-12
SET ZZ BASE T-IK-WSM 17°/30° (White Scanmarker and screw)	ZBCD9033	10			10-12
SET ZZ BASE T-IK-WSM 3.5/4.0/4.5/5.5 (White Scanmarker and screw)	ZBCD9031	15			10-12
SET ZZ BASE T-IK-ZZMUA-L15 3.5/4.0/4.5/5.5 (Zirkonzahn Multi Unit Abutment non hex L15 and	BSAH9031L15	25	15		10-12
SET ZZ BASE T-IK-ZZMUA-L25 3.5/4.0/4.5/5.5 (Zirkonzahn Multi Unit Abutment non hex L25 and	BSAH9031L25	25	15		10-12
SET ZZ BASE T-IK-ZZMUA-L35 3.5/4.0/4.5/5.5 (Zirkonzahn Multi Unit Abutment non hex L35 and	BSAH9031L35	25	15		10-12
SET ZZ BASE T-IK-ZZMUA-L45 3.5/4.0/4.5/5.5 (Zirkonzahn Multi Unit Abutment non hex L45 and	BSAH9031L45	25	15		10-12
ZZ BASE T-IK-ASCSC 3.5/4.0/4.5/5.5 (Angled screw-channel Screw)	BSKE9001	15			10-12
ZZ BASE T-IK-ICSC 0° (Impression coping Screw)	BSFE9002	15			10-12
ZZ BASE T-IK-ICSC 17°/30° (Impression coping Screw)	BSFE9003	15			10-12
ZZ BASE T-IK-ICSC 3.5/4.0/4.5/5.5 (Impression coping Screw)	BSFE9001	15			10-12
ZZ BASE T-IK-MI 0° (Model implant)	MIAG9002				10-12
ZZ BASE T-IK-MI 17°/30° (Model implant)	MIAG9003				10-12
ZZ BASE T-IK-MI 3.5/4.0/4.5/5.5 (Model implant)	MIAG9001				10-12
ZZ BASE T-IK-SC 0° (Screw)	BSAE9002	10			10-12
ZZ BASE T-IK-SC 17°/30° (Screw)	BSAE9003	10			10-12
ZZ BASE T-IK-SC 3.5/4.0/4.5/5.5 (Screw)	BSAE9001	15			10-12
ZZ BASE T-IK-WMSM-L12 0° (White Metal Scanmarker L12 with Screw)	ZBED9032L12	15			10-12
ZZ BASE T-IK-WMSM-L12 17°/30° (White Metal Scanmarker L12 with Screw)	ZBED9033L12	15			10-12
ZZ BASE T-IK-WMSM-L12 3.5/4.0/4.5/5.5 (White Metal Scanmarker L12 with Screw)	ZBED9031L12	15			10-12
ZZ BASE T-IK-ZSC 0° (Zirconia Screw)	BSAE9032	10			10-12
ZZ BASE T-IK-ZSC 17°/30° (Zirconia Screw)	BSAE9033	10			10-12
ZZ BASE T-IK-ZSC 3.5/4.0/4.5/5.5 (Zirconia Screw)	BSAE9031	15			10-12

IMPLANT DIRECT™ LEGACY™

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE ID-LE-ABCH 3.0 (Raw-Abutment® customizable hex and screw)	BSAF5791	30			10-12
SET ZZ BASE ID-LE-ABCH 3.5 (Raw-Abutment® customizable hex and screw)	BSAF5782	30			10-12
SET ZZ BASE ID-LE-ABCH 4.5 (Raw-Abutment® customizable hex and screw)	BSAF5783	30			10-12
SET ZZ BASE ID-LE-ABCH 5.7 (Raw-Abutment® customizable hex and screw)	BSAF5784	30			10-12
SET ZZ BASE ID-LE-ABCH-D14 3.0 (Raw-Abutment® customizable hex D14 and screw)	BSAF5791D14	30			10-12
SET ZZ BASE ID-LE-ABCH-D14 3.5 (Raw-Abutment® customizable hex D14 and screw)	BSAF5782D14	30			10-12
SET ZZ BASE ID-LE-ABCH-D14 4.5 (Raw-Abutment® customizable hex D14 and screw)	BSAF5783D14	30			10-12
SET ZZ BASE ID-LE-ABCH-D14 5.7 (Raw-Abutment® customizable hex D14 and screw)	BSAF5784D14	30			10-12
SET ZZ BASE ID-LE-CCTB 3.0 (Conical cemented titan base non hex and screw)	BSAD5731	30			10-12
SET ZZ BASE ID-LE-CCTB 3.5 (Conical cemented titan base non hex and screw)	BSAD5732	30			10-12
SET ZZ BASE ID-LE-CCTB 4.5 (Conical cemented titan base non hex and screw)	BSAD5733	30			10-12
SET ZZ BASE ID-LE-CCTB 5.7 (Conical cemented titan base non hex and screw)	BSAD5734	30			10-12
SET ZZ BASE ID-LE-CCTB-K85 3.0 (Conical cemented titan base non hex K85 and screw)	BSAD5731K85	30			10-12
SET ZZ BASE ID-LE-CCTB-K85 3.5 (Conical cemented titan base non hex K85 and screw)	BSAD5732K85	30			10-12
SET ZZ BASE ID-LE-CCTB-K85 4.5 (Conical cemented titan base non hex K85 and screw)	BSAD5733K85	30			10-12
SET ZZ BASE ID-LE-CCTB-K85 5.7 (Conical cemented titan base non hex K85 and screw)	BSAD5734K85	30			10-12
SET ZZ BASE ID-LE-ICOTH 3.0 (Impression coping open tray hex and Screw and Driver)	BSAQ5741	15			10-12
SET ZZ BASE ID-LE-ICOTH 3.5 (Impression coping open tray hex and Screw and Driver)	BSAQ5742	15			10-12
SET ZZ BASE ID-LE-ICOTH 4.5 (Impression coping open tray hex and Screw and Driver)	BSAQ5743	15			10-12
SET ZZ BASE ID-LE-ICOTH 5.7 (Impression coping open tray hex and Screw and Driver)	BSAQ5744	15			10-12
SET ZZ BASE ID-LE-NTB 3.0 (Narrow titan base non hex and screw)	BSAP5731	30			10-12
SET ZZ BASE ID-LE-NTB 3.5 (Narrow titan base non hex and screw)	BSAP5732	30			10-12
SET ZZ BASE ID-LE-NTBH 3.0 (Narrow titan base hex and screw)	BSAP5781	30			10-12
SET ZZ BASE ID-LE-NTBH 3.5 (Narrow titan base hex and screw)	BSAP5782	30			10-12
SET ZZ BASE ID-LE-PCTBH 3.0 (Parallel cemented titan base hex and screw)	BSAB5781	30			10-12
SET ZZ BASE ID-LE-PCTBH 3.5 (Parallel cemented titan base hex and screw)	BSAB5782	30			10-12
SET ZZ BASE ID-LE-PCTBH 4.5 (Parallel cemented titan base hex and screw)	BSAB5783	30			10-12
SET ZZ BASE ID-LE-PCTBH 5.7 (Parallel cemented titan base hex and screw)	BSAB5784	30			10-12
SET ZZ BASE ID-LE-PCTBH-K85 3.0 (Parallel cemented titan base hex K85 and screw)	BSAB5781K85	30			10-12

IMPLANT DIRECT™ LEGACY™

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
SET ZZ BASE ID-LE-PCTBH-K85 3.5 (Parallel cemented titan base hex K85 and screw)	BSAB5782K85	30			10-12
SET ZZ BASE ID-LE-PCTBH-K85 4.5 (Parallel cemented titan base hex K85 and screw)	BSAB5783K85	30			10-12
SET ZZ BASE ID-LE-PCTBH-K85 5.7 (Parallel cemented titan base hex K85 and screw)	BSAB5784K85	30			10-12
SET ZZ BASE ID-LE-SA 3.0 (Scanalog and screw)	ZBDE5731	15			10-12
SET ZZ BASE ID-LE-SA 3.5 (Scanalog and screw)	ZBDE5732	15			10-12
SET ZZ BASE ID-LE-SA 4.5 (Scanalog and screw)	ZBDE5733	15			10-12
SET ZZ BASE ID-LE-SA 5.7 (Scanalog and screw)	ZBDE5734	15			10-12
SET ZZ BASE ID-LE-SM 3.0 (Scanmarker and screw)	ZBAD5731	30			10-12
SET ZZ BASE ID-LE-SM 3.5 (Scanmarker and screw)	ZBAD5732	30			10-12
SET ZZ BASE ID-LE-SM 4.5 (Scanmarker and screw)	ZBAD5733	30			10-12
SET ZZ BASE ID-LE-SM 5.7 (Scanmarker and screw)	ZBAD5734	30			10-12
SET ZZ BASE ID-LE-WSM 3.0 (White Scanmarker and screw)	ZBCD5731	30			10-12
SET ZZ BASE ID-LE-WSM 3.5 (White Scanmarker and screw)	ZBCD5732	30			10-12
SET ZZ BASE ID-LE-WSM 4.5 (White Scanmarker and screw)	ZBCD5733	30			10-12
SET ZZ BASE ID-LE-WSM 5.7 (White Scanmarker and screw)	ZBCD5734	30			10-12
ZZ BASE ID-LE-ASCSC 3.0 (Angled screw-channel Screw)	BSKE5701	20			10-12
ZZ BASE ID-LE-ASCSC 3.5 (Angled screw-channel Screw)	BSKE5702	20			10-12
ZZ BASE ID-LE-ICSC 3.0 (Impression coping Screw)	BSFE5701	15			10-12
ZZ BASE ID-LE-ICSC 3.5 (Impression coping Screw)	BSFE5702	15			10-12
ZZ BASE ID-LE-ICSC 4.5 (Impression coping Screw)	BSFE5703	15			10-12
ZZ BASE ID-LE-ICSC 5.7 (Impression coping Screw)	BSFE5704	15			10-12
ZZ BASE ID-LE-SC 3.0 (Screw)	BSAE5701	30			10-12
ZZ BASE ID-LE-SC 3.5 (Screw)	BSAE5702	30			10-12
ZZ BASE ID-LE-SC 4.5 (Screw)	BSAE5703	30			10-12
ZZ BASE ID-LE-SC 5.7 (Screw)	BSAE5704	30			10-12
ZZ BASE ID-LE-WMSM-L12 3.0 (White Metal Scanmarker L12 with Screw)	ZBED5731L12	15			10-12
ZZ BASE ID-LE-WMSM-L12 3.5 (White Metal Scanmarker L12 with Screw)	ZBED5732L12	15			10-12
ZZ BASE ID-LE-WMSM-L12 4.5 (White Metal Scanmarker L12 with Screw)	ZBED5733L12	15			10-12
ZZ BASE ID-LE-WMSM-L12 5.7 (White Metal Scanmarker L12 with Screw)	ZBED5734L12	15			10-12

IMPLANT DIRECT™ LEGACY™

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
ZZ BASE ID-LE-ZSC 3.0 (Zirconia screw)	BSAE5731	30			10-12
ZZ BASE ID-LE-ZSC 3.5 (Zirconia screw)	BSAE5732	30			10-12
ZZ BASE ID-LE-ZSC 4.5 (Zirconia screw)	BSAE5733	30			10-12
ZZ BASE ID-LE-ZSC 5.7 (Zirconia screw)	BSAE5734	30			10-12

IMPLANT DIRECT™ OVERDENTURE ABUTMENT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE ID-ODA-CCTB (Conical cemented titan base non hex and screw)	BSAD4431	30			10-12
SET ZZ BASE ID-ODA-ICOT (Impression coping open tray non hex and Screw and Driver)	BSAQ4461	15			10-12
SET ZZ BASE ID-ODA-SA (Scananalog and screw)	ZBDE4431	15			10-12
SET ZZ BASE ID-ODA-SM (Scanmarker and screw)	ZBAD4431	30			10-12
SET ZZ BASE ID-ODA-WSM (White Scanmarker and screw)	ZBCD4431	30			10-12
ZZ BASE ID-ODA-ICSC (Impression coping Screw)	BSFE4401	15			10-12
ZZ BASE ID-ODA-SC (Screw)	BSAE4401	30			10-12
ZZ BASE ID-ODA-WMSM-L12 (White Metal Scanmarker L12 with Screw)	ZBED4431L12	15			10-12
ZZ BASE ID-ODA-ZSC (Zirconia Screw)	BSAE4431	30			10-12

INTRA-LOCK® INTERNATIONAL CONIC ABUTMENT SYSTEM

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
SET ZZ BASE IL-CA-CCTB 17° (Conical cemented titan base non hex and screw)	BSAD11532	30			10-12
SET ZZ BASE IL-CA-CCTB CS/30° (Conical cemented titan base non hex and screw)	BSAD11531	30			10-12
SET ZZ BASE IL-CA-ICOT 17° (Impression coping open tray non hex and Screw and Driver)	BSAQ11562	15			10-12
SET ZZ BASE IL-CA-ICOT CS/30° (Impression coping open tray non hex and Screw and Driver)	BSAQ11561	15			10-12
SET ZZ BASE IL-CA-SM 17° (Scanmarker and screw)	ZBAD11532	30			10-12
SET ZZ BASE IL-CA-SM CS/30° (Scanmarker and screw)	ZBAD11531	30			10-12
ZZ BASE IL-CA-ICSC 17° (Impression coping Screw)	BSFE11502	15			10-12
ZZ BASE IL-CA-ICSC CS/30° (Impression coping Screw)	BSFE11501	15			10-12
ZZ BASE IL-CA-MI 17° (Model implant)	MIAG11502				10-12
ZZ BASE IL-CA-MI CS/30° (Model implant)	MIAG11501				10-12
ZZ BASE IL-CA-SC 17° (Screw)	BSAE11502	30			10-12
ZZ BASE IL-CA-SC CS/30° (Screw)	BSAE11501	30			10-12
ZZ BASE IL-CA-WMSM-L12 17° (White Metal Scanmarker L12 with Screw)	ZBED11532L12	15			10-12
ZZ BASE IL-CA-WMSM-L12 CS/30° (White Metal Scanmarker L12 with Screw)	ZBED11531L12	15			10-12
ZZ BASE IL-CA-ZSC 17° (Zirconia Screw)	BSAE11532	30			10-12
ZZ BASE IL-CA-ZSC CS/30° (Zirconia Screw)	BSAE11531	30			10-12

INTRA-LOCK® INTERNATIONAL FLATONE®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
SET ZZ BASE IL-FO-CCTB (Conical cemented titan base non hex and screw)	BSAD11432	30			10-12
SET ZZ BASE IL-FO-CCTB N (Conical cemented titan base non hex and screw)	BSAD11431	30			10-12
SET ZZ BASE IL-FO-ICOT (Impression coping open tray non hex and Screw and Driver)	BSAQ11462	15			10-12
SET ZZ BASE IL-FO-ICOT N (Impression coping open tray non hex and Screw and Driver)	BSAQ11461	15			10-12
SET ZZ BASE IL-FO-SM (Scanmarker and screw)	ZBAD11432	30			10-12
SET ZZ BASE IL-FO-SM N (Scanmarker and screw)	ZBAD11431	30			10-12
ZZ BASE IL-FO-ICSC (Impression coping Screw)	BSFE11402	15			10-12
ZZ BASE IL-FO-ICSC N (Impression coping Screw)	BSFE11401	15			10-12
ZZ BASE IL-FO-SC (Screw)	BSAE11402	30			10-12
ZZ BASE IL-FO-SC N (Screw)	BSAE11401	30			10-12
ZZ BASE IL-FO-WMSM-L12 (White Metal Scanmarker L12 with Screw)	ZBED11432L12	15			10-12
ZZ BASE IL-FO-WMSM-L12 N (White Metal Scanmarker L12 with Screw)	ZBED11431L12	15			10-12
ZZ BASE IL-FO-ZSC (Zirconia Screw)	BSAE11432	30			10-12
ZZ BASE IL-FO-ZSC N (Zirconia Screw)	BSAE11431	30			10-12

INTRA-LOCK® INTERNATIONAL INTERNAL IMPLANTS

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE IL-II-ABCH Sq (Raw-Abutment® customizable hex and screw)	BSAF8191	30			10-12
SET ZZ BASE IL-II-ABCH St (Raw-Abutment® customizable hex and screw)	BSAF8192	30			10-12
SET ZZ BASE IL-II-ABCH Wd (Raw-Abutment® customizable hex and screw)	BSAF8193	30			10-12
SET ZZ BASE IL-II-CCTB Sq (Conical cemented titan base non hex and screw)	BSAD8131	30			10-12
SET ZZ BASE IL-II-CCTB St (Conical cemented titan base non hex and screw)	BSAD8132	30			10-12
SET ZZ BASE IL-II-CCTB Wd (Conical cemented titan base non hex and screw)	BSAD8133	30			10-12
SET ZZ BASE IL-II-PCTBH Sq (Parallel cemented titan base hex and screw)	BSAB8181	30			10-12
SET ZZ BASE IL-II-PCTBH St (Parallel cemented titan base hex and screw)	BSAB8182	30			10-12
SET ZZ BASE IL-II-PCTBH Wd (Parallel cemented titan base hex and screw)	BSAB8183	30			10-12
SET ZZ BASE IL-II-SA Sq (Scananalog and screw)	ZBDE8131	15			10-12
SET ZZ BASE IL-II-SA St (Scananalog and screw)	ZBDE8132	15			10-12
SET ZZ BASE IL-II-SA Wd (Scananalog and screw)	ZBDE8133	15			10-12
SET ZZ BASE IL-II-SM Sq (Scanmarker and screw)	ZBAD8131	30			10-12
SET ZZ BASE IL-II-SM St (Scanmarker and screw)	ZBAD8132	30			10-12
SET ZZ BASE IL-II-SM Wd (Scanmarker and screw)	ZBAD8133	30			10-12
SET ZZ BASE IL-II-WSM Sq (White Scanmarker and screw)	ZBCD8131	30			10-12
SET ZZ BASE IL-II-WSM St (White Scanmarker and screw)	ZBCD8132	30			10-12
SET ZZ BASE IL-II-WSM Wd (White Scanmarker and screw)	ZBCD8133	30			10-12
ZZ BASE IL-II-SC Sq (Screw)	BSAE8101	30			10-12
ZZ BASE IL-II-SC St (Screw)	BSAE8102	30			10-12
ZZ BASE IL-II-SC Wd (Screw)	BSAE8103	30			10-12
ZZ BASE IL-II-ZSC Sq (Zirconia Screw)	BSAE8131	30			10-12
ZZ BASE IL-II-ZSC St (Zirconia Screw)	BSAE8132	30			10-12
ZZ BASE IL-II-ZSC Wd (Zirconia Screw)	BSAE8133	30			10-12

INTRA-LOCK® INTERNATIONAL UNIHEX™

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE IL-UH-ABCH 3.4 (Raw-Abutment® customizable hex and screw)	BSAF7991	30			10-12
SET ZZ BASE IL-UH-ABCH 4.0 (Raw-Abutment® customizable hex and screw)	BSAF7992	30			10-12
SET ZZ BASE IL-UH-ABCH 4.75 (Raw-Abutment® customizable hex and screw)	BSAF7993	30			10-12
SET ZZ BASE IL-UH-CCTB 3.4 (Conical cemented titan base non hex and screw)	BSAD7931	30			10-12
SET ZZ BASE IL-UH-CCTB 4.0 (Conical cemented titan base non hex and screw)	BSAD7932	30			10-12
SET ZZ BASE IL-UH-CCTB 4.75 (Conical cemented titan base non hex and screw)	BSAD7933	30			10-12
SET ZZ BASE IL-UH-ICOTH 3.4 (Impression coping open tray hex and Screw and Driver)	BSAQ7941	15			10-12
SET ZZ BASE IL-UH-ICOTH 4.0 (Impression coping open tray hex and Screw and Driver)	BSAQ7942	15			10-12
SET ZZ BASE IL-UH-ICOTH 4.75 (Impression coping open tray hex and Screw and Driver)	BSAQ7943	15			10-12
SET ZZ BASE IL-UH-PCTBH 3.4 (Parallel cemented titan base hex and screw)	BSAB7981	30			10-12
SET ZZ BASE IL-UH-PCTBH 4.0 (Parallel cemented titan base hex and screw)	BSAB7982	30			10-12
SET ZZ BASE IL-UH-PCTBH 4.75 (Parallel cemented titan base hex and screw)	BSAB7983	30			10-12
SET ZZ BASE IL-UH-SM 3.4 (Scanmarker and screw)	ZBAD7931	30			10-12
SET ZZ BASE IL-UH-SM 4.0 (Scanmarker and screw)	ZBAD7932	30			10-12
SET ZZ BASE IL-UH-SM 4.75 (Scanmarker and screw)	ZBAD7933	30			10-12
SET ZZ BASE IL-UH-WSM 3.4 (White Scanmarker and screw)	ZBCD7931	30			10-12
SET ZZ BASE IL-UH-WSM 4.0 (White Scanmarker and screw)	ZBCD7932	30			10-12
SET ZZ BASE IL-UH-WSM 4.75 (White Scanmarker and screw)	ZBCD7933	30			10-12
ZZ BASE IL-UH-ICSC 3.4 (Impression coping Screw)	BSFE7901	15			10-12
ZZ BASE IL-UH-ICSC 4.0 (Impression coping Screw)	BSFE7902	15			10-12
ZZ BASE IL-UH-ICSC 4.75 (Impression coping Screw)	BSFE7903	15			10-12
ZZ BASE IL-UH-SC 3.4 (Screw)	BSAE7901	30			10-12
ZZ BASE IL-UH-SC 4.0 (Screw)	BSAE7902	30			10-12
ZZ BASE IL-UH-SC 4.75 (Screw)	BSAE7903	30			10-12
ZZ BASE IL-UH-WMSM-L12 3.4 (White Metal Scanmarker L12 with Screw)	ZBED7931L12	15			10-12
ZZ BASE IL-UH-WMSM-L12 4.0 (White Metal Scanmarker L12 with Screw)	ZBED7932L12	15			10-12
ZZ BASE IL-UH-WMSM-L12 4.75 (White Metal Scanmarker L12 with Screw)	ZBED7933L12	15			10-12
ZZ BASE IL-UH-ZSC 3.4 (Zirconia Screw)	BSAE7931	30			10-12
ZZ BASE IL-UH-ZSC 4.0 (Zirconia Screw)	BSAE7932	30			10-12

INTRA-LOCK® INTERNATIONAL UNIHEX™

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE IL-UH-ZSC 4.75 (Zirconia Screw)	BSAE7933	30			10-12

K3PRO® MINI KONUS/STANDARD

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE K3-ABCH Standard 3 Pro (Raw-Abutment® customizable hex and screw)	BSAF6382	25-30			10-12
SET ZZ BASE K3-ASCTBH-K80 Mini 2 Pro (Angled screw-channel titan base hex K80 and screw)	BSAT6381K80	15			10-12
SET ZZ BASE K3-ASCTBH-K80 Standard 3 Pro (Angled screw-channel titan base hex K80 and screw)	BSAT6382K80	20			10-12
SET ZZ BASE K3-ASCTB-K80 Mini 2 Pro (Angled screw-channel titan base non hex K80 and screw)	BSAT6331K80	15			10-12
SET ZZ BASE K3-ASCTB-K80 Standard 3 Pro (Angled screw-channel titan base non hex K80 and screw)	BSAT6332K80	20			10-12
SET ZZ BASE K3-CCTB Mini 2 Pro (Conical cemented titan base non hex and screw)	BSAD6331	15-20			10-12
SET ZZ BASE K3-CCTB Standard 3 Pro (Conical cemented titan base non hex and screw)	BSAD6332	25-30			10-12
SET ZZ BASE K3-CCTB-L15 Mini 2 Pro (Conical cemented titan base non hex L15 and screw)	BSAD6331L15	15-20			10-12
SET ZZ BASE K3-CCTB-L15 Standard 3 Pro (Conical cemented titan base non hex L15 and screw)	BSAD6332L15	25-30			10-12
SET ZZ BASE K3-ICOTH Mini 2 Pro (Impression coping open tray hex and Screw and Driver)	BSAQ6341	15			10-12
SET ZZ BASE K3-ICOTH Standard 3 Pro (Impression coping open tray hex and Screw and Driver)	BSAQ6342	15			10-12
SET ZZ BASE K3-PCTBH Mini 2 Pro (Parallel cemented titan base hex and screw)	BSAB6381	15-20			10-12
SET ZZ BASE K3-PCTBH Standard 3 Pro (Parallel cemented titan base hex and screw)	BSAB6382	25-30			10-12
SET ZZ BASE K3-PCTBH-K85 Mini 2 Pro (Parallel cemented titan base hex K85 and screw)	BSAB6381K85	15-20			10-12
SET ZZ BASE K3-PCTBH-K85 Standard 3 Pro (Parallel cemented titan base hex K85 and screw)	BSAB6382K85	25-30			10-12
SET ZZ BASE K3-PCTBH-L15 Mini 2 Pro (Parallel cemented titan base hex L15 and screw)	BSAB6391L15	15-20			10-12
SET ZZ BASE K3-PCTBH-L15 Standard 3 Pro (Parallel cemented titan base hex L15 and screw)	BSAB6392L15	25-30			10-12
SET ZZ BASE K3-SA Mini 2 Pro (Scananalog and screw)	ZBDE6331	15			10-12
SET ZZ BASE K3-SA Standard 3 Pro (Scananalog and screw)	ZBDE6332	15			10-12
SET ZZ BASE K3-SM Mini 2 Pro (Scanmarker and screw)	ZBAD6331	15-20			10-12
SET ZZ BASE K3-SM Standard 3 Pro (Scanmarker and screw)	ZBAD6332	25-30			10-12
SET ZZ BASE K3-WSM Mini 2 Pro (White Scanmarker and screw)	ZBCD6331	15-20			10-12
SET ZZ BASE K3-WSM Standard 3 Pro (White Scanmarker and screw)	ZBCD6332	25-30			10-12
ZZ BASE K3-ASCSC Mini 2 Pro (Angled screw-channel Screw)	BSKE6301	15			10-12
ZZ BASE K3-ASCSC Standard 3 Pro (Angled screw-channel Screw)	BSKE6302	15			10-12
ZZ BASE K3-ICSC Mini 2 Pro (Impression coping Screw)	BSFE6301	15			10-12
ZZ BASE K3-ICSC Standard 3 Pro (Impression coping Screw)	BSFE6302	15			10-12
ZZ BASE K3-SC Mini 2 Pro (Screw)	BSAE6301	15-20			10-12
ZZ BASE K3-SC Standard 3 Pro (Screw)	BSAE6302	25-30			10-12

K3PRO® MINI KONUS/STANDARD

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
ZZ BASE K3-WMSM-L12 Mini 2 Pro (White Metal Scanmarker L12 with Screw)	ZBED6331L12	15			10-12
ZZ BASE K3-WMSM-L12 Standard 3 Pro (White Metal Scanmarker L12 with Screw)	ZBED6332L12	15			10-12
ZZ BASE K3-ZSC Mini 2 Pro (Zirconia Screw)	BSAE6331	15-20			10-12
ZZ BASE K3-ZSC Standard 3 Pro (Zirconia Screw)	BSAE6332	25-30			10-12

KLOCKNER® ESSENTIAL® CONE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE K-EC-ABCH 3.5/4.0/4.5 (Raw-Abutment® customizable hex and screw)	BSAF2881	30			10-12
SET ZZ BASE K-EC-ABCH 4.8 (Raw-Abutment® customizable hex and screw)	BSAF2882	30			10-12
SET ZZ BASE K-EC-CCTB 3.5/4.0/4.5 (Conical cemented titan base non hex and screw)	BSAD2831	30			10-12
SET ZZ BASE K-EC-CCTB 4.8 (Conical cemented titan base non hex and screw)	BSAD2832	30			10-12
SET ZZ BASE K-EC-CCTB-K85 3.5/4.0/4.5 (Conical cemented titan base non hex K85 and screw)	BSAD2831K85	30			10-12
SET ZZ BASE K-EC-CCTB-K85 4.8 (Conical cemented titan base non hex K85 and screw)	BSAD2832K85	30			10-12
SET ZZ BASE K-EC-ICOTH 3.5/4.0/4.5 (Impression coping open tray hex and Screw and Driver)	BSAQ2841	15			10-12
SET ZZ BASE K-EC-ICOTH 4.8 (Impression coping open tray hex and Screw and Driver)	BSAQ2842	15			10-12
SET ZZ BASE K-EC-PCTBH 3.5/4.0/4.5 (Parallel cemented titan base hex and screw)	BSAB2881	30			10-12
SET ZZ BASE K-EC-PCTBH 4.8 (Parallel cemented titan base hex and screw)	BSAB2882	30			10-12
SET ZZ BASE K-EC-PCTBH-K85 3.5/4.0/4.5 (Parallel cemented titan base hex K85 and screw)	BSAB2881K85	30			10-12
SET ZZ BASE K-EC-PCTBH-K85 4.8 (Parallel cemented titan base hex K85 and screw)	BSAB2882K85	30			10-12
SET ZZ BASE K-EC-SA 3.5/4.0/4.5 (Scananalog and screw)	ZBDE2831	15			10-12
SET ZZ BASE K-EC-SA 4.8 (Scananalog and screw)	ZBDE2832	15			10-12
SET ZZ BASE K-EC-SM 3.5/4.0/4.5 (Scanmarker and screw)	ZBAD2831	30			10-12
SET ZZ BASE K-EC-SM 4.8 (Scanmarker and screw)	ZBAD2832	30			10-12
ZZ BASE K-EC-ASCSC 3.5/4.0/4.5 (Angled screw-channel Screw)	BSKE2801	20			10-12
ZZ BASE K-EC-ICSC 3.5/4.0/4.5 (Impression coping Screw)	BSFE2801	15			10-12
ZZ BASE K-EC-ICSC 4.8 (Impression coping Screw)	BSFE2802	15			10-12
ZZ BASE K-EC-SC 3.5/4.0/4.5 (Screw)	BSAE2801	30			10-12
ZZ BASE K-EC-SC 4.8 (Screw)	BSAE2802	30			10-12
ZZ BASE K-EC-WMSM-L12 3.5/4.0/4.5 (White Metal Scanmarker L12 with Screw)	ZBED2831L12	15			10-12
ZZ BASE K-EC-WMSM-L12 4.8 (White Metal Scanmarker L12 with Screw)	ZBED2832L12	15			10-12
ZZ BASE K-EC-ZSC 3.5/4.0/4.5 (Zirconia screw)	BSAE2831	30			10-12
ZZ BASE K-EC-ZSC 4.8 (Zirconia screw)	BSAE2832	30			10-12

KLOCKNER® NK2/SK2

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE K-NK2/SK2-ABCH 3.3/3.8/4.2 (Raw-Abutment® customizable hex and screw)	BSAF2991	25			10-12
SET ZZ BASE K-NK2/SK2-ABCH 4.7 (Raw-Abutment® customizable hex and screw)	BSAF2992	25			10-12
SET ZZ BASE K-NK2/SK2-CCTB 3.3/3.8/4.2 (Conical cemented titan base non hex and screw)	BSAD2931	25			10-12
SET ZZ BASE K-NK2/SK2-CCTB 4.7 (Conical cemented titan base non hex and screw)	BSAD2932	25			10-12
SET ZZ BASE K-NK2/SK2-ICOTH 3.3/3.8/4.2 (Impression coping open tray hex and Screw and Dr	BSAQ2941	15			10-12
SET ZZ BASE K-NK2/SK2-ICOTH 4.7 (Impression coping open tray hex and Screw and Driver)	BSAQ2942	15			10-12
SET ZZ BASE K-NK2/SK2-NTB 3.3/3.8/4.2 (Narrow titan base non hex and screw)	BSAP2931	25			10-12
SET ZZ BASE K-NK2/SK2-NTBH 3.3/3.8/4.2 (Narrow titan base hex and screw)	BSAP2981	25			10-12
SET ZZ BASE K-NK2/SK2-PCTBH 3.3/3.8/4.2 (Parallel cemented titan base hex and screw)	BSAB2981	25			10-12
SET ZZ BASE K-NK2/SK2-PCTBH 4.7 (Parallel cemented titan base hex and screw)	BSAB2982	25			10-12
SET ZZ BASE K-NK2/SK2-SM 3.3/3.8/4.2 (Scanmarker and screw)	ZBAD2931	25			10-12
SET ZZ BASE K-NK2/SK2-SM 4.7 (Scanmarker and screw)	ZBAD2932	25			10-12
ZZ BASE K-NK2/SK2-ICSC 3.3/3.8/4.2 (Impression coping Screw)	BSFE2901	15			10-12
ZZ BASE K-NK2/SK2-ICSC 4.7 (Impression coping Screw)	BSFE2902	15			10-12
ZZ BASE K-NK2/SK2-SC 3.3/3.8/4.2 (Screw)	BSAE2901	25			10-12
ZZ BASE K-NK2/SK2-SC 4.7 (Screw)	BSAE2902	25			10-12
ZZ BASE K-NK2/SK2-WMSM-L12 3.3/3.8/4.2 (White Metal Scanmarker L12 with Screw)	ZBED2931L12	15			10-12
ZZ BASE K-NK2/SK2-WMSM-L12 4.7 (White Metal Scanmarker L12 with Screw)	ZBED2932L12	15			10-12
ZZ BASE K-NK2/SK2-ZSC 3.3/3.8/4.2 (Zirconia screw)	BSAE2931	25			10-12
ZZ BASE K-NK2/SK2-ZSC 4.7 (Zirconia screw)	BSAE2932	25			10-12

MEDENTIKA® MEDENTIBASE®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE MD-MB-CCTB (Conical cemented titan base non hex and screw)	BSAD8631	15			10-12
SET ZZ BASE MD-MB-ICOT (Impression coping open tray non hex and Screw and Driver)	BSAQ8661	15			10-12
SET ZZ BASE MD-MB-SA (Scanalog and screw)	ZBDE8631	15			10-12
SET ZZ BASE MD-MB-SM (Scanmarker and screw)	ZBAD8631	15			10-12
SET ZZ BASE MD-MB-WSM (White Scanmarker and screw)	ZBCD8631	15			10-12
ZZ BASE MD-MB-ICSC (Impression coping Screw)	BSFE8601	15			10-12
ZZ BASE MD-MB-MI (Model implant)	MIAG8601				10-12
ZZ BASE MD-MB-SC (Screw)	BSAE8601	15			10-12
ZZ BASE MD-MB-WMSM-L12 (White Metal Scanmarker L12 with Screw)	ZBED8631L12	15			10-12
ZZ BASE MD-MB-ZSC (Zirconia Screw)	BSAE8631	15			10-12

MEDENTIKA® M-IMPLANT / MICROCONE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE MD-MI-ABCH 3.0 (Raw-Abutment® customizable hex and screw)	BSAF8781	15			10-12
SET ZZ BASE MD-MI-ABCH 3.5-5.0 (Raw-Abutment® customizable hex and screw)	BSAF8782	25			10-12
SET ZZ BASE MD-MI-CCTB 3.0 (Conical cemented titan base non hex and screw)	BSAD8731	15			10-12
SET ZZ BASE MD-MI-CCTB 3.5-5.0 (Conical cemented titan base non hex and screw)	BSAD8732	25			10-12
SET ZZ BASE MD-MI-CCTB-K85 3.0 (Conical cemented titan base non hex K85 and screw)	BSAD8731K85	15			10-12
SET ZZ BASE MD-MI-CCTB-K85 3.5-5.0 (Conical cemented titan base non hex K85 and screw)	BSAD8732K85	25			10-12
SET ZZ BASE MD-MI-ICOTH 3.0 (Impression coping open tray hex and Screw and Driver)	BSAQ8741	15			10-12
SET ZZ BASE MD-MI-ICOTH 3.5-5.0 (Impression coping open tray hex and Screw and Driver)	BSAQ8742	15			10-12
SET ZZ BASE MD-MI-NTB 3.0 (Narrow titan base non hex and screw)	BSAP8731	15			10-12
SET ZZ BASE MD-MI-NTB 3.5-5.0 (Narrow titan base non hex and screw)	BSAP8732	25			10-12
SET ZZ BASE MD-MI-NTBH 3.0 (Narrow titan base hex and screw)	BSAP8781	15			10-12
SET ZZ BASE MD-MI-NTBH 3.5-5.0 (Narrow titan base hex and screw)	BSAP8782	25			10-12
SET ZZ BASE MD-MI-PCTBH 3.0 (Parallel cemented titan base hex and screw)	BSAB8781	25			10-12
SET ZZ BASE MD-MI-PCTBH 3.5-5.0 (Parallel cemented titan base hex and screw)	BSAB8782	25			10-12
SET ZZ BASE MD-MI-PCTBH-K85 3.0 (Parallel cemented titan base hex K85 and screw)	BSAB8781K85	15			10-12
SET ZZ BASE MD-MI-PCTBH-K85 3.5-5.0 (Parallel cemented titan base hex K85 and screw)	BSAB8782K85	25			10-12
SET ZZ BASE MD-MI-SA 3.0 (Scananalog and screw)	ZBDE8731	15			10-12
SET ZZ BASE MD-MI-SA 3.5-5.0 (Scananalog and screw)	ZBDE8732	15			10-12
SET ZZ BASE MD-MI-SM 3.0 (Scanmarker and screw)	ZBAD8731	15			10-12
SET ZZ BASE MD-MI-SM 3.5-5.0 (Scanmarker and screw)	ZBAD8732	25			10-12
SET ZZ BASE MD-MI-WSM 3.0 (White Scanmarker and screw)	ZBCD8731	15			10-12
SET ZZ BASE MD-MI-WSM 3.5-5.0 (White Scanmarker and screw)	ZBCD8732	25			10-12
SET ZZ BASE MD-MI-ZZMUA-L15 3.0 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH8731L15	25	15		10-12
SET ZZ BASE MD-MI-ZZMUA-L15 3.5-5.0 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH8732L15	25	25		10-12
SET ZZ BASE MD-MI-ZZMUA-L25 3.0 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH8731L25	25	15		10-12
SET ZZ BASE MD-MI-ZZMUA-L25 3.5-5.0 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH8732L25	25	25		10-12
SET ZZ BASE MD-MI-ZZMUA-L35 3.0 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH8731L35	25	15		10-12
SET ZZ BASE MD-MI-ZZMUA-L35 3.5-5.0 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH8732L35	25	25		10-12
SET ZZ BASE MD-MI-ZZMUA-L45 3.0 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH8731L45	25	15		10-12

MEDENTIKA® M-IMPLANT / MICROCONE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
SET ZZ BASE MD-MI-ZZMUA-L45 3.5-5.0 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH8732L45	25	25		10-12
ZZ BASE MD-MI-ASCSC 3.0 (Angled screw-channel Screw)	BSKE8701	15			10-12
ZZ BASE MD-MI-ASCSC 3.5-5.0 (Angled screw-channel Screw)	BSKE8702	20			10-12
ZZ BASE MD-MI-ICSC 3.0 (Impression coping Screw)	BSFE8701	15			10-12
ZZ BASE MD-MI-ICSC 3.5-5.0 (Impression coping Screw)	BSFE8702	15			10-12
ZZ BASE MD-MI-MI 3.0 (Model implant)	MIAG8701				10-12
ZZ BASE MD-MI-MI 3.5-5.0 (Model implant)	MIAG8702				10-12
ZZ BASE MD-MI-SC 3.0 (Screw)	BSAE8701	15			10-12
ZZ BASE MD-MI-SC 3.5-5.0 (Screw)	BSAE8702	25			10-12
ZZ BASE MD-MI-WMSM-L12 3.0 (White Metal Scanmarker L12 with Screw)	ZBED8731L12	15			10-12
ZZ BASE MD-MI-WMSM-L12 3.5-5.0 (White Metal Scanmarker L12 with Screw)	ZBED8732L12	15			10-12
ZZ BASE MD-MI-ZSC 3.0 (Zirconia Screw)	BSAE8731	15			10-12
ZZ BASE MD-MI-ZSC 3.5-5.0 (Zirconia Screw)	BSAE8732	25			10-12

MEDENTIS MEDICAL ICX®-TEMLANT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE I-T-SM 3.45-4.80 (Scanmarker and screw)	ZBAD1931	25			10-12
SET ZZ BASE I-T-ABCH 3.45-4.8 (Raw-Abutment® customizable hex and screw)	BSAF1981	25			10-12
SET ZZ BASE I-T-ABCH-D14 3.45-4.8 (Raw-Abutment® customizable hex D14 and screw)	BSAF1981D14	25			10-12
SET ZZ BASE I-T-ASCTBH-K80 3.45-4.8 (Angled screw-channel titan base hex K80 and screw)	BSAT1981K80	20			10-12
SET ZZ BASE I-T-ASCTB-K80 3.45-4.8 (Angled screw-channel titan base non hex K80 and screw)	BSAT1931K80	20			10-12
SET ZZ BASE I-T-CCTB 3.45-4.8 (Conical cemented titan base non hex and screw)	BSAD1931	25			10-12
SET ZZ BASE I-T-CCTB-K85 3.45-4.80 (Conical cemented titan base non hex K85 and screw)	BSAD1931K85	25			10-12
SET ZZ BASE I-T-CCTB-L20 3.45-4.8 (Conical cemented titan base non hex L20 and screw)	BSAD1931L20	25			10-12
SET ZZ BASE I-T-CCTB-L30 3.45-4.8 (Conical cemented titan base non hex L30 and screw)	BSAD1931L30	25			10-12
SET ZZ BASE I-T-CCTB-L40 3.45-4.8 (Conical cemented titan base non hex L40 and screw)	BSAD1931L40	25			10-12
SET ZZ BASE I-T-ICOT 3.45-4.8 (Impression coping open tray non hex and Screw and Driver)	BSAQ1961	15			10-12
SET ZZ BASE I-T-ICOTH 3.45-4.80 (Impression coping open tray hex and Screw and Driver)	BSAQ1941	15			10-12
SET ZZ BASE I-T-NTB 3.45-4.8 (Narrow titan base non hex and screw)	BSAP1931	25			10-12
SET ZZ BASE I-T-NTBH 3.45-4.8 (Narrow titan base hex and screw)	BSAP1981	25			10-12
SET ZZ BASE I-T-PCTBH 3.45-4.80 (Parallel cemented titan base hex and screw)	BSAB1981	25			10-12
SET ZZ BASE I-T-PCTBH-K85 3.45-4.80 (Parallel cemented titan base hex K85 and screw)	BSAB1981K85	25			10-12
SET ZZ BASE I-T-PCTBH-L20 3.45-4.8 (Parallel cemented titan base hex L20 and screw)	BSAB1981L20	25			10-12
SET ZZ BASE I-T-PCTBH-L30 3.45-4.8 (Parallel cemented titan base hex L30 and screw)	BSAB1981L30	25			10-12
SET ZZ BASE I-T-PCTBH-L40 3.45-4.8 (Parallel cemented titan base hex L40 and screw)	BSAB1981L40	25			10-12
SET ZZ BASE I-T-SA 3.45-4.8 (Scananalog and screw)	ZBDE1931	15			10-12
SET ZZ BASE I-T-ZZMUA-L15 3.45-4.8 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH1931L15	25	25		10-12
SET ZZ BASE I-T-ZZMUA-L25 3.45-4.8 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH1931L25	25	25		10-12
SET ZZ BASE I-T-ZZMUA-L35 3.45-4.8 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH1931L35	25	25		10-12
SET ZZ BASE I-T-ZZMUA-L45 3.45-4.8 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH1931L45	25	25		10-12
ZZ BASE I-T-ICSC 3.45-4.80 (Impression coping Screw)	BSFE1901	15			10-12
ZZ BASE I-T-ASCSC 3.45-4.8 (Angled screw-channel Screw)	BSKE1901	20			10-12
ZZ BASE I-T-LBSC 3.45-4.80 (Lab Screw)	BSHE1901	25			10-12
ZZ BASE I-T-MI 3.45-4.80 (Model implant)	MIAG1901				10-12
ZZ BASE I-T-SC 3.45-4.80 (Screw)	BSAE1901	25			10-12

MEDENTIS MEDICAL ICX®-TEMPLANT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE I-T-WMSM-L10 3.45-4.8 (White Metal Scanmarker L10 with Screw)	ZBED1931L10	15			10-12
ZZ BASE I-T-WMSM-L12 3.45-4.8 (White Metal Scanmarker L12 with Screw)	ZBED1931L12	15			10-12
ZZ BASE I-T-WMSM-L15 3.45-4.8 (White Metal Scanmarker L15 with Screw)	ZBED1931L15	15			10-12
ZZ BASE I-T-ZSC 3.45-4.80 (Zirconia screw)	BSAE1931	25			10-12
ZZ BASE I-T-ZSC-G 3.45-4.8 (Zirconia Screw Gold)	BSCE1931	25			10-12

MEDENTIS MEDICAL ICX-MULTI®-KONZEPT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE I-MK-CCTB (Conical cemented titan base non hex and screw)	BSAD10031	27			10-12
SET ZZ BASE I-MK-ICOT (Impression coping open tray non hex and Screw and Driver)	BSAQ10061	15			10-12
SET ZZ BASE I-MK-SA (Scananalog and screw)	ZBDE10031	15			10-12
SET ZZ BASE I-MK-SM (Scanmarker and screw)	ZBAD10031	27			10-12
SET ZZ BASE I-MK-WSM (White Scanmarker and screw)	ZBCD10031	27			10-12
ZZ BASE I-MK-ICSC (Impression coping Screw)	BSFE10001	15			10-12
ZZ BASE I-MK-MI (Model implant)	MIAG10001				10-12
ZZ BASE I-MK-SC (Screw)	BSAE10001	27			10-12
ZZ BASE I-MK-WMSM-L12 (White Metal Scanmarker L12 with Screw)	ZBED10031L12	15			10-12
ZZ BASE I-MK-ZSC (Zirconia screw)	BSAE10031	27			10-12

MEGAGEN ANYONE® INTERNAL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE ME-AO-ABCH R/W (Raw-Abutment® customizable hex and screw)	BSAF10382	35			10-12
SET ZZ BASE ME-AO-ABCH S (Raw-Abutment® customizable hex and screw)	BSAF10381	35			10-12
SET ZZ BASE ME-AO-CCTB R/W (Conical cemented titan base non hex and screw)	BSAD10332	35			10-12
SET ZZ BASE ME-AO-CCTB-K85 R/W (Conical cemented titan base non hex K85 and screw)	BSAD10332K85	35			10-12
SET ZZ BASE ME-AO-CCTB-K85 S (Conical cemented titan base non hex K85 and screw)	BSAD10331K85	35			10-12
SET ZZ BASE ME-AO-CCTB-L15 S (Conical cemented titan base non hex L15 and screw)	BSAD10331L15	35			10-12
SET ZZ BASE ME-AO-ICOTH R/W (Impression coping open tray hex and Screw and Driver)	BSAQ10342	15			10-12
SET ZZ BASE ME-AO-ICOTH S (Impression coping open tray hex and Screw and Driver)	BSAQ10341	15			10-12
SET ZZ BASE ME-AO-NTB R/W (Narrow titan base non hex and screw)	BSAP10332	35			10-12
SET ZZ BASE ME-AO-NTB S (Narrow titan base non hex and screw)	BSAP10331	35			10-12
SET ZZ BASE ME-AO-NTBH R/W (Narrow titan base hex and screw)	BSAP10382	35			10-12
SET ZZ BASE ME-AO-NTBH S (Narrow titan base hex and screw)	BSAP10381	35			10-12
SET ZZ BASE ME-AO-PCTBH R/W (Parallel cemented titan base hex and screw)	BSAB10382	35			10-12
SET ZZ BASE ME-AO-PCTBH-K85 R/W (Parallel cemented titan base hex K85 and screw)	BSAB10382K85	35			10-12
SET ZZ BASE ME-AO-PCTBH-K85 S (Parallel cemented titan base hex K85 and screw)	BSAB10381K85	35			10-12
SET ZZ BASE ME-AO-PCTBH-L15 S (Parallel cemented titan base hex L15 and screw)	BSAB10391L15	35			10-12
SET ZZ BASE ME-AO-SA R/W (Scananalog and screw)	ZBDE10332	15			10-12
SET ZZ BASE ME-AO-SA S (Scananalog and screw)	ZBDE10331	15			10-12
SET ZZ BASE ME-AO-SM R/W (Scanmarker and screw)	ZBAD10332	35			10-12
SET ZZ BASE ME-AO-SM S (Scanmarker and screw)	ZBAD10331	35			10-12
ZZ BASE ME-AO-ASCSC R/W (Angled screw-channel Screw)	BSKE10302	20			10-12
ZZ BASE ME-AO-ASCSC S (Angled screw-channel Screw)	BSKE10301	20			10-12
ZZ BASE ME-AO-ICSC R/W (Impression coping Screw)	BSFE10302	15			10-12
ZZ BASE ME-AO-ICSC S (Impression coping Screw)	BSFE10301	15			10-12
ZZ BASE ME-AO-MI R/W (Model implant)	MIAG10302				10-12
ZZ BASE ME-AO-MI S (Model implant)	MIAG10301				10-12
ZZ BASE ME-AO-SC R/W (Screw)	BSAE10302	35			10-12
ZZ BASE ME-AO-SC S (Screw)	BSAE10301	35			10-12
ZZ BASE ME-AO-WMSM-L12 R/W (White Metal Scanmarker L12 with Screw)	ZBED10332L12	15			10-12

MEGAGEN ANYONE® INTERNAL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
ZZ BASE ME-AO-WMSM-L12 S (White Metal Scanmarker L12 with Screw)	ZBED10331L12	15			10-12
ZZ BASE ME-AO-ZSC R/W (Zirconia screw)	BSAE10332	35			10-12
ZZ BASE ME-AO-ZSC S (Zirconia screw)	BSAE10331	35			10-12

MEGAGEN ANYRIDGE®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE ME-AR-ABCH 3.5-8.0 (Raw-Abutment® customizable hex and screw)	BSAF6881	30			10-12
SET ZZ BASE ME-AR-ABCH-D14 3.5-8.0 (Raw-Abutment® customizable hex D14 and screw)	BSAF6881D14	30			10-12
SET ZZ BASE ME-AR-ASCTBH-K80 3.5-8.0 (Angled screw-channel titan base hex K80 and screw)	BSAT6881K80	20			10-12
SET ZZ BASE ME-AR-ASCTB-K80 3.5-8.0 (Angled screw-channel titan base non hex K80 and scr	BSAT6831K80	20			10-12
SET ZZ BASE ME-AR-CCTB 3.5-8.0 (Conical cemented titan base non hex and screw)	BSAD6831	30			10-12
SET ZZ BASE ME-AR-CCTB-K85 3.5-8.0 (Conical cemented titan base non hex K85 and screw)	BSAD6831K85	30			10-12
SET ZZ BASE ME-AR-CCTB-L20 3.5-8.0 (Conical cemented titan base non hex L20 and screw)	BSAD6831L20	30			10-12
SET ZZ BASE ME-AR-CCTB-L30 3.5-8.0 (Conical cemented titan base non hex L30 and screw)	BSAD6831L30	30			10-12
SET ZZ BASE ME-AR-CCTB-L40 3.5-8.0 (Conical cemented titan base non hex L40 and screw)	BSAD6831L40	30			10-12
SET ZZ BASE ME-AR-ICOT 3.5-8.0 (Impression coping open tray non hex and Screw and Driver)	BSAQ6861	15			10-12
SET ZZ BASE ME-AR-ICOTH 3.5-8.0 (Impression coping open tray hex and Screw and Driver)	BSAQ6841	15			10-12
SET ZZ BASE ME-AR-NTB 3.5-8.0 (Narrow titan base non hex and screw)	BSAP6831	30			10-12
SET ZZ BASE ME-AR-NTBH 3.5-8.0 (Narrow titan base hex and screw)	BSAP6881	30			10-12
SET ZZ BASE ME-AR-PCTBH 3.5-8.0 (Parallel cemented titan base hex and screw)	BSAB6881	30			10-12
SET ZZ BASE ME-AR-PCTBH-K85 3.5-8.0 (Parallel cemented titan base hex K85 and screw)	BSAB6881K85	30			10-12
SET ZZ BASE ME-AR-PCTBH-L20 3.5-8.0 (Parallel cemented titan base hex L20 and screw)	BSAB6891L20	30			10-12
SET ZZ BASE ME-AR-PCTBH-L30 3.5-8.0 (Parallel cemented titan base hex L30 and screw)	BSAB6891L30	30			10-12
SET ZZ BASE ME-AR-PCTBH-L40 3.5-8.0 (Parallel cemented titan base hex L40 and screw)	BSAB6891L40	30			10-12
SET ZZ BASE ME-AR-SA 3.5-8.0 (Scananalog and screw)	ZBDE6831	15			10-12
SET ZZ BASE ME-AR-SM 3.5-8.0 (Scanmarker and screw)	ZBAD6831	30			10-12
SET ZZ BASE ME-AR-ZZMUA-L15 3.5-8.0 (Zirkonzahn Multi Unit Abutment non hex L15 and Scre	BSAH6831L15	25	25-35		10-12
SET ZZ BASE ME-AR-ZZMUA-L25 3.5-8.0 (Zirkonzahn Multi Unit Abutment non hex L25 and Scre	BSAH6831L25	25	25-35		10-12
SET ZZ BASE ME-AR-ZZMUA-L35 3.5-8.0 (Zirkonzahn Multi Unit Abutment non hex L35 and Scre	BSAH6831L35	25	25-35		10-12
SET ZZ BASE ME-AR-ZZMUA-L45 3.5-8.0 (Zirkonzahn Multi Unit Abutment non hex L45 and Scre	BSAH6831L45	25	25-35		10-12
ZZ BASE ME-AR-ASCSC 3.5-8.0 (Angled screw-channel Screw)	BSKE6801	20			10-12
ZZ BASE ME-AR-ICSC 3.5-8.0 (Impression coping Screw)	BSFE6801	15			10-12
ZZ BASE ME-AR-MI 3.5-8.0 (Model implant)	MIAG6801				10-12
ZZ BASE ME-AR-SC 3.5-8.0 (Screw)	BSAE6801	30			10-12
ZZ BASE ME-AR-WMSM-L12 3.5-8.0 (White Metal Scanmarker L12 with Screw)	ZBED6831L12	15			10-12

MEGAGEN ANYRIDGE®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE ME-AR-ZSC 3.5-8.0 (Zirconia Screw)	BSAE6831	25-30			10-12

MEGAGEN EXFEEL® EXTERNAL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE ME-EX-ABCH 3.5 (Raw-Abutment® customizable hex and screw)	BSAF13591	25			10-12
SET ZZ BASE ME-EX-ABCH 4.1 (Raw-Abutment® customizable hex and screw)	BSAF13592	35			10-12
SET ZZ BASE ME-EX-ABCH 5.1 (Raw-Abutment® customizable hex and screw)	BSAF13593	35			10-12
SET ZZ BASE ME-EX-ABCH-D14 3.5 (Raw-Abutment® customizable hex D14 and screw)	BSAF13591D14	25			10-12
SET ZZ BASE ME-EX-ABCH-D14 4.1 (Raw-Abutment® customizable hex D14 and screw)	BSAF13592D14	35			10-12
SET ZZ BASE ME-EX-ABCH-D14 5.1 (Raw-Abutment® customizable hex D14 and screw)	BSAF13593D14	35			10-12
SET ZZ BASE ME-EX-CCTB 3.5 (Conical cemented titan base non hex and screw)	BSAD13531	25			10-12
SET ZZ BASE ME-EX-CCTB 4.1 (Conical cemented titan base non hex and screw)	BSAD13532	35			10-12
SET ZZ BASE ME-EX-CCTB 5.1 (Conical cemented titan base non hex and screw)	BSAD13533	35			10-12
SET ZZ BASE ME-EX-CCTB-K85 3.5 (Conical cemented titan base non hex K85 and screw)	BSAD13531K85	25			10-12
SET ZZ BASE ME-EX-CCTB-K85 4.1 (Conical cemented titan base non hex K85 and screw)	BSAD13532K85	35			10-12
SET ZZ BASE ME-EX-CCTB-K85 5.1 (Conical cemented titan base non hex K85 and screw)	BSAD13533K85	35			10-12
SET ZZ BASE ME-EX-CCTB-L10 3.5 (Conical cemented titan base non hex L10 and screw)	BSAD13531L10	25			10-12
SET ZZ BASE ME-EX-CCTB-L10 4.1 (Conical cemented titan base non hex L10 and screw)	BSAD13532L10	35			10-12
SET ZZ BASE ME-EX-CCTB-L10 5.1 (Conical cemented titan base non hex L10 and screw)	BSAD13533L10	35			10-12
SET ZZ BASE ME-EX-CCTB-L10-G 4.1 (Conical cemented titan base non hex L10 and screw Gold)	BSAD13532L10	35			10-12
SET ZZ BASE ME-EX-CCTB-L20 3.5 (Conical cemented titan base non hex L20 and screw)	BSAD13531L20	25			10-12
SET ZZ BASE ME-EX-CCTB-L20 4.1 (Conical cemented titan base non hex L20 and screw)	BSAD13532L20	35			10-12
SET ZZ BASE ME-EX-CCTB-L20 5.1 (Conical cemented titan base non hex L20 and screw)	BSAD13533L20	35			10-12
SET ZZ BASE ME-EX-CCTB-L30 3.5 (Conical cemented titan base non hex L30 and screw)	BSAD13531L30	25			10-12
SET ZZ BASE ME-EX-CCTB-L30 4.1 (Conical cemented titan base non hex L30 and screw)	BSAD13532L30	35			10-12
SET ZZ BASE ME-EX-CCTB-L30 5.1 (Conical cemented titan base non hex L30 and screw)	BSAD13533L30	35			10-12
SET ZZ BASE ME-EX-ICOTH 3.5 (Impression coping open tray hex and Screw and Driver)	BSAQ13541	15			10-12
SET ZZ BASE ME-EX-ICOTH 4.1 (Impression coping open tray hex and Screw and Driver)	BSAQ13542	15			10-12
SET ZZ BASE ME-EX-ICOTH 5.1 (Impression coping open tray hex and Screw and Driver)	BSAQ13543	15			10-12
SET ZZ BASE ME-EX-NTB 3.5 (Narrow titan base non hex and screw)	BSAP13531	25			10-12
SET ZZ BASE ME-EX-NTB 4.1 (Narrow titan base non hex and screw)	BSAP13532	35			10-12
SET ZZ BASE ME-EX-NTBH 3.5 (Narrow titan base hex and screw)	BSAP13581	25			10-12
SET ZZ BASE ME-EX-NTBH 4.1 (Narrow titan base hex and screw)	BSAP13582	35			10-12

All information is subject to change. Errors and omissions excepted. Version: 21.09.2023

MEGAGEN EXFEEL® EXTERNAL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE ME-EX-PCTBH 3.5 (Parallel cemented titan base hex and screw)	BSAB13581	25			10-12
SET ZZ BASE ME-EX-PCTBH 4.1 (Parallel cemented titan base hex and screw)	BSAB13582	35			10-12
SET ZZ BASE ME-EX-PCTBH 5.1 (Parallel cemented titan base hex and screw)	BSAB13583	35			10-12
SET ZZ BASE ME-EX-PCTBH-K85 3.5 (Parallel cemented titan base hex K85 and screw)	BSAB13581K85	25			10-12
SET ZZ BASE ME-EX-PCTBH-K85 4.1 (Parallel cemented titan base hex K85 and screw)	BSAB13582K85	35			10-12
SET ZZ BASE ME-EX-PCTBH-K85 5.1 (Parallel cemented titan base hex K85 and screw)	BSAB13583K85	35			10-12
SET ZZ BASE ME-EX-PCTBH-L10 3.5 (Parallel cemented titan base hex L10 and screw)	BSAB13591L10	25			10-12
SET ZZ BASE ME-EX-PCTBH-L10 4.1 (Parallel cemented titan base hex L10 and screw)	BSAB13592L10	35			10-12
SET ZZ BASE ME-EX-PCTBH-L10 5.1 (Parallel cemented titan base hex L10 and screw)	BSAB13593L10	35			10-12
SET ZZ BASE ME-EX-PCTBH-L10-G 4.1 (Parallel cemented titan base hex L10 and screw Gold)	BSCB13592L10	35			10-12
SET ZZ BASE ME-EX-PCTBH-L20 3.5 (Parallel cemented titan base hex L20 and screw)	BSAB13591L20	25			10-12
SET ZZ BASE ME-EX-PCTBH-L20 4.1 (Parallel cemented titan base hex L20 and screw)	BSAB13592L20	35			10-12
SET ZZ BASE ME-EX-PCTBH-L20 5.1 (Parallel cemented titan base hex L20 and screw)	BSAB13593L20	35			10-12
SET ZZ BASE ME-EX-PCTBH-L30 3.5 (Parallel cemented titan base hex L30 and screw)	BSAB13591L30	25			10-12
SET ZZ BASE ME-EX-PCTBH-L30 4.1 (Parallel cemented titan base hex L30 and screw)	BSAB13592L30	35			10-12
SET ZZ BASE ME-EX-PCTBH-L30 5.1 (Parallel cemented titan base hex L30 and screw)	BSAB13593L30	35			10-12
SET ZZ BASE ME-EX-PCTBH-L40 3.5 (Parallel cemented titan base hex L40 and screw)	BSAB13591L40	25			10-12
SET ZZ BASE ME-EX-PCTBH-L40 4.1 (Parallel cemented titan base hex L40 and screw)	BSAB13592L40	35			10-12
SET ZZ BASE ME-EX-SA 3.5 (Scananalog and screw)	ZBDE13531	15			10-12
SET ZZ BASE ME-EX-SA 4.1 (Scananalog and screw)	ZBDE13532	15			10-12
SET ZZ BASE ME-EX-SA 5.1 (Scananalog and screw)	ZBDE13533	15			10-12
SET ZZ BASE ME-EX-SM 3.5 (Scanmarker and screw)	ZBAD13531	25			10-12
SET ZZ BASE ME-EX-SM 4.1 (Scanmarker and screw)	ZBAD13532	35			10-12
SET ZZ BASE ME-EX-SM 5.1 (Scanmarker and screw)	ZBAD13533	35			10-12
SET ZZ BASE ME-EX-WSM 3.5 (White Scanmarker and screw)	ZBCD13531	25			10-12
SET ZZ BASE ME-EX-WSM 4.1 (White Scanmarker and screw)	ZBCD13532	35			10-12
SET ZZ BASE ME-EX-WSM 5.1 (White Scanmarker and screw)	ZBCD13533	35			10-12
SET ZZ BASE ME-EX-ZZMUA17 3.5 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and impl)	BSAI13531	20			10-12
SET ZZ BASE ME-EX-ZZMUA17 4.1 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and impl)	BSAI13532	30			10-12

MEGAGEN EXFEEL® EXTERNAL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE ME-EX-ZZMUA17 5.1 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and impl)	BSAI13533	35			10-12
SET ZZ BASE ME-EX-ZZMUA17-L10 3.5 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10)	BSAI13531L10	20			10-12
SET ZZ BASE ME-EX-ZZMUA17-L10 4.1 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10)	BSAI13532L10	30			10-12
SET ZZ BASE ME-EX-ZZMUA17-L10 5.1 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10)	BSAI13533L10	35			10-12
SET ZZ BASE ME-EX-ZZMUA-L15 3.5 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH13531L15	25	25		10-12
SET ZZ BASE ME-EX-ZZMUA-L15 4.1 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH13532L15	25	35		10-12
SET ZZ BASE ME-EX-ZZMUA-L15 5.1 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH13533L15	25	35		10-12
SET ZZ BASE ME-EX-ZZMUA-L25 3.5 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH13531L25	25	25		10-12
SET ZZ BASE ME-EX-ZZMUA-L25 4.1 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH13532L25	25	35		10-12
SET ZZ BASE ME-EX-ZZMUA-L25 5.1 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH13533L25	25	35		10-12
SET ZZ BASE ME-EX-ZZMUA-L35 3.5 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH13531L35	25	25		10-12
SET ZZ BASE ME-EX-ZZMUA-L35 4.1 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH13532L35	25	35		10-12
SET ZZ BASE ME-EX-ZZMUA-L35 5.1 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH13533L35	25	35		10-12
SET ZZ BASE ME-EX-ZZMUA-L45 3.5 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH13531L45	25	25		10-12
SET ZZ BASE ME-EX-ZZMUA-L45 4.1 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH13532L45	25	35		10-12
SET ZZ BASE ME-EX-ZZMUA-L45 5.1 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH13533L45	25	35		10-12
ZZ BASE ME-EX-BSC 3.5 (Black Screw)	BSIE13501	25			10-12
ZZ BASE ME-EX-BSC 4.1 (Black Screw)	BSIE13502	35			10-12
ZZ BASE ME-EX-BSC 5.1 (Black Screw)	BSIE13503	35			10-12
ZZ BASE ME-EX-BZSC 3.5 (Black Zirconia screw)	BSIE13531	25			10-12
ZZ BASE ME-EX-BZSC 4.1 (Black Zirconia screw)	BSIE13532	35			10-12
ZZ BASE ME-EX-BZSC 5.1 (Black Zirconia screw)	BSIE13533	35			10-12
ZZ BASE ME-EX-ICSC 3.5 (Impression coping Screw)	BSFE13501	15			10-12
ZZ BASE ME-EX-ICSC 4.1 (Impression coping Screw)	BSFE13502	15			10-12
ZZ BASE ME-EX-ICSC 5.1 (Impression coping Screw)	BSFE13503	15			10-12
ZZ BASE ME-EX-LBSC 3.5 (Lab Screw)	BSHE13501	25			10-12
ZZ BASE ME-EX-LBSC 4.1 (Lab Screw)	BSHE13502	35			10-12
ZZ BASE ME-EX-LBSC 5.1 (Lab Screw)	BSHE13503	35			10-12
ZZ BASE ME-EX-LBZSC 3.5 (Lab Zirconia screw)	BSHE13531	25			10-12

MEGAGEN EXFEEL® EXTERNAL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
ZZ BASE ME-EX-LBZSC 4.1 (Lab Zirconia screw)	BSHE13532	35			10-12
ZZ BASE ME-EX-LBZSC 5.1 (Lab Zirconia screw)	BSHE13533	35			10-12
ZZ BASE ME-EX-MI 3.5 (Model implant)	MIAG13501				10-12
ZZ BASE ME-EX-MI 4.1 (Model implant)	MIAG13502				10-12
ZZ BASE ME-EX-MI 5.1 (Model implant)	MIAG13503				10-12
ZZ BASE ME-EX-SC 3.5 (Screw)	BSAE13501	25			10-12
ZZ BASE ME-EX-SC 4.1 (Screw)	BSAE13502	35			10-12
ZZ BASE ME-EX-SC 5.1 (Screw)	BSAE13503	35			10-12
ZZ BASE ME-EX-SC-G 4.1 (Screw Gold)	BSCE13502	35			10-12
ZZ BASE ME-EX-WMSM-L12 3.5 (White Metal Scanmarker L12 with Screw)	ZBED13531L12	15			10-12
ZZ BASE ME-EX-WMSM-L12 4.1 (White Metal Scanmarker L12 with Screw)	ZBED13532L12	15			10-12
ZZ BASE ME-EX-WMSM-L12 5.1 (White Metal Scanmarker L12 with Screw)	ZBED13533L12	15			10-12
ZZ BASE ME-EX-ZSC 3.5 (Zirconia screw)	BSAE13531	25			10-12
ZZ BASE ME-EX-ZSC 4.1 (Zirconia screw)	BSAE13532	35			10-12
ZZ BASE ME-EX-ZSC 5.1 (Zirconia screw)	BSAE13533	35			10-12
ZZ BASE ME-EX-ZSC-G 3.5 (Zirconia screw Gold)	BSCE13531	25			10-12
ZZ BASE ME-EX-ZSC-G 4.1 (Zirconia screw Gold)	BSCE13532	35			10-12
ZZ BASE ME-EX-ZZMUASC 3.5 (Zirkonzahn Multi Unit Abutment Screw)	BSEE13501	20			10-12
ZZ BASE ME-EX-ZZMUASC 4.1 (Zirkonzahn Multi Unit Abutment Screw)	BSEE13502	30			10-12
ZZ BASE ME-EX-ZZMUASC 5.1 (Zirkonzahn Multi Unit Abutment Screw)	BSEE13503	35			10-12

MEGAGEN EXFEEL® INTERNAL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE ME-IN-ABCH 4.8 (Raw-Abutment® customizable hex and screw)	BSAF13681	35			10-12
SET ZZ BASE ME-IN-ABCH-D14 4.8 (Raw-Abutment® customizable hex D14 and screw)	BSAF13681D14	35			10-12
SET ZZ BASE ME-IN-CCTB 4.8 (Conical cemented titan base non hex and screw)	BSAD13631	35			10-12
SET ZZ BASE ME-IN-CCTB-K85 4.8 (Conical cemented titan base non hex K85 and screw)	BSAD13631K85	35			10-12
SET ZZ BASE ME-IN-CCTB-L15 4.8 (Conical cemented titan base non hex L15 and screw)	BSAD13631L15	35			10-12
SET ZZ BASE ME-IN-PCTBH 4.8 (Parallel cemented titan base hex and screw)	BSAB13671	35			10-12
SET ZZ BASE ME-IN-PCTBH 4.8 (Parallel cemented titan base hex and screw)	BSAB13681	35			10-12
SET ZZ BASE ME-IN-PCTBH-K85 4.8 (Parallel cemented titan base hex K85 and screw)	BSAB13681K85	35			10-12
SET ZZ BASE ME-IN-PCTBH-L15 4.8 (Parallel cemented titan base hex L15 and screw)	BSAB13691L15	35			10-12
SET ZZ BASE ME-IN-SA 4.8 (Scananalog and screw)	ZBDE13631	15			10-12
SET ZZ BASE ME-IN-SM 4.8 (Scanmarker and screw)	ZBAD13631	35			10-12
SET ZZ BASE ME-IN-WSM 4.8 (White Scanmarker and screw)	ZBCD13631	35			10-12
SET ZZ BASE ME-IN-ZZMUA17 4.8 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and imple)	BSAI13631	30			10-12
SET ZZ BASE ME-IN-ZZMUA17 4.8 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and imple)	BSAI13641	30			10-12
SET ZZ BASE ME-IN-ZZMUA-L15 4.8 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH13631L15	25	35		10-12
SET ZZ BASE ME-IN-ZZMUA-L25 4.8 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH13631L25	25	35		10-12
SET ZZ BASE ME-IN-ZZMUA-L35 4.8 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH13631L35	25	35		10-12
SET ZZ BASE ME-IN-ZZMUA-L45 4.8 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH13631L45	25	35		10-12
ZZ BASE ME-IN-BLZSC 4.8 (Black Long Zirconia Screw)	BSIE13681	35			10-12
ZZ BASE ME-IN-BSC 4.8 (Black Screw)	BSIE13601	35			10-12
ZZ BASE ME-IN-BZSC 4.8 (Black Zirconia screw)	BSIE13631	35			10-12
ZZ BASE ME-IN-LBLZSC 4.8 (Lab Long Zirconia Screw)	BSHE13681	35			10-12
ZZ BASE ME-IN-LBSC 4.8 (Lab Screw)	BSHE13601	35			10-12
ZZ BASE ME-IN-LBZSC 4.8 (Lab Zirconia screw)	BSHE13631	35			10-12
ZZ BASE ME-IN-LZSC 4.8 (Long zirconia screw)	BSAE13681	35			10-12
ZZ BASE ME-IN-LZSC-G 4.8 (Long Zirconia Screw Gold)	BSCE13681	35			10-12
ZZ BASE ME-IN-MI 4.8 (Model implant)	MIAG13601				10-12
ZZ BASE ME-IN-SC 4.8 (Screw)	BSAE13601	35			10-12
ZZ BASE ME-IN-ZSC 4.8 (Zirconia screw)	BSAE13631	35			10-12

MEGAGEN EXFEEL® INTERNAL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE ME-IN-ZSC-G 4.8 (Zirconia Screw Gold)	BSCE13631	35			10-12
ZZ BASE ME-IN-ZZMUASC 4.8 (Zirkonzahn Multi Unit Abutment Screw)	BSEE13601	30			10-12

MIS® C1

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE MS-C-ABCH Narrow (Raw-Abutment® customizable hex and screw)	BSAF5881	30			10-12
SET ZZ BASE MS-C-ABCH Standard (Raw-Abutment® customizable hex and screw)	BSAF5882	30			10-12
SET ZZ BASE MS-C-ABCH Wide (Raw-Abutment® customizable hex and screw)	BSAF5883	30			10-12
SET ZZ BASE MS-C-ABCH-D14 Narrow (Raw-Abutment® customizable hex D14 and screw)	BSAF5881D14	30			10-12
SET ZZ BASE MS-C-ABCH-D14 Standard (Raw-Abutment® customizable hex D14 and screw)	BSAF5882D14	30			10-12
SET ZZ BASE MS-C-ABCH-D14 Wide (Raw-Abutment® customizable hex D14 and screw)	BSAF5883D14	30			10-12
SET ZZ BASE MS-C-ASCTBH-K80 Narrow (Angled screw-channel titan base hex K80 and screw)	BSAT5881K80	20			10-12
SET ZZ BASE MS-C-ASCTBH-K80 Standard (Angled screw-channel titan base hex K80 and screw)	BSAT5882K80	20			10-12
SET ZZ BASE MS-C-ASCTB-K80 Narrow (Angled screw-channel titan base non hex K80 and screw)	BSAT5831K80	20			10-12
SET ZZ BASE MS-C-ASCTB-K80 Standard (Angled screw-channel titan base non hex K80 and screw)	BSAT5832K80	20			10-12
SET ZZ BASE MS-C-CCTB Narrow (Conical cemented titan base non hex and screw)	BSAD5831	30			10-12
SET ZZ BASE MS-C-CCTB Standard (Conical cemented titan base non hex and screw)	BSAD5832	30			10-12
SET ZZ BASE MS-C-CCTB Wide (Conical cemented titan base non hex and screw)	BSAD5833	30			10-12
SET ZZ BASE MS-C-CCTB-K85 Narrow (Conical cemented titan base non hex K85 and screw)	BSAD5831K85	30			10-12
SET ZZ BASE MS-C-CCTB-K85 Standard (Conical cemented titan base non hex K85 and screw)	BSAD5832K85	30			10-12
SET ZZ BASE MS-C-CCTB-K85 Wide (Conical cemented titan base non hex K85 and screw)	BSAD5833K85	30			10-12
SET ZZ BASE MS-C-CCTB-L20 Narrow (Conical cemented titan base non hex L20 and screw)	BSAD5831L20	30			10-12
SET ZZ BASE MS-C-CCTB-L20 Standard (Conical cemented titan base non hex L20 and screw)	BSAD5832L20	30			10-12
SET ZZ BASE MS-C-CCTB-L20 Wide (Conical cemented titan base non hex L20 and screw)	BSAD5833L20	30			10-12
SET ZZ BASE MS-C-CCTB-L30 Narrow (Conical cemented titan base non hex L30 and screw)	BSAD5831L30	30			10-12
SET ZZ BASE MS-C-CCTB-L30 Standard (Conical cemented titan base non hex L30 and screw)	BSAD5832L30	30			10-12
SET ZZ BASE MS-C-CCTB-L30 Wide (Conical cemented titan base non hex L30 and screw)	BSAD5833L30	30			10-12
SET ZZ BASE MS-C-CCTB-L40 Narrow (Conical cemented titan base non hex L40 and screw)	BSAD5831L40	30			10-12
SET ZZ BASE MS-C-CCTB-L40 Standard (Conical cemented titan base non hex L40 and screw)	BSAD5832L40	30			10-12
SET ZZ BASE MS-C-CCTB-L40 Wide (Conical cemented titan base non hex L40 and screw)	BSAD5833L40	30			10-12
SET ZZ BASE MS-C-ICOTH Narrow (Impression coping open tray hex and Screw and Driver)	BSAQ5841	15			10-12
SET ZZ BASE MS-C-ICOTH Standard (Impression coping open tray hex and Screw and Driver)	BSAQ5842	15			10-12
SET ZZ BASE MS-C-ICOTH Wide (Impression coping open tray hex and Screw and Driver)	BSAQ5843	15			10-12
SET ZZ BASE MS-C-NTB Narrow (Narrow titan base non hex and screw)	BSAP5831	30			10-12

All information is subject to change. Errors and omissions excepted. Version: 21.09.2023

MIS® C1

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE MS-C-NTB Standard (Narrow titan base non hex and screw)	BSAP5832	30			10-12
SET ZZ BASE MS-C-NTBH Narrow (Narrow titan base hex and screw)	BSAP5881	30			10-12
SET ZZ BASE MS-C-NTBH Standard (Narrow titan base hex and screw)	BSAP5882	30			10-12
SET ZZ BASE MS-C-PCTBH Narrow (Parallel cemented titan base hex and screw)	BSAB5881	30			10-12
SET ZZ BASE MS-C-PCTBH Standard (Parallel cemented titan base hex and screw)	BSAB5882	30			10-12
SET ZZ BASE MS-C-PCTBH Wide (Parallel cemented titan base hex and screw)	BSAB5883	30			10-12
SET ZZ BASE MS-C-PCTBH-K85 Narrow (Parallel cemented titan base hex K85 and screw)	BSAB5881K85	30			10-12
SET ZZ BASE MS-C-PCTBH-K85 Standard (Parallel cemented titan base hex K85 and screw)	BSAB5882K85	30			10-12
SET ZZ BASE MS-C-PCTBH-K85 Wide (Parallel cemented titan base hex K85 and screw)	BSAB5883K85	30			10-12
SET ZZ BASE MS-C-PCTBH-L20 Narrow (Parallel cemented titan base hex L20 and screw)	BSAB5881L20	30			10-12
SET ZZ BASE MS-C-PCTBH-L20 Standard (Parallel cemented titan base hex L20 and screw)	BSAB5882L20	30			10-12
SET ZZ BASE MS-C-PCTBH-L20 Wide (Parallel cemented titan base hex L20 and screw)	BSAB5883L20	30			10-12
SET ZZ BASE MS-C-PCTBH-L30 Narrow (Parallel cemented titan base hex L30 and screw)	BSAB5881L30	30			10-12
SET ZZ BASE MS-C-PCTBH-L30 Standard (Parallel cemented titan base hex L30 and screw)	BSAB5882L30	30			10-12
SET ZZ BASE MS-C-PCTBH-L30 Wide (Parallel cemented titan base hex L30 and screw)	BSAB5883L30	30			10-12
SET ZZ BASE MS-C-PCTBH-L40 Narrow (Parallel cemented titan base hex L40 and screw)	BSAB5881L40	30			10-12
SET ZZ BASE MS-C-PCTBH-L40 Standard (Parallel cemented titan base hex L40 and screw)	BSAB5882L40	30			10-12
SET ZZ BASE MS-C-PCTBH-L40 Wide (Parallel cemented titan base hex L40 and screw)	BSAB5883L40	30			10-12
SET ZZ BASE MS-C-SA Narrow (Scananalog and screw)	ZBDE5831	15			10-12
SET ZZ BASE MS-C-SA Standard (Scananalog and screw)	ZBDE5832	15			10-12
SET ZZ BASE MS-C-SA Wide (Scananalog and screw)	ZBDE5833	15			10-12
SET ZZ BASE MS-C-SM Narrow (Scanmarker and screw)	ZBAD5831	30			10-12
SET ZZ BASE MS-C-SM Standard (Scanmarker and screw)	ZBAD5832	30			10-12
SET ZZ BASE MS-C-SM Wide (Scanmarker and screw)	ZBAD5833	30			10-12
SET ZZ BASE MS-C-WSM Narrow (White Scanmarker and screw)	ZBCD5831	30			10-12
SET ZZ BASE MS-C-WSM Standard (White Scanmarker and screw)	ZBCD5832	30			10-12
SET ZZ BASE MS-C-WSM Wide (White Scanmarker and screw)	ZBCD5833	30			10-12
ZZ BASE MS-C-ASCSC Narrow (Angled screw-channel Screw)	BSKE5801	20			10-12
ZZ BASE MS-C-ASCSC Standard (Angled screw-channel Screw)	BSKE5802	20			10-12

MIS® C1

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
ZZ BASE MS-C-ICSC Narrow (Impression coping Screw)	BSFE5801	15			10-12
ZZ BASE MS-C-ICSC Standard (Impression coping Screw)	BSFE5802	15			10-12
ZZ BASE MS-C-ICSC Wide (Impression coping Screw)	BSFE5803	15			10-12
ZZ BASE MS-C-MI Narrow (Model implant)	MIAG5801				10-12
ZZ BASE MS-C-MI Standard (Model implant)	MIAG5802				10-12
ZZ BASE MS-C-MI Wide (Model implant)	MIAG5803				10-12
ZZ BASE MS-C-SC Narrow (Screw)	BSAE5801	30			10-12
ZZ BASE MS-C-SC Standard (Screw)	BSAE5802	30			10-12
ZZ BASE MS-C-SC Wide (Screw)	BSAE5803	30			10-12
ZZ BASE MS-C-WMSM-L12 Narrow (White Metal Scanmarker L12 with Screw)	ZBED5831L12	15			10-12
ZZ BASE MS-C-WMSM-L12 Standard (White Metal Scanmarker L12 with Screw)	ZBED5832L12	15			10-12
ZZ BASE MS-C-WMSM-L12 Wide (White Metal Scanmarker L12 with Screw)	ZBED5833L12	15			10-12
ZZ BASE MS-C-ZSC Narrow (Zirconia screw)	BSAE5831	30			10-12
ZZ BASE MS-C-ZSC Standard (Zirconia screw)	BSAE5832	30			10-12
ZZ BASE MS-C-ZSC Wide (Zirconia screw)	BSAE5833	30			10-12

MIS® MULTI UNIT ABUTMENT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE MS-MUA-SM (Scanmarker and screw)	ZBAD1731	20			10-12
SET ZZ BASE MS-MUA-CCTB (Conical cemented titan base non hex and screw)	BSAD1731	20			10-12
SET ZZ BASE MS-MUA-CCTB-K85 (Conical cemented titan base non hex K85 and screw)	BSAD1731K85	20			10-12
SET ZZ BASE MS-MUA-ICOT (Impression coping open tray non hex and Screw and Driver)	BSAQ1761	15			10-12
SET ZZ BASE MS-MUA-SA (Scananalog and screw)	ZBDE1731	15			10-12
SET ZZ BASE MS-MUA-WSM (White Scanmarker and screw)	ZBCD1731	20			10-12
ZZ BASE MS-MUA-ICSC (Impression coping Screw)	BSFE1701	15			10-12
ZZ BASE MS-MUA-SC (Screw)	BSAE1701	20			10-12
ZZ BASE MS-MUA-SC-G (Screw Gold)	BSCE1701	20			10-12
ZZ BASE MS-MUA-WMSM-L12 (White Metal Scanmarker L12 with Screw)	ZBED1731L12	15			10-12
ZZ BASE MS-MUA-ZSC (Zirconia screw)	BSAE1731	20			10-12
ZZ BASE MS-MUA-ZSC-G (Zirconia Screw Gold)	BSCE1731	20			10-12

MIS® MULTI UNIT SYSTEM

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE MS-MUS-CCTB (Conical cemented titan base non hex and screw)	BSAD5931	30			10-12
SET ZZ BASE MS-MUS-ICOT (Impression coping open tray non hex and Screw and Driver)	BSAQ5961	15			10-12
SET ZZ BASE MS-MUS-SA (Scananalog and screw)	ZBDE5931	15			10-12
SET ZZ BASE MS-MUS-SM (Scanmarker and screw)	ZBAD5931	30			10-12
SET ZZ BASE MS-MUS-WSM (White Scanmarker and screw)	ZBCD5931	30			10-12
ZZ BASE MS-MUS-ICSC (Impression coping Screw)	BSFE5901	15			10-12
ZZ BASE MS-MUS-MI (Model implant)	MIAG5901				10-12
ZZ BASE MS-MUS-SC (Screw)	BSAE5901	30			10-12
ZZ BASE MS-MUS-WMSM-L12 (White Metal Scanmarker L12 with Screw)	ZBED5931L12	15			10-12
ZZ BASE MS-MUS-ZSC (Zirconia screw)	BSAE5931	25			10-12

MIS® SEVEN

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE MS-S-SM Standard (Scanmarker and screw)	ZBAD1631	30			10-12
SET ZZ BASE MS-S-SM Wide (Scanmarker and screw)	ZBAD1632	30			10-12
SET ZZ BASE MS-S-ABCH Narrow (Raw-Abutment® customizable hex and screw)	BSAF1683	30			10-12
SET ZZ BASE MS-S-ABCH Standard (Raw-Abutment® customizable hex and screw)	BSAF1681	30			10-12
SET ZZ BASE MS-S-ABCH Wide (Raw-Abutment® customizable hex and screw)	BSAF1682	30			10-12
SET ZZ BASE MS-S-ABCH-D14 Narrow (Raw-Abutment® customizable hex D14 and screw)	BSAF1683D14	30			10-12
SET ZZ BASE MS-S-ABCH-D14 Standard (Raw-Abutment® customizable hex D14 and screw)	BSAF1681D14	30			10-12
SET ZZ BASE MS-S-ABCH-D14 Wide (Raw-Abutment® customizable hex D14 and screw)	BSAF1682D14	30			10-12
SET ZZ BASE MS-S-ASCTBH-K80 Narrow (Angled screw-channel titan base hex K80 and screw)	BSAT1683K80	20			10-12
SET ZZ BASE MS-S-ASCTBH-K80 Standard (Angled screw-channel titan base hex K80 and screw)	BSAT1681K80	20			10-12
SET ZZ BASE MS-S-ASCTBH-K80 Wide (Angled screw-channel titan base hex K80 and screw)	BSAT1682K80	20			10-12
SET ZZ BASE MS-S-ASCTB-K80 Narrow (Angled screw-channel titan base non hex K80 and screw)	BSAT1633K80	20			10-12
SET ZZ BASE MS-S-ASCTB-K80 Standard (Angled screw-channel titan base non hex K80 and screw)	BSAT1631K80	20			10-12
SET ZZ BASE MS-S-ASCTB-K80 Wide (Angled screw-channel titan base non hex K80 and screw)	BSAT1632K80	20			10-12
SET ZZ BASE MS-S-CCTB Narrow (Conical cemented titan base non hex and screw)	BSAD1633	30			10-12
SET ZZ BASE MS-S-CCTB Standard (Conical cemented titan base non hex and screw)	BSAD1631	30			10-12
SET ZZ BASE MS-S-CCTB Wide (Conical cemented titan base non hex and screw)	BSAD1632	30			10-12
SET ZZ BASE MS-S-CCTB-G Narrow (Conical cemented titan base non hex and screw Gold)	BSCD1633	30			10-12
SET ZZ BASE MS-S-CCTB-G Standard (Conical cemented titan base non hex and screw Gold)	BSCD1631	30			10-12
SET ZZ BASE MS-S-CCTB-G Wide (Conical cemented titan base non hex and screw Gold)	BSCD1632	30			10-12
SET ZZ BASE MS-S-CCTB-K85 Narrow (Conical cemented titan base non hex K85 and screw)	BSAD1633K85	30			10-12
SET ZZ BASE MS-S-CCTB-K85 Standard (Conical cemented titan base non hex K85 and screw)	BSAD1631K85	30			10-12
SET ZZ BASE MS-S-CCTB-K85 Wide (Conical cemented titan base non hex K85 and screw)	BSAD1632K85	30			10-12
SET ZZ BASE MS-S-CCTB-L10 Narrow (Conical cemented titan base non hex L10 and screw)	BSAD1633L10	30			10-12
SET ZZ BASE MS-S-CCTB-L10 Standard (Conical cemented titan base non hex L10 and screw)	BSAD1631L10	30			10-12
SET ZZ BASE MS-S-CCTB-L10 Wide (Conical cemented titan base non hex L10 and screw)	BSAD1632L10	30			10-12
SET ZZ BASE MS-S-CCTB-L10-G Narrow (Conical cemented titan base non hex L10 and screw Gold)	BSCD1633L10	30			10-12
SET ZZ BASE MS-S-CCTB-L10-G Standard (Conical cemented titan base non hex L10 and screw Gold)	BSCD1631L10	30			10-12
SET ZZ BASE MS-S-CCTB-L10-G Wide (Conical cemented titan base non hex L10 and screw Gold)	BSCD1632L10	30			10-12

MIS® SEVEN

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE MS-S-CCTB-L20 Narrow (Conical cemented titan base non hex L20 and screw)	BSAD1633L20	30			10-12
SET ZZ BASE MS-S-CCTB-L20 Standard (Conical cemented titan base non hex L20 and screw)	BSAD1631L20	30			10-12
SET ZZ BASE MS-S-CCTB-L20 Wide (Conical cemented titan base non hex L20 and screw)	BSAD1632L20	30			10-12
SET ZZ BASE MS-S-CCTB-L20-G Narrow (Conical cemented titan base non hex L20 and screw G)	BSCD1633L20	30			10-12
SET ZZ BASE MS-S-CCTB-L20-G Standard (Conical cemented titan base non hex L20 and screw)	BSCD1631L20	30			10-12
SET ZZ BASE MS-S-CCTB-L20-G Wide (Conical cemented titan base non hex L20 and screw Gol)	BSCD1632L20	30			10-12
SET ZZ BASE MS-S-CCTB-L30 Narrow (Conical cemented titan base non hex L30 and screw)	BSAD1633L30	30			10-12
SET ZZ BASE MS-S-CCTB-L30 Standard (Conical cemented titan base non hex L30 and screw)	BSAD1631L30	30			10-12
SET ZZ BASE MS-S-CCTB-L30 Wide (Conical cemented titan base non hex L30 and screw)	BSAD1632L30	30			10-12
SET ZZ BASE MS-S-CCTB-L30-G Narrow (Conical cemented titan base non hex L30 and screw G)	BSCD1633L30	30			10-12
SET ZZ BASE MS-S-CCTB-L30-G Standard (Conical cemented titan base non hex L30 and screw)	BSCD1631L30	30			10-12
SET ZZ BASE MS-S-CCTB-L30-G Wide (Conical cemented titan base non hex L30 and screw Gol)	BSCD1632L30	30			10-12
SET ZZ BASE MS-S-CCTB-L40 Narrow (Conical cemented titan base non hex L40 and screw)	BSAD1633L40	30			10-12
SET ZZ BASE MS-S-CCTB-L40 Standard (Conical cemented titan base non hex L40 and screw)	BSAD1631L40	30			10-12
SET ZZ BASE MS-S-CCTB-L40 Wide (Conical cemented titan base non hex L40 and screw)	BSAD1632L40	30			10-12
SET ZZ BASE MS-S-CCTB-L40-G Narrow (Conical cemented titan base non hex L40 and screw G)	BSCD1633L40	30			10-12
SET ZZ BASE MS-S-CCTB-L40-G Standard (Conical cemented titan base non hex L40 and screw)	BSCD1631L40	30			10-12
SET ZZ BASE MS-S-CCTB-L40-G Wide (Conical cemented titan base non hex L40 and screw Gol)	BSCD1632L40	30			10-12
SET ZZ BASE MS-S-ICOT Narrow (Impression coping open tray non hex and Screw and Driver)	BSAQ1663	15			10-12
SET ZZ BASE MS-S-ICOT Standard (Impression coping open tray non hex and Screw and Driver)	BSAQ1661	15			10-12
SET ZZ BASE MS-S-ICOT Wide (Impression coping open tray non hex and Screw and Driver)	BSAQ1662	15			10-12
SET ZZ BASE MS-S-ICOTH Narrow (Impression coping open tray hex and Screw and Driver)	BSAQ1643	15			10-12
SET ZZ BASE MS-S-ICOTH Standard (Impression coping open tray hex and Screw and Driver)	BSAQ1641	15			10-12
SET ZZ BASE MS-S-ICOTH Wide (Impression coping open tray hex and Screw and Driver)	BSAQ1642	15			10-12
SET ZZ BASE MS-S-NTB NARROW (Narrow titan base non hex and screw)	BSAP1633	30			10-12
SET ZZ BASE MS-S-NTB STANDARD (Narrow titan base non hex and screw)	BSAP1631	30			10-12
SET ZZ BASE MS-S-NTBH NARROW (Narrow titan base hex and screw)	BSAP1683	30			10-12
SET ZZ BASE MS-S-NTBH STANDARD (Narrow titan base hex and screw)	BSAP1681	30			10-12
SET ZZ BASE MS-S-PCTBH Narrow (Parallel cemented titan base hex and screw)	BSAB1683	30			10-12

MIS® SEVEN

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE MS-S-PCTBH Standard (Parallel cemented titan base hex and screw)	BSAB1681	30			10-12
SET ZZ BASE MS-S-PCTBH Wide (Parallel cemented titan base hex and screw)	BSAB1682	30			10-12
SET ZZ BASE MS-S-PCTBH-G Narrow (Parallel cemented titan base hex and screw Gold)	BSCB1683	30			10-12
SET ZZ BASE MS-S-PCTBH-G Wide (Parallel cemented titan base hex and screw Gold)	BSCB1682	30			10-12
SET ZZ BASE MS-S-PCTBH-K85 Narrow (Parallel cemented titan base hex K85 and screw)	BSAB1683K85	30			10-12
SET ZZ BASE MS-S-PCTBH-K85 Standard (Parallel cemented titan base hex K85 and screw)	BSAB1681K85	30			10-12
SET ZZ BASE MS-S-PCTBH-K85 Wide (Parallel cemented titan base hex K85 and screw)	BSAB1682K85	30			10-12
SET ZZ BASE MS-S-PCTBH-L10 Narrow (Parallel cemented titan base hex L10 and screw)	BSAB1693L10	30			10-12
SET ZZ BASE MS-S-PCTBH-L10 Standard (Parallel cemented titan base hex L10 and screw)	BSAB1691L10	30			10-12
SET ZZ BASE MS-S-PCTBH-L10 Wide (Parallel cemented titan base hex L10 and screw)	BSAB1692L10	30			10-12
SET ZZ BASE MS-S-PCTBH-L10-G Narrow (Parallel cemented titan base hex L10 and screw Gold)	BSCB1693L10	30			10-12
SET ZZ BASE MS-S-PCTBH-L10-G Standard (Parallel cemented titan base hex L10 and screw Gold)	BSCB1691L10	30			10-12
SET ZZ BASE MS-S-PCTBH-L10-G Wide (Parallel cemented titan base hex L10 and screw Gold)	BSCB1692L10	30			10-12
SET ZZ BASE MS-S-PCTBH-L20 Narrow (Parallel cemented titan base hex L20 and screw)	BSAB1693L20	30			10-12
SET ZZ BASE MS-S-PCTBH-L20 Standard (Parallel cemented titan base hex L20 and screw)	BSAB1691L20	30			10-12
SET ZZ BASE MS-S-PCTBH-L20 Wide (Parallel cemented titan base hex L20 and screw)	BSAB1692L20	30			10-12
SET ZZ BASE MS-S-PCTBH-L20-G Narrow (Parallel cemented titan base hex L20 and screw Gold)	BSCB1693L20	30			10-12
SET ZZ BASE MS-S-PCTBH-L20-G Standard (Parallel cemented titan base hex L20 and screw Gold)	BSCB1691L20	30			10-12
SET ZZ BASE MS-S-PCTBH-L20-G Wide (Parallel cemented titan base hex L20 and screw Gold)	BSCB1692L20	30			10-12
SET ZZ BASE MS-S-PCTBH-L30 Narrow (Parallel cemented titan base hex L30 and screw)	BSAB1693L30	30			10-12
SET ZZ BASE MS-S-PCTBH-L30 Standard (Parallel cemented titan base hex L30 and screw)	BSAB1691L30	30			10-12
SET ZZ BASE MS-S-PCTBH-L30 Wide (Parallel cemented titan base hex L30 and screw)	BSAB1692L30	30			10-12
SET ZZ BASE MS-S-PCTBH-L30-G Narrow (Parallel cemented titan base hex L30 and screw Gold)	BSCB1693L30	30			10-12
SET ZZ BASE MS-S-PCTBH-L30-G Standard (Parallel cemented titan base hex L30 and screw Gold)	BSCB1691L30	30			10-12
SET ZZ BASE MS-S-PCTBH-L30-G Wide (Parallel cemented titan base hex L30 and screw Gold)	BSCB1692L30	30			10-12
SET ZZ BASE MS-S-PCTBH-L40 Narrow (Parallel cemented titan base hex L40 and screw)	BSAB1693L40	30			10-12
SET ZZ BASE MS-S-PCTBH-L40 Standard (Parallel cemented titan base hex L40 and screw)	BSAB1691L40	30			10-12
SET ZZ BASE MS-S-PCTBH-L40 Wide (Parallel cemented titan base hex L40 and screw)	BSAB1692L40	30			10-12
SET ZZ BASE MS-S-PCTBH-L40-G Narrow (Parallel cemented titan base hex L40 and screw Gold)	BSCB1693L40	30			10-12

MIS® SEVEN

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE MS-S-PCTBH-L40-G Standard (Parallel cemented titan base hex L40 and screw Gold)	BSCB1691L40	30			10-12
SET ZZ BASE MS-S-PCTBH-L40-G Wide (Parallel cemented titan base hex L40 and screw Gold)	BSCB1692L40	30			10-12
SET ZZ BASE MS-S-SA Narrow (Scanalog and screw)	ZBDE1633	15			10-12
SET ZZ BASE MS-S-SA Standard (Scanalog and screw)	ZBDE1631	15			10-12
SET ZZ BASE MS-S-SA Wide (Scanalog and screw)	ZBDE1632	15			10-12
SET ZZ BASE MS-S-SM Narrow (Scanmarker and screw)	ZBAD1633	30			10-12
SET ZZ BASE MS-S-WSM Narrow (White Scanmarker and screw)	ZBCD1633	30			10-12
SET ZZ BASE MS-S-WSM Standard (White Scanmarker and screw)	ZBCD1631	30			10-12
SET ZZ BASE MS-S-WSM Wide (White Scanmarker and screw)	ZBCD1632	30			10-12
SET ZZ BASE MS-S-ZZMUA17 Narrow Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and in)	BSAI1633	20			10-12
SET ZZ BASE MS-S-ZZMUA17 Narrow Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and in)	BSAI1643	20			10-12
SET ZZ BASE MS-S-ZZMUA17 Standard Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and in)	BSAI1631	25			10-12
SET ZZ BASE MS-S-ZZMUA17 Standard Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and in)	BSAI1641	25			10-12
SET ZZ BASE MS-S-ZZMUA17 Wide Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and imp)	BSAI1632	25			10-12
SET ZZ BASE MS-S-ZZMUA17 Wide Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and imp)	BSAI1642	25			10-12
SET ZZ BASE MS-S-ZZMUA17-L10 Narrow Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L)	BSAI1633L10	20			10-12
SET ZZ BASE MS-S-ZZMUA17-L10 Narrow Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 L)	BSAI1643L10	20			10-12
SET ZZ BASE MS-S-ZZMUA17-L10 Standard Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L)	BSAI1631L10	25			10-12
SET ZZ BASE MS-S-ZZMUA17-L10 Standard Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 L)	BSAI1641L10	25			10-12
SET ZZ BASE MS-S-ZZMUA17-L10 Wide Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10)	BSAI1632L10	25			10-12
SET ZZ BASE MS-S-ZZMUA17-L10 Wide Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 L10)	BSAI1642L10	25			10-12
SET ZZ BASE MS-S-ZZMUA-L05 Wide (Zirkonzahn Multi Unit Abutment non hex L05 and Screw)	BSAH1632L05	25	30		10-12
SET ZZ BASE MS-S-ZZMUA-L15 Narrow (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH1633L15	25	30		10-12
SET ZZ BASE MS-S-ZZMUA-L15 Standard (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH1631L15	25	30		10-12
SET ZZ BASE MS-S-ZZMUA-L15 Wide (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH1632L15	25	30		10-12
SET ZZ BASE MS-S-ZZMUA-L25 Narrow (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH1633L25	25	30		10-12
SET ZZ BASE MS-S-ZZMUA-L25 Standard (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH1631L25	25	30		10-12
SET ZZ BASE MS-S-ZZMUA-L25 Wide (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH1632L25	25	30		10-12
SET ZZ BASE MS-S-ZZMUA-L35 Narrow (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH1633L35	25	30		10-12

All information is subject to change. Errors and omissions excepted. Version: 21.09.2023

MIS® SEVEN

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE MS-S-ZZMUA-L35 Standard (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH1631L35	25	30		10-12
SET ZZ BASE MS-S-ZZMUA-L35 Wide (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH1632L35	25	30		10-12
SET ZZ BASE MS-S-ZZMUA-L45 Narrow (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH1633L45	25	30		10-12
SET ZZ BASE MS-S-ZZMUA-L45 Standard (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH1631L45	25	30		10-12
SET ZZ BASE MS-S-ZZMUA-L45 Wide (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH1632L45	25	30		10-12
ZZ BASE MS-S-ASCSC Narrow (Angled screw-channel Screw)	BSKE1603	20			10-12
ZZ BASE MS-S-ASCSC Standard (Angled screw-channel Screw)	BSKE1601	20			10-12
ZZ BASE MS-S-ASCSC Wide (Angled screw-channel Screw)	BSKE1602	20			10-12
ZZ BASE MS-S-BSC Narrow (Black Screw)	BSIE1603	30			10-12
ZZ BASE MS-S-BSC Standard (Black Screw)	BSIE1601	30			10-12
ZZ BASE MS-S-BSC Wide (Black Screw)	BSIE1602	30			10-12
ZZ BASE MS-S-BZSC Narrow (Black Zirconia screw)	BSIE1633	30			10-12
ZZ BASE MS-S-BZSC Standard (Black Zirconia screw)	BSIE1631	30			10-12
ZZ BASE MS-S-BZSC Wide (Black Zirconia screw)	BSIE1632	30			10-12
ZZ BASE MS-S-ICSC Narrow (Impression coping Screw)	BSFE1603	15			10-12
ZZ BASE MS-S-ICSC Standard (Impression coping Screw)	BSFE1601	15			10-12
ZZ BASE MS-S-ICSC Wide (Impression coping Screw)	BSFE1602	15			10-12
ZZ BASE MS-S-LBSC Narrow (Lab Screw)	BSHE1603	30			10-12
ZZ BASE MS-S-LBSC Standard (Lab Screw)	BSHE1601	30			10-12
ZZ BASE MS-S-LBSC Wide (Lab Screw)	BSHE1602	30			10-12
ZZ BASE MS-S-LBZSC Narrow (Lab Zirconia screw)	BSHE1633	30			10-12
ZZ BASE MS-S-LBZSC Standard (Lab Zirconia screw)	BSHE1631	30			10-12
ZZ BASE MS-S-LBZSC Wide (Lab Zirconia screw)	BSHE1632	30			10-12
ZZ BASE MS-S-MI Narrow (Model implant)	MIAG1603				10-12
ZZ BASE MS-S-MI Standard (Model implant)	MIAG1601				10-12
ZZ BASE MS-S-MI Wide (Model implant)	MIAG1602				10-12
ZZ BASE MS-S-SC Narrow (Screw)	BSAE1603	30			10-12
ZZ BASE MS-S-SC Standard (Screw)	BSAE1601	30			10-12
ZZ BASE MS-S-SC Wide (Screw)	BSAE1602	30			10-12

MIS® SEVEN

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
ZZ BASE MS-S-SC-G Narrow (Screw Gold)	BSCE1603	30			10-12
ZZ BASE MS-S-SC-G Standard (Screw Gold)	BSCE1601	30			10-12
ZZ BASE MS-S-SC-G Wide (Screw Gold)	BSCE1602	30			10-12
ZZ BASE MS-S-WMSM-L12 Narrow (White Metal Scanmarker L12 with Screw)	ZBED1633L12	15			10-12
ZZ BASE MS-S-WMSM-L12 Standard (White Metal Scanmarker L12 with Screw)	ZBED1631L12	15			10-12
ZZ BASE MS-S-WMSM-L12 Wide (White Metal Scanmarker L12 with Screw)	ZBED1632L12	15			10-12
ZZ BASE MS-S-ZSC Narrow (Zirconia screw)	BSAE1633	30			10-12
ZZ BASE MS-S-ZSC Standard (Zirconia screw)	BSAE1631	30			10-12
ZZ BASE MS-S-ZSC Wide (Zirconia screw)	BSAE1632	30			10-12
ZZ BASE MS-S-ZSC-G Narrow (Zirconia Screw Gold)	BSCE1633	30			10-12
ZZ BASE MS-S-ZSC-G Standard (Zirconia Screw Gold)	BSCE1631	30			10-12
ZZ BASE MS-S-ZSC-G Wide (Zirconia Screw Gold)	BSCE1632	30			10-12
ZZ BASE MS-S-ZZMUASC Narrow (Zirkonzahn Multi Unit Abutment Screw)	BSEE1603	20			10-12
ZZ BASE MS-S-ZZMUASC Standard (Zirkonzahn Multi Unit Abutment Screw)	BSEE1601	25			10-12
ZZ BASE MS-S-ZZMUASC Wide (Zirkonzahn Multi Unit Abutment Screw)	BSEE1602	25			10-12

MIS® V3					
Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE MS-V3-ABCH Narrow (Raw-Abutment® customizable hex and screw)	BSAF9381	30			10-12
SET ZZ BASE MS-V3-ABCH Standard (Raw-Abutment® customizable hex and screw)	BSAF9382	30			10-12
SET ZZ BASE MS-V3-ABCH-D14 Standard (Raw-Abutment® customizable hex D14 and screw)	BSAF9382D14	30			10-12
SET ZZ BASE MS-V3-CCTB Narrow (Conical cemented titan base non hex and screw)	BSAD9331	30			10-12
SET ZZ BASE MS-V3-CCTB Standard (Conical cemented titan base non hex and screw)	BSAD9332	30			10-12
SET ZZ BASE MS-V3-CCTB-K85 Narrow (Conical cemented titan base non hex K85 and screw)	BSAD9331K85	30			10-12
SET ZZ BASE MS-V3-CCTB-K85 Standard (Conical cemented titan base non hex K85 and screw)	BSAD9332K85	30			10-12
SET ZZ BASE MS-V3-ICOTH Narrow (Impression coping open tray hex and Screw and Driver)	BSAQ9341	15			10-12
SET ZZ BASE MS-V3-ICOTH Standard (Impression coping open tray hex and Screw and Driver)	BSAQ9342	15			10-12
SET ZZ BASE MS-V3-NTB Standard (Narrow titan base non hex and screw)	BSAP9332	30			10-12
SET ZZ BASE MS-V3-NTBH Standard (Narrow titan base hex and screw)	BSAP9382	30			10-12
SET ZZ BASE MS-V3-PCTBH Narrow (Parallel cemented titan base hex and screw)	BSAB9381	30			10-12
SET ZZ BASE MS-V3-PCTBH Standard (Parallel cemented titan base hex and screw)	BSAB9382	30			10-12
SET ZZ BASE MS-V3-PCTBH-K85 Narrow (Parallel cemented titan base hex K85 and screw)	BSAB9381K85	30			10-12
SET ZZ BASE MS-V3-PCTBH-K85 Standard (Parallel cemented titan base hex K85 and screw)	BSAB9382K85	30			10-12
SET ZZ BASE MS-V3-SA Narrow (Scananalog and screw)	ZBDE9331	15			10-12
SET ZZ BASE MS-V3-SA Standard (Scananalog and screw)	ZBDE9332	15			10-12
SET ZZ BASE MS-V3-SM Narrow (Scanmarker and screw)	ZBAD9331	30			10-12
SET ZZ BASE MS-V3-SM Standard (Scanmarker and screw)	ZBAD9332	30			10-12
SET ZZ BASE MS-V3-WSM Narrow (White Scanmarker and screw)	ZBCD9331	30			10-12
SET ZZ BASE MS-V3-WSM Standard (White Scanmarker and screw)	ZBCD9332	30			10-12
ZZ BASE MS-V3-ASCSC Narrow (Angled screw-channel Screw)	BSKE9301	20			10-12
ZZ BASE MS-V3-ICSC Narrow (Impression coping Screw)	BSFE9301	15			10-12
ZZ BASE MS-V3-ICSC Standard (Impression coping Screw)	BSFE9302	15			10-12
ZZ BASE MS-V3-MI Narrow (Model implant)	MIAG9301				10-12
ZZ BASE MS-V3-MI Standard (Model implant)	MIAG9302				10-12
ZZ BASE MS-V3-SC Narrow (Screw)	BSAE9301	30			10-12
ZZ BASE MS-V3-SC Standard (Screw)	BSAE9302	30			10-12
ZZ BASE MS-V3-WMSM-L12 Standard (White Metal Scanmarker L12 with Screw)	ZBED9332L12	15			10-12

MIS® V3

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE MS-V3-ZSC Narrow (Zirconia Screw)	BSAE9331	30			10-12
ZZ BASE MS-V3-ZSC Standard (Zirconia screw)	BSAE9332	30			10-12

MOZO-GRAU® TICARE® TAPERED SCREW®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE MG-TS-CCTB 3.4 10° (Conical cemented titan base non hex and screw)	BSAD5131	20			10-12
SET ZZ BASE MG-TS-CCTB 4.1 10° (Conical cemented titan base non hex and screw)	BSAD5132	20			10-12
SET ZZ BASE MG-TS-CCTB 4.1 30° (Conical cemented titan base non hex and screw)	BSAD5133	20			10-12
SET ZZ BASE MG-TS-CCTB 5.0 10° (Conical cemented titan base non hex and screw)	BSAD5134	20			10-12
SET ZZ BASE MG-TS-CCTB 5.0 30° (Conical cemented titan base non hex and screw)	BSAD5135	20			10-12
SET ZZ BASE MG-TS-PCTBH 4.1 10° (Parallel cemented titan base hex and screw)	BSAB5182	20			10-12
SET ZZ BASE MG-TS-PCTBH 5.0 10° (Parallel cemented titan base hex and screw)	BSAB5184	20			10-12
SET ZZ BASE MG-TS-SM 3.4 10° (Scanmarker and screw)	ZBAD5131	20			10-12
SET ZZ BASE MG-TS-SM 4.1 10° (Scanmarker and screw)	ZBAD5132	20			10-12
SET ZZ BASE MG-TS-SM 4.1 30° (Scanmarker and screw)	ZBAD5133	20			10-12
SET ZZ BASE MG-TS-SM 5.0 10° (Scanmarker and screw)	ZBAD5134	20			10-12
SET ZZ BASE MG-TS-SM 5.0 30° (Scanmarker and screw)	ZBAD5135	20			10-12
ZZ BASE MG-TS-SC 3.4 10° (Screw)	BSAE5101	20			10-12
ZZ BASE MG-TS-SC 4.1 10° (Screw)	BSAE5102	20			10-12
ZZ BASE MG-TS-SC 4.1 30° (Screw)	BSAE5103	20			10-12
ZZ BASE MG-TS-SC 5.0 10° (Screw)	BSAE5104	20			10-12
ZZ BASE MG-TS-SC 5.0 30° (Screw)	BSAE5105	20			10-12
ZZ BASE MG-TS-ZSC 3.4 10° (Zirconia screw)	BSAE5131	20			10-12
ZZ BASE MG-TS-ZSC 4.1 10° (Zirconia screw)	BSAE5132	20			10-12
ZZ BASE MG-TS-ZSC 4.1 30° (Zirconia screw)	BSAE5133	20			10-12
ZZ BASE MG-TS-ZSC 5.0 10° (Zirconia screw)	BSAE5134	20			10-12
ZZ BASE MG-TS-ZSC 5.0 30° (Zirconia screw)	BSAE5135	20			10-12

NEOBIOTECH IS

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE NT-IS-ABCH 4.3 (Raw-Abutment® customizable hex and screw)	BSAF7381	30			10-12
SET ZZ BASE NT-IS-ASCTBH-K80 4.3 (Angled screw-channel titan base hex K80 and screw)	BSAT7381K80	20			10-12
SET ZZ BASE NT-IS-ASCTB-K80 4.3 (Angled screw-channel titan base non hex K80 and screw)	BSAT7331K80	20			10-12
SET ZZ BASE NT-IS-CCTB 4.3 (Conical cemented titan base non hex and screw)	BSAD7331	30			10-12
SET ZZ BASE NT-IS-CCTB-K85 4.3 (Conical cemented titan base non hex K85 and screw)	BSAD7331K85	30			10-12
SET ZZ BASE NT-IS-ICOTH 4.3 (Impression coping open tray hex and Screw and Driver)	BSAQ7341	15			10-12
SET ZZ BASE NT-IS-NTB 4.3 (Narrow titan base non hex and screw)	BSAP7331	30			10-12
SET ZZ BASE NT-IS-NTBH 4.3 (Narrow titan base hex and screw)	BSAP7381	30			10-12
SET ZZ BASE NT-IS-PCTBH 4.3 (Parallel cemented titan base hex and screw)	BSAB7381	30			10-12
SET ZZ BASE NT-IS-PCTBH-K85 4.3 (Parallel cemented titan base hex K85 and screw)	BSAB7381K85	30			10-12
SET ZZ BASE NT-IS-SA 4.3 (Scanalog and screw)	ZBDE7331	15			10-12
SET ZZ BASE NT-IS-SM 4.3 (Scanmarker and screw)	ZBAD7331	30			10-12
SET ZZ BASE NT-IS-WSM 4.3 (White Scanmarker and screw)	ZBCD7331	30			10-12
ZZ BASE NB-IS-ZSC 4.3 (Zirconia Screw)	BSAE7331	30			10-12
ZZ BASE NT-IS-ASCSC 4.3 (Angled screw-channel Screw)	BSKE7301	20			10-12
ZZ BASE NT-IS-ICSC 4.3 (Impression coping Screw)	BSFE7301	15			10-12
ZZ BASE NT-IS-SC 4.3 (Screw)	BSAE7301	30			10-12
ZZ BASE NT-IS-WMSM-L12 4.3 (White Metal Scanmarker L12 with Screw)	ZBED7331L12	15			10-12
ZZ BASE NT-IS-WMSM-L15 4.3 (White Metal Scanmarker L15 with Screw)	ZBED7331L15	15			10-12

NEOSS® PROACTIVE®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE NE-PA-ABCH 3.25 narrow (Raw-Abutment® customizable hex and screw)	BSAF3991	32			10-12
SET ZZ BASE NE-PA-ABCH 3.5/4.0/4.5/5.0/5.5 tapered (Raw-Abutment® customizable hex and screw)	BSAF3992	32			10-12
SET ZZ BASE NE-PA-ASCTBH-K80 3.25 narrow (Angled screw-channel titan base hex K80 and screw)	BSAT3981K80	20			10-12
SET ZZ BASE NE-PA-ASCTBH-K80 3.5/4.0/4.5/5.0/5.5 tapered (Angled screw-channel titan base hex and screw)	BSAT3982K80	20			10-12
SET ZZ BASE NE-PA-ASCTB-K80 3.25 narrow (Angled screw-channel titan base non hex K80 and screw)	BSAT3931K80	20			10-12
SET ZZ BASE NE-PA-ASCTB-K80 3.5/4.0/4.5/5.0/5.5 tapered (Angled screw-channel titan base non hex and screw)	BSAT3932K80	20			10-12
SET ZZ BASE NE-PA-CCTB 3.25 narrow (Conical cemented titan base non hex and screw)	BSAD3931	32			10-12
SET ZZ BASE NE-PA-CCTB 3.5/4.0/4.5/5.0/5.5 tapered (Conical cemented titan base non hex and screw)	BSAD3932	32			10-12
SET ZZ BASE NE-PA-CCTB-K85 3.25 narrow (Conical cemented titan base non hex K85 and screw)	BSAD3931K85	32			10-12
SET ZZ BASE NE-PA-CCTB-K85 3.5/4.0/4.5/5.0/5.5 tapered (Conical cemented titan base non hex and screw)	BSAD3932K85	32			10-12
SET ZZ BASE NE-PA-CCTB-L10 3.25 narrow (Conical cemented titan base non hex L10 and screw)	BSAD3931L10	32			10-12
SET ZZ BASE NE-PA-CCTB-L10 3.5/4.0/4.5/5.0/5.5 tapered (Conical cemented titan base non hex and screw)	BSAD3932L10	32			10-12
SET ZZ BASE NE-PA-CCTB-L20 3.25 narrow (Conical cemented titan base non hex L20 and screw)	BSAD3931L20	32			10-12
SET ZZ BASE NE-PA-CCTB-L20 3.5/4.0/4.5/5.0/5.5 tapered (Conical cemented titan base non hex and screw)	BSAD3932L20	32			10-12
SET ZZ BASE NE-PA-CCTB-L30 3.25 narrow (Conical cemented titan base non hex L30 and screw)	BSAD3931L30	32			10-12
SET ZZ BASE NE-PA-CCTB-L30 3.5/4.0/4.5/5.0/5.5 tapered (Conical cemented titan base non hex and screw)	BSAD3932L30	32			10-12
SET ZZ BASE NE-PA-CCTB-L40 3.25 narrow (Conical cemented titan base non hex L40 and screw)	BSAD3931L40	32			10-12
SET ZZ BASE NE-PA-CCTB-L40 3.5/4.0/4.5/5.0/5.5 tapered (Conical cemented titan base non hex and screw)	BSAD3932L40	32			10-12
SET ZZ BASE NE-PA-ICOT 3.25 narrow (Impression coping open tray non hex and Screw and Drive)	BSAQ3961	15			10-12
SET ZZ BASE NE-PA-ICOT 3.5/4.0/4.5/5.0/5.5 tapered (Impression coping open tray non hex and screw)	BSAQ3962	15			10-12
SET ZZ BASE NE-PA-ICOTH 3.25 Narrow (Impression coping open tray hex and Screw and Drive)	BSAQ3941	15			10-12
SET ZZ BASE NE-PA-ICOTH 3.5/4.0/4.5/5.0/5.5 Tapered (Impression coping open tray hex and screw)	BSAQ3942	15			10-12
SET ZZ BASE NE-PA-NTB 3.25 narrow (Narrow titan base non hex and screw)	BSAP3931	32			10-12
SET ZZ BASE NE-PA-NTBH 3.25 narrow (Narrow titan base hex and screw)	BSAP3981	32			10-12
SET ZZ BASE NE-PA-PCTBH 3.25 narrow (Parallel cemented titan base hex and screw)	BSAB3981	32			10-12
SET ZZ BASE NE-PA-PCTBH 3.5/4.0/4.5/5.0/5.5 tapered (Parallel cemented titan base hex and screw)	BSAB3982	32			10-12
SET ZZ BASE NE-PA-PCTBH-K85 3.25 narrow (Parallel cemented titan base hex K85 and screw)	BSAB3981K85	32			10-12
SET ZZ BASE NE-PA-PCTBH-K85 3.5/4.0/4.5/5.0/5.5 tapered (Parallel cemented titan base hex and screw)	BSAB3982K85	32			10-12
SET ZZ BASE NE-PA-PCTBH-L10 3.25 narrow (Parallel cemented titan base hex L10 and screw)	BSAB3981L10	32			10-12

NEOSS® PROACTIVE®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE NE-PA-PCTBH-L10 3.5/4.0/4.5/5.0/5.5 tapered (Parallel cemented titan base hex L	BSAB3982L10	32			10-12
SET ZZ BASE NE-PA-PCTBH-L20 3.25 narrow (Parallel cemented titan base hex L20 and screw)	BSAB3981L20	32			10-12
SET ZZ BASE NE-PA-PCTBH-L20 3.5/4.0/4.5/5.0/5.5 tapered (Parallel cemented titan base hex L	BSAB3982L20	32			10-12
SET ZZ BASE NE-PA-PCTBH-L30 3.25 narrow (Parallel cemented titan base hex L30 and screw)	BSAB3981L30	32			10-12
SET ZZ BASE NE-PA-PCTBH-L30 3.5/4.0/4.5/5.0/5.5 tapered (Parallel cemented titan base hex L	BSAB3982L30	32			10-12
SET ZZ BASE NE-PA-PCTBH-L40 3.25 narrow (Parallel cemented titan base hex L40 and screw)	BSAB3981L40	32			10-12
SET ZZ BASE NE-PA-PCTBH-L40 3.5/4.0/4.5/5.0/5.5 tapered (Parallel cemented titan base hex L	BSAB3982L40	32			10-12
SET ZZ BASE NE-PA-SA 3.25 narrow (Scananalog and screw)	ZBDE3931	15			10-12
SET ZZ BASE NE-PA-SA 3.5/4.0/4.5/5.0/5.5 tapered (Scananalog and screw)	ZBDE3932	15			10-12
SET ZZ BASE NE-PA-SM 3.25 narrow (Scanmarker and screw)	ZBAD3931	32			10-12
SET ZZ BASE NE-PA-SM 3.5/4.0/4.5/5.0/5.5 tapered (Scanmarker and screw)	ZBAD3932	32			10-12
SET ZZ BASE NE-PA-ZZMUA-L15 3.25 narrow (Zirkonzahn Multi Unit Abutment non hex L15 and	BSAH3931L15	25	32		10-12
SET ZZ BASE NE-PA-ZZMUA-L15 3.5/4.0/4.5/5.0/5.5 tapered (Zirkonzahn Multi Unit Abutment noi	BSAH3932L15	25	32		10-12
SET ZZ BASE NE-PA-ZZMUA-L25 3.25 narrow (Zirkonzahn Multi Unit Abutment non hex L25 and	BSAH3931L25	25	32		10-12
SET ZZ BASE NE-PA-ZZMUA-L25 3.5/4.0/4.5/5.0/5.5 tapered (Zirkonzahn Multi Unit Abutment noi	BSAH3932L25	25	32		10-12
SET ZZ BASE NE-PA-ZZMUA-L35 3.25 narrow (Zirkonzahn Multi Unit Abutment non hex L35 and	BSAH3931L35	25	32		10-12
SET ZZ BASE NE-PA-ZZMUA-L35 3.5/4.0/4.5/5.0/5.5 tapered (Zirkonzahn Multi Unit Abutment noi	BSAH3932L35	25	32		10-12
SET ZZ BASE NE-PA-ZZMUA-L45 3.25 narrow (Zirkonzahn Multi Unit Abutment non hex L45 and	BSAH3931L45	25	32		10-12
SET ZZ BASE NE-PA-ZZMUA-L45 3.5/4.0/4.5/5.0/5.5 tapered (Zirkonzahn Multi Unit Abutment noi	BSAH3932L45	25	32		10-12
ZZ BASE NE-PA-ASCSC 3.25 narrow (Angled screw-channel Screw)	BSKE3901	20			10-12
ZZ BASE NE-PA-ASCSC 3.5/4.0/4.5/5.0/5.5 tapered (Angled screw-channel Screw)	BSKE3902	20			10-12
ZZ BASE NE-PA-ICSC 3.25 Narrow (Impression coping Screw)	BSFE3901	15			10-12
ZZ BASE NE-PA-ICSC 3.5/4.0/4.5/5.0/5.5 Tapered (Impression coping Screw)	BSFE3902	15			10-12
ZZ BASE NE-PA-MI 3.25 narrow (Model implant)	MIAG3901				10-12
ZZ BASE NE-PA-MI 3.5/4.0/4.5/5.0/5.5 tapered (Model implant)	MIAG3902				10-12
ZZ BASE NE-PA-SC 3.25 narrow (Screw)	BSAE3901	32			10-12
ZZ BASE NE-PA-SC 3.5/4.0/4.5/5.0/5.5 tapered (Screw)	BSAE3902	32			10-12
ZZ BASE NE-PA-WMSM-L12 3.25 (White Metal Scanmarker L12 with Screw)	ZBED3931L12	15			10-12
ZZ BASE NE-PA-WMSM-L12 3.5/4.0/4.5/5.0/5.5 (White Metal Scanmarker L12 with Screw)	ZBED3932L12	15			10-12

NEOSS® PROACTIVE®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
ZZ BASE NE-PA-ZSC 3.25 narrow (Zirconia screw)	BSAE3931	32			10-12
ZZ BASE NE-PA-ZSC 3.5/4.0/4.5/5.0/5.5 tapered (Zirconia screw)	BSAE3932	32			10-12

NOBEL BIOCARE® BRÄNEMARK SYSTEM® MKIII/SHORTY/GROOVY® / NOBE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE NB-B-ABCH 3.5 (Raw-Abutment® customizable hex and screw)	BSAF0291	25			10-12
SET ZZ BASE NB-B-ABCH 4.1 (Raw-Abutment® customizable hex and screw)	BSAF0292	35			10-12
SET ZZ BASE NB-B-ABCH 5.1 (Raw-Abutment® customizable hex and screw)	BSAF0293	35			10-12
SET ZZ BASE NB-B-ABCH-D14 3.5 (Raw-Abutment® customizable hex D14 and screw)	BSAF0291D14	25			10-12
SET ZZ BASE NB-B-ABCH-D14 4.1 (Raw-Abutment® customizable hex D14 and screw)	BSAF0292D14	35			10-12
SET ZZ BASE NB-B-ABCH-D14 5.1 (Raw-Abutment® customizable hex D14 and screw)	BSAF0293D14	35			10-12
SET ZZ BASE NB-B-ASCTBH-K80 3.5 (Angled screw-channel titan base hex K80 and screw)	BSAT0281K80	20			10-12
SET ZZ BASE NB-B-ASCTBH-K80 4.1 (Angled screw-channel titan base hex K80 and screw)	BSAT0282K80	20			10-12
SET ZZ BASE NB-B-ASCTB-K80 3.5 (Angled screw-channel titan base non hex K80 and screw)	BSAT0231K80	20			10-12
SET ZZ BASE NB-B-ASCTB-K80 4.1 (Angled screw-channel titan base non hex K80 and screw)	BSAT0232K80	20			10-12
SET ZZ BASE NB-B-CCTB 3.5 (Conical cemented titan base non hex and screw)	BSAD0231	25			10-12
SET ZZ BASE NB-B-CCTB 4.1 (Conical cemented titan base non hex and screw)	BSAD0232	35			10-12
SET ZZ BASE NB-B-CCTB 5.1 (Conical cemented titan base non hex and screw)	BSAD0233	35			10-12
SET ZZ BASE NB-B-CCTB-G 3.5 (Conical cemented titan base non hex and screw Gold)	BSCD0231	25			10-12
SET ZZ BASE NB-B-CCTB-G 3.5 (Conical cemented titan base non hex and screw Gold)	BSCD0241	25			10-12
SET ZZ BASE NB-B-CCTB-G 4.1 (Conical cemented titan base non hex and screw Gold)	BSCD0232	35			10-12
SET ZZ BASE NB-B-CCTB-G 4.1 (Conical cemented titan base non hex and screw Gold)	BSCD0242	35			10-12
SET ZZ BASE NB-B-CCTB-G 5.1 (Conical cemented titan base non hex and screw Gold)	BSCD0233	35			10-12
SET ZZ BASE NB-B-CCTB-G 5.1 (Conical cemented titan base non hex and screw Gold)	BSCD0243	35			10-12
SET ZZ BASE NB-B-CCTB-K85 3.5 (Conical cemented titan base non hex K85 and screw)	BSAD0231K85	25			10-12
SET ZZ BASE NB-B-CCTB-K85 4.1 (Conical cemented titan base non hex K85 and screw)	BSAD0232K85	35			10-12
SET ZZ BASE NB-B-CCTB-K85 5.1 (Conical cemented titan base non hex K85 and screw)	BSAD0233K85	35			10-12
SET ZZ BASE NB-B-CCTB-L10 3.5 (Conical cemented titan base non hex L10 and screw)	BSAD0231L10	25			10-12
SET ZZ BASE NB-B-CCTB-L10 4.1 (Conical cemented titan base non hex L10 and screw)	BSAD0232L10	35			10-12
SET ZZ BASE NB-B-CCTB-L10 5.1 (Conical cemented titan base non hex L10 and screw)	BSAD0233L10	35			10-12
SET ZZ BASE NB-B-CCTB-L10-G 3.5 (Conical cemented titan base non hex L10 and screw Gold)	BSCD0231L10	25			10-12
SET ZZ BASE NB-B-CCTB-L10-G 4.1 (Conical cemented titan base non hex L10 and screw Gold)	BSCD0232L10	35			10-12
SET ZZ BASE NB-B-CCTB-L20 3.5 (Conical cemented titan base non hex L20 and screw)	BSAD0231L20	25			10-12
SET ZZ BASE NB-B-CCTB-L20 4.1 (Conical cemented titan base non hex L20 and screw)	BSAD0232L20	35			10-12

NOBEL BIOCARE® BRÄNEMARK SYSTEM® MKIII/SHORTY/GROOVY® / NOBE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE NB-B-CCTB-L20 5.1 (Conical cemented titan base non hex L20 and screw)	BSAD0233L20	35			10-12
SET ZZ BASE NB-B-CCTB-L20-G 3.5 (Conical cemented titan base non hex L20 and screw Gold)	BSCD0231L20	25			10-12
SET ZZ BASE NB-B-CCTB-L20-G 4.1 (Conical cemented titan base non hex L20 and screw Gold)	BSCD0232L20	35			10-12
SET ZZ BASE NB-B-CCTB-L30 3.5 (Conical cemented titan base non hex L30 and screw)	BSAD0231L30	25			10-12
SET ZZ BASE NB-B-CCTB-L30 4.1 (Conical cemented titan base non hex L30 and screw)	BSAD0232L30	35			10-12
SET ZZ BASE NB-B-CCTB-L30 5.1 (Conical cemented titan base non hex L30 and screw)	BSAD0233L30	35			10-12
SET ZZ BASE NB-B-CCTB-L30-G 4.1 (Conical cemented titan base non hex L30 and screw Gold)	BSCD0232L30	35			10-12
SET ZZ BASE NB-B-ICOT 3.5 (Impression coping open tray non hex and Screw and Driver)	BSAQ0261	15			10-12
SET ZZ BASE NB-B-ICOT 4.1 (Impression coping open tray non hex and Screw and Driver)	BSAQ0262	15			10-12
SET ZZ BASE NB-B-ICOTH 3.5 (Impression coping open tray hex and Screw and Driver)	BSAQ0241	15			10-12
SET ZZ BASE NB-B-ICOTH 4.1 (Impression coping open tray hex and Screw and Driver)	BSAQ0242	15			10-12
SET ZZ BASE NB-B-ICOTH 5.1 (Impression coping open tray hex and Screw and Driver)	BSAQ0243	15			10-12
SET ZZ BASE NB-B-NTB 3.5 (Narrow titan base non hex and screw)	BSAP0231	25			10-12
SET ZZ BASE NB-B-NTB 4.1 (Narrow titan base non hex and screw)	BSAP0232	35			10-12
SET ZZ BASE NB-B-NTBH 3.5 (Narrow titan base hex and screw)	BSAP0281	25			10-12
SET ZZ BASE NB-B-NTBH 4.1 (Narrow titan base hex and screw)	BSAP0282	35			10-12
SET ZZ BASE NB-B-PCTBH 3.5 (Parallel cemented titan base hex and screw)	BSAB0281	25			10-12
SET ZZ BASE NB-B-PCTBH 3.5 (Parallel cemented titan base hex and screw)	BSAB0291	25			10-12
SET ZZ BASE NB-B-PCTBH 4.1 (Parallel cemented titan base hex and screw)	BSAB0282	35			10-12
SET ZZ BASE NB-B-PCTBH 4.1 (Parallel cemented titan base hex and screw)	BSAB0292	35			10-12
SET ZZ BASE NB-B-PCTBH 5.1 (Parallel cemented titan base hex and screw)	BSAB0283	35			10-12
SET ZZ BASE NB-B-PCTBH 5.1 (Parallel cemented titan base hex and screw)	BSAB0293	35			10-12
SET ZZ BASE NB-B-PCTBH-G 3.5 (Parallel cemented titan base hex and screw Gold)	BSCB0291	25			10-12
SET ZZ BASE NB-B-PCTBH-G 4.1 (Parallel cemented titan base hex and screw Gold)	BSCB0292	35			10-12
SET ZZ BASE NB-B-PCTBH-G 4.1 (Parallel cemented titan base hex and screw Gold)	BSCB0282	35			10-12
SET ZZ BASE NB-B-PCTBH-G 5.1 (Parallel cemented titan base hex and screw Gold)	BSCB0293	35			10-12
SET ZZ BASE NB-B-PCTBH-G 5.1 (Parallel cemented titan base hex and screw Gold)	BSCB0283	35			10-12
SET ZZ BASE NB-B-PCTBH-K85 3.5 (Parallel cemented titan base hex K85 and screw)	BSAB0281K85	25			10-12
SET ZZ BASE NB-B-PCTBH-K85 4.1 (Parallel cemented titan base hex K85 and screw)	BSAB0282K85	35			10-12

NOBEL BIOCARE® BRÄNEMARK SYSTEM® MKIII/SHORTY/GROOVY® / NOBE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE NB-B-PCTBH-K85 5.1 (Parallel cemented titan base hex K85 and screw)	BSAB0283K85	35			10-12
SET ZZ BASE NB-B-PCTBH-L10 3.5 (Parallel cemented titan base hex L10 and screw)	BSAB0291L10	25			10-12
SET ZZ BASE NB-B-PCTBH-L10 4.1 (Parallel cemented titan base hex L10 and screw)	BSAB0292L10	35			10-12
SET ZZ BASE NB-B-PCTBH-L10 5.1 (Parallel cemented titan base hex L10 and screw)	BSAB0293L10	35			10-12
SET ZZ BASE NB-B-PCTBH-L10-G 4.1 (Parallel cemented titan base hex L10 and screw Gold)	BSCB0292L10	35			10-12
SET ZZ BASE NB-B-PCTBH-L10-G 5.1 (Parallel cemented titan base hex L10 and screw Gold)	BSCB0293L10	35			10-12
SET ZZ BASE NB-B-PCTBH-L20 3.5 (Parallel cemented titan base hex L20 and screw)	BSAB0291L20	25			10-12
SET ZZ BASE NB-B-PCTBH-L20 4.1 (Parallel cemented titan base hex L20 and screw)	BSAB0292L20	35			10-12
SET ZZ BASE NB-B-PCTBH-L20 5.1 (Parallel cemented titan base hex L20 and screw)	BSAB0293L20	35			10-12
SET ZZ BASE NB-B-PCTBH-L20-G 4.1 (Parallel cemented titan base hex L20 and screw Gold)	BSCB0292L20	35			10-12
SET ZZ BASE NB-B-PCTBH-L20-G 5.1 (Parallel cemented titan base hex L20 and screw Gold)	BSCB0293L20	35			10-12
SET ZZ BASE NB-B-PCTBH-L30 3.5 (Parallel cemented titan base hex L30 and screw)	BSAB0291L30	25			10-12
SET ZZ BASE NB-B-PCTBH-L30 4.1 (Parallel cemented titan base hex L30 and screw)	BSAB0292L30	35			10-12
SET ZZ BASE NB-B-PCTBH-L30 5.1 (Parallel cemented titan base hex L30 and screw)	BSAB0293L30	35			10-12
SET ZZ BASE NB-B-PCTBH-L30-G 4.1 (Parallel cemented titan base hex L30 and screw Gold)	BSCB0292L30	35			10-12
SET ZZ BASE NB-B-PCTBH-L40 3.5 (Parallel cemented titan base hex L40 and screw)	BSAB0291L40	25			10-12
SET ZZ BASE NB-B-PCTBH-L40 4.1 (Parallel cemented titan base hex L40 and screw)	BSAB0292L40	35			10-12
SET ZZ BASE NB-B-PCTBH-L40-G 4.1 (Parallel cemented titan base hex L40 and screw Gold)	BSCB0292L40	35			10-12
SET ZZ BASE NB-B-SA 3.5 (Scanalog and screw)	ZBDE0231	15			10-12
SET ZZ BASE NB-B-SA 4.1 (Scanalog and screw)	ZBDE0232	15			10-12
SET ZZ BASE NB-B-SA 5.1 (Scanalog and screw)	ZBDE0233	15			10-12
SET ZZ BASE NB-B-SM 3.5 (Scanmarker and screw)	ZBAD0231	25			10-12
SET ZZ BASE NB-B-SM 4.1 (Scanmarker and screw)	ZBAD0232	35			10-12
SET ZZ BASE NB-B-SM 5.1 (Scanmarker and screw)	ZBAD0233	35			10-12
SET ZZ BASE NB-B-WSM 4.1 (White Scanmarker and screw)	ZBCD0232	35			10-12
SET ZZ BASE NB-B-WSM 5.1 (White Scanmarker and screw)	ZBCD0233	35			10-12
SET ZZ BASE NB-B-ZZMUA17 3.5 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and implar	BSAI0231	20			10-12
SET ZZ BASE NB-B-ZZMUA17 4.1 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and implar	BSAI0232	30			10-12
SET ZZ BASE NB-B-ZZMUA17 5.1 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and implar	BSAI0233	35			10-12

NOBEL BIOCARE® BRÄNEMARK SYSTEM® MKIII/SHORTY/GROOVY® / NOBE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE NB-B-ZZMUA17-L10 3.5 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10 ai	BSAI0231L10	20			10-12
SET ZZ BASE NB-B-ZZMUA17-L10 4.1 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10 ai	BSAI0232L10	30			10-12
SET ZZ BASE NB-B-ZZMUA17-L10 5.1 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10 ai	BSAI0233L10	35			10-12
SET ZZ BASE NB-B-ZZMUA-L15 3.5 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH0231L15	25	25		10-12
SET ZZ BASE NB-B-ZZMUA-L15 4.1 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH0232L15	25	35		10-12
SET ZZ BASE NB-B-ZZMUA-L15 5.1 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH0233L15	25	35		10-12
SET ZZ BASE NB-B-ZZMUA-L25 3.5 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH0231L25	25	25		10-12
SET ZZ BASE NB-B-ZZMUA-L25 4.1 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH0232L25	25	35		10-12
SET ZZ BASE NB-B-ZZMUA-L25 5.1 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH0233L25	25	35		10-12
SET ZZ BASE NB-B-ZZMUA-L35 3.5 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH0231L35	25	25		10-12
SET ZZ BASE NB-B-ZZMUA-L35 4.1 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH0232L35	25	35		10-12
SET ZZ BASE NB-B-ZZMUA-L35 5.1 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH0233L35	25	35		10-12
SET ZZ BASE NB-B-ZZMUA-L45 3.5 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH0231L45	25	25		10-12
SET ZZ BASE NB-B-ZZMUA-L45 4.1 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH0232L45	25	35		10-12
SET ZZ BASE NB-B-ZZMUA-L45 5.1 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH0233L45	25	35		10-12
ZZ BASE NB-B-ASCSC 3.5 (Angled screw-channel Screw)	BSKE0201	20			10-12
ZZ BASE NB-B-ASCSC 4.1 (Angled screw-channel Screw)	BSKE0202	20			10-12
ZZ BASE NB-B-BSC 3.5 (Black Screw)	BSIE0201	25			10-12
ZZ BASE NB-B-BSC 4.1 (Black Screw)	BSIE0202	35			10-12
ZZ BASE NB-B-BSC 5.1 (Black Screw)	BSIE0203	35			10-12
ZZ BASE NB-B-BZSC 3.5 (Black Zirconia screw)	BSIE0231	25			10-12
ZZ BASE NB-B-BZSC 4.1 (Black Zirconia screw)	BSIE0232	35			10-12
ZZ BASE NB-B-BZSC 5.1 (Black Zirconia screw)	BSIE0233	35			10-12
ZZ BASE NB-B-ICSC 3.5 (Impression coping Screw)	BSFE0201	15			10-12
ZZ BASE NB-B-ICSC 4.1 (Impression coping Screw)	BSFE0202	15			10-12
ZZ BASE NB-B-ICSC 5.1 (Impression coping Screw)	BSFE0203	15			10-12
ZZ BASE NB-B-LBSC 3.5 (Lab Screw)	BSHE0201	25			10-12
ZZ BASE NB-B-LBSC 4.1 (Lab Screw)	BSHE0202	35			10-12
ZZ BASE NB-B-LBSC 5.1 (Lab Screw)	BSHE0203	35			10-12

NOBEL BIOCARE® BRÄNEMARK SYSTEM® MKIII/SHORTY/GROOVY® / NOBE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE NB-B-LBZSC 3.5 (Lab Zirconia screw)	BSHE0231	25			10-12
ZZ BASE NB-B-LBZSC 4.1 (Lab Zirconia screw)	BSHE0232	35			10-12
ZZ BASE NB-B-LBZSC 5.1 (Lab Zirconia screw)	BSHE0233	35			10-12
ZZ BASE NB-B-MI 3.5 (Model implant)	MIAG0201				10-12
ZZ BASE NB-B-MI 4.1 (Model implant)	MIAG0202				10-12
ZZ BASE NB-B-MI 5.1 (Model implant)	MIAG0203				10-12
ZZ BASE NB-B-SC 3.5 (Screw)	BSAE0201	25			10-12
ZZ BASE NB-B-SC 4.1 (Screw)	BSAE0202	35			10-12
ZZ BASE NB-B-SC 5.1 (Screw)	BSAE0203	35			10-12
ZZ BASE NB-B-SC-G 3.5 (Screw Gold)	BSCE0201	25			10-12
ZZ BASE NB-B-SC-G 4.1 (Screw Gold)	BSCE0202	35			10-12
ZZ BASE NB-B-SC-G 5.1 (Screw Gold)	BSCE0203	35			10-12
ZZ BASE NB-B-WMSM-L10 3.5 (White Metal Scanmarker L10 with Screw)	ZBED0231L10	15			10-12
ZZ BASE NB-B-WMSM-L10 4.1 (White Metal Scanmarker L10 with Screw)	ZBED0232L10	15			10-12
ZZ BASE NB-B-WMSM-L10 5.1 (White Metal Scanmarker L10 with Screw)	ZBED0233L10	15			10-12
ZZ BASE NB-B-WMSM-L12 3.5 (White Metal Scanmarker L12 with Screw)	ZBED0231L12	15			10-12
ZZ BASE NB-B-WMSM-L12 4.1 (White Metal Scanmarker L12 with Screw)	ZBED0232L12	15			10-12
ZZ BASE NB-B-WMSM-L12 5.1 (White Metal Scanmarker L12 with Screw)	ZBED0233L12	15			10-12
ZZ BASE NB-B-WMSM-L15 3.5 (White Metal Scanmarker L15 with Screw)	ZBED0231L15	15			10-12
ZZ BASE NB-B-WMSM-L15 4.1 (White Metal Scanmarker L15 with Screw)	ZBED0232L15	15			10-12
ZZ BASE NB-B-WMSM-L15 5.1 (White Metal Scanmarker L15 with Screw)	ZBED0233L15	15			10-12
ZZ BASE NB-B-ZSC 3.5 (Zirconia screw)	BSAE0231	25			10-12
ZZ BASE NB-B-ZSC 4.1 (Zirconia screw)	BSAE0232	35			10-12
ZZ BASE NB-B-ZSC 5.1 (Zirconia screw)	BSAE0233	35			10-12
ZZ BASE NB-B-ZSC-G 3.5 (Zirconia screw Gold)	BSCE0231	25			10-12
ZZ BASE NB-B-ZSC-G 4.1 (Zirconia screw Gold)	BSCE0232	35			10-12
ZZ BASE NB-B-ZSC-G 5.1 (Zirconia screw Gold)	BSCE0233	35			10-12
ZZ BASE NB-B-ZZMUASC 3.5 (Zirkonzahn Multi Unit Abutment Screw)	BSEE0201	20			10-12
ZZ BASE NB-B-ZZMUASC 4.1 (Zirkonzahn Multi Unit Abutment Screw)	BSEE0202	30			10-12

NOBEL BIOCARE® BRÄNEMARK SYSTEM® MKIII/SHORTY/GROOVY® / NOBE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE NB-B-ZZMUASC 5.1 (Zirkonzahn Multi Unit Abutment Screw)	BSEE0203	35			10-12

NOBEL BIOCARE® MULTI-UNIT ABUTMENT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE NB-MUA-CCTB NP/RP (Conical cemented titan base non hex and screw)	BSAD0441	15			10-12
SET ZZ BASE NB-MUA-CCTB NP/RP (Conical cemented titan base non hex and screw)	BSAD0431	15			10-12
SET ZZ BASE NB-MUA-CCTB WP (Conical cemented titan base non hex and screw)	BSAD0432	35			10-12
SET ZZ BASE NB-MUA-CCTB-G WP (Conical cemented titan base non hex and screw Gold)	BSCD0432	35			10-12
SET ZZ BASE NB-MUA-CCTB-K85 NP/RP (Conical cemented titan base non hex K85 and screw)	BSAD0431K85	15			10-12
SET ZZ BASE NB-MUA-ICOT NP/RP (Impression coping open tray non hex and Screw and Driver)	BSAQ0461	15			10-12
SET ZZ BASE NB-MUA-ICOT WP (Impression coping open tray non hex and Screw and Driver)	BSAQ0462	15			10-12
SET ZZ BASE NB-MUA-SA NP/RP (Scanalog and screw)	ZBDE0431	15			10-12
SET ZZ BASE NB-MUA-SA WP (Scanalog and screw)	ZBDE0432	15			10-12
SET ZZ BASE NB-MUA-SM NP/RP (Scanmarker and screw)	ZBAD0431	15			10-12
SET ZZ BASE NB-MUA-SM WP (Scanmarker and screw)	ZBAD0432	35			10-12
SET ZZ BASE NB-MUA-WSM WP (White Scanmarker and screw)	ZBCD0432	35			10-12
ZZ BASE NB-MUA-BSC NP/RP (Black Screw)	BSIE0401	15			10-12
ZZ BASE NB-MUA-BSC WP (Black Screw)	BSIE0402	35			10-12
ZZ BASE NB-MUA-BZSC NP/RP (Black Zirconia screw)	BSIE0431	15			10-12
ZZ BASE NB-MUA-BZSC WP (Black Zirconia screw)	BSIE0432	35			10-12
ZZ BASE NB-MUA-ICSC NP/RP (Impression coping Screw)	BSFE0401	15			10-12
ZZ BASE NB-MUA-ICSC WP (Impression coping Screw)	BSFE0402	15			10-12
ZZ BASE NB-MUA-LBSC NP/RP (Lab Screw)	BSHE0401	15			10-12
ZZ BASE NB-MUA-LBSC WP (Lab Screw)	BSHE0402	35			10-12
ZZ BASE NB-MUA-LBZSC NP/RP (Lab Zirconia screw)	BSHE0431	15			10-12
ZZ BASE NB-MUA-LBZSC WP (Lab Zirconia screw)	BSHE0432	35			10-12
ZZ BASE NB-MUA-MI NP/RP (Model implant)	MIAG0401				10-12
ZZ BASE NB-MUA-MI WP (Model implant)	MIAG0402				10-12
ZZ BASE NB-MUA-SC NP/RP (Screw)	BSAE0401	15			10-12
ZZ BASE NB-MUA-SC WP (Screw)	BSAE0402	35			10-12
ZZ BASE NB-MUA-SC-G NP/RP (Screw Gold)	BSCE0401	15			10-12
ZZ BASE NB-MUA-SC-G WP (Screw Gold)	BSCE0402	35			10-12
ZZ BASE NB-MUA-WMSM-L10 NP/RP (White Metal Scanmarker L10 with Screw)	ZBED0431L10	15			10-12

NOBEL BIOCARE® MULTI-UNIT ABUTMENT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE NB-MUA-WMSM-L10 WP (White Metal Scanmarker L10 with Screw)	ZBED0432L10	15			10-12
ZZ BASE NB-MUA-WMSM-L12 NP/RP (White Metal Scanmarker L12 with Screw)	ZBED0431L12	15			10-12
ZZ BASE NB-MUA-WMSM-L12 WP (White Metal Scanmarker L12 with Screw)	ZBED0432L12	15			10-12
ZZ BASE NB-MUA-WMSM-L15 NP/RP (White Metal Scanmarker L15 with Screw)	ZBED0431L15	15			10-12
ZZ BASE NB-MUA-WMSM-L15 WP (White Metal Scanmarker L15 with Screw)	ZBED0432L15	15			10-12
ZZ BASE NB-MUA-ZSC NP/RP (Zirconia screw)	BSAE0431	15			10-12
ZZ BASE NB-MUA-ZSC WP (Zirconia screw)	BSAE0432	35			10-12
ZZ BASE NB-MUA-ZSC-G NP/RP (Zirconia screw Gold)	BSCE0431	15			10-12
ZZ BASE NB-MUA-ZSC-G WP (Zirconia screw Gold)	BSCE0432	35			10-12

NOBEL BIO CARE® NOBEL ACTIVE®/NOBEL REPLACE® CC/NOBEL PARALLEL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
SET ZZ BASE NB-A-ABCH CC 3.0 (Raw-Abutment® customizable hex and screw)	BSAF2483	15			10-12
SET ZZ BASE NB-A-ABCH NP (Raw-Abutment® customizable hex and screw)	BSAF2481	25			10-12
SET ZZ BASE NB-A-ABCH RP (Raw-Abutment® customizable hex and screw)	BSAF2482	35			10-12
SET ZZ BASE NB-A-ABCH WP (Raw-Abutment® customizable hex and screw)	BSAF2484	35			10-12
SET ZZ BASE NB-A-ABCH-D14 CC 3.0 (Raw-Abutment® customizable hex D14 and screw)	BSAF2483D14	15			10-12
SET ZZ BASE NB-A-ABCH-D14 NP (Raw-Abutment® customizable hex D14 and screw)	BSAF2481D14	25			10-12
SET ZZ BASE NB-A-ABCH-D14 RP (Raw-Abutment® customizable hex D14 and screw)	BSAF2482D14	35			10-12
SET ZZ BASE NB-A-ABCH-D14 WP (Raw-Abutment® customizable hex D14 and screw)	BSAF2484D14	35			10-12
SET ZZ BASE NB-A-ABH-L15 NP (Abutment hex L15 and screw)	BSAC2431L15	25			10-12
SET ZZ BASE NB-A-ABH-L15 RP (Abutment hex L15 and screw)	BSAC2432L15	35			10-12
SET ZZ BASE NB-A-ABH-L15 WP (Abutment hex L15 and screw)	BSAC2434L15	35			10-12
SET ZZ BASE NB-A-ABH-L30 NP (Abutment hex L30 and screw)	BSAC2431L30	25			10-12
SET ZZ BASE NB-A-ABH-L30 RP (Abutment hex L30 and screw)	BSAC2432L30	35			10-12
SET ZZ BASE NB-A-ABH-L30 WP (Abutment hex L30 and screw)	BSAC2434L30	35			10-12
SET ZZ BASE NB-A-ASCTBH-K80 CC 3.0 (Angled screw-channel titan base hex K80 and screw)	BSAT2483K80	15			10-12
SET ZZ BASE NB-A-ASCTBH-K80 NP (Angled screw-channel titan base hex K80 and screw)	BSAT2481K80	20			10-12
SET ZZ BASE NB-A-ASCTBH-K80 RP (Angled screw-channel titan base hex K80 and screw)	BSAT2482K80	20			10-12
SET ZZ BASE NB-A-ASCTBH-K80 WP (Angled screw-channel titan base hex K80 and screw)	BSAT2484K80	20			10-12
SET ZZ BASE NB-A-ASCTB-K80 CC 3.0 (Angled screw-channel titan base non hex K80 and screw)	BSAT2433K80	15			10-12
SET ZZ BASE NB-A-ASCTB-K80 NP (Angled screw-channel titan base non hex K80 and screw)	BSAT2431K80	20			10-12
SET ZZ BASE NB-A-ASCTB-K80 RP (Angled screw-channel titan base non hex K80 and screw)	BSAT2432K80	20			10-12
SET ZZ BASE NB-A-ASCTB-K80 WP (Angled screw-channel titan base non hex K80 and screw)	BSAT2434K80	20			10-12
SET ZZ BASE NB-A-CCTB CC 3.0 (Conical cemented titan base non hex and screw)	BSAD2433	15			10-12
SET ZZ BASE NB-A-CCTB NP (Conical cemented titan base non hex and screw)	BSAD2431	25			10-12
SET ZZ BASE NB-A-CCTB RP (Conical cemented titan base non hex and screw)	BSAD2432	35			10-12
SET ZZ BASE NB-A-CCTB WP (Conical cemented titan base non hex and screw)	BSAD2434	35			10-12
SET ZZ BASE NB-A-CCTB-G CC 3.0 (Conical cemented titan base non hex and screw Gold)	BSCD2433	15			10-12
SET ZZ BASE NB-A-CCTB-G RP (Conical cemented titan base non hex and screw Gold)	BSCD2432	35			10-12
SET ZZ BASE NB-A-CCTB-G WP (Conical cemented titan base non hex and screw Gold)	BSCD2434	35			10-12

All information is subject to change. Errors and omissions excepted. Version: 21.09.2023

NOBEL BIO CARE® NOBEL ACTIVE® NOBEL REPLACE® CC/NOBEL PARALLEL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
SET ZZ BASE NB-A-CCTB-K85 CC 3.0 (Conical cemented titan base non hex K85 and screw)	BSAD2433K85	15			10-12
SET ZZ BASE NB-A-CCTB-K85 NP (Conical cemented titan base non hex K85 and screw)	BSAD2431K85	25			10-12
SET ZZ BASE NB-A-CCTB-K85 RP (Conical cemented titan base non hex K85 and screw)	BSAD2432K85	35			10-12
SET ZZ BASE NB-A-CCTB-K85 WP (Conical cemented titan base non hex K85 and screw)	BSAD2434K85	35			10-12
SET ZZ BASE NB-A-CCTB-L20 CC 3.0 (Conical cemented titan base non hex L20 and screw)	BSAD2433L20	15			10-12
SET ZZ BASE NB-A-CCTB-L20 NP (Conical cemented titan base non hex L20 and screw)	BSAD2431L20	25			10-12
SET ZZ BASE NB-A-CCTB-L20 RP (Conical cemented titan base non hex L20 and screw)	BSAD2432L20	35			10-12
SET ZZ BASE NB-A-CCTB-L20 WP (Conical cemented titan base non hex L20 and screw)	BSAD2434L20	35			10-12
SET ZZ BASE NB-A-CCTB-L20-G CC 3.0 (Conical cemented titan base non hex L20 and screw Gold)	BSCD2433L20	15			10-12
SET ZZ BASE NB-A-CCTB-L20-G NP (Conical cemented titan base non hex L20 and screw Gold)	BSCD2431L20	25			10-12
SET ZZ BASE NB-A-CCTB-L20-G RP (Conical cemented titan base non hex L20 and screw Gold)	BSCD2432L20	35			10-12
SET ZZ BASE NB-A-CCTB-L20-K85 NP (Conical cemented titan base non hex L20 K85 and screw)	BSAD2431L20K85	25			10-12
SET ZZ BASE NB-A-CCTB-L20-K85 RP (Conical cemented titan base non hex L20 K85 and screw)	BSAD2432L20K85	35			10-12
SET ZZ BASE NB-A-CCTB-L30 CC 3.0 (Conical cemented titan base non hex L30 and screw)	BSAD2433L30	15			10-12
SET ZZ BASE NB-A-CCTB-L30 NP (Conical cemented titan base non hex L30 and screw)	BSAD2431L30	25			10-12
SET ZZ BASE NB-A-CCTB-L30 RP (Conical cemented titan base non hex L30 and screw)	BSAD2432L30	35			10-12
SET ZZ BASE NB-A-CCTB-L30 WP (Conical cemented titan base non hex L30 and screw)	BSAD2434L30	35			10-12
SET ZZ BASE NB-A-CCTB-L30-G CC 3.0 (Conical cemented titan base non hex L30 and screw Gold)	BSCD2433L30	15			10-12
SET ZZ BASE NB-A-CCTB-L30-G NP (Conical cemented titan base non hex L30 and screw Gold)	BSCD2431L30	25			10-12
SET ZZ BASE NB-A-CCTB-L30-G RP (Conical cemented titan base non hex L30 and screw Gold)	BSCD2432L30	35			10-12
SET ZZ BASE NB-A-CCTB-L30-K85 NP (Conical cemented titan base non hex L30 K85 and screw)	BSAD2431L30K85	25			10-12
SET ZZ BASE NB-A-CCTB-L30-K85 RP (Conical cemented titan base non hex L30 K85 and screw)	BSAD2432L30K85	35			10-12
SET ZZ BASE NB-A-CCTB-L40 NP (Conical cemented titan base non hex L40 and screw)	BSAD2431L40	25			10-12
SET ZZ BASE NB-A-CCTB-L40 RP (Conical cemented titan base non hex L40 and screw)	BSAD2432L40	35			10-12
SET ZZ BASE NB-A-CCTB-L40 WP (Conical cemented titan base non hex L40 and screw)	BSAD2434L40	35			10-12
SET ZZ BASE NB-A-CCTB-L40-G NP (Conical cemented titan base non hex L40 and screw Gold)	BSCD2431L40	25			10-12
SET ZZ BASE NB-A-CCTB-L40-G RP (Conical cemented titan base non hex L40 and screw Gold)	BSCD2432L40	35			10-12
SET ZZ BASE NB-A-CCTB-L40-K85 NP (Conical cemented titan base non hex L40 K85 and screw)	BSAD2431L40K85	25			10-12
SET ZZ BASE NB-A-ICOT CC 3.0 (Impression coping open tray non hex and Screw and Driver)	BSAQ2463	15			10-12

NOBEL BIOCARE® NOBELACTIVE®/NOBELREPLACE® CC/NOBEL PARALLEL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE NB-A-ICOT NP (Impression coping open tray non hex and Screw and Driver)	BSAQ2461	15			10-12
SET ZZ BASE NB-A-ICOT RP (Impression coping open tray non hex and Screw and Driver)	BSAQ2462	15			10-12
SET ZZ BASE NB-A-ICOT WP (Impression coping open tray non hex and Screw and Driver)	BSAQ2464	15			10-12
SET ZZ BASE NB-A-ICOTH CC 3.0 (Impression coping open tray hex and Screw and Driver)	BSAQ2443	15			10-12
SET ZZ BASE NB-A-ICOTH NP (Impression coping open tray hex and Screw and Driver)	BSAQ2441	15			10-12
SET ZZ BASE NB-A-ICOTH RP (Impression coping open tray hex and Screw and Driver)	BSAQ2442	15			10-12
SET ZZ BASE NB-A-ICOTH WP (Impression coping open tray hex and Screw and Driver)	BSAQ2444	15			10-12
SET ZZ BASE NB-A-NTB CC 3.0 (Narrow titan base non hex and screw)	BSAP2433	15			10-12
SET ZZ BASE NB-A-NTB NP (Narrow titan base non hex and screw)	BSAP2431	25			10-12
SET ZZ BASE NB-A-NTB RP (Narrow titan base non hex and screw)	BSAP2432	35			10-12
SET ZZ BASE NB-A-NTBH CC 3.0 (Narrow titan base hex and screw)	BSAP2483	15			10-12
SET ZZ BASE NB-A-NTBH NP (Narrow titan base hex and screw)	BSAP2481	25			10-12
SET ZZ BASE NB-A-NTBH RP (Narrow titan base hex and screw)	BSAP2482	35			10-12
SET ZZ BASE NB-A-PCTBH CC 3.0 (Parallel cemented titan base hex and screw)	BSAB2483	15			10-12
SET ZZ BASE NB-A-PCTBH NP (Parallel cemented titan base hex and screw)	BSAB2481	25			10-12
SET ZZ BASE NB-A-PCTBH RP (Parallel cemented titan base hex and screw)	BSAB2482	35			10-12
SET ZZ BASE NB-A-PCTBH WP (Parallel cemented titan base hex and screw)	BSAB2484	35			10-12
SET ZZ BASE NB-A-PCTBH-G CC 3.0 (Parallel cemented titan base hex and screw Gold)	BSCB2483	15			10-12
SET ZZ BASE NB-A-PCTBH-G NP (Parallel cemented titan base hex and screw Gold)	BSCB2481	25			10-12
SET ZZ BASE NB-A-PCTBH-G RP (Parallel cemented titan base hex and screw Gold)	BSCB2482	35			10-12
SET ZZ BASE NB-A-PCTBH-G WP (Parallel cemented titan base hex and screw Gold)	BSCB2484	35			10-12
SET ZZ BASE NB-A-PCTBH-K85 CC 3.0 (Parallel cemented titan base hex K85 and screw)	BSAB2483K85	15			10-12
SET ZZ BASE NB-A-PCTBH-K85 NP (Parallel cemented titan base hex K85 and screw)	BSAB2481K85	25			10-12
SET ZZ BASE NB-A-PCTBH-K85 RP (Parallel cemented titan base hex K85 and screw)	BSAB2482K85	35			10-12
SET ZZ BASE NB-A-PCTBH-K85 WP (Parallel cemented titan base hex K85 and screw)	BSAB2484K85	35			10-12
SET ZZ BASE NB-A-PCTBH-L10-K85 WP (Parallel cemented titan base hex L10 K85 and screw)	BSAB2484L10K85	35			10-12
SET ZZ BASE NB-A-PCTBH-L20 CC 3.0 (Parallel cemented titan base hex L20 and screw)	BSAB2493L20	15			10-12
SET ZZ BASE NB-A-PCTBH-L20 NP (Parallel cemented titan base hex L20 and screw)	BSAB2491L20	25			10-12
SET ZZ BASE NB-A-PCTBH-L20 RP (Parallel cemented titan base hex L20 and screw)	BSAB2492L20	35			10-12

NOBEL BIOCARE® NOBELACTIVE®/NOBELREPLACE® CC/NOBEL PARALLEL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE NB-A-PCTBH-L20 WP (Parallel cemented titan base hex L20 and screw)	BSAB2494L20	35			10-12
SET ZZ BASE NB-A-PCTBH-L20-G CC 3.0 (Parallel cemented titan base hex L20 and screw Gold)	BSCB2493L20	15			10-12
SET ZZ BASE NB-A-PCTBH-L20-K85 CC 3.0 (Parallel cemented titan base hex L20 K85 and screw)	BSAB2483L20K85	15			10-12
SET ZZ BASE NB-A-PCTBH-L20-K85 NP (Parallel cemented titan base hex L20 K85 and screw)	BSAB2481L20K85	25			10-12
SET ZZ BASE NB-A-PCTBH-L20-K85 RP (Parallel cemented titan base hex L20 K85 and screw)	BSAB2482L20K85	35			10-12
SET ZZ BASE NB-A-PCTBH-L20-K85 WP (Parallel cemented titan base hex L20 K85 and screw)	BSAB2484L20K85	35			10-12
SET ZZ BASE NB-A-PCTBH-L30 CC 3.0 (Parallel cemented titan base hex L30 and screw)	BSAB2493L30	15			10-12
SET ZZ BASE NB-A-PCTBH-L30 NP (Parallel cemented titan base hex L30 and screw)	BSAB2491L30	25			10-12
SET ZZ BASE NB-A-PCTBH-L30 RP (Parallel cemented titan base hex L30 and screw)	BSAB2492L30	35			10-12
SET ZZ BASE NB-A-PCTBH-L30 WP (Parallel cemented titan base hex L30 and screw)	BSAB2494L30	35			10-12
SET ZZ BASE NB-A-PCTBH-L30-G CC 3.0 (Parallel cemented titan base hex L30 and screw Gold)	BSCB2493L30	15			10-12
SET ZZ BASE NB-A-PCTBH-L30-G NP (Parallel cemented titan base hex L30 and screw Gold)	BSCB2491L30	25			10-12
SET ZZ BASE NB-A-PCTBH-L30-G RP (Parallel cemented titan base hex L30 and screw Gold)	BSCB2492L30	35			10-12
SET ZZ BASE NB-A-PCTBH-L30-K85 CC 3.0 (Parallel cemented titan base hex L30 K85 and screw)	BSAB2483L30K85	15			10-12
SET ZZ BASE NB-A-PCTBH-L30-K85 NP (Parallel cemented titan base hex L30 K85 and screw)	BSAB2481L30K85	25			10-12
SET ZZ BASE NB-A-PCTBH-L30-K85 RP (Parallel cemented titan base hex L30 K85 and screw)	BSAB2482L30K85	35			10-12
SET ZZ BASE NB-A-PCTBH-L30-K85 WP (Parallel cemented titan base hex L30 K85 and screw)	BSAB2484L30K85	35			10-12
SET ZZ BASE NB-A-PCTBH-L40 NP (Parallel cemented titan base hex L40 and screw)	BSAB2491L40	25			10-12
SET ZZ BASE NB-A-PCTBH-L40 RP (Parallel cemented titan base hex L40 and screw)	BSAB2492L40	35			10-12
SET ZZ BASE NB-A-PCTBH-L40 WP (Parallel cemented titan base hex L40 and screw)	BSAB2494L40	35			10-12
SET ZZ BASE NB-A-PCTBH-L40-G NP (Parallel cemented titan base hex L40 and screw Gold)	BSCB2491L40	25			10-12
SET ZZ BASE NB-A-PCTBH-L40-G RP (Parallel cemented titan base hex L40 and screw Gold)	BSCB2492L40	35			10-12
SET ZZ BASE NB-A-PCTBH-L40-K85 CC 3.0 (Parallel cemented titan base hex L40 K85 and screw)	BSAB2483L40K85	35			10-12
SET ZZ BASE NB-A-PCTBH-L40-K85 RP (Parallel cemented titan base hex L40 K85 and screw)	BSAB2482L40K85	35			10-12
SET ZZ BASE NB-A-SA CC 3.0 (Scanalog and screw)	ZBDE2433	15			10-12
SET ZZ BASE NB-A-SA NP (Scanalog and screw)	ZBDE2431	15			10-12
SET ZZ BASE NB-A-SA RP (Scanalog and screw)	ZBDE2432	15			10-12
SET ZZ BASE NB-A-SA WP (Scanalog and screw)	ZBDE2434	15			10-12
SET ZZ BASE NB-A-SM CC 3.0 (Scanmarker and screw)	ZBAD2433	15			10-12

NOBEL BIO CARE® NOBEL ACTIVE®/NOBEL REPLACE® CC/NOBEL PARALLEL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE NB-A-SM NP (Scanmarker and screw)	ZBAD2431	25			10-12
SET ZZ BASE NB-A-SM RP (Scanmarker and screw)	ZBAD2432	35			10-12
SET ZZ BASE NB-A-SM WP (Scanmarker and screw)	ZBAD2434	35			10-12
SET ZZ BASE NB-A-WSM CC 3.0 (White Scanmarker and screw)	ZBCD2433	15			10-12
SET ZZ BASE NB-A-WSM RP (White Scanmarker and screw)	ZBCD2432	35			10-12
SET ZZ BASE NB-A-WSM WP (White Scanmarker and screw)	ZBCD2434	35			10-12
SET ZZ BASE NB-A-ZZMUA17 CC 3.0 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and im	BSAI2433	15			10-12
SET ZZ BASE NB-A-ZZMUA17 CC 3.0 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and im	BSAI2443	15			10-12
SET ZZ BASE NB-A-ZZMUA17 NP Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and implar	BSAI2431	20			10-12
SET ZZ BASE NB-A-ZZMUA17 NP Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and implar	BSAI2441	20			10-12
SET ZZ BASE NB-A-ZZMUA17 RP Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and implar	BSAI2432	30			10-12
SET ZZ BASE NB-A-ZZMUA17 RP Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and implar	BSAI2442	30			10-12
SET ZZ BASE NB-A-ZZMUA17 WP Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and impla	BSAI2434	30			10-12
SET ZZ BASE NB-A-ZZMUA17 WP Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and impla	BSAI2444	30			10-12
SET ZZ BASE NB-A-ZZMUA17-L10 CC 3.0 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L1	BSAI2433L10	15			10-12
SET ZZ BASE NB-A-ZZMUA17-L10 CC 3.0 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 L1	BSAI2443L10	15			10-12
SET ZZ BASE NB-A-ZZMUA17-L10 NP Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10 ai	BSAI2431L10	20			10-12
SET ZZ BASE NB-A-ZZMUA17-L10 NP Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 L10 ai	BSAI2441L10	20			10-12
SET ZZ BASE NB-A-ZZMUA17-L10 RP Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10 ai	BSAI2432L10	30			10-12
SET ZZ BASE NB-A-ZZMUA17-L10 RP Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 L10 ai	BSAI2442L10	30			10-12
SET ZZ BASE NB-A-ZZMUA17-L10 WP Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10 a	BSAI2434L10	30			10-12
SET ZZ BASE NB-A-ZZMUA17-L10 WP Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 L10 a	BSAI2444L10	30			10-12
SET ZZ BASE NB-A-ZZMUA-L15 CC 3.0 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw	BSAH2433L15	25	15		10-12
SET ZZ BASE NB-A-ZZMUA-L15 NP (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH2431L15	25	25		10-12
SET ZZ BASE NB-A-ZZMUA-L15 RP (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH2432L15	25	35		10-12
SET ZZ BASE NB-A-ZZMUA-L15 WP (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH2434L15	25	35		10-12
SET ZZ BASE NB-A-ZZMUA-L25 CC 3.0 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw	BSAH2433L25	25	15		10-12
SET ZZ BASE NB-A-ZZMUA-L25 NP (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH2431L25	25	25		10-12
SET ZZ BASE NB-A-ZZMUA-L25 RP (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH2432L25	25	35		10-12

NOBEL BIO CARE® NOBEL ACTIVE®/NOBEL REPLACE® CC/NOBEL PARALLEL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
SET ZZ BASE NB-A-ZZMUA-L25 WP (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH2434L25	25	35		10-12
SET ZZ BASE NB-A-ZZMUA-L35 CC 3.0 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH2433L35	25	15		10-12
SET ZZ BASE NB-A-ZZMUA-L35 NP (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH2431L35	25	25		10-12
SET ZZ BASE NB-A-ZZMUA-L35 RP (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH2432L35	25	35		10-12
SET ZZ BASE NB-A-ZZMUA-L35 WP (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH2434L35	25	35		10-12
SET ZZ BASE NB-A-ZZMUA-L45 CC 3.0 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH2433L45	25	15		10-12
SET ZZ BASE NB-A-ZZMUA-L45 NP (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH2431L45	25	25		10-12
SET ZZ BASE NB-A-ZZMUA-L45 RP (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH2432L45	25	35		10-12
SET ZZ BASE NB-A-ZZMUA-L45 WP (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH2434L45	25	35		10-12
ZZ BASE NB-A-ASCSC CC 3.0 (Angled screw-channel Screw)	BSKE2403	15			10-12
ZZ BASE NB-A-ASCSC NP (Angled screw-channel Screw)	BSKE2401	20			10-12
ZZ BASE NB-A-ASCSC RP (Angled screw-channel Screw)	BSKE2402	20			10-12
ZZ BASE NB-A-ASCSC WP (Angled screw-channel Screw)	BSKE2404	20			10-12
ZZ BASE NB-A-BSC CC 3.0 (Black Screw)	BSIE2403	15			10-12
ZZ BASE NB-A-BSC NP (Black Screw)	BSIE2401	25			10-12
ZZ BASE NB-A-BSC RP (Black Screw)	BSIE2402	35			10-12
ZZ BASE NB-A-BSC WP (Black Screw)	BSIE2404	35			10-12
ZZ BASE NB-A-BZSC CC 3.0 (Black Zirconia screw)	BSIE2433	15			10-12
ZZ BASE NB-A-BZSC NP (Black Zirconia screw)	BSIE2431	25			10-12
ZZ BASE NB-A-BZSC RP (Black Zirconia screw)	BSIE2432	35			10-12
ZZ BASE NB-A-BZSC WP (Black Zirconia screw)	BSIE2434	35			10-12
ZZ BASE NB-A-ICSC CC 3.0 (Impression coping Screw)	BSFE2403	15			10-12
ZZ BASE NB-A-ICSC NP (Impression coping Screw)	BSFE2401	15			10-12
ZZ BASE NB-A-ICSC RP (Impression coping Screw)	BSFE2402	15			10-12
ZZ BASE NB-A-ICSC WP (Impression coping Screw)	BSFE2404	15			10-12
ZZ BASE NB-A-LBSC CC 3.0 (Lab Screw)	BSHE2403	15			10-12
ZZ BASE NB-A-LBSC NP (Lab Screw)	BSHE2401	25			10-12
ZZ BASE NB-A-LBSC RP (Lab Screw)	BSHE2402	35			10-12
ZZ BASE NB-A-LBSC WP (Lab Screw)	BSHE2404	35			10-12

NOBEL BIO CARE® NOBEL ACTIVE®/NOBEL REPLACE® CC/NOBEL PARALLEL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
ZZ BASE NB-A-LBZSC CC 3.0 (Lab Zirconia screw)	BSHE2433	15			10-12
ZZ BASE NB-A-LBZSC NP (Lab Zirconia screw)	BSHE2431	25			10-12
ZZ BASE NB-A-LBZSC RP (Lab Zirconia screw)	BSHE2432	35			10-12
ZZ BASE NB-A-LBZSC WP (Lab Zirconia screw)	BSHE2434	35			10-12
ZZ BASE NB-A-MI CC 3.0 (Model implant)	MIAG2403				10-12
ZZ BASE NB-A-MI NP (Model implant)	MIAG2401				10-12
ZZ BASE NB-A-MI RP (Model implant)	MIAG2402				10-12
ZZ BASE NB-A-MI WP (Model implant)	MIAG2404				10-12
ZZ BASE NB-A-SC CC 3.0 (Screw)	BSAE2403	15			10-12
ZZ BASE NB-A-SC NP (Screw)	BSAE2401	25			10-12
ZZ BASE NB-A-SC RP (Screw)	BSAE2402	35			10-12
ZZ BASE NB-A-SC WP (Screw)	BSAE2404	35			10-12
ZZ BASE NB-A-SC-G CC 3.0 (Screw Gold)	BSCE2403	15			10-12
ZZ BASE NB-A-SC-G WP (Screw Gold)	BSCE2404	35			10-12
ZZ BASE NB-A-WMSM-L10 CC 3.0 (White Metal Scanmarker L10 with Screw)	ZBED2433L10	15			10-12
ZZ BASE NB-A-WMSM-L10 NP (White Metal Scanmarker L10 with Screw)	ZBED2431L10	15			10-12
ZZ BASE NB-A-WMSM-L10 RP (White Metal Scanmarker L10 with Screw)	ZBED2432L10	15			10-12
ZZ BASE NB-A-WMSM-L10 WP (White Metal Scanmarker L10 with Screw)	ZBED2434L10	15			10-12
ZZ BASE NB-A-WMSM-L12 CC 3.0 (White Metal Scanmarker L12 with Screw)	ZBED2433L12	15			10-12
ZZ BASE NB-A-WMSM-L12 NP (White Metal Scanmarker L12 with Screw)	ZBED2431L12	15			10-12
ZZ BASE NB-A-WMSM-L12 RP (White Metal Scanmarker L12 with Screw)	ZBED2432L12	15			10-12
ZZ BASE NB-A-WMSM-L12 WP (White Metal Scanmarker L12 with Screw)	ZBED2434L12	15			10-12
ZZ BASE NB-A-WMSM-L15 CC 3.0 (White Metal Scanmarker L15 with Screw)	ZBED2433L15	15			10-12
ZZ BASE NB-A-WMSM-L15 NP (White Metal Scanmarker L15 with Screw)	ZBED2431L15	15			10-12
ZZ BASE NB-A-WMSM-L15 RP (White Metal Scanmarker L15 with Screw)	ZBED2432L15	15			10-12
ZZ BASE NB-A-WMSM-L15 WP (White Metal Scanmarker L15 with Screw)	ZBED2434L15	15			10-12
ZZ BASE NB-A-ZSC CC 3.0 (Zirconia screw)	BSAE2433	15			10-12
ZZ BASE NB-A-ZSC NP (Zirconia screw)	BSAE2431	25			10-12
ZZ BASE NB-A-ZSC RP (Zirconia screw)	BSAE2432	35			10-12

NOBEL BIO CARE® NOBEL ACTIVE®/NOBEL REPLACE® CC/NOBEL PARALLEL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
ZZ BASE NB-A-ZSC WP (Zirconia screw)	BSAE2434	35			10-12
ZZ BASE NB-A-ZSC-G CC 3.0 (Zirconia Screw Gold)	BSCE2433	15			10-12
ZZ BASE NB-A-ZSC-G NP (Zirconia Screw Gold)	BSCE2431	25			10-12
ZZ BASE NB-A-ZSC-G RP (Zirconia Screw Gold)	BSCE2432	35			10-12
ZZ BASE NB-A-ZSC-G WP (Zirconia Screw Gold)	BSCE2434	35			10-12
ZZ BASE NB-A-ZZMUASC CC 3.0 (Zirkonzahn Multi Unit Abutment Screw)	BSEE2403	15			10-12
ZZ BASE NB-A-ZZMUASC NP (Zirkonzahn Multi Unit Abutment Screw)	BSEE2401	20			10-12
ZZ BASE NB-A-ZZMUASC RP (Zirkonzahn Multi Unit Abutment Screw)	BSEE2402	30			10-12
ZZ BASE NB-A-ZZMUASC WP (Zirkonzahn Multi Unit Abutment Screw)	BSEE2404	30			10-12

NOBEL BIO CARE® NOBEL REPLACE®/REPLACE SELECT TAPERED/NOBEL SPE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
SET ZZ BASE NB-R-ABCH 3.5 (Raw-Abutment® customizable hex and screw)	BSAF0391	35			10-12
SET ZZ BASE NB-R-ABCH 4.3 (Raw-Abutment® customizable hex and screw)	BSAF0392	35-45			10-12
SET ZZ BASE NB-R-ABCH 5.0 (Raw-Abutment® customizable hex and screw)	BSAF0393	35-45			10-12
SET ZZ BASE NB-R-ABCH 6.0 (Raw-Abutment® customizable hex and screw)	BSAF0394	35-45			10-12
SET ZZ BASE NB-R-ABCH-D14 3.5 (Raw-Abutment® customizable hex D14 and screw)	BSAF0391D14	35			10-12
SET ZZ BASE NB-R-ABCH-D14 4.3 (Raw-Abutment® customizable hex D14 and screw)	BSAF0392D14	35-45			10-12
SET ZZ BASE NB-R-ABCH-D14 5.0 (Raw-Abutment® customizable hex D14 and screw)	BSAF0393D14	35-45			10-12
SET ZZ BASE NB-R-ABCH-D14 6.0 (Raw-Abutment® customizable hex D14 and screw)	BSAF0394D14	35-45			10-12
SET ZZ BASE NB-R-ASCTBH-K80 3.5 (Angled screw-channel titan base hex K80 and screw)	BSAT0381K80	20			10-12
SET ZZ BASE NB-R-ASCTBH-K80 3.5 Type 2 (Angled screw-channel titan base hex K80 Type 2 and screw)	BSAT0391K80	20			10-12
SET ZZ BASE NB-R-ASCTBH-K80 4.3 (Angled screw-channel titan base hex K80 and screw)	BSAT0382K80	20			10-12
SET ZZ BASE NB-R-ASCTBH-K80 4.3 Type 2 (Angled screw-channel titan base hex K80 Type 2 and screw)	BSAT0392K80	20			10-12
SET ZZ BASE NB-R-ASCTB-K80 3.5 (Angled screw-channel titan base non hex K80 and screw)	BSAT0331K80	20			10-12
SET ZZ BASE NB-R-ASCTB-K80 4.3 (Angled screw-channel titan base non hex K80 and screw)	BSAT0332K80	20			10-12
SET ZZ BASE NB-R-CCTB 3.5 (Conical cemented titan base non hex and screw)	BSAD0331	35			10-12
SET ZZ BASE NB-R-CCTB 4.3 (Conical cemented titan base non hex and screw)	BSAD0332	35-45			10-12
SET ZZ BASE NB-R-CCTB 5.0 (Conical cemented titan base non hex and screw)	BSAD0333	35-45			10-12
SET ZZ BASE NB-R-CCTB 6.0 (Conical cemented titan base non hex and screw)	BSAD0334	35-45			10-12
SET ZZ BASE NB-R-CCTB-G 3.5 (Conical cemented titan base non hex and screw Gold)	BSCD0331	35			10-12
SET ZZ BASE NB-R-CCTB-G 4.3 (Conical cemented titan base non hex and screw Gold)	BSCD0332	35-45			10-12
SET ZZ BASE NB-R-CCTB-G 5.0 (Conical cemented titan base non hex and screw Gold)	BSCD0333	35-45			10-12
SET ZZ BASE NB-R-CCTB-G 6.0 (Conical cemented titan base non hex and screw Gold)	BSCD0334	35-45			10-12
SET ZZ BASE NB-R-CCTB-K85 3.5 (Conical cemented titan base non hex K85 and screw)	BSAD0331K85	35			10-12
SET ZZ BASE NB-R-CCTB-K85 4.3 (Conical cemented titan base non hex K85 and screw)	BSAD0332K85	35-45			10-12
SET ZZ BASE NB-R-CCTB-K85 5.0 (Conical cemented titan base non hex K85 and screw)	BSAD0333K85	35-45			10-12
SET ZZ BASE NB-R-CCTB-K85 6.0 (Conical cemented titan base non hex K85 and screw)	BSAD0334K85	35-45			10-12
SET ZZ BASE NB-R-CCTB-L10 3.5 (Conical cemented titan base non hex L10 and screw)	BSAD0331L10	35			10-12
SET ZZ BASE NB-R-CCTB-L10 4.3 (Conical cemented titan base non hex L10 and screw)	BSAD0332L10	35-45			10-12
SET ZZ BASE NB-R-CCTB-L10 5.0 (Conical cemented titan base non hex L10 and screw)	BSAD0333L10	35-45			10-12

NOBEL BIO CARE® NOBEL REPLACE®/REPLACE SELECT TAPERED/NOBEL SPE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
SET ZZ BASE NB-R-CCTB-L10 6.0 (Conical cemented titan base non hex L10 and screw)	BSAD0334L10	35-45			10-12
SET ZZ BASE NB-R-CCTB-L10-G 3.5 (Conical cemented titan base non hex L10 and screw Gold)	BSCD0331L10	35			10-12
SET ZZ BASE NB-R-CCTB-L10-G 4.3 (Conical cemented titan base non hex L10 and screw Gold)	BSCD0332L10	35-45			10-12
SET ZZ BASE NB-R-CCTB-L10-G 5.0 (Conical cemented titan base non hex L10 and screw Gold)	BSCD0333L10	35-45			10-12
SET ZZ BASE NB-R-CCTB-L10-G 6.0 (Conical cemented titan base non hex L10 and screw Gold)	BSCD0334L10	35-45			10-12
SET ZZ BASE NB-R-CCTB-L20 3.5 (Conical cemented titan base non hex L20 and screw)	BSAD0331L20	35			10-12
SET ZZ BASE NB-R-CCTB-L20 4.3 (Conical cemented titan base non hex L20 and screw)	BSAD0332L20	35-45			10-12
SET ZZ BASE NB-R-CCTB-L20 5.0 (Conical cemented titan base non hex L20 and screw)	BSAD0333L20	35-45			10-12
SET ZZ BASE NB-R-CCTB-L20 6.0 (Conical cemented titan base non hex L20 and screw)	BSAD0334L20	35-45			10-12
SET ZZ BASE NB-R-CCTB-L20-G 3.5 (Conical cemented titan base non hex L20 and screw Gold)	BSCD0331L20	35			10-12
SET ZZ BASE NB-R-CCTB-L20-G 4.3 (Conical cemented titan base non hex L20 and screw Gold)	BSCD0332L20	35-45			10-12
SET ZZ BASE NB-R-CCTB-L20-G 5.0 (Conical cemented titan base non hex L20 and screw Gold)	BSCD0333L20	35-45			10-12
SET ZZ BASE NB-R-CCTB-L30 3.5 (Conical cemented titan base non hex L30 and screw)	BSAD0331L30	35			10-12
SET ZZ BASE NB-R-CCTB-L30 4.3 (Conical cemented titan base non hex L30 and screw)	BSAD0332L30	35-45			10-12
SET ZZ BASE NB-R-CCTB-L30 5.0 (Conical cemented titan base non hex L30 and screw)	BSAD0333L30	35-45			10-12
SET ZZ BASE NB-R-CCTB-L30 6.0 (Conical cemented titan base non hex L30 and screw)	BSAD0334L30	35-45			10-12
SET ZZ BASE NB-R-CCTB-L30-G 3.5 (Conical cemented titan base non hex L30 and screw Gold)	BSCD0331L30	35			10-12
SET ZZ BASE NB-R-CCTB-L30-G 3.5 (Conical cemented titan base non hex L30 and screw Gold)	BSCD0334L20	35-45			10-12
SET ZZ BASE NB-R-CCTB-L30-G 4.3 (Conical cemented titan base non hex L30 and screw Gold)	BSCD0332L30	35-45			10-12
SET ZZ BASE NB-R-CCTB-L30-G 5.0 (Conical cemented titan base non hex L30 and screw Gold)	BSCD0333L30	35-45			10-12
SET ZZ BASE NB-R-CCTB-L30-G 6.0 (Conical cemented titan base non hex L30 and screw Gold)	BSCD0334L30	35-45			10-12
SET ZZ BASE NB-R-ICOT 3.5 (Impression coping open tray non hex and Screw and Driver)	BSAQ0361	15			10-12
SET ZZ BASE NB-R-ICOT 4.3 (Impression coping open tray non hex and Screw and Driver)	BSAQ0362	15			10-12
SET ZZ BASE NB-R-ICOT 6.0 (Impression coping open tray non hex and Screw and Driver)	BSAQ0364	15			10-12
SET ZZ BASE NB-R-ICOTH 3.5 (Impression coping open tray hex and Screw and Driver)	BSAQ0341	15			10-12
SET ZZ BASE NB-R-ICOTH 4.3 (Impression coping open tray hex and Screw and Driver)	BSAQ0342	15			10-12
SET ZZ BASE NB-R-ICOTH 5.0 (Impression coping open tray hex and Screw and Driver)	BSAQ0343	15			10-12
SET ZZ BASE NB-R-ICOTH 6.0 (Impression coping open tray hex and Screw and Driver)	BSAQ0344	15			10-12
SET ZZ BASE NB-R-NTB 3.5 (Narrow titan base non hex and screw)	BSAP0331	35			10-12

All information is subject to change. Errors and omissions excepted. Version: 21.09.2023

NOBEL BIO CARE® NOBEL REPLACE®/REPLACE SELECT TAPERED/NOBEL SPE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
SET ZZ BASE NB-R-NTB 4.3 (Narrow titan base non hex and screw)	BSAP0332	35-45			10-12
SET ZZ BASE NB-R-NTBH 3.5 (Narrow titan base hex and screw)	BSAP0381	35			10-12
SET ZZ BASE NB-R-NTBH 4.3 (Narrow titan base hex and screw)	BSAP0382	35-45			10-12
SET ZZ BASE NB-R-PCTBH 3.5 (Parallel cemented titan base hex and screw)	BSAB0381	35			10-12
SET ZZ BASE NB-R-PCTBH 3.5 (Parallel cemented titan base hex and screw)	BSAB0391	35			10-12
SET ZZ BASE NB-R-PCTBH 4.3 (Parallel cemented titan base hex and screw)	BSAB0392	35-45			10-12
SET ZZ BASE NB-R-PCTBH 5.0 (Parallel cemented titan base hex and screw)	BSAB0383	35-45			10-12
SET ZZ BASE NB-R-PCTBH 6.0 (Parallel cemented titan base hex and screw)	BSAB0384	35-45			10-12
SET ZZ BASE NB-R-PCTBH-G 3.5 (Parallel cemented titan base hex and screw Gold)	BSCB0381	35			10-12
SET ZZ BASE NB-R-PCTBH-G 4.3 (Parallel cemented titan base hex and screw Gold)	BSCB0382	35-45			10-12
SET ZZ BASE NB-R-PCTBH-G 4.3 (Parallel cemented titan base hex and screw Gold)	BSCB0392	35-45			10-12
SET ZZ BASE NB-R-PCTBH-G 5.0 (Parallel cemented titan base hex and screw Gold)	BSCB0383	35-45			10-12
SET ZZ BASE NB-R-PCTBH-G 6.0 (Parallel cemented titan base hex and screw Gold)	BSCB0384	35-45			10-12
SET ZZ BASE NB-R-PCTBH-K85 3.5 (Parallel cemented titan base hex K85 and screw)	BSAB0381K85	35			10-12
SET ZZ BASE NB-R-PCTBH-K85 4.3 (Parallel cemented titan base hex K85 and screw)	BSAB0382K85	35-45			10-12
SET ZZ BASE NB-R-PCTBH-K85 5.0 (Parallel cemented titan base hex K85 and screw)	BSAB0383K85	35-45			10-12
SET ZZ BASE NB-R-PCTBH-K85 6.0 (Parallel cemented titan base hex K85 and screw)	BSAB0384K85	35-45			10-12
SET ZZ BASE NB-R-PCTBH-L10 3.5 (Parallel cemented titan base hex L10 and screw)	BSAB0391L10	35			10-12
SET ZZ BASE NB-R-PCTBH-L10 4.3 (Parallel cemented titan base hex L10 and screw)	BSAB0392L10	35-45			10-12
SET ZZ BASE NB-R-PCTBH-L10 5.0 (Parallel cemented titan base hex L10 and screw)	BSAB0393L10	35-45			10-12
SET ZZ BASE NB-R-PCTBH-L10 6.0 (Parallel cemented titan base hex L10 and screw)	BSAB0394L10	35-45			10-12
SET ZZ BASE NB-R-PCTBH-L10-G 3.5 (Parallel cemented titan base hex L10 and screw Gold)	BSCB0391L10	35			10-12
SET ZZ BASE NB-R-PCTBH-L10-G 4.3 (Parallel cemented titan base hex L10 and screw Gold)	BSCB0392L10	35-45			10-12
SET ZZ BASE NB-R-PCTBH-L10-G 5.0 (Parallel cemented titan base hex L10 and screw Gold)	BSCB0393L10	35-45			10-12
SET ZZ BASE NB-R-PCTBH-L10-G 6.0 (Parallel cemented titan base hex L10 and screw Gold)	BSCB0394L10	35-45			10-12
SET ZZ BASE NB-R-PCTBH-L20 3.5 (Parallel cemented titan base hex L20 and screw)	BSAB0391L20	35			10-12
SET ZZ BASE NB-R-PCTBH-L20 4.3 (Parallel cemented titan base hex L20 and screw)	BSAB0392L20	35-45			10-12
SET ZZ BASE NB-R-PCTBH-L20 5.0 (Parallel cemented titan base hex L20 and screw)	BSAB0393L20	35-45			10-12
SET ZZ BASE NB-R-PCTBH-L20 6.0 (Parallel cemented titan base hex L20 and screw)	BSAB0394L20	35-45			10-12

NOBEL BIO CARE® NOBEL REPLACE®/REPLACE SELECT TAPERED/NOBEL SPE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
SET ZZ BASE NB-R-PCTBH-L20-G 3.5 (Parallel cemented titan base hex L20 and screw Gold)	BSCB0391L20	35			10-12
SET ZZ BASE NB-R-PCTBH-L20-G 4.3 (Parallel cemented titan base hex L20 and screw Gold)	BSCB0392L20	35-45			10-12
SET ZZ BASE NB-R-PCTBH-L20-G 5.0 (Parallel cemented titan base hex L20 and screw Gold)	BSCB0393L20	35-45			10-12
SET ZZ BASE NB-R-PCTBH-L20-G 6.0 (Parallel cemented titan base hex L20 and screw Gold)	BSCB0394L20	35-45			10-12
SET ZZ BASE NB-R-PCTBH-L30 3.5 (Parallel cemented titan base hex L30 and screw)	BSAB0391L30	35			10-12
SET ZZ BASE NB-R-PCTBH-L30 4.3 (Parallel cemented titan base hex L30 and screw)	BSAB0392L30	35-45			10-12
SET ZZ BASE NB-R-PCTBH-L30 5.0 (Parallel cemented titan base hex L30 and screw)	BSAB0393L30	35-45			10-12
SET ZZ BASE NB-R-PCTBH-L30 6.0 (Parallel cemented titan base hex L30 and screw)	BSAB0394L30	35-45			10-12
SET ZZ BASE NB-R-PCTBH-L30-G 3.5 (Parallel cemented titan base hex L30 and screw Gold)	BSCB0391L30	35			10-12
SET ZZ BASE NB-R-PCTBH-L30-G 4.3 (Parallel cemented titan base hex L30 and screw Gold)	BSCB0392L30	35-45			10-12
SET ZZ BASE NB-R-PCTBH-L30-G 5.0 (Parallel cemented titan base hex L30 and screw Gold)	BSCB0393L30	35-45			10-12
SET ZZ BASE NB-R-PCTBH-L30-G 6.0 (Parallel cemented titan base hex L30 and screw Gold)	BSCB0394L30	35-45			10-12
SET ZZ BASE NB-R-PCTBH-L40 3.5 (Parallel cemented titan base hex L40 and screw)	BSAB0391L40	35			10-12
SET ZZ BASE NB-R-PCTBH-L40 4.3 (Parallel cemented titan base hex L40 and screw)	BSAB0392L40	35-45			10-12
SET ZZ BASE NB-R-PCTBH-L40 5.0 (Parallel cemented titan base hex L40 and screw)	BSAB0393L40	35-45			10-12
SET ZZ BASE NB-R-PCTBH-L40-G 3.5 (Parallel cemented titan base hex L40 and screw Gold)	BSCB0391L40	35			10-12
SET ZZ BASE NB-R-PCTBH-L40-G 4.3 (Parallel cemented titan base hex L40 and screw Gold)	BSCB0392L40	35-45			10-12
SET ZZ BASE NB-R-PCTBH-L40-G 5.0 (Parallel cemented titan base hex L40 and screw Gold)	BSCB0393L40	35-45			10-12
SET ZZ BASE NB-R-SA 3.5 (Scanalog and screw)	ZBDE0331	15			10-12
SET ZZ BASE NB-R-SA 4.3 (Scanalog and screw)	ZBDE0332	15			10-12
SET ZZ BASE NB-R-SA 5.0 (Scanalog and screw)	ZBDE0333	15			10-12
SET ZZ BASE NB-R-SA 6.0 (Scanalog and screw)	ZBDE0334	15			10-12
SET ZZ BASE NB-R-SM 3.5 (Scanmarker and screw)	ZBAD0331	35			10-12
SET ZZ BASE NB-R-SM 4.3 (Scanmarker and screw)	ZBAD0332	35-45			10-12
SET ZZ BASE NB-R-SM 5.0 (Scanmarker and screw)	ZBAD0333	35-45			10-12
SET ZZ BASE NB-R-SM 6.0 (Scanmarker and screw)	ZBAD0334	35-45			10-12
SET ZZ BASE NB-R-WSM 3.5 (White Scanmarker and screw)	ZBCD0331	35			10-12
SET ZZ BASE NB-R-WSM 4.3 (White Scanmarker and screw)	ZBCD0332	35-45			10-12
SET ZZ BASE NB-R-WSM 5.0 (White Scanmarker and screw)	ZBCD0333	35-45			10-12

All information is subject to change. Errors and omissions excepted. Version: 21.09.2023

NOBEL BIOCARE® NOBELREPLACE®/REPLACE SELECT TAPERED/NOBELSPE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE NB-R-WSM 6.0 (White Scanmarker and screw)	ZBCD0334	35-45			10-12
SET ZZ BASE NB-R-ZZMUA17 3.5 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and implar	BSAI0331	25			10-12
SET ZZ BASE NB-R-ZZMUA17 3.5 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and implar	BSAI0341	25			10-12
SET ZZ BASE NB-R-ZZMUA17 4.3 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and implar	BSAI0332	25			10-12
SET ZZ BASE NB-R-ZZMUA17 4.3 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and implar	BSAI0342	25			10-12
SET ZZ BASE NB-R-ZZMUA17 5.0 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and implar	BSAI0333	30			10-12
SET ZZ BASE NB-R-ZZMUA17 5.0 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and implar	BSAI0343	30			10-12
SET ZZ BASE NB-R-ZZMUA17-L10 3.5 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10 a	BSAI0331L10	25			10-12
SET ZZ BASE NB-R-ZZMUA17-L10 3.5 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 L10 a	BSAI0341L10	25			10-12
SET ZZ BASE NB-R-ZZMUA17-L10 4.3 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10 a	BSAI0332L10	25			10-12
SET ZZ BASE NB-R-ZZMUA17-L10 4.3 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 L10 a	BSAI0342L10	25			10-12
SET ZZ BASE NB-R-ZZMUA17-L10 5.0 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10 a	BSAI0333L10	30			10-12
SET ZZ BASE NB-R-ZZMUA17-L10 5.0 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 L10 a	BSAI0343L10	30			10-12
SET ZZ BASE NB-R-ZZMUA-L05 4.3 (Zirkonzahn Multi Unit Abutment non hex L05 and Screw)	BSAH0332L05	25	35		10-12
SET ZZ BASE NB-R-ZZMUA-L05 5.0 (Zirkonzahn Multi Unit Abutment non hex L05 and Screw)	BSAH0333L05	25	35		10-12
SET ZZ BASE NB-R-ZZMUA-L15 3.5 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH0331L15	25	35		10-12
SET ZZ BASE NB-R-ZZMUA-L15 4.3 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH0332L15	25	35		10-12
SET ZZ BASE NB-R-ZZMUA-L15 5.0 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH0333L15	25	35		10-12
SET ZZ BASE NB-R-ZZMUA-L15 6.0 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH0334L15	25	35		10-12
SET ZZ BASE NB-R-ZZMUA-L25 3.5 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH0331L25	25	35		10-12
SET ZZ BASE NB-R-ZZMUA-L25 4.3 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH0332L25	25	35		10-12
SET ZZ BASE NB-R-ZZMUA-L25 5.0 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH0333L25	25	35		10-12
SET ZZ BASE NB-R-ZZMUA-L25 6.0 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH0334L25	25	35		10-12
SET ZZ BASE NB-R-ZZMUA-L35 3.5 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH0331L35	25	35		10-12
SET ZZ BASE NB-R-ZZMUA-L35 4.3 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH0332L35	25	35		10-12
SET ZZ BASE NB-R-ZZMUA-L35 5.0 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH0333L35	25	35		10-12
SET ZZ BASE NB-R-ZZMUA-L35 6.0 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH0334L35	25	35		10-12
SET ZZ BASE NB-R-ZZMUA-L45 3.5 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH0331L45	25	35		10-12
SET ZZ BASE NB-R-ZZMUA-L45 4.3 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH0332L45	25	35		10-12

NOBEL BIO CARE® NOBEL REPLACE®/REPLACE SELECT TAPERED/NOBEL SPE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
SET ZZ BASE NB-R-ZZMUA-L45 5.0 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH0333L45	25	35		10-12
SET ZZ BASE NB-R-ZZMUA-L45 6.0 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH0334L45	25	35		10-12
ZZ BASE NB-R-ASCSC 3.5 (Angled screw-channel Screw)	BSKE0301	20			10-12
ZZ BASE NB-R-ASCSC 4.3 (Angled screw-channel Screw)	BSKE0302	20			10-12
ZZ BASE NB-R-BSC 3.5 (Black Screw)	BSIE0301	35			10-12
ZZ BASE NB-R-BSC 4.3 (Black Screw)	BSIE0302	35-45			10-12
ZZ BASE NB-R-BSC 5.0 (Black Screw)	BSIE0303	35-45			10-12
ZZ BASE NB-R-BSC 6.0 (Black Screw)	BSIE0304	35-45			10-12
ZZ BASE NB-R-BZSC 3.5 (Black Zirconia screw)	BSIE0331	35			10-12
ZZ BASE NB-R-BZSC 4.3 (Black Zirconia screw)	BSIE0332	35-45			10-12
ZZ BASE NB-R-BZSC 5.0 (Black Zirconia screw)	BSIE0333	35-45			10-12
ZZ BASE NB-R-BZSC 6.0 (Black Zirconia screw)	BSIE0334	35-45			10-12
ZZ BASE NB-R-ICSC 3.5 (Impression coping Screw)	BSFE0301	15			10-12
ZZ BASE NB-R-ICSC 4.3 (Impression coping Screw)	BSFE0302	15			10-12
ZZ BASE NB-R-ICSC 5.0 (Impression coping Screw)	BSFE0303	15			10-12
ZZ BASE NB-R-ICSC 6.0 (Impression coping Screw)	BSFE0304	15			10-12
ZZ BASE NB-R-LBSC 3.5 (Lab Screw)	BSHE0301	35			10-12
ZZ BASE NB-R-LBSC 4.3 (Lab Screw)	BSHE0302	35-45			10-12
ZZ BASE NB-R-LBSC 5.0 (Lab Screw)	BSHE0303	35-45			10-12
ZZ BASE NB-R-LBSC 6.0 (Lab Screw)	BSHE0304	35-45			10-12
ZZ BASE NB-R-LBZSC 3.5 (Lab Zirconia screw)	BSHE0331	35			10-12
ZZ BASE NB-R-LBZSC 4.3 (Lab Zirconia screw)	BSHE0332	35-45			10-12
ZZ BASE NB-R-LBZSC 5.0 (Lab Zirconia screw)	BSHE0333	35-45			10-12
ZZ BASE NB-R-LBZSC 6.0 (Lab Zirconia screw)	BSHE0334	35-45			10-12
ZZ BASE NB-R-MI 3.5 (Model implant)	MIAG0301				10-12
ZZ BASE NB-R-MI 4.3 (Model implant)	MIAG0302				10-12
ZZ BASE NB-R-MI 5.0 (Model implant)	MIAG0303				10-12
ZZ BASE NB-R-MI 6.0 (Model implant)	MIAG0304				10-12
ZZ BASE NB-R-SC 3.5 (Screw)	BSAE0301	35			10-12

NOBEL BIOCARE® NOBELREPLACE®/REPLACE SELECT TAPERED/NOBELSPE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE NB-R-SC 4.3 (Screw)	BSAE0302	35-45			10-12
ZZ BASE NB-R-SC 5.0 (Screw)	BSAE0303	35-45			10-12
ZZ BASE NB-R-SC 6.0 (Screw)	BSAE0304	35-45			10-12
ZZ BASE NB-R-SC-G 3.5 (Screw Gold)	BSCE0301	35			10-12
ZZ BASE NB-R-SC-G 4.3 (Screw Gold)	BSCE0302	35-45			10-12
ZZ BASE NB-R-SC-G 5.0 (Screw Gold)	BSCE0303	35-45			10-12
ZZ BASE NB-R-SC-G 6.0 (Screw Gold)	BSCE0304	35-45			10-12
ZZ BASE NB-R-WMSM-L10 3.5 (White Metal Scanmarker L10 with Screw)	ZBED0331L10	15			10-12
ZZ BASE NB-R-WMSM-L10 4.3 (White Metal Scanmarker L10 with Screw)	ZBED0332L10	15			10-12
ZZ BASE NB-R-WMSM-L10 5.0 (White Metal Scanmarker L10 with Screw)	ZBED0333L10	15			10-12
ZZ BASE NB-R-WMSM-L10 6.0 (White Metal Scanmarker L10 with Screw)	ZBED0334L10	15			10-12
ZZ BASE NB-R-WMSM-L12 3.5 (White Metal Scanmarker L12 with Screw)	ZBED0331L12	15			10-12
ZZ BASE NB-R-WMSM-L12 4.3 (White Metal Scanmarker L12 with Screw)	ZBED0332L12	15			10-12
ZZ BASE NB-R-WMSM-L12 5.0 (White Metal Scanmarker L12 with Screw)	ZBED0333L12	15			10-12
ZZ BASE NB-R-WMSM-L12 6.0 (White Metal Scanmarker L12 with Screw)	ZBED0334L12	15			10-12
ZZ BASE NB-R-WMSM-L15 3.5 (White Metal Scanmarker L15 with Screw)	ZBED0331L15	15			10-12
ZZ BASE NB-R-WMSM-L15 4.3 (White Metal Scanmarker L15 with Screw)	ZBED0332L15	15			10-12
ZZ BASE NB-R-WMSM-L15 5.0 (White Metal Scanmarker L15 with Screw)	ZBED0333L15	15			10-12
ZZ BASE NB-R-WMSM-L15 6.0 (White Metal Scanmarker L15 with Screw)	ZBED0334L15	15			10-12
ZZ BASE NB-R-ZSC 3.5 (Zirconia screw)	BSAE0331	35			10-12
ZZ BASE NB-R-ZSC 4.3 (Zirconia screw)	BSAE0332	35-45			10-12
ZZ BASE NB-R-ZSC 5.0 (Zirconia screw)	BSAE0333	35-45			10-12
ZZ BASE NB-R-ZSC 6.0 (Zirconia screw)	BSAE0334	35-45			10-12
ZZ BASE NB-R-ZSC-G 3.5 (Zirconia screw Gold)	BSCE0331	35			10-12
ZZ BASE NB-R-ZSC-G 4.3 (Zirconia screw Gold)	BSCE0332	35-45			10-12
ZZ BASE NB-R-ZSC-G 5.0 (Zirconia screw Gold)	BSCE0333	35-45			10-12
ZZ BASE NB-R-ZSC-G 6.0 (Zirconia screw Gold)	BSCE0334	35-45			10-12
ZZ BASE NB-R-ZZMUASC 3.5 (Zirkonzahn Multi Unit Abutment Screw)	BSEE0301	25			10-12
ZZ BASE NB-R-ZZMUASC 4.3 (Zirkonzahn Multi Unit Abutment Screw)	BSEE0302	25			10-12

NOBEL BIOCARE® NOBELREPLACE®/REPLACE SELECT TAPERED/NOBELSPE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE NB-R-ZZMUASC 5.0 (Zirkonzahn Multi Unit Abutment Screw)	BSEE0303	30			10-12

NOBEL BIO CARE® NOBELZYGOMA

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE NB-Z-ABCH B4.1 (Raw-Abutment® customizable hex and screw)	BSAF9991	35			10-12
SET ZZ BASE NB-Z-CCTB B4.1 (Conical cemented titan base non hex and screw)	BSAD9931	35			10-12
SET ZZ BASE NB-Z-NTB B4.1 (Narrow titan base non hex and screw)	BSAP9931	35			10-12
SET ZZ BASE NB-Z-NTBH B4.1 (Narrow titan base hex and screw)	BSAP9981	35			10-12
SET ZZ BASE NB-Z-PCTBH B4.1 (Parallel cemented titan base hex and screw)	BSAB9981	35			10-12
SET ZZ BASE NB-Z-SM B4.1 (Scanmarker and screw)	ZBAD9931	35			10-12
SET ZZ BASE NB-Z-WSM B4.1 (White Scanmarker and screw)	ZBCD9931	35			10-12
ZZ BASE NB-Z-MI B4.1 (Model implant)	MIAG9901				10-12
ZZ BASE NB-Z-SC B4.1 (Screw)	BSAE9901	35			10-12
ZZ BASE NB-Z-ZSC B4.1 (Zirconia Screw)	BSAE9931	35			10-12

NOBEL BIO CARE® ON1™

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE NB-ON-CCTB NP (Conical cemented titan base non hex and screw)	BSAD12631	35			10-12
SET ZZ BASE NB-ON-CCTB RP (Conical cemented titan base non hex and screw)	BSAD12632	35			10-12
SET ZZ BASE NB-ON-CCTB WP (Conical cemented titan base non hex and screw)	BSAD12633	35			10-12
SET ZZ BASE NB-ON-ICOT NP (Impression coping open tray non hex and Screw and Driver)	BSAQ12641	15			10-12
SET ZZ BASE NB-ON-ICOT RP (Impression coping open tray non hex and Screw and Driver)	BSAQ12642	15			10-12
SET ZZ BASE NB-ON-ICOT WP (Impression coping open tray non hex and Screw and Driver)	BSAQ12643	15			10-12
SET ZZ BASE NB-ON-PCTBH NP (Parallel cemented titan base hex and screw)	BSAB12681	35			10-12
SET ZZ BASE NB-ON-PCTBH RP (Parallel cemented titan base hex and screw)	BSAB12682	35			10-12
SET ZZ BASE NB-ON-PCTBH RP (Parallel cemented titan base hex and screw)	BSAB12692	35			10-12
SET ZZ BASE NB-ON-PCTBH WP (Parallel cemented titan base hex and screw)	BSAB12683	35			10-12
SET ZZ BASE NB-ON-SA NP (Scananalog and screw)	ZBDE12631	15			10-12
SET ZZ BASE NB-ON-SA RP (Scananalog and screw)	ZBDE12632	15			10-12
SET ZZ BASE NB-ON-SA WP (Scananalog and screw)	ZBDE12633	15			10-12
SET ZZ BASE NB-ON-SM NP (Scanmarker and screw)	ZBAD12631	35			10-12
SET ZZ BASE NB-ON-SM RP (Scanmarker and screw)	ZBAD12632	35			10-12
SET ZZ BASE NB-ON-SM WP (Scanmarker and screw)	ZBAD12633	35			10-12
SET ZZ BASE NB-ON-WSM NP (White Scanmarker and screw)	ZBCD12631	35			10-12
SET ZZ BASE NB-ON-WSM RP (White Scanmarker and screw)	ZBCD12632	35			10-12
SET ZZ BASE NB-ON-WSM WP (White Scanmarker and screw)	ZBCD12633	35			10-12
ZZ BASE NB-ON-ICSC NP (Impression coping Screw)	BSFE12601	15			10-12
ZZ BASE NB-ON-ICSC RP (Impression coping Screw)	BSFE12602	15			10-12
ZZ BASE NB-ON-ICSC WP (Impression coping Screw)	BSFE12603	15			10-12
ZZ BASE NB-ON-MI NP (Model implant)	MIAG12601				10-12
ZZ BASE NB-ON-MI RP (Model implant)	MIAG12602				10-12
ZZ BASE NB-ON-SC NP (Screw)	BSAE12601	35			10-12
ZZ BASE NB-ON-SC RP (Screw)	BSAE12602	35			10-12
ZZ BASE NB-ON-SC WP (Screw)	BSAE12603	35			10-12
ZZ BASE NB-ON-WMSM-L12 NP (White Metal Scanmarker L12 with Screw)	ZBED12631L12	15			10-12
ZZ BASE NB-ON-WMSM-L12 RP (White Metal Scanmarker L12 with Screw)	ZBED12632L12	15			10-12

All information is subject to change. Errors and omissions excepted. Version: 21.09.2023

NOBEL BIOCARE® ON1™

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
ZZ BASE NB-ON-WMSM-L12 WP (White Metal Scanmarker L12 with Screw)	ZBED12633L12	15			10-12
ZZ BASE NB-ON-ZSC NP (Zirconia Screw)	BSAE12631	35			10-12
ZZ BASE NB-ON-ZSC RP (Zirconia Screw)	BSAE12632	35			10-12
ZZ BASE NB-ON-ZSC WP (Zirconia Screw)	BSAE12633	35			10-12

OSSTEM IMPLANT CONVERTIBLE ABUTMENT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE OS-GS-CA-CCTB 4.0 (Conical cemented titan base non hex and screw)	BSAD5531	30			10-12
SET ZZ BASE OS-GS-CA-CCTB 5.0 (Conical cemented titan base non hex and screw)	BSAD5532	30			10-12
SET ZZ BASE OS-GS-CA-CCTB 6.0 (Conical cemented titan base non hex and screw)	BSAD5533	30			10-12
SET ZZ BASE OS-GS-CA-ICOT 4.0 (Impression coping open tray non hex and Screw and Driver)	BSAQ5541	15			10-12
SET ZZ BASE OS-GS-CA-ICOT 5.0 (Impression coping open tray non hex and Screw and Driver)	BSAQ5542	15			10-12
SET ZZ BASE OS-GS-CA-ICOT 6.0 (Impression coping open tray non hex and Screw and Driver)	BSAQ5543	15			10-12
SET ZZ BASE OS-GS-CA-PCTBH 4.0 (Parallel cemented titan base hex and screw)	BSAB5581	30			10-12
SET ZZ BASE OS-GS-CA-PCTBH 5.0 (Parallel cemented titan base hex and screw)	BSAB5582	30			10-12
SET ZZ BASE OS-GS-CA-PCTBH 6.0 (Parallel cemented titan base hex and screw)	BSAB5583	30			10-12
SET ZZ BASE OS-GS-CA-SM 4.0 (Scanmarker and screw)	ZBAD5531	30			10-12
SET ZZ BASE OS-GS-CA-SM 5.0 (Scanmarker and screw)	ZBAD5532	30			10-12
SET ZZ BASE OS-GS-CA-SM 6.0 (Scanmarker and screw)	ZBAD5533	30			10-12
SET ZZ BASE OS-GS-SA 4.0 (Scananalog and screw)	ZBDE5531	15			10-12
SET ZZ BASE OS-GS-SA 5.0 (Scananalog and screw)	ZBDE5532	15			10-12
SET ZZ BASE OS-GS-SA 6.0 (Scananalog and screw)	ZBDE5533	15			10-12
ZZ BASE OS-GS-CA-ICSC 4.0 (Impression coping Screw)	BSFE5501	15			10-12
ZZ BASE OS-GS-CA-ICSC 5.0 (Impression coping Screw)	BSFE5502	15			10-12
ZZ BASE OS-GS-CA-ICSC 6.0 (Impression coping Screw)	BSFE5503	15			10-12
ZZ BASE OS-GS-CA-SC 4.0 (Screw)	BSAE5501	30			10-12
ZZ BASE OS-GS-CA-SC 5.0 (Screw)	BSAE5502	30			10-12
ZZ BASE OS-GS-CA-SC 6.0 (Screw)	BSAE5503	30			10-12
ZZ BASE OS-GS-CA-WMSM-L12 4.0 (White Metal Scanmarker L12 with Screw)	ZBED5531L12	15			10-12
ZZ BASE OS-GS-CA-WMSM-L12 5.0 (White Metal Scanmarker L12 with Screw)	ZBED5532L12	15			10-12
ZZ BASE OS-GS-CA-WMSM-L12 6.0 (White Metal Scanmarker L12 with Screw)	ZBED5533L12	15			10-12
ZZ BASE OS-GS-CA-ZSC 4.0 (Zirconia screw)	BSAE5531	30			10-12
ZZ BASE OS-GS-CA-ZSC 5.0 (Zirconia screw)	BSAE5532	30			10-12
ZZ BASE OS-GS-CA-ZSC 6.0 (Zirconia screw)	BSAE5533	30			10-12

OSSTEM IMPLANT GS/TS

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE OS-GS-ABC 4.0/4.5/5.0 (Raw-Abutment® customizable non hex and screw)	BSAF4231	30			10-12
SET ZZ BASE OS-GS-ABCH 3.5 (Raw-Abutment® customizable hex and screw)	BSAF4282	20			10-12
SET ZZ BASE OS-GS-ABCH 4.0/4.5/5.0 (Raw-Abutment® customizable hex and screw)	BSAF4281	30			10-12
SET ZZ BASE OS-GS-ABCH-D14 3.5 (Raw-Abutment® customizable hex D14 and screw)	BSAF4282D14	20			10-12
SET ZZ BASE OS-GS-ABCH-D14 4.0/4.5/5.0 (Raw-Abutment® customizable hex D14 and screw)	BSAF4281D14	30			10-12
SET ZZ BASE OS-GS-ASCTBH-K80 3.5 (Angled screw-channel titan base hex K80 and screw)	BSAT4282K80	20			10-12
SET ZZ BASE OS-GS-ASCTBH-K80 4.0/4.5/5.0 (Angled screw-channel titan base hex K80 and screw)	BSAT4281K80	20			10-12
SET ZZ BASE OS-GS-ASCTB-K80 3.5 (Angled screw-channel titan base non hex K80 and screw)	BSAT4232K80	20			10-12
SET ZZ BASE OS-GS-ASCTB-K80 4.0/4.5/5.0 (Angled screw-channel titan base non hex K80 and screw)	BSAT4231K80	20			10-12
SET ZZ BASE OS-GS-CCTB 3.5 (Conical cemented titan base non hex and screw)	BSAD4242	20			10-12
SET ZZ BASE OS-GS-CCTB 3.5 (Conical cemented titan base non hex and screw)	BSAD4232	20			10-12
SET ZZ BASE OS-GS-CCTB 4.0/4.5/5.0 (Conical cemented titan base non hex and screw)	BSAD4241	30			10-12
SET ZZ BASE OS-GS-CCTB-K85 3.5 (Conical cemented titan base non hex K85 and screw)	BSAD4232K85	20			10-12
SET ZZ BASE OS-GS-CCTB-K85 4.0/4.5/5.0 (Conical cemented titan base non hex K85 and screw)	BSAD4231K85	30			10-12
SET ZZ BASE OS-GS-CCTB-L10 3.5 (Conical cemented titan base non hex L10 and screw)	BSAD4232L10	20			10-12
SET ZZ BASE OS-GS-CCTB-L10 4.0/4.5/5.0 (Conical cemented titan base non hex L10 and screw)	BSAD4231L10	30			10-12
SET ZZ BASE OS-GS-CCTB-L20 3.5 (Conical cemented titan base non hex L20 and screw)	BSAD4232L20	20			10-12
SET ZZ BASE OS-GS-CCTB-L20 4.0/4.5/5.0 (Conical cemented titan base non hex L20 and screw)	BSAD4231L20	30			10-12
SET ZZ BASE OS-GS-CCTB-L30 3.5 (Conical cemented titan base non hex L30 and screw)	BSAD4232L30	20			10-12
SET ZZ BASE OS-GS-CCTB-L30 4.0/4.5/5.0 (Conical cemented titan base non hex L30 and screw)	BSAD4231L30	30			10-12
SET ZZ BASE OS-GS-CCTB-L40 3.5 (Conical cemented titan base non hex L40 and screw)	BSAD4232L40	20			10-12
SET ZZ BASE OS-GS-CCTB-L40 4.0/4.5/5.0 (Conical cemented titan base non hex L40 and screw)	BSAD4231L40	30			10-12
SET ZZ BASE OS-GS-ICOT 3.5 (Impression coping open tray non hex and Screw and Driver)	BSAQ4262	15			10-12
SET ZZ BASE OS-GS-ICOT 4.0/4.5/5.0 (Impression coping open tray non hex and Screw and Driver)	BSAQ4261	15			10-12
SET ZZ BASE OS-GS-ICOTH 3.5 (Impression coping open tray hex and Screw and Driver)	BSAQ4242	15			10-12
SET ZZ BASE OS-GS-ICOTH 4.0/4.5/5.0 (Impression coping open tray hex and Screw and Driver)	BSAQ4241	15			10-12
SET ZZ BASE OS-GS-NTB 3.5 (Narrow titan base non hex and screw)	BSAP4232	20			10-12
SET ZZ BASE OS-GS-NTB 4.0/4.5/5.0 (Narrow titan base non hex and screw)	BSAP4231	30			10-12
SET ZZ BASE OS-GS-NTBH 3.5 (Narrow titan base hex and screw)	BSAP4282	20			10-12

OSSTEM IMPLANT GS/TS

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE OS-GS-NTBH 4.0/4.5/5.0 (Narrow titan base hex and screw)	BSAP4281	30			10-12
SET ZZ BASE OS-GS-PCTBH 3.5 (Parallel cemented titan base hex and screw)	BSAB4282	20			10-12
SET ZZ BASE OS-GS-PCTBH 3.5 (Parallel cemented titan base hex and screw)	BSAB4292	20			10-12
SET ZZ BASE OS-GS-PCTBH 4.0/4.5/5.0 (Parallel cemented titan base hex and screw)	BSAB4281	30			10-12
SET ZZ BASE OS-GS-PCTBH 4.0/4.5/5.0 (Parallel cemented titan base hex and screw)	BSAB4291	30			10-12
SET ZZ BASE OS-GS-PCTBH-K85 3.5 (Parallel cemented titan base hex K85 and screw)	BSAB4282K85	20			10-12
SET ZZ BASE OS-GS-PCTBH-K85 4.0/4.5/5.0 (Parallel cemented titan base hex K85 and screw)	BSAB4281K85	30			10-12
SET ZZ BASE OS-GS-PCTBH-L10 3.5 (Parallel cemented titan base hex L10 and screw)	BSAB4282L10	20			10-12
SET ZZ BASE OS-GS-PCTBH-L10 3.5 (Parallel cemented titan base hex L10 and screw)	BSAB4292L10	20			10-12
SET ZZ BASE OS-GS-PCTBH-L10 4.0/4.5/5.0 (Parallel cemented titan base hex L10 and screw)	BSAB4291L10	30			10-12
SET ZZ BASE OS-GS-PCTBH-L10 4.0/4.5/5.0 (Parallel cemented titan base hex L10 and screw)	BSAB4281L10	30			10-12
SET ZZ BASE OS-GS-PCTBH-L20 3.5 (Parallel cemented titan base hex L20 and screw)	BSAB4292L20	20			10-12
SET ZZ BASE OS-GS-PCTBH-L20 4.0/4.5/5.0 (Parallel cemented titan base hex L20 and screw)	BSAB4291L20	30			10-12
SET ZZ BASE OS-GS-PCTBH-L30 3.5 (Parallel cemented titan base hex L30 and screw)	BSAB4292L30	20			10-12
SET ZZ BASE OS-GS-PCTBH-L30 4.0/4.5/5.0 (Parallel cemented titan base hex L30 and screw)	BSAB4291L30	30			10-12
SET ZZ BASE OS-GS-PCTBH-L40 3.5 (Parallel cemented titan base hex L40 and screw)	BSAB4292L40	20			10-12
SET ZZ BASE OS-GS-PCTBH-L40 4.0/4.5/5.0 (Parallel cemented titan base hex L40 and screw)	BSAB4291L40	30			10-12
SET ZZ BASE OS-GS-SA 3.5 (Scanalog and screw)	ZBDE4232	15			10-12
SET ZZ BASE OS-GS-SA 4.0/4.5/5.0 (Scanalog and screw)	ZBDE4231	15			10-12
SET ZZ BASE OS-GS-SM 3.5 (Scanmarker and screw)	ZBAD4232	20			10-12
SET ZZ BASE OS-GS-SM 4.0/4.5/5.0 (Scanmarker and screw)	ZBAD4231	30			10-12
SET ZZ BASE OS-GS-WSM 3.5 (White Scanmarker and screw)	ZBCD4232	20			10-12
SET ZZ BASE OS-GS-WSM 4.0/4.5/5.0 (White Scanmarker and screw)	ZBCD4231	30			10-12
SET ZZ BASE OS-GS-ZZMUA-L15 3.5 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH4232L15	25	20		10-12
SET ZZ BASE OS-GS-ZZMUA-L15 4.0/4.5/5.0 (Zirkonzahn Multi Unit Abutment non hex L15 and S	BSAH4231L15	25	30		10-12
SET ZZ BASE OS-GS-ZZMUA-L25 3.5 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH4232L25	25	20		10-12
SET ZZ BASE OS-GS-ZZMUA-L25 4.0/4.5/5.0 (Zirkonzahn Multi Unit Abutment non hex L25 and S	BSAH4231L25	25	30		10-12
SET ZZ BASE OS-GS-ZZMUA-L35 3.5 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH4232L35	25	20		10-12
SET ZZ BASE OS-GS-ZZMUA-L35 4.0/4.5/5.0 (Zirkonzahn Multi Unit Abutment non hex L35 and S	BSAH4231L35	25	30		10-12

All information is subject to change. Errors and omissions excepted. Version: 21.09.2023

OSSTEM IMPLANT GS/TS

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
SET ZZ BASE OS-GS-ZZMUA-L45 3.5 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH4232L45	25	20		10-12
SET ZZ BASE OS-GS-ZZMUA-L45 4.0/4.5/5.0 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH4231L45	25	30		10-12
ZZ BASE OS-GS-ASCSC 3.5 (Angled screw-channel Screw)	BSKE4202	20			10-12
ZZ BASE OS-GS-ASCSC 4.0/4.5/5.0 (Angled screw-channel Screw)	BSKE4201	20			10-12
ZZ BASE OS-GS-ICSC 3.5 (Impression coping Screw)	BSFE4202	15			10-12
ZZ BASE OS-GS-ICSC 4.0/4.5/5.0 (Impression coping Screw)	BSFE4201	15			10-12
ZZ BASE OS-GS-MI 3.5 (Model implant)	MIAG4202				10-12
ZZ BASE OS-GS-MI 4.0/4.5/5.0 (Model implant)	MIAG4201				10-12
ZZ BASE OS-GS-SC 3.5 (Screw)	BSAE4202	20			10-12
ZZ BASE OS-GS-SC 4.0/4.5/5.0 (Screw)	BSAE4201	30			10-12
ZZ BASE OS-GS-WMSM-L12 3.5 (White Metal Scanmarker L12 with Screw)	ZBED4232L12	15			10-12
ZZ BASE OS-GS-WMSM-L12 4.0/4.5/5.0 (White Metal Scanmarker L12 with Screw)	ZBED4231L12	15			10-12
ZZ BASE OS-GS-ZSC 3.5 (Zirconia screw)	BSAE4232	20			10-12
ZZ BASE OS-GS-ZSC 4.0/4.5/5.0 (Zirconia screw)	BSAE4231	30			10-12

OSSTEM IMPLANT GS/TS MULTI ABUTMENT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE OS-GS-CCTB MA (Conical cemented titan base non hex and screw)	BSAD9131	20			10-12
SET ZZ BASE OS-GS-CCTB-K85 MA (Conical cemented titan base non hex K85 and screw)	BSAD9131K85	20			10-12
SET ZZ BASE OS-GS-ICOT MA (Impression coping open tray non hex and Screw and Driver)	BSAQ9161	15			10-12
SET ZZ BASE OS-GS-SA MA (Scanalog and screw)	ZBDE9131	15			10-12
SET ZZ BASE OS-GS-SM MA (Scanmarker and screw)	ZBAD9131	20			10-12
SET ZZ BASE OS-GS-WSM MA (White Scanmarker and screw)	ZBCD9131	20			10-12
ZZ BASE OS-GS-ICSC MA (Impression coping Screw)	BSFE9101	15			10-12
ZZ BASE OS-GS-MA (Model implant)	MIAG9101				10-12
ZZ BASE OS-GS-SC MA (Screw)	BSAE9101	20			10-12
ZZ BASE OS-GS-WMSM-L12 MA (White Metal Scanmarker L12 with Screw)	ZBED9131L12	15			10-12
ZZ BASE OS-GS-ZSC MA (Zirconia Screw)	BSAE9131	20			10-12

OSSTEM IMPLANT US

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE OS-US-ABCH 3.5 (Raw-Abutment® customizable hex and screw)	BSAF4391	30			10-12
SET ZZ BASE OS-US-ABCH 4.1 (Raw-Abutment® customizable hex and screw)	BSAF4392	30			10-12
SET ZZ BASE OS-US-ABCH 5.1 (Raw-Abutment® customizable hex and screw)	BSAF4393	30			10-12
SET ZZ BASE OS-US-ABCH-D14 3.5 (Raw-Abutment® customizable hex D14 and screw)	BSAF4391D14	30			10-12
SET ZZ BASE OS-US-ABCH-D14 4.1 (Raw-Abutment® customizable hex D14 and screw)	BSAF4392D14	30			10-12
SET ZZ BASE OS-US-ABCH-D14 5.1 (Raw-Abutment® customizable hex D14 and screw)	BSAF4393D14	30			10-12
SET ZZ BASE OS-US-CCTB 3.5 (Conical cemented titan base non hex and screw)	BSAD4331	30			10-12
SET ZZ BASE OS-US-CCTB 4.1 (Conical cemented titan base non hex and screw)	BSAD4332	30			10-12
SET ZZ BASE OS-US-CCTB 5.1 (Conical cemented titan base non hex and screw)	BSAD4333	30			10-12
SET ZZ BASE OS-US-CCTB-L10 3.5 (Conical cemented titan base non hex L10 and screw)	BSAD4331L10	30			10-12
SET ZZ BASE OS-US-CCTB-L10 4.1 (Conical cemented titan base non hex L10 and screw)	BSAD4332L10	30			10-12
SET ZZ BASE OS-US-ICOTH 3.5 (Impression coping open tray hex and Screw and Driver)	BSAQ4341	15			10-12
SET ZZ BASE OS-US-ICOTH 4.1 (Impression coping open tray hex and Screw and Driver)	BSAQ4342	15			10-12
SET ZZ BASE OS-US-ICOTH 5.1 (Impression coping open tray hex and Screw and Driver)	BSAQ4343	15			10-12
SET ZZ BASE OS-US-NTB 3.5 (Narrow titan base non hex and screw)	BSAP4331	30			10-12
SET ZZ BASE OS-US-NTB 4.1 (Narrow titan base non hex and screw)	BSAP4332	30			10-12
SET ZZ BASE OS-US-NTBH 3.5 (Narrow titan base hex and screw)	BSAP4381	30			10-12
SET ZZ BASE OS-US-NTBH 4.1 (Narrow titan base hex and screw)	BSAP4382	30			10-12
SET ZZ BASE OS-US-PCTBH 3.5 (Parallel cemented titan base hex and screw)	BSAB4381	30			10-12
SET ZZ BASE OS-US-PCTBH 4.1 (Parallel cemented titan base hex and screw)	BSAB4382	30			10-12
SET ZZ BASE OS-US-PCTBH 5.1 (Parallel cemented titan base hex and screw)	BSAB4383	30			10-12
SET ZZ BASE OS-US-PCTBH-L10 3.5 (Parallel cemented titan base hex L10 and screw)	BSAB4391L10	30			10-12
SET ZZ BASE OS-US-PCTBH-L10 4.1 (Parallel cemented titan base hex L10 and screw)	BSAB4392L10	30			10-12
SET ZZ BASE OS-US-SA 3.5 (Scananalog and screw)	ZBDE4331	15			10-12
SET ZZ BASE OS-US-SA 4.1 (Scananalog and screw)	ZBDE4332	15			10-12
SET ZZ BASE OS-US-SA 5.1 (Scananalog and screw)	ZBDE4333	15			10-12
SET ZZ BASE OS-US-SM 3.5 (Scanmarker and screw)	ZBAD4331	30			10-12
SET ZZ BASE OS-US-SM 4.1 (Scanmarker and screw)	ZBAD4332	30			10-12
SET ZZ BASE OS-US-SM 5.1 (Scanmarker and screw)	ZBAD4333	30			10-12

OSSTEM IMPLANT US

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
SET ZZ BASE OS-US-WSM 3.5 (White Scanmarker and screw)	ZBCD4331	30			10-12
SET ZZ BASE OS-US-WSM 4.1 (White Scanmarker and screw)	ZBCD4332	30			10-12
SET ZZ BASE OS-US-WSM 5.1 (White Scanmarker and screw)	ZBCD4333	30			10-12
ZZ BASE OS-US-ICSC 3.5 (Impression coping Screw)	BSFE4301	15			10-12
ZZ BASE OS-US-ICSC 4.1 (Impression coping Screw)	BSFE4302	15			10-12
ZZ BASE OS-US-ICSC 5.1 (Impression coping Screw)	BSFE4303	15			10-12
ZZ BASE OS-US-SC 3.5 (Screw)	BSAE4301	30			10-12
ZZ BASE OS-US-SC 4.1 (Screw)	BSAE4302	30			10-12
ZZ BASE OS-US-SC 5.1 (Screw)	BSAE4303	30			10-12
ZZ BASE OS-US-WMSM-L12 3.5 (White Metal Scanmarker L12 with Screw)	ZBED4331L12	15			10-12
ZZ BASE OS-US-WMSM-L12 4.1 (White Metal Scanmarker L12 with Screw)	ZBED4332L12	15			10-12
ZZ BASE OS-US-WMSM-L12 5.1 (White Metal Scanmarker L12 with Screw)	ZBED4333L12	15			10-12
ZZ BASE OS-US-ZSC 3.5 (Zirconia screw)	BSAE4331	30			10-12
ZZ BASE OS-US-ZSC 4.1 (Zirconia screw)	BSAE4332	30			10-12
ZZ BASE OS-US-ZSC 5.1 (Zirconia screw)	BSAE4333	30			10-12

PALTOP® CONICAL ACTIVE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE PA-CA-ABCH 3.25/3.75/4.2/5.0 (Raw-Abutment® customizable hex and screw)	BSAF11081	25			10-12
SET ZZ BASE PA-CA-CCTB 3.25/3.75/4.2/5.0 (Conical cemented titan base non hex and screw)	BSAD11031	25			10-12
SET ZZ BASE PA-CA-NTB 3.25/3.75/4.2/5.0 (Narrow titan base non hex and screw)	BSAP11031	25			10-12
SET ZZ BASE PA-CA-NTBH 3.25/3.75/4.2/5.0 (Narrow titan base hex and screw)	BSAP11081	25			10-12
SET ZZ BASE PA-CA-PCTBH 3.25/3.75/4.2/5.0 (Parallel cemented titan base hex and screw)	BSAB11081	25			10-12
SET ZZ BASE PA-CA-SM 3.25/3.75/4.2/5.0 (Scanmarker and screw)	ZBAD11031	25			10-12
ZZ BASE PA-CA-MI 3.25/3.75/4.2/5.0 (Model implant)	MIAG11001				10-12
ZZ BASE PA-CA-SC 3.25/3.75/4.2/5.0 (Screw)	BSAE11001	25			10-12
ZZ BASE PA-CA-ZSC 3.25/3.75/4.2/5.0 (Zirconia screw)	BSAE11031	25			10-12

PALTOP® INTERNAL HEX CONNECTION

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE PA-IH CCTB 3.25 (Conical cemented titan base non hex and screw)	BSAD10931	25-35			10-12
SET ZZ BASE PA-IH CCTB 3.75/4.2/5.0 (Conical cemented titan base non hex and screw)	BSAD10932	30-35			10-12
SET ZZ BASE PA-IH CCTB 6.0 (Conical cemented titan base non hex and screw)	BSAD10933	30-35			10-12
SET ZZ BASE PA-IH-ABCH 3.25 (Raw-Abutment® customizable hex and screw)	BSAF10981	25-35			10-12
SET ZZ BASE PA-IH-ABCH 3.75/4.2/5.0 (Raw-Abutment® customizable hex and screw)	BSAF10982	30-35			10-12
SET ZZ BASE PA-IH-ABCH 6.0 (Raw-Abutment® customizable hex and screw)	BSAF10983	30-35			10-12
SET ZZ BASE PA-IH-NTB 3.25 (Narrow titan base non hex and screw)	BSAP10931	25-35			10-12
SET ZZ BASE PA-IH-NTBH 3.25 (Narrow titan base hex and screw)	BSAP10981	25-35			10-12
SET ZZ BASE PA-IH-PCTBH 3.25 (Parallel cemented titan base hex and screw)	BSAB10981	25-35			10-12
SET ZZ BASE PA-IH-PCTBH 3.75/4.2/5.0 (Parallel cemented titan base hex and screw)	BSAB10982	30-35			10-12
SET ZZ BASE PA-IH-PCTBH 6.0 (Parallel cemented titan base hex and screw)	BSAB10983	30-35			10-12
SET ZZ BASE PA-IH-SA 3.25 (Scananalog and screw)	ZBDE10931	15			10-12
SET ZZ BASE PA-IH-SA 3.75/4.2/5.0 (Scananalog and screw)	ZBDE10932	15			10-12
SET ZZ BASE PA-IH-SM 3.25 (Scanmarker and screw)	ZBAD10931	25-35			10-12
SET ZZ BASE PA-IH-SM 3.75/4.2/5.0 (Scanmarker and screw)	ZBAD10932	30-35			10-12
SET ZZ BASE PA-IH-SM 6.0 (Scanmarker and screw)	ZBAD10933	30-35			10-12
ZZ BASE PA-IH-ASCSC 3.25 (Angled screw-channel Screw)	BSKE10901	20			10-12
ZZ BASE PA-IH-ASCSC 3.75/4.2/5.0 (Angled screw-channel Screw)	BSKE10902	20			10-12
ZZ BASE PA-IH-MI 3.25 (Model implant)	MIAG10901				10-12
ZZ BASE PA-IH-MI 3.75/4.2/5.0 (Model implant)	MIAG10902				10-12
ZZ BASE PA-IH-SC 3.25 (Screw)	BSAE10901	25-35			10-12
ZZ BASE PA-IH-SC 3.75/4.2/5.0 (Screw)	BSAE10902	30-35			10-12
ZZ BASE PA-IH-SC 6.0 (Screw)	BSAE10903	30-35			10-12
ZZ BASE PA-IH-ZSC 3.25 (Zirconia screw)	BSAE10931	25-35			10-12
ZZ BASE PA-IH-ZSC 3.75/4.2/5.0 (Zirconia screw)	BSAE10932	30-35			10-12
ZZ BASE PA-IH-ZSC 6.0 (Zirconia screw)	BSAE10933	30-35			10-12

PALTOP® MULTI-UNIT ABUTMENT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE PA-MU-CCTB (Conical cemented titan base non hex and screw)	BSAD11131	30-35			10-12
SET ZZ BASE PA-MU-ICOT (Impression coping open tray non hex and Screw and Driver)	BSAQ11161	15			10-12
SET ZZ BASE PA-MU-SA (Scanalog and screw)	ZBDE11131	15			10-12
SET ZZ BASE PA-MU-SM (Scanmarker and screw)	ZBAD11131	30-35			10-12
ZZ BASE PA-MU-ICSC (Impression coping Screw)	BSFE11101	15			10-12
ZZ BASE PA-MU-MI (Model implant)	MIAG11101				10-12
ZZ BASE PA-MU-SC (Screw)	BSAE11101	30-35			10-12
ZZ BASE PA-MU-WMSM-L12 (White Metal Scanmarker L12 with Screw)	ZBED11131L12	15			10-12
ZZ BASE PA-MU-ZSC (Zirconia screw)	BSAE11131	30-35			10-12

PALTOP® SINGLE-UNIT ABUTMENT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE PA-SU-CCTB (Conical cemented titan base non hex and screw)	BSAD11231	30-35			10-12
SET ZZ BASE PA-SU-PCTBH (Parallel cemented titan base hex and screw)	BSAB11281	30-35			10-12
SET ZZ BASE PA-SU-SM (Scanmarker and screw)	ZBAD11231	30-35			10-12
ZZ BASE PA-SU-SC (Screw)	BSAE11201	30-35			10-12
ZZ BASE PA-SU-ZSC (Zirconia screw)	BSAE11231	30-35			10-12

PHIBO® TSA® ADVANCE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE P-TS-ABCH S3 (Raw-Abutment® customizable hex and screw)	BSAF8491	35			10-12
SET ZZ BASE P-TS-ABCH S4 (Raw-Abutment® customizable hex and screw)	BSAF8492	35			10-12
SET ZZ BASE P-TS-ABCH S5 (Raw-Abutment® customizable hex and screw)	BSAF8493	35			10-12
SET ZZ BASE P-TS-CCTB S3 (Conical cemented titan base non hex and screw)	BSAD8431	35			10-12
SET ZZ BASE P-TS-CCTB S4 (Conical cemented titan base non hex and screw)	BSAD8432	35			10-12
SET ZZ BASE P-TS-CCTB S5 (Conical cemented titan base non hex and screw)	BSAD8433	35			10-12
SET ZZ BASE P-TS-CCTB-K85 S3 (Conical cemented titan base non hex K85 and screw)	BSAD8431K85	35			10-12
SET ZZ BASE P-TS-CCTB-K85 S4 (Conical cemented titan base non hex K85 and screw)	BSAD8432K85	35			10-12
SET ZZ BASE P-TS-ICOTH S4 (Impression coping open tray hex and Screw and Driver)	BSAQ8442	15			10-12
SET ZZ BASE P-TS-ICOTH S5 (Impression coping open tray hex and Screw and Driver)	BSAQ8443	15			10-12
SET ZZ BASE P-TS-PCTBH S3 (Parallel cemented titan base hex and screw)	BSAB8481	35			10-12
SET ZZ BASE P-TS-PCTBH S4 (Parallel cemented titan base hex and screw)	BSAB8482	35			10-12
SET ZZ BASE P-TS-PCTBH S5 (Parallel cemented titan base hex and screw)	BSAB8483	35			10-12
SET ZZ BASE P-TS-SA S3 (Scanalog and screw)	ZBDE8431	15			10-12
SET ZZ BASE P-TS-SA S4 (Scanalog and screw)	ZBDE8432	15			10-12
SET ZZ BASE P-TS-SA S5 (Scanalog and screw)	ZBDE8433	15			10-12
SET ZZ BASE P-TS-SM S3 (Scanmarker and screw)	ZBAD8431	35			10-12
SET ZZ BASE P-TS-SM S4 (Scanmarker and screw)	ZBAD8432	35			10-12
SET ZZ BASE P-TS-SM S5 (Scanmarker and screw)	ZBAD8433	35			10-12
SET ZZ BASE P-TS-WSM S3 (White Scanmarker and screw)	ZBCD8431	35			10-12
SET ZZ BASE P-TS-WSM S4 (White Scanmarker and screw)	ZBCD8432	35			10-12
SET ZZ BASE P-TS-WSM S5 (White Scanmarker and screw)	ZBCD8433	35			10-12
ZZ BASE P-TS-ICSC S4 (Impression coping Screw)	BSFE8402	15			10-12
ZZ BASE P-TS-ICSC S5 (Impression coping Screw)	BSFE8403	15			10-12
ZZ BASE P-TS-SC S3 (Screw)	BSAE8401	35			10-12
ZZ BASE P-TS-SC S4 (Screw)	BSAE8402	35			10-12
ZZ BASE P-TS-SC S5 (Screw)	BSAE8403	35			10-12
ZZ BASE P-TS-WMSM-L12 S4 (White Metal Scanmarker L12 with Screw)	ZBED8432L12	15			10-12
ZZ BASE P-TS-WMSM-L12 S5 (White Metal Scanmarker L12 with Screw)	ZBED8433L12	15			10-12

PHIBO® TSA® ADVANCE

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
ZZ BASE P-TS-ZSC S3 (Zirconia Screw)	BSAE8431	35			10-12
ZZ BASE P-TS-ZSC S4 (Zirconia Screw)	BSAE8432	35			10-12
ZZ BASE P-TS-ZSC S5 (Zirconia Screw)	BSAE8433	35			10-12

PRE-MILLED ABUTMENT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
Straumann® NC Pre-milled Abut. Blank, Meden. Ø11.5, TAN	ZZ-022-0028	35	0		10-12
Straumann® NC Pre-milled Abut. Blank, Meden. Ø15.8, TAN	ZZ-022-0029	35	0		10-12
Straumann® RC Pre-milled Abut. Blank, Meden. Ø11.5, TAN	ZZ-022-0030	35	0		10-12
Straumann® RC Pre-milled Abut. Blank, Meden. Ø15.8, TAN	ZZ-022-0031	35	0		10-12
Straumann® RN Pre-milled Abut. Blank, Meden. Ø11.5, TAN	ZZ-022-0034	35	0		10-12
Straumann® RN Pre-milled Abut. Blank, Meden. Ø15.8, TAN	ZZ-022-0035	35	0		10-12
Straumann® WN Pre-milled Abut. Blank, Meden. Ø11.5, TAN	ZZ-022-0036	35	0		10-12
Straumann® WN Pre-milled Abut. Blank, Meden. Ø15.8, TAN	ZZ-022-0037	35	0		10-12

SGS DENTAL CONICAL PLATFORM

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE SG-CP-ABCH 2.1 Narrow (Raw-Abutment® customizable hex and screw)	BSAF12381	25			10-12
SET ZZ BASE SG-CP-ABCH 3.5 (Raw-Abutment® customizable hex and screw)	BSAF12382	25			10-12
SET ZZ BASE SG-CP-CCTB 2.1 Narrow (Conical cemented titan base non hex and screw)	BSAD12331	25			10-12
SET ZZ BASE SG-CP-CCTB 3.5 (Conical cemented titan base non hex and screw)	BSAD12332	25			10-12
SET ZZ BASE SG-CP-NTB 2.1 Narrow (Narrow titan base non hex and screw)	BSAP12331	25			10-12
SET ZZ BASE SG-CP-NTB 3.5 (Narrow titan base non hex and screw)	BSAP12332	25			10-12
SET ZZ BASE SG-CP-NTBH 2.1 Narrow (Narrow titan base hex and screw)	BSAP12381	25			10-12
SET ZZ BASE SG-CP-NTBH 3.5 (Narrow titan base hex and screw)	BSAP12382	25			10-12
SET ZZ BASE SG-CP-PCTBH 2.1 Narrow (Parallel cemented titan base hex and screw)	BSAB12381	25			10-12
SET ZZ BASE SG-CP-PCTBH 3.5 (Parallel cemented titan base hex and screw)	BSAB12382	25			10-12
SET ZZ BASE SG-CP-SM 2.1 Narrow (Scanmarker and screw)	ZBAD12331	25			10-12
SET ZZ BASE SG-CP-SM 3.5 (Scanmarker and screw)	ZBAD12332	25			10-12
ZZ BASE SG-CP-SC 2.1 Narrow (Screw)	BSAE12301	25			10-12
ZZ BASE SG-CP-SC 3.5 (Screw)	BSAE12302	25			10-12
ZZ BASE SG-CP-ZSC 2.1 Narrow (Zirconia Screw)	BSAE12331	25			10-12
ZZ BASE SG-CP-ZSC 3.5 (Zirconia Screw)	BSAE12332	25			10-12

SGS DENTAL INTERNAL HEXAGON

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE SG-IH-ABCH 3.75 (Raw-Abutment® customizable hex and screw)	BSAF12881	25			10-12
SET ZZ BASE SG-IH-CCTB 3.75 (Conical cemented titan base non hex and screw)	BSAD12831	25			10-12
SET ZZ BASE SG-IH-ICOTH 3.75 (Impression coping open tray hex and Screw and Driver)	BSAQ12841	15			10-12
SET ZZ BASE SG-IH-PCTBH 3.75 (Parallel cemented titan base hex and screw)	BSAB12881	25			10-12
SET ZZ BASE SG-IH-SM 3.75 (Scanmarker and screw)	ZBAD12831	25			10-12
ZZ BASE SG-IH-ICSC 3.75 (Impression coping Screw)	BSFE12801	15			10-12
ZZ BASE SG-IH-SC 3.75 (Screw)	BSAE12801	25			10-12
ZZ BASE SG-IH-WMSM-L12 3.75 (White Metal Scanmarker L12 with Screw)	ZBED12831L12	15			10-12
ZZ BASE SG-IH-ZSC 3.75 (Zirconia Screw)	BSAE12831	25			10-12

SGS DENTAL THE-ONE MULTI-UNIT ABUTMENTS

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE SG-MU-CCTB S16/S17 (Conical cemented titan base non hex and screw)	BSAD12432	25			10-12
SET ZZ BASE SG-MU-CCTB S16S/S17S (Conical cemented titan base non hex and screw)	BSAD12431	25			10-12
SET ZZ BASE SG-MU-ICOT S16/S17 (Impression coping open tray non hex and Screw and Driver)	BSAQ12462	15			10-12
SET ZZ BASE SG-MU-ICOT S16S/S17S (Impression coping open tray non hex and Screw and Driver)	BSAQ12461	15			10-12
SET ZZ BASE SG-MU-SM S16/S17 (Scanmarker and screw)	ZBAD12432	25			10-12
SET ZZ BASE SG-MU-SM S16S/S17S (Scanmarker and screw)	ZBAD12431	25			10-12
ZZ BASE SG-MU-ICSC S16/S17 (Impression coping Screw)	BSFE12402	15			10-12
ZZ BASE SG-MU-ICSC S16S/S17S (Impression coping Screw)	BSFE12401	15			10-12
ZZ BASE SG-MU-SC S16/S17 (Screw)	BSAE12402	25			10-12
ZZ BASE SG-MU-SC S16S/S17S (Screw)	BSAE12401	25			10-12
ZZ BASE SG-MU-WMSM-L12 S16/S17 (White Metal Scanmarker L12 with Screw)	ZBED12432L12	15			10-12
ZZ BASE SG-MU-WMSM-L12 S16S/S17S (White Metal Scanmarker L12 with Screw)	ZBED12431L12	15			10-12
ZZ BASE SG-MU-ZSC S16/S17 (Zirconia screw)	BSAE12432	25			10-12
ZZ BASE SG-MU-ZSC S16S/S17S (Zirconia screw)	BSAE12431	25			10-12

SIC® INVENT SICACE®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE S-S-ABCH 3.3 (Raw-Abutment® customizable hex and screw)	BSAF1291	20			10-12
SET ZZ BASE S-S-ABCH 4.2 (Raw-Abutment® customizable hex and screw)	BSAF1292	20			10-12
SET ZZ BASE S-S-ASCTBH-K80 3.3 (Angled screw-channel titan base hex K80 and screw)	BSAT1281K80	20			10-12
SET ZZ BASE S-S-ASCTBH-K80 4.2 (Angled screw-channel titan base hex K80 and screw)	BSAT1282K80	20			10-12
SET ZZ BASE S-S-ASCTB-K80 3.3 (Angled screw-channel titan base non hex K80 and screw)	BSAT1231K80	20			10-12
SET ZZ BASE S-S-ASCTB-K80 4.2 (Angled screw-channel titan base non hex K80 and screw)	BSAT1232K80	20			10-12
SET ZZ BASE S-S-CCTB 3.3 (Conical cemented titan base non hex and screw)	BSAD1231	20			10-12
SET ZZ BASE S-S-CCTB 4.2 (Conical cemented titan base non hex and screw)	BSAD1232	20			10-12
SET ZZ BASE S-S-CCTB-K85 3.3 (Conical cemented titan base non hex K85 and screw)	BSAD1231K85	20			10-12
SET ZZ BASE S-S-CCTB-K85 4.2 (Conical cemented titan base non hex K85 and screw)	BSAD1232K85	20			10-12
SET ZZ BASE S-S-CCTB-L10 3.3 (Conical cemented titan base non hex L10 and screw)	BSAD1231L10	20			10-12
SET ZZ BASE S-S-CCTB-L10 4.2 (Conical cemented titan base non hex L10 and screw)	BSAD1232L10	20			10-12
SET ZZ BASE S-S-CCTB-L20 3.3 (Conical cemented titan base non hex L20 and screw)	BSAD1231L20	20			10-12
SET ZZ BASE S-S-CCTB-L20 4.2 (Conical cemented titan base non hex L20 and screw)	BSAD1232L20	20			10-12
SET ZZ BASE S-S-CCTB-L30 3.3 (Conical cemented titan base non hex L30 and screw)	BSAD1231L30	20			10-12
SET ZZ BASE S-S-CCTB-L30 4.2 (Conical cemented titan base non hex L30 and screw)	BSAD1232L30	20			10-12
SET ZZ BASE S-S-CCTB-L40 3.3 (Conical cemented titan base non hex L40 and screw)	BSAD1231L40	20			10-12
SET ZZ BASE S-S-CCTB-L40 4.2 (Conical cemented titan base non hex L40 and screw)	BSAD1232L40	20			10-12
SET ZZ BASE S-S-ICOT 3.3 (Impression coping open tray non hex and Screw and Driver)	BSAQ1261	15			10-12
SET ZZ BASE S-S-ICOT 4.2 (Impression coping open tray non hex and Screw and Driver)	BSAQ1262	15			10-12
SET ZZ BASE S-S-ICOETH 3.3 (Impression coping open tray hex and Screw and Driver)	BSAQ1241	15			10-12
SET ZZ BASE S-S-ICOETH 4.2 (Impression coping open tray hex and Screw and Driver)	BSAQ1242	15			10-12
SET ZZ BASE S-S-NTB 3.3 (Narrow titan base non hex and screw)	BSAP1231	20			10-12
SET ZZ BASE S-S-NTBH 3.3 (Narrow titan base hex and screw)	BSAP1281	20			10-12
SET ZZ BASE S-S-PCTBH 3.3 (Parallel cemented titan base hex and screw)	BSAB1281	20			10-12
SET ZZ BASE S-S-PCTBH 4.2 (Parallel cemented titan base hex and screw)	BSAB1282	20			10-12
SET ZZ BASE S-S-PCTBH-K85 3.3 (Parallel cemented titan base hex K85 and screw)	BSAB1281K85	20			10-12
SET ZZ BASE S-S-PCTBH-K85 4.2 (Parallel cemented titan base hex K85 and screw)	BSAB1282K85	20			10-12
SET ZZ BASE S-S-PCTBH-L10 3.3 (Parallel cemented titan base hex L10 and screw)	BSAB1281L10	20			10-12

All information is subject to change. Errors and omissions excepted. Version: 21.09.2023

SIC® INVENT SICACE®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE S-S-PCTBH-L10 4.2 (Parallel cemented titan base hex L10 and screw)	BSAB1282L10	20			10-12
SET ZZ BASE S-S-PCTBH-L20 3.3 (Parallel cemented titan base hex L20 and screw)	BSAB1281L20	20			10-12
SET ZZ BASE S-S-PCTBH-L20 4.2 (Parallel cemented titan base hex L20 and screw)	BSAB1282L20	20			10-12
SET ZZ BASE S-S-PCTBH-L30 3.3 (Parallel cemented titan base hex L30 and screw)	BSAB1281L30	20			10-12
SET ZZ BASE S-S-PCTBH-L30 4.2 (Parallel cemented titan base hex L30 and screw)	BSAB1282L30	20			10-12
SET ZZ BASE S-S-PCTBH-L40 3.3 (Parallel cemented titan base hex L40 and screw)	BSAB1281L40	20			10-12
SET ZZ BASE S-S-PCTBH-L40 4.2 (Parallel cemented titan base hex L40 and screw)	BSAB1282L40	20			10-12
SET ZZ BASE S-S-SA 3.3 (Scananalog and screw)	ZBDE1231	15			10-12
SET ZZ BASE S-S-SA 4.2 (Scananalog and screw)	ZBDE1232	15			10-12
SET ZZ BASE S-S-SM 3.3 (Scanmarker and screw)	ZBAD1231	20			10-12
SET ZZ BASE S-S-SM 4.2 (Scanmarker and screw)	ZBAD1232	20			10-12
SET ZZ BASE S-S-WSM 3.3 (White Scanmarker and screw)	ZBCD1231	20			10-12
SET ZZ BASE S-S-WSM 4.2 (White Scanmarker and screw)	ZBCD1232	20			10-12
SET ZZ BASE S-S-ZZMUA-L15 3.3 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH1231L15	25	20		10-12
SET ZZ BASE S-S-ZZMUA-L15 4.2 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH1232L15	25	20		10-12
SET ZZ BASE S-S-ZZMUA-L25 3.3 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH1231L25	25	20		10-12
SET ZZ BASE S-S-ZZMUA-L25 4.2 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH1232L25	25	20		10-12
SET ZZ BASE S-S-ZZMUA-L35 3.3 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH1231L35	25	20		10-12
SET ZZ BASE S-S-ZZMUA-L35 4.2 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH1232L35	25	20		10-12
SET ZZ BASE S-S-ZZMUA-L45 3.3 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH1231L45	25	20		10-12
SET ZZ BASE S-S-ZZMUA-L45 4.2 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH1232L45	25	20		10-12
ZZ BASE S-S-ASCSC 3.3 (Angled screw-channel Screw)	BSKE1201	20			10-12
ZZ BASE S-S-ASCSC 4.2 (Angled screw-channel Screw)	BSKE1202	20			10-12
ZZ BASE S-S-ICSC 3.3 (Impression coping Screw)	BSFE1201	15			10-12
ZZ BASE S-S-ICSC 4.2 (Impression coping Screw)	BSFE1202	15			10-12
ZZ BASE S-S-SC 3.3 (Screw)	BSAE1201	20			10-12
ZZ BASE S-S-SC 4.2 (Screw)	BSAE1202	20			10-12
ZZ BASE S-S-SC-G 3.3 (Screw Gold)	BSCE1201	20			10-12
ZZ BASE S-S-SC-G 4.2 (Screw Gold)	BSCE1202	20			10-12

SIC® INVENT SICACE®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE S-S-WMSM-L12 3.3 (White Metal Scanmarker L12 with Screw)	ZBED1231L12	15			10-12
ZZ BASE S-S-WMSM-L12 4.2 (White Metal Scanmarker L12 with Screw)	ZBED1232L12	15			10-12
ZZ BASE S-S-ZSC 3.3 (Zirconia screw)	BSAE1231	20			10-12
ZZ BASE S-S-ZSC 4.2 (Zirconia screw)	BSAE1232	20			10-12
ZZ BASE S-S-ZSC-G 3.3 (Zirconia Screw Gold)	BSCE1231	20			10-12
ZZ BASE S-S-ZSC-G 4.2 (Zirconia Screw Gold)	BSCE1232	20			10-12

SOUTHERN IMPLANTS® DEEP CONICAL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE SI-DC-ABCH 3.0 (Raw-Abutment® customizable hex and screw)	BSAF9781	15			10-12
SET ZZ BASE SI-DC-ABCH 3.5/4.0 (Raw-Abutment® customizable hex and screw)	BSAF9782	20			10-12
SET ZZ BASE SI-DC-ABCH 5.0 (Raw-Abutment® customizable hex and screw)	BSAF9783	25-32			10-12
SET ZZ BASE SI-DC-CCTB 3.0 (Conical cemented titan base non hex and screw)	BSAD9731	15			10-12
SET ZZ BASE SI-DC-CCTB 3.5/4.0 (Conical cemented titan base non hex and screw)	BSAD9732	20			10-12
SET ZZ BASE SI-DC-CCTB 5.0 (Conical cemented titan base non hex and screw)	BSAD9733	25-32			10-12
SET ZZ BASE SI-DC-CCTB-K85 3.0 (Conical cemented titan base non hex K85 and screw)	BSAD9731K85	15			10-12
SET ZZ BASE SI-DC-CCTB-K85 3.5/4.0 (Conical cemented titan base non hex K85 and screw)	BSAD9732K85	20			10-12
SET ZZ BASE SI-DC-CCTB-K85 5.0 (Conical cemented titan base non hex K85 and screw)	BSAD9733K85	25-32			10-12
SET ZZ BASE SI-DC-PCTBH 3.0 (Parallel cemented titan base hex and screw)	BSAB9781	15			10-12
SET ZZ BASE SI-DC-PCTBH 3.5/4.0 (Parallel cemented titan base hex and screw)	BSAB9782	20			10-12
SET ZZ BASE SI-DC-PCTBH 5.0 (Parallel cemented titan base hex and screw)	BSAB9783	25-32			10-12
SET ZZ BASE SI-DC-PCTBH-K85 3.0 (Parallel cemented titan base hex K85 and screw)	BSAB9781K85	15			10-12
SET ZZ BASE SI-DC-PCTBH-K85 3.5/4.0 (Parallel cemented titan base hex K85 and screw)	BSAB9782K85	20			10-12
SET ZZ BASE SI-DC-PCTBH-K85 5.0 (Parallel cemented titan base hex K85 and screw)	BSAB9783K85	25-32			10-12
SET ZZ BASE SI-DC-SA 3.0 (Scananalog and screw)	ZBDE9731	15			10-12
SET ZZ BASE SI-DC-SA 3.5/4.0 (Scananalog and screw)	ZBDE9732	15			10-12
SET ZZ BASE SI-DC-SA 5.0 (Scananalog and screw)	ZBDE9733	15			10-12
SET ZZ BASE SI-DC-SM 3.0 (Scanmarker and screw)	ZBAD9731	15			10-12
SET ZZ BASE SI-DC-SM 3.5/4.0 (Scanmarker and screw)	ZBAD9732	20			10-12
SET ZZ BASE SI-DC-SM 5.0 (Scanmarker and screw)	ZBAD9733	25-32			10-12
SET ZZ BASE SI-DC-WSM 3.0 (White Scanmarker and screw)	ZBCD9731	15			10-12
SET ZZ BASE SI-DC-WSM 3.5/4.0 (White Scanmarker and screw)	ZBCD9732	20			10-12
SET ZZ BASE SI-DC-WSM 5.0 (White Scanmarker and screw)	ZBCD9733	25-32			10-12
ZZ BASE SI-DC-ASCSC 3.0 (Angled Screw Channel Screw)	BSKE9701	15			10-12
ZZ BASE SI-DC-MI 3.0 (Model implant)	MIAG9701				10-12
ZZ BASE SI-DC-MI 3.5/4.0 (Model implant)	MIAG9702				10-12
ZZ BASE SI-DC-MI 5.0 (Model implant)	MIAG9703				10-12
ZZ BASE SI-DC-SC 3.0 (Screw)	BSAE9701	15			10-12

SOUTHERN IMPLANTS® DEEP CONICAL

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE SI-DC-SC 3.5/4.0 (Screw)	BSAE9702	20			10-12
ZZ BASE SI-DC-SC 5.0 (Screw)	BSAE9703	25-32			10-12
ZZ BASE SI-DC-WMSM-L12 3.0 (White Metal Scanmarker L12 with Screw)	ZBED9731L12	15			10-12
ZZ BASE SI-DC-WMSM-L12 3.5/4.0 (White Metal Scanmarker L12 with Screw)	ZBED9732L12	15			10-12
ZZ BASE SI-DC-ZSC 3.0 (Zirconia Screw)	BSAE9731	15			10-12
ZZ BASE SI-DC-ZSC 3.5/4.0 (Zirconia Screw)	BSAE9732	20			10-12
ZZ BASE SI-DC-ZSC 5.0 (Zirconia Screw))	BSAE9733	25-32			10-12

SOUTHERN IMPLANTS® EXTERNAL HEX

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE SI-EH-ABCH 3.25 (Raw-Abutment® customizable hex and screw)	BSAF13091	32			10-12
SET ZZ BASE SI-EH-ABCH 4.0 (Raw-Abutment® customizable hex and screw)	BSAF13092	32			10-12
SET ZZ BASE SI-EH-CCTB 3.25 (Conical cemented titan base non hex and screw)	BSAD13031	32			10-12
SET ZZ BASE SI-EH-CCTB 4.0 (Conical cemented titan base non hex and screw)	BSAD13032	32			10-12
SET ZZ BASE SI-EH-CCTB-K85 3.25 (Conical cemented titan base non hex K85 and screw)	BSAD13031K85	32			10-12
SET ZZ BASE SI-EH-CCTB-K85 4.0 (Conical cemented titan base non hex K85 and screw)	BSAD13032K85	32			10-12
SET ZZ BASE SI-EH-CCTB-L10 3.25 (Conical cemented titan base non hex L10 and screw)	BSAD13031L10	32			10-12
SET ZZ BASE SI-EH-CCTB-L20 3.25 (Conical cemented titan base non hex L20 and screw)	BSAD13031L20	32			10-12
SET ZZ BASE SI-EH-CCTB-L30 3.25 (Conical cemented titan base non hex L30 and screw)	BSAD13031L30	32			10-12
SET ZZ BASE SI-EH-CCTB-L40 3.25 (Conical cemented titan base non hex L40 and screw)	BSAD13031L40	32			10-12
SET ZZ BASE SI-EH-ICOTH 3.25 (Impression coping open tray hex and Screw and Driver)	BSAQ13041	15			10-12
SET ZZ BASE SI-EH-ICOTH 4.0 (Impression coping open tray hex and Screw and Driver)	BSAQ13042	15			10-12
SET ZZ BASE SI-EH-NTB 3.25 (Narrow titan base non hex and screw)	BSAP13031	32			10-12
SET ZZ BASE SI-EH-NTBH 3.25 (Narrow titan base hex and screw)	BSAP13081	32			10-12
SET ZZ BASE SI-EH-PCTBH 3.25 (Parallel cemented titan base hex and screw)	BSAB13081	32			10-12
SET ZZ BASE SI-EH-PCTBH 4.0 (Parallel cemented titan base hex and screw)	BSAB13082	32			10-12
SET ZZ BASE SI-EH-PCTBH-L10 3.25 (Parallel cemented titan base hex L10 and screw)	BSAB13091L10	32			10-12
SET ZZ BASE SI-EH-PCTBH-L20 3.25 (Parallel cemented titan base hex L20 and screw)	BSAB13091L20	32			10-12
SET ZZ BASE SI-EH-PCTBH-L30 3.25 (Parallel cemented titan base hex L30 and screw)	BSAB13091L30	32			10-12
SET ZZ BASE SI-EH-PCTBH-L40 3.25 (Parallel cemented titan base hex L40 and screw)	BSAB13091L40	32			10-12
SET ZZ BASE SI-EH-SA 3.25 (Scananalog and screw)	ZBDE13031	15			10-12
SET ZZ BASE SI-EH-SA 4.0 (Scananalog and screw)	ZBDE13032	15			10-12
SET ZZ BASE SI-EH-SM 3.25 (Scanmarker and screw)	ZBAD13031	32			10-12
SET ZZ BASE SI-EH-SM 4.0 (Scanmarker and screw)	ZBAD13032	32			10-12
ZZ BASE SI-EH-ICSC 3.25 (Impression coping Screw)	BSFE13001	15			10-12
ZZ BASE SI-EH-ICSC 4.0 (Impression coping Screw)	BSFE13002	15			10-12
ZZ BASE SI-EH-MI 3.25 (Model implant)	MIAG13001				10-12
ZZ BASE SI-EH-MI 4.0 (Model implant)	MIAG13002				10-12
ZZ BASE SI-EH-SC 3.25 (Screw)	BSAE13001	32			10-12

SOUTHERN IMPLANTS® EXTERNAL HEX

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE SI-EH-SC 4.0 (Screw)	BSAE13002	32			10-12
ZZ BASE SI-EH-WMSM-L12 3.25 (White Metal Scanmarker L12 with Screw)	ZBED13031L12	15			10-12
ZZ BASE SI-EH-WMSM-L12 4.0 (White Metal Scanmarker L12 with Screw)	ZBED13032L12	15			10-12
ZZ BASE SI-EH-ZSC 3.25 (Zirconia Screw)	BSAE13031	32			10-12
ZZ BASE SI-EH-ZSC 4.0 (Zirconia Screw)	BSAE13032	32			10-12

SOUTHERN IMPLANTS® INTERNAL HEX

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE SI-IH-ABCH 3.7/4.2/5.0 (Raw-Abutment® customizable hex and screw)	BSAF9881	32			10-12
SET ZZ BASE SI-IH-CCTB 3.7/4.2/5.0 (Conical cemented titan base non hex and screw)	BSAD9831	32			10-12
SET ZZ BASE SI-IH-ICOTH 3.7/4.2/5.0 (Impression coping open tray hex and Screw and Driver)	BSAQ9841	15			10-12
SET ZZ BASE SI-IH-PCTBH 3.7/4.2/5.0 (Parallel cemented titan base hex and screw)	BSAB9881	32			10-12
SET ZZ BASE SI-IH-SM 3.7/4.2/5.0 (Scanmarker and screw)	ZBAD9831	32			10-12
SET ZZ BASE SI-IH-WSM 3.7/4.2/5.0 (White Scanmarker and screw)	ZBCD9831	32			10-12
ZZ BASE SI-IH-ICSC 3.7/4.2/5.0 (Impression coping Screw)	BSFE9801	15			10-12
ZZ BASE SI-IH-MI 3.7/4.2/5.0 (Model implant)	MIAG9801				10-12
ZZ BASE SI-IH-SC 3.7/4.2/5.0 (Screw)	BSAE9801	32			10-12
ZZ BASE SI-IH-WMSM-L12 3.7/4.2/5.0 (White Metal Scanmarker L12 with Screw)	ZBED9831L12	15			10-12
ZZ BASE SI-IH-ZSC 3.7/4.2/5.0 (Zirconia Screw)	BSAE9831	32			10-12

SOUTHERN IMPLANTS® IT CONNECTION

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
SET ZZ BASE SI-IC-ABCH 4.8 (Raw-Abutment® customizable hex and screw)	BSAF9681	25			10-12
SET ZZ BASE SI-IC-ABCH 6.5 (Raw-Abutment® customizable hex and screw)	BSAF9682	25			10-12
SET ZZ BASE SI-IC-CCTB 4.8 (Conical cemented titan base non hex and screw)	BSAD9631	25			10-12
SET ZZ BASE SI-IC-CCTB 6.5 (Conical cemented titan base non hex and screw)	BSAD9632	25			10-12
SET ZZ BASE SI-IC-ICOTH 4.8 (Impression coping open tray hex and Screw and Driver)	BSAQ9641	15			10-12
SET ZZ BASE SI-IC-ICOTH 6.5 (Impression coping open tray hex and Screw and Driver)	BSAQ9642	15			10-12
SET ZZ BASE SI-IC-PCTBH 4.8 (Parallel cemented titan base hex and screw)	BSAB9681	25			10-12
SET ZZ BASE SI-IC-PCTBH 6.5 (Parallel cemented titan base hex and screw)	BSAB9682	25			10-12
SET ZZ BASE SI-IC-SM 4.8 (Scanmarker and screw)	ZBAD9631	20			10-12
SET ZZ BASE SI-IC-SM 6.5 (Scanmarker and screw)	ZBAD9632	20			10-12
SET ZZ BASE SI-IC-WSM 4.8 (White Scanmarker and screw)	ZBCD9631	25			10-12
SET ZZ BASE SI-IC-WSM 6.5 (White Scanmarker and screw)	ZBCD9632	25			10-12
ZZ BASE SI-IC-ICSC 4.8 (Impression coping Screw)	BSFE9601	15			10-12
ZZ BASE SI-IC-ICSC 6.5 (Impression coping Screw)	BSFE9602	15			10-12
ZZ BASE SI-IC-MI 4.8 (Model implant)	MIAG9601				10-12
ZZ BASE SI-IC-SC 4.8 (Screw)	BSAE9601	25			10-12
ZZ BASE SI-IC-SC 6.5 (Screw)	BSAE9602	25			10-12
ZZ BASE SI-IC-WMSM-L12 4.8 (White Metal Scanmarker L12 with Screw)	ZBED9631L12	15			10-12
ZZ BASE SI-IC-WMSM-L12 6.5 (White Metal Scanmarker L12 with Screw)	ZBED9632L12	15			10-12
ZZ BASE SI-IC-ZSC 4.8 (Zirconia Screw)	BSAE9631	25			10-12
ZZ BASE SI-IC-ZSC 6.5 (Zirconia Screw)	BSAE9632	25			10-12

STRAUMANN® BONE LEVEL®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
Drill Guide for Straumann® Bone Level® implants RC 4.1 (for Rescue Kit)	WRAC0552				10-12
ScanAnalog for Straumann® Bone Level® implants RC 4.1 (for Rescue Kit Training)	WRAC0532				10-12
SET ZZ BASE ST-BL-ABCH NC 3.3 (Raw-Abutment® customizable hex and screw)	BSAF0581	35			10-12
SET ZZ BASE ST-BL-ABCH RC 4.1 (Raw-Abutment® customizable hex and screw)	BSAF0582	35			10-12
SET ZZ BASE ST-BL-ABCH SC (Raw-Abutment® customizable hex and screw)	BSAF0583	25			10-12
SET ZZ BASE ST-BL-ABCH-D14 NC 3.3 (Raw-Abutment® customizable hex D14 and screw)	BSAF0581D14	35			10-12
SET ZZ BASE ST-BL-ABCH-D14 RC 4.1 (Raw-Abutment® customizable hex D14 and screw)	BSAF0582D14	35			10-12
SET ZZ BASE ST-BL-ASCTBH-K80 NC 3.3 (Angled screw-channel titan base hex K80 and screw)	BSAT0581K80	20			10-12
SET ZZ BASE ST-BL-ASCTBH-K80 NC 3.3 Type 2 (Angled screw-channel titan base hex K80 Typ	BSAT0591K80	20			10-12
SET ZZ BASE ST-BL-ASCTBH-K80 RC 4.1 (Angled screw-channel titan base hex K80 and screw)	BSAT0582K80	20			10-12
SET ZZ BASE ST-BL-ASCTBH-K80 RC 4.1 Type 2 (Angled screw-channel titan base hex K80 Typ	BSAT0592K80	20			10-12
SET ZZ BASE ST-BL-ASCTBH-K80 SC (Angled screw-channel titan base hex K80 and screw)	BSAT0583K80	20			10-12
SET ZZ BASE ST-BL-ASCTBH-K80 SC Type 2 (Angled screw-channel titan base hex K80 Type 2	BSAT0593K80	20			10-12
SET ZZ BASE ST-BL-ASCTB-K80 NC 3.3 (Angled screw-channel titan base non hex K80 and scre	BSAT0531K80	20			10-12
SET ZZ BASE ST-BL-ASCTB-K80 RC 4.1 (Angled screw-channel titan base non hex K80 and scre	BSAT0532K80	20			10-12
SET ZZ BASE ST-BL-ASCTB-K80 SC (Angled screw-channel titan base non hex K80 and screw)	BSAT0533K80	20			10-12
SET ZZ BASE ST-BL-CCTB NC 3.3 (Conical cemented titan base non hex and screw)	BSAD0531	35			10-12
SET ZZ BASE ST-BL-CCTB RC 4.1 (Conical cemented titan base non hex and screw)	BSAD0532	35			10-12
SET ZZ BASE ST-BL-CCTB SC (Conical cemented titan base non hex and screw)	BSAD0533	25			10-12
SET ZZ BASE ST-BL-CCTB-G NC 3.3 (Conical cemented titan base non hex and screw Gold)	BSCD0531	35			10-12
SET ZZ BASE ST-BL-CCTB-G SC (Conical cemented titan base non hex and screw Gold)	BSCD0533	25			10-12
SET ZZ BASE ST-BL-CCTB-K85 NC 3.3 (Conical cemented titan base non hex K85 and screw)	BSAD0531K85	35			10-12
SET ZZ BASE ST-BL-CCTB-K85 RC 4.1 (Conical cemented titan base non hex K85 and screw)	BSAD0532K85	35			10-12
SET ZZ BASE ST-BL-CCTB-K85 SC (Conical cemented titan base non hex K85 and screw)	BSAD0533K85	25			10-12
SET ZZ BASE ST-BL-CCTB-L10 NC 3.3 (Conical cemented titan base non hex L10 and screw)	BSAD0531L10	35			10-12
SET ZZ BASE ST-BL-CCTB-L10 RC 4.1 (Conical cemented titan base non hex L10 and screw)	BSAD0532L10	35			10-12
SET ZZ BASE ST-BL-CCTB-L10-G NC 3.3 (Conical cemented titan base non hex L10 and screw C	BSCD0531L10	35			10-12
SET ZZ BASE ST-BL-CCTB-L10-G RC 4.1 (Conical cemented titan base non hex L10 and screw C	BSCD0532L10	35			10-12
SET ZZ BASE ST-BL-CCTB-L20 NC 3.3 (Conical cemented titan base non hex L20 and screw)	BSAD0531L20	35			10-12

STRAUMANN® BONE LEVEL®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE ST-BL-CCTB-L20 RC 4.1 (Conical cemented titan base non hex L20 and screw)	BSAD0532L20	35			10-12
SET ZZ BASE ST-BL-CCTB-L20-G NC 3.3 (Conical cemented titan base non hex L20 and screw C)	BSCD0531L20	35			10-12
SET ZZ BASE ST-BL-CCTB-L20-G RC 4.1 (Conical cemented titan base non hex L20 and screw C)	BSCD0532L20	35			10-12
SET ZZ BASE ST-BL-CCTB-L30 NC 3.3 (Conical cemented titan base non hex L30 and screw)	BSAD0531L30	35			10-12
SET ZZ BASE ST-BL-CCTB-L30 RC 4.1 (Conical cemented titan base non hex L30 and screw)	BSAD0532L30	35			10-12
SET ZZ BASE ST-BL-CCTB-L30-G NC 3.3 (Conical cemented titan base non hex L30 and screw C)	BSCD0531L30	35			10-12
SET ZZ BASE ST-BL-CCTB-L30-G RC 4.1 (Conical cemented titan base non hex L30 and screw C)	BSCD0532L30	35			10-12
SET ZZ BASE ST-BL-ICOT NC 3.3 (Impression coping open tray non hex and Screw and Driver)	BSAQ0561	15			10-12
SET ZZ BASE ST-BL-ICOT RC 4.1 (Impression coping open tray non hex and Screw and Driver)	BSAQ0562	15			10-12
SET ZZ BASE ST-BL-ICOT SC (Impression coping open tray non hex and Screw and Driver)	BSAQ0563	15			10-12
SET ZZ BASE ST-BL-ICOTH NC 3.3 (Impression coping open tray hex and Screw and Driver)	BSAQ0541	15			10-12
SET ZZ BASE ST-BL-ICOTH RC 4.1 (Impression coping open tray hex and Screw and Driver)	BSAQ0542	15			10-12
SET ZZ BASE ST-BL-ICOTH SC (Impression coping open tray hex and Screw and Driver)	BSAQ0543	15			10-12
SET ZZ BASE ST-BL-NTB NC 3.3 (Narrow titan base non hex and screw)	BSAP0531	35			10-12
SET ZZ BASE ST-BL-NTB RC 4.1 (Narrow titan base non hex and screw)	BSAP0532	35			10-12
SET ZZ BASE ST-BL-NTB SC (Narrow titan base non hex and screw)	BSAP0533	25			10-12
SET ZZ BASE ST-BL-NTBH NC 3.3 (Narrow titan base hex and screw)	BSAP0581	35			10-12
SET ZZ BASE ST-BL-NTBH RC 4.1 (Narrow titan base hex and screw)	BSAP0582	35			10-12
SET ZZ BASE ST-BL-NTBH SC (Narrow titan base hex and screw)	BSAP0583	25			10-12
SET ZZ BASE ST-BL-PCTBH NC 3.3 (Parallel cemented titan base hex and screw)	BSAB0591	35			10-12
SET ZZ BASE ST-BL-PCTBH RC 4.1 (Parallel cemented titan base hex and screw)	BSAB0592	35			10-12
SET ZZ BASE ST-BL-PCTBH SC (Parallel cemented titan base hex and screw)	BSAB0583	25			10-12
SET ZZ BASE ST-BL-PCTBH-G NC 3.3 (Parallel cemented titan base hex and screw Gold)	BSCB0591	35			10-12
SET ZZ BASE ST-BL-PCTBH-G SC (Parallel cemented titan base hex and screw Gold)	BSCB0583	25			10-12
SET ZZ BASE ST-BL-PCTBH-K85 NC 3.3 (Parallel cemented titan base hex K85 and screw)	BSAB0581K85	35			10-12
SET ZZ BASE ST-BL-PCTBH-K85 RC 4.1 (Parallel cemented titan base hex K85 and screw)	BSAB0582K85	35			10-12
SET ZZ BASE ST-BL-PCTBH-L20 NC 3.3 (Parallel cemented titan base hex L20 and screw)	BSAB0591L20	35			10-12
SET ZZ BASE ST-BL-PCTBH-L20 RC 4.1 (Parallel cemented titan base hex L20 and screw)	BSAB0592L20	35			10-12
SET ZZ BASE ST-BL-PCTBH-L20-G NC 3.3 (Parallel cemented titan base hex L20 and screw Gold)	BSCB0591L20	35			10-12

STRAUMANN® BONE LEVEL®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE ST-BL-PCTBH-L20-G RC 4.1 (Parallel cemented titan base hex L20 and screw Gold)	BSCB0592L20	35			10-12
SET ZZ BASE ST-BL-PCTBH-L30 NC 3.3 (Parallel cemented titan base hex L30 and screw)	BSAB0591L30	35			10-12
SET ZZ BASE ST-BL-PCTBH-L30 RC 4.1 (Parallel cemented titan base hex L30 and screw)	BSAB0592L30	35			10-12
SET ZZ BASE ST-BL-PCTBH-L30-G NC 3.3 (Parallel cemented titan base hex L30 and screw Gold)	BSCB0591L30	35			10-12
SET ZZ BASE ST-BL-PCTBH-L30-G RC 4.1 (Parallel cemented titan base hex L30 and screw Gold)	BSCB0592L30	35			10-12
SET ZZ BASE ST-BL-PCTBH-L40 NC 3.3 (Parallel cemented titan base hex L40 and screw)	BSAB0591L40	35			10-12
SET ZZ BASE ST-BL-PCTBH-L40 RC 4.1 (Parallel cemented titan base hex L40 and screw)	BSAB0592L40	35			10-12
SET ZZ BASE ST-BL-PCTBH-L40-G NC 3.3 (Parallel cemented titan base hex L40 and screw Gold)	BSCB0591L40	35			10-12
SET ZZ BASE ST-BL-PCTBH-L40-G RC 4.1 (Parallel cemented titan base hex L40 and screw Gold)	BSCB0592L40	35			10-12
SET ZZ BASE ST-BL-SA NC 3.3 (Scananalog and screw)	ZBDE0531	15			10-12
SET ZZ BASE ST-BL-SA RC 4.1 (Scananalog and screw)	ZBDE0532	15			10-12
SET ZZ BASE ST-BL-SA SC (Scananalog and screw)	ZBDE0533	15			10-12
SET ZZ BASE ST-BL-SM NC 3.3 (Scanmarker and screw)	ZBAD0531	35			10-12
SET ZZ BASE ST-BL-SM RC 4.1 (Scanmarker and screw)	ZBAD0532	35			10-12
SET ZZ BASE ST-BL-SM SC (Scanmarker and screw)	ZBAD0533	25			10-12
SET ZZ BASE ST-BL-ZZMUA17 NC 3.3 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and implant)	BSAI0531	20			10-12
SET ZZ BASE ST-BL-ZZMUA17 NC 3.3 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and implant)	BSAI0541	20			10-12
SET ZZ BASE ST-BL-ZZMUA17 RC 4.1 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and implant)	BSAI0532	20			10-12
SET ZZ BASE ST-BL-ZZMUA17 RC 4.1 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and implant)	BSAI0542	20			10-12
SET ZZ BASE ST-BL-ZZMUA17 SC Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and implant)	BSAI0533	15			10-12
SET ZZ BASE ST-BL-ZZMUA17 SC Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and implant)	BSAI0543	15			10-12
SET ZZ BASE ST-BL-ZZMUA17-L10 NC 3.3 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L)	BSAI0531L10	20			10-12
SET ZZ BASE ST-BL-ZZMUA17-L10 NC 3.3 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 L)	BSAI0541L10	20			10-12
SET ZZ BASE ST-BL-ZZMUA17-L10 RC 4.1 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L)	BSAI0532L10	20			10-12
SET ZZ BASE ST-BL-ZZMUA17-L10 RC 4.1 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 L)	BSAI0542L10	20			10-12
SET ZZ BASE ST-BL-ZZMUA17-L10 SC Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10)	BSAI0533L10	15			10-12
SET ZZ BASE ST-BL-ZZMUA17-L10 SC Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 L10)	BSAI0543L10	15			10-12
SET ZZ BASE ST-BL-ZZMUA30 SC Type 1 (Zirkonzahn Multi Unit Abutment 30° Type 1 and implant)	BSAK0533	15			10-12
SET ZZ BASE ST-BL-ZZMUA30 SC Type 2 (Zirkonzahn Multi Unit Abutment 30° Type 2 and implant)	BSAK0543	15			10-12

STRAUMANN® BONE LEVEL®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE ST-BL-ZZMUA30-L10 SC Type 1 (Zirkonzahn Multi Unit Abutment 30° Type 1 L10)	BSAK0533L10	15			10-12
SET ZZ BASE ST-BL-ZZMUA30-L10 SC Type 2 (Zirkonzahn Multi Unit Abutment 30° Type 2 L10)	BSAK0543L10	15			10-12
SET ZZ BASE ST-BL-ZZMUA-L15 NC 3.3 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH0531L15	25	30		10-12
SET ZZ BASE ST-BL-ZZMUA-L15 RC 4.1 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH0532L15	25	30		10-12
SET ZZ BASE ST-BL-ZZMUA-L15 SC (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH0533L15	25	25		10-12
SET ZZ BASE ST-BL-ZZMUA-L25 NC 3.3 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH0531L25	25	30		10-12
SET ZZ BASE ST-BL-ZZMUA-L25 RC 4.1 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH0532L25	25	30		10-12
SET ZZ BASE ST-BL-ZZMUA-L25 SC (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH0533L25	25	25		10-12
SET ZZ BASE ST-BL-ZZMUA-L35 NC 3.3 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH0531L35	25	30		10-12
SET ZZ BASE ST-BL-ZZMUA-L35 RC 4.1 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH0532L35	25	30		10-12
SET ZZ BASE ST-BL-ZZMUA-L35 SC (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH0533L35	25	25		10-12
SET ZZ BASE ST-BL-ZZMUA-L45 NC 3.3 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH0531L45	25	30		10-12
SET ZZ BASE ST-BL-ZZMUA-L45 RC 4.1 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH0532L45	25	30		10-12
SET ZZ BASE ST-BL-ZZMUA-L45 SC (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH0533L45	25	25		10-12
ZZ BASE ST-BL-ASCSC NC 3.3 (Angled screw-channel Screw)	BSKE0501	20			10-12
ZZ BASE ST-BL-ASCSC RC 4.1 (Angled screw-channel Screw)	BSKE0502	20			10-12
ZZ BASE ST-BL-ASCSC SC (Angled screw-channel Screw)	BSKE0503	20			10-12
ZZ BASE ST-BL-BSC NC 3.3 (Black Screw)	BSIE0501	35			10-12
ZZ BASE ST-BL-BSC RC 4.1 (Black Screw)	BSIE0502	35			10-12
ZZ BASE ST-BL-BSC SC (Black Screw)	BSIE0503	25			10-12
ZZ BASE ST-BL-BZSC NC 3.3 (Black Zirconia screw)	BSIE0531	35			10-12
ZZ BASE ST-BL-BZSC RC 4.1 (Black Zirconia screw)	BSIE0532	35			10-12
ZZ BASE ST-BL-BZSC SC (Black Zirconia screw)	BSIE0533	25			10-12
ZZ BASE ST-BL-ICSC NC 3.3 (Impression coping Screw)	BSFE0501	15			10-12
ZZ BASE ST-BL-ICSC RC 4.1 (Impression coping Screw)	BSFE0502	15			10-12
ZZ BASE ST-BL-ICSC SC (Impression coping Screw)	BSFE0503	15			10-12
ZZ BASE ST-BL-LBSC NC 3.3 (Lab Screw)	BSHE0501	35			10-12
ZZ BASE ST-BL-LBSC RC 4.1 (Lab Screw)	BSHE0502	35			10-12
ZZ BASE ST-BL-LBSC SC (Lab Screw)	BSHE0503	25			10-12

STRAUMANN® BONE LEVEL®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE ST-BL-LBZSC NC 3.3 (Lab Zirconia screw)	BSHE0531	35			10-12
ZZ BASE ST-BL-LBZSC RC 4.1 (Lab Zirconia screw)	BSHE0532	35			10-12
ZZ BASE ST-BL-LBZSC SC (Lab Zirconia screw)	BSHE0533	25			10-12
ZZ BASE ST-BL-MI NC 3.3 (Model implant)	MIAG0501				10-12
ZZ BASE ST-BL-MI RC 4.1 (Model implant)	MIAG0502				10-12
ZZ BASE ST-BL-MI SC (Model implant)	MIAG0503				10-12
ZZ BASE ST-BL-SC NC 3.3 (Screw)	BSAE0501	35			10-12
ZZ BASE ST-BL-SC RC 4.1 (Screw)	BSAE0502	35			10-12
ZZ BASE ST-BL-SC SC (Screw)	BSAE0503	25			10-12
ZZ BASE ST-BL-SC-G SC (Screw Gold)	BSCE0503	25			10-12
ZZ BASE ST-BL-WMSM-L10 NC 3.3 (White Metal Scanmarker L10 with Screw)	ZBED0531L10	15			10-12
ZZ BASE ST-BL-WMSM-L10 RC 4.1 (White Metal Scanmarker L10 with Screw)	ZBED0532L10	15			10-12
ZZ BASE ST-BL-WMSM-L10 SC (White Metal Scanmarker L10 with Screw)	ZBED0533L10	15			10-12
ZZ BASE ST-BL-WMSM-L12 NC 3.3 (White Metal Scanmarker L12 with Screw)	ZBED0531L12	15			10-12
ZZ BASE ST-BL-WMSM-L12 RC 4.1 (White Metal Scanmarker L12 with Screw)	ZBED0532L12	15			10-12
ZZ BASE ST-BL-WMSM-L12 SC (White Metal Scanmarker L12 with Screw)	ZBED0533L12	15			10-12
ZZ BASE ST-BL-WMSM-L15 NC 3.3 (White Metal Scanmarker L15 with Screw)	ZBED0531L15	15			10-12
ZZ BASE ST-BL-WMSM-L15 RC 4.1 (White Metal Scanmarker L15 with Screw)	ZBED0532L15	15			10-12
ZZ BASE ST-BL-WMSM-L15 SC (White Metal Scanmarker L15 with Screw)	ZBED0533L15	15			10-12
ZZ BASE ST-BL-ZSC NC 3.3 (Zirconia screw)	BSAE0531	35			10-12
ZZ BASE ST-BL-ZSC RC 4.1 (Zirconia screw)	BSAE0532	35			10-12
ZZ BASE ST-BL-ZSC SC (Zirconia screw)	BSAE0533	25			10-12
ZZ BASE ST-BL-ZSC-G NC 3.3 (Zirconia Screw Gold)	BSCE0531	35			10-12
ZZ BASE ST-BL-ZSC-G RC 4.1 (Zirconia Screw Gold)	BSCE0532	35			10-12
ZZ BASE ST-BL-ZSC-G SC (Zirconia Screw Gold)	BSCE0533	25			10-12
ZZ BASE ST-BL-ZZMUASC NC 3.3 (Zirkonzahn Multi Unit Abutment Screw)	BSEE0501	20			10-12
ZZ BASE ST-BL-ZZMUASC RC 4.1 (Zirkonzahn Multi Unit Abutment Screw)	BSEE0502	20			10-12
ZZ BASE ST-BL-ZZMUASC SC (Zirkonzahn Multi Unit Abutment Screw)	BSEE0503	15			10-12

STRAUMANN® MULTIBASE ABUTMENT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE ST-MA-CCTB 3.5 (Conical cemented titan base non hex and screw)	BSAD5331	15			10-12
SET ZZ BASE ST-MA-CCTB 4.5 (Conical cemented titan base non hex and screw)	BSAD5332	15			10-12
SET ZZ BASE ST-MA-CCTB-K85 3.5 (Conical cemented titan base non hex K85 and screw)	BSAD5331K85	15			10-12
SET ZZ BASE ST-MA-CCTB-K85 4.5 (Conical cemented titan base non hex K85 and screw)	BSAD5332K85	15			10-12
SET ZZ BASE ST-MA-ICOT 3.5 (Impression coping open tray non hex and Screw and Driver)	BSAQ5361	15			10-12
SET ZZ BASE ST-MA-ICOT 4.5 (Impression coping open tray non hex and Screw and Driver)	BSAQ5362	15			10-12
SET ZZ BASE ST-MA-SA 3.5 (Scananalog and screw)	ZBDE5331	15			10-12
SET ZZ BASE ST-MA-SA 4.5 (Scananalog and screw)	ZBDE5332	15			10-12
SET ZZ BASE ST-MA-SM 3.5 (Scanmarker and screw)	ZBAD5331	15			10-12
SET ZZ BASE ST-MA-SM 4.5 (Scanmarker and screw)	ZBAD5332	15			10-12
SET ZZ BASE ST-MA-WSM 3.5 (White Scanmarker and screw)	ZBCD5331	15			10-12
SET ZZ BASE ST-MA-WSM 4.5 (White Scanmarker and screw)	ZBCD5332	15			10-12
ZZ BASE ST-MA-ICSC 3.5 (Impression coping Screw)	BSFE5301	15			10-12
ZZ BASE ST-MA-ICSC 4.5 (Impression coping Screw)	BSFE5302	15			10-12
ZZ BASE ST-MA-SC 3.5 (Screw)	BSAE5301	15			10-12
ZZ BASE ST-MA-SC 4.5 (Screw)	BSAE5302	15			10-12
ZZ BASE ST-MA-WMSM-L12 3.5 (White Metal Scanmarker L12 with Screw)	ZBED5331L12	15			10-12
ZZ BASE ST-MA-WMSM-L12 4.5 (White Metal Scanmarker L12 with Screw)	ZBED5332L12	15			10-12
ZZ BASE ST-MA-ZSC 3.5 (Zirconia screw)	BSAE5331	15			10-12
ZZ BASE ST-MA-ZSC 4.5 (Zirconia screw)	BSAE5332	15			10-12

STRAUMANN® SCREW-RETAINED ABUTMENT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE ST-SA-CCTB 3.5 (Conical cemented titan base non hex and screw)	BSAD8031	15			10-12
SET ZZ BASE ST-SA-CCTB 3.5 (Conical cemented titan base non hex and screw)	BSAD8041	15			10-12
SET ZZ BASE ST-SA-CCTB 4.6 (Conical cemented titan base non hex and screw)	BSAD8032	15			10-12
SET ZZ BASE ST-SA-CCTB 4.6 (Conical cemented titan base non hex and screw)	BSAD8042	15			10-12
SET ZZ BASE ST-SA-CCTB-K85 3.5 (Conical cemented titan base non hex K85 and screw)	BSAD8031K85	15			10-12
SET ZZ BASE ST-SA-CCTB-K85 4.6 (Conical cemented titan base non hex K85 and screw)	BSAD8042K85	15			10-12
SET ZZ BASE ST-SA-ICOT 3.5 (Impression coping open tray non hex and Screw and Driver)	BSAQ8061	15			10-12
SET ZZ BASE ST-SA-ICOT 4.6 (Impression coping open tray non hex and Screw and Driver)	BSAQ8062	15			10-12
SET ZZ BASE ST-SA-SA 3.5 (Scananalog and screw)	ZBDE8031	15			10-12
SET ZZ BASE ST-SA-SA 4.6 (Scananalog and screw)	ZBDE8032	15			10-12
SET ZZ BASE ST-SA-SM 3.5 (Scanmarker and screw)	ZBAD8031	15			10-12
SET ZZ BASE ST-SA-SM 4.6 (Scanmarker and screw)	ZBAD8032	15			10-12
SET ZZ BASE ST-SA-WSM 3.5 (White Scanmarker and screw)	ZBCD8031	15			10-12
ZZ BASE ST-SA-ASCSC 3.5 (Angled screw-channel Screw)	BSKE8001	15			10-12
ZZ BASE ST-SA-ASCSC 4.6 (Angled screw-channel Screw)	BSKE8002	15			10-12
ZZ BASE ST-SA-ICSC 3.5 (Impression coping Screw)	BSFE8001	15			10-12
ZZ BASE ST-SA-ICSC 4.6 (Impression coping Screw)	BSFE8002	15			10-12
ZZ BASE ST-SA-SC 3.5 (Screw)	BSAE8001	15			10-12
ZZ BASE ST-SA-SC 4.6 (Screw)	BSAE8002	15			10-12
ZZ BASE ST-SA-WMSM-L12 3.5 (White Metal Scanmarker L12 with Screw)	ZBED8031L12	15			10-12
ZZ BASE ST-SA-WMSM-L12 4.6 (White Metal Scanmarker L12 with Screw)	ZBED8032L12	15			10-12
ZZ BASE ST-SA-ZSC 3.5 (Zirconia Screw)	BSAE8031	15			10-12
ZZ BASE ST-SA-ZSC 4.6 (Zirconia Screw)	BSAE8032	15			10-12

STRAUMANN® TISSUE LEVEL (STANDARD PLUS NARROW NECK CROSSFIT®)

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE ST-SO-ABCH NNC 3.5 (Raw-Abutment® customizable hex and screw)	BSAF0684	35			10-12
SET ZZ BASE ST-SO-ABCH RN 4.8 (Raw-Abutment® customizable hex and screw)	BSAF0682	35			10-12
SET ZZ BASE ST-SO-ABCH WN 6.5 (Raw-Abutment® customizable hex and screw)	BSAF0683	35			10-12
SET ZZ BASE ST-SO-ABCH-D14 NNC 3.5 (Raw-Abutment® customizable hex D14 and screw)	BSAF0684D14	35			10-12
SET ZZ BASE ST-SO-ABCH-D14 RN 4.8 (Raw-Abutment® customizable hex D14 and screw)	BSAF0682D14	35			10-12
SET ZZ BASE ST-SO-ABCH-D14 WN 6.5 (Raw-Abutment® customizable hex D14 and screw)	BSAF0683D14	35			10-12
SET ZZ BASE ST-SO-ASCTBH-K80 NNC 3.5 (Angled screw-channel titan base hex K80 and screw)	BSAT0684K80	20			10-12
SET ZZ BASE ST-SO-ASCTBH-K80 RN 4.8 (Angled screw-channel titan base hex K80 and screw)	BSAT0682K80	20			10-12
SET ZZ BASE ST-SO-ASCTBH-K80 WN 6.5 (Angled screw-channel titan base hex K80 and screw)	BSAT0683K80	20			10-12
SET ZZ BASE ST-SO-ASCTB-K80 NNC 3.5 (Angled screw-channel titan base non hex K80 and screw)	BSAT0634K80	20			10-12
SET ZZ BASE ST-SO-ASCTB-K80 RN 4.8 (Angled screw-channel titan base non hex K80 and screw)	BSAT0632K80	20			10-12
SET ZZ BASE ST-SO-ASCTB-K80 WN 6.5 (Angled screw-channel titan base non hex K80 and screw)	BSAT0633K80	20			10-12
SET ZZ BASE ST-SO-CCTB NN 3.5 (Conical cemented titan base non hex and screw)	BSAD0631	35			10-12
SET ZZ BASE ST-SO-CCTB NNC 3.5 (Conical cemented titan base non hex and screw)	BSAD0634	35			10-12
SET ZZ BASE ST-SO-CCTB RN 4.8 (Conical cemented titan base non hex and screw)	BSAD0622	35			10-12
SET ZZ BASE ST-SO-CCTB RN 4.8 (Conical cemented titan base non hex and screw)	BSAD0642	35			10-12
SET ZZ BASE ST-SO-CCTB RN 4.8 (Conical cemented titan base non hex and screw)	BSAD0632	35			10-12
SET ZZ BASE ST-SO-CCTB WN 6.5 (Conical cemented titan base non hex and screw)	BSAD0643	35			10-12
SET ZZ BASE ST-SO-CCTB WN 6.5 (Conical cemented titan base non hex and screw)	BSAD0633	35			10-12
SET ZZ BASE ST-SO-CCTB-G NN 3.5 (Conical cemented titan base non hex and screw Gold)	BSCD0631	35			10-12
SET ZZ BASE ST-SO-CCTB-G NNC 3.5 (Conical cemented titan base non hex and screw Gold)	BSCD0634	35			10-12
SET ZZ BASE ST-SO-CCTB-G RN 4.8 (Conical cemented titan base non hex and screw Gold)	BSCD0642	35			10-12
SET ZZ BASE ST-SO-CCTB-G WN 6.5 (Conical cemented titan base non hex and screw Gold)	BSCD0633	35			10-12
SET ZZ BASE ST-SO-CCTB-G WN 6.5 (Conical cemented titan base non hex and screw Gold)	BSCD0643	35			10-12
SET ZZ BASE ST-SO-CCTB-K85 NNC 3.5 (Conical cemented titan base non hex K85 and screw)	BSAD0634K85	35			10-12
SET ZZ BASE ST-SO-CCTB-K85 RN 4.8 (Conical cemented titan base non hex K85 and screw)	BSAD0632K85	35			10-12
SET ZZ BASE ST-SO-CCTB-K85 WN 6.5 (Conical cemented titan base non hex K85 and screw)	BSAD0633K85	35			10-12
SET ZZ BASE ST-SO-CCTB-K85 WN 6.5 (Conical cemented titan base non hex K85 and screw)	BSAD0643K85	35			10-12
SET ZZ BASE ST-SO-CCTB-L15 RN 4.8 (Conical cemented titan base non hex L15 and screw)	BSAD0632L15	35			10-12

STRAUMANN® TISSUE LEVEL (STANDARD PLUS NARROW NECK CROSSFIT®)

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
SET ZZ BASE ST-SO-CCTB-L15 WN 6.5 (Conical cemented titan base non hex L15 and screw)	BSAD0633L15	35			10-12
SET ZZ BASE ST-SO-CCTB-L15-G RN 4.8 (Conical cemented titan base non hex L15 and screw)	BSCD0632L15	35			10-12
SET ZZ BASE ST-SO-CCTB-L15-G WN 6.5 (Conical cemented titan base non hex L15 and screw)	BSCD0633L15	35			10-12
SET ZZ BASE ST-SO-CCTB-L20 NNC 3.5 (Conical cemented titan base non hex L20 and screw)	BSAD0634L20	35			10-12
SET ZZ BASE ST-SO-CCTB-L20-G NNC 3.5 (Conical cemented titan base non hex L20 and screw)	BSCD0634L20	35			10-12
SET ZZ BASE ST-SO-CCTB-L30 NNC 3.5 (Conical cemented titan base non hex L30 and screw)	BSAD0634L30	35			10-12
SET ZZ BASE ST-SO-CCTB-L30-G NNC 3.5 (Conical cemented titan base non hex L30 and screw)	BSCD0634L30	35			10-12
SET ZZ BASE ST-SO-ICOT NNC 3.5 (Impression coping open tray non hex and Screw and Driver)	BSAQ0664	15			10-12
SET ZZ BASE ST-SO-ICOT RN 4.8 (Impression coping open tray non hex and Screw and Driver)	BSAQ0662	15			10-12
SET ZZ BASE ST-SO-ICOT WN 6.5 (Impression coping open tray non hex and Screw and Driver)	BSAQ0663	15			10-12
SET ZZ BASE ST-SO-ICOTH NNC 3.5 (Impression coping open tray hex and Screw and Driver)	BSAQ0644	15			10-12
SET ZZ BASE ST-SO-ICOTH RN 4.8 (Impression coping open tray hex and Screw and Driver)	BSAQ0642	15			10-12
SET ZZ BASE ST-SO-ICOTH WN 6.5 (Impression coping open tray hex and Screw and Driver)	BSAQ0643	15			10-12
SET ZZ BASE ST-SO-NTB NNC 3.5 (Narrow titan base non hex and screw)	BSAP0634	35			10-12
SET ZZ BASE ST-SO-NTBH NNC 3.5 (Narrow titan base hex and screw)	BSAP0684	35			10-12
SET ZZ BASE ST-SO-PCTBH NN 3.5 (Parallel cemented titan base hex and screw)	BSAB0681	35			10-12
SET ZZ BASE ST-SO-PCTBH NNC 3.5 (Parallel cemented titan base hex and screw)	BSAB0684	35			10-12
SET ZZ BASE ST-SO-PCTBH RN 4.8 (Parallel cemented titan base hex and screw)	BSAB0682	35			10-12
SET ZZ BASE ST-SO-PCTBH RN 4.8 (Parallel cemented titan base hex and screw)	BSAB0622	35			10-12
SET ZZ BASE ST-SO-PCTBH RN 4.8 (Parallel cemented titan base hex and screw)	BSAB0672	35			10-12
SET ZZ BASE ST-SO-PCTBH RN 4.8 (Parallel cemented titan base hex and screw)	BSAB0692	35			10-12
SET ZZ BASE ST-SO-PCTBH WN 6.5 (Parallel cemented titan base hex and screw)	BSAB0683	35			10-12
SET ZZ BASE ST-SO-PCTBH WN 6.5 (Parallel cemented titan base hex and screw)	BSAB0693	35			10-12
SET ZZ BASE ST-SO-PCTBH-G NNC 3.5 (Parallel cemented titan base hex and screw Gold)	BSCB0684	35			10-12
SET ZZ BASE ST-SO-PCTBH-G RN 4.8 (Parallel cemented titan base hex and screw Gold)	BSCB0692	35			10-12
SET ZZ BASE ST-SO-PCTBH-G WN 6.5 (Parallel cemented titan base hex and screw Gold)	BSCB0683	35			10-12
SET ZZ BASE ST-SO-PCTBH-G WN 6.5 (Parallel cemented titan base hex and screw Gold)	BSCB0693	35			10-12
SET ZZ BASE ST-SO-PCTBH-K85 NNC 3.5 (Parallel cemented titan base hex K85 and screw)	BSAB0684K85	35			10-12
SET ZZ BASE ST-SO-PCTBH-K85 RN 4.8 (Parallel cemented titan base hex K85 and screw)	BSAB0682K85	35			10-12

STRAUMANN® TISSUE LEVEL (STANDARD PLUS NARROW NECK CROSSFIT®)

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE ST-SO-PCTBH-K85 WN 6.5 (Parallel cemented titan base hex K85 and screw)	BSAB0683K85	35			10-12
SET ZZ BASE ST-SO-PCTBH-K85 WN 6.5 (Parallel cemented titan base hex K85 and screw)	BSAB0693K85	35			10-12
SET ZZ BASE ST-SO-PCTBH-L10 NNC 3.5 (Parallel cemented titan base hex L10 and screw)	BSAB0694L10	35			10-12
SET ZZ BASE ST-SO-PCTBH-L10-G NNC 3.5 (Parallel cemented titan base hex L10 and screw G)	BSCB0694L10	35			10-12
SET ZZ BASE ST-SO-PCTBH-L15 RN 4.8 (Parallel cemented titan base hex L15 and screw)	BSAB0692L15	35			10-12
SET ZZ BASE ST-SO-PCTBH-L15 WN 6.5 (Parallel cemented titan base hex L15 and screw)	BSAB0693L15	35			10-12
SET ZZ BASE ST-SO-PCTBH-L15-G RN 4.8 (Parallel cemented titan base hex L15 and screw Gol)	BSCB0692L15	35			10-12
SET ZZ BASE ST-SO-PCTBH-L20 NNC 3.5 (Parallel cemented titan base hex L20 and screw)	BSAB0694L20	35			10-12
SET ZZ BASE ST-SO-PCTBH-L20-G NNC 3.5 (Parallel cemented titan base hex L20 and screw G)	BSCB0694L20	35			10-12
SET ZZ BASE ST-SO-PCTBH-L30 NNC 3.5 (Parallel cemented titan base hex L30 and screw)	BSAB0694L30	35			10-12
SET ZZ BASE ST-SO-PCTBH-L30-G NNC 3.5 (Parallel cemented titan base hex L30 and screw G)	BSCB0694L30	35			10-12
SET ZZ BASE ST-SO-SA NNC 3.5 (Scananalog and screw)	ZBDE0634	15			10-12
SET ZZ BASE ST-SO-SA RN 4.8 (Scananalog and screw)	ZBDE0632	15			10-12
SET ZZ BASE ST-SO-SA WN 6.5 (Scananalog and screw)	ZBDE0633	15			10-12
SET ZZ BASE ST-SO-SM NNC 3.5 (Scanmarker and screw)	ZBAD0634	35			10-12
SET ZZ BASE ST-SO-SM RN 4.8 (Scanmarker and screw)	ZBAD0632	35			10-12
SET ZZ BASE ST-SO-SM RN 4.8 (Scanmarker with Synocta and screw)	ZBAD0642	35			10-12
SET ZZ BASE ST-SO-SM WN 6.5 (Scanmarker and screw)	ZBAD0633	35			10-12
SET ZZ BASE ST-SO-SM WN 6.5 (Scanmarker with Synocta and screw)	ZBAD0643	35			10-12
SET ZZ BASE ST-SO-WSM NNC 3.5 (White Scanmarker and screw)	ZBCD0634	35			10-12
SET ZZ BASE ST-SO-ZZMUA17 NNC 3.5 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and i)	BSAI0634	20			10-12
SET ZZ BASE ST-SO-ZZMUA17 NNC 3.5 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and i)	BSAI0644	20			10-12
SET ZZ BASE ST-SO-ZZMUA17 RN 4.8 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and i)	BSAI0632	30			10-12
SET ZZ BASE ST-SO-ZZMUA17 RN 4.8 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and i)	BSAI0642	30			10-12
SET ZZ BASE ST-SO-ZZMUA17-L10 NNC 3.5 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and i)	BSAI0634L10	20			10-12
SET ZZ BASE ST-SO-ZZMUA17-L10 NNC 3.5 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and i)	BSAI0644L10	20			10-12
SET ZZ BASE ST-SO-ZZMUA-L05 WN 6.5 (Zirkonzahn Multi Unit Abutment non hex L05 and Scre)	BSAH0633L05	25	35		10-12
SET ZZ BASE ST-SO-ZZMUA-L15 NNC 3.5 (Zirkonzahn Multi Unit Abutment non hex L15 and Scr)	BSAH0634L15	25	35		10-12
SET ZZ BASE ST-SO-ZZMUA-L15 RN 4.8 (Zirkonzahn Multi Unit Abutment non hex L15 and Scre)	BSAH0632L15	25	35		10-12

STRAUMANN® TISSUE LEVEL (STANDARD PLUS NARROW NECK CROSSFIT®)

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE ST-SO-ZZMUA-L15 WN 6.5 (Zirkonzahn Multi Unit Abutment non hex L15 and Scre	BSAH0633L15	25	35		10-12
SET ZZ BASE ST-SO-ZZMUA-L25 NNC 3.5 (Zirkonzahn Multi Unit Abutment non hex L25 and Scr	BSAH0634L25	25	35		10-12
SET ZZ BASE ST-SO-ZZMUA-L25 RN 4.8 (Zirkonzahn Multi Unit Abutment non hex L25 and Scre	BSAH0632L25	25	35		10-12
SET ZZ BASE ST-SO-ZZMUA-L35 NNC 3.5 (Zirkonzahn Multi Unit Abutment non hex L35 and Scr	BSAH0634L35	25	35		10-12
SET ZZ BASE ST-SO-ZZMUA-L35 RN 4.8 (Zirkonzahn Multi Unit Abutment non hex L35 and Scre	BSAH0632L35	25	35		10-12
SET ZZ BASE ST-SO-ZZMUA-L45 NNC 3.5 (Zirkonzahn Multi Unit Abutment non hex L45 and Scr	BSAH0634L45	25	35		10-12
SET ZZ BASE ST-SO-ZZMUA-L45 RN 4.8 (Zirkonzahn Multi Unit Abutment non hex L45 and Scre	BSAH0632L45	25	35		10-12
ZZ BASE ST-SO-ASCSC NNC 3.5 (Angled screw-channel Screw)	BSKE0604	20			10-12
ZZ BASE ST-SO-ASCSC RN 4.8 (Angled screw-channel Screw)	BSKE0602	20			10-12
ZZ BASE ST-SO-ASCSC WN 6.5 (Angled screw-channel Screw)	BSKE0603	20			10-12
ZZ BASE ST-SO-BLZSC RN 4.8 (Black Long Zirconia Screw)	BSIE0682	35			10-12
ZZ BASE ST-SO-BLZSC WN 6.5 (Black Long Zirconia Screw)	BSIE0683	35			10-12
ZZ BASE ST-SO-BSC NNC 3.5 (Black Screw)	BSIE0604	35			10-12
ZZ BASE ST-SO-BSC RN 4.8 (Black Screw)	BSIE0602	35			10-12
ZZ BASE ST-SO-BSC WN 6.5 (Black Screw)	BSIE0603	35			10-12
ZZ BASE ST-SO-BZSC NNC 3.5 (Black Zirconia screw)	BSIE0634	35			10-12
ZZ BASE ST-SO-BZSC RN 4.8 (Black Zirconia screw)	BSIE0632	35			10-12
ZZ BASE ST-SO-BZSC WN 6.5 (Black Zirconia screw)	BSIE0633	35			10-12
ZZ BASE ST-SO-ICSC NNC 3.5 (Impression coping Screw)	BSFE0604	15			10-12
ZZ BASE ST-SO-ICSC RN 4.8 (Impression coping Screw)	BSFE0602	15			10-12
ZZ BASE ST-SO-ICSC WN 6.5 (Impression coping Screw)	BSFE0603	15			10-12
ZZ BASE ST-SO-LBLZSC RN 4.8 (Lab Long Zirconia Screw)	BSHE0682	35			10-12
ZZ BASE ST-SO-LBLZSC WN 6.5 (Lab Long Zirconia Screw)	BSHE0683	35			10-12
ZZ BASE ST-SO-LBSC NNC 3.5 (Lab Screw)	BSHE0604	35			10-12
ZZ BASE ST-SO-LBSC RN 4.8 (Lab Screw)	BSHE0602	35			10-12
ZZ BASE ST-SO-LBSC WN 6.5 (Lab Screw)	BSHE0603	35			10-12
ZZ BASE ST-SO-LBZSC NNC 3.5 (Lab Zirconia screw)	BSHE0634	35			10-12
ZZ BASE ST-SO-LBZSC RN 4.8 (Lab Zirconia screw)	BSHE0632	35			10-12
ZZ BASE ST-SO-LBZSC WN 6.5 (Lab Zirconia screw)	BSHE0633	35			10-12

STRAUMANN® TISSUE LEVEL (STANDARD PLUS NARROW NECK CROSSFIT®)

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE ST-SO-LZSC RN 4.8 (Long zirconia screw)	BSAE0682	35			10-12
ZZ BASE ST-SO-LZSC WN 6.5 (Long zirconia screw)	BSAE0683	35			10-12
ZZ BASE ST-SO-LZSC-G RN 4.8 (Long Zirconia Screw Gold)	BSCE0682	35			10-12
ZZ BASE ST-SO-LZSC-G WN 6.5 (Long Zirconia Screw Gold)	BSCE0683	35			10-12
ZZ BASE ST-SO-MI NNC 3.5 (Model implant)	MIAG0604				10-12
ZZ BASE ST-SO-MI RN 4.8 (Model implant)	MIAG0602				10-12
ZZ BASE ST-SO-MI WN 6.5 (Model implant)	MIAG0603				10-12
ZZ BASE ST-SO-SC NN 3.5 (Screw)	BSAE0601	35			10-12
ZZ BASE ST-SO-SC NNC 3.5 (Screw)	BSAE0604	35			10-12
ZZ BASE ST-SO-SC RN 4.8 (Screw)	BSAE0602	35			10-12
ZZ BASE ST-SO-SC WN 6.5 (Screw)	BSAE0603	35			10-12
ZZ BASE ST-SO-SC-G NN 3.5 (Screw Gold)	BSCE0601	35			10-12
ZZ BASE ST-SO-SC-G NNC 3.5 (Screw Gold)	BSCE0604	35			10-12
ZZ BASE ST-SO-WMSM-L10 NN 3.5 (White Metal Scanmarker L10 with Screw)	ZBED0631L10	15			10-12
ZZ BASE ST-SO-WMSM-L10 NNC 3.5 (White Metal Scanmarker L10 with Screw)	ZBED0634L10	15			10-12
ZZ BASE ST-SO-WMSM-L10 RN 4.8 (White Metal Scanmarker L10 with Screw)	ZBED0632L10	15			10-12
ZZ BASE ST-SO-WMSM-L10 WN 6.5 (White Metal Scanmarker L10 with Screw)	ZBED0633L10	15			10-12
ZZ BASE ST-SO-WMSM-L12 NN 3.5 (White Metal Scanmarker L12 with Screw)	ZBED0631L12	15			10-12
ZZ BASE ST-SO-WMSM-L12 NNC 3.5 (White Metal Scanmarker L12 with Screw)	ZBED0634L12	15			10-12
ZZ BASE ST-SO-WMSM-L12 RN 4.8 (White Metal Scanmarker L12 with Screw)	ZBED0632L12	15			10-12
ZZ BASE ST-SO-WMSM-L12 WN 6.5 (White Metal Scanmarker L12 with Screw)	ZBED0633L12	15			10-12
ZZ BASE ST-SO-WMSM-L15 NN 3.5 (White Metal Scanmarker L15 with Screw)	ZBED0631L15	15			10-12
ZZ BASE ST-SO-WMSM-L15 NNC 3.5 (White Metal Scanmarker L15 with Screw)	ZBED0634L15	15			10-12
ZZ BASE ST-SO-WMSM-L15 RN 4.8 (White Metal Scanmarker L15 with Screw)	ZBED0632L15	15			10-12
ZZ BASE ST-SO-WMSM-L15 WN 6.5 (White Metal Scanmarker L15 with Screw)	ZBED0633L15	15			10-12
ZZ BASE ST-SO-ZSC NN 3.5 (Zirconia screw)	BSAE0631	35			10-12
ZZ BASE ST-SO-ZSC NNC 3.5 (Zirconia screw)	BSAE0634	35			10-12
ZZ BASE ST-SO-ZSC RN 4.8 (Zirconia screw)	BSAE0632	35			10-12
ZZ BASE ST-SO-ZSC WN 6.5 (Zirconia screw)	BSAE0633	35			10-12

STRAUMANN® TISSUE LEVEL (STANDARD PLUS NARROW NECK CROSSFIT®)

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE ST-SO-ZSC-G NN 3.5 (Zirconia Screw Gold)	BSCE0631	35			10-12
ZZ BASE ST-SO-ZSC-G RN 4.8 (Zirconia Screw Gold)	BSCE0632	35			10-12
ZZ BASE ST-SO-ZSC-G WN 6.5 (Zirconia Screw Gold)	BSCE0633	35			10-12
ZZ BASE ST-SO-ZZMUASC NNC 3.5 (Zirkonzahn Multi Unit Abutment Screw)	BSEE0604	20			10-12
ZZ BASE ST-SO-ZZMUASC RN 4.8 (Zirkonzahn Multi Unit Abutment Screw)	BSEE0602	30			10-12

STRAUMANN® TLX

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
ZZ BASE ST-TX-MI NT (Model implant)	MIAG16401				10-12
ZZ BASE ST-TX-MI RT (Model implant)	MIAG16402				10-12

SWEDEN & MARTINA OUTLINK²

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE S&M-O-SM 3.30 (Scanmarker and screw)	ZBAD1331	30			10-12
SET ZZ BASE S&M-O-SM 4.10 (Scanmarker and screw)	ZBAD1332	35			10-12
SET ZZ BASE S&M-O-ASCTBH-K80 3.3 (Angled screw-channel titan base hex K80 and screw)	BSAT1381K80	20			10-12
SET ZZ BASE S&M-O-ASCTBH-K80 4.1 (Angled screw-channel titan base hex K80 and screw)	BSAT1382K80	20			10-12
SET ZZ BASE S&M-O-ASCTB-K80 3.3 (Angled screw-channel titan base non hex K80 and screw)	BSAT1331K80	20			10-12
SET ZZ BASE S&M-O-ASCTB-K80 4.1 (Angled screw-channel titan base non hex K80 and screw)	BSAT1332K80	20			10-12
SET ZZ BASE S&M-O-CCTB 3.30 (Conical cemented titan base non hex and screw)	BSAD1331	30			10-12
SET ZZ BASE S&M-O-CCTB 4.10 (Conical cemented titan base non hex and screw)	BSAD1332	35			10-12
SET ZZ BASE S&M-O-CCTB 5.0 (Conical cemented titan base non hex and screw)	BSAD1333	35			10-12
SET ZZ BASE S&M-O-CCTB-K85 3.3 (Conical cemented titan base non hex K85 and screw)	BSAD1331K85	30			10-12
SET ZZ BASE S&M-O-CCTB-K85 4.1 (Conical cemented titan base non hex K85 and screw)	BSAD1332K85	35			10-12
SET ZZ BASE S&M-O-CCTB-K85 5.0 (Conical cemented titan base non hex K85 and screw)	BSAD1333K85	35			10-12
SET ZZ BASE S&M-O-ICOTH 3.3 (Impression coping open tray hex and Screw and Driver)	BSAQ1341	15			10-12
SET ZZ BASE S&M-O-ICOTH 4.1 (Impression coping open tray hex and Screw and Driver)	BSAQ1342	15			10-12
SET ZZ BASE S&M-O-ICOTH 5.0 (Impression coping open tray hex and Screw and Driver)	BSAQ1343	15			10-12
SET ZZ BASE S&M-O-NTB 3.3 (Narrow titan base non hex and screw)	BSAP1331	30			10-12
SET ZZ BASE S&M-O-NTB 4.1 (Narrow titan base non hex and screw)	BSAP1332	35			10-12
SET ZZ BASE S&M-O-NTBH 3.3 (Narrow titan base hex and screw)	BSAP1381	30			10-12
SET ZZ BASE S&M-O-NTBH 4.1 (Narrow titan base hex and screw)	BSAP1382	35			10-12
SET ZZ BASE S&M-O-PCTBH 3.3 (Parallel cemented titan base hex and screw)	BSAB1381	30			10-12
SET ZZ BASE S&M-O-PCTBH 4.1 (Parallel cemented titan base hex and screw)	BSAB1382	35			10-12
SET ZZ BASE S&M-O-PCTBH 4.1 SP (Parallel cemented titan base hex and screw)	BSAB1384	35			10-12
SET ZZ BASE S&M-O-PCTBH 5.0 (Parallel cemented titan base hex and screw)	BSAB1383	35			10-12
SET ZZ BASE S&M-O-PCTBH-K85 3.3 (Parallel cemented titan base hex K85 and screw)	BSAB1381K85	30			10-12
SET ZZ BASE S&M-O-PCTBH-K85 4.1 (Parallel cemented titan base hex K85 and screw)	BSAB1382K85	35			10-12
SET ZZ BASE S&M-O-PCTBH-K85 5.0 (Parallel cemented titan base hex K85 and screw)	BSAB1383K85	35			10-12
SET ZZ BASE S&M-O-SA 3.3 (Scananalog and screw)	ZBDE1331	15			10-12
SET ZZ BASE S&M-O-SA 4.1 (Scananalog and screw)	ZBDE1332	15			10-12
SET ZZ BASE S&M-O-SA 5.0 (Scananalog and screw)	ZBDE1333	15			10-12

SWEDEN & MARTINA OUTLINK²

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE S&M-O-SM 5.0 (Scanmarker and screw)	ZBAD1333	35			10-12
SET ZZ BASE S&M-O-WSM 3.3 (White Scanmarker and screw)	ZBCD1331	30			10-12
SET ZZ BASE S&M-O-WSM 4.1 (White Scanmarker and screw)	ZBCD1332	35			10-12
SET ZZ BASE S&M-O-WSM 5.0 (White Scanmarker and screw)	ZBCD1333	35			10-12
ZZ BASE S&M-O-SC 3.30 (Screw)	BSAE1301	30			10-12
ZZ BASE S&M-O-SC 4.10 (Screw)	BSAE1302	35			10-12
ZZ BASE S&M-O-SC 4.10 SP (Screw)	BSAE1304	35			10-12
ZZ BASE S&M-O-ZSC 3.30 (Zirconia screw)	BSAE1331	30			10-12
ZZ BASE S&M-O-ZSC 4.10 (Zirconia screw)	BSAE1332	35			10-12
ZZ BASE S&M-O-ZSC 4.10 SP (Zirconia screw)	BSAE1334	35			10-12
ZZ BASE S&M-O-ASCSC 3.3 (Angled screw-channel Screw)	BSKE1301	20			10-12
ZZ BASE S&M-O-ASCSC 4.1 (Angled screw-channel Screw)	BSKE1302	20			10-12
ZZ BASE S&M-O-ICSC 3.3 (Impression coping Screw)	BSFE1301	15			10-12
ZZ BASE S&M-O-ICSC 4.1 (Impression coping Screw)	BSFE1302	15			10-12
ZZ BASE S&M-O-ICSC 5.0 (Impression coping Screw)	BSFE1303	15			10-12
ZZ BASE S&M-O-SC 5.0 (Screw)	BSAE1303	35			10-12
ZZ BASE S&M-O-SC-G 3.3 (Screw Gold)	BSCE1301	30			10-12
ZZ BASE S&M-O-SC-G 4.1 (Screw Gold)	BSCE1302	35			10-12
ZZ BASE S&M-O-SC-G 5.0 (Screw Gold)	BSCE1303	35			10-12
ZZ BASE S&M-O-WMSM-L12 3.3 (White Metal Scanmarker L12 with Screw)	ZBED1331L12	15			10-12
ZZ BASE S&M-O-WMSM-L12 4.1 (White Metal Scanmarker L12 with Screw)	ZBED1332L12	15			10-12
ZZ BASE S&M-O-WMSM-L12 5.0 (White Metal Scanmarker L12 with Screw)	ZBED1333L12	15			10-12
ZZ BASE S&M-O-ZSC 5.0 (Zirconia screw)	BSAE1333	35			10-12
ZZ BASE S&M-O-ZSC-G 3.3 (Zirconia Screw Gold)	BSCE1331	30			10-12

SWEDEN & MARTINA P.A.D® MULTI UNIT ABUTMENT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE S&M-MUA-SM (Scanmarker and screw)	ZBAD1431	15			10-12
SET ZZ BASE S&M-MUA-CCTB (Conical cemented titan base non hex and screw)	BSAD1431	15			10-12
SET ZZ BASE S&M-MUA-ICOT (Impression coping open tray non hex and Screw and Driver)	BSAQ1461	15			10-12
SET ZZ BASE S&M-MUA-SA (Scananalog and screw)	ZBDE1431	15			10-12
ZZ BASE S&M-MUA-ICSC (Impression coping Screw)	BSFE1401	15			10-12
ZZ BASE S&M-MUA-SC (Screw)	BSAE1401	15			10-12
ZZ BASE S&M-MUA-SC-G (Screw Gold)	BSCE1401	15			10-12
ZZ BASE S&M-MUA-WMSM-L12 (White Metal Scanmarker L12 with Screw)	ZBED1431L12	15			10-12
ZZ BASE S&M-MUA-ZSC (Zirconia screw)	BSAE1431	15			10-12

SWEDEN & MARTINA PRAMA

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
SET ZZ BASE S&M-PR-ABCH 3.8/4.25/5.0 (Raw-Abutment® customizable hex and screw)	BSAF12291	20-25			10-12
SET ZZ BASE S&M-PR-ASCTBH-K80 3.8/4.25/5.0 (Angled screw-channel titan base hex K80 and	BSAT12281K80	20			10-12
SET ZZ BASE S&M-PR-ASCTB-K80 3.8/4.25/5.0 (Angled screw-channel titan base non hex K80 a	BSAT12231K80	20			10-12
SET ZZ BASE S&M-PR-CCTB 3.8/4.25/5.0 (Conical cemented titan base non hex and screw)	BSAD12231	20-25			10-12
SET ZZ BASE S&M-PR-CCTB-K85 3.8/4.25/5.0 (Conical cemented titan base non hex K85 and sc	BSAD12231K85	20-25			10-12
SET ZZ BASE S&M-PR-ICOTH 3.8/4.25/5.0 (Impression coping open tray hex and Screw and Driv	BSAQ12241	15			10-12
SET ZZ BASE S&M-PR-PCTBH 3.8/4.25/5.0 (Parallel cemented titan base hex and screw)	BSAB12281	20-25			10-12
SET ZZ BASE S&M-PR-PCTBH-K85 3.8/4.25/5.0 (Parallel cemented titan base hex K85 and screw	BSAB12281K85	20-25			10-12
SET ZZ BASE S&M-PR-SA 3.8/4.25/5.0 (Scanalog and screw)	ZBDE12231	15			10-12
SET ZZ BASE S&M-PR-SM 3.8/4.25/5.0 (Scanmarker and screw)	ZBAD12231	20-25			10-12
ZZ BASE S&M-PR-ASCSC 3.8/4.25/5.0 (Angled screw-channel Screw)	BSKE12201	20			10-12
ZZ BASE S&M-PR-ICSC 3.8/4.25/5.0 (Impression coping Screw)	BSFE12201	15			10-12
ZZ BASE S&M-PR-MI 3.8/4.25/5.0 (Model implant)	MIAG12201				10-12
ZZ BASE S&M-PR-SC 3.8/4.25/5.0 (Screw)	BSAE12201	20-25			10-12
ZZ BASE S&M-PR-WMSM-L12 3.8/4.25/5.0 (White Metal Scanmarker L12 with Screw)	ZBED12231L12	15			10-12
ZZ BASE S&M-PR-ZSC 3.8/4.25/5.0 (Zirconia screw)	BSAE12231	20-25			10-12

SWEDEN & MARTINA PREMIUM KOHNO®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE S&M-K-SM 4.25 (Scanmarker and screw)	ZBAD1533	35			10-12
SET ZZ BASE S&M-K-ABCH 3.30 (Raw-Abutment® customizable hex and screw)	BSAF1591	30			10-12
SET ZZ BASE S&M-K-ABCH 3.80 (Raw-Abutment® customizable hex and screw)	BSAF1592	30			10-12
SET ZZ BASE S&M-K-ABCH 4.25 (Raw-Abutment® customizable hex and screw)	BSAF1593	35			10-12
SET ZZ BASE S&M-K-ABCH 5.0 (Raw-Abutment® customizable hex and screw)	BSAF1594	35			10-12
SET ZZ BASE S&M-K-ABCH 6.0 (Raw-Abutment® customizable hex and screw)	BSAF1595	35			10-12
SET ZZ BASE S&M-K-ASCTBH-K80 3.3 (Angled screw-channel titan base hex K80 and screw)	BSAT1581K80	20			10-12
SET ZZ BASE S&M-K-ASCTBH-K80 3.8 (Angled screw-channel titan base hex K80 and screw)	BSAT1582K80	20			10-12
SET ZZ BASE S&M-K-ASCTBH-K80 4.25 (Angled screw-channel titan base hex K80 and screw)	BSAT1583K80	20			10-12
SET ZZ BASE S&M-K-ASCTB-K80 3.3 (Angled screw-channel titan base non hex K80 and screw)	BSAT1531K80	20			10-12
SET ZZ BASE S&M-K-ASCTB-K80 3.8 (Angled screw-channel titan base non hex K80 and screw)	BSAT1532K80	20			10-12
SET ZZ BASE S&M-K-ASCTB-K80 4.25 (Angled screw-channel titan base non hex K80 and screw)	BSAT1533K80	20			10-12
SET ZZ BASE S&M-K-CCTB 3.30 (Conical cemented titan base non hex and screw)	BSAD1531	30			10-12
SET ZZ BASE S&M-K-CCTB 3.80 (Conical cemented titan base non hex and screw)	BSAD1532	30			10-12
SET ZZ BASE S&M-K-CCTB 4.25 (Conical cemented titan base non hex and screw)	BSAD1533	35			10-12
SET ZZ BASE S&M-K-CCTB 5.0 (Conical cemented titan base non hex and screw)	BSAD1534	35			10-12
SET ZZ BASE S&M-K-CCTB 6.0 (Conical cemented titan base non hex and screw)	BSAD1535	35			10-12
SET ZZ BASE S&M-K-CCTB-K85 3.3 (Conical cemented titan base non hex K85 and screw)	BSAD1531K85	30			10-12
SET ZZ BASE S&M-K-CCTB-K85 3.8 (Conical cemented titan base non hex K85 and screw)	BSAD1532K85	30			10-12
SET ZZ BASE S&M-K-CCTB-K85 4.25 (Conical cemented titan base non hex K85 and screw)	BSAD1533K85	35			10-12
SET ZZ BASE S&M-K-CCTB-K85 5.0 (Conical cemented titan base non hex K85 and screw)	BSAD1534K85	35			10-12
SET ZZ BASE S&M-K-CCTB-L10 6.0 (Conical cemented titan base non hex L10 and screw)	BSAD1535L10	35			10-12
SET ZZ BASE S&M-K-CCTB-L20 3.3 (Conical cemented titan base non hex L20 and screw)	BSAD1531L20	30			10-12
SET ZZ BASE S&M-K-CCTB-L20 3.8 (Conical cemented titan base non hex L20 and screw)	BSAD1532L20	30			10-12
SET ZZ BASE S&M-K-CCTB-L20 4.25 (Conical cemented titan base non hex L20 and screw)	BSAD1533L20	35			10-12
SET ZZ BASE S&M-K-CCTB-L20 5.0 (Conical cemented titan base non hex L20 and screw)	BSAD1534L20	35			10-12
SET ZZ BASE S&M-K-CCTB-L20 6.0 (Conical cemented titan base non hex L20 and screw)	BSAD1535L20	35			10-12
SET ZZ BASE S&M-K-CCTB-L30 3.3 (Conical cemented titan base non hex L30 and screw)	BSAD1531L30	30			10-12
SET ZZ BASE S&M-K-CCTB-L30 3.8 (Conical cemented titan base non hex L30 and screw)	BSAD1532L30	30			10-12

SWEDEN & MARTINA PREMIUM KOHNO®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE S&M-K-CCTB-L30 4.25 (Conical cemented titan base non hex L30 and screw)	BSAD1533L30	35			10-12
SET ZZ BASE S&M-K-CCTB-L30 5.0 (Conical cemented titan base non hex L30 and screw)	BSAD1534L30	35			10-12
SET ZZ BASE S&M-K-CCTB-L30 6.0 (Conical cemented titan base non hex L30 and screw)	BSAD1535L30	35			10-12
SET ZZ BASE S&M-K-CCTB-L40 3.3 (Conical cemented titan base non hex L40 and screw)	BSAD1531L40	30			10-12
SET ZZ BASE S&M-K-CCTB-L40 3.8 (Conical cemented titan base non hex L40 and screw)	BSAD1532L40	30			10-12
SET ZZ BASE S&M-K-CCTB-L40 4.25 (Conical cemented titan base non hex L40 and screw)	BSAD1533L40	35			10-12
SET ZZ BASE S&M-K-CCTB-L40 5.0 (Conical cemented titan base non hex L40 and screw)	BSAD1534L40	35			10-12
SET ZZ BASE S&M-K-CCTB-L40 6.0 (Conical cemented titan base non hex L40 and screw)	BSAD1535L40	35			10-12
SET ZZ BASE S&M-K-ICOT 3.3 (Impression coping open tray non hex and Screw and Driver)	BSAQ1561	15			10-12
SET ZZ BASE S&M-K-ICOT 3.8 (Impression coping open tray non hex and Screw and Driver)	BSAQ1562	15			10-12
SET ZZ BASE S&M-K-ICOT 4.25 (Impression coping open tray non hex and Screw and Driver)	BSAQ1563	15			10-12
SET ZZ BASE S&M-K-ICOT 5.0 (Impression coping open tray non hex and Screw and Driver)	BSAQ1564	15			10-12
SET ZZ BASE S&M-K-ICOT 6.0 (Impression coping open tray non hex and Screw and Driver)	BSAQ1565	15			10-12
SET ZZ BASE S&M-K-ICOTH 3.3 (Impression coping open tray hex and Screw and Driver)	BSAQ1541	15			10-12
SET ZZ BASE S&M-K-ICOTH 3.8 (Impression coping open tray hex and Screw and Driver)	BSAQ1542	15			10-12
SET ZZ BASE S&M-K-ICOTH 4.25 (Impression coping open tray hex and Screw and Driver)	BSAQ1543	15			10-12
SET ZZ BASE S&M-K-ICOTH 5.0 (Impression coping open tray hex and Screw and Driver)	BSAQ1544	15			10-12
SET ZZ BASE S&M-K-ICOTH 6.0 (Impression coping open tray hex and Screw and Driver)	BSAQ1545	15			10-12
SET ZZ BASE S&M-K-NTB 3.3 (Narrow titan base non hex and screw)	BSAP1531	30			10-12
SET ZZ BASE S&M-K-NTBH 3.3 (Narrow titan base hex and screw)	BSAP1581	30			10-12
SET ZZ BASE S&M-K-PCTBH 3.30 (Parallel cemented titan base hex and screw)	BSAB1581	30			10-12
SET ZZ BASE S&M-K-PCTBH 3.80 (Parallel cemented titan base hex and screw)	BSAB1582	30			10-12
SET ZZ BASE S&M-K-PCTBH 4.25 (Parallel cemented titan base hex and screw)	BSAB1583	35			10-12
SET ZZ BASE S&M-K-PCTBH 5.0 (Parallel cemented titan base hex and screw)	BSAB1584	35			10-12
SET ZZ BASE S&M-K-PCTBH 6.0 (Parallel cemented titan base hex and screw)	BSAB1585	35			10-12
SET ZZ BASE S&M-K-PCTBH-K85 3.3 (Parallel cemented titan base hex K85 and screw)	BSAB1581K85	30			10-12
SET ZZ BASE S&M-K-PCTBH-K85 3.8 (Parallel cemented titan base hex K85 and screw)	BSAB1582K85	30			10-12
SET ZZ BASE S&M-K-PCTBH-K85 4.25 (Parallel cemented titan base hex K85 and screw)	BSAB1583K85	35			10-12
SET ZZ BASE S&M-K-PCTBH-K85 5.0 (Parallel cemented titan base hex K85 and screw)	BSAB1584K85	35			10-12

All information is subject to change. Errors and omissions excepted. Version: 21.09.2023

SWEDEN & MARTINA PREMIUM KOHNO®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE S&M-K-PCTBH-L10 6.0 (Parallel cemented titan base hex L10 and screw)	BSAB1585L10	35			10-12
SET ZZ BASE S&M-K-PCTBH-L20 3.3 (Parallel cemented titan base hex L20 and screw)	BSAB1581L20	30			10-12
SET ZZ BASE S&M-K-PCTBH-L20 3.8 (Parallel cemented titan base hex L20 and screw)	BSAB1582L20	30			10-12
SET ZZ BASE S&M-K-PCTBH-L20 4.25 (Parallel cemented titan base hex L20 and screw)	BSAB1583L20	35			10-12
SET ZZ BASE S&M-K-PCTBH-L20 5.0 (Parallel cemented titan base hex L20 and screw)	BSAB1584L20	35			10-12
SET ZZ BASE S&M-K-PCTBH-L20 6.0 (Parallel cemented titan base hex L20 and screw)	BSAB1585L20	35			10-12
SET ZZ BASE S&M-K-PCTBH-L30 3.3 (Parallel cemented titan base hex L30 and screw)	BSAB1581L30	30			10-12
SET ZZ BASE S&M-K-PCTBH-L30 3.8 (Parallel cemented titan base hex L30 and screw)	BSAB1582L30	30			10-12
SET ZZ BASE S&M-K-PCTBH-L30 4.25 (Parallel cemented titan base hex L30 and screw)	BSAB1583L30	35			10-12
SET ZZ BASE S&M-K-PCTBH-L30 5.0 (Parallel cemented titan base hex L30 and screw)	BSAB1584L30	35			10-12
SET ZZ BASE S&M-K-PCTBH-L30 6.0 (Parallel cemented titan base hex L30 and screw)	BSAB1585L30	35			10-12
SET ZZ BASE S&M-K-PCTBH-L40 3.3 (Parallel cemented titan base hex L40 and screw)	BSAB1581L40	30			10-12
SET ZZ BASE S&M-K-PCTBH-L40 3.8 (Parallel cemented titan base hex L40 and screw)	BSAB1582L40	30			10-12
SET ZZ BASE S&M-K-PCTBH-L40 4.25 (Parallel cemented titan base hex L40 and screw)	BSAB1583L40	35			10-12
SET ZZ BASE S&M-K-PCTBH-L40 5.0 (Parallel cemented titan base hex L40 and screw)	BSAB1584L40	35			10-12
SET ZZ BASE S&M-K-PCTBH-L40 6.0 (Parallel cemented titan base hex L40 and screw)	BSAB1585L40	35			10-12
SET ZZ BASE S&M-K-SA 3.3 (Scanalog and screw)	ZBDE1531	15			10-12
SET ZZ BASE S&M-K-SA 3.8 (Scanalog and screw)	ZBDE1532	15			10-12
SET ZZ BASE S&M-K-SA 4.25 (Scanalog and screw)	ZBDE1533	15			10-12
SET ZZ BASE S&M-K-SA 5.0 (Scanalog and screw)	ZBDE1534	15			10-12
SET ZZ BASE S&M-K-SM 3.3 (Scanmarker and screw)	ZBAD1531	30			10-12
SET ZZ BASE S&M-K-SM 3.8 (Scanmarker and screw)	ZBAD1532	30			10-12
SET ZZ BASE S&M-K-SM 5.0 (Scanmarker and screw)	ZBAD1534	35			10-12
SET ZZ BASE S&M-K-SM 6.0 (Scanmarker and screw)	ZBAD1535	35			10-12
SET ZZ BASE S&M-K-WSM 3.30 (White Scanmarker and screw)	ZBCD1531	30			10-12
SET ZZ BASE S&M-K-WSM 3.80 (White Scanmarker and screw)	ZBCD1532	30			10-12
SET ZZ BASE S&M-K-WSM 4.25 (White Scanmarker and screw)	ZBCD1533	35			10-12
SET ZZ BASE S&M-K-WSM 5.0 (White Scanmarker and screw)	ZBCD1534	35			10-12
SET ZZ BASE S&M-K-WSM 6.0 (White Scanmarker and screw)	ZBCD1535	35			10-12

SWEDEN & MARTINA PREMIUM KOHNO®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE S&M-K-SC 3.30 (Screw)	BSAE1501	30			10-12
ZZ BASE S&M-K-SC 3.80 (Screw)	BSAE1502	30			10-12
ZZ BASE S&M-K-SC 4.25 (Screw)	BSAE1503	35			10-12
ZZ BASE S&M-K-ZSC 3.30 (Zirconia screw)	BSAE1531	30			10-12
ZZ BASE S&M-K-ZSC 3.80 (Zirconia screw)	BSAE1532	30			10-12
ZZ BASE S&M-K-ZSC 4.25 (Zirconia screw)	BSAE1533	35			10-12
ZZ BASE S&M-K-ASCSC 3.3 (Angled screw-channel Screw)	BSKE1501	20			10-12
ZZ BASE S&M-K-ASCSC 3.8 (Angled screw-channel Screw)	BSKE1502	20			10-12
ZZ BASE S&M-K-ASCSC 4.25 (Angled screw-channel Screw)	BSKE1503	20			10-12
ZZ BASE S&M-K-ICSC 3.3 (Impression coping Screw)	BSFE1501	15			10-12
ZZ BASE S&M-K-ICSC 3.8 (Impression coping Screw)	BSFE1502	15			10-12
ZZ BASE S&M-K-ICSC 4.25 (Impression coping Screw)	BSFE1503	15			10-12
ZZ BASE S&M-K-ICSC 5.0 (Impression coping Screw)	BSFE1504	15			10-12
ZZ BASE S&M-K-ICSC 6.0 (Impression coping Screw)	BSFE1505	15			10-12
ZZ BASE S&M-K-SC 5.0 (Screw)	BSAE1504	35			10-12
ZZ BASE S&M-K-SC 6.0 (Screw)	BSAE1505	35			10-12
ZZ BASE S&M-K-SC-G 3.3 (Screw Gold)	BSCE1501	30			10-12
ZZ BASE S&M-K-SC-G 3.8 (Screw Gold)	BSCE1502	30			10-12
ZZ BASE S&M-K-SC-G 4.25 (Screw Gold)	BSCE1503	35			10-12
ZZ BASE S&M-K-SC-G 5.0 (Screw Gold)	BSCE1504	35			10-12
ZZ BASE S&M-K-SC-G 6.0 (Screw Gold)	BSCE1505	35			10-12
ZZ BASE S&M-K-WMSM-L10 3.8 (White Metal Scanmarker L10 with Screw)	ZBED1532L10	15			10-12
ZZ BASE S&M-K-WMSM-L12 3.3 (White Metal Scanmarker L12 with Screw)	ZBED1531L12	15			10-12
ZZ BASE S&M-K-WMSM-L12 3.8 (White Metal Scanmarker L12 with Screw)	ZBED1532L12	15			10-12
ZZ BASE S&M-K-WMSM-L12 4.25 (White Metal Scanmarker L12 with Screw)	ZBED1533L12	15			10-12
ZZ BASE S&M-K-WMSM-L12 5.0 (White Metal Scanmarker L12 with Screw)	ZBED1534L12	15			10-12
ZZ BASE S&M-K-WMSM-L12 6.0 (White Metal Scanmarker L12 with Screw)	ZBED1535L12	15			10-12
ZZ BASE S&M-K-ZSC 5.0 (Zirconia screw)	BSAE1534	35			10-12
ZZ BASE S&M-K-ZSC 6.0 (Zirconia screw)	BSAE1535	35			10-12

SWEDEN & MARTINA PREMIUM KOHNO®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE S&M-K-ZSC-G 3.3 (Zirconia Screw Gold)	BSCE1531	30			10-12
ZZ BASE S&M-K-ZSC-G 3.8 (Zirconia Screw Gold)	BSCE1532	30			10-12
ZZ BASE S&M-K-ZSC-G 4.25 (Zirconia Screw Gold)	BSCE1533	35			10-12
ZZ BASE S&M-K-ZSC-G 5.0 (Zirconia Screw Gold)	BSCE1534	35			10-12
ZZ BASE S&M-K-ZSC-G 6.0 (Zirconia Screw Gold)	BSCE1535	35			10-12

THOMMEN MEDICAL SPI®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE T-SPI-ABCH 3.5 (Raw-Abutment® customizable hex and screw)	BSAF3491	15			10-12
SET ZZ BASE T-SPI-ABCH 4.0 (Raw-Abutment® customizable hex and screw)	BSAF3492	25			10-12
SET ZZ BASE T-SPI-ABCH 4.5 (Raw-Abutment® customizable hex and screw)	BSAF3493	25			10-12
SET ZZ BASE T-SPI-ABCH 5.0 (Raw-Abutment® customizable hex and screw)	BSAF3494	25			10-12
SET ZZ BASE T-SPI-ABCH 6.0 (Raw-Abutment® customizable hex and screw)	BSAF3495	25			10-12
SET ZZ BASE T-SPI-ABCH-D14 3.5 (Raw-Abutment® customizable hex D14 and screw)	BSAF3491D14	15			10-12
SET ZZ BASE T-SPI-ABCH-D14 4.0 (Raw-Abutment® customizable hex D14 and screw)	BSAF3492D14	25			10-12
SET ZZ BASE T-SPI-ASCTBH-K80 3.5 (Angled screw-channel titan base hex K80 and screw)	BSAT3481K80	15			10-12
SET ZZ BASE T-SPI-ASCTBH-K80 4.0 (Angled screw-channel titan base hex K80 and screw)	BSAT3482K80	20			10-12
SET ZZ BASE T-SPI-ASCTBH-K80 4.5 (Angled screw-channel titan base hex K80 and screw)	BSAT3483K80	20			10-12
SET ZZ BASE T-SPI-ASCTB-K80 3.5 (Angled screw-channel titan base non hex K80 and screw)	BSAT3431K80	15			10-12
SET ZZ BASE T-SPI-ASCTB-K80 4.0 (Angled screw-channel titan base non hex K80 and screw)	BSAT3432K80	20			10-12
SET ZZ BASE T-SPI-ASCTB-K80 4.5 (Angled screw-channel titan base non hex K80 and screw)	BSAT3433K80	20			10-12
SET ZZ BASE T-SPI-CCTB 3.5 (Conical cemented titan base non hex and screw)	BSAD3431	15			10-12
SET ZZ BASE T-SPI-CCTB 4.0 (Conical cemented titan base non hex and screw)	BSAD3432	25			10-12
SET ZZ BASE T-SPI-CCTB 4.5 (Conical cemented titan base non hex and screw)	BSAD3433	25			10-12
SET ZZ BASE T-SPI-CCTB 5.0 (Conical cemented titan base non hex and screw)	BSAD3434	25			10-12
SET ZZ BASE T-SPI-CCTB 6.0 (Conical cemented titan base non hex and screw)	BSAD3435	25			10-12
SET ZZ BASE T-SPI-CCTB-K85 3.5 (Conical cemented titan base non hex K85 and screw)	BSAD3431K85	15			10-12
SET ZZ BASE T-SPI-CCTB-K85 4.0 (Conical cemented titan base non hex K85 and screw)	BSAD3432K85	25			10-12
SET ZZ BASE T-SPI-CCTB-K85 4.5 (Conical cemented titan base non hex K85 and screw)	BSAD3433K85	25			10-12
SET ZZ BASE T-SPI-CCTB-K85 5.0 (Conical cemented titan base non hex K85 and screw)	BSAD3434K85	25			10-12
SET ZZ BASE T-SPI-CCTB-K85 6.0 (Conical cemented titan base non hex K85 and screw)	BSAD3435K85	25			10-12
SET ZZ BASE T-SPI-CCTB-L20 3.5 (Conical cemented titan base non hex L20 and screw)	BSAD3431L20	15			10-12
SET ZZ BASE T-SPI-CCTB-L20 4.0 (Conical cemented titan base non hex L20 and screw)	BSAD3432L20	25			10-12
SET ZZ BASE T-SPI-CCTB-L20 4.5 (Conical cemented titan base non hex L20 and screw)	BSAD3433L20	25			10-12
SET ZZ BASE T-SPI-CCTB-L20 5.0 (Conical cemented titan base non hex L20 and screw)	BSAD3434L20	25			10-12
SET ZZ BASE T-SPI-CCTB-L20 6.0 (Conical cemented titan base non hex L20 and screw)	BSAD3435L20	25			10-12
SET ZZ BASE T-SPI-CCTB-L30 3.5 (Conical cemented titan base non hex L30 and screw)	BSAD3431L30	15			10-12

THOMMEN MEDICAL SPI®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE T-SPI-CCTB-L30 4.0 (Conical cemented titan base non hex L30 and screw)	BSAD3432L30	25			10-12
SET ZZ BASE T-SPI-CCTB-L30 4.5 (Conical cemented titan base non hex L30 and screw)	BSAD3433L30	25			10-12
SET ZZ BASE T-SPI-CCTB-L30 5.0 (Conical cemented titan base non hex L30 and screw)	BSAD3434L30	25			10-12
SET ZZ BASE T-SPI-CCTB-L30 6.0 (Conical cemented titan base non hex L30 and screw)	BSAD3435L30	25			10-12
SET ZZ BASE T-SPI-CCTB-L40 3.5 (Conical cemented titan base non hex L40 and screw)	BSAD3431L40	15			10-12
SET ZZ BASE T-SPI-CCTB-L40 4.0 (Conical cemented titan base non hex L40 and screw)	BSAD3432L40	25			10-12
SET ZZ BASE T-SPI-CCTB-L40 4.5 (Conical cemented titan base non hex L40 and screw)	BSAD3433L40	25			10-12
SET ZZ BASE T-SPI-CCTB-L40 5.0 (Conical cemented titan base non hex L40 and screw)	BSAD3434L40	25			10-12
SET ZZ BASE T-SPI-CCTB-L40 6.0 (Conical cemented titan base non hex L40 and screw)	BSAD3435L40	25			10-12
SET ZZ BASE T-SPI-ICOTH 3.5 (Impression coping open tray hex and Screw and Driver)	BSAQ3441	15			10-12
SET ZZ BASE T-SPI-ICOTH 4.0 (Impression coping open tray hex and Screw and Driver)	BSAQ3442	15			10-12
SET ZZ BASE T-SPI-ICOTH 4.5 (Impression coping open tray hex and Screw and Driver)	BSAQ3443	15			10-12
SET ZZ BASE T-SPI-ICOTH 5.0 (Impression coping open tray hex and Screw and Driver)	BSAQ3444	15			10-12
SET ZZ BASE T-SPI-ICOTH 6.0 (Impression coping open tray hex and Screw and Driver)	BSAQ3445	15			10-12
SET ZZ BASE T-SPI-PCTBH 3.5 (Parallel cemented titan base hex and screw)	BSAB3481	15			10-12
SET ZZ BASE T-SPI-PCTBH 4.0 (Parallel cemented titan base hex and screw)	BSAB3482	25			10-12
SET ZZ BASE T-SPI-PCTBH 4.5 (Parallel cemented titan base hex and screw)	BSAB3483	25			10-12
SET ZZ BASE T-SPI-PCTBH 5.0 (Parallel cemented titan base hex and screw)	BSAB3484	25			10-12
SET ZZ BASE T-SPI-PCTBH 6.0 (Parallel cemented titan base hex and screw)	BSAB3485	25			10-12
SET ZZ BASE T-SPI-PCTBH-K85 3.5 (Parallel cemented titan base hex K85 and screw)	BSAB3481K85	15			10-12
SET ZZ BASE T-SPI-PCTBH-K85 4.0 (Parallel cemented titan base hex K85 and screw)	BSAB3482K85	25			10-12
SET ZZ BASE T-SPI-PCTBH-K85 4.5 (Parallel cemented titan base hex K85 and screw)	BSAB3483K85	25			10-12
SET ZZ BASE T-SPI-PCTBH-K85 5.0 (Parallel cemented titan base hex K85 and screw)	BSAB3484K85	25			10-12
SET ZZ BASE T-SPI-PCTBH-K85 6.0 (Parallel cemented titan base hex K85 and screw)	BSAB3485K85	25			10-12
SET ZZ BASE T-SPI-PCTBH-L20 3.5 (Parallel cemented titan base hex L20 and screw)	BSAB3491L20	15			10-12
SET ZZ BASE T-SPI-PCTBH-L20 4.0 (Parallel cemented titan base hex L20 and screw)	BSAB3492L20	25			10-12
SET ZZ BASE T-SPI-PCTBH-L20 4.5 (Parallel cemented titan base hex L20 and screw)	BSAB3493L20	25			10-12
SET ZZ BASE T-SPI-PCTBH-L20 5.0 (Parallel cemented titan base hex L20 and screw)	BSAB3494L20	25			10-12
SET ZZ BASE T-SPI-PCTBH-L20 6.0 (Parallel cemented titan base hex L20 and screw)	BSAB3495L20	25			10-12

THOMMEN MEDICAL SPI®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE T-SPI-PCTBH-L30 3.5 (Parallel cemented titan base hex L30 and screw)	BSAB3491L30	15			10-12
SET ZZ BASE T-SPI-PCTBH-L30 4.0 (Parallel cemented titan base hex L30 and screw)	BSAB3492L30	25			10-12
SET ZZ BASE T-SPI-PCTBH-L30 4.5 (Parallel cemented titan base hex L30 and screw)	BSAB3493L30	25			10-12
SET ZZ BASE T-SPI-PCTBH-L30 5.0 (Parallel cemented titan base hex L30 and screw)	BSAB3494L30	25			10-12
SET ZZ BASE T-SPI-PCTBH-L30 6.0 (Parallel cemented titan base hex L30 and screw)	BSAB3495L30	25			10-12
SET ZZ BASE T-SPI-PCTBH-L40 3.5 (Parallel cemented titan base hex L40 and screw)	BSAB3491L40	15			10-12
SET ZZ BASE T-SPI-PCTBH-L40 4.0 (Parallel cemented titan base hex L40 and screw)	BSAB3492L40	25			10-12
SET ZZ BASE T-SPI-PCTBH-L40 4.5 (Parallel cemented titan base hex L40 and screw)	BSAB3493L40	25			10-12
SET ZZ BASE T-SPI-PCTBH-L40 5.0 (Parallel cemented titan base hex L40 and screw)	BSAB3494L40	25			10-12
SET ZZ BASE T-SPI-PCTBH-L40 6.0 (Parallel cemented titan base hex L40 and screw)	BSAB3495L40	25			10-12
SET ZZ BASE T-SPI-SA 3.5 (Scanalog and screw)	ZBDE3431	15			10-12
SET ZZ BASE T-SPI-SA 4.0 (Scanalog and screw)	ZBDE3432	15			10-12
SET ZZ BASE T-SPI-SA 4.5 (Scanalog and screw)	ZBDE3433	15			10-12
SET ZZ BASE T-SPI-SA 5.0 (Scanalog and screw)	ZBDE3434	15			10-12
SET ZZ BASE T-SPI-SA 6.0 (Scanalog and screw)	ZBDE3435	15			10-12
SET ZZ BASE T-SPI-SM 3.5 (Scanmarker and screw)	ZBAD3431	15			10-12
SET ZZ BASE T-SPI-SM 4.0 (Scanmarker and screw)	ZBAD3432	25			10-12
SET ZZ BASE T-SPI-SM 4.5 (Scanmarker and screw)	ZBAD3433	25			10-12
SET ZZ BASE T-SPI-SM 5.0 (Scanmarker and screw)	ZBAD3434	25			10-12
SET ZZ BASE T-SPI-SM 6.0 (Scanmarker and screw)	ZBAD3435	25			10-12
SET ZZ BASE T-SPI-WSM 3.5 (White Scanmarker and screw)	ZBCD3431	15			10-12
SET ZZ BASE T-SPI-WSM 4.0 (White Scanmarker and screw)	ZBCD3432	25			10-12
SET ZZ BASE T-SPI-WSM 4.5 (White Scanmarker and screw)	ZBCD3433	25			10-12
SET ZZ BASE T-SPI-WSM 5.0 (White Scanmarker and screw)	ZBCD3434	25			10-12
SET ZZ BASE T-SPI-WSM 6.0 (White Scanmarker and screw)	ZBCD3435	25			10-12
ZZ BASE T-SPI-ASCSC 3.5 (Angled screw-channel Screw)	BSKE3401	15			10-12
ZZ BASE T-SPI-ASCSC 4.0 (Angled screw-channel Screw)	BSKE3402	20			10-12
ZZ BASE T-SPI-ASCSC 4.5 (Angled screw-channel Screw)	BSKE3403	20			10-12
ZZ BASE T-SPI-ICSC 3.5 (Impression coping Screw)	BSFE3401	15			10-12

THOMMEN MEDICAL SPI®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE T-SPI-ICSC 4.0 (Impression coping Screw)	BSFE3402	15			10-12
ZZ BASE T-SPI-ICSC 4.5 (Impression coping Screw)	BSFE3403	15			10-12
ZZ BASE T-SPI-ICSC 5.0 (Impression coping Screw)	BSFE3404	15			10-12
ZZ BASE T-SPI-ICSC 6.0 (Impression coping Screw)	BSFE3405	15			10-12
ZZ BASE T-SPI-MI 3.5 (Model implant)	MIAG3401				10-12
ZZ BASE T-SPI-MI 4.0 (Model implant)	MIAG3402				10-12
ZZ BASE T-SPI-MI 4.5 (Model implant)	MIAG3403				10-12
ZZ BASE T-SPI-MI 5.0 (Model implant)	MIAG3404				10-12
ZZ BASE T-SPI-MI 6.0 (Model implant)	MIAG3405				10-12
ZZ BASE T-SPI-SC 3.5 (Screw)	BSAE3401	15			10-12
ZZ BASE T-SPI-SC 4.0 (Screw)	BSAE3402	25			10-12
ZZ BASE T-SPI-SC 4.5 (Screw)	BSAE3403	25			10-12
ZZ BASE T-SPI-SC 5.0 (Screw)	BSAE3404	25			10-12
ZZ BASE T-SPI-SC 6.0 (Screw)	BSAE3405	25			10-12
ZZ BASE T-SPI-WMSM-L12 3.5 (White Metal Scanmarker L12 with Screw)	ZBED3431L12	15			10-12
ZZ BASE T-SPI-WMSM-L12 4.0 (White Metal Scanmarker L12 with Screw)	ZBED3432L12	15			10-12
ZZ BASE T-SPI-WMSM-L12 4.5 (White Metal Scanmarker L12 with Screw)	ZBED3433L12	15			10-12
ZZ BASE T-SPI-WMSM-L12 5.0 (White Metal Scanmarker L12 with Screw)	ZBED3434L12	15			10-12
ZZ BASE T-SPI-WMSM-L12 6.0 (White Metal Scanmarker L12 with Screw)	ZBED3435L12	15			10-12
ZZ BASE T-SPI-ZSC 3.5 (Zirconia screw)	BSAE3431	25			10-12
ZZ BASE T-SPI-ZSC 4.0 (Zirconia screw)	BSAE3432	25			10-12
ZZ BASE T-SPI-ZSC 4.5 (Zirconia screw)	BSAE3433	25			10-12
ZZ BASE T-SPI-ZSC 5.0 (Zirconia screw)	BSAE3434	25			10-12
ZZ BASE T-SPI-ZSC 6.0 (Zirconia screw)	BSAE3435	25			10-12

THOMMEN MEDICAL SPI® VARIOMULTI

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE T-VM-CCTB 4.8 (Conical cemented titan base non hex and screw)	BSAD6131	25			10-12
SET ZZ BASE T-VM-SA 4.8 (Scananalog and screw)	ZBDE6131	15			10-12
SET ZZ BASE T-VM-SM 4.8 (Scanmarker and screw)	ZBAD6131	25			10-12
SET ZZ BASE T-VM-WSM 4.8 (White Scanmarker and screw)	ZBCD6131	25			10-12
ZZ BASE T-VM-SC 4.8 (Screw)	BSAE6101	25			10-12
ZZ BASE T-VM-ZSC 4.8 (Zirconia screw)	BSAE6131	25			10-12

VARIOBASE® FOR CROWN

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
Straumann® Basal Screw NC	ZZ-025-2900	35			10-12
Straumann® Basal Screw NNC	ZZ-048-313	35			10-12
Straumann® Basal Screw RC	ZZ-025-4900	35			10-12
Straumann® Basal Screw Synocta®	ZZ-048-356	35			10-12
Straumann® NC Variobase® Abutment for crown AH 3.5 with screw	ZZ-025-2921	35			10-12
Straumann® NC Variobase® for crown AH 5.5 with screw	ZZ-022-0027	35			10-12
Straumann® NC Variobase® for Crown incl. screw, Ø 3.8mm, AH 3.5mm, GH 2mm, TAN	ZZ-022-0102	35			10-12
Straumann® NC Variobase® for Crown incl. screw, Ø 3.9mm, AH 3.5mm, GH 3mm, TAN	ZZ-022-0104	35			10-12
Straumann® NC Variobase® for Crown incl. screw, Ø 3.9mm, AH 5.5mm, GH 2mm, TAN	ZZ-022-0106	35			10-12
Straumann® NC Variobase® for Crown incl. screw, Ø 3.9mm, AH 5.5mm, GH 3mm, TAN	ZZ-022-0108	35			10-12
Straumann® NNC Variobase® Abutment for crown AH 3.75 with screw	ZZ-048-709	35			10-12
Straumann® NNC Variobase® for crown AH 5.75 with screw	ZZ-022-0021	35			10-12
Straumann® RC Variobase® Abutment for crown AH 3.5 with screw	ZZ-025-4921	35			10-12
Straumann® RC Variobase® for crown AH 5.5 with screw	ZZ-022-0026	35			10-12
Straumann® RC Variobase® for Crown incl. screw, Ø 4.5mm, AH 3.5mm, GH 2mm, TAN	ZZ-022-0103	35			10-12
Straumann® RC Variobase® for Crown incl. screw, Ø 4.5mm, AH 3.5mm, GH 3mm, TAN	ZZ-022-0105	35			10-12
Straumann® RC Variobase® for Crown incl. screw, Ø 4.5mm, AH 5.5mm, GH 2mm, TAN	ZZ-022-0107	35			10-12
Straumann® RC Variobase® for Crown incl. screw, Ø 4.5mm, AH 5.5mm, GH 3mm, TAN	ZZ-022-0109	35			10-12
Straumann® RN Variobase® Abutment for crown AH 4.0 with screw	ZZ-048-710	35			10-12
Straumann® RN Variobase® for crown AH 6 with screw	ZZ-022-0022	35			10-12
Straumann® WN Variobase® Abutment for crown AH 4.5 with screw	ZZ-048-711	35			10-12
Straumann® WN Variobase® for crown AH 6.5 with screw	ZZ-022-0023	35			10-12

WARANTEC ONEPLANT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
SET ZZ BASE WT-OI-ABCH 3.6-5.3 (Raw-Abutment® customizable hex and screw)	BSAF7481	30			10-12
ZZ BASE WT-OI-SC 3.6-5.3 (Screw)	BSAE7401	30			10-12

ZIMMER DENTAL® TAPERED SCREW-VENT®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE Z-TS-ABCH 3.5 (Raw-Abutment® customizable hex and screw)	BSAF1181	30			10-12
SET ZZ BASE Z-TS-ABCH 4.5 (Raw-Abutment® customizable hex and screw)	BSAF1182	30			10-12
SET ZZ BASE Z-TS-ABCH 5.7 (Raw-Abutment® customizable hex and screw)	BSAF1183	30			10-12
SET ZZ BASE Z-TS-ABCH-D14 3.5 (Raw-Abutment® customizable hex D14 and screw)	BSAF1181D14	30			10-12
SET ZZ BASE Z-TS-ABCH-D14 4.5 (Raw-Abutment® customizable hex D14 and screw)	BSAF1182D14	30			10-12
SET ZZ BASE Z-TS-ABCH-D14 5.7 (Raw-Abutment® customizable hex D14 and screw)	BSAF1183D14	30			10-12
SET ZZ BASE Z-TS-ASCTBH-K80 3.5 (Angled screw-channel titan base hex K80 and screw)	BSAT1181K80	20			10-12
SET ZZ BASE Z-TS-ASCTBH-K80 4.5 (Angled screw-channel titan base hex K80 and screw)	BSAT1182K80	20			10-12
SET ZZ BASE Z-TS-ASCTB-K80 3.5 (Angled screw-channel titan base non hex K80 and screw)	BSAT1131K80	20			10-12
SET ZZ BASE Z-TS-ASCTB-K80 4.5 (Angled screw-channel titan base non hex K80 and screw)	BSAT1132K80	20			10-12
SET ZZ BASE Z-TS-CCTB 3.5 (Conical cemented titan base non hex and screw)	BSAD1131	30			10-12
SET ZZ BASE Z-TS-CCTB 4.5 (Conical cemented titan base non hex and screw)	BSAD1132	30			10-12
SET ZZ BASE Z-TS-CCTB 5.7 (Conical cemented titan base non hex and screw)	BSAD1133	30			10-12
SET ZZ BASE Z-TS-CCTB-G 3.5 (Conical cemented titan base non hex and screw Gold)	BSCD1131	30			10-12
SET ZZ BASE Z-TS-CCTB-K85 3.5 (Conical cemented titan base non hex K85 and screw)	BSAD1131K85	30			10-12
SET ZZ BASE Z-TS-CCTB-K85 4.5 (Conical cemented titan base non hex K85 and screw)	BSAD1132K85	30			10-12
SET ZZ BASE Z-TS-CCTB-K85 5.7 (Conical cemented titan base non hex K85 and screw)	BSAD1133K85	30			10-12
SET ZZ BASE Z-TS-CCTB-L10 3.5 (Conical cemented titan base non hex L10 and screw)	BSAD1131L10	30			10-12
SET ZZ BASE Z-TS-CCTB-L10 4.5 (Conical cemented titan base non hex L10 and screw)	BSAD1132L10	30			10-12
SET ZZ BASE Z-TS-CCTB-L10 5.7 (Conical cemented titan base non hex L10 and screw)	BSAD1133L10	30			10-12
SET ZZ BASE Z-TS-CCTB-L10-G 3.5 (Conical cemented titan base non hex L10 and screw Gold)	BSCD1131L10	30			10-12
SET ZZ BASE Z-TS-CCTB-L20 3.5 (Conical cemented titan base non hex L20 and screw)	BSAD1131L20	30			10-12
SET ZZ BASE Z-TS-CCTB-L20 4.5 (Conical cemented titan base non hex L20 and screw)	BSAD1132L20	30			10-12
SET ZZ BASE Z-TS-CCTB-L20 5.7 (Conical cemented titan base non hex L20 and screw)	BSAD1133L20	30			10-12
SET ZZ BASE Z-TS-CCTB-L20-G 3.5 (Conical cemented titan base non hex L20 and screw Gold)	BSCD1131L20	30			10-12
SET ZZ BASE Z-TS-CCTB-L30 3.5 (Conical cemented titan base non hex L30 and screw)	BSAD1131L30	30			10-12
SET ZZ BASE Z-TS-CCTB-L30 4.5 (Conical cemented titan base non hex L30 and screw)	BSAD1132L30	30			10-12
SET ZZ BASE Z-TS-CCTB-L30 5.7 (Conical cemented titan base non hex L30 and screw)	BSAD1133L30	30			10-12
SET ZZ BASE Z-TS-CCTB-L30-G 3.5 (Conical cemented titan base non hex L30 and screw Gold)	BSCD1131L30	30			10-12

ZIMMER DENTAL® TAPERED SCREW-VENT®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE Z-TS-CCTB-L40 3.5 (Conical cemented titan base non hex L40 and screw)	BSAD1131L40	30			10-12
SET ZZ BASE Z-TS-CCTB-L40 4.5 (Conical cemented titan base non hex L40 and screw)	BSAD1132L40	30			10-12
SET ZZ BASE Z-TS-CCTB-L40 5.7 (Conical cemented titan base non hex L40 and screw)	BSAD1133L40	30			10-12
SET ZZ BASE Z-TS-ICOT 3.5 (Impression coping open tray non hex and Screw and Driver)	BSAQ1161	15			10-12
SET ZZ BASE Z-TS-ICOT 4.5 (Impression coping open tray non hex and Screw and Driver)	BSAQ1162	15			10-12
SET ZZ BASE Z-TS-ICOT 5.7 (Impression coping open tray non hex and Screw and Driver)	BSAQ1163	15			10-12
SET ZZ BASE Z-TS-ICOTH 3.5 (Impression coping open tray hex and Screw and Driver)	BSAQ1141	15			10-12
SET ZZ BASE Z-TS-ICOTH 4.5 (Impression coping open tray hex and Screw and Driver)	BSAQ1142	15			10-12
SET ZZ BASE Z-TS-ICOTH 5.7 (Impression coping open tray hex and Screw and Driver)	BSAQ1143	15			10-12
SET ZZ BASE Z-TS-NTB 3.5 (Narrow titan base non hex and screw)	BSAP1131	30			10-12
SET ZZ BASE Z-TS-NTBH 3.5 (Narrow titan base hex and screw)	BSAP1181	30			10-12
SET ZZ BASE Z-TS-PCTBH 3.5 (Parallel cemented titan base hex and screw)	BSAB1181	30			10-12
SET ZZ BASE Z-TS-PCTBH 4.5 (Parallel cemented titan base hex and screw)	BSAB1182	30			10-12
SET ZZ BASE Z-TS-PCTBH 5.7 (Parallel cemented titan base hex and screw)	BSAB1183	30			10-12
SET ZZ BASE Z-TS-PCTBH-G 3.5 (Parallel cemented titan base hex and screw Gold)	BSCB1181	30			10-12
SET ZZ BASE Z-TS-PCTBH-G 4.5 (Parallel cemented titan base hex and screw Gold)	BSCB1182	30			10-12
SET ZZ BASE Z-TS-PCTBH-K85 3.5 (Parallel cemented titan base hex K85 and screw)	BSAB1181K85	30			10-12
SET ZZ BASE Z-TS-PCTBH-K85 4.5 (Parallel cemented titan base hex K85 and screw)	BSAB1182K85	30			10-12
SET ZZ BASE Z-TS-PCTBH-K85 5.7 (Parallel cemented titan base hex K85 and screw)	BSAB1183K85	30			10-12
SET ZZ BASE Z-TS-PCTBH-L10 3.5 (Parallel cemented titan base hex L10 and screw)	BSAB1191L10	30			10-12
SET ZZ BASE Z-TS-PCTBH-L10 4.5 (Parallel cemented titan base hex L10 and screw)	BSAB1192L10	30			10-12
SET ZZ BASE Z-TS-PCTBH-L10 5.7 (Parallel cemented titan base hex L10 and screw)	BSAB1193L10	30			10-12
SET ZZ BASE Z-TS-PCTBH-L10-G 3.5 (Parallel cemented titan base hex L10 and screw Gold)	BSCB1191L10	30			10-12
SET ZZ BASE Z-TS-PCTBH-L20 3.5 (Parallel cemented titan base hex L20 and screw)	BSAB1191L20	30			10-12
SET ZZ BASE Z-TS-PCTBH-L20 4.5 (Parallel cemented titan base hex L20 and screw)	BSAB1192L20	30			10-12
SET ZZ BASE Z-TS-PCTBH-L20 5.7 (Parallel cemented titan base hex L20 and screw)	BSAB1193L20	30			10-12
SET ZZ BASE Z-TS-PCTBH-L20-G 3.5 (Parallel cemented titan base hex L20 and screw Gold)	BSCB1191L20	30			10-12
SET ZZ BASE Z-TS-PCTBH-L30 3.5 (Parallel cemented titan base hex L30 and screw)	BSAB1191L30	30			10-12
SET ZZ BASE Z-TS-PCTBH-L30 4.5 (Parallel cemented titan base hex L30 and screw)	BSAB1192L30	30			10-12

ZIMMER DENTAL® TAPERED SCREW-VENT®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE Z-TS-PCTBH-L30 5.7 (Parallel cemented titan base hex L30 and screw)	BSAB1193L30	30			10-12
SET ZZ BASE Z-TS-PCTBH-L30-G 3.5 (Parallel cemented titan base hex L30 and screw Gold)	BSCB1191L30	30			10-12
SET ZZ BASE Z-TS-PCTBH-L40 3.5 (Parallel cemented titan base hex L40 and screw)	BSAB1191L40	30			10-12
SET ZZ BASE Z-TS-PCTBH-L40 4.5 (Parallel cemented titan base hex L40 and screw)	BSAB1192L40	30			10-12
SET ZZ BASE Z-TS-PCTBH-L40 5.7 (Parallel cemented titan base hex L40 and screw)	BSAB1193L40	30			10-12
SET ZZ BASE Z-TS-PCTBH-L40-G 3.5 (Parallel cemented titan base hex L40 and screw Gold)	BSCB1191L40	30			10-12
SET ZZ BASE Z-TS-SA 3.5 (Scanalog and screw)	ZBDE1131	15			10-12
SET ZZ BASE Z-TS-SA 4.5 (Scanalog and screw)	ZBDE1132	15			10-12
SET ZZ BASE Z-TS-SA 5.7 (Scanalog and screw)	ZBDE1133	15			10-12
SET ZZ BASE Z-TS-SM 3.5 (Scanmarker and screw)	ZBAD1131	30			10-12
SET ZZ BASE Z-TS-SM 4.5 (Scanmarker and screw)	ZBAD1132	30			10-12
SET ZZ BASE Z-TS-SM 5.7 (Scanmarker and screw)	ZBAD1133	30			10-12
SET ZZ BASE Z-TS-WSM 5.7 (White Scanmarker and screw)	ZBCD1133	30			10-12
SET ZZ BASE Z-TS-ZZMUA17 3.5 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and implan	BSAI1131	25			10-12
SET ZZ BASE Z-TS-ZZMUA17 3.5 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and implan	BSAI1141	25			10-12
SET ZZ BASE Z-TS-ZZMUA17 4.5 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and implan	BSAI1132	25			10-12
SET ZZ BASE Z-TS-ZZMUA17 4.5 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and implan	BSAI1142	25			10-12
SET ZZ BASE Z-TS-ZZMUA17 5.7 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 and implan	BSAI1133	25			10-12
SET ZZ BASE Z-TS-ZZMUA17 5.7 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 and implan	BSAI1143	25			10-12
SET ZZ BASE Z-TS-ZZMUA17-L10 3.5 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10 ar	BSAI1131L10	25			10-12
SET ZZ BASE Z-TS-ZZMUA17-L10 3.5 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 L10 ar	BSAI1141L10	25			10-12
SET ZZ BASE Z-TS-ZZMUA17-L10 4.5 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10 ar	BSAI1132L10	25			10-12
SET ZZ BASE Z-TS-ZZMUA17-L10 4.5 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 L10 ar	BSAI1142L10	25			10-12
SET ZZ BASE Z-TS-ZZMUA17-L10 5.7 Type 1 (Zirkonzahn Multi Unit Abutment 17° Type 1 L10 ar	BSAI1133L10	25			10-12
SET ZZ BASE Z-TS-ZZMUA17-L10 5.7 Type 2 (Zirkonzahn Multi Unit Abutment 17° Type 2 L10 ar	BSAI1143L10	25			10-12
SET ZZ BASE Z-TS-ZZMUA30 3.5 Type 1 (Zirkonzahn Multi Unit Abutment 30° Type 1 and implan	BSAK1131	25			10-12
SET ZZ BASE Z-TS-ZZMUA30 3.5 Type 2 (Zirkonzahn Multi Unit Abutment 30° Type 2 and implan	BSAK1141	25			10-12
SET ZZ BASE Z-TS-ZZMUA30 4.5 Type 1 (Zirkonzahn Multi Unit Abutment 30° Type 1 and implan	BSAK1132	25			10-12
SET ZZ BASE Z-TS-ZZMUA30 4.5 Type 2 (Zirkonzahn Multi Unit Abutment 30° Type 2 and implan	BSAK1142	25			10-12

ZIMMER DENTAL® TAPERED SCREW-VENT®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE Z-TS-ZZMUA30 5.7 Type 1 (Zirkonzahn Multi Unit Abutment 30° Type 1 and implant)	BSAK1133	25			10-12
SET ZZ BASE Z-TS-ZZMUA30 5.7 Type 2 (Zirkonzahn Multi Unit Abutment 30° Type 2 and implant)	BSAK1143	25			10-12
SET ZZ BASE Z-TS-ZZMUA30-L10 3.5 Type 1 (Zirkonzahn Multi Unit Abutment 30° Type 1 L10 ar)	BSAK1131L10	25			10-12
SET ZZ BASE Z-TS-ZZMUA30-L10 3.5 Type 2 (Zirkonzahn Multi Unit Abutment 30° Type 2 L10 ar)	BSAK1141L10	25			10-12
SET ZZ BASE Z-TS-ZZMUA30-L10 4.5 Type 1 (Zirkonzahn Multi Unit Abutment 30° Type 1 L10 ar)	BSAK1132L10	25			10-12
SET ZZ BASE Z-TS-ZZMUA30-L10 4.5 Type 2 (Zirkonzahn Multi Unit Abutment 30° Type 2 L10 ar)	BSAK1142L10	25			10-12
SET ZZ BASE Z-TS-ZZMUA30-L10 5.7 Type 1 (Zirkonzahn Multi Unit Abutment 30° Type 1 L10 ar)	BSAK1133L10	25			10-12
SET ZZ BASE Z-TS-ZZMUA30-L10 5.7 Type 2 (Zirkonzahn Multi Unit Abutment 30° Type 2 L10 ar)	BSAK1143L10	25			10-12
SET ZZ BASE Z-TS-ZZMUA-L05 3.5 (Zirkonzahn Multi Unit Abutment non hex L05 and Screw)	BSAH1131L05	25	30		10-12
SET ZZ BASE Z-TS-ZZMUA-L15 3.5 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH1131L15	25	30		10-12
SET ZZ BASE Z-TS-ZZMUA-L15 4.5 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH1132L15	25	30		10-12
SET ZZ BASE Z-TS-ZZMUA-L15 5.7 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH1133L15	25	30		10-12
SET ZZ BASE Z-TS-ZZMUA-L25 3.5 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH1131L25	25	30		10-12
SET ZZ BASE Z-TS-ZZMUA-L25 4.5 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH1132L25	25	30		10-12
SET ZZ BASE Z-TS-ZZMUA-L25 5.7 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH1133L25	25	30		10-12
SET ZZ BASE Z-TS-ZZMUA-L35 3.5 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH1131L35	25	30		10-12
SET ZZ BASE Z-TS-ZZMUA-L35 4.5 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH1132L35	25	30		10-12
SET ZZ BASE Z-TS-ZZMUA-L35 5.7 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH1133L35	25	30		10-12
SET ZZ BASE Z-TS-ZZMUA-L45 3.5 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH1131L45	25	30		10-12
SET ZZ BASE Z-TS-ZZMUA-L45 4.5 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH1132L45	25	30		10-12
SET ZZ BASE Z-TS-ZZMUA-L45 5.7 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH1133L45	25	30		10-12
ZZ BASE Z-TS-ASCSC 3.5 (Angled screw-channel Screw)	BSKE1101	20			10-12
ZZ BASE Z-TS-ASCSC 4.5 (Angled screw-channel Screw)	BSKE1102	20			10-12
ZZ BASE Z-TS-BSC 3.5 (Black Screw)	BSIE1101	30			10-12
ZZ BASE Z-TS-BSC 4.5 (Black Screw)	BSIE1102	30			10-12
ZZ BASE Z-TS-BSC 5.7 (Black Screw)	BSIE1103	30			10-12
ZZ BASE Z-TS-BZSC 3.5 (Black Zirconia screw)	BSIE1131	30			10-12
ZZ BASE Z-TS-BZSC 4.5 (Black Zirconia screw)	BSIE1132	30			10-12
ZZ BASE Z-TS-BZSC 5.7 (Black Zirconia screw)	BSIE1133	30			10-12

ZIMMER DENTAL® TAPERED SCREW-VENT®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE Z-TS-ICSC 3.5 (Impression coping Screw)	BSFE1101	15			10-12
ZZ BASE Z-TS-ICSC 4.5 (Impression coping Screw)	BSFE1102	15			10-12
ZZ BASE Z-TS-ICSC 5.7 (Impression coping Screw)	BSFE1103	15			10-12
ZZ BASE Z-TS-LBSC 3.5 (Lab Screw)	BSHE1101	30			10-12
ZZ BASE Z-TS-LBSC 4.5 (Lab Screw)	BSHE1102	30			10-12
ZZ BASE Z-TS-LBSC 5.7 (Lab Screw)	BSHE1103	30			10-12
ZZ BASE Z-TS-LBZSC 3.5 (Lab Zirconia screw)	BSHE1131	30			10-12
ZZ BASE Z-TS-LBZSC 4.5 (Lab Zirconia screw)	BSHE1132	30			10-12
ZZ BASE Z-TS-LBZSC 5.7 (Lab Zirconia screw)	BSHE1133	30			10-12
ZZ BASE Z-TS-MI 3.5 (Model implant)	MIAG1101				10-12
ZZ BASE Z-TS-MI 4.5 (Model implant)	MIAG1102				10-12
ZZ BASE Z-TS-MI 5.7 (Model implant)	MIAG1103				10-12
ZZ BASE Z-TS-SC 3.5 (Screw)	BSAE1101	30			10-12
ZZ BASE Z-TS-SC 4.5 (Screw)	BSAE1102	30			10-12
ZZ BASE Z-TS-SC 5.7 (Screw)	BSAE1103	30			10-12
ZZ BASE Z-TS-SC-G 3.5 (Screw Gold)	BSCE1101	30			10-12
ZZ BASE Z-TS-SC-G 4.5 (Screw Gold)	BSCE1102	30			10-12
ZZ BASE Z-TS-SC-G 5.7 (Screw Gold)	BSCE1103	30			10-12
ZZ BASE Z-TS-WMSM-L10 3.5 (White Metal Scanmarker L10 with Screw)	ZBED1131L10	15			10-12
ZZ BASE Z-TS-WMSM-L10 4.5 (White Metal Scanmarker L10 with Screw)	ZBED1132L10	15			10-12
ZZ BASE Z-TS-WMSM-L10 5.7 (White Metal Scanmarker L10 with Screw)	ZBED1133L10	15			10-12
ZZ BASE Z-TS-WMSM-L12 3.5 (White Metal Scanmarker L12 with Screw)	ZBED1131L12	15			10-12
ZZ BASE Z-TS-WMSM-L12 4.5 (White Metal Scanmarker L12 with Screw)	ZBED1132L12	15			10-12
ZZ BASE Z-TS-WMSM-L12 5.7 (White Metal Scanmarker L12 with Screw)	ZBED1133L12	15			10-12
ZZ BASE Z-TS-WMSM-L15 3.5 (White Metal Scanmarker L15 with Screw)	ZBED1131L15	15			10-12
ZZ BASE Z-TS-WMSM-L15 4.5 (White Metal Scanmarker L15 with Screw)	ZBED1132L15	15			10-12
ZZ BASE Z-TS-WMSM-L15 5.7 (White Metal Scanmarker L15 with Screw)	ZBED1133L15	15			10-12
ZZ BASE Z-TS-ZSC 3.5 (Zirconia screw)	BSAE1131	30			10-12
ZZ BASE Z-TS-ZSC 4.5 (Zirconia screw)	BSAE1132	30			10-12

ZIMMER DENTAL® TAPERED SCREW-VENT®

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE Z-TS-ZSC 5.7 (Zirconia screw)	BSAE1133	30			10-12
ZZ BASE Z-TS-ZSC-G 3.5 (Zirconia Screw Gold)	BSCE1131	20			10-12
ZZ BASE Z-TS-ZSC-G 4.5 (Zirconia Screw Gold)	BSCE1132	20			10-12
ZZ BASE Z-TS-ZSC-G 5.7 (Zirconia Screw Gold)	BSCE1133	20			10-12
ZZ BASE Z-TS-ZZMUASC 3.5 (Zirkonzahn Multi Unit Abutment Screw)	BSEE1101	25			10-12
ZZ BASE Z-TS-ZZMUASC 4.5 (Zirkonzahn Multi Unit Abutment Screw)	BSEE1102	25			10-12
ZZ BASE Z-TS-ZZMUASC 5.7 (Zirkonzahn Multi Unit Abutment Screw)	BSEE1103	25			10-12

ZIMMER DENTAL® TAPERED SCREW-VENT® MULTI UNIT ABUTMENT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE Z-TS-MUA-CCTB (Conical cemented titan base non hex and screw)	BSAD6031	30			10-12
SET ZZ BASE Z-TS-MUA-CCTB-K85 (Conical cemented titan base non hex K85 and screw)	BSAD6031K85	30			10-12
SET ZZ BASE Z-TS-MUA-ICOT (Impression coping open tray non hex and Screw and Driver)	BSAQ6061	15			10-12
SET ZZ BASE Z-TS-MUA-SM (Scanmarker and screw)	ZBAD6031	30			10-12
SET ZZ BASE Z-TS-MUA-WSM (White Scanmarker and screw)	ZBCD6031	30			10-12
SET ZZ BASE Z-TS-SA MUA (Scananalog and screw)	ZBDE6031	15			10-12
ZZ BASE Z-TS-MUA-ICSC (Impression coping Screw)	BSFE6001	15			10-12
ZZ BASE Z-TS-MUA-SC (Screw)	BSAE6001	30			10-12
ZZ BASE Z-TS-MUA-WMSM-L12 (White Metal Scanmarker L12 with Screw)	ZBED6031L12	15			10-12
ZZ BASE Z-TS-MUA-ZSC (Zirconia screw)	BSAE6031	30			10-12

ZIRKONZAHN ACCESSORIES

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
Universal Screwdriver Set	SIAK6011				10-12
White Metal Scanmarker DSM L6	ZBFD0101L06	15			10-12
ZZ BASE ZZ-A-IC OT-Driver (Impression coping open tray Driver)	BSGE0001				10-12
ZZ BASE ZZ-A-SA SC (Scanalog Screw)	BSDE0001	15			10-12

ZIRKONZAHN LOC-CONNECTOR

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin temporaries
ZZ-LO-BAR 1.8 (ZZ LOC-Connector Bar)	BSAR10432	25			10-12
ZZ-LO-BAR 2.5 (ZZ LOC-Connector Bar)	BSAR10433	25			10-12
ZZ-LO-BAR 3.86 (ZZ LOC-Connector Bar)	BSAR10431	25			10-12

ZIRKONZAHN MULTI UNIT ABUTMENT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
SET ZZ BASE ZZ-MUA-CCTB (Conical cemented titan base non hex and screw)	BSAD9231	25			10-12
SET ZZ BASE ZZ-MUA-CCTB-K85 (Conical cemented titan base non hex K85 and screw)	BSAD9241K85	25			10-12
SET ZZ BASE ZZ-MUA-CCTB-L10 (Conical cemented titan base non hex L10 and screw)	BSAD9231L10	25			10-12
SET ZZ BASE ZZ-MUA-IC (Impression coping non hex and impression coping screw)	BSAQ9231	10			10-12
SET ZZ BASE ZZ-MUA-ICH-OP (Impression coping hex and impression coping screw)	BSAQ9233	10			10-12
SET ZZ BASE ZZ-MUA-ICH-SP (Impression coping hex and impression coping screw)	BSAQ9232	10			10-12
SET ZZ BASE ZZ-MUA-ICOT (Impression coping open tray non hex and Screw and Driver)	BSAQ9261	15			10-12
SET ZZ BASE ZZ-MUA-ICOTH OP (Impression coping open tray hex and Screw and Driver)	BSAQ9243	15			10-12
SET ZZ BASE ZZ-MUA-ICOTH SP (Impression coping open tray hex and Screw and Driver)	BSAQ9242	15			10-12
SET ZZ BASE ZZ-MUA-NTBH-OP (Narrow titan base hex and screw)	BSAP9283	25			10-12
SET ZZ BASE ZZ-MUA-NTBH-SP (Narrow titan base hex and screw)	BSAP9282	25			10-12
SET ZZ BASE ZZ-MUA-SA (Scanalog and screw)	ZBDE9231	15			10-12
SET ZZ BASE ZZ-MUA-SM (Scanmarker and screw)	ZBAD9231	25			10-12
SET ZZ BASE ZZ-MUA-SM OP (Scanmarker and screw)	ZBAD9233	25			10-12
SET ZZ BASE ZZ-MUA-SM SP (Scanmarker and screw)	ZBAD9232	25			10-12
SET ZZ BASE ZZ-MUA-WSM (White scanmarker and screw)	ZBCD9231	25			10-12
SET ZZ BASE ZZ-MUA-WSM OP (White Scanmarker and screw)	ZBCD9233	25			10-12
SET ZZ BASE ZZ-MUA-WSM SP (White Scanmarker and screw)	ZBCD9232	25			10-12
ZZ BASE ZZ-MUA-BSC (Black Screw)	BSIE9201	25			10-12
ZZ BASE ZZ-MUA-BZSC (Black Zirconia screw)	BSIE9231	25			10-12
ZZ BASE ZZ-MUA-CSC (Cover Screw)	BSAS9201	10			10-12
ZZ BASE ZZ-MUA-CSC-GO (Cover Screw Gold)	BSAS9201GO	10			10-12
ZZ BASE ZZ-MUA-CSC-PI (Cover Screw Pink)	BSAS9201PI	10			10-12
ZZ BASE ZZ-MUA-ICSC (Impression coping Screw)	BSFE9201	15			10-12
ZZ BASE ZZ-MUA-ICSC OP (Impression coping Screw)	BSFE9203	15			10-12
ZZ BASE ZZ-MUA-ICSC SP (Impression coping Screw)	BSFE9202	15			10-12
ZZ BASE ZZ-MUA-LBSC (Lab Screw)	BSHE9201	25			10-12
ZZ BASE ZZ-MUA-LBZSC (Lab Zirconia screw)	BSHE9231	25			10-12
ZZ BASE ZZ-MUA-MI (Model implant)	MIAG9201				10-12

ZIRKONZAHN MULTI UNIT ABUTMENT

Item description	Item number	Recommended torque in Ncm Screw/Zirconia Screw	Recommended torque in Ncm Multi Unit Abutment	Recommended torque in Ncm LOC-Connector	Recommended torque in Ncm composite resin tempraries
ZZ BASE ZZ-MUA-SC (Screw)	BSAE9201	25			10-12
ZZ BASE ZZ-MUA-WMSM-L10 (White Metal Scanmarker L10 with Screw)	ZBED9231L10	15			10-12
ZZ BASE ZZ-MUA-WMSM-L12 (White Metal Scanmarker L12 with Screw)	ZBED9231L12	15			10-12
ZZ BASE ZZ-MUA-WMSM-L12 OP (White Metal Scanmarker L12 with Screw)	ZBED9233L12	15			10-12
ZZ BASE ZZ-MUA-WMSM-L12 SP (White Metal Scanmarker L12 with Screw)	ZBED9232L12	15			10-12
ZZ BASE ZZ-MUA-ZSC (Zirconia Screw)	BSAE9231	25			10-12
ZZ-MUA-ZZLO-L25 4.8 (ZZ LOC-Connector L25)	BSAM9231L25			25	10-12